

JAN 24 1997

K964694

**BOEHRINGER
MANNHEIM
CORPORATION**

510(k) Summary

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

1. Submitter name, address, contact Boehringer Mannheim Corporation
2400 Bisso Lane
P.O. Box 4117
Concord, CA 94524-4117
(510) 674 - 0690 extension 8240
FAX 510 687 - 1850

Contact Person: Betsy Soares-Maddox

Date Prepared: November 19, 1996

2. Device name Proprietary name: Elecsys® LH Assay

Common name: Electrochemiluminescence assay for the determination of human luteinizing hormone (LH).

Classification name: System, Test, Human Luteinizing Hormone

3. Predicate device We claim substantial equivalence to the Enzymun® LH Assay (K900799).

4. Device Description Sandwich principle. Total duration of assay: 18 minutes.
•1st incubation (9min.): 20 µL of sample, a biotinylated monoclonal LH-specific antibody (75 µL) and a monoclonal LH-specific antibody labeled with a ruthenium complex (75 µL)** react to form a sandwich complex.
•2nd incubation (9min.): after addition of streptavidin-coated microparticles (30 µL), the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
**Tris(2,2'-bipyridyl)ruthenium(II) complex (Ru(bpy)²⁺)

Continued on next page

510(k) Summary, Continued

**4.
Device
Description**

- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier (0.4 second read frame).
 - Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent bar code.
-

**5.
Intended use**

Immunoassay for the in vitro quantitative determination of human luteinizing hormone in human serum and plasma.

**6.
Comparison
to predicate
device**

The Boehringer Mannheim Elecsys® LH Assay is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed Enzymun® LH Assay (K900779).

The following table compares the Elecsys® LH Assay with the predicate device, Enzymun® LH Assay. Specific data on the performance of the test have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device is provided in attachment 6.

Similarities:

- Intended Use: Immunoassay for the in vitro quantitative determination of human luteinizing hormone (LH)
 - Sample type: Serum and plasma
 - Antibody: Same pair of monoclonal mouse anti-LH antibodies
 - Solid phase binding principle: Streptavidin/Biotin
-

Continued on next page

510(k) Summary, Continued

6. Comparison to predicate device cont.

Differences:

Feature	Elecsys® LH	Enzymun-Test® LH
Assay Standardization	Enzmun® LH	WHO # 80/552
Detection method	Electrochemiluminescence	ELISA/1-step sandwich assay using streptavidin technology
Instrument required	Elecsys® 2010	ES 300
Calibration Stability	A calibration is recommended every 7 days if kits is not consumed; 4 weeks with same reagent lot if reagent is consumed within 7 days.	Full calibration required every 2 weeks. One-point calibration required every run.

Performance Characteristics:

Feature	Elecsys® LH			Enzymun-Test® LH		
Precision	Modified NCCLS (mIU/mL):			Modified NCCLS (mIU/mL):		
Level	<u>Low</u>	<u>Mid</u>	<u>High</u>	<u>Low</u>	<u>Mid</u>	<u>High</u>
N	60	60	60	120	120	120
Within-Run Mean	0.54	9.38	50.72	3.6	13.5	59.5
%CV	1.82	1.13	0.81	2.9	3.8	1.5
Total Mean	0.54	9.38	50.72	3.6	13.5	59.5
%CV	5.17	1.97	1.99	4.4	4.7	3.9
Lower Detection Limit	0.10 mIU/mL			0.50 mIU/mL		

Continued on next page

510(k) Summary, Continued

6. Comparison to predicate device, (cont.)

Performance Characteristics:

Feature	Elecsys® LH	Enzymun-Test® LH																								
Linearity	0.1 - 200 mIU/mL (with a deviation from a linear line of ±10%)	0.5 - 150 mIU/mL (with a deviation from a linear line of ±10%)																								
Method Comparison	<p>Vs Enzymun-Test® LH</p> <p><u>Least Squares</u></p> <p>y = 1.00x - 0.199</p> <p>r = 0.993</p> <p>SEE = 1.141</p> <p>N = 166</p> <p><u>Passing/Bablok</u></p> <p>y = 0.964x + 0.040</p> <p>r = 0.993</p> <p>SEE = 0.456</p> <p>N = 166</p>	<p>Vs Enzymun-Test® LH</p> <p><u>Least Squares</u></p> <p>y = 0.93x + 0.42</p> <p>r = 0.953</p> <p>SEE = 5.079</p> <p>N = 62</p>																								
Interfering substances	No interference at:	No interference at:																								
Bilirubin	25.0 mg/dL	64.5 mg/dL																								
Hemoglobin	1 g/dL	1 g/dL																								
Lipemia	1500 mg/dL	1250 mg/dL																								
Biotin	30 ng/mL	40 ng/mL																								
Rheumatoid Factor	1500 U/mL	no interference																								
Specificity	<table><tr><th>Level tested</th><th>% Cross-reactivity</th></tr><tr><td>FSH</td><td>300 mIU/mL < 0.1</td></tr><tr><td>HCG</td><td>600 IU/mL < 0.1</td></tr><tr><td>TSH</td><td>300 µIU/mL < 0.1</td></tr><tr><td>HGH</td><td>600 µIU/mL < 0.1</td></tr><tr><td>HPL</td><td>13.80 pmol/mL < 0.1</td></tr></table>	Level tested	% Cross-reactivity	FSH	300 mIU/mL < 0.1	HCG	600 IU/mL < 0.1	TSH	300 µIU/mL < 0.1	HGH	600 µIU/mL < 0.1	HPL	13.80 pmol/mL < 0.1	<table><tr><th>Level tested</th><th>% Cross-reactivity</th></tr><tr><td>200 mU/mL</td><td>0.00</td></tr><tr><td>200 U/mL</td><td>0.00</td></tr><tr><td>100 µU/mL</td><td>0.00</td></tr><tr><td>200 ng/mL</td><td>0.00</td></tr><tr><td>---</td><td>---</td></tr></table>	Level tested	% Cross-reactivity	200 mU/mL	0.00	200 U/mL	0.00	100 µU/mL	0.00	200 ng/mL	0.00	---	---
Level tested	% Cross-reactivity																									
FSH	300 mIU/mL < 0.1																									
HCG	600 IU/mL < 0.1																									
TSH	300 µIU/mL < 0.1																									
HGH	600 µIU/mL < 0.1																									
HPL	13.80 pmol/mL < 0.1																									
Level tested	% Cross-reactivity																									
200 mU/mL	0.00																									
200 U/mL	0.00																									
100 µU/mL	0.00																									
200 ng/mL	0.00																									
---	---																									