

MAR 21 2002

## 510(k) Summary

K020635

### 1. Company Identification

Eastman Kodak Company  
343 State Street  
Rochester, NY 14650

Establishment Registration: 1317267

### 2. Contact Person

Susan Pate  
Regulatory Affairs Associate  
(716)-781-6314

### 3. 510(k) Summary Preparation Date

2/26/2002

### 4. Device Name

Trade Name: KODAK DirectView CR 800 System  
KODAK DirectView CR 900 System

Common Name: Storage phosphor reader with software modification

### 5. Device Classification

Class II

### 6. Indications for Use

The KODAK DirectView CR 800 System and the KODAK DirectView CR 900 System are compact laser scanners capable of reading the latent image formed on a storage phosphor imaging plate and producing a digital image for projection radiology applications.

### 7. Description of Device

The KODAK DirectView CR 800 is a single cassette computed radiography (CR) system while the KODAK DirectView CR 900 is an autoloading CR system. Each CR system contains a storage phosphor reader, CR cassette, HeNe laser, QA workstation with image processing software, and a DICOM network connection.

These CR systems are capable of reading the latent image formed on a storage phosphor imaging plate in a CR cassette and producing a digital image. The system scans the phosphor plate, communicates the image data to a host (e.g. soft copy review workstation or laser printer) and erases the residual image on the plate. An internal imaging plate eraser removes the residual image. The plate is then returned to the cassette, ready for the next exposure.

### **Software Modification**

Anti-scatter grids are commonly used in radiography to reduce scattered x-rays and improve image contrast and signal-to-noise ratio. The new software modification uses a grid detection and suppression algorithm to reduce line artifacts caused by the use of anti-scatter grids. These grid artifacts, when present, can be distracting to radiologists when reviewing CR images, particularly in soft copy. The user may select the amount of grid suppression to be applied within a range of 0-15 through preference settings at the operator console of the CR system.

All other features and capabilities of the CR 800 and the CR 900 system remain unchanged by this software modification.

## **8. Substantial Equivalence**

The analytical characteristics of the grid detection and suppression algorithm were evaluated for a wide variety of grids commonly in use in the medical community. The algorithm was effective at reducing grid induced artifacts.

An evaluation was performed, with hard copy CR images, to study the impact of grid suppression on diagnostic image quality. The effect of grid suppression on images with and without grids was assessed across several grid frequencies, orientations and a variety of radiographic exams (e.g. chest, abdomen, spine, shoulder, pelvis). On a scale of 1 to 9 with 7 being acceptable for interpretation, ten radiologists rated the exams as acceptable for interpretation or better when grid suppression was applied. Average scores ranged from 7.2 with no grid suppression to 7.7 with maximum grid suppression applied. The differences in rating for with and without grid suppression were not clinically significant.

The KODAK DirectView CR 800 and KODAK DirectView CR 900 with grid suppression are substantially equivalent to the previously cleared device [K923544].



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 21 2002

Ms. Susan Pate  
Regulatory Affairs Associate  
Eastman Kodak Company  
343 State Street  
ROCHESTER NY 14650

Re: K020635  
Trade/Device Name: Kodak DirectView CR 800  
and CR 900 Systems  
Regulation Number: 21 CFR 892.1630  
Regulation Name: Electrostatic x-ray imaging system  
Regulatory Class: II  
Product Code: 90 MQB  
Dated: February 26, 2002  
Received: February 27, 2002

Dear Ms. Pate:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Attachment 9 - Indications for Use

510(k) Number (if known): ~~K020635~~ K020635

Device Name: KODAK DirectView CR 800 and KODAK DirectView CR 900 Systems

Indications of Use: The KODAK DirectView CR 800 System and the KODAK Direct View CR 900 System are compact laser scanners capable of reading the latent image formed on a storage phosphor imaging plate and producing a digital image for projection radiography applications.

---

Concurrence of CDRH, Office of Device Evaluation

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter

  
(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K020635