

OCT - 7 2003

K031413

**510(k) Summary**  
**for the VigorMist COMPRESSOR NEBULIZER**  
*(per 21CFR807.92)*

**1. SUBMITTER**

Vega Technologies Inc.  
11F-13, No.100  
Chang Chun Road  
Taipei City  
Taiwan 104

Contact Person: Joseph Lu  
Telephone: +886-2541-6996  
Fax: +886-2521-3803  
E-mail: v6996@ms24.hinet.net  
Date Prepared: January 03, 2003

**2. DEVICE NAME**

Proprietary Name: **VigorMist** Compressor Nebulizer  
(Model: CN-01W)  
K031413  
Common/Usual Name: Compressor Nebulizer  
Classification Name: Nebulizer (Direct patient interface)  
Regulation Number: 21 CFR, 868.5630

**3. PREDICATE DEVICE**

- OMRON Compressor Nebulizer (K914836)
- VEGA Ultrasonic Nebulizer (k002831)

#### **4. DEVICE DESCRIPTION**

The VigorMist™ Compressor is a small air compressor designed to provide sufficient air pressure and flow to power a hand held nebulizer. It measures 11 7/8" x 7 1/8" x 4 1/2" and weights 4.6 lbs.

The device has a thermal protector that will automatically shut off the device when overheated. The operating components are located internally. The compressor, some minor wiring and exhaust/intake tubing are located inside.. External components include a switch, filter with housing, AC cord and cover or accessories compartment. The accessories coming with the compressor include a nebulizer, Air tube and mouthpiece.

#### **5. INTENDED USE**

The VigorMist™ Compressor is intended to be used with a compatible pneumatic nebulizer (specifically the WestMed VixOne) to convert certain inhalable drugs into an aerosol form for inhalation by a patient for the treatment of asthma, COPD, and other respiratory ailments. The device is intended for the home care market. It is intended for use with a single adult, pediatric, or infant patient.

#### **6. PERFORMANCE TESTING**

Testing provided in this premarket notification includes biocompatibility, standards conformity.. The VigorMist™ successfully completed testing and demonstrates that the product fulfills performance specifications.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

OCT - 7 2003

Vega Technologies, Incorporated  
C/O Mr. Tzu – Wei Li  
Responsible Third Party Official  
Center for Measurement Standards of Industrial Technology Research Institute  
Bldg 16, 321 Kuang Fu Road Sec. 2  
Hsinchu  
TAIWAN 30042, R.O.C.

Re: K031413  
Trade/Device Name: Vigormist Compressor Nebulizer (Model: CN-01W)  
Regulation Number: 868.6250  
Regulation Name: Portable Air Compressor  
Regulatory Class: II  
Product Code: BTL, CAF  
Dated: September 18, 2003  
Received: September 23, 2003

Dear Mr. Li:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4646. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: K031413

Device Name: **VIGORMIST Compressor Nebulizer** (model: CN-01W)

Indications For Use:

The **VigorMist Compressor** is intended to be used with a compatible pneumatic nebulizer (specifically the WestMed VixOne) to convert certain inhalable drugs into an aerosol form for inhalation by a patient for the treatment of asthma, COPD, and other respiratory ailments. The device is intended for the home care market. It is intended for use with a single adult, pediatric, or infant patient.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

510(k) Number: K031413

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)