

MAR 24 2004

K 040051

1/3

## 510(k) Summary

The following 510(k) summary has been prepared pursuant to requirements specified in 21CFR§807.92(a).

### 807.92(a)(1)

#### **Submitter Information**

Name: Eagle Technology, Inc.  
Address: 1-2-23 Hirosaka, Kanazawa, Ishikawa, Japan 920-1331  
Phone: +81-76-263-3636  
Facsimile: +81-76-294-6745  
Contact Person: Hideyuki Fujihira  
Date: 1/5/2004

### 807.92(a)(2)

Trade Name: MEGvision EQ1000C Series  
Common Name: Magnetoencephalographic (MEG) Device  
Classification Name(s): Electroencephalograph  
Classification Number: 21CFR882.1400 GWQ

### 807.92(a)(3)

#### **Predicate Device(s)**

Device Name	510(k) #	Manufacturer
CTF "Whole-Cortex MEG System"	K971329	CTF Systems, Inc.
Neuromag-122	K962764	Neuromag Ltd.

3 pages

Eagle Technology, Inc.

807.92 (a)(4)**Device Description**

The MEGvision integrates up to 320 dc-SQUID axial gradiometer with PCs and data acquisition software in order to measure the magnetic signals generated by the intercellular dendritic currents. These detectors positioned in a helmet shaped array give the user the ability to record the electrical activity of the entire surface of the brain cortex simultaneously without having to move the position of the sensor array.

807.92(a)(5)**Intended Use(s)**

The MEGvision non-invasively measures the magnetoencephalographic (MEG) signals produced by the electrical activities by the tissue activities of the brain. These signals, position, direction, and sensitivity of the sensors are acquired and displayed, and may be interpreted by trained clinicians to help localize these active areas. The locations may be correlated to anatomical structure of the brain.

807.92(a)(6)**Technological Characteristics****TABLE 1.1 Comparison to the Predicate Devices**

	<b>Eagle Technology, Inc. MEGvision</b>	<b>CTF Systems, Inc. "Whole-Cortex MEG System" (K971329)</b>	<b>Neuromag Ltd. "Neuromag-122" (K962764)</b>
No. of SQUID detectors/ channels for MEG data:	64 to 320	64 to 200	122
Operating Principle	Superconducting flux transformer coupled with dc-SQUID driven by digitally controlled analog flux locked loop circuit	Superconducting flux transformer coupled with dc-SQUID controlled by digital flux-locked loop	Superconducting flux transformer coupled with dc-SQUID controlled by analog flux-locked loop
No. of auxiliary channels for other types of data	166	88	166
Gradiometer:	1 axial first order gradiometer per location	1 axial first order gradiometer per location	2 orthogonal planar first order gradiometers per location
Intersensor spacing	20mm to 25mm (160 sensor configuration)	32 mm (150 sensor configuration)	43-44 mm
Gradiometer placement	64 to 320 location distributed across the helmet shaped lower tip of a dewar	64 to 200 locations distributed across the helmet shaped lower tip of a dewar (optional Caucasian or Oriental head shape)	61 locations distributed across the helmet shaped lower tip of a dewar.

	<b>Eagle Technology, Inc. MEGvision</b>	<b>CTF Systems, Inc. "Whole-Cortex MEG System" (K971329)</b>	<b>Neuromag Ltd. "Neuromag-122" (K962764)</b>
Cryogen used:	Liquid helium	Liquid helium	Liquid helium
Coverage	One acquisition to cover entire head	One acquisition to cover entire head	One acquisition to cover entire head
Gantry	Floor mounted fixed gantry.	Floor mounted, standard gantry is fixed. Optional gantry tilts to 90 degrees	Floor mounted, standard gantry tilts up to 30 degrees. Optional gantry tilts to 45 degrees
Patient Position	Lying on back	Seated, or lying on back with optional bed	Seated or supine. Optional chair insert for children
Head Position Indicator	Included	Included	Available
Computer	Personal Computer with Windows	HP workstation with UNIX environment	HP workstation with UNIX environment
Networking Capabilities	Ethernet connections to other network system available	Ethernet connections to other workstations included	Ethernet connections to other workstations available
Magnetically Shielded Room Accessories	Interior DC lights, video camera and two-way intercom for patients	Interior DC lights, video camera and monitor and two-way intercom for monitoring patients	Video monitor and two-way intercom for monitoring patients
Intended Use	The <b>MEGvision</b> is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissue in the brain. When interpreted by a trained clinician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.	The <b>CTF "Whole-Cortex MEG System"</b> , is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissue in the brain. When interpreted by a trained clinician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.	The <b>Neuromag-122 system</b> is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissues in the brain. When interpreted by a trained clinician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 24 2004

Eagle Technology, Inc.  
c/o Mr. Lewis Fisher  
Eagle Technology North America, LLC  
25 Bisbee Court, Suite B  
Santa Fe, New Mexico 87508

Re: K040051  
Trade/Device Name: MEGvision EQ 1000C Series  
Regulation Number: 21 CFR 882.1400  
Regulation Name: Electroencephalograph  
Regulatory Class: II  
Product Code: GWQ  
Dated: January 5, 2004  
Received: January 12, 2004

Dear Mr. Fisher:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Lew Fisher

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

*Miriam C. Provost*  
for Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

**510(k) Number (if known):** K040051

**Device Name:** MEGvision EQ1000C Series

**Indications For Use:** The MEGvision is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissue in the brain. When interpreted by a trained technician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.

**Prescription Use** XX  
(Part 21 CFR 801 Subpart D)

**AND/OR**

**Over-The-Counter Use** \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost  
(Division Sign-Off)  
Division of General, Restorative,  
and Neurological Devices

Page 8

**510(k) Number** K040051