

**Model 7250 Jewel® AF
AF Only Clinical Evaluation**

Introduction

Marshall Stanton, MD

**Circulatory System Devices Panel
Presentation**



Dec 5, 2000

confidential

0263 '00 DEC -8 19:54

System Description

Model 7250 Jewel AF



Model 9454 Patient Activator

2

confidential

System Description

Model 7250 Jewel® AF

- Studied /FDA Approved in VT/AT population June 2000
- Comparable to other Medtronic ventricular ICDs
- Includes features intended for treatment and prevention of atrial tachyarrhythmias
 - Atrial Tachyarrhythmia Termination
 - Antitachycardia Pacing (ATP)
 - High Frequency Burst (HFB)
 - Atrial Shock – automatic, or patient activated
 - Atrial Tachyarrhythmia Prevention
 - Atrial Rate Stabilization (ARS)
 - Switch Back Delay ("Cool Down")

3

confidential

System Description

- Model 9464 Patient Activator
 - Battery-powered, radio-frequency hand-held device
 - Used to self-activate shocks for atrial tachyarrhythmias
- Lead Model 6937A
 - 9 French, unipolar high voltage lead designed for placement in the CS or SVC
 - Similar to 6937 lead with slight differences

4

confidential

Model 7250 Jewel® AF Proposed Indications For Use

The Jewel AF system is intended to provide pacing, cardioversion, and defibrillation for treatment of patients with:

- Symptomatic, drug refractory atrial tachyarrhythmias and/or
- Life threatening ventricular tachyarrhythmias

5

confidential

Role of Device-Based Therapy For Atrial Tachyarrhythmias

- Part of an overall treatment strategy
- Symptomatic and drug refractory patients
- Estimated at 5% of total AF population
- For patients who need more control of their therapy
- For patients with frequent hospital visits for cardioversions

6

confidential

Presenters

- Michael Gold, MD, PhD Study Results
- David Schwartzman, MD Case Presentations
- David Newman, MD Quality of Life

- Medtronic, Inc Study Questions
- Dennis Connolly, MD Linda Johnson, PhD
- Scott Brown, PhD Roy Martin, DVM
- Reece Holbrook, BSEE Rahul Mehra, PhD
- Nirav Sheth, PhD

7

confidential

Model 7250 Jewel® AF AF-Only Study

Study Results

Michael Gold, MD, PhD

Circulatory System Devices Panel
Dec. 5, 2000

8

confidential

Study Purpose

- Demonstrate the safety and efficacy of the Model 7250 Jewel AF® in patients suffering from symptomatic, drug refractory atrial tachyarrhythmias *without* ventricular ICD indications.

9

confidential

Study Design

- Multi-center IDE Study
- Prospective follow up study for safety and efficacy of treatment therapies
- Randomized, crossover study for evaluation of prevention therapies
 - 3 months ON, 3 months OFF

10

confidential

Study Inclusion Criteria

- Patient must have experienced 2 atrial fibrillation/flutter episodes within last 3 months
 - 1 episode with ECG documentation
- Episodes must be symptomatic
- Patient must be drug refractory or intolerant
 - patient failed ≥ 1 antiarrhythmic drug
- Patient must be in sinus rhythm at time of implant
 - if post-cardioversion SR ≥ 1 hour

11

confidential

Study Exclusion Criteria

- Chronic AF (< 1hr SR post cardioversion)
- History of uncontrolled angina
- History of sustained VT/VF (>30 seconds of hemodynamically unstable VT/VF)
- NYHA Class IV heart failure
- Cardiac surgery within 1 month
- Atrial thrombus detected within last 6 months prior to implant
- CVA within 1 year
- Unwilling to give informed consent
- Mechanical tricuspid heart valve
- Life expectancy < 1 year

12

confidential

Study Experience

- 144 implanted / 146 enrolled
 - From 11/21/1997 – 11/12/1999
 - United States (107), Europe(33), Canada (6)
- Follow-up
 - Database cutoff: May 31, 2000
 - Mean (SD): 12.7 (± 6.1) months
 - Cumulative patient months: 1835
 - Range: 0.1 – 25.9 months

13

confidential

Patient Characteristics

• Male	71%	• AF Type	
• Mean Age	62 \pm 12 yrs	– Paroxysmal AF	35%
• CAD	32%	– Incessant AF*	65%
• MI	19%	• Primary Indication	
• CHF	29%	– AF	74%
• Mean EF	51 \pm 18 %	– AF/Aflutter	23%
• % EF \leq 40%	31%	– Aflutter	3%
• LA Size	46 \pm 8 mm		

*Incessant AF (aka: Persistent AF) is defined as AF that does not spontaneously convert to SR, but can be cardioverted.

14

confidential

Primary Objectives

- Safety
 - Estimate the relative risk (RR) of system/procedure-related complications for the Model 7250 vs the Model 7219D
 - The 95% upper confidence bound will be < 3.0
- Efficacy
 - Estimate the efficacy of atrial tachyarrhythmia termination therapies of the Model 7250
 - The 95% lower confidence bound will be $> 75\%$

15

confidential

Comparison to Historical Control (Model 7219D ICD) - Analysis Methods

- Comparison to 7219D
 - primary safety endpoint: complication-free survival
 - Secondary objective: survival from all-cause mortality
- Adjusted for differences in baseline patient characteristics in the two populations
 - Multi-variate Cox Proportional Hazards Regression Model

16

confidential

Comparison to Historical Control (Model 7219D ICD) Analysis Method

Method:

- Step 1. Compare two populations with univariate analyses (p-value 0.10) to find significant factors
- Step 2. Determine which factors were predictive of the outcome variable with univariate analyses (p-value 0.10 or RR of at least 100 ± 20%) to find significant predictors
- Step 3. Variables identified in Steps 1 & 2 were put into a "best fit" model using a step-wise approach (p-value 0.05 or RR of at least 100 ± 20%)
- Variables in the final model for complication-free survival were:
 - region (US, non-US), gender, CAD, hypertension, cardiomyopathy, NYHA Class, CABG, SMVT
 - Variables in the final model for survival from all-cause mortality were:
 - gender, CAD, MI, hypertension, cardiomyopathy, CHF, NYHA Class, CABG, valve surgery, AFib, AFlutter, SMVT

17

confidential

Safety Primary Endpoint Results

- Relative Risk (RR) of system/procedure related complications for the Model 7250 vs the Model 7219D
 - 1.31 (95% LCB 0.76, 95% UCB: 2.25)
 - Safety Objective Met: 95% UCB 2.25 < 3.0

18

confidential

System/Procedure Related Complications

Adverse Event	Number of Events	Number of Patients (%)
Lead Dislodgements	11	10(6.8%)
Atrial Fibrillation	3	3(2.1%)
Hematoma	2	2(1.4%)
Infection	2	2(1.4%)
Allergic Reaction	1	1(0.7%)
Anxiety	1	1(0.7%)
Patient Unable To Tolerate Therapy	1	1(0.7%)
Inappropriate Detection	1	1(0.7%)
Lead Failure	1	1(0.7%)
Pacemaker Syndrome	1	1(0.7%)
Skin Irritation	1	1(0.7%)
Undersensing	1	1(0.7%)
Total	26	23

19

confidential

Complication-free Survival

- System/procedure related complication-free survival @ 6-months


7250 AF Only Number at risk: 108	86.6% (95%CI: 79.8, 91.3)
7219D Number at risk: 173	91.6% (95% CI: 88.8, 93.8)

20

confidential

AT/AF Rhythm Classification

Programmable Detection Zones



21

confidential

Efficacy: Success/Failure of Termination Therapy

- Device Classification of Termination
 - Success is defined as 5 consecutive sinus beats within 3 minutes of therapy delivery without redetection.
- Reported Therapy Efficacy
 - Crude proportion is the number of episodes terminated by a specific therapy divided by the number of episodes treated with that type of therapy.
 - What is the probability that a randomly selected episode will be terminated?
 - Generalized Estimating Equation Adjusted Proportion (GEE) is the average of the average efficacy for each patient; corrects for multiple episodes in individual patients.
 - What is the probability that a randomly selected episode from a randomly selected patients will be terminated?

22

confidential

Efficacy Primary Endpoint Results

- **91% efficacy for atrial tachyarrhythmia termination therapies**
 - Evaluates all episodes with at least 1 shock in the therapy sequence
 - 1092 of 1200 AT/AF episodes in 107 patients terminated
 - Crude Efficacy Proportion: 91.0%
 - GEE Efficacy Proportion: 85.9% (95% LCB:81.7%, 95% UCB 89.2%)
- **Efficacy Objective Met: LCB 81.7% > 75%**

23

confidential

Secondary Objectives & Additional Analyses

24

confidential

Efficacy – Number of Shocks per Episode

- Number of Atrial Shocks Per Episode
 - 1200 AT/AF episodes in 107 patients
 - Mean # shocks per episode = 1.19
 - 1 shock 86.1%
 - 2 shocks 10.1%
 - 3 shocks 2.1%
 - 4 shocks 0.7%
 - 5 shocks 0.1%
 - 6 shocks 0.2%
 - 10 shocks 0.1%

25

confidential

Efficacy Of Atrial Shock Therapy For AF

- 92.4% Atrial Shock Therapy Efficacy for AF episodes only
 - 1868 total AF episodes in 102 patients
 - 723 in 85 patients treated with shock
 - 668 successfully terminated
 - Crude efficacy proportion: 92.4%
 - GEE proportion: 88.4% (95% CI: 83.3%, 92.1%)

26

confidential



Efficacy – Patient Activated Shocks

- 92.8% Patient Activated Shock Therapy Efficacy for AT/AF
 - Evaluates all episodes with at least 1 patient activated shock in the therapy sequence
 - 519 of 559 AT/AF episodes in 67 patients terminated
 - Crude Efficacy Proportion: 92.8%
 - GEE Efficacy Proportion: 89.1% (95% CI: 83.8%, 92.8%)

27

confidential

Pacing Efficacy			
	Episodes (patients)	Efficacy Crude Proportion	Efficacy GEE
ATP & HFB for AT/AF	1560/4466 (109 patients)	34.9%	28%
ATP for AT	1049/2720 (89 patients)	38.6%	32.1%
HFB for AT/AF	449/2964 (99 patients)	15.1%	12.2%
ATP & HFB for AT	1212/2896 (92 patients)	41.9%	35.5%
HFB for AT	163/1394 (74 patients)	11.7%	10.6%
HFB for AF	286/1570 (86 patients)	18.2%	14.1%




Mortality

- Estimate the relative risk (RR) of death as compared to the 7219D study
 - Adjusted RR for the 7250 is 0.51
 - 0.51(95% LCB 0.12, 95%UCB 2.17)

31

confidential

Comparison of 7250 AF-Only & 7219D Survival Curves



Survival From All-Cause Mortality

- Comparison of 6-month Kaplan-Meier all-cause survival rates

7250 AF-Only	98.6%
Number at risk: 122	(95% CI: 94.5, 99.6)
7219D	96.4%
Number at risk: 192	(95% CI: 93.7, 98.0)

33

confidential

Positive Predictive Value of Atrial Detection

- 98.8% PPV for detection of atrial tachyarrhythmias
 - 4913 total spontaneous atrial episodes detected by the device
 - 4859 atrial episodes appropriately detected (per investigator classification)
 - Crude Proportion PPV: 98.8%
 - GEE Proportion: 98.6% (95%CI: 96.0%, 99.5%)

34

confidential

Effects of Prevention Therapies

- Evaluation of Atrial Rate Stabilization and Switch Back Delay Pacing Features
 - 75 patients completed the randomized 3 months ON vs OFF protocol
 - n = 38 ON-OFF
 - n = 37 OFF-ON
 - No significant difference in frequency of episodes (p=0.72)

35

confidential

Atrial DFTs

- Mean Atrial DFTs
 - two-tiered step-up protocol
 - At implant (n = 86)
 - 6.8 joules \pm 4.8 joules

36

confidential

Atrial Shock Induced VT/VF

- No incidence of atrial shock induced VT/VF
– 0% (95% CI: 0.0%, 0.3%)

37

confidential

Incidence of Spontaneous VT/VF


- 11 (7.6%) patients had 67 spontaneous, appropriately detected ventricular episodes
 - 16 episodes of spontaneous "VF" (CL < 300 ms)
 - 51 episodes of VT (CL 310-360 ms)
- VT/VF therapies
 - 57 episodes in 6 patients successfully treated with ATP or shock
 - 10 episodes in 5 patients terminated spontaneously after 9-220 sec.

38

confidential

Atrial Shock Programming

Physicians consistently programmed shocks ON.



39

confidential

Use of 9464 Patient Activator

- 67 patients used the 9464 Patient Activator to treat 559 episodes
- 90.5% of episodes treated with 1 shock
- 92.8% (519) of episodes were successfully terminated (89.1% GEE)
- Over 70% of patients had patient activated shock programmed "ON" at last contact


40

confidential

Use of 9464 Patient Activator

Patients consistently used the Patient Activator over time; in total, treating 46.6% (GEE) of all long episodes.

Use of Patient Activated Shocks For Episodes > 30 Minutes



41

confidential

6937A Lead Patient Experience

- 55% (79/144) AF Only and 7% (35/530) VT/AT patients implanted with the 6937A lead
- 111 of 114 implanted in the CS
- Mean Follow-up 12.8 months \pm 6.2
- 6937A Lead related AEs: 3 atrial lead dislodgements; 1 subclavian vein thrombosis
- 3 month complication-free survival
 - 97.3% (95% LCB 93.2%, 95% UCB xx.x%)
- Atrial DFTs (n = 63)
 - Mean (SD) 6.2 joules \pm 4.6

42

confidential

Conclusion - Risk vs. Benefit

Risks – Similar to ventricular defibrillators

Benefits – New tool for AF treatment strategy

- Highly effective Atrial Defibrillator
- Maintenance of Sinus Rhythm
- Pacing Therapies offer incremental efficacy with low risk
- Patient Control of therapy
- Enhanced Monitoring for chronic, non-static disease
- Improved Quality of Life
- Physician and Patient Acceptance
- Ventricular Arrhythmia Protection

43

confidential

Conclusion

The 7250 Jewel AF Only Clinical
Evaluation shows
the Jewel AF Model 7250
System is
safe and effective
in patients with
Atrial Tachyarrhythmias.

44

confidential
