



**Figure 1. Testing Algorithm for Small Minipools
(Currently Applies only to Whole Blood)**



¹ Additional Testing Options: Dilute and retest Individual Donations with same NAT
Perform Alternate NAT on individual donations

² Further Testing: Perform Alternate NAT on individual donations
and/or retest for serologic markers after passage of time

**Figure 2. Testing Algorithm Involving Retests on Sub-Pools
(Currently Applies to Both Source Plasma and Whole Blood)**




¹ There can be several levels of deconstruction using original or freshly pooled Subpools.

² Options for Additional Testing

- a) If the reactive deconstruction Subpool(s) were archived from the construction of the Master Pool:
 - repool & retest smallest reactive Subpool with same NAT; If new pool is still reactive, then do alternate NAT on Individual Donations, or dilute individual donations and retest with same NAT, or
 - dilute individual donations and retest with same NAT, or
 - perform alternate NAT on individual donations.
- b) If the reactive deconstruction Subpools were freshly made (i.e. pooling events independent of the original Master Pool):
 - dilute individual donations and retest with same NAT, or
 - perform alternate NAT on individual donations.

Figure 3. Sub-Algorithm for Retesting Involving Sub-Pools



¹ Options for Additional Testing:

Dilute Subpools and retest
Retest Subpools with alternate NAT