

FOOD AND DRUG ADMINISTRATION
 CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

CIRCULATORY SYSTEM DEVICES ADVISORY PANEL

6531 '01 JUL 30 10:14

MEETING

+ + + + +

TUESDAY,

JULY 10, 2001

+ + + + +

**This transcript has not
 been edited and FDA
 makes no representation
 regarding its accuracy**

The panel met at 8:00 a.m. in Salon A of the
 Gaithersburg Marriott Washingtonian Center, 9751
 Washingtonian Boulevard, Gaithersburg, Maryland, Dr.
 Julie Swain, Acting Chairperson, presiding.

PRESENT:

JULIE SWAIN, M.D.	Acting Chairperson
SALIM AZIZ, M.D.	Voting Member
ROBERT M. DACEY	Consumer Representative
MICHAEL DOMANSKI	Consultant
MARK HAIGNEY, M.D.	Consultant
TED KAPTCHUK, O.M.D.	Consultant
WARREN K. LASKEY, M.D.	Voting Member
MICHAEL C. MORTON	Industry Representative
ILEANA PINA, M.D.	Consultant
JANET T. WITTES, Ph.D.	Voting Member
MEGAN MOYNAHAN	Executive Secretary

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
 1323 RHODE ISLAND AVE., N.W.
 WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

A-G-E-N-D-A

Call to Order, by Julie Swain, M.D., Acting Chairperson . . .	4
Introductory Remarks, by Bernie Statland, M.D., Director, Office of Device Evaluation	9
Open Public Hearing	11
Sponsor Presentation:	
Guidant Corporation P010012, Contak CD, and EasyTrak Lead System:	
Introduction and Background, Dale DeVries, VP Clinical and Regulatory Affairs, Guidant	11
Safety and Efficacy, Device Description, Methods and Results, Patrick Yong, Principal Clinical Research Associate, Guidant	21
Tests and Clinical Evaluation, Michael Higgenbottom, M.D., Cardiopulmonary Exercise Core Lab	57
FDA Presentation:	
Helen S. Barold, M.D.	69
Gerry Gray, M.D.	76
Open Committee Discussion and Recommendations	84
Open Public Hearing	222
Voting	223
Lunch	234

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

Sponsor Presentation:

Medtronic Corporation P010015, Medtronic InSync
Atrial Synchronous Biventricular Pacing Device
and Attain Lead System:

Introduction and Background, Marshall
Stanton, M.D., Medical Director,
Medtronic 236

Design/Methodology, William T. Abraham,
M.D., University of Kentucky 239

Safety Results, Anne B. Curtis, M.D.,
University of Florida 257

Efficacy Results and Conclusions,
Dr. Abraham. 269

FDA Presentation:

Mitchell J. Shein 288
Helen S. Barold, M.D. 292

Open Committee Discussion and Recommendations . 302

Voting 380

Open Public Hearing 397

Adjourn 398

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

P-R-O-C-E-E-D-I-N-G-S

(8:02 a.m.)

DR. SWAIN: Would everyone please sign in.

On the outside, there's some sign-up sheets. Welcome this morning. I'm Julie Swain. I'm a cardiovascular surgeon on the faculty at Harvard Medical School, but this year working at NASA on space station research.

And what I'd like to do first is to call this meeting to order, and it's a meeting of the Circulatory System Devices Panel, and this morning we're going to deal with the Guidant Contak CD and EasyTrak Lead System for the treatment of congestive heart failure.

There's a slight change in the printed agenda that you have for this afternoon, in that we'll have a break after the FDA presentation, which will be approximately 3:45 p.m.

First of all, I'd like to have our panel members and the people sitting at the table introduce themselves and their institution and their area of specialty.

MR. DILLARD: Jim Dillard. I'm the Director

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of the Division of Cardiovascular and Respiratory
2 Devices at the Food and Drug Administration.

3 DR. DOMANSKI: I'm Mike Domanski. I'm a
4 cardiologist at the National Heart, Lung, and Blood
5 Institute.

6 DR. LASKEY: Warren Laskey. I'm an
7 interventional cardiologist at the University of
8 Maryland.

9 DR. PINA: Ileana Pina, Heart Failure
10 Transplantation, Case Western Reserve in Cleveland.

11 DR. HAIGNEY: Bart Haigney. I'm Director of
12 Cardiology at Uniformed Services in Bethesda.

13 DR. KRUCOFF: Mitch Krucoff. I'm an
14 interventional cardiologist at Duke Medical Center and
15 the Director of Interventional Devices, Clinical
16 Trials at the Duke Clinical Research Institute.

17 DR. WITTES: I'm Janet Wittes, a
18 statistician at Statistics Collaborative in D.C.

19 MS. MOYNAHAN: I'm Megan Moynahan. I'm the
20 Executive Secretary of the Circulatory System Devices
21 Panel.

22 DR. AZIZ: I'm Salim Aziz, a cardiovascular

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 surgeon, University of Colorado, Denver.

2 DR. KAPTCHUK: Ted Kaptchuk, Assistant
3 Professor of Medicine, Harvard Medical School.

4 MR. MORTON: I'm Michael Morton. I'm with
5 W.L. Gore and Associates. I'm the industry
6 representative.

7 MR. DACEY: I'm Robert Dacey, consumer
8 representative from Boulder, Colorado.

9 DR. SWAIN: Thank you. Ms. Moynahan, read
10 the conflict of interest.

11 MS. MOYNAHAN: The following announcement
12 addresses conflict of interest issues associated with
13 this meeting and is made part of the record to
14 preclude even the appearance of an impropriety.

15 To determine if any conflict existed, the
16 agency reviewed the submitted agenda for this meeting
17 and all financial interests reported by the committee
18 participants.

19 The conflict of interest statute prohibits
20 special government employees from participating in
21 matters that could affect their or their employer's
22 financial interest.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The agency has determined, however, that the
2 participation of certain members and consultants, the
3 need for whose services outweighs the potential
4 conflict of interest involved, is in the best interest
5 of the government.

6 We would like to note for the record that
7 the agency took into consideration matters regarding
8 Dr. Salim Aziz, Warren Laskey, Mitchell Krucoff, Mark
9 Haigney and Ileana Pina.

10 Each of these panelists reported interests
11 in firms at issue, but in matters that are concluded
12 or not related to today's agenda. The agency has
13 determined, therefore, that they may participate fully
14 in all discussions.

15 The agency also would like to note that, due
16 to the regulations governing covered relationships,
17 the panel chair, Dr. Cynthia Tracey, will not
18 participate in today's deliberations.

19 In the event that the discussion involves
20 any other products or firms not already on the agenda,
21 for which an FDA participant has a financial interest,
22 the participant should excuse him or herself from such

1 involvement and the exclusion will be noted for the
2 record.

3 With respect to all other participants, we
4 ask in the interest of fairness that all persons
5 making statements or presentations disclose any
6 current or previous financial involvement with any
7 firm whose products they may wish to comment upon.

8 DR. SWAIN: And just to remind you that
9 anyone speaking in the public session or for the
10 companies, before you speak, please introduce
11 yourself, your name and position, whether you're an
12 owner or own stock in a company, whether you're an
13 employee of a company that's related to these devices,
14 or whether you're paid for your transportation here or
15 research grants, or things of that sort.

16 The first part of the meeting is an open
17 public -- excuse me. One more. Voting status statement
18 by Ms. Moynahan.

19 MS. MOYNAHAN: "Pursuant to the authority
20 granted under the Medical Devices Advisory Committee
21 Charter, dated October 27, 1990, and as amended August
22 18, 1999, I appoint the following individuals as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 voting members of the Circulatory System Devices Panel
2 for this meeting on July 10, 2001: Mitchell Krucoff,
3 Michael Domanski, Julie Swain, Ted Kaptchuk, Ileana
4 Pina, and Mark Haigney. In addition, I appoint Dr.
5 Julie Swain to serve as panel chair for the duration
6 of this meeting.

7 For the record, Dr. Pina is a consultant to
8 the Cardiovascular and Renal Drugs Advisory Committee
9 of the Center for Drug Evaluation and Research, and
10 the other individuals are consultants to this panel.

11 They are all special government employees
12 and have undergone the customary conflict of interest
13 review and have reviewed the material to be considered
14 at this meeting." And it's signed by Dr. David W.
15 Feigal, Director, Center for Devices and Radiological
16 Health.

17 DR. SWAIN: Thank you. Now the first part
18 of the meeting is an open public hearing and --

19 MS. MOYNAHAN: Sorry. I have one other
20 thing. We'll have a few introductory remarks by Dr.
21 Bernie Statland. Dr. Statland is the Director of
22 FDA's Office of Device Evaluation.

1 DR. STATLAND: Good morning. My name is
2 Bernie Statland. I'm the Director of the Office of
3 Device Evaluation. My job is very easy. The people
4 around me do all the work, and I have the opportunity
5 of greeting you.

6 It's a very exciting day for all of us and,
7 hopefully, by the end of the day we'll have learned a
8 lot, we'll have gained a lot, and we'll move forward.

9 I would like to really extend three thank
10 yous as we start off the day. I would first of all
11 like to thank the panel members for coming here,
12 giving up of your time, your expertise and your
13 participation. Without all of you, we would not be
14 able to move forward.

15 And, second, I would like to thank the
16 individuals within the FDA, within the division headed
17 by Jim Dillard, the Division of Cardiovascular and
18 Respiratory Diseases, DCRD, for the tremendous amount
19 of effort that they have placed.

20 But last but not least, I would like to
21 thank the companies, industry, that really have
22 invested their intellectual capital, their financial

1 capital and their initiative to bring forward devices
2 that, hopefully, today and also other devices we'll
3 hear about that will make a difference.

4 And without eating up any more of your
5 valuable time, I welcome you. I hope it's a productive
6 meeting, a meeting where there will be a lot of give
7 and take. We'll learn something from it. And thank
8 you all very much.

9 DR. SWAIN: Thank you. Okay. Finally, we'll
10 get to the open public hearing part, and there were no
11 prior requests to speak. Is there anyone in the
12 audience who wishes to address the panel on this
13 morning's topic?

14 If not, we will close the open public
15 hearing part and we'll start with the sponsor's
16 presentation and, again, remind you to introduce
17 yourself, and your position, and any conflict. And
18 this lasts approximately one hour.

19 MR. DeVRIES: Good morning. My name is Dale
20 DeVries. I'm Vice President of Clinical and Regulatory
21 Affairs for Guidant Corporation. I do own Guidant
22 stock.

1 It's my pleasure to be addressing this panel
2 and starting the discussion related to the first-ever
3 PMA device that may be considered for approval, where
4 we have cardiac resynchronization therapy in
5 combination with an ICD.

6 I would also like to express my thanks for
7 all the effort that went forward in bringing this to
8 the panel today, in particular, to the FDA, the
9 reviewers and all the staff at the FDA; to the panel
10 members in preparation by review of the materials that
11 were prepared and sent to you and considered for
12 review today; to the clinicians and consultants who
13 will be speaking on our behalf related to this trial;
14 to my Guidant associates; to all the general public
15 and other interested parties that are here today.

16 Why are we here today? Obviously, it's to
17 review the existing evidence for the safety and
18 effectiveness for cardiac resynchronization therapy
19 when combined with an ICD that's already proven.

20 We'd like to confirm that there's a patient
21 population that clearly benefits from this therapy and
22 to make this important therapy available for the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 management of heart failure patients.

2 I'll do a few brief introductory comments.
3 Then we'll turn the presentation over to Pat Yong, a
4 Principle Clinical Research Associate for Guidant. Pat
5 will go through device description, methods and
6 results of our trial.

7 He'll start the presentation by going over
8 the safety and efficacy results for the all-patient
9 population. Then he'll take a few moments to describe
10 the process and rationale that we used in going
11 through a subpatient population. Then Pat will review
12 the safety and effectiveness of the second population.

13 In addition to that, we'll have a
14 presentation by Dr. Higginbotham. He's an expert in
15 exercise testing from Duke University. He will go
16 through some of the tests that we used to evaluate
17 functional status and quality of life and then make a
18 few comments about the clinical evaluation and
19 importance of these measurements.

20 Obviously, when you go through a new device
21 and a new technology and you try to make it available
22 to the general public, there are a lot of things that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 may occur during the course of this journey.

2 One is, as we continue to learn, our basic
3 knowledge level will grow. In addition to that, the
4 evidence will continue to build related to the therapy
5 and devices that we're considering.

6 Also, new tools might become available to
7 physicians along the way. In addition to that, for
8 corporations such as ours, technologies will improve.
9 We'll incorporate the new evidence that we've gathered
10 and include that in the new devices as we bring them
11 forward, and this is occurring at an accelerated rate.

12 In addition to that, the clinical practices
13 that physicians use in their health care of heart
14 failure patients will change.

15 One of the things I wanted to do was just
16 make a few comments about the activities that we had
17 related to the study. First of all, we wanted to do
18 an overview of the safety of this study.

19 Second, the effectiveness of the therapy,
20 and in addition to that, we wanted to consider the
21 benefit and risk associated with this device and this
22 therapy, in particular, related to the Contak CD.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Safety. The device and the procedures are
2 safe, and that is based on the information that we're
3 bringing forward. In addition to that, we needed to
4 make sure that we were not creating any additional
5 harm for the patients who may receive this device.

6 Effectiveness. The trial did not achieve
7 its primary endpoints for effectiveness. In fact, I
8 wish I were standing in front of you here today with
9 a nice package all tied up with a ribbon in a neat
10 bow. It's not true.

11 The Contak -- the CRT trial that we brought
12 forward did not achieve clinical significance.
13 However, there is a significant amount of information
14 collected in this trial related to the clinical
15 benefit, where we have reasonable assurances of the
16 benefit that will be received by the patients.

17 In addition to that, we want to spend some
18 time discussing the rationale and identification of a
19 product of a group of patients that clearly benefit
20 from this therapy.

21 Cardiac resynchronization therapy is
22 effective in patients with moderate to severe heart

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 failure; that is, clinically meaningful results.

2 Benefit versus risk. These patients are
3 already indicated for an IC. This was one of the
4 challenges that we have inside our corporation as we
5 reviewed the results. It's important to remember that
6 these patients are already exposed to the risk of
7 getting an implantable device, in particular, an
8 implantable defibrillator.

9 In addition to that, the benefits of CRT
10 therapy outweigh the risk of placing a left side
11 coronary venous lead.

12 I'd like to make a few comments about the
13 overall process and study chronology for the study.
14 I know some of you may have thought that the clinical
15 section of the panel package that you received was
16 fairly involved.

17 First of all, we started with the Ventak CHF
18 study, the original design. That was a procedure that
19 required opening of the chest to allow placement of
20 the leads. The enrollment was fairly slow at the early
21 phases of the trial.

22 We received feedback from physicians that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the incremental risk of placing the lead was greater
2 than they wanted to have for some of the patients that
3 were very, very sick. We had received this
4 information early on when we were considering heart
5 failure devices.

6 Our response to this was to develop an
7 EasyTrak lead system for placement of the leads on the
8 left side of the heart, without having such an
9 invasive process.

10 In addition to that, we needed to modify the
11 generator so that the connection system for the lead
12 could be managed. This generator was the most recent
13 version of the ICD that we had available on the
14 market.

15 We continued enrollment and enrollment
16 actually accelerated. We completed enrollment, which
17 we refer to here as Phase 1. We enrolled
18 approximately 250 patients.

19 On the very same day, ironically, that we
20 completed enrollment in this first study, we received
21 information from the FDA that the requirements for all
22 sponsors of CRT trials had changed. T h e y

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 wanted to make sure that we had a minimum of six
2 months of continuous data on the control patient
3 population. In addition to that, they wanted to have
4 six months of continuous data for the arm of the study
5 where we'd have the therapy included for the patients.

6 We managed to modify the trial, change the
7 endpoints accordingly, and continue to enroll in the
8 trial. Enrollment was completed at the end of 2000.
9 The PMA was summarized based on a January cut-off and
10 submitted to the FDA in February of 2001.

11 In addition to that which is customary for
12 new devices and new therapies of this sort, we were
13 asked to update the panel information related to the
14 most recent cut-off related to this patient population
15 so that you would have that information for
16 consideration in your panel pack.

17 We're here today at the panel meeting. Now
18 several other things occurred along the course of this
19 trial. New drugs became available to physicians in
20 the management of the heart failure population.

21 In addition to that, heart failure
22 physicians and primary care physicians explained to us

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 one of their major goals was to keep heart failure
2 patients out of the hospital. They have managed to do
3 this. They have changed their patient management
4 program for many of the population that have heart
5 failure.

6 This study was no exception. There were
7 multiple changes along the way. These changes were
8 incorporated into trial and the resulting data that we
9 have for you today.

10 So, in summary, there is a subgroup analysis
11 that we have for you to consider. There's strong
12 clinical evidence related to the performance of this
13 device and this therapy.

14 The results are consistent with other
15 trials, and the benefits are incremental to
16 contemporary heart failure therapy treatment and also
17 to patients who would also receive an ICD. There's a
18 strong case for approval.

19 We did not do this trial alone. We had 47
20 centers involved in the trial and all of the support
21 staff at those centers.

22 We have several physicians and consultants

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 with us here today. You can see that we have
2 identified them by name, their clinical discipline,
3 their study involvement, and their title affiliation.

4 As each one of these people speak, we'll ask
5 them to disclose their affiliation, because they're
6 very difficult to read.

7 First of all, Dr. Boehmer. Dr. Boehmer is
8 a heart failure specialist and cardiologist. He's a
9 member of our events committee and a principal
10 investigator in this trial.

11 Dr. Foster is a cardiologist and director of
12 our echo lab. Dr. Steve Higgins is an
13 electrophysiologist and principal investigator in this
14 trial. Dr. Higgins is also one of the largest
15 enrollers and has a large body of information related
16 to the use of this device and the lead system.

17 Dr. Larntz is an independent
18 biostatistician. He was a statistical consultant for
19 this trial, coaching and counseling us related to the
20 activities as to how we summarized our information.

21 In addition to that, Dr. Larntz was also
22 involved in some of the actual calculations and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 statistical reports that we completed.

2 Dr. Mester is a cardiologist who has a lot
3 of experience in vascular intervention. He's a
4 principal investigator in this trial and has one of
5 the highest success rates related to the lead implant
6 and device use.

7 Dr. Mester was also instrumental in the
8 development of some of the training programs related
9 to this new therapy. We also have Dr. Saxon. Dr. Saxon
10 is an electrophysiologist, a principal investigator in
11 this trial, and a consultant in charge of our core
12 lab.

13 We have several representatives from Guidant
14 available for responding to your questions today.
15 This is the group of individuals that will be
16 responding to most of the questions you have today.
17 It's obvious that the benefit of this product is very
18 important and the clinical outcome is what most of
19 judgement related to this product will be today. With
20 that, I'd like to turn the presentation over to Pat
21 Yong.

22 MR. YONG: Hello. My name is Patrick Yong

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 and I'm an employee of the sponsor. I'm going to
2 start out with a therapy and device description, then
3 move to the description of the stent design, and then
4 finally to the endpoints and the associated results
5 that support the indications we seek.

6 Patients who were enrolled in the Ventak
7 CHF/ Contak CD study were all characterized as having
8 dilated cardiomyopathy. Dilated cardiomyopathy is
9 associated with ventricular remodeling and fibroid
10 ingrowth as the heart enlarges.

11 As this fibrotic ingrowth invades the
12 heart's natural conduction system, it may lead to a
13 second important characteristic of this patient
14 population, that of an intraventricular conduction
15 delay, which is indicated by a white curex on the ECG.

16 This intraventricular conduction delay could
17 lead to an asynchronous ventricular contraction with
18 a loss of pumping effectiveness.

19 Cardiac resynchronization therapy or CRT is
20 intended to restore ventricular synchrony by pacing
21 both ventricles simultaneously.

22 Unlike pharmacologic therapies that are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 inotropic in nature, CRT has been shown to improve
2 both ventricular function and increase the ventricular
3 efficiency of the heart.

4 In order to achieve biventricular
5 stimulation, we have to be able to pace the left side
6 of the heart. And we do this by taking advantage of
7 the coronary venous vasculature, which surrounds the
8 surface of the heart. This gives us ready access to
9 the left ventricle.

10 By placing a guide wire into the desired
11 location in the coronary venous vasculature, a lead
12 can be advanced, using the over-the-wire technique,
13 similar to that used in interventional cardiology, to
14 get the lead placed in its desired location.

15 On the left is the investigational EasyTrak
16 lead. It's unipolar, has passive fixation, and has
17 started moving. It rides over the guide wire to its
18 final destination.

19 On the right is the family of guide
20 catheters that are used by the electrophysiologist to
21 help cannulate the ostium of the coronary sinus. It
22 also serves as a conduit through which the lead is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 placed into the heart.

2 The second mechanical component of the
3 system is the Contak CD heart failure device. It's
4 based on the commercially available Ventak AD3 ICD.
5 and, therefore, has all the standard features
6 associated with it: detecting arrhythmias, delivering
7 shock for defibrillation and delivering anti-
8 tachycardia pacing to treat monoformic V-tach. In
9 addition to the EasyTrak lead, this device uses
10 commercially available right atrial and left
11 ventricular cardiofibrillation leads. Because this
12 device delivers both CRT and has ICD capability, we
13 use the acronym CRTD to describe it.

14 In terms of how the device works, how we get
15 biventricular sensing and biventricular stimulation,
16 this diagram shows the header. In addition to the two
17 normal ports that are used to connect the
18 defibrillation leads and connect the right ventricular
19 sensing and right atrial sensing leads, an additional
20 port has been added to accommodate the EasyTrak lead.

21 This biventricular output is hardwired in
22 the header itself, so the left ventricular output and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the right ventricular output are tied together.
2 Therefore, we get both biventricular stimulation and
3 biventricular sensing.

4 In this particular device, the outputs for
5 the left ventricular and right ventricular channels
6 are not independently programmable.

7 As a consequence of making the header this
8 way, we optimize our biventricular and tachycardia
9 pacing through ICD.

10 The device had to undergo a number of
11 stringent tests before we could get to clinical
12 trials, starting out with design verification testing,
13 or DVT. For the Contak CD volt generator, or VG, it
14 consisted of electrical and mechanical testing,
15 testing of battery capacity, electromagnetic
16 capability, as well as a software DVT for both the
17 pulse generator and the program application which is
18 used for programming the device.

19 There was also similar device verification
20 testing that took place with the EasyTrak lead. The
21 tests includes the axial load, electrical resistance,
22 insulation integrity, pacing impedance, and fatigue

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 resistance.

2 We also tested the LV-1 connector, which was
3 used in these track leads, and its compatibility with
4 PG header. And finally, we had to consider how the
5 lead worked with all the implant accessories that are
6 used.

7 Once we considered the pace maker -- the PG
8 load and the lethal load, we now have to consider them
9 as a system. We had system design validation,
10 including a systems feature test that simulated use
11 under real-world conditions. We also had to look at
12 safety risk analysis, including hazard analysis,
13 reliability and friction analysis, and component
14 qualification testing.

15 We also had to consider the biocompatibility
16 evaluation, the PG, the lead and the accessories. All
17 this culminated in the animal studies, which were
18 tested in a feline model.

19 Now we'll turn to the study design. The
20 study was designed to demonstrate the safety and
21 effectiveness of CRTD in the population study. The
22 major criteria in the study included a VT/VF

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 indication for ICD implantation. Therefore, every
2 patient enrolled in the study is going to undergo a
3 device implant.

4 Furthermore, patients had to have
5 symptomatic heart failure, which would be New York
6 Heart Class II through IV, while on heart failure drug
7 therapy.

8 Patients had to have left ventricular
9 dysfunction and an interventricular conduction delay
10 with a measured QRS of at least 120 milliseconds.
11 Furthermore, patients had to be in sinus rhythm with
12 no indication for a bradycardiac pacemaker.

13 In terms of the study's scope, this study
14 was conducted at 47 investigational centers, with 581
15 patients enrolled.

16 Fourteen patients did not undergo an implant
17 procedure. Sixty-six patients did undergo the implant
18 procedure, but did not receive the investigational
19 system, leaving the 501 patients implanted with the
20 investigational system.

21 As Dale DeVries pointed out, this study was
22 conducted in two phases. Phase 1 was the original

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 study design. When the study was originally designed,
2 it was built to answer the question: Does CRT improve
3 chronic functional status?

4 The original design was that of a
5 randomized, double-blind, cross-over. We chose a
6 cross-over design because our initial look was at
7 patients who would be doing a thoracotomy, and we
8 wanted to use the most efficient design possible to
9 minimize the number of patients exposed.

10 To achieve double blinding, we had to turn
11 to a team. We had one investigator, the
12 electrophysiologist, who would be responsible for
13 programming the device and he obviously would know if
14 the patient was programmed.

15 But the second individual, who would be the
16 heart failure specialist, would be responsible for
17 following up the patient and that individual was
18 blinded as to the pacing mode. Finally, the patient
19 themselves are blinded as well.

20 The primary outcome is measured in terms of
21 functional status. The patients were given an
22 exercise test using a modified Naughton protocol on a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 treadmill. This protocol was performed at the
2 randomization visit, again in three months, and again
3 at six months. Peak VO_2 would be used as the primary
4 endpoint. The original sample size required for Phase
5 1 of 248 patients, were enrolled.

6 The diagram at the bottom shows how the
7 study was laid out. What makes the Ventak CHF/Contak
8 CD study unique in its design is the pressing need
9 that patients had for an ICD.

10 All these patients have VT/VF, and it's
11 important to get a defibrillator in them as soon as
12 possible. So it wasn't always possible to insure that
13 the patients were on the right medications or have
14 adequate doses.

15 Therefore, what we did was to implant the
16 device and give the physicians a one-month period of
17 time. This one-month period of time would allow the
18 patient to recover from the surgery, it allowed them
19 to recover from any cardiac arrest or any arrhythmia
20 they had.

21 But it also gives physicians the opportunity
22 to adjust medications as need be before the randomized

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 study is carried out.

2 Patients were then randomized to get either
3 three months of CRT, followed by three months or no
4 CRT; or, alternately, three months of no CRT followed
5 by CRT.

6 At the end of this intense study phase,
7 patients were continued to be followed at three-month
8 intervals for device evaluation and so forth to ensure
9 the device is working properly long term.

10 After we enrolled our original sample size
11 of 248 patients, there was a dialogue between Guidant
12 and the FDA concerning the design of the study that
13 was looking at CRT.

14 We then modified our study and went forward
15 with Phase 2. By looking at Phase 2, we looked at six
16 months of continuous data, rather than a cross-over
17 design with a new primary hypothesis, that of
18 determining: Does CRTD slow heart failure
19 progression?

20 We still retained the elements of the
21 original Phase 1, that we would still continue to look
22 at whether or not CRTD improves chronic functional

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 status.

2 So the features of the new study design are
3 as follows. Instead of a cross-over design with three
4 months duration each, we have a parallel design with
5 six months.

6 We also now consider morbidity and mortality
7 as the primary outcome. Again, we retained the
8 elements of the original design; that is, patients
9 were implanted. We had no CRT for the first month to
10 give physicians an opportunity to follow the patients
11 and adjust the medications, and then the patients were
12 randomized through the six months of CRT or six months
13 of no CRT. Again, at the end of the six-month period,
14 patients would continue to be followed.

15 In terms of how data were integrated between
16 the two phases, the patients enrolled in Phase 1 would
17 contribute data from the first three-month period.
18 Patients enrolled in Phase 2 would contribute data
19 throughout the entire six-month period.

20 Based upon the sample size available at the
21 time of analysis, we have roughly half the patients in
22 Phase 1, half the patients in Phase 2 for a follow-up

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 time four and a half months during the intensive
2 therapy phase.

3 In terms of study organization, because it
4 was important to consider the effect on functional
5 status, we went through a core laboratory and they
6 were used to evaluate all data from exercise testing.

7 We also used the services of a external
8 events committee. It consisted of four
9 electrocardiologists, with three heart failure
10 specialists and one electrophysiologist serving.

11 It was their task to review and adjudicate
12 all deaths and all hospitalizations that took place
13 during the course of the study. Furthermore, this
14 committee was blind to the randomized therapy while
15 they made their deliberations.

16 We also had an independent statistician who
17 provided statistical recommendations and helped
18 perform the covariate analyses.

19 The patient demographics of the Ventak
20 CHF/Contak CD study are very similar to that of a
21 standard defibrillator population with a high
22 proportion of coronary artery disease.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 What's notable is that 58 percent of the
2 patients are New York Heart Class III, and a left
3 bundle branch block was the most predominant sort of
4 interventricular conduction delay.

5 We also find that patients were very well
6 medicated at the time they were treated in the study.
7 We determined that there were no clinically
8 significant differences between the CRT and the
9 control groups at the time of enrollment.

10 Now we'll turn to the endpoints and the
11 study results. We consider our results in three
12 distinct ways. First of all, we have the EasyTrak
13 lead. We'll look at it and its performance.

14 Then we'll look at what happens when we take
15 the EasyTrak lead and combine it with CRTD to see how
16 well it works as a system.

17 Finally, we'll consider CRT and how well it
18 works as a therapy.

19 Starting with the EasyTrak lead, EasyTrak
20 lead safety was determined on the rate of the lead-
21 related adverse event rates. Effectiveness was judged
22 in terms of lead performance, that is, pacing, sensing

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 and impedance. We also considered the implant success
2 rate.

3 Shown on the left is the lead-related
4 adverse event rate of approximately 12 percent. The
5 95 percent confidence -- 95 confidence interval is
6 well within the acceptance boundary. Therefore, we do
7 meet the safety endpoint.

8 Shown on the right are the three most common
9 types of adverse effects associated with the EasyTrak
10 lead. The first one is elevated left ventricular
11 thresholds, which were seen in 29 patients, or 6.5
12 percent overall.

13 Twenty-five of these 29 were resolved with
14 reprogramming -- I'm sorry. With repositioning, or 86
15 percent.

16 There were four patients in whom given the
17 lead could not be repositioned or the investigator
18 elected not to try it. That happened in 14 percent in
19 which the therapy was abandoned.

20 The second most common type was that of
21 double counting. Because of the biventricular sensing
22 it's possible to sense the QRS twice, leading to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 inappropriate therapy.

2 We saw that happen in ten patients, or in
3 2.2 percent. It resulted in lead reprogramming in
4 five, resulted in lead revision in four, and one
5 patient expired before lead revision could be carried
6 out. However, that death was unrelated to the device
7 or the therapy.

8 The third most common was that of coronary
9 venous trauma, which was reported in ten patients, or
10 2.2 percent.

11 In all ten of these situations in the study,
12 no intervention was necessary and no cardiac tamponade
13 was noted. Furthermore, there was no short- or long-
14 term sequelae resulting from these coronary venous
15 traumas. Therefore, the adverse event rate is within
16 our safety standards.

17 Turning now to effectiveness, what we first
18 have to consider is getting the EasyTrak lead into the
19 patient cannot be measured with the implant success
20 rate.

21 Overall, we were able to implant the device
22 in 87 percent of the patients in whom it was tested.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 What we consider the number one reason for
2 not being able to implant the device was either
3 inability to locate the ostium of the coronary sinus
4 or inability to cannulate it and get a stable
5 position. That happened most commonly. When we take
6 that into account, we were able to implant it in 91
7 percent.

8 We also considered the impact of a learning
9 curve. Shown on the right is the patient population
10 divided with quartiles based upon investigator
11 experience. We find that over time with increased
12 investigator experience, the implant's success rate
13 rises over time. Furthermore, the ability to find and
14 cannulate the coronary sinus os also improves as well.
15 So investigator experience improving, we get to 91
16 percent in class success rate.

17 Also of importance is that of procedure
18 time. The skin-to-skin time to place the entire
19 system for the first quartile was about three-and-a-
20 half hours. But as time goes on and investigators get
21 more experienced, the mean time is now reduced to two
22 hours by the time you get to the fourth quartile.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Let's now consider the effectiveness of the
2 lead once it's in place. Starting out with left
3 ventricular thresholds, the 95-percent tolerance
4 interval was well within the boundaries set forth at
5 the outset of the study.

6 Our threshold was about 1.8 to 1.9 volts
7 throughout the study and was remarkably stable over
8 time. The second effective set-point was that of
9 lead impedance. This is the biventricular lead
10 impedance and represents the parallel combination of
11 the left ventricular and right ventricular leads.

12 The 95-percent confidence interval, again,
13 is well above the standard set forth at the start of
14 the study. Similar to the pacing thresholds, the lead
15 impedance is, again, remarkably stable over time.

16 The third effective set-point was that the
17 biventricular R-wave amplitude, or the ability to
18 sense the rhythm. Very similar to the other
19 endpoints, it was also stable over time with a mean
20 value of about ten millivolts. The mean value was,
21 again, well above the acceptance value.

22 So to consider the performance and safety of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the EasyTrak lead, all safety and effectiveness
2 endpoints have been met.

3 Next we consider how the EasyTrak lead and
4 CRTD, when used as a system, interact. In terms of
5 safety for the implanted system, safety was measured
6 in terms of the severe device-related adverse event
7 rate and in terms of outcome mortality.

8 For effectiveness, we considered a
9 combination of CRT and the EasyTrak lead to see how it
10 affected ICD performance. The two regimens used here
11 were detection time and a success rate of anti-
12 tachycardia pacing or ATP for the termination of
13 monomorphic V-tach.

14 Starting out with the safety endpoints,
15 both the severe device-related adverse event rate and
16 the outcome mortality rate were both well within the
17 standards set forth at the outset of the study.
18 Therefore, we meet the safety standard so we've
19 combined the EasyTrak and CRTD together.

20 When we consider the performance of the
21 combination, we started out looking at the induced
22 ventricular fibrillation detection time. Because we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 have two-lead sensing, we wanted to make sure that the
2 combination of two-lead sensing did not do anything to
3 degrade the ability of the device to detect induced VF
4 and, clearly, we see that that doesn't happen.

5 The Contak CD average detection time was 2.2
6 seconds. That compares very favorably to that of a
7 standard fibrillator, the Ventak AV1, of 2.0 seconds.

8 For ATP conversion efficacy, again, we have
9 biventricular ATP. We first considered that of
10 induced MVT, which is either tested at the time of
11 implant, or it could be deferred.

12 This was tested experimentally in 44
13 patients. The conversion rate with CRTD was 64
14 percent. While this was less than we anticipated, it
15 was similar to that of published studies about
16 terminal testing of ATP by using the right side, which
17 rated 59 to 80 percent.

18 We also tested the ability of the device to
19 treat spontaneous MVT. In 196 patients, an empiric ATP
20 scheme was used. The conversion rate of these episodes
21 was 88 percent. Again, very similar to that of
22 published studies for right ventricular ATP between 89

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 and 92 percent.

2 Therefore, ventricular fibrillation sensing
3 and the treatment of ATP with biventricular ATP are
4 similar to that for conventional ICDs.

5 To summarize the system's safety, for the
6 EasyTrak lead, we demonstrated safety and
7 effectiveness of the device. We also find that if we
8 combine the EasyTrak lead as part of a system, that
9 the system still remained safe and IC performance
10 remains robust.

11 We also had additional experience as well in
12 that there are two separate studies that can look at
13 how well the EasyTrak lead system worked.

14 One is the continuation of the Ventak
15 CHF/Contak CD study beyond the therapy phase. That's
16 near completion. We also have the European Registry,
17 which has been completed.

18 Both of these studies provided prospective
19 longitudinal assessment of the Contak CD EasyTrak
20 system and CD monitoring. The focus of these two
21 studies are complementary.

22 From the Ventak CHF/Contak CD study, we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 continue with our analysis on morbidity and mortality.
2 From the European Registry, we get additional
3 information about functional status.

4 Between the US and the European studies, we
5 have over 1,500 patients enrolled. In the United
6 States, we have over 500 patients enrolled at 47
7 centers. The mean follow-up time is 16 months.

8 We have over 100 patients now followed for
9 over two years, with three and a half years being the
10 maximum available to us now. The cumulative
11 experience is over 7900 patient-months.

12 When you consider our European experience,
13 we have enrolled 1,000 patients at 248 centers. The
14 mean follow-up there is four months with a maximum
15 follow-up of 20 months, with an additional 4400
16 patient-months of experience.

17 We'll now consider the third investigational
18 component of the study, that of CRT effectiveness. We
19 had complementary endpoints that were used to evaluate
20 the effectiveness of CRT and they represent
21 complementary perspectives.

22 These endpoints were prospectively powered

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 at 80 percent with a five percent alpha level. The
2 first consideration is that of the progression of
3 heart failure for our primary endpoint. T h i s
4 analysis considered events and the time it took to
5 reach those events.

6 The three components consisted of all-cause
7 mortality, hospitalization for heart failure, which
8 lasted at least 23 hours and VT/VF events which
9 required device intervention.

10 The complementary perspective without a
11 secondary one, was looking at the functional status.
12 The endpoints specified in the investigational plan
13 were peak VO_2 . This was measured with an assisted,
14 limited exercise test performed on a treadmill and
15 quality of life.

16 Quality of life was determined by the
17 Minnesota Living with Heart Failure questionnaire. We
18 also formed ancillary analyses that would help
19 complement the functional status endpoints. We
20 looked at six-minute walk, V_E/VCO_2 slope and changes
21 in the New York Heart Class.

22 We also provided covariate analysis. The

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 purpose of the covariate analysis was defined as any
2 factors other than treatment that could affect the
3 outcome.

4 As stipulated in the protocol, clinically
5 relevant variables were selected before processing any
6 analysis for the primary and secondary endpoints at
7 the conclusion of the study.

8 We utilized the services of physicians who
9 formed our Contak CD events committee. They provided
10 five variables, based on their clinical experience,
11 that were associated by them with the progression of
12 heart failure.

13 The five clinical variables provided were
14 that of New York Heart Class, bundle branch
15 morphology, etiology, whether ischemic or non-
16 ischemic, left ventricular ejection fraction and the
17 QRS width.

18 Once we have these covariants identified, we
19 were then able to proceed with the primary and the
20 secondary analyses.

21 In the primary analysis, we look at the
22 time-to-event analysis and the progression of heart

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 failure and we were then able to look at the
2 longitudinal analysis of functional status over time.

3 Starting with progression of heart failure,
4 on the Y-axis is the result of reduction with CRT with
5 respect to the control group with no CRT.

6 We find that a 56 percent relative reduction
7 in mortality, 25 percent relative reduction in heart
8 failure hospitalization and a 13 percent relative
9 reduction in VT/VF events.

10 The composite endpoint was a 19 percent
11 reduction overall, which was not statistically
12 significant. However, every component of the index was
13 consistent and in a direction favorable to CRT with no
14 clinical evidence of harm.

15 When we considered the functional status
16 endpoints, after six months of CRT in this patient
17 population, we saw .7 milligram per kilogram per
18 minute improvement in peak VO_2 , which approached but
19 did not achieve statistical significance.

20 With quality of life, we saw improvement in
21 both groups, though we did not detect any
22 statistically significant difference in quality of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 life, between CRT and no CRT.

2 Turning now to the ancillary analysis, for
3 a six-minute walk, we saw a 22-meter improvement in
4 the CRT group after six months, compared to the
5 control group. Similar to peak VO_2 , it approached,
6 but did not achieve statistical significance.

7 For the V_E/VCO_2 slope, again, after six
8 months of CRT, we did not detect any differences.

9 The final ancillary analysis was that of New
10 York Heart Class. After six months we see that 68
11 percent of the patients are New York Heart Class I or
12 II, compared to 81 percent of patients after six
13 months of CRT. Again, this approached but did not
14 achieve statistical significance.

15 So to summarize our CRT effectiveness,
16 statistical significance was not reached for the
17 primary or the secondary endpoints. However, from the
18 event analysis, we do see a positive directional
19 effect of CRT upon all individual components of the
20 index -- mortality, heart failure hospitalization and
21 VT/VF events.

22 When we consider functional status, we also

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 see a modest trend toward clinical improvement with
2 CRT in the indices of peak VO₂, six-minute walk and
3 New York Heart Class.

4 So at the conclusion of the study, we sort
5 of stepped back and reflected on what we had learned
6 from the study and what we could determine.

7 First of all, let's consider what we've
8 learned about the EasyTrak lead and the CRTD EasyTrak
9 system.

10 We know the system is safe. The lead is
11 safe; the combination of lead with the conventional
12 ICD is also safe.

13 We also know the devices perform as
14 designed, that the EasyTrak lead can be placed with
15 high confidence in a decent amount of time, and
16 electrical performance is stable as well. We also
17 know that, when used as a system, that IC performance
18 is not compromised.

19 In terms of CT, when you consider safety, we
20 don't see any evidence that CRT is associated with
21 clinical harm in this patient population. In terms of
22 effectiveness of CRT, we see a positive directional

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 effect on the components of heart failure progression
2 as well as trends to improvement in functional status.

3 We also see a high degree of physician and
4 patient preference for the therapy. When we polled
5 patients in the last follow-up conducted, 97 percent
6 of the patients were programmed to CRTD at the last
7 follow-up visit.

8 We also have to consider the results of the
9 covariate analysis. The covariates provided by the
10 Heart Failure Events Committee were able to identify
11 for us a patient population with advanced heart
12 failure. After looking at this patient population,
13 that gives us encouragement to continue further looks.

14 The greatest improvements in peak VO_2 and
15 quality of life were associated with the severity of
16 baseline heart failure. That is, the more sick the
17 patient was, the more likely you are to see
18 improvement.

19 We also see emerging data from other studies
20 as well. Published studies were conducted while the
21 study was in progress, and the published evidence of
22 CRT in a similar patient population also showed

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 results.

2 You can look at the French Pilot Study,
3 which looked at Heart Class III/IV patients; the
4 InSync Study, which looked at Heart Class III/IV;
5 MUSTIC, which enrolled Class III; and PATH-CHF, which
6 included patients with New York Heart Class III/IV
7 heart failure. So positive results were reported in
8 these other studies with a similar patient population.

9 This is what we know, but then we have to
10 stop and consider what we don't know, and we have to
11 walk through some questions. We first of all have to
12 consider the risk/benefit ratio of CRT in the study
13 populace with New York Heart Class III/IV heart
14 failure.

15 First of all, we have to ask ourselves in
16 terms of risk: Is it possible that CRTD in a
17 specialized patient population could cause harm in a
18 patient with advanced heart failure? We don't know if
19 the lead or system safety looks different for these
20 patients than for the general population. We also
21 don't know the effect of heart failure progression in
22 these patients as well.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 On the flip side, consider benefit. Is it
2 possible these patients may have a greater benefit in
3 terms of heart failure progression and in terms of a
4 functional status. We didn't know if it is possible
5 that some or all functional status variables may show
6 improvement.

7 And, finally, what is the magnitude of
8 improvement, compared with other heart failure studies
9 and that of other reported CRTD studies.

10 With that in mind, we decided to proceed
11 with further analysis. We wanted to come to the
12 impact of CRT on a subgroup of patients with advanced
13 heart failure. We sought the Council for Independent
14 Statisticians to advise the process, that first of
15 all, multiple hypothesis testing must be kept to a
16 minimum to avoid detecting any spurious results.

17 It was also decided to use the same five
18 clinical variables obtained from the independent
19 through -- excuse me -- an independent group of
20 physicians and use those to do a separate analysis;
21 and, finally, without any advance knowledge of what
22 the impact of its covariates are, that the covariates

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 should be discrete. They should also be stable. We
2 should consider covariates that are stable over time
3 and are not in flux while the patient is in the study.

4 Baseline characteristics should be similar
5 for the control versus CRT group and, importantly, p-
6 values should be interpreted with care from this type
7 of analysis.

8 Our observations should include a
9 satisfactory sample size, and also importantly, that
10 clinical merit is a reasonable consideration in
11 assessing these findings. Therefore, with this
12 plan and these ground rules in place, we then
13 proceeded to further analysis.

14 Something else you have to take into
15 consideration was to recognize that the New York Heart
16 Class can change, that patients were enrolled in a
17 class first, and then physicians had the opportunity
18 to adjust medications. This was necessary because the
19 VT/VF events were negative.

20 At the time of enrollment in the study,
21 patients were predominantly New York Heart Class III.
22 But before any intervention took place and while

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 physicians had the opportunity to follow their
2 patients we see a shift towards higher heart classes.
3 So a number of patients are now in New York Heart
4 Class I or II who weren't there at time of enrollment.

5 It's the patients who remain in Class III or
6 IV after the one-month waiting period that constitute
7 the core of the advances heart failure group.

8 So let's consider the covariate analysis.
9 Looking at our pre-specified endpoints of peak VO_2 and
10 quality of life, we found that New York Heart Class
11 III or IV patients had a statistically significant
12 relationship with CRT.

13 We also saw that there was decreasing left
14 ventricular ejection fraction and widening QRS. We
15 also saw interactions with peak VO_2 ; however, not with
16 quality of life.

17 When we consider etiology and bundle branch
18 morphology, we found no relationships between these
19 covariates and the outcome.

20 Therefore, it's the New York Heart Class
21 III/IV patients at the time of randomization who were
22 found to be the only covariate that had a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 statistically significant relationship with both
2 endpoints.

3 Furthermore, the subgroup of patients is
4 consistent with that of other published studies, such
5 as the French Pilot Study, InSync, MUSTIC. They all
6 are considering a very similar patient population.

7 Let's look at the results we get from this
8 analysis. Starting out with the sample size, roughly
9 half the patients who enrolled in the study were
10 advanced heart failure group and the sample size we
11 obtained is similar to that we originally estimated
12 for Phase 1 study, using parallel arms.

13 If we separate out these patients, and
14 reconsider the endpoints for lead effectiveness, lead
15 safety, system effectiveness and system safety, we
16 find that we still meet all the safety and
17 effectiveness endpoints for the subgroup.

18 So let's go back now and reconsider the
19 question of heart failure. In this patient
20 population, again, similar to that of the original
21 all-patient population, we see that everything is
22 moving in the right direction, it's consistent and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 favorable towards CRT with respect to no CRT, and an
2 overall 25 percent reduction in progression of heart
3 failure. Again, no clinical evidence that CRT is
4 harmful in this patient population.

5 We now walk through the functional status
6 endpoints. The patients with advance heart failure
7 after six months of CRT showed an impressive 2.1 mL
8 per kilogram improvement over the control group, which
9 pretty much stayed flat.

10 Quality of life in this patient population
11 improved by nearly 11 points over the control group
12 after six months of CRT.

13 For the six-minute walk distance, we found
14 that a patient with CRT had a 48-meter improvement
15 over those patients who were randomized to the control
16 group. The V_E/VCO_2 slope was also favorable towards
17 CRT with a 3.7-meter improvement.

18 We finally looked at the shift in New York
19 Heart Class. We find that after six months of no CRT,
20 42 percent patients are now in the New York Heart
21 Class I or II. After six months of CRT, 72 percent of
22 patients are now near New York Heart Class I or II.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The results that we achieved are consistent
2 with clinical expectations. We went to the literature
3 to try to determine what represents a clinically
4 meaningful change to V-tach. If we start with peak
5 VO_2 , that numbers between one and two mLs per kilogram
6 per minute have been cited as the clinically relevant
7 improvements for patient population, and we were able
8 to achieve a 2.1 mL per kilogram per minute
9 improvement.

10 The designers of the quality of life
11 questionnaire designed it so that a five-point
12 improvement would be clinically meaningful, and we
13 were able to achieve an 11-point improvement.

14 In studies involving the six-minute walk
15 test, it was true that you need to see an improvement
16 of at least 45 meters to be clinically meaningful and
17 in this patient population we saw a 14-meter
18 improvement.

19 In studies with a V_E/VCO_2 slope, grades
20 between minus three and minus 12 have been cited as
21 clinically important improvements, and here we see a
22 minus 3.7 unit improvement.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Based upon other published studies of heart
2 failure, based upon other CRT studies that have been
3 published, we are very consistent with what's already
4 there.

5 So, finally, to summarize the relevance of
6 the observed CRT effects in this patient population.
7 First of all, the data we get are concordant. We have
8 identified patients with advanced symptomatic heart
9 failure who've been found to benefit with CRT. The
10 data are also consistent.

11 First of all, when we look at all of the
12 clinical variables, peak VO_2 , quality of life, six-
13 minute walk, V_E/VCO_2 slope or New York Heart Class, we
14 see clinically meaningful changes in all of them, all
15 of them in the direction that favors CRT.

16 Furthermore, if you look at our sample size,
17 that's very close to the sample size needed to show
18 the original functional status improvements in Phase
19 1 of the study.

20 It's also important to notice that the
21 magnitude of changes we see are indeed clinically
22 relevant. They match that of other heart failure

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 studies and are consistent with external studies which
2 also consider CRT.

3 Finally, there's the effect of the meta-
4 analysis. If we look at the individual components, we
5 saw a positive directional effect of CRT upon all the
6 individual components, again, with no clinical
7 evidence of harm.

8 So, to summarize, for the Ventak CHF/Contak
9 CD biventricular study, we believe we've provided
10 reasonable assurance that the Contak CD EasyTrak
11 system is safe and effective in the indicated patient
12 population.

13 Let's first consider safety. That we either
14 look at the lead or the system or CRT, it was found to
15 be safe for the entire patient study and it was found
16 to be safe in the advanced heart rate group.

17 This group is going to need to find its
18 patients with New York Heart Class III/IV, while on
19 heart failure drug therapy, left ventricular
20 dysfunction, defined as a LVEF less than 35 percent,
21 of a wide QRS, at least 120 milliseconds, and who are
22 also indicated for an ICD. Physicians were to take in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 all patients studied, as well as the advanced heart
2 failure group.

3 We consider effectiveness. But the
4 effectiveness of the devices themselves the EasyTrak
5 lead and the system are effective for the entire
6 patient population and remain effective in the
7 advanced heart failure group.

8 Finally, in terms of CRT and its
9 effectiveness, that we've demonstrated clinical
10 assurance that CRT is effective in a patient
11 population with advanced heart failure.

12 The clinical data before you supports the
13 proposed labeling that we seek. With that I'll now
14 turn the floor over to Dr. Michael Higgenbotham, who
15 will comment on the clinical relevance.

16 DR. HIGGENBOTHAM: Thanks, Pat. Panel
17 members, ladies and gentlemen. I must announce that I
18 have operated as a core laboratory for cardio-
19 pulmonary exercise testing in the Contak CD study and
20 have no other interest in Guidant.

21 When you look at the functional status of
22 heart failure patients, we sort of look at the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 functional incapacity as leading to several different
2 problems in patients with heart failure.

3 The primary problem, of course, is that
4 these patients are unable to achieve certain levels of
5 peak exercise and certain levels of sustainable
6 exercise.

7 But there are some secondary problems that
8 occur in a functional incapacity as well. They, of
9 course, are the symptoms, the unpleasant symptoms of
10 shortness of breath, and fatigue, and anxiety that
11 accompany attempts to scope certain levels of physical
12 activity.

13 And last but not least are the impacts that
14 the functional incapacity has on interactions with
15 other people, that lead to the very important
16 financial and social impacts on the patient's quality
17 of life.

18 So any worthy assessment of functional
19 capacity in heart failure patients, of course, has to
20 embrace those three sort of domains. The primary
21 problem of functional incapacity and the consequences
22 of it.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The Contak CD study achieved that. It
2 achieved an appropriate selection of influence, which
3 were these four measured here. The peak oxygen
4 uptake, six-minute walk, quality of life questionnaire
5 and the New York Heart Classification.

6 The peak oxygen uptake, of course, defined
7 the level of peak activity this certain individual
8 could obtain. The six-minute walk described another
9 element of exercise, which is a sustained ability of
10 exercise over a reasonable period of time.

11 And on the quality of life questionnaire
12 that was used, assessed the patient's impression of
13 the types of things that could be comfortably achieved
14 in the physical domain in this questionnaire, and also
15 the affect that that had on the patient's interaction
16 with other individuals, the psychosocial consequences
17 of functional incapacity and improvements in
18 functional capacity.

19 Finally, NYHA class, which is the other
20 important component of this sort of collage of this
21 collection of functional evaluations looked at the
22 patient's impression of how his incapacity affected

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 his day-to-day activities.

2 Now with those four, although they're all
3 essential, peak oxygen uptake is somewhat the most
4 robust. It's the most reliable measure of functional
5 status because it's the only one of all of those
6 measurements that objectively measures cardiac
7 reserve.

8 Peak oxygen uptake is a pretty good non-
9 invasive estimate, in fact, of exercise cardiac
10 output. And I know we like oxygen uptake because it's
11 independent of the protocol we used and of methods.

12 It doesn't matter whether there are minor
13 departures from the protocol. It doesn't matter what
14 instrument is used to measure exercise tolerance, and
15 we get pretty much the same answer when we used
16 maximal oxygen uptake.

17 And to those of us interested in exercise
18 physiology, it also gives us a common language, a
19 common currency with which to communicate.

20 An exercise time of five minutes means
21 nothing to us. An exercise oxygen uptake of eight ml
22 per kilograms per minute means everything to us.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Maximal oxygen uptake is important also
2 because it's reproducible without a great amount of
3 vulnerability to placebo or training effects.

4 This is -- it distinguishes maximal oxygen
5 uptake from measurements like maximal workload,
6 maximal exercise time, which are tremendously
7 susceptible to differences in motivation, and
8 differences in mechanical efficiency, which lead to
9 progressive increase in exercise time, as we have well
10 learned from a multitude of studies looking at
11 pharmacologic interventions.

12 Finally, peak oxygen uptake is the endpoint
13 that doesn't have to apologize to anybody. It is a
14 primary measure of quality of life and an objective
15 one that need not necessarily correlate with anything
16 else.

17 It's very much the gold standard of the
18 objective element of quality of life. One of the
19 objections, of course, to maximum oxygen consumption
20 is that day-to-day life is not a maximal event.

21 Surely, maximal exercise tolerance is not
22 the determinative of a patient's ability -- a person's

1 ability to perform day-to-day activity.

2 That's true, except when you're a heart
3 failure patient. It's true that in normal individuals
4 a difference of 20 percent, even 50 percent of maximal
5 oxygen uptake doesn't alter at all the extent to which
6 you can carry out day-to-day activities.

7 But this illustration shows very well that
8 the opposite prevails in patients who are impaired. If
9 we take 20 cc's per kilogram per minute broadly as the
10 low-end of normal, we see the impact and the
11 relationship.

12 And this is a figure that I modified after
13 Norman Jones' illustration. The greater the maximal
14 oxygen uptake is reduced, the more impact it has on
15 day-to-day activities.

16 In fact, it's not true that day-to-day life
17 for an impaired heart failure patient doesn't get into
18 the maximal oxygen uptake domain.

19 You can see that the sicker the patients are
20 toward the left hand side of this curve, the steeper
21 the relationship is.

22 And it makes sense that smaller increments

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 or decrements in physical capacity should have more
2 profound influences on day-to-day activity the more
3 people are impaired.

4 The second point made by this illustration
5 is the two cc's of oxygen uptake is a lot. If you
6 look at the difference between eight cc's per kilogram
7 per minute, ten and 12, you're looking at completely
8 different situations in terms of the independence of
9 the patient's life.

10 And two cc's makes a difference between a
11 patient that's stuck at home and one that can be taken
12 out to see friends or to go to the mall with some
13 assistance from friends or relatives.

14 Two cc's between the maximal oxygen uptake
15 of ten and 12 gets them on the phone to their friends
16 and relatives saying that they'll be going out by
17 themselves. It's sort of the mark of independence when
18 you get up to that maximal oxygen uptake.

19 They are the direct associations that shows
20 you the magnitude that oxygen uptake has on the day-
21 to-day quality of living.

22 Not only in the Contak CD study, where there

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 were these large changes, two cc's per kilogram per
2 minute seen in the advanced heart failure group, but
3 everything was concordant.

4 Not only as Pat showed you earlier, there
5 were significant changes and concordant changes in the
6 quality of life measurement, the six-minute walk
7 measurement and this V_E/VCO_2 ratio. Increments seem
8 to fall into place. Concordant moving in the right
9 direction, but also according to the right quantity.

10 It's a little hard for you to think about a
11 quality of life measurement of 10.9, but just consider
12 for a moment what the physical equivalent of a 48
13 meter change is in a six-minute walk.

14 If you have a six-minute race between two
15 heart failure patients and up to six minutes, one ends
16 up 48 meters ahead of the next, you don't have to
17 correlate that with very much to understand that that
18 is a profound change in physical performance.

19 The V_E/VCO_2 ratio is a metabolical
20 surrogate. It's a metabolic measurement that measures
21 the excessive ventilation that occurs in heart failure
22 patients out of proportion to the primary driver of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 ventilation during exercise, namely, CO_2 .

2 And what it describes is an increase in the
3 dead space, a mismatch between ventilation and
4 perfusion of the lungs that happens because of the
5 inadequate cardiac output that you're seeing in heart
6 failure. So it measures the patient's ability to
7 distribute blood to this particular organ.

8 And the problem with V_E/V_{CO_2} , of course, is
9 that it's a surrogate measure. The patient doesn't say
10 I feel worse today because I my V_E/V_{CO_2} slope is a
11 little bit higher. And that is a little bit of a
12 problem in interpretation.

13 But there are two major advantages that make
14 us keep on wanting to measure these measurements.
15 Number one, they're mechanistic. That's so important
16 for us in exercise physiology to know that the changes
17 in exercise tolerance that we see have a mechanistic
18 basis. It's something that you can't shake.

19 Now the second thing, talking of lack of
20 shakeability, is the V_E/V_{CO_2} ratios. Very reproducible
21 and it's not dependent on motivation. Most of the
22 data that we see in this V_E/V_{CO_2} slope is from

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 submaximal domain of exercise, something you can't
2 affect.

3 So it adds very much to the robustness of
4 our confidence that in this subgroup of patients,
5 something really was going on that was physiological.

6 Another thing I'd like to emphasize is that
7 we hardly every see this kind of concordance or
8 magnitude of exercise responses in pharmacologic
9 studies.

10 I have personally not seen data over a six-
11 month period where pharmacologic intervention safely
12 improves exercise tolerance and gives such a beautiful
13 concordance in all of the estimates of functional
14 capacity.

15 It's not unique though because in other CRT
16 studies, interestingly, show very similar data. I
17 show here the three controlled trials that have been
18 done with CRT and published fairly recently.

19 And without dwelling too much on statistical
20 significance, because occasional -- in this grid have
21 not achieved .05 statistical level of significance.
22 But all of them have come very close. Most of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 them have been significant and the qualitative and
2 quantitative concordance in those three studies is
3 remarkable.

4 When you think of the spotty or completely
5 absent improvement in exercise tolerance in
6 pharmacologic studies, three out of three isn't bad.

7 So I conclude that each of the endpoints
8 selected for this study were good ones. They were
9 totally appropriate. None of them was redundant. They
10 weren't repetitive. They looked at different elements
11 of functional incapacity.

12 They were complimentary, in that every one
13 of them, not only reinforced the validity of the other
14 in this subgroup, but added new information as to the
15 clinical relevance of the findings.

16 And, of course, the important thing for us
17 to consider is that CRT seemed to benefit this whole
18 profile of objective and subjective measurements.

19 The changes that we were seeing in the
20 functional status were positive, very internally
21 consistent, implying that there is reality in these
22 measurements.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The other thing that tends to reinforce my
2 confidence in these objective and subjective
3 measurements was that they compared very favorably
4 with those observed in other heart failure studies,
5 favorably in terms of pharmacologic studies, but very
6 concordant with what was seen, what has been seen in
7 other mechanical intervention studies.

8 We probably, more importantly, or as
9 importantly to these findings being real, was that
10 they were clinically important. The changes in
11 magnitude of these meant that the types of clinical
12 changes that you see, some metabolic measurements,
13 some easier, like the six-minute walk test to
14 interpret, had real clinical meaning for these heart
15 failure patients.

16 So I conclude that this probably represents
17 some sort of break through in this long struggle that
18 we had to identify an intervention in heart failure.
19 But over a very prolonged period of time, and it seems
20 to get better as time goes on, we see a safe and
21 effective improvement in functional class for these
22 heart failure patients, and that, I think is a major

1 advance in heart failure therapy. Thanks a lot.

2 DR. SWAIN: Mr. DeVries. Does the sponsor
3 have any further comments to make at this time?
4 Great. Thank you very much and thank you for staying
5 on time. Excellent presentations. The next will be the
6 FDA presentation, Dr. Barold.

7 DR. BAROLD: Good morning. This PNA was
8 submitted in what we call a modular form, meaning that
9 the manufacturing information, the device description
10 and pre-clinical and non-clinical laboratories studies
11 were submitted prior to the clinical data, were avidly
12 evaluated by the FDA and these modules have been
13 subsequently closed.

14 Today we'll be presenting the information
15 from the clinical module, which is the last module to
16 be evaluated. Next slide, please.

17 The review team at the FDA was quite
18 extensive, and I'd like to express my appreciation for
19 all of the support that they helped to put this
20 presentation together. Today you will be hearing from
21 myself and also from Dr. Gerry Gray in regards to some
22 of the statistical analysis. Next slide, please.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The preclinical testing was by the sponsor
2 and I would just again add that the generator of these
3 tests results have all been evaluated by the FDA and
4 they met the appropriate standards for testing. Next
5 slide, please.

6 Just as a basic device description to remind
7 you, this is a full functional ICD and dual chamber
8 pacemaker. It has a unique feature of biventricular
9 pacing capabilities. This is achieved by tying both
10 the right ventricular and left ventricular leads
11 together and they receive both simultaneous sensing
12 and pacing capabilities.

13 The EasyTrak left ventricular lead, which
14 will also be evaluated today, was placed into the
15 coronary vena system. Next slide, please.

16 The sponsor went through a complete study
17 method, the phase one and phase two. Today he will be
18 presenting the data from the phase two part of the
19 study and just to remind you very briefly, that every
20 patient received an implant.

21 These were patients that were indicated for
22 an ICD, so every single patient got the same implant.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Then after a 30-day waiting period, the
2 patients were then randomized to either the cardiac
3 resynchronization therapy or pacing on or pacing off.

4 And then after the six-month point, the
5 investigators were allowed to turn the cardiac
6 resynchronization therapy on. And as you can see, the
7 majority of the investigators did. Next slide,
8 please.

9 I'd like to read you the indications for use
10 of this device as defined by the sponsor. They are in
11 patients who have advanced symptomatic heart failure,
12 defined as New York Heart Association Class III and
13 four, including left ventricular dysfunction in an
14 ejection fraction less than or equal to 35 percent, a
15 wide QRS complex defined as a QRS greater than or
16 equal to 120 milliseconds while on heart failure drug
17 therapy and patients who are at high risk of sudden
18 cardiac death due to ventricular arrhythmias.

19 The primary study endpoint to this study
20 were to show a slowing of the progression of heart
21 failure as defined as a composite of all-cause death,
22 heart failure related hospitalizations and ventricular

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 tachycardia and ventricular fibrillation requiring
2 device therapy.

3 A modified endpoint, as was explained by the
4 sponsor, also included adverse events for heart
5 failure, and this would typically be patients that
6 presented to the emergency room for diuretic therapy,
7 but were not hospitalized.

8 The study was powered to detect a 25 percent
9 reduction in the event rates, and I'd like to remind
10 you that the control event rates were assumed to be a
11 15 percent death rate, a 30 percent rate for
12 hospitalization for congestive heart failure and a 20
13 percent rate of ventricular tachycardia ventricular
14 fibrillation.

15 These numbers were obtained from the Precise
16 Study, which was a randomized placebo controlled study
17 of carvedilol. Based on this, the sponsor calculated
18 a sample size of 308 patients. Next slide, please.

19 Additional study endpoints included an
20 improvement in the functional status, as measured by
21 the peak VO_2 , the V_E/VCO_2 slope and the six-minute hall
22 walk.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 They also measured the quality of life using
2 the Minnesota Living with Heart Failure quality of
3 life questionnaire, did a change in New York Heart
4 Association function class, the ATP conversion
5 efficacy, VT detection time, severe device related
6 adverse events and operative mortality. Next slide,
7 please.

8 They also looked at the appropriate study
9 endpoints to evaluate lead efficacy and safety. Next
10 slide, please.

11 The inclusion criteria the sponsor did talk
12 about. I would just like to remind you that these are
13 patients that are indicated for an ICD and, in
14 addition, they had to have symptomatic heart failure,
15 despite what was defined as optimal drug therapy,
16 although not a specific New York Heart Association
17 class was required to be included in this study.

18 The patients also had to have an ejection
19 fraction less than 35 percent and a QRS duration
20 greater than or equal to 120 milliseconds. Next slide,
21 please.

22 The inclusion criteria are all listed here

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 and the major exclusion criteria were that they were
2 not allowed to have a general indication for permanent
3 pacing. Next slide, please.

4 The analysis performed, which our
5 statistician will detail later, the primary endpoint
6 used a rank-based method in survival analysis. The
7 secondary endpoints using analysis of covariants. Next
8 slide, please.

9 Patient accountability. There were 581
10 patients enrolled in this study. 501 were actually
11 implanted and 490 randomized. There were 248 patients
12 in the treatment group and 253 in the control.

13 In the advanced heart failure subgroup,
14 which we will discuss later, there are 120 patients in
15 the therapy group and 116 in the control group. Next
16 slide, please.

17 These are the baseline characteristics of
18 all the patients which the sponsor has detailed. I
19 would just like to point out the New York Heart
20 Association class at the time of implant, or at the
21 time of enrollment, in which approximately one third
22 of the patients were in Class II. There were no major

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 differences between the therapy on or off groups.
2 Next slide, please.

3 At this point, I'd like to mention the
4 change in New York Heart Association class from
5 baseline, when the device was implanted, to 30 days
6 later, when the patient was actually randomized in the
7 study.

8 As you can see, there was a shift in the New
9 York Heart Association class from the time of implant
10 to randomization. This shift, however, was fairly
11 equal between the two treatment groups.

12 And I would also like to point out that the
13 advanced heart failure subgroup, which will be
14 discussed later, consists of this group of patients
15 here that are in Class III and IV at the time of
16 randomization. Next slide, please.

17 These just outline the baseline
18 characteristics of this advanced heart failure
19 subgroup. They are very similar to the patient's
20 characteristics associated with all patients, with the
21 exception of the New York Heart Association class.

22 Now these are the baseline characteristics

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 at the time of implant. And as you can see, a certain
2 percentage of those patients were in Class II heart
3 failure and then advanced to Class III heart failure
4 at the 30-day point. Next slide, please.

5 At this point, I'd like to turn the
6 presentation over to our statistician, Dr. Gerry Gray.

7 DR. GRAY: Good morning. My name is Gerry
8 Gray. I was the statistical reviewer for this
9 submission. I just want to talk a little bit about
10 some of the statistical issues regarding subgroup
11 analyses and the AHF subgroup

12 There were, actually, two different subgroup
13 analyses that have been submitted for this PMA. The
14 first one, on the original round, was a subgroup
15 called non-right bundle branch block, NYHA III/IV at
16 enrollment. There are 290 patients in that group.

17 In the next round we have the subgroup that
18 we're going to talk about now, the advanced heart
19 failure subgroup, which consisted of patients who are
20 NYHA class III/IV at randomization. And, again, that
21 was one month after enrollment and implantation.
22 There were 165 patients in common between those two

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 groups. Next slide.

2 There were five covariants that were defined
3 or selected by the Heart Failure Events Committee for
4 the analyses. NYHA class, bundle branch block
5 morphology, ischemic or not, QRS width, left
6 ventricular ejection fraction. And the NYHA class is
7 the one that's being used to define the AHF subgroup
8 here.

9 Again, I want you to recall that there was
10 a one-month waiting -- a one-month transition period
11 or waiting period between the time of implantation and
12 the time of randomization. Next slide.

13 Now this graphic shows you, it traces out
14 the NYHA class of all the patients in the study from
15 the time of enrollment to the time of randomization
16 one month later, and then at the six month follow up
17 time.

18 And what this is -- there's a lot of stuff
19 in here. But what you need to note mostly are these
20 highlighted bars right there and right there that show
21 you the amount of switching that went on between the
22 enrollment time, when patients were implanted with a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 device, and the randomization time, when the CRT
2 therapy was turned on or not. Next slide, please.

3 The reason we look at that is because we
4 want the covariants -- if you're doing a statistical
5 analysis using covariants, you don't want them to be
6 dependent on the treatment.

7 In other words, they should be something
8 that was either measured before the treatment began or
9 they are something that can't possibly be affected by
10 the treatment, like age or gender.

11 And in this case we're talking about using
12 NYHA class one month after implantation. And it
13 certainly was before the CRT therapy was either turned
14 on or not, but when you look at that graph previously,
15 you wonder was that affected by the implantation.

16 Is the implantation in the following of the
17 patients doing something to them that is actually --
18 might be part of treatment. I don't know.

19 A lot of that could just be due to something
20 called regression to the mean. In other words, when
21 patients were enrolled they were sicker and so one
22 month later, they happen to -- a large portion or some

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 portion of them happen to get better.

2 And the problem here is that if the
3 covariants are effected by treatment, then it's very
4 difficult to interpret the results, because if you use
5 a covariant that is affected by the treatment and
6 adjust for it, then you're sort of adjusting for
7 treatment to some -- it's very difficult to interpret
8 what comes out in the end. Next slide, please.

9 Okay. Leaving all that aside for now, for
10 the NYHA subgroup, if you're going to do -- use a
11 subgroup analysis to make a judgement, there are a lot
12 of caveats that you have to keep in mind when you're
13 interpreting those results.

14 And I think the major one is that in general
15 there's a tendency to over-interpret the significance
16 of what you see.

17 And the main reason for that is because when
18 you analyze multiple subgroups, when you look through
19 and you use NYHA class or QRS width, or bundle branch
20 block morphology, or whatever, to define a subgroup,
21 the chance of finding one that is somehow, and this is
22 in quotes, "significant" is very high.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 And the P values -- statistically, the P
2 values you get out from that are very difficult to
3 interpret.

4 In this case we have five endpoints.
5 Treating all the endpoints the same, we have five
6 endpoints and we have five different covariants that
7 you can use.

8 And you can use those combinations to create
9 many, many different subgroups. So from a statistical
10 point of view, you need to proceed with some caution
11 here. Next slide, please.

12 Now from a statistical point of view now is
13 there a justification for looking at this AHS subgroup
14 by itself? And what I have here are five different
15 criteria that you might use to judge a subgroup
16 analysis.

17 First of all, is it prospectively defined?
18 Was there enough information about that subgroup
19 before the study even began that we knew in advance we
20 were going to look at that subgroup? And the answer
21 here is no.

22 NYHA class was defined as a covariant, but

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 it was not seen as any more important than any of the
2 others, and it wasn't used to define any subgroup.

3 Secondly, is there a plausible biological
4 explanation for the results that you see when you
5 examine that subgroup. And in this case I think the
6 answer is yes, based on the other information.

7 Is the subgroup analysis -- is it an
8 analysis of the primary endpoint. Again, it was
9 prospectively defined. In this case, I think no, based
10 on the endpoints that we actually see that were
11 significant.

12 Is there a treatment effect in the overall
13 analysis? In other words, is there an effect of the
14 treatment in the overall analysis that you see that
15 something going on and you decide we're going to look
16 and try to understand where this is happening, and the
17 answer in this case is no.

18 And, finally, is there some interaction of
19 the treatment with the variable that we've used to
20 define the subgroup that would lead us to use that to
21 focus in on that subgroup? And, again, the answer in
22 this case is no.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Now there's another bullet that got left off
2 of this slide, which is there some independent outside
3 evidence that would lead you to this subgroup? And
4 that's not mine -- I won't comment on that. Next
5 slide, please.

6 So leaving aside now the question of whether
7 or not we should have examined this subgroup or not,
8 now that we are looking at the AHF subgroup, how are
9 we going to judge the statistical significance of what
10 we see?

11 And the question here is really what is a
12 significant P value? And the problem is that
13 statistically what you would consider a significant P
14 value gets smaller with each subgroup considered and
15 with each analysis that you do. It's really an
16 adjustment for a multiplicity effect.

17 In this case, because we're looking at --
18 for the AHF subgroup, we're focusing in on the
19 secondary endpoints, we really didn't define
20 prospectively the analysis we were going to use. That
21 wasn't agreed on in the IDE.

22 And in this case there are a bunch of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 different analyses that you might have used that were
2 all probably reasonable. You might be doing a single
3 time point, versus some repeated measures, et cetera,
4 et cetera. There's a long list of things that you
5 might have done.

6 And probably even more important than that
7 is a question of how many subgroups you considered
8 before you arrived at the one that you're focusing in
9 on.

10 In this case, I need to point out that
11 exploratory analyses count. In other words, if you get
12 your data and do -- make some graphs and do a few
13 simple tests for a bunch of different covariants and
14 based on that exclude some of them from consideration,
15 you've really, in effect, tested those.

16 And whether or not you do a formal test at
17 the end, doesn't really matter because you've looked
18 at the data and used them to make that judgement. How
19 much adjustment is necessary, because there's no way
20 to really go back and understand how many subgroup
21 analyses were even considered.

22 So just to summarize what I said, the use of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the NYHA class at the time of randomization one month
2 after implementation raises some questions about its
3 proper use as a covariant.

4 The justification from the statistical point
5 of view, the justification for considering the AHA
6 subgroup based on the data in this study is weak.

7 And finally, if you do go ahead and look at
8 that subgroup, there really is a very difficult
9 problem of how to interpret the P values that you get
10 out to judge what is statistically significant or not.

11 So with that, I'm going to turn it back over
12 to Helen to continue the presentation.

13 DR. SWAIN: Let's turn now to the actual
14 data. This is the composite endpoint for all patients
15 in this study. You can see it's broken down into death
16 from any cause, hospitalization, the adverse events
17 for heart failure and recurrent VTVF.

18 Now, again, the original endpoint did not
19 include the adverse event for heart failure, but I
20 will be talking about the modified composite endpoint.

21 You can see that there is a 23 percent
22 reduction in the modified composite endpoint with a P

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 value of 0.12 and the exact numbers are listed here.

2 I just would like to remind you that the
3 predicted control rates were quite different than what
4 was actually observed. The observed rates were much
5 lower than what was actually predicted. Next slide,
6 please.

7 Here are -- it's just a graphical
8 representation of that data that was put together by
9 Dr. Gray for each of the four parts of the composite
10 endpoint. Next slide, please.

11 And here's the composite endpoint for the
12 advanced heart failure subgroup. You can see in this
13 case that the event rates are a little bit closer in
14 the control group to the predicted rates but, again,
15 don't quite match that. There is a 29 percent overall
16 reduction in -- seen in the advanced heart failure
17 subgroup with a P value of 0.11. Next slide, please.

18 I'd like to examine the mortality that was
19 seen in this study. In both groups there was 11
20 patients who died in the treatment group of all cause
21 mortality, broken down into seven patients with a
22 cardiac cause and four patients who died of pump

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 failure.

2 In the no pacing group, there was 16
3 patients who died of all cause mortality, four
4 patients died of cardiac causes and nine patients died
5 of pump failure. These are the P values for the two
6 groups. Next slide, please.

7 From a functional standpoint, here's the
8 data for the peak VO_2 and the all patient group for
9 the three month and six month. And you can see the
10 individual numbers listed here with the associated P -
11 values. Next slide, please.

12 This is a graphical representation of this
13 data, which shows the individual variability in this
14 data and the control -- and the treatment group at
15 baseline three months and six months. There are more
16 patients, obviously, in the baseline to three months
17 because those were the phase one patients.

18 And you can see the colored lines are the
19 mean values for this. But you can see there's a
20 tremendous amount of variability in the data. Next
21 slide, please.

22 If you can, again, take out the advanced

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 heart failure subgroup and you compare from baseline
2 to the six months, there was more of an improvement
3 seen in the treatment group over the non-treatment
4 group. Next slide, please.

5 As far as the V_E/VCO_2 slope, both groups
6 showed an improvement in the all patient group.
7 However, there was no significant difference between
8 the groups. Next slide, please.

9 And here is the -- again, the individual
10 data which, again, illustrates the high variability
11 seen in all patients. Next slide, please.

12 Here's the data for the change in quality of
13 life at three months and six months. You can see that
14 both groups saw an improvement, both treatment and
15 control, and you can see the incremental change at six
16 months between the treatment and control with a P
17 value of 0.44 Next slide, please.

18 Again, just to point out the incredible
19 variability seen in the all patient groups here, and
20 with the colored lines showing the means. Next slide,
21 please.

22 Here is the advanced heart failure subgroup

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 data. And you can see -- well, there's a little bit
2 of a problem with that P value there. You can see the
3 data at three months and six months, and then
4 incremental change seen between the two groups. Next
5 slide, please.

6 And, again, here's the data for the six-
7 minute hall walk at three and six months. You do see
8 an improvement in both groups at the three and six-
9 month points. The incremental change between the two
10 groups and the associated P value. Next slide, please.

11 And, again, the data shows a tremendous
12 amount of variability with the mean value shown in
13 color. Next slide, please. And here is the advanced
14 heart failure subgroup. You can see that the
15 difference at six months between the two groups is 48
16 meters with an associated P value.

17 Next slide, please. This slide shows the
18 change in New York Heart Association functional class
19 in the all patient. These are patients that changed -
20 - decreased by two or more, one, no change. And you
21 see there isn't an improvement in both the treatment
22 and the non-treatment group. Next slide, please.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 And here's that associated change in the
2 advance heart failure subgroup and you can see in the
3 treatment group the numbers are quite low for this
4 analysis, but there did appear to be more patients
5 improved one or two classes in the treatment group.
6 Next slide, please.

7 I'd also like to point out the ATP
8 conversion efficacy. Remember that this device's
9 primary effect is an implantable cardioverter
10 defibrillator, so we need to assure that it is able to
11 deliver that therapy.

12 They tested the conversion rate, or they
13 looked at the conversion rate in the EP lab and the
14 spontaneous conversion rate. I do think that the
15 spontaneous conversion rate is what we should be
16 looking at, and they showed an 88 percent success
17 rate, and these numbers were similar in their advanced
18 heart failure subgroup. Next slide, please.

19 Also very important to look at is the
20 ventricular fibrillation detection time. And there was
21 no significant increase in the detection time with the
22 addition of the left ventricular lead and, again, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 data was for the advanced heart failure subgroup. Next
2 slide.

3 Now, let's take a look at the severe device-
4 related adverse events and operative mortality
5 associated with this. The sponsor hypothesized this
6 rate to be approximately 20 percent, and the actual
7 rate was 1.2 percent.

8 During the study there were five generator
9 failures, four of which required new implants. One was
10 noted prior to implant.

11 The operative mortality was 3.4 percent for
12 the thoracotomy procedure and 2.0 percent for the
13 transthoracic procedure. In the peri-operative
14 mortality there were 12 patients that died in the
15 peri-operative period at a rate of 2.1 percent. Next
16 slide, please.

17 This slide just illustrates the
18 hospitalizations for heart failure in the all patients
19 and the advanced heart failure subgroup.

20 You can see that there were 48
21 hospitalizations for congestive heart failure in the
22 treatment group and 48 hospitalizations for heart

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 failure in the no treatment group.

2 And in the subgroup there were 36
3 hospitalizations for heart failure in the treatment
4 group and 32 in the no treatment group. Next slide,
5 please.

6 Let's take a look at the EasyTrak lead
7 safety. There were 72 or 13.9 percent of patients had
8 a lead-related adverse event, the most common of which
9 was left ventricular lead dislodgement and there were
10 29 of those seen, or 6.5 percent.

11 Some of the more serious complications that
12 we saw included five cases of coronary sinus
13 perforations and one guide wire fracture that was
14 subsequently removed by snare.

15 Sixty-nine patients could not have the Easy
16 Lead Trak placed and the majority of those were due to
17 problems located in cannulating the coronary sinus.
18 Next slide, please.

19 One of the more important adverse events
20 that we identified were incidents of coronary sinus
21 trauma. Remember, that's what differentiates this
22 lead implant from normal ICD's since you are putting

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 a lead in the coronary sinus. So it's important that
2 we look at that particular part.

3 And we asked the sponsor to provide us data
4 from all of their EasyTrak lead implants. They are
5 using this in a number of studies, so we felt it was
6 necessary to combine all of those leads.

7 So that's why you'll see in all of the
8 EasyTrak lead implants and the information we have is
9 1,374 lead implants. There were 39 cases of coronary
10 sinus trauma.

11 Out of those there were 20 dissections, 17
12 perforations, and two episodes of tamponade. 36 of
13 these, or 92 percent resolved without any
14 intervention.

15 Two of these caused subsequent death of the
16 patients and there were possibly another two that may
17 have been related to coronary sinus trauma. Next
18 slide, please.

19 The sponsor did a very nice job of going
20 through the EasyTrak lead results. But here's the
21 exact data. The adverse events associated with the
22 left ventricular lead was 10.9 percent and there was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 a procedurally related adverse event of 13.9.

2 You can see that the implant success rate
3 was 86.7 percent and the leads performed as the
4 sponsor expected them to do. Next slide, please.

5 So as far as the EasyTrak lead points, the
6 pacing thresholds stabilized after one month, which is
7 to be expected with these leads, in that they're
8 sensing an impedance endpoints. Next slide, please.

9 So in summary, the sponsor has met all the
10 preclinical and manufacturing requirements for this
11 device. They met their safety and lead performance
12 endpoints. They did not satisfy the effectiveness
13 endpoints when evaluating the all patients.

14 However, in the advance heart failure
15 subgroup, there does appear to be more improvement
16 with the cardiac resynchronization therapy in most of
17 the functional implants.

18 At this point, would you like to discuss the
19 PMA and then I'll read the questions, or would you
20 like me to read them now?

21 DR. SWAIN: Actually, why don't you just
22 quickly go over the questions to remind our panel

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 members that we need to have them address the
2 questions.

3 MS. BAROLD: Okay. For the record then, I
4 will go ahead and read the exact wording of the
5 questions.

6 Question one, in the clinical report,
7 section 4.3.1.4 identifies the adverse events,
8 complications and observations for the system as a
9 whole and each individual component, including the
10 EasyTrak lead system.

11 Question 1A, the rate of coronary sinus
12 trauma observed in this study with the EasyTrak lead
13 was three to four percent. Please discuss safety
14 issues associated with the implantation of a third
15 lead in the coronary venous system and comment on
16 whether the data in the PMA supports the safety of the
17 lead system for the proposed indication.

18 Additionally, please discuss the clinical
19 importance of the overall adverse events,
20 complications, observations and comment on whether the
21 data provides reasonable assurance of safety of this
22 device system.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Question 2. The primary effectiveness
2 endpoints in this study were composite of all cause
3 mortality, hospitalizations for heart failure and
4 episodes of ventricular tachycardia, ventricular
5 fibrillation requiring therapy. The secondary
6 endpoints for peak VO_2 , V_E/VCO_2 slope, six-minute hall
7 walk and quality of life questionnaire.

8 We ask you to please discuss the clinical
9 relevance of the effectiveness endpoints for the
10 patient population. The clinical study was designed
11 with six months of follow up. Please comment on
12 whether this point is adequate to provide a reasonable
13 estimate of device safety and effectiveness.

14 A subgroup analysis performed on those
15 patients with Class III and IV heart failure showed a
16 more favorable outcome in the secondary endpoints.
17 Please discuss whether the data in the PMA provides a
18 reasonable assurance of safety and effectiveness in
19 this group.

20 Question 3. The control group saw
21 improvements in their functional status, quality of
22 life, six-minute hall walk and New York Heart

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Association functional class.

2 Please comment on this improvement in the
3 control group as it relates to the improvement in the
4 treatment group and if this relationship changes with
5 the subgroup analysis of patients with advanced heart
6 failure.

7 Additionally, please comment on the clinical
8 relevance that this finding has on the observed
9 effectiveness of cardia resynchronization therapy.

10 Question 4. Please discuss whether the data
11 in the PMA provide reasonable assurance of
12 effectiveness for this device in the patient
13 population study?

14 Question 5. One aspect of the premarket
15 evaluation of a new market is the review of its
16 labelling. The labelling must indicate which patients
17 are appropriate for the treatment, identify potential
18 adverse events with the use of this device and explain
19 how the product should be used to maximize benefits
20 and minimize adverse events.

21 If you recommend approval of this PMA,
22 please address the following questions regarding the

1 product labelling, as found in section three. In the
2 indications portion of the labelling it states, "This
3 device is indicated for patients with advanced
4 symptomatic heart failure, as defined by New York
5 Heart Association Class III and IV, including left
6 ventricular dysfunction of wide QS complex while on
7 heart failure drug therapy and have current
8 indications for an ICD.

9 Based on the data provided, is this
10 indication supported by the data provided? Please
11 comment on whether the indication statement identifies
12 the appropriate patient population for treatment with
13 this device.

14 Also, please comment on the operator
15 instructions as to whether they adequately describe
16 how the device should be used to maximize the benefits
17 and minimize adverse events. Please comment and
18 provide any other recommendations or comments
19 regarding the labelling of this device.

20 Question 6. Please identify and discuss the
21 items that you believe should be contained in a
22 physician's training program for this device. For

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 example, please comment on whether training should be
2 required for proper placements of the EasyTrak lead
3 system.

4 Question 7. Based on the clinical data
5 provided in the panel packet, do you believe that
6 additional clinical follow up or post-market studies
7 are necessary to evaluate the long-term effects of
8 biventricular pacing on heart failure.

9 If so, how would you design a study,
10 including study design, duration, sample size, patient
11 characteristics, for example, is a Q restoration of
12 120 milliseconds long enough to suggest significant
13 ventricular desynchrony and what other measurements
14 could be substituted, and what other additional
15 potential endpoints should be looked at? That
16 concludes the FDA presentation.

17 DR. SWAIN: Thank you very much. I have to
18 comment that that was an excellent review by the FDA
19 reviewers and the work that was put in. And the panel
20 package is very well done. I still remember a decade
21 or ago when we were given, I think, five to six feet
22 of data, taller than me, to review. So I compliment

1 the FDA.

2 And what we're going to do is we're going to
3 start having questions from our primary reviewer, Dr.
4 Domanski, for approximately 15 minutes.

5 The way I run the panel is about 15 minutes
6 from the primary, ten minutes from each of the other
7 reviewers and then we just keep going in circles till
8 everybody has all of the questions asked and answered
9 that they wish.

10 So we'll start with Mike and then we'll take
11 a break when Mike's finished for 15 minutes, and then
12 reconvene. Mike.

13 DR. DOMANSKI: Well, I'm going to -- I do
14 want to also compliment the FDA group that put this
15 together. You know, you really did a beautiful job.
16 It's very nicely done and I think it tells a story
17 almost by itself that strikes me in looking at it as
18 pretty straightforward.

19 I think I have some real concerns about this
20 study, and I'm going to cut early to the chase on it.
21 And what I would like to is I am very concerned that
22 you have not demonstrated effectiveness of this device

1 with the studies you've presented.

2 But I don't want to be wrong about that,
3 because a tremendous amount of effort and money and so
4 forth goes into these activities, and it's important
5 that one be smart about putting on the market things
6 that are useful to the public.

7 But I do think that we -- I think that
8 within this study you need to demonstrate
9 effectiveness and safety, of course, and it has to be
10 within this study.

11 We've had a -- this brief preamble is worth
12 it to me anyway. I think that we've seen -- I've seen
13 over many years now on this panel devices sometimes
14 come to us that turned out to be quite useful as time
15 went by, but which were not presented in a way that
16 made that clear in the application.

17 And I'm concerned that that may be the case
18 here. I'm not convinced that it's useful, but I think
19 it may well be a useful maneuver.

20 So what I'd like to do is track through your
21 effectiveness endpoints, using the FDA's workup and
22 give you an opportunity to respond in a way that makes

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com