

sgg

AT

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

CIRCULATORY SYSTEM DEVICES PANEL

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

Monday, September 9, 2002

10:30 a.m.

Hilton Washington DC North
629 Perry Parkway
Gaithersburg, Maryland

PARTICIPANTS

Warren K. Laskey, M.D., Acting Chairperson

Elisa Harvey, DVM, Interim Executive Secretary
Geretta Wood, Executive Secretary

VOTING MEMBERS:

Salim Aziz, M.D.
Julie A. Freischlag, M.D.

CONSULTANTS:

Bruce A. Perler, M.D.
Anthony Comerota, M.D.
Christopher J. White, M.D.
Kenneth E. Najarian, M.D.
Anne C. Roberts, M.D.
Michael J. Pentecost, M.D.
Kent R. Bailey, Ph.D.
Anton Sidawy, M.D.
Ileana Pina, M.D.

CONSUMER REPRESENTATIVE:

Robert A. Dacey

INDUSTRY REPRESENTATIVE:

Andrew K. Balo

FOOD AND DRUG ADMINISTRATION:

Bram Zuckerman, M.D.
Doyle Gantt
Dorothy Abel
Paul Chandeysson, M.D.

C O N T E N T S

Call to Order	4
Open Public Hearing:	
Rodney A. White, M.D.	10
W.L. Gore & Associates, Inc. Presentation, P0200044, Excluder Bifurcated Endoprosthesis:	
Introduction, Jon Sininger	19
Product and Study Overview, David Williams	23
Abdominal Aortic Aneurysm Background, David C. Brewster, M.D., Harvard Medical School and Massachusetts General Hospital	28
Trial Design and Study Management, David C. Naftel, M.D., University of Alabama	37
Pivotal Study Clinical Results, Jon S. Matsumura, M.D., Northwestern University Medical School	41
FDA Presentation, A. Doyle Gantt	71
Open Committee Discussion	87
Open Public Hearing:	
Takao Ohki, M.D.	246
Mark Fellingner, M.D.	247
Roy K. Greenberg, M.D.	248
Rodney A. White, M.D.	250
Final Comments from the FDA, Dorothy Abel	253
Final Comments from the Sponsors	255
Recommendations and Vote	256

P R O C E E D I N G S

DR. LASKEY: It is time for us to come to order. The topic to be discussed today is the premarket application for the Gore bifurcated endoprosthesis, P020040. I would like to have the executive secretary read the conflict of interest statement now.

DR. HARVEY: The following announcement addresses conflict of interest issues associated with this meeting and is made part of the record to preclude even the appearance of an impropriety.

To determine if any conflict existed, the agency reviewed the submitted agenda for this meeting and all financial interests reported by the committee participants. The conflict of interest statutes prohibit special government employees from participating in matters that could affect their or their employers' financial interests. The agency has determined, however, that the participation of certain members and consultants, the need for whose services outweighs the potential conflict of interest involved, is in the best interest of the government.

Therefore, a waiver has been granted for Dr. Bruce Perler for his interest in a firm that

1 could potentially be affected by the panel's
2 recommendations. The waiver involves a grant to
3 his employer for a competitor study in which he is
4 not involved in data generation or analysis, and
5 for which funding is less than \$100,000 per year.
6 Copies of this waiver may be obtained from the
7 agency's Freedom of Information Office, Room 12A-15
8 of the Parklawn Building.

9 We would like to note for the record that
10 the agency took into consideration other matters
11 regarding Drs. Julie Freischlag, Kenneth Najarian,
12 Anne Roberts and Michael Pentecost. Each of these
13 panelists reported interest in firms at issue but
14 in matters that are not related to today's agenda.
15 The agency has determined, therefore, that they may
16 participate fully in all discussions.

17 The agency also would like to note that
18 due to the regulations governing covered
19 relationships, the panel chair, Dr. Cynthia Tracey
20 will not participate in today's deliberations.

21 In the event that the discussions involve
22 any other products or firms not already on the
23 agenda for which an FDA participant has a financial
24 interest, the participant should excuse him or
25 herself from such involvement and the exclusion

1 will be noted for the record.

2 With respect to all other participants, we
3 ask in the interest of fairness that all persons
4 making statements or presentations disclose any
5 current or previous financial involvement with any
6 firm whose products they may wish to comment upon.

7 DR. LASKEY: Thank you. I would like to
8 have the panel members now introduce themselves,
9 beginning to my right, please.

10 MR. BALO: Andy Balo, DexCom, Inc.,
11 industry rep.

12 DR. AZIZ: Salim Aziz, clinical associate
13 professor, University of Colorado.

14 DR. COMEROTA: Anthony Comerota, vascular
15 surgeon, director of the Jobst Vascular Center and
16 professor at University of Michigan, Ann Arbor.

17 DR. PENTECOST: Michael Pentecost,
18 professor and chairman of radiology at Georgetown.

19 DR. BAILEY: Kent Bailey, biostatistician,
20 Mayo Clinic.

21 MS. WOOD: Geretta Wood, executive
22 secretary.

23 DR. HARVEY: Elisa Harvey, interim
24 executive secretary for this panel meeting.

25 DR. LASKEY: Warren Laskey. I am an

1 interventional cardiologist from Baltimore.

2 DR. SIDAWY: Tony Sidawy. I am chief of
3 surgery at the VA Medical Center here, and
4 professor of surgery at George Washington and
5 Georgetown Universities.

6 DR. FREISCHLAG: Julie Freischlag. I am a
7 vascular surgeon and chief of vascular surgery at
8 UCLA.

9 DR. NAJARIAN: Kenneth Najarian,
10 interventional radiologist and professor of
11 radiology at the University of Vermont.

12 DR. ROBERTS: Anne Roberts, interventional
13 radiologist and professor and chief of vascular and
14 interventional radiology at UC San Diego.

15 DR. PERLER: Bruce Perler, chief of
16 vascular surgery at Johns Hopkins.

17 DR. WHITE: I am Chris White. I am an
18 interventional cardiologist, and I am from the
19 Ochsner Clinic in New Orleans.

20 MR. DACEY: Robert Dacey, consumer
21 representative, Boulder County, Colorado.

22 DR. ZUCKERMAN: Bram Zuckerman, director,
23 Division of Cardiovascular Devices, Food and Drug
24 Administration.

25 DR. LASKEY: Thank you, all. Elisa, would

1 you read the voting status please?

2 DR. HARVEY: Yes. This is an appointment
3 to temporary voting status. Pursuant to the
4 authority granted under the Medical Devices
5 Advisory Committee Charter, dated October 27, 1990
6 and as amended August 18, 1999, I appoint the
7 following individuals as voting members of the
8 Circulatory System Devices Panel for this meeting,
9 on September 9th, 2002: Anthony Comerota,
10 Christopher White, Kenneth Najarian, Anne Roberts,
11 Michael Pentecost, Bruce Perler, Kent Bailey and
12 Anton Sidawy.

13 For the record, these people are special
14 government employees and are consultants to this
15 panel under the Medical Devices Advisory Committee.
16 They have undergone the customary conflict of
17 interest review and have reviewed the material to
18 be considered at this meeting. In addition, I
19 appoint Dr. Warren Laskey to serve as panel chair
20 for the duration of this meeting.

21 It is signed by Dr. David Feigal, Director
22 for the Center of Devices and Radiological Health,
23 August 30th, 2002.

24 In addition, I have another voting status
25 to read: Pursuant to the authority granted under

1 the Medical Devices Advisory Committee Charter of
2 the Center for Devices and Radiological Health,
3 dated October 22, 1990 and as amended August 18th,
4 1999, I appoint the following individual as a
5 voting member of the Circulatory System Devices
6 Panel for the meeting on September 9th, 2002,
7 Ileana Pina, M.D.

8 For the record, Dr. Pina is a consultant
9 to the Cardiovascular and Renal Drugs Advisory
10 Committee of the Center for Drug Evaluation and
11 Research. She is a special government employee who
12 has undergone the customary conflict of interest
13 review and has reviewed the material to be
14 considered at this meeting. It is signed by
15 William Hubbard, Senior Associate Commissioner for
16 Policy and Planning, on behalf of Linda Skladany,
17 Senior Associate Commissioner for External
18 Relations, September 2, 2002.

19 DR. LASKEY: Thank you. At this point I
20 would like to open this portion of the meeting, the
21 public hearing, and to ask the audience if there is
22 anyone who wishes to address the panel on the day's
23 topic preferably. Dr. Rodney White had sent a
24 letter requesting time before the panel. Is Dr.
25 White in the audience?

1 **Open Public Hearing**

2 DR. RODNEY WHITE: Yes. My name is Rodney
3 White. I am a vascular surgeon from Harbor UCLA.
4 I am a member and chairman of the Lifeline Registry
5 Committee, which is the topic for today; secretary
6 of the Society for Vascular Surgery. I think my
7 greatest conflict is that I am a clinical vascular
8 surgeon who makes my living implanting these
9 devices and showing up at meetings like this, and I
10 have been the PI or co-investigator in many of the
11 clinical trials that are currently under
12 evaluation.

13 [Slide]

14 What I wanted to speak to you briefly
15 about this morning is the Lifeline Registry. This
16 is a project that was initiated back in 1998, prior
17 to the approval of any of the endoluminal graft
18 devices. At that particular time there was
19 interest by not only the manufacturers and the
20 clinical investigators but by the various agencies,
21 and particularly the FDA, to look at issues that
22 may be developed and related to endoluminal
23 vascular grafts.

24 [Slide]

25 The Lifeline Registry goals were to do two

1 things, to provide a longitudinal observational
2 database where endoluminal graft performance could
3 be defined and evaluated and then, secondarily,
4 over time, as will become apparent during today's
5 talk, that there are surveillance issues that need
6 to be addressed to follow these patients
7 appropriately, and the attempt, again, from this
8 multifacet aspect was to develop those as issues
9 arise and make clinical tools that could easily
10 follow these patients available.

11 [Slide]

12 The Lifeline Registry has other aspects to
13 it. The web page which is, again, supported by the
14 industrial partners and by the SVS and AVS, has not
15 only patient information but updates data on a six
16 monthly basis from the Registry. This is published
17 not only in the Journal of Vascular Surgery but is
18 also updated periodically so that the data is
19 available to everybody.

20 [Slide]

21 As I have mentioned before, the mission is
22 to provide longitudinal consecutive data, and with
23 reference to the panel meetings all of the
24 manufacturers that have received approval and the
25 data set that will be presented today will become

1 part of the Registry. There is an option for the
2 manufacturers to be able to use the Registry for
3 storage of their postmarket surveillance data, and
4 that makes then an easily available, and because of
5 the high level of compliance that goes with these
6 studies, reliable database to make these long-term
7 observations.

8 There are also initiatives to work with
9 the new VA cooperative study which was recently
10 approved, and with our Canadian collaborators who
11 have similar studies, so that the attempt would be
12 on a voluntary basis to have a very large registry
13 of data at a high compliance level to address
14 issues as they develop.

15 [Slide]

16 The key stakeholders then are the
17 clinicians in the societies, the Lifeline
18 Foundation which is the funding arm of the SVS that
19 our industrial advisory committee is made up of.
20 All of the major manufacturers are participating in
21 this, and the federal agencies, including NIH, FDA
22 and CMS, have ex officio seats so that, hopefully,
23 all of the relevant individual people are there.

24 [Slide]

25 With regard to surveillance, this is a big

1 clinical problem. We are relying on CT imaging.
2 There is a major amount of data that is required,
3 and the issues involved are at many levels, all the
4 way from patient compliance to storage and
5 efficient cost and effective relay of this to a
6 site that could be accepted.

7 [Slide]

8 I list here the major manufacturers who
9 have throughout this project supported this.
10 Again, this effort was initiated in 1998 prior to
11 any endoluminal graft being implanted. So, one of
12 my comments would be that from all respects this is
13 an effort by industry and by the agencies to be
14 very proactive in terms of following these
15 patients.

16 [Slide]

17 I show you an example of a patient who has
18 had a device for six years. The issues related to
19 imaging and what happens to that aneurysm over time
20 are particularly relevant, and it is a new scenario
21 we have not dealt with before. There are many
22 measurements involved in this; what happens to the
23 aneurysm and the fixation sites? Are there leaks?
24 Are the grafts patent? Things that we have never
25 dealt with before clinically.

1 [Slide]

2 Through an interactive system that is
3 developed to be able to collect this data from
4 initial paper forms and now through an electronic
5 format, the attempt is to make available not only
6 to the manufacturers and the agencies through the
7 PMA data sets that are stored in the Registry, but
8 also as a clinical tool that could be developed for
9 surveillance mechanisms.

10 [Slide]

11 I had mentioned to you that in the Journal
12 of Vascular Surgery there is a publication every
13 six months. This is just one of the tables from
14 the March issue of this last year, which shows that
15 there are 1600 patients now. We are also looking
16 at control patients that become available and as
17 the data set grows, this now is the highest
18 compliance and the largest volume of patients that
19 is available with follow-ups in the three- to-six
20 year range. So, this is becoming quite mature and
21 able to address many of the issues.

22 [Slide]

23 All the data is stored on a web site
24 through the New England Research Institute, which
25 is our administrative arm.

1 [Slide]

2 It has a secured site that enables us to
3 do that, and through a series of tables, and I
4 won't go through it, there is data on each patient
5 relevant to measurements.

6 [Slide]

7 And then corresponding imaging is stored
8 so that it is readily available, able to be
9 analyzed in retrospect.

10 [Slide]

11 To summarize this, what I would say is
12 that there are two papers I would refer you to, one
13 in the Journal of Vascular Surgery, one that
14 globally describes the Registry and how it
15 operates, and the second one, which is the first
16 data report that was published in June--this will
17 be the format again, every six months all the data
18 published on the web site, accessible and for
19 relevant questions from anyone that could be
20 addressed. These are unauthored publications in an
21 attempt to make this readily available.

22 [Slide]

23 If anyone is interested in contacting the
24 New England Research Institute, which is the
25 administrative arm, they can supply more

1 information. Thank you for the opportunity to
2 present this.

3 DR. LASKEY: Thank you. Are there any
4 questions from the panel members?

5 DR. COMEROTA: Rod, that is a nice review.
6 You mentioned that patients who were so treated
7 with endografts would become part of the Registry,
8 but then you also went on to say that it is
9 voluntary. Could you just clarify that please?

10 DR. RODNEY WHITE: Yes, there are actually
11 two parts to the Registry. The first part, which
12 we call Part A, is similar to what you are going to
13 hear today, the data set that is submitted for the
14 PMA submissions. For Medtronic and Guidant,
15 following their approval in 1999, they have
16 continued to update their data. So, that is the
17 Registry data set. After today, will be the
18 submission of the Gore data and with subsequent
19 instances we will, hopefully, be able to get all
20 the manufacturers' data so that there is voluntary
21 compliance of their submissions, although they
22 can--and if the way I say this isn't correct,
23 please, the agency will tell me, but they have
24 offered to the manufacturers that they can use the
25 Registry to store their data. When they store the

1 data they can download it by the automated system
2 and use that as part of their annual reports. In
3 that way, it is then made a very high compliance
4 level data set across the industry, available to
5 everybody to view the data, to be able to see how
6 that works.

7 So, it is voluntary and each manufacturer
8 can do this on their own but, again, I would
9 emphasize there has been an effort across the
10 industry to do this in collaboration before it even
11 became an issue.

12 Secondly, the clinical tool that is
13 available is available to individual practitioners
14 and that, obviously, would be voluntary but we
15 offer that. We ask them to consent their patients
16 according to the IRB regulations so that we can
17 follow those folks over time. There is also even a
18 new ability that has been worked out for
19 investigator IDEs to capture that data by the same
20 mechanism. So, it is voluntary but has turned out
21 to be a good repository and, hopefully, a way to
22 work issues through all the relevant parties to
23 solve any problems that come up.

24 DR. LASKEY: Thank you. This is a very
25 important area. Dr. Zuckerman?

1 DR. ZUCKERMAN: Yes, Dr. White has
2 indicated one possible mechanism or pathway by
3 which manufacturers, after device approval, have
4 periodically updated the agency with required data.
5 This isn't to say, however, that this is the only
6 way it can be done. From the agency's perspective,
7 with the other two mentioned manufacturers we were
8 interested in periodic update reports. What you
9 have heard here is one mechanism for generating
10 such data.

11 DR. LASKEY: I just have one final
12 question, if I might. Maybe I missed it on the
13 slide, but the support for this Registry derives
14 from?

15 DR. RODNEY WHITE: There are two sources
16 of funding. The first is an ongoing commitment
17 from the Society for Vascular Surgery and American
18 Society for Vascular Surgery, which is from the
19 Lifeline Foundation itself, the funding arm. The
20 major funding comes from what we call the
21 industrial advisory committee, which is constituted
22 of each of the major manufacturers that make these
23 devices. They all, to a company, have on an annual
24 basis now, for four years, supported that effort.
25 So, the finances are between the academic

1 societies, the Lifeline, and the manufacturers.

2 DR. LASKEY: Thank you very much.

3 DR. RODNEY WHITE: Thank you.

4 DR. LASKEY: Are there any other members
5 of the audience requesting time?

6 [No response]

7 Again, thank you, Dr. White. I would like
8 to close this portion of the open public hearing
9 and move on. I would like to move to the sponsor's
10 presentation at this point. I just want to remind
11 people, we are shooting for a twelve o'clock break
12 for lunch, to stay on schedule. Dr. Harvey?

13 DR. HARVEY: Please remember to introduce
14 yourself when you come to the podium to speak, and
15 to state your conflict of interest and also to use
16 the mike whenever you are asked any questions and
17 need to come forward.

18 **Sponsor Presentation**

19 **Introduction**

20 MR. SININGER: Good morning.

21 [Slide]

22 I am John Sininger, with W.L. Gore &
23 Associates. I am responsible for Gore's Medical
24 Products Division. W.L. Gore is the sponsor for
25 this premarket application for the Excluder

1 Bifurcated Endoprosthesis. Because this is such a
2 mouthful, we will be referring to this frequently
3 in the presentation as EBE. So, as you hear that
4 through the presentation you will know what we are
5 referring to.

6 [Slide]

7 Gore is a 44-year old high technology
8 company engaged in the development, manufacturing
9 and sales of a broad range of high technology
10 products. Gore's history and our reputation in all
11 the markets that we serve is that we provide only
12 the highest quality and highest performance
13 products in all of the markets that we serve. Gore
14 provides a broad range of products, from
15 sophisticated aerospace applications to Gortex
16 fabrics, which most people know Gore by, to
17 microfiltration, industrial filtration products
18 and, of course, Gore's Medical Products Division
19 which is presenting this presentation today.

20 [Slide]

21 Gore Medical has been in the business of
22 developing, making and selling products for almost
23 25 years, over 25 years. In that period of time
24 Gore has developed many products that really serve
25 many different patient populations, and Gore was

1 actually a real pioneer in some of the early work
2 for products to repair peripheral vasculature.

3 [Slide]

4 In this period of time there have been
5 over seven and a half million clinical implants of
6 Gore medical products worldwide. The significance
7 of this is that this clinical background offers us
8 a clear understanding of the need for a safe and
9 effective treatment for abdominal aortic aneurysm
10 repair.

11 I have had the privilege of being with
12 Gore and Gore's Medical Products Division for 25 of
13 those years. I recall being part of some of the
14 very early development work for not only our
15 peripheral vascular grafts but also our first
16 aortic graft. It was a real pleasure and a
17 privilege to be able to provide a product that made
18 such a significant difference to patients with this
19 life-threatening AAA disease. It is even more
20 significant for me today, 25 years later, to be
21 standing here, introducing our presenters with a
22 product that we believe makes an even greater
23 difference to the patients who receive these
24 products. We have seen these products make a real
25 difference in the lives of patients who are

1 receiving these devices.

2 It is for this reason that I am very
3 excited about the opportunity to be here to
4 introduce our presenters who will be presenting
5 data which we believe supports the primary safety
6 and efficacy of the device. We hope your review of
7 the data substantiates that conclusion.

8 [Slide]

9 Our agenda today is that Dave Williams,
10 who is a Gore associate, will be presenting an
11 overview of the device and study overview. Dr.
12 David Brewster, from Harvard Medical School, will
13 be providing a background to abdominal aortic
14 aneurysm repair. Dr. David Naftel will be talking
15 about trial design and trial management. Finally,
16 Dr. Jon Matsumura, who is also our principal
17 investigator, will be presenting the clinical
18 results.

19 In addition to these individuals, there
20 are a number of clinical investigators who are
21 here, as well as a number of other Gore associates
22 who are here to answer any questions you may have.

23 Finally, I would like to sincerely thank
24 all of you, the FDA, for all your time and
25 consideration and efforts in reviewing all this

1 data, and look forward to a lively discussion and
2 your conclusions. Thank you. Dave Williams?

3 **Product and Study Overview**

4 MR. WILLIAMS: Thank you, John.

5 [Slide]

6 Good morning and thank you, ladies and
7 gentlemen of the panel for the opportunity to
8 present today.

9 [Slide]

10 I will do a quick overview of the Excluder
11 Bifurcated Endoprosthesis, or EBE, device
12 description including its deployment; briefly
13 summarize the preclinical evaluation; and then
14 provide a brief overview of the clinical evaluation
15 experience.

16 [Slide]

17 The device design is bifurcated, modular
18 and it has a fully supported self-expanding nitinol
19 stent which basically supports an EPTFE or PTFE
20 vascular graft on the blood contact surface.

21 These are the various components of the
22 modular device. You can see the two primary
23 components. We have the trunk ipsolateral leg. We
24 have the contralateral leg. We have an aortic
25 extender and an iliac extender.

1 There is a unique feature in this device
2 in that the outer nitinol stent is attached to the
3 underlying PTFE graft material in a sutureless
4 fashion. It uses fluoropolymer films to bond the
5 stent to the underlying graft. At the proximal end
6 or the aortic trunk end of the device you can see
7 that there are anchors for active fixation, and
8 there is also an external sealing cuff to aid in
9 hemostasis relative to the proximal application of
10 the device.

11 The device also relies on oversizing as
12 part of the fundamental design performance in that
13 you choose device sizes based on the patient's
14 healthy anatomy in both the proximal aortic neck,
15 infrarenal neck, as well as the common or external
16 iliac vessels. So, the device is relying on both
17 active and passive fixation and oversizing to
18 create fixation and hemostatic seal to exclude the
19 aneurysm.

20 [Slide]

21 This is a picture of the device's two main
22 components as they would be assembled in situ. So,
23 you have the trunk ipsolateral leg component with
24 the contralateral leg component overlapping or
25 docking into that primary component.

1 [Slide]

2 The device is loaded or constrained down
3 onto a delivery catheter. The trunk ipsilateral
4 leg component is constrained or loaded onto an 18
5 French delivery profile catheter. You can see here
6 that there is a PTFE sleeve or corset that is used
7 to hold the device on the catheter in the
8 constrained position. That corset or constraining
9 sleeve is laced in place with a single PTFE fiber
10 which we refer to as the deployment line. That
11 line runs the length of the delivery catheter and
12 exits here, in the hub end or the operator end of
13 the catheter, and is connected to this deployment
14 knob. In the hub you can also see the Y valve or
15 the Touhy-Bourst which contains the guidewire
16 lumen, a fleshing port in addition to the
17 deployment knob.

18 [Slide]

19 This is an image of the contralateral leg
20 component partially deployed. This gives you a
21 feel for the constraining sleeve as it is being
22 unlaced through the deployment line retraction,
23 allowing the nitinol stent to self-expand,
24 deploying the device into position.

25 [Slide]

1 Next we will see a computer animation
2 briefly of the positioning and deployment of the
3 device in the infrarenal aortic anatomy. The
4 anatomy is accessed by retrograde guidewire
5 cannulation. Over the guidewire an 18 French
6 vascular access sheath is placed into the
7 infrarenal aortic anatomy. With the sheath in
8 place, the trunk ipsolateral component can then be
9 tracked into a proximate position. The vascular
10 introducer sheath is withdrawn to expose the
11 device. Then, under fluoroscopic visualization the
12 proximal end markers and the contralateral and
13 ipsolateral orientation markers can be located to
14 properly position the direction of the deployment
15 of the two limb components.

16 Once the device has been fine-tuned or
17 positioned relative to the lowest renal artery and
18 to proper lateral orientation of the limbs, the
19 deployment knob is pulled; the corset opens and
20 allows the self-expanding stent to actuate the
21 device deployment.

22 Trunk ballooning is recommended at this
23 point to further optimize the dilatation of the
24 device. Contralateral access is gained via
25 guidewire and central lumen position is confirmed.

1 Then, a 12 French sheath is placed inside the
2 contralateral leg hole. The contralateral leg is
3 delivered and deployed in similar fashion.
4 Adjunctive ballooning in the top and bottom of the
5 device is recommended, and adjunctive aortic
6 extenders or iliac extenders may be placed to
7 further optimize the procedure.

8 [Slide]

9 The preclinical evaluation summary is that
10 all evaluations, including toxicology,
11 biocompatibility, in vitro and in vivo tests
12 demonstrate that the EBE system meets the
13 functional requirements for aortic endovascular
14 devices.

15 [Slide]

16 A quick overview of the various EBE
17 clinical studies is listed here. The device that
18 you are considering today is our first generation
19 Excluder device which went through both European
20 and U.S. feasibility trials starting in late '97
21 and ending in mid-'98.

22 The trial data under consideration today
23 is the pivotal trial which effectively started to
24 enroll patients in December of '98 and stopped
25 enrolling patients in January of 2000. There was a

1 continued access portion to this pivotal trial and,
2 subsequent to that, we have continued to study, in
3 an IDE format, in the United States a second
4 generation EBE device and those studies are
5 ongoing.

6 [Slide]

7 The worldwide clinical experience with the
8 EBE was initiated in Europe in September of 1997.
9 The EBE has been commercially available outside the
10 United States since 1998, and with the continuing
11 U.S. clinical trial evaluations, combined with the
12 rest of the world commercial use experience, we now
13 have exceeded 4400 implants or in excess of 10,000
14 individual component pieces.

15 [Slide]

16 Refocusing back onto the pivotal study
17 under consideration today, the data includes events
18 through February 29th of this year. Although the
19 primary and secondary hypotheses for the study were
20 evaluated to a 12-month endpoint, the protocol
21 amendments and patient consents allowed for patient
22 follow-up out to five years of this particular
23 patient cohort. Thank you for your attention.

24 **Abdominal Aortic Aneurysm Background**

25 DR. BREWSTER: Mr. Chairman, distinguished

1 panel members, good morning.

2 [Slide]

3 My name is David Brewster. I am a
4 clinical professor of surgery at Harvard Medical
5 School and at the Massachusetts General Hospital,
6 in Boston where I also serve as director of
7 endovascular surgery.

8 W.L. Gore has paid my expenses to be here
9 with you today, but I have no financial interests
10 in the device or the company, nor in the outcome of
11 this meeting today.

12 I am a board-certified vascular surgeon
13 with over 25 years of experience in the management
14 of patients with abdominal aortic aneurysms.
15 During this time, I have repaired approximately
16 1500 aneurysms by conventional open surgical graft
17 repair, and during the last eight years I have had
18 considerable experience with the alternative mode
19 of therapy being considered today. During this
20 time I have treated nearly 500 patients with
21 endovascular grafts, employing a wide variety of
22 different devices and participating as principal
23 investigator in five FDA clinical Phase II trials
24 including, of course, the Gore EBE trial being
25 presented this morning.

1 In the next few minutes I will review some
2 facts regarding aortic aneurysm in the hope that
3 this information will serve as a backdrop or
4 yardstick, if you will, by which to evaluate the
5 results of the EBE clinical trial and, therefore,
6 aid you in your deliberations.

7 Apologies in advance to those panel
8 members already familiar with this material, but my
9 goal is to ensure that all panelists are acquainted
10 with basic facts regarding epidemiology and natural
11 history of aneurysms, as well as the anticipated
12 outcome of traditional open surgical repair.

13 [Slide]

14 Aortic aneurysms, no doubt, represent an
15 important public health problem. Approximately
16 200,000 new cases are diagnosed in the U.S. each
17 year and 50,000 procedures approximately are
18 performed per year for aneurysm repair. The
19 principal goal of these procedures is to prevent
20 aneurysm rupture, which is the 13th leading cause
21 of death in the U.S. and 10th leading cause if one
22 considers only men over the age of 65 which, of
23 course, is the most common patient cohort.

24 [Slide]

25 It is well recognized that aneurysms are

1 being encountered more frequently in contemporary
2 practice. This is due to both better diagnosis and
3 recognition by a variety of imaging techniques, as
4 well as an apparent true increase in prevalence.
5 This latter phenomenon appears largely attributable
6 to the well-documented aging of our population.
7 The occurrence of an aneurysm, as well as the
8 rupture risk, are known to increase sharply with
9 age. As a consequence of such trends, several
10 projections indicate a substantial increase in the
11 number of patients with aneurysms who will require
12 treatment in the next several decades, many of them
13 likely elderly and with co-morbidities that make
14 them at increased risk for standard open repair.

15 [Slide]

16 The expected natural history of an
17 aneurysm is gradual expansion leading to eventual
18 rupture unless this process is interrupted first by
19 death of the host or from another cause. The goal
20 of treatment, therefore, becomes treatment to
21 prevent rupture in susceptible individuals.
22 Decision-making currently lacks true scientific
23 precision and, rather, represents a reasoned
24 balancing of estimated risks of rupture versus
25 repair, and an individualized approach to each

1 patient is emphasized based upon his or her own
2 age, co-morbid conditions and, very importantly,
3 treatment preferences. Considerable clinical
4 judgment remains vitally important in such
5 decision-making.

6 [Slide]

7 Although a number of factors contribute to
8 rupture risk, it is widely accepted that aneurysm
9 size, as measured by maximal diameter, is the most
10 important determinant of rupture risk. While it is
11 well established that truly small aneurysms have a
12 low risk of rupture, this risk begins to sharply
13 increase after the aneurysm reaches 4.5-5.0 cm in
14 size.

15 While the fairly wide range of estimated
16 rupture risk, indicated here from literature
17 review, indicates the considerable differences
18 reported in the literature from one series to
19 another, a recent meta-analysis indicates an annual
20 rupture risk of approximately 10 percent per year
21 for aneurysms in the size range typical of those
22 treated in the EBE clinical trial, as will be
23 presented shortly.

24 [Slide]

25 Left untreated, aneurysms in the size

1 range relevant to most clinical decision-making can
2 be expected to enlarge approximately 10 percent per
3 year in maximal diameter. Hence, for aneurysms
4 typical of the EBE trial, enlargement of
5 approximately 0.5 cm or 5 mm per year would be
6 anticipated.

7 [Slide]

8 During the past five decades standard open
9 operative repair has been well established as a
10 very effective and durable method of repair.
11 Despite this, however, considerable room for
12 improvement in the outcomes of therapy continue to
13 exist. Although many individual referral-based
14 reports from institutions of excellence suggest
15 mortality rates of open repair well below five
16 percent, many recent population-based series from
17 large statewide or national databases reveal a
18 real-world mortality more in the range of five to
19 ten percent even in current practice.

20 In addition, all vascular surgeons
21 recognize the substantial morbidity and
22 complication rates of this extensive standard
23 repair. A rate that is often substantially higher
24 in elderly patients are those with associated
25 co-morbidities. Patients who are often a typical

1 cohort are a sizeable percentage of those requiring
2 treatment. Even in the best of circumstances,
3 recovering from open repair takes many months and,
4 indeed, several recent quality of life studies
5 indicate that a significant number of older
6 patients never quite regain their preoperative
7 baseline functional status.

8 For all of these reasons, many high risk
9 patients are often currently denied open surgical
10 repair and left with the fear and concern that
11 rupture may unpredictably occur at any time.

12 [Slide]

13 The goals of endovascular repair are to
14 achieve a repair quite similar to that of open
15 graft insertion, but to accomplish this in a manner
16 which is less invasive by working within the
17 vascular system and using small incisions and
18 rather minimal anesthesia. The endovascular device
19 is placed within the aneurysm sac, and with secure
20 anchoring and fixation above and below the aneurysm
21 in relatively normal and healthy arterial segments
22 exclusion of the weakened portion of the aorta from
23 arterial circulation and pressure is achieved,
24 thus, eliminating the danger of rupture.

25 [Slide]

1 I would like to conclude with a brief
2 review of a concept which is unique to endovascular
3 aneurysm repair, that of endoleak. Endoleak
4 denotes continued perfusion of the aneurysm sac as
5 detected by one of several imaging modalities.
6 Endoleaks have been classified by the source or
7 cause of such failure to totally exclude the
8 aneurysm from the circulation. Type IV leaks refer
9 to those with transgraft seepage which may occur as
10 an early phenomenon in devices constructed of
11 porous fabrics. These are usually self-limited and
12 of little to no clinical importance. In contrast,
13 type I leaks are those due to failure to achieve a
14 hemostatic seal at either the proximal or distal
15 attachment zones, while type III leaks refer to
16 leakage or continued perfusion of the sac caused by
17 defects in the graft material itself or leakage at
18 junction points of modular devices.

19 Because both of these types, type I and
20 type III leaks transmit full arterial pressure to
21 the aneurysm sac, both are generally accepted as
22 potentially dangerous and an indication for further
23 intervention.

24 In contrast, type II leaks are caused by
25 reversed flow in normal arterial branches which may

1 remain patent. When the aneurysm is excluded by
2 the endoluminal device and pressure falls to low
3 levels within the sac normal antigrade flow may
4 reverse and lead to continued perfusion via lumbar
5 arteries or patent inferior mesenteric vessel.

6 Unlike type I and type III leaks, however,
7 the clinical significance of the common type II
8 endoleaks is much more uncertain as many seal
9 spontaneously at later intervals and undesirable
10 clinical outcomes, such as further growth and
11 aneurysm rupture, are very infrequent.

12 [Slide]

13 In summary, pertinent facts to remember as
14 we hear the results of the EBE clinical trial are
15 that aneurysms of the 5-6 cm size range, typical of
16 those aneurysms treated in the trial, enlarge on an
17 average rate of approximately 10 percent or half a
18 centimeter a year and carry an annual rupture risk
19 of approximately 10 percent per year. Obviously,
20 the therapy seeks to alter and improve on these
21 natural history behaviors.

22 Although open surgical repair remains a
23 very effective and durable treatment, morbidity and
24 mortality risks remain substantial and other
25 limitations exist. In appropriate patients

1 endovascular repair offers a safe and effective
2 alternative with many potential advantages. Thank
3 you.

4 **Trial Design and Study Management**

5 DR. NAFTEL: My name is David Naftel, and
6 I appreciate this opportunity to speak to you.

7 [Slide]

8 I am a consultant for W.L. Gore and I have
9 no financial interest except fee for service. I am
10 a professor of biostatistics and professor of
11 surgery at the University of Alabama in Birmingham.
12 I will be discussing the trial design and study
13 management this morning.

14 [Slide]

15 The first indication for use, the Excluder
16 Bifurcated Endoprosthesis is intended to exclude
17 the aneurysm from blood circulation in patients
18 diagnosed with infrarenal AAA disease and who have
19 appropriate anatomy. It is this indication for use
20 that drives the corpus of this study.

21 [Slide]

22 So, the purpose of the study that we will
23 discuss this morning is to determine the efficacy
24 and the safety of the EBE for the primary treatment
25 of infrarenal AAA.

1 [Slide]

2 There are two main hypotheses. The
3 primary safety hypothesis is that subjects treated
4 with the EBE have a proportion of major adverse
5 events that is less than subjects treated with open
6 surgical repair as evaluated through 12 months. A
7 major adverse event is defined as any one of the
8 following: First, requires therapy and short
9 hospitalization; or requires major therapy and
10 unplanned increase in level of care and prolonged
11 hospitalization; or permanent adverse sequelae or
12 death.

13 [Slide]

14 The primary efficacy hypothesis is that
15 the EBE is an effective treatment method to exclude
16 the aneurysm from blood circulation when used in
17 the primary treatment of infrarenal AAA as
18 evaluated at 12 months. Efficacy is defined as all
19 of the following: Absence of endoleaks with or
20 without treatment; absence of aneurysm enlargement,
21 defined as greater than or equal to 5 mm; and
22 absence of major device efficacy complications.

23 [Slide]

24 There are also secondary hypotheses. That
25 is, compared with the control subjects, the EBE

1 subjects will have shorter stay in the intensive
2 care unit; shorter hospital length of stay; and
3 they will return to normal activities faster.

4 [Slide]

5 Here is the study design. It was a
6 multicenter, prospective, intent-to-treat design.
7 It is non-randomized but there are concurrent open
8 surgical controls. The hypotheses are all focused
9 on a 12-month duration. There is an independent
10 core lab at Cleveland Clinic Foundation. There is
11 a clinical events committee to review the major and
12 minor adverse events, and then a data safety
13 monitoring board to continually monitor the safety
14 of the study.

15 [Slide]

16 The primary safety and efficacy endpoints
17 that were focused on in designing the study were,
18 first of all, a 15 percent difference in major
19 adverse events between the two groups at one year
20 and at least an 80 percent primary efficacy success
21 also at one year. The 15 percent difference was
22 used in the calculation of the sample size and it
23 produced a minimum number of available subjects at
24 one year to be 78 control patients and 156 EBE
25 patients. This was based on a two-sided comparison

1 with an alpha level of 0.05 and 80 percent power.
2 We used a ratio of two EBE to one control subject.

3 [Slide]

4 Standard statistical methods were used.
5 Multivariable analyses were used to produce
6 risk-adjusted comparisons of the two groups. These
7 included both logistic regression and Cox
8 proportional hazard. For the time-related events
9 we used Kaplan-Meier for time to death and also
10 time to first major adverse event. The Nelson
11 method was used to produce cumulative adverse
12 events across time on a per patient basis. Other
13 standard methods were used to compare the two
14 groups.

15 [Slide]

16 A number of inclusion and exclusion
17 criteria were used, and here I will focus only on
18 the anatomic criteria. For the control group only
19 there had to be planned or expected use of
20 infrarenal clamp. For the EBE group there had to
21 be proximal aortic neck length greater than or
22 equal to 15 mm; a proximal aortic neck angulation
23 less than 60 degrees; and no significant thrombus
24 at the arterial implant site.

25 [Slide]

1 Here are the follow-up requirements in the
2 study that were adhered to. Contrast-enhanced CT
3 was performed at one month in the EBE group; at
4 three months if an endoleak had been found at one
5 month; and then also a CT at six months and at one
6 year and annually. The one-year and annual CTs
7 were performed also for the control subjects.

8 Abdominal x-rays were performed at
9 discharge in the EBE group and again at six months
10 and annually. Bilateral ankle brachial index and
11 physical exams were conducted in both groups at
12 discharge, one month, six months and 12 months. In
13 addition, the EBE group had a physical exam at
14 three months if an endoleak was found at one month.

15 [Slide]

16 A variety of centers, 19 centers across
17 the country including academic, non-academic and
18 community hospitals and a variety of specialists
19 for vascular disease were included in the study.
20 Thank you.

21 Pivotal Study Clinical Results

22 DR. MATSUMURA: Good morning.

23 [Slide]

24 My name is Jon Matsumura. I am a paid
25 consultant for W.L. Gore. I am also a

1 board-certified vascular surgeon and have
2 concentrated most of my professional academic
3 interest in endovascular therapy of aortic
4 aneurysms. I am grateful for the privilege to work
5 as the PI with the 19-site investigators in concert
6 with the trial sponsor. It is also my pleasure to
7 present the pivotal study data to you this morning.

8 [Slide]

9 Just to reiterate, the 12-month follow-up
10 which I will be presenting you first had a data
11 cut-off point in June of 2001. This included 260
12 patients in the EBE and 101 patients in the control
13 group.

14 [Slide]

15 Let's go right to the pre-procedure
16 results. These are clinical characteristics which
17 were found to be significantly different between
18 the two groups. Specifically, the EBE group had an
19 average age that was three years older than the
20 control group although there was a wide range. The
21 EBE group had a higher proportion of men compared
22 to the control group, and the EBE group had a lower
23 proportion of patients with symptomatic abdominal
24 aortic aneurysm.

25 There were many other clinical

1 characteristics that were evaluated to determine
2 the comparability of the two treatment groups.
3 This included a medical history of coronary-artery
4 disease, arrhythmia, valvular heart disease, CHF,
5 stroke, history of inflammatory aneurysm or family
6 history of aneurysms or other aneurysms, history of
7 peripheral arterial disease or prior vascular
8 interventions, and none of these were different
9 between the two groups.

10 [Slide]

11 Additional clinical characteristics that
12 were compared include long-term steroid use,
13 history of thrombotic event, emphysema, smoking
14 history, renal failure or paralysis, erectile
15 dysfunction in men, hepatic dysfunction, bleeding
16 disorder and history of cancer. There were no
17 differences between the two groups in any of these
18 clinical characteristics.

19 [Slide]

20 We also used many of the
21 Society-determined risk factor score systems such
22 as the ASA by the anesthesiologists and the New
23 York Heart Association, which had no differences
24 between the two groups. The SVS joint societies
25 risk factor score system is compared here. I will

1 point out that in the hyperlipidemia subcategory
2 there was a difference and the EBE group had more
3 hyperlipidemia than the control group, although
4 there was no difference in the composite SVS risk
5 score index between the two groups. Based on these
6 comparisons, we figured that the two treatment
7 groups are comparable.

8 [Slide]

9 We also looked at many anatomic variables
10 as well as disease states for the arterial anatomy.
11 Shown here are six that were significantly
12 different between the two groups. The EBE group
13 had an average aneurysm size 3 mm smaller than the
14 control group. The proximal aortic neck was
15 longer, narrower and had less angulation in the EBE
16 group compared to the control group. The left and
17 right common iliac arteries were smaller in the EBE
18 group compared to the control group. These five
19 differences would be expected given the protocol
20 requirements for endovascular aneurysm repair.

21 I won't go through all of them but we
22 studied 40 additional pretreatment anatomic and
23 disease variables that are in the PMA and there are
24 no differences between the treatment groups in
25 these other variables.

1 [Slide]

2 This figure examines in more detail the
3 differences in aneurysm size. Although there was a
4 difference in mean aneurysm size between the EBE
5 and the control group, you can see that a wide
6 range of aneurysm sizes were treated in both
7 groups. In addition, the majority of aneurysms
8 were over 5 cm in size in both treatment groups.

9 [Slide]

10 Let's get to what the immediate procedure
11 results are. In terms of EBE deployment evaluation
12 at the initial procedure, there were 235 patients
13 enrolled in the EBE group. All of the patients
14 received one trunk ipsolateral leg and one
15 contralateral leg and 100 percent deployment
16 success. In addition to those components, a third
17 or 32 percent of the patients had either one or
18 more aortic extenders, one or more iliac extenders
19 or one or more both aortic and iliac extenders as
20 part of their initial treatment, all of which were
21 deployed successfully.

22 [Slide]

23 Some of the immediate procedure results
24 showed improved outcomes in the EBE and the control
25 group and are shown here. The mean procedure time

1 was shorter in the EBE group. The mean blood loss
2 was less compared to the controls, and the chance
3 that you would require homologous blood transfusion
4 was less in the EBE group compared to control.

5 [Slide]

6 We are going to go into the 12-month data,
7 and before I show the results of the actual
8 comparisons I want to show you the accountability.
9 At one month these are the number of patients
10 available who had not died, withdrawn or been lost
11 to follow-up. At 12 months there are 81 controls
12 and 215 EBE patients available for follow-up. You
13 can see that we had over 90 percent compliance with
14 follow-up clinical visits at the time points for
15 each of the two treatment groups.

16 [Slide]

17 To remind you, our primary safety
18 hypothesis was designed to evaluate major adverse
19 events through the 12-month time point.

20 [Slide]

21 This is a breakdown, first, of the major
22 adverse events in the two groups. If you look at
23 any major adverse event, there were 57 percent of
24 the control patients who had one and 14 percent of
25 the EBE group, and this was a significant reduction

1 in major adverse events in the first 30 days.

2 When further breaking this down into organ
3 systems or subgroups of major adverse events, this
4 reduction of major adverse events was in several
5 categories. In bleeding and pulmonary there was a
6 12-fold reduction. In cardiac there was a 4-fold
7 reduction. In bowel, an 8-fold reduction and in
8 vascular a 6-fold reduction in complications in the
9 EBE group compared to the control group.

10 [Slide]

11 When you look at the major adverse events
12 that occur between the 30-day time point and 12
13 months, 25 percent of the control group experienced
14 one and 27 percent in the EBE group, which are not
15 different rates. It is important to note that on
16 the clinical events committee we considered
17 interventions in the EBE group that were performed
18 for endoleak or aneurysm enlargement where the
19 patients stayed a day in the hospital, such as a
20 coil embolization, as major adverse events in this
21 group. If you look at the other categories broken
22 down by organ system, there were no differences in
23 the rate of adverse events after the 30-day time
24 point to 12 months between EBE and the control
25 group.

1 [Slide]

2 This figure puts those two data sets
3 together and I am going to take some time to go
4 through this. On the X axis is the months after
5 procedure, 0, 12 or 14. On the Y axis is
6 cumulative major adverse events as a rate per
7 patient. The yellow curve is the EBE group; the
8 white curve is the control group.

9 You can see that in the first 30 days
10 there is a marked increase in the rate of adverse
11 events in the control group compared to the EBE
12 group. But after that time point these curves are
13 relatively parallel and the ongoing rate of major
14 adverse events after the first month is similar in
15 the two groups.

16 Numerically, this can be seen here in this
17 table. If you are in the control group you had an
18 average of 1.2 major adverse events per patient in
19 the first month. If you were in the EBE group you
20 had an average of one event per five patients. At
21 12 months there continued to be a difference. You
22 had a chance of 1.8 adverse events per patient in
23 the control group and 0.9 adverse events in the EBE
24 group. These are significantly different by a
25 Nelson test.

1 [Slide]

2 Another way of comparing major adverse
3 events is not just to look at how many
4 complications you have on a per patient basis, but
5 what is the chance that you will have any major
6 adverse event. I think there might be concern that
7 some patients have many adverse events so we also
8 want to look at the proportion or the chance that
9 you will have any or the first major adverse event.

10 This is a Kaplan-Meier depiction of that.
11 Again, on the X axis is the time and on the Y axis
12 freedom from first major adverse event. The yellow
13 again is the EBE and the white is the control. At
14 one month, again, there is a marked reduction in
15 the chance of having even one major adverse event,
16 from 86 percent in the EBE group to 43 percent in
17 the control group. At 12 months the freedom from
18 major adverse events continues to be different.
19 There is a 62 percent chance in the EBE group of
20 never having had a major adverse event and a 35
21 percent chance in the control group, and these are
22 significant by log rank.

23 [Slide]

24 In addition to these univariable analyses,
25 we conducted the multivariable analysis. This is

1 logistic regression looking at independent risk
2 factors for early major adverse events and also to
3 give you a risk-adjusted estimate of the risk
4 associated with the treatment group.

5 Four risk factors were identified other
6 than treatment group. A history of myocardial
7 infarction, history of thrombotic event, and SVS
8 pulmonary risk score of one or greater, and a lower
9 platelet count were independent risk factors for
10 early major adverse events. More importantly, we
11 determined that the control group was a strong and
12 independent risk factor for early major adverse
13 events with a 12-fold odds ratio.

14 [Slide]

15 We did a similar multivariable analysis
16 looking at late major adverse events. Using the
17 Cox model, these five risk factors were found to be
18 independent risk factors for late major adverse
19 events: an older age, smaller body mass index, a
20 history of prior vascular intervention, a history
21 of symptomatic aneurysm and an increased proximal
22 neck angle.

23 In this model we forced in-treatment group
24 to see if it could predict late major adverse
25 events and the EBE treatment group was not an

1 independent risk factor for late major adverse
2 events.

3 [Slide]

4 We also looked at survival. Although the
5 study wasn't powered to detect this it is of
6 obvious interest. The survival is similar in the
7 EBE and control groups. At 12 months the control
8 group survival is 94 percent and 92 percent in the
9 EBE group. These are not different.

10 [Slide]

11 We did a multivariable analysis to look at
12 survival as well. These five risk factors were
13 found to be independent risk factors for mortality
14 in the Cox model. Again, an SVS pulmonary risk
15 score of one or greater, a history of erectile
16 dysfunction if you are a man, a lower platelet
17 count, a lower initial ankle brachial index, and a
18 larger difference between your maximum aneurysm
19 diameter and proximal aortic diameter upon entry
20 into the study. Again, we forced in-treatment
21 group to see if there was any effect from treatment
22 group allocation, and treatment group is not an
23 independent risk factor for mortality.

24 [Slide]

25 Recently, reporting standards and the

1 evaluation of endovascular repair have been
2 published out of the joint societies, and they
3 defined a primary outcome measure of endovascular
4 repair as aneurysm-related deaths. These would be
5 defined as deaths due to aneurysm rupture, or
6 related to the primary procedure, or a secondary
7 procedure such as a later open surgical conversion.

8 We took a cautious interpretation of what
9 "related" means and said that it is any death
10 within 30 days of a primary or secondary procedure
11 or during the same hospitalization. With this
12 definition, we calculated aneurysm-related survival
13 in the two groups. The survival curves are
14 similar. There is a 98 percent aneurysm-related
15 survival at one year in the control group and the
16 EBE group.

17 [Slide]

18 In summary of our safety analysis,
19 compared with open surgical repair for the primary
20 treatment of aneurysm, the data demonstrate that
21 EBE is safe. There is a lower rate of major
22 adverse events, similar overall survival and
23 similar aneurysm-related survival and there were no
24 device-related deaths.

25 [Slide]

1 Just to remind you, we are going to switch
2 to the efficacy evaluation. The hypothesis was
3 based on evaluation exclusion from the blood
4 circulation at 12 months, and it had those three
5 components of endoleak, aneurysm size increase and
6 device efficacy complications. I will go through
7 each of them.

8 [Slide]

9 These are the endoleak rates from our core
10 lab. The one month is in yellow; the 12 months is
11 in grey. The majority of the patients did not have
12 an endoleak. Of those who had an endoleak, most of
13 the endoleaks were of the type II variety. There
14 were some type I endoleaks at a lower frequency and
15 some endoleaks at a lower frequency of
16 indeterminate origin. There were no type III or
17 type IV leaks found by the core lab.

18 [Slide]

19 Aneurysm diameter size change evaluation
20 by the core lab is shown in this bar graph.
21 Fourteen percent of the patients had an aneurysm
22 size decrease of 5 mm or more; 79 percent of the
23 patients had no change in their aneurysm size; 7
24 percent of the patients had an aneurysm size
25 increase of 5 mm or more.

1 [Slide]

2 In addressing aneurysm size increase and
3 endoleaks, it is important to note that during the
4 course of this study the investigators met several
5 times, at least annually in our investigator
6 meetings, and we formed collectively a treatment
7 guideline set, and this is it. We felt that
8 aneurysms with type I endoleaks, type III endoleaks
9 and enlargement regardless of endoleak status
10 should be intensively studied and considered for
11 catheter-based reintervention or conversion to open
12 repair.

13 It is important that this consideration
14 include the local investigator and the attending
15 physician's assessment. As Dr. Brewster mentioned,
16 there is still significant judgment used in the
17 treatment of these patients and it includes the
18 individual patient's co-morbidities, life
19 expectancy and, of course, the patient's own
20 personal choices. I would add that the
21 endovascular treatment of aortic aneurysms is
22 really an evolving process and these guidelines may
23 change in the future.

24 [Slide]

25 These are the reinterventions that were

1 conducted during the first year in the pivotal
2 study. There were 15 patients who had 17
3 reinterventions. Of the 15 patients, the majority
4 of interventions were done for an endoleak but
5 there was one patient who had a ligation performed,
6 who had both an endoleak and aneurysm size
7 increase, a ligation of a hypergastric artery. The
8 other 16 procedures in the 14 patients were all
9 catheter-based embolization procedures or other
10 embolization procedures.

11 [Slide]

12 This is the third component of the
13 efficacy hypothesis, the major device efficacy
14 complications. I will go through this table slowly
15 as well. There were no patients who had access
16 failure. As mentioned, there was 100 percent
17 deployment success. There were no intraoperative
18 or early conversions in this group.

19 There were two patients who had occlusion
20 of a branch artery, one early and one late. One
21 was a hypergastric artery that was inadvertently
22 occluded and led to lasting butt claudication. It
23 was determined to be a major DEC. There was
24 another patient who developed occlusion of the left
25 renal artery and, at five weeks, underwent an

1 ilioarenal bypass for that problem.

2 There were no patients found to have lumen
3 obstruction or extrusion or erosion. No patients
4 were found by the sites to have prosthetic material
5 fatigue, and we will have more on that later.
6 There was one patient who had prosthesis migration
7 that required therapy with aortic cuffs. This
8 patient had a main trunk component placed at the
9 procedure. After deployment, it was noted later to
10 move down a couple of centimeters caudally and an
11 additional aortic extender was placed to treat that
12 patient.

13 There were no patients with prosthesis
14 realignment, and in the pivotal study there were no
15 patients with aneurysm rupture. This, aneurysm
16 rupture, is an important device efficacy
17 complication, and I will point out that there is
18 one rupture in the U.S. feasibility study; two
19 ruptures in the European experience. The details
20 are in your panel pack and perhaps we will explore
21 those in the Q&A as well.

22 [Slide]

23 In summary of the device performance,
24 there was 100 percent patency. There were no limb
25 occlusions or clinical adverse events related to

1 device patency. There was 100 percent freedom from
2 aneurysm rupture.

3 [Slide]

4 On this slide I am going to start from the
5 bottom. If you take those three efficacy
6 complications, you have 27 with endoleak at 12
7 months; 13 with aneurysm enlargement; and 3 that I
8 just described with the major DECs, device efficacy
9 complications. Because these 43 complications
10 occurred in 38 patients, they are overlapping.
11 That is how you get to 38 patients, and we used the
12 denominator cautiously of 196 patients who had
13 12-month CT core lab data available, giving a
14 primary efficacy success of 80.6 percent with these
15 95 percent confidence intervals. I will reiterate
16 that these were based on core lab assessment for
17 endoleak and aneurysm enlargement.

18 [Slide]

19 Our assessment of the efficacy data is
20 that the EBE is an effective treatment method to
21 exclude the aneurysm from the blood circulation.

22 [Slide]

23 The core lab looked at other imaging
24 findings as well in its review of CTs and abdominal
25 x-rays to evaluate device integrity, device patency

1 and trunk and component migration.

2 [Slide]

3 In terms of device integrity, the core lab
4 identified a fracture in a discharge film of one
5 patient, at a rate of 0.4 percent for that
6 interval. No fractures were subsequently
7 identified in the 12 months. I will discuss this
8 fracture in the next slide. Again, because this is
9 an important issue, device integrity or fractures,
10 I will point out that there is one fracture in the
11 second generation trial and there is a fracture
12 identified in the European experience. Hopefully,
13 we will discuss those in more detail during the
14 Q&A.

15 [Slide]

16 In terms of the pivotal study, this
17 patient had the fracture visible on the discharge
18 film, which was the only x-ray performed in that
19 patient. There were no clinical consequences. At
20 12 months a CT scan was performed and evaluated
21 both by the site and the core lab, and no endoleak,
22 no aneurysm enlargement, and no migration was
23 identified. Unfortunately, the patient was
24 diagnosed with inoperable liver cancer in the
25 second year and died of this. No autopsy was

1 performed and the device is unavailable for
2 analysis.

3 [Slide]

4 Some of the other core lab imaging
5 findings--a small percentage of patients were
6 identified by core lab review with device lumen
7 narrowing; trunk migration of 10 mm or more
8 relative to the arterial landmarks; and component
9 migration of 10 mm or more relative to other
10 components. In these patients there were no type I
11 or type III endoleaks; no aneurysm enlargement; and
12 no vascular adverse events or reinterventions.

13 [Slide]

14 We also had these secondary hypotheses
15 which basically deal with how did the patients
16 recover from the procedure.

17 [Slide]

18 These are the data on the secondary
19 outcomes. The EBE patients had a 10-fold reduction
20 in length of ICU stay; a 5-fold reduction in mean
21 length of hospital stay; a reduction in time to
22 ambulation; and also a reduction in the time to
23 return to normal activities as reported by the
24 patients themselves.

25 [Slide]

1 When you look at the pivotal study results
2 and the 12-month endpoint compared to open surgery,
3 the EBE is safe. We had 100 percent of the devices
4 successfully deployed and patent. There was faster
5 recovery; a striking reduction in major adverse
6 events; similar survival both overall and aneurysm
7 related. We had clinically effective aneurysm
8 exclusion. There were no conversions in 12 months
9 and no aneurysm ruptures.

10 [Slide]

11 In addition to the 12-month data showing
12 the safety and effectiveness, we continued to
13 follow these patients in the clinical trial, and I
14 have the privilege also to present to you the
15 24-month data today to answer the question are the
16 12-month study results sustained.

17 [Slide]

18 This was rigorous and diligent follow-up
19 at two years. Again, of the patients available for
20 follow-up, we had over 90 percent compliance with
21 the clinical visits in both treatment groups.

22 [Slide]

23 The survival curve going out to two years
24 and 93 percent of the patients are still alive in
25 that control group, 87 percent in the EBE. There

1 is no significant difference in these two.

2 [Slide]

3 Aneurysm-related survival, 98 percent in
4 the control group, 98 percent in the EBE group at
5 24 months, no significant difference in
6 aneurysm-related survival.

7 [Slide]

8 Endoleak results by the core lab are very
9 similar. The majority of the patients do not have
10 an endoleak. Of those who do have an endoleak--the
11 24-month data, by the way, is in grey; the 12-month
12 in yellow. The majority of the endoleaks are the
13 branch variety type. There is a small frequency of
14 patients who have type I endoleak and endoleak of
15 undetermined source. There were no type III or
16 type IV endoleaks.

17 [Slide]

18 Aneurysm size change at 24 months--again,
19 the 1-year data from 1-12 months is in gold or
20 yellow; the grey is 1-24 months, and 19 percent of
21 patients at 2 years have aneurysm size decrease of
22 5 mm or more. The majority of patients have no
23 change in aneurysm size and 14 percent of the
24 patients have aneurysm size increase of 5 mm or
25 more.

1 [Slide]

2 There were 11 patients who had 12
3 reinterventions in the second year. Nine of those
4 patients had catheter-based embolizations but I
5 want to talk about the three who had late
6 conversion to open repair. Two of these patients
7 had aneurysm enlargement with no endoleak
8 identified on preoperative imaging. One case had
9 an endoleak and aneurysm enlargement, and this
10 patient declined to have a catheter-based
11 embolization. All three were converted to open
12 surgical repair and were discharged. However,
13 there is one unfortunate patient who was readmitted
14 and actually died 24 days following the conversion
15 of endocarditis. There were no signs of graft
16 infection at the initial procedure with negative
17 cultures. Many of the details on these patients
18 are in your panel pack.

19 [Slide]

20 The other findings that the core lab is
21 looking at are for integrity, lumen narrowing and
22 migration, as shown here. In the 24-month
23 follow-up no other wire-form fractures have been
24 identified. A small percentage of patients have
25 radiographic evidence of device lumen narrowing,

1 trunk migration or component migration but, again,
2 there were no clinical consequences in any of these
3 patients. There were no type I or type III
4 endoleaks or aneurysm enlargement seen by the core
5 lab, and no clinical vascular adverse events or
6 reinterventions.

7 [Slide]

8 Again, this is our cumulative major
9 adverse events rates now extended out to 24 months
10 on the X axis. Again, the Y axis is the rate of
11 adverse events per patient. You have seen the left
12 half of this graph before. When you follow out the
13 cumulative major adverse events, they continue to
14 run essentially parallel. There is a persistent
15 difference, 1.9 cumulative major adverse events per
16 patient at two years in the control group compared
17 to 1.1 adverse events per patient in the EBE group.
18 Again, this is significantly different by the
19 Nelson.

20 [Slide]

21 So, to answer the question about the
22 24-month data, what does it show? Are the 12-month
23 results sustained? Yes, the 24-month data
24 substantiates the findings of the 12-month data.

25 [Slide]

1 In conclusion, in this presentation of the
2 pivotal study data, we feel that the EBE is safe
3 and effective for treatment of abdominal aortic
4 aneurysms and 100 percent of the devices were
5 successfully deployed and patent. There was faster
6 recovery. There is a striking and persistent
7 reduction in major adverse events. There is
8 similar survival both overall and aneurysm related.
9 There is clinically effective aneurysm exclusion
10 with rare conversions and no aneurysm ruptures.
11 Thank you for your attention.

12 DR. LASKEY: I would like to thank and
13 applaud this morning's presenters. You not only
14 stayed on time but you are a tad early. Therefore,
15 you have the privilege of responding to some early
16 questioning from the panel.

17 [Laughter]

18 So, let's just take five minutes. Are
19 there any burning questions from any of the panel
20 members before we break for lunch?

21 DR. BAILEY: Could I just ask a quick one?
22 I understand that the enrollment was parallel
23 groups, all surgical candidates and then, according
24 to their anatomy, they were divided into the two
25 assigned treatments. I think I saw that the

1 statistical plan was to have a two to one ratio
2 and, indeed, that seems to be close to what
3 happened. Was this by chance or was there some
4 mechanism to actually achieve this ratio? Or, is
5 that just the natural ratio that comes in the door?

6 DR. MATSUMURA: I guess I will answer that
7 question. The ratio is relatively two to one. The
8 clinical criteria for enrollment were very similar.
9 We didn't go through all those because they are in
10 your pack and in the protocol and are similar to
11 other trials. The only differences were where we
12 showed the infrarenal anatomy had to meet the
13 Excluder specifications, the EBE specifications in
14 the test group, and in the control group they had
15 to meet the criteria that an infrarenal clamp was
16 planned.

17 I think that if you look across sites,
18 there is a little bit of variation in the ratio of
19 two to one. I don't know if it is in the panel
20 pack but it is in the PMA, two to one. But sites
21 were told ahead of time that that was our
22 enrollment goal of two to one, and I suspect that
23 as they were seeing patients, you know, that they
24 had that in consideration.

25 I remember this meeting with the

1 investigators, and we had asked them not to be
2 enrolling in other trials during their enrollment
3 for this trial and to put all the patients who
4 qualify in. So, I think it is rather fortuitous
5 that that came out that way.

6 DR. BREWSTER: I think I would just add
7 from a real-world clinical perspective in terms of
8 the clinician interacting with the patient, once a
9 site had enrolled an adequate number of control
10 patients in the trial, I think the natural tendency
11 of an investigator or center would be to not
12 necessarily continue to enroll control patients
13 because there is a certain follow-up burden, and so
14 forth. So, once we felt at a particular site that
15 an adequate number of controls had been enrolled, I
16 think we probably ceased to enroll controls. That
17 probably fosters the difference as well.

18 DR. LASKEY: Tony?

19 DR. COMEROTA: Jon, that was a very
20 complete description of the results. My question
21 is not burning but one of curiosity. You mentioned
22 that proximal neck angle or increased proximal neck
23 angle was an independent risk factor and I am
24 presuming that that applies for the control group
25 as well the EBE group. Now, the protocol design,

1 of course, excluded proximal neck angle of greater
2 than 60 degrees in the EBE group. Is this new and
3 interesting information that we can carry away that
4 a neck angle increases risk for an operation in
5 patients with abdominal aortic aneurysms?

6 DR. MATSUMURA: I think it is new
7 information. I hope you don't carry it away
8 because one of our investigators has plans to
9 analyze that. But before this study was conducted
10 and analyzed we couldn't find any significant
11 anatomic predictors of risk in the literature, and
12 we did a fairly extensive search, which is in the
13 executive summary of the PMA and maybe in the panel
14 pack. But many people have looked at clinical risk
15 factors and, therefore, those are the ones that we
16 really wanted to stratify. So, in our analysis,
17 with Dr. Naftel's help, we did conduct this
18 analysis really for risk adjustment and we wanted
19 to include all the data that we had available, and
20 we had extensive data on anatomy that was really
21 very impressive in its detail and completeness. Of
22 all the anatomic variables we tried to throw in the
23 model, I think proximal neck angle only made it for
24 late adverse events, not early and not survival.

25 But I think it is interesting. I can

1 remember that David called me up and said, well,
2 why is this and what does a clinician think of
3 this? I guess I am not going to write that paper;
4 another sub-investigator is going to do that, but I
5 suspect that it has something to do with more neck
6 angulation being probably a marker for more
7 advanced disease and perhaps either surgeons treat
8 those patients because they have a larger aneurysm,
9 or maybe it is just a marker that goes with
10 something else about advanced disease. But I can't
11 imagine that the neck angle itself makes it harder.
12 The answer to your other question, is it applicable
13 to control and EBE, I believe it is. It is for
14 both groups.

15 DR. NAFTEL: I will just say that for all
16 the models, in addition to analyzing all the
17 patients together, in each case we applied the
18 models to just the control and then just to the EBE
19 to make sure there is no interaction and the
20 results were consistent, and they were in each
21 case.

22 DR. LASKEY: Ileana, one more question and
23 then we will break for lunch.

24 DR. PINA: Just out of probably sheer
25 ignorance, what do you do with anticoagulation? I

1 notice that some of the events that are labeled as
2 stats for other issues came out of CBAs and
3 peripheral embolization. How do you handle the
4 anticoagulation?

5 DR. BREWSTER: There was a clinical
6 adverse event committee that negotiated or
7 considered adverse events, identified by Dr.
8 Matsumura, the study primary investigator, which
9 included all identifiable adverse events. The
10 purpose was to better classify these, more
11 accurately classify these in order to clarify
12 reporting such as we have had this morning.

13 The initial study also had a rather large
14 category of so-called "other" events. Another
15 purpose of this adverse events committee, which met
16 fairly often and included the primary investigator
17 of the study itself, at least two site
18 investigators and a member of the data safety
19 monitoring board, reclassified these "other" events
20 into appropriate categories, again, to better
21 clarify reporting.

22 DR. PINA: Did you leave anticoagulation
23 up to the investigators or did you have a set
24 protocol? In other words, did the patients have to
25 be on Coumadin for X number of days? Did they have

1 to be on aspirin? I mean, all these people have
2 some sort of vascular disease.

3 DR. BREWSTER: There was no protocol
4 requirement in terms of postoperative
5 anticoagulation. That was left to whatever the
6 standard practice of the investigator might be. I
7 don't think any patients though were electively
8 anticoagulated in terms of Coumadin, for instance.
9 Many of them, no doubt, were put on aspirin. The
10 protocol, in terms of perioperative management, was
11 similar to standard open repair in that all
12 patients were advised to undergo perioperative
13 heparinization at the time of implant.

14 DR. LASKEY: Gentlemen, thank you. That
15 was a very articulate presentation. Thank you for
16 staying on time, and we will see you again at one
17 o'clock. I would like to adjourn for lunch.

18 [Whereupon, at 12:00 noon the proceedings
19 were recessed for lunch, to resume at 1:00 p.m.]
20

- - -

1 A F T E R N O O N P R O C E E D I N G S

2 DR. LASKEY: It is shortly after 1:00. We
3 should resume. Let's reopen the session with the
4 FDA presentation.

5 **FDA Presentation**

6 MR. GANTT: Good afternoon.

7 [Slide]

8 My name is Doyle Gantt. I am a senior
9 biomedical engineer reviewer and one of the lead
10 reviewers on this PMA application. Dorothy Abel is
11 the other lead reviewer on this application.

12 [Slide]

13 My presentation will include the
14 following, an introduction of the review team at
15 FDA; a summary of the FDA review; and the questions
16 for panel consideration. We had an opportunity to
17 see the sponsor's presentation prior to their
18 presentation this morning, and it accurately
19 summarizes the data reviewed by the agency so these
20 data will not be repeated in this presentation.

21 [Slide]

22 As with most PMAs like the one being
23 discussed today, review of the documents involves a
24 large number of reviewers that have provided
25 reviews in their areas of expertise. Included were

1 clinical, statistical, in vivo, in vitro--

2 [Slide]

3 --as well as biocompatibility, packaging,
4 sterilization, bioresearch monitoring,
5 manufacturing, QSR regulation and patient labeling,
6 and I would like to acknowledge all those
7 individuals who contributed to the review of this
8 application.

9 [Slide]

10 I would now like to begin with a summary
11 of the FDA review of the application.

12 [Slide]

13 First let's start with the preclinical.
14 FDA reviews of the biocompatibility, in vivo animal
15 studies, manufacturing and sterilization
16 information, including packaging and shelf-life,
17 have been completed and there are no issues
18 regarding these areas for the panel to discuss.

19 [Slide]

20 FDA review also included an assessment of
21 the device integrity and there are a number of
22 factors that I think we need to consider when
23 looking at this issue. First of all, as with other
24 stents used in the vascular system, endovascular
25 grafts may be subject to conditions which may

1 result in loss of device integrity.

2 [Slide]

3 Another factor, depending upon the
4 location and type of breach in integrity, there may
5 or may not be an immediate or eventual clinical
6 consequence. Another factor which must be
7 considered in review of this issue is the
8 difficulty in identifying and confirming the
9 presence of structural failures in vivo. The
10 sponsor didn't talk much about this in this
11 morning's presentation, but in review of the
12 failure analyses on this subject it became quite
13 clear that these things are very difficult to view
14 using standard x-ray techniques.

15 [Slide]

16 Prior to sending out the panel packs,
17 there were two reports of wire-form fractures
18 identified by the core laboratory, one at discharge
19 in a patient enrolled in the Phase II study, and
20 the other at 12 months in a patient enrolled in the
21 ongoing second generation device study. As was
22 mentioned by the sponsor this morning in their
23 presentation, a second generation device study has
24 been initiated to obtain data for a similar device.
25 Although this device is not the subject of this

1 PMA, the device is comparable enough from a
2 structural standpoint that we feel it is important
3 to consider this as part of the review of this
4 device as well.

5 [Slide]

6 Upon learning of these reports, the
7 sponsor did conduct a failure analysis and they
8 have communicated those findings to us. There have
9 been no adverse effects associated with either of
10 the two reports and there is not any conclusive
11 evidence to verify the presence or absence of the
12 fractures. As I mentioned earlier, they are very
13 difficult to visualize using x-ray.

14 [Slide]

15 Both of these reported fractures were
16 identified in the main body of the graft, not in a
17 seal zone or point of attachment to the aorta,
18 another factor that we believe is important in
19 considering the significance of the integrity
20 issue. The FDA review of the failure analysis of
21 these two reports has been completed, with no
22 additional information being requested of the
23 sponsor.

24 [Slide]

25 Finally, the sponsor has recently reported

1 a fracture in an explanted device. The fracture
2 was also located in the main body in the bifurcated
3 region of the device and there is very limited
4 information available at this time about this
5 particular report.

6 [Slide]

7 Now I would like to switch gears a little
8 bit and just go over a brief summary of the
9 clinical review that was conducted by FDA. This is
10 just an overview slide of the clinical study. As
11 was mentioned earlier, the pivotal study provided
12 primary safety and effectiveness data. As you
13 heard in this morning's presentations, this was a
14 non-randomized study with concurrent open surgical
15 control, consisting of patients who were not
16 eligible for treatment with endovascular graft due
17 to anatomical restrictions.

18 [Slide]

19 Some of the notable issues that we
20 addressed during the review of the clinical data
21 included the appropriateness of the non-randomized
22 study design; difficulty in enrolling patients,
23 primarily because this is a male dominated disease;
24 the number of, reasons for, and outcome of patients
25 converted to open surgical repair; clarification of

1 the rate of major adverse events after one month;
2 and clarification of the number of type I and III
3 endoleaks and aneurysm enlargements.

4 [Slide]

5 In summary, all the FDA requests for
6 additional information have been satisfied, and the
7 review team has identified the following questions
8 that we would like the panel to consider during
9 their discussion of this application.

10 [Slide]

11 Question number 1, the primary safety
12 endpoint of the clinical study was the rate of
13 major complications as evaluated through 12 months.
14 Additionally, data were presented for individual
15 adverse effects, analyses were provided for risk
16 factors associated with adverse events, and causes
17 of death are provided. A summary of the 24-month
18 results is also included. Please comment on
19 whether the results of the clinical study provide
20 reasonable assurance of safety in the intended
21 population.

22 [Slide]

23 Question number 2, the primary
24 effectiveness endpoint of the clinical study was
25 exclusion of the infrarenal abdominal aortic

1 aneurysm from the blood circulation defined by
2 absence of aneurysm enlargement and endoleaks, as
3 evaluated through 12 months. Additionally, data
4 regarding potential problems associated with
5 endovascular treatment, for example migration,
6 aneurysm enlargement, endoleaks, ruptures,
7 conversion, device integrity, are presented. A
8 summary of the 24-month results is also included.
9 Please comment on whether the results of the
10 clinical study provide reasonable assurance of
11 effectiveness in the intended population.

12 [Slide]

13 Number 3, the core laboratory has reported
14 two cases of wire-form fractures, one identified at
15 discharge in a patient enrolled in the pivotal
16 clinical study, and the other at 12 months in a
17 patient enrolled in the ongoing second generation
18 device study. There have been no adverse events
19 associated with either report and there is not
20 conclusive evidence to verify the presence or
21 absence of the fractures. Both reported fractures
22 were identified in the main body of the graft, not
23 in a seal zone or point of attachment to the aorta.

24 [Slide]

25 As a continuation, after the panel packs

1 were sent to the panel, the sponsor reported a
2 wire-form fracture which was recently identified
3 during the sponsor's analysis of a device explanted
4 in Germany. Details concerning the length of
5 implantation, implanting physician identity, and
6 device lot and serial numbers remain unavailable.
7 Based on the sponsor's analysis, it appears that
8 the fracture, which was also located in the main
9 body of the graft in the crotch of the bifurcation,
10 did not result in any clinical complications and
11 the ends of the wire did not appear to be
12 protruding through the device material or the
13 surrounding tissue. Please comment on the
14 significance of these observations.

15 [Slide]

16 One aspect of the premarket evaluation of
17 a new product is the review of its labeling. The
18 labeling must indicate which patients are
19 appropriate for treatment, identify potential
20 adverse events with the use of the device, and
21 explain how the product should be used to maximize
22 clinical benefit and minimize adverse events.

23 [Slide]

24 If the panel recommends approval for this
25 device, then we would like the panel to address the

1 following questions concerning the label.

2 Does the indication for use, as stated
3 below, adequately define the patient population
4 studied, and for which the device will be marketed?

5 The Excluder Endoprosthesis is intended to
6 exclude the aneurysm from the blood circulation in
7 patients diagnosed with infrarenal AAA disease who
8 have appropriate anatomy.

9 As a point of reference, we included an
10 addendum to the panel questions that were sent out
11 in the panel packs. That addendum includes the
12 indications for use statement for each of the
13 currently approved endovascular devices used in the
14 treatment of AAA. For convenience, I have a series
15 of slides that we can project during the panel
16 discussion to make that discussion a little bit
17 easier.

18 [Slide]

19 The second question related to the label,
20 based on the clinical investigation experience, are
21 there any additional warnings, precautions, or
22 contraindications that you think should be
23 included, either specific to this device or from a
24 generic standpoint for endovascular grafts?

25 Again, I have a series of slides, which

1 was an addendum which was included in the panel
2 pack, that describes the proposed label that the
3 company has provided to us concerning the warnings,
4 precautions and contraindications and we can
5 project those if there is a need during the panel
6 discussion to see the proposed label.

7 [Slide]

8 The third question related to the label,
9 please comment on whether the instructions for use
10 adequately describe how the device is to be
11 delivered.

12 [Slide]

13 Finally, do you have any other comments on
14 the label?

15 [Slide]

16 Question 5 is please comment on the
17 adequacy of the proposed physician training plan,
18 as described in the panel package.

19 [Slide]

20 Finally, the sponsor is proposing to
21 conduct a post-approval study on the patients
22 enrolled in the pivotal clinical study, that is, it
23 started with 235 test patients and 99 controls.
24 Five-year follow-up on all patients who are alive
25 and not withdrawn from the study will be obtained

1 in accordance with the clinical protocol approved
2 under the IDE for this device. Please comment on
3 the acceptability of this plan, as briefly
4 described in the panel package. As one final note
5 on that matter, this is very consistent with the
6 five-year post-approval studies being conducted by
7 the other approved devices that are on the market.

8 With that, I will end my presentation and
9 open it up for questions if there are questions of
10 me, or if you would like to get started with the
11 panel discussion that could happen as well.

12 DR. LASKEY: Does anybody have any
13 questions for the presenter at this point? Dr.
14 Pina?

15 DR. PINA: Thank you for your
16 presentation. In your review you have a paragraph
17 about the adjunctive procedures that were needed
18 during the implementation. How does that compare
19 to other similar devices on the market as far as
20 percentage of adjunctive procedures that are needed
21 at the time of implantation?

22 MR. GANTT: I am not sure of the response
23 to that question. I might ask one of the other
24 reviewers. Paul? Our clinical reviewer, Paul
25 Chandeysson.

1 DR. CHANDEYSSON: My name is Paul
2 Chandeysson. The rate of adjunctive procedures is
3 relatively low for this type of device, seven
4 percent.

5 DR. WHITE: Could I ask you another
6 question before you leave? I was interested, and I
7 am confused when I read the panel data, when you
8 talked about the number of audited core laboratory
9 images, specifically you talked about 155 paired CT
10 images for aneurysm growth. Are you familiar with
11 that part? Where did you get that information? I
12 couldn't find that kind of audited information in
13 the PMA. Where is the number of exams that were
14 actually done to look at aneurysm growth? Is that
15 somewhere in the PMA?

16 DR. CHANDEYSSON: I thought that was
17 somewhere because that is where those numbers come
18 from. It is possible it is an incorrect number but
19 I thought that the number of paired CT studies was
20 there.

21 DR. WHITE: The reason it is important is
22 that the denominator becomes the frequency of the
23 growth. So, your number of 155 paired studies is
24 the only 155 I can find in the submission. Maybe
25 the sponsor can help or the core laboratory can

1 help with that because that denominator number is
2 going to end up being crucial in deciding what was
3 the percentage of aneurysms that grew.

4 MR. GANTT: I might add something here.
5 Keep in mind that you have an annotated version of
6 the application. We have sent you a condensed
7 version of the entire submission as the panel pack.
8 I don't know if the sponsor wants to comment
9 further.

10 DR. LASKEY: Well, not at this point; we
11 will get to that when we have them come to the
12 table, but if there are questions for you
13 specifically from the panel. Tony?

14 DR. COMEROTA: Will you address the
15 statistical analysis either on safety or efficacy?

16 MR. GANTT: We didn't bring the
17 statistician with us for that part of review, but
18 if there are some general questions about the
19 statistical review we would be happy to respond to
20 that.

21 DR. WHITE: Well, most of us on the panel
22 are not statisticians and the data as presented,
23 from a clinician's perspective, looks reasonably
24 compelling. Yet, there were some questions raised
25 by the statistician, not regarding safety, of

1 course, but regarding efficacy and I just wanted to
2 have that addressed if it were possible.

3 MR. GANTT: I believe the only thing that
4 came up that was somewhat controversial in nature
5 during the review, the statistical review, was one
6 of the primary effectiveness endpoints and whether
7 or not we would be able to allow a particular claim
8 regarding the effectiveness of the device. Paul,
9 do you have any further information about that and
10 how we decided to resolve that?

11 DR. CHANDEYSSON: I think the issue was
12 about whether the point value of the effectiveness,
13 which was something above 80 percent, was
14 sufficient or whether the lower 95 percent
15 confidence interval would have to be considered,
16 and that was below the projected 80 percent. That
17 was the issue.

18 DR. LASKEY: If I am not mistaken, there
19 was also a very important point made about
20 surrogates that is worth emphasizing either now or
21 later when we get to it. I must say just as a
22 point of procedure, Dr. Zuckerman, we usually have
23 a short little précis presented by the FDA
24 statistician. We are just not having that today,
25 and there are a few items of contention that it

1 would be worthwhile having--was it Gerry Gray? I
2 forget who did this.

3 DR. ZUCKERMAN: You know, those points are
4 noted, but Dr. Bailey is here to help us out.

5 DR. BAILEY: Yes, but I don't know anymore
6 about what Dr. Kamer wrote. So, in his absence we
7 will have to ask you guys.

8 DR. LASKEY: Well, that may fall out of
9 the discussion. One more question, Ileana, and
10 then we will move on to the primary reviewers.

11 DR. PINA: I am a little concerned about
12 the deaths. I have been through each and every one
13 of them that you listed in the packet here, and
14 some of them that are listed as being pneumonia or
15 sepsis are actually the result of a cardiac event
16 and I counted several sudden deaths that were not
17 classified as sudden deaths but if I looked at the
18 history and I were sitting on an adjudication
19 committee, they would be sudden death for whatever
20 the etiology.

21 So, I am a little concerned about the
22 cardiovascular risk here. I mean, these are
23 patients who have extensive vascular disease and I
24 am not entirely surprised, but I think they should
25 be called what they are. Pneumonia is secondary to

1 the patient having been admitted with an arrest.
2 There are also several CVAs, which is the reason I
3 was asking about the anticoagulation protocol
4 because if we are going to sit here and make
5 recommendations and there are CVAs involved, and
6 these patients have cruddy aortas and manipulation
7 in there is going to, you know, let loose some
8 stuff, I am very, very concerned about that.

9 There are obviously the cancers and all
10 those that are way, way out, but some of these
11 occurred within a month or two, the first
12 event--unstable angina; there are some myocardial
13 infarctions and I counted three or four sudden
14 deaths. There are a couple of endocarditis. There
15 is a pericarditis that actually sounds more like
16 endocarditis than pericarditis. And these are
17 things that we need to think about if we approve
18 and when we are making the recommendations.

19 DR. ZUCKERMAN: I would like to make one
20 correction. Although Mr. Gary Kamer isn't here to
21 help us interpret the FDA's statistical review, Dr.
22 Gerry Gray will be here this afternoon. Dr. Gray
23 is our team leader for cardiovascular stats and he
24 can answer questions that are brought up by the
25 review done by Mr. Kamer.

1 DR. LASKEY: Thank you. Let's move on to
2 have the panel present their discussion. The
3 reviewers are Drs. Najarian and Comerota. Why
4 don't we begin with Dr. Comerota? May we have the
5 sponsor and their representatives step forward,
6 please?

7 **Open Committee Discussion**

8 *DR. COMEROTA: Well, thank you very much.
9 I will begin by congratulating the presenters for
10 very elegant presentations and a review of the
11 data, and also thank the reviewers for the FDA for
12 a very complete summary, at least from my
13 perspective.

14 In terms of the background of abdominal
15 aortic aneurysms, I think perhaps Dr. Brewster's
16 look at the risk of rupture was slightly
17 pessimistic in terms of rupture, at least from the
18 smaller sizes of the aneurysm. But I think what is
19 also true is that despite attention to this entity
20 the death rate from ruptured aneurysms over the
21 last two to three decades has not diminished
22 despite improvement in patient care when it is
23 offered. Certainly, an operation is the standard
24 by which all treatment modalities are to be judged.

25 We all recognize the advances in

1 technology and the zeal for endovascular repair of
2 abdominal aortic aneurysms, as well as the pressure
3 from the lay public to have this offered to them.
4 So, certainly we can understand the reason for the
5 design of the protocol being non-randomized, which
6 certainly is a criticism of the protocol by any who
7 might view it. Nonetheless, it certainly is
8 understandable.

9 The device description and its delivery
10 and the technique was well summarized, and the
11 manufacturer recommends oversizing of somewhere
12 between 10-21 percent of the graft to the aortic
13 attachment, and somewhere between 5 and 26 percent
14 for the iliac attachments.

15 The proposed indication was reviewed and
16 the endoprosthesis, the EBE, was recommended for
17 patients with appropriate anatomy. Just to
18 redefine that, appropriate anatomy is an aortic
19 neck of 1.5 cm or more, an angle of the aortic neck
20 of 60 degrees or less and, of course, ilio-femoral
21 morphology compatible with access and the delivery
22 system and, of course, no thrombus at the aortic or
23 the iliac implantation sites which might compromise
24 the seal of the graft to the iliac artery
25 interface.

1 The manufacturers presented their
2 feasibility study which was performed in 30
3 patients, 28 men and two women. Then, following
4 the feasibility study they proceeded to their
5 controlled clinical trial. During follow-up of the
6 feasibility study endoleaks were detected in 21
7 percent of the patients at three months; 25 percent
8 at 6 months; and 20 percent of the patients at 12
9 months. An increase in size of the aneurysms in
10 those patients in the feasibility study by 5 mm or
11 more was observed in 17 percent by the 12-month
12 follow-up, and that is by way of background of
13 course.

14 One patient ruptured the aneurysm, a
15 77-year old woman who was essentially treated in
16 violation of the study protocol in that her aortic
17 neck angle measured 90 degrees at the time of
18 implantation of the endograft. She had an endoleak
19 identified. The recommendations were that she be
20 converted to open repair. She refused any
21 conversion to open repair and subsequently
22 ruptured, I believe, three years later.

23 Another interesting patient is a 75-year
24 old gentleman who demonstrated aneurysm growth at
25 36 months follow-up. The patient was converted to

1 open repair. There were elevated pressures
2 measured in the aneurysm sac, however, no endoleak
3 was identified at the time of open conversion. An
4 interesting observation, in my opinion, was that
5 there was serous fluid which was drained from the
6 sac. The graft was intact. Gross and fluoroscopic
7 examination of the explanted device did not reveal
8 any perforation or fracture or device failure.

9 The last death was a 75-year old gentleman
10 hospitalized with evidence of sepsis two months
11 after the endograft was placed. Blood cultures
12 were positive for Staph. aureus. The infected
13 graft was removed successfully. The patient
14 withdrew from the study during follow-up.

15 The control of the clinical trial was
16 performed and the results were reviewed with us
17 this morning in elegant fashion. The pivotal trial
18 was a concurrently controlled clinical trial, not
19 randomized. The patients met the inclusion and
20 exclusion criteria, and those being suitable
21 candidates for open repair of their aneurysm with
22 intended placement of the aortic clamp in an
23 infrarenal location.

24 I will address the underlying assumption
25 that the complications of surgery were related to

1 the medical condition of the patient rather than
2 the anatomy of the aortic aneurysm. It, therefore,
3 appeared that a non-randomized control group such
4 as this would offer reasonably valid comparison for
5 the test group. I do challenge that underlying
6 assumption that the anatomy of the aneurysm is not
7 important in terms of an associated risk factor
8 because, as it turned out, 11 percent of the
9 control patients had suprarenal clamping of their
10 abdominal aortic aneurysm during repair and I think
11 most of us would agree that clamping of the aorta
12 above the renal arteries would be associated with a
13 high complication rate than infrarenal clamping in
14 most centers.

15 There was a required sample size that was
16 calculated upon the assumption that there would be
17 a 10 percent complication rate in the test group
18 and 25 percent complication rate in the control
19 group. Then efficacy measures were calculated
20 based upon the sample size.

21 That comes into play in terms of the
22 efficacy analysis that was performed by the
23 statistician and the conclusions from that
24 statistical efficacy analysis, which I guess we
25 should read into the record for completeness sake.

1 In terms of the risk of the surgical
2 patients, I would just indicate that in my opinion
3 I think the surgical patients in this trial were at
4 increased risk compared to the non-operated
5 patients, and that more patients were symptomatic
6 in the surgical group. There were more females in
7 the surgical group. The anatomic considerations
8 were as I reviewed, with 11 percent of them having
9 suprarenal clamping. This morning, in Dr.
10 Matsumura's presentation, we learned that an
11 increasing angle of the aortic neck was probably a
12 long-term risk factor and, of course, by definition
13 we had a greater angle in the surgical patients
14 than the endograft patients. So, the assumption
15 that anatomy is not an important risk consideration
16 I believe is not particularly valid.

17 In the study, I think we have to
18 compliment the investigators both endovascularly
19 and surgically on achieving an exceptionally low
20 operative mortality rate, one percent 30-day
21 mortality rate in those patients undergoing
22 endovascular repair and, as you saw, a zero percent
23 30-day mortality rate in the surgical group. But I
24 don't think we, as surgeons, would agree that a
25 zero percent mortality rate at 30 days means a zero

1 percent operative mortality. If we are critical of
2 ourselves, we realize that there is a two percent
3 operative mortality. Two patients died after the
4 30-day window but they did not survive the
5 hospitalization for the aneurysm repair and died as
6 a direct cause of complications that were
7 experienced during their operation. That is an
8 important consideration, of course.

9 In terms of the safety data, the principal
10 safety analysis looks very favorable. The safety
11 data were a comparison of the number of patient
12 deaths, as well as other adverse events. I
13 mentioned the 30-day mortality. The early adverse
14 events commonly observed in the control group are
15 calculated as 57 percent, and in the excluded group
16 as 14 percent which, of course, was highly
17 statistically significant. Interestingly and
18 surprisingly to me, there were no open conversions
19 reported before 24 months, another remarkable
20 observation.

21 Three patients had conversions after the
22 24-month time period due to aneurysm enlargement,
23 and no leak was found at the time of the open
24 conversion in those patients. The observation of
25 serous drainage or serous fluid within the aneurysm

1 sac in at least two of those three open conversions
2 I thought was an interesting observation and raises
3 a future question.

4 Follow-up CT scans regarding trunk
5 migration, regarding component migration
6 demonstrated an exceptionally low rate of true
7 migration and component migration.

8 There were significantly fewer major
9 adverse event rates in the Excluder group compared
10 to the control group. The specific events that
11 were individually significantly reduced are
12 bleeding, pulmonary complications, cardiac
13 complications and gastrointestinal complications,
14 as were reviewed this morning.

15 The overall death rates I think we need to
16 be cognizant of because during not only the
17 one-year but during the entire follow-up period
18 there was a 14 percent death rate in the endograft
19 group compared to about a 17 percent death rate in
20 the control group, and the overwhelming, if not all
21 of those deaths, were not directly related to an
22 aneurysm cause, aneurysm etiology or intervention
23 for their aneurysm.

24 In terms of effectiveness of the Excluder
25 Endoprosthesis in the management of aortic aneurysm

1 patients, there was 100 percent delivery rate. In
2 68 percent of the patients only the trunk
3 ipsolateral and contralateral limb components of
4 the device were required. The aortic extender was
5 used in seven percent and one or more iliac
6 extenders were used in 23 percent. Only three
7 percent of the patients required both an aortic and
8 iliac extender prosthesis.

9 We heard the results of the core
10 laboratory reports regarding the presence of
11 endoleaks. The number of endoleaks was relatively
12 small, especially compared to other devices
13 currently available. The number of type I
14 endoleaks were preciously small by my observation.

15 In terms of secondary outcomes, it appears
16 that secondary outcomes demonstrate significant
17 benefit in the EBE group compared to the control
18 group, this being reduction in blood loss; the
19 reduction in the transfusion requirements;
20 significantly more rapid procedure time and
21 decreased length of ICU stay, reduced hospital
22 stay; and quicker time to ambulation and recovery
23 to patients' normal activities.

24 These observations are not in isolation.
25 The European experience was reported. At the time

1 of the report that we received 234 patients were
2 treated with the Excluder device and were entered
3 into the EuroStar Registry from 33 centers. There
4 were no conversions to surgery and there was no
5 operative mortality in the patients entered into
6 the EuroStar Registry. There was one potential
7 aneurysm rupture several months following the
8 deployment of the Excluder graft, and the reported
9 rate of endoleak at 12 months in the EuroStar
10 Registry was 11 percent.

11 So, from my perspective, it appears that,
12 indeed, the EBE device met the requirements of
13 safety. In terms of efficacy, from a clinician's
14 perspective, it appears that it met the
15 requirements for efficacy. However, as I reviewed
16 the statistician's report, the FDA statistician's
17 report, it did not meet the statistical
18 requirements for efficacy based upon the a priori
19 effectiveness goal set by the manufacturer of a
20 rate of at least 80 percent. As I read the
21 statistical analysis, that has to do with the
22 confidence interval being somewhere between 77
23 percent and 95 percent rather than 80 and 95
24 percent. And, I am going to leave the statistical
25 argument up to the rest of the panel, not must

1 myself, but I thought it important to read that
2 into the record.

3 I would also make the point that while
4 conversion to open surgery has been low, the
5 majority of those patients who were converted with
6 this particular device demonstrated an intact graft
7 with no endoleak but with clear or serous fluid
8 within the sac. I wonder if this is a unique
9 property of the PTFE itself in terms of either
10 allowing some serous drainage or promoting that
11 type of response from the surrounding tissues, and
12 perhaps this is something that we can address as a
13 panel.

14 In terms of summary and my conclusions, I
15 believe the sponsor of the Exclude Bifurcated
16 Endoprosthesis has reported their data in a rather
17 complete fashion including non-randomized
18 controlled clinical trial of the Excluder versus
19 conventional open surgical repair. It demonstrated
20 significantly lower morbidity than conventional
21 surgery for infrarenal abdominal aortic aneurysms.
22 The mortality was very low in both groups and not
23 different. The clinical utility endpoints such as
24 blood loss, blood transfusions and those that I
25 have summarized were significantly lower in the

1 Excluder Endoprosthesis group. So the device is
2 safe and the rate of successful implantation is
3 enviably high.

4 The effectiveness of the Excluder is
5 measured by subsequent aneurysm rupture. That is
6 very high. Only one patient suffered a ruptured
7 aneurysm subsequent to attempted endograft
8 placement but not in the controlled trial. That
9 patient, as I mentioned, refused conversion. So,
10 compared to other devices on the market the
11 endograft treatment of abdominal aortic aneurysm
12 with the Excluder appears to offer excellent safety
13 and effectiveness with good durability.

14 I would also raise just one final question
15 in terms of current data regarding the management
16 of patients with aneurysms less than 5.5 cm. Of
17 course, we are aware that two randomized trials
18 have been published since the initiation of this
19 trial, demonstrating that elective intervention of
20 the "small" aneurysm demonstrated no benefit
21 compared to careful surveillance. That may become
22 an issue for subsequent management of all types of
23 patients with aneurysms in the future, or may
24 become an issue in terms of future trials looking
25 at less invasive methods of management of patients

1 with smaller aneurysms compared to the natural
2 history of those patients undergoing careful
3 surveillance.

4 Mr. Chairman, that is my report. Thank
5 you.

6 DR. LASKEY: Tony, do you have any
7 specific questions you wanted to ask of the
8 presenters today?

9 DR. COMEROTA: Well, one is the
10 observation of that serous fluid within the sac in
11 those who were converted that demonstrated no
12 evidence of endoleak although the aneurysm was
13 enlarging and there was increased pressure within
14 the sac.

15 The other question, and I suspect that it
16 may be a bit unfair to pose it, is in terms of are
17 there going to be recommendations of management
18 based upon size? And, will there be different
19 considerations of intervention for an
20 endoprosthesis compared to standard open repair for
21 the patient with the smaller, i.e., less than 5.5
22 cm abdominal aortic aneurysm?

23 DR. MATSUMURA: Dr. Comerota, I think I
24 will start with a response to the first question
25 regarding the serous observation. I don't recall

1 the specific question but I think you are asking
2 what do we think about that. You talked about four
3 patients with conversion, and I just want to
4 clarify that one of those was in the feasibility
5 study and did not have an endoleak visible and had
6 aneurysm enlargement, and the conversion was about
7 three years later.

8 Three of the conversions that you spoke
9 about were in the pivotal study. One of those
10 patients did have a type II endoleak and refused a
11 coil embolization. Presumably, that may be related
12 to the growth. That patient, upon conversion, did
13 well.

14 DR. COMEROTA: At the time of the
15 conversion, however, the demonstration of the
16 endoleak--correct me if I am in error, the
17 demonstration of the endoleak was temporally
18 removed from the conversion by quite a period of
19 time, and at the time of conversion there was no
20 demonstration of endoleak. Is that correct?

21 DR. BREWSTER: Could you say that again,
22 Dr. Comerota?

23 DR. COMEROTA: If we are talking about the
24 same patient, I believe that there was a patient
25 who had an endoleak demonstrated early during the