

3836 '00 MAR 13 AM 106

APPROVAL ORDER

OOM-0809

AAV /



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Sherri K. Lakota
Manager, Regulatory Affairs
Alcon Laboratories
6201 South Freeway
Fort Worth, Texas 76134

FEB 24 2000

Re: P990023
Cellugel® Ophthalmic Viscosurgical Device (OVD)
Filed: May 7, 1999
Amended: June 1, July 23, August 3 and August 13, September 3 and September 9, November 23, December 10 and December 14, 1999; February 14, February 16 and February 17, 2000.

Dear Ms. Lakota:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has completed its review of your premarket approval application (PMA) for the Cellugel® Ophthalmic Viscosurgical Device (OVD). This device is indicated for use during surgery in the anterior segment of the eye. It is designed to create and maintain space, to protect the corneal endothelium and other intraocular tissues and to manipulate tissues during surgery. It may also be used to coat intraocular lenses and other instruments during cataract extraction and IOL insertion. We are pleased to inform you that the PMA is approved subject to the conditions described below and in the "Conditions of Approval" (enclosed). You may begin commercial distribution of the device upon receipt of this letter.

The sale, distribution, and use of this device are restricted to prescription use in accordance with 21 CFR 801.109 within the meaning of section 520(e) of the Federal Food, Drug, and Cosmetic Act (the act) under the authority of section 515(d)(1)(B)(ii) of the act. FDA has also determined that, to ensure the safe and effective use of the device, the device is further restricted within the meaning of section 520(e) under the authority of section 515(d)(1)(B)(ii) insofar as the sale, distribution, and use must not violate sections 502(q) and (r) of the act.

Expiration dating for this device has been established and approved at one year. The device should be stored at room temperature (15° – 30°C) and should be protected from light and freezing. This is to advise you that the protocol you used to establish this expiration dating is considered an approved protocol for the purpose of extending the expiration dating as provided by 21 CFR 814.39(a)(7).]

CDRH will notify the public of its decision to approve your PMA by making available a summary of the safety and effectiveness data upon which the approval is based.

The information can be found on the FDA CDRH Internet HomePage located at <http://www.fda.gov/cdrh/pmapage.html>. Written requests for this information can also be made to the Dockets Management Branch, (HFA-305), Food and Drug Administration, 5630 Fishers Lane, Rm. 1061, Rockville, MD 20852. The written request should include the PMA number or docket number. Within 30 days from the date that this information is placed on the Internet, any interested person may seek review of this decision by requesting an opportunity for administrative review, either through a hearing or review by an independent advisory committee, under section 515(g) of the Federal Food, Drug, and Cosmetic Act (the act).

Failure to comply with the conditions of approval invalidates this approval order. Commercial distribution of a device that is not in compliance with these conditions is a violation of the act.

You are reminded that, as soon as possible and before commercial distribution of your device, you must submit an amendment to this PMA submission with copies of all approved labeling in final printed form. As part of our reengineering effort, the Office of Device Evaluation is piloting a new process for review of final printed labeling. The labeling will not routinely be reviewed by FDA staff when PMA applicants include with their submission of the final printed labeling a cover letter stating that the final printed labeling is identical to the labeling approved in draft form. If the final printed labeling is not identical, any changes from the final draft labeling should be highlighted and explained in the amendment. Please see the CDRH Pilot for Review of Final Printed Labeling document at <http://www.fda.gov/cdrh/pmat/pilotpmat.html> for further details.


All required documents should be submitted in triplicate, unless otherwise specified, to the address below and should reference the above PMA number to facilitate processing.

PMA Document Mail Center (HFZ-401)
Center for Devices and Radiological Health
Food and Drug Administration
9200 Corporate Blvd.
Rockville, Maryland 20850

Page 3 -Ms. Sherri K. Lakota

If you have any questions concerning this approval order, please contact Mr. Lawrence J. Romanell at (301) 594-2053.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Philip J. Phillips". The signature is fluid and cursive, with the first name "Philip" being more prominent than the last name "Phillips".

Philip J. Phillips
Deputy Director for Science and
Regulatory Policy
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Issued: 3-4-98

CONDITIONS OF APPROVAL

APPROVED LABELING. As soon as possible, and before commercial distribution of your device, submit three copies of an amendment to this PMA submission with copies of all approved labeling in final printed form to the PMA Document Mail Center (HFZ-401), Center for Devices and Radiological Health, Food and Drug Administration (FDA), 9200 Corporate Blvd., Rockville, Maryland 20850.

ADVERTISEMENT. No advertisement or other descriptive printed material issued by the applicant or private label distributor with respect to this device shall recommend or imply that the device may be used for any use that is not included in the FDA approved labeling for the device. If the FDA approval order has restricted the sale, distribution and use of the device to prescription use in accordance with 21 CFR 801.109 and specified that this restriction is being imposed in accordance with the provisions of section 520(e) of the act under the authority of section 515(d)(1)(B)(ii) of the act, all advertisements and other descriptive printed material issued by the applicant or distributor with respect to the device shall include a brief statement of the intended uses of the device and relevant warnings, precautions, side effects and contraindications.

PREMARKET APPROVAL APPLICATION (PMA) SUPPLEMENT. Before making any change affecting the safety or effectiveness of the device, submit a PMA supplement for review and approval by FDA unless the change is of a type for which a "Special PMA Supplement-Changes Being Effectuated" is permitted under 21 CFR 814.39(d) or an alternate submission is permitted in accordance with 21 CFR 814.39(e). A PMA supplement or alternate submission shall comply with applicable requirements under 21 CFR 814.39 of the final rule for Premarket Approval of Medical Devices.

All situations which require a PMA supplement cannot be briefly summarized, please consult the PMA regulation for further guidance. The guidance provided below is only for several key instances.

A PMA supplement must be submitted when unanticipated adverse effects, increases in the incidence of anticipated adverse effects, or device failures necessitate a labeling, manufacturing, or device modification.

A PMA supplement must be submitted if the device is to be modified and the modified device should be subjected to animal or laboratory or clinical testing designed to determine if the modified device remains safe and effective.

A "Special PMA Supplement - Changes Being Effectuated" is limited to the labeling, quality control and manufacturing process changes specified under 21 CFR 814.39(d)(2). It allows for the addition of, but not the replacement of previously approved, quality control specifications and test methods. These changes may be implemented before FDA approval upon acknowledgment by FDA that the submission is being processed as a "Special PMA Supplement - Changes Being Effectuated." This acknowledgment is in addition to that issued by the PMA Document Mail Center for all PMA supplements submitted. This procedure is not applicable to changes in device design, composition, specifications, circuitry, software or energy source.

Alternate submissions permitted under 21 CFR 814.39(e) apply to changes that otherwise require approval of a PMA supplement before implementation of the change and include the use of a 30-day PMA supplement or annual postapproval report. FDA must have previously indicated in an advisory opinion to the affected industry or in correspondence with the applicant that the alternate submission is permitted for the change. Before such can occur, FDA and the PMA applicant(s) involved must agree upon any needed testing protocol, test results, reporting format, information to be reported, and the alternate submission to be used.

POSTAPPROVAL REPORTS. Continued approval of this PMA is contingent upon the submission of postapproval reports required under 21 CFR 814.84 at intervals of 1 year from the date of approval of the original PMA. Postapproval reports for supplements approved under the original PMA, if applicable, are to be included in the next and subsequent annual reports for the original PMA unless specified otherwise in the approval order for the PMA supplement. Two copies identified as "Annual Report" and bearing the applicable PMA reference number are to be submitted to the PMA Document Mail Center (HFZ-401), Center for Devices and Radiological Health, Food and Drug Administration, 9200 Corporate Blvd., Rockville, Maryland 20850. The postapproval report shall indicate the beginning and ending date of the period covered by the report and shall include the following information required by 21 CFR 814.84:

(1) Identification of changes described in 21 CFR 814.39(a) and changes required to be reported to FDA under 21 CFR 814.39(b).

(2) Bibliography and summary of the following information not previously submitted as part of the PMA and that is known to or reasonably should be known to the applicant:

(a) unpublished reports of data from any clinical investigations or nonclinical laboratory studies involving the device or related devices ("related" devices include devices which are the same or substantially similar to the applicant's device); and

(b) reports in the scientific literature concerning the device.

If, after reviewing the bibliography and summary, FDA concludes that agency review of one or more of the above reports is required, the applicant shall submit two copies of each identified report when so notified by FDA.

ADVERSE REACTION AND DEVICE DEFECT REPORTING. As provided by 21 CFR 814.82(a)(9), FDA has determined that in order to provide continued reasonable assurance of the safety and effectiveness of the device, the applicant shall submit 3 copies of a written report identified, as applicable, as an "Adverse Reaction Report" or "Device Defect Report" to the PMA Document Mail Center (HFZ-401), Center for Devices and Radiological Health, Food and Drug Administration, 9200 Corporate Blvd., Rockville, Maryland 20850 within 10 days after the applicant receives or has knowledge of information concerning:

(1) A mix-up of the device or its labeling with another article.

(2) Any adverse reaction, side effect, injury, toxicity, or sensitivity reaction that is attributable to the device and

(a) has not been addressed by the device's labeling or

(b) has been addressed by the device's labeling, but is occurring with unexpected severity or frequency.

(3) Any significant chemical, physical or other change or deterioration in the device or any failure of the device to meet the specifications established in the approved PMA that could not cause or contribute to death or serious injury but are not correctable by adjustments or other maintenance procedures described in the approved labeling. The report shall include a discussion of the applicant's assessment of the change, deterioration or failure and any proposed or implemented corrective action by the applicant. When such events are correctable by adjustments or other maintenance procedures described in the approved labeling, all such events known to the applicant shall be included in the Annual Report described under "Postapproval Reports" above unless specified otherwise in the conditions of approval to this PMA. This postapproval report shall appropriately categorize these events and include the number of reported and otherwise known instances of each category during the reporting period. Additional information regarding the events discussed above shall be submitted by the applicant when determined by FDA to be necessary to provide continued reasonable assurance of the safety and effectiveness of the device for its intended use.

REPORTING UNDER THE MEDICAL DEVICE REPORTING (MDR) REGULATION. The Medical Device Reporting (MDR) Regulation became effective on December 13, 1984. This regulation was replaced by the reporting requirements of the Safe Medical Devices Act of 1990 which became effective July 31, 1996 and requires that all manufacturers and importers of medical devices, including in vitro diagnostic devices, report to the FDA whenever they receive or otherwise become aware of information, from any source, that reasonably suggests that a device marketed by the manufacturer or importer:

(1) May have caused or contributed to a death or serious injury; or

(2) Has malfunctioned and such device or similar device marketed by the manufacturer or importer would be likely to cause or contribute to a death or serious injury if the malfunction were to recur.

The same events subject to reporting under the MDR Regulation may also be subject to the above "Adverse Reaction and Device Defect Reporting" requirements in the "Conditions of Approval" for this PMA. FDA has determined that such duplicative reporting is unnecessary. Whenever an event involving a device is subject to reporting under both the MDR Regulation and the "Conditions of Approval" for a PMA, the manufacturer shall submit the appropriate reports required by the MDR Regulation within the time frames as identified in 21 CFR 803.10(c) using FDA Form 3500A, i.e., 30 days after becoming aware of a reportable death, serious injury, or malfunction as described in 21 CFR 803.50 and 21 CFR 803.52 and 5 days after becoming aware that a reportable MDR event requires remedial action to prevent an unreasonable risk of substantial harm to the public health. The manufacturer is responsible for submitting a baseline report on FDA Form 3417 for a device when the device model is first reported under 21 CFR 803.50. This baseline report is to include the PMA reference number. Any written report and its envelope is to be specifically identified, e.g., "Manufacturer Report," "5-Day Report," "Baseline Report," etc.

Any written report is to be submitted to:

Food and Drug Administration
Center for Devices and Radiological Health
Medical Device Reporting
PO Box 3002
Rockville, Maryland 20847-3002

Copies of the MDR Regulation (FOD # 336&1336) and FDA publications entitled "An Overview of the Medical Device Reporting Regulation" (FOD # 509) and "Medical Device Reporting for Manufacturers" (FOD #987) are available on the CDRH WWW Home Page. They are also available through CDRH's Fact-On-Demand (F-O-D) at 800-899-0381. Written requests for information can be made by sending a facsimile to CDRH's Division of Small Manufacturers Assistance (DSMA) at 301-443-8818.

SUMMARY OF SAFETY AND
EFFECTIVENESS DATA (SSED)

SUMMARY OF SAFETY AND EFFECTIVENESS

I. General Information

- A. Device Generic Name: Intraocular Fluid (Hydroxypropyl methylcellulose)
- B. Device Trade Name: CELLUGEL®
Ophthalmic Viscosurgical Device
- C. Applicant's Name/Address: ALCON LABORATORIES, INC.
6201 South Freeway
Fort Worth, Texas 76134
- D. Premarket Approval (PMA) Application Number: P990023
- E. Date of Panel Recommendation: N/A
- F. Dates of Good Manufacturing Practice (GMP) Inspection:

07/02/1999
- Manufacturing Site: Alcon-Couvreur SA
Rijksweg 14
B-2870 Puurs, Belgium
- G. Date of Notice of Approval to the Applicant: FEB 24 2000

II. Indications for Use

CELLUGEL® is indicated for use during surgery in the anterior segment of the eye. CELLUGEL® is designed to create and maintain space, to protect the corneal endothelium and other intraocular tissues and to manipulate tissues during surgery. CELLUGEL® may also be used to coat intraocular lenses and instruments during cataract extraction and intraocular lens insertion.

III. Contraindications

At present, there are no known contraindications to the use of CELLUGEL® ophthalmic viscosurgical device (OVD) when used as recommended.

IV. Precautions

Precautions are limited to those normally associated with the surgical procedure being performed. As with all ophthalmic viscosurgical devices, a transient rise in intraocular pressure (IOP) in the early postoperative period has been reported in some cases. It is therefore recommended that CELLUGEL® be removed by the anterior chamber by thorough irrigation and/or aspiration at the end of surgery to minimize post-operative IOP increases. IOP should be monitored postsurgically and appropriate therapy instituted if significant increases occur. In addition to the above, the following precautions should be observed:

- Do not reuse cannulas.
- Use only if material is clear.
- Avoid trapping air bubbles within CELLUGEL® before injection.
- This product contains dry natural rubber.

V. Device Description

CELLUGEL® Ophthalmic Viscosurgical Device is a sterile, nonpyrogenic, noninflammatory viscoelastic solution of highly purified non-proteinaceous 2% hydroxypropyl methylcellulose (HPMC) with an average molecular weight of 300,000 daltons dissolved in an isotonic, physiological buffer.

Each milliliter of CELLUGEL® contains 2% HPMC, 0.525% sodium chloride, 0.075% potassium chloride, 0.048% calcium chloride, 0.03% magnesium chloride, 0.39% sodium acetate, 0.17% sodium citrate, and water for injection q.s.

The osmolarity of CELLUGEL® is 315 ± 35 mOsm/kg, the pH 7.2 ± 0.4 , and the viscosity $30,000 \pm 10,000$ cps (at 0.2 sec^{-1} , 25°C).

VI. Alternate Practices and Procedures

Prior to ophthalmic viscosurgical devices, air, gases, or irrigating solutions were utilized as anterior chamber maintainers and surgical aids. Numerous other ophthalmic viscosurgical devices exist today and have been on the market since 1983.

VII. Marketing History

CELLUGEL® has been marketed and sold internationally between the years 1991-1996 by Vision Biology, Inc., (VBI) in countries whose Ministries of Health have approved the sale. Alcon Laboratories, Inc. purchased the device from VBI in 1996.

CELLUGEL® has been CE Marked by Alcon Laboratories, Inc., under the Medical Device Directive in February 1999. Product has been marketed in the EU bearing the CE Mark beginning in May 1999. More than 130,000 syringes have been marketed since

1991. CELLUGEL® has not been removed from the market in any countries for reasons related to safety or effectiveness.

Table 1

VBI Approvals to Market CELLUGEL

COUNTRY	APPROVAL DATE
Sweden	8/26/88
Hong Kong	8/26/88
The Netherlands	8/26/88
Chile	8/26/88
Switzerland	8/26/88
Portugal	8/26/88
Peru	9/22/88
Denmark	1989
Andorra	1989

VIII. Potential Adverse Effects of the Device on Health

In Clinical Studies C-96-48 and C-98-22, adverse events were reported in patients receiving CELLUGEL® and in patients receiving the control substance (a commercially available sodium hyaluronate viscoelastic that has been on the market for at least five years). In the two clinical studies, a total of 348 patients received CELLUGEL® and a total of 344 patients received the control OVD. Adverse events occurring at a frequency $\geq 1\%$ are presented in Table 2. Adverse events occurring at a rate of $< 1\%$ are listed in the text following Table 2.

No patients were discontinued from C-98-22 and no patients were discontinued from C-96-48 due to a device-related adverse event.

Table 2

Ophthalmological Adverse Events Occurring at a Rate $\geq 1\%$

Observation	OVD	C-96-48 ^a		C-98-22 ^b	
		N	%	N	%
External Slit-lamp Observations ^c	Cellugel	110	55.3	3	2.0
	Control	87	44.2	1	0.7
Posterior Capsule Haze	Cellugel	94	47.2	13	8.7
	Control	87	44.2	13	8.8
Intraocular Slit-lamp Observations ^d	Cellugel	70	35.2		
	Control	80	40.8		
Macular Degeneration	Cellugel	43	17.1	17	11.4
	Control	34	17.3	20	13.6
Lid Observations ^e	Cellugel	34	17.1	2	1.3
	Control	35	17.9	1	0.7
Posterior Segment Observations ^f	Cellugel	26	13.1	2	1.3
	Control	24	12.2	4	2.7
Nd: YAG Posterior Capsulotomy	Cellugel	20	10.1		
	Control	11	5.6		
Dry Eye ^g	Cellugel	11	5.5	5	3.4
	Control	8	4.1	2	1.4
Iris Atrophy	Cellugel	10	5.0		
	Control	6	3.1		
Macular Edema	Cellugel	9	4.5	1	0.7
	Control	8	4.1	3	2.0
Secondary Glaucoma	Cellugel	6	3.0		
	Control	5	2.6		
Hyphema	Cellugel	5	2.5		
	Control	2	1.0		
IOL Repositioning or Exchange	Cellugel	4	2.0		
	Control	2	1.0		
IOP > 40 mmHg	Cellugel	3	1.5	6	4.0
	Control	2	1.0	8	5.4
Vitreous in the Anterior Chamber	Cellugel	2	1.0	1	0.7
	Control	4	2.0		
Endothelial Damage	Cellugel	2	1.0	3	2.0
	Control			3	2.0
Cells (AC Cells \geq grade 3)	Cellugel			2	1.3
	Control				

Corneal Edema (\geq grade 3)	Cellugel			1	0.7
	Control			2	1.4
Nd: YAG Anterior Synechiolysis	Cellugel	2	1.0		
	Control	1	0.5		
Retina Procedure	Cellugel	2	1.0		
	Control	4	2.0		
Lid Procedure	Cellugel	2	1.0		
	Control	3	1.5		
Conjunctival Cyst/Filament Removal	Cellugel	2	1.0		
	Control				
Subjective Complaints ^h	Cellugel			6	4.0
	Control			7	4.8

^a Clinical Study C-96-48 - Cellugel (N=199); Control (N=196 with one patient not returning for follow-up).

^b Clinical Study C-98-22 - Cellugel (N=149); Control (N=147).

^c Includes conjunctival injection, conjunctival hemorrhage, superficial punctate, keratitis, ecchymosis, arcus senilis, conjunctival chemosis, pinguecula, subconjunctival hemorrhage, hyperemia, conjunctival gape and corneal abrasion.

^d Includes corneal folds, Descemets folds, endothelial folds, striae, guttata, trace endothelial changes, cortical remnants, endothelial pigment, endothelial debris and microcystic corneal edema.

^e Includes blepharitis, dermatochalasis, lid edema, ptosis, collarettes, and chalazion.

^f Includes posterior capsular folds/wrinkling, retinal pigment epithelial changes and posterior vitreous detachment.

^g Includes poor tear film.

^h Includes foreign body sensation, ocular pain and diplopia.

Other ophthalmic adverse events considered unrelated to the use of the OVD and occurring among patients at a rate of <1% included: eye discomfort, IOL membrane, puritus, retinal hemorrhage, blurred vision, IOL repositioning with vitrectomy, removal of residual lens cortex and foreign body removal.

IX. Summary of Preclinical Studies

An extensive battery of toxicity studies have been performed with CELLUGEL Ophthalmic Viscosurgical Device to evaluate the safety of this material as an adjunctive device for use during intraocular surgery.

Toxicology testing was conducted in accordance with ISO 10993 and the draft ISO Viscoelastic standard (ISO/WD 15798.2). All tests were conducted in compliance with Good Laboratory Practices (21 CFR 58) regulations.

No evidence of cytotoxicity, hemolysis, sensitization, mutagenic potential, or ocular irritation was found in any test performed on CELLUGEL. The results of these studies are summarized:

Table 3
Toxicological Studies

Study Type	Study Title	Species	Dose	# Animals Per Dose	Outcome
Cytotoxicity	<i>In vitro</i> Cytotoxicity Study (Agar Overlay Method) in the L929 Mouse Fibroblast Cell Line	L929 Mouse Fibroblast Cell	0.1 mL	N/A	Noncytotoxic
	<i>In vitro</i> Cytotoxicity Study (MEM Elution) in the L929 Mouse Fibroblast Cell Line	L929 Mouse Fibroblast Cell	25% Test solution	N/A	Noncytotoxic
Mutagenicity	Ames Mutagenicity	<i>Salmonella typhimurium</i>	0.1 mL	N/A	No mutagenic effects
	<i>E. coli</i> Plate Incorporation Mutagenicity Assay	<i>E. coli</i>	100 mg/mL	N/A	No mutagenic effects
Single Dose Toxicity	Acute Intraperitoneal Toxicity in Mice	Mouse (non-Swiss Albino CFI derived)	6 mL/kg	10	6 mL/kg, not a toxic dose
	Acute Oral Toxicity in the Rat	Rat (Sprague Dawley)	5 g/kg	10	5 g/kg, not a toxic dose
Immuno-genicity	Dermal Sensitization Study (Maximization Method) in Guinea Pig	Guinea pig	Induction 0.1 ml intra-dermal Induction	10	No immuno-genic effects

Study Type	Study Title	Species	Dose	# Animals Per Dose	Outcome
			and challenge 0.3 mL topical		
	Systemic Antigenicity in the Guinea Pig	Guinea pig	10 mL/kg of a 25% test solution	6	No immunogenic effects
Hemolysis	<i>In vitro</i> Hemolysis Study (direct contact)	N/A	0.2 mL blood added to a 20% test solution	N/A	Nonhemolytic

Study Type	Study Title	Species	Dose	# Animals Per Dose	Outcome
Local Tolerance	Primary Skin Irritation Test in the rabbit	Rabbits (NZW)	0.5 mL	6	Nonirritating
	Intraocular Irritation Study In the Rabbit (with Tonometry and Specular Photography)	Rabbits (NZW)	0.15 mL anterior chamber	6	Nonirritating
	Intraocular Irritation Evaluation of CELLUGEL in rabbits	Rabbits (NZW)	0.1 mL anterior chamber injection	8	Nonirritating
	Intraocular Irritation Evaluation of CELLUGEL in Primates	Cynomolgous monkey	0.1 mL anterior chamber injection	4	Nonirritating
	Intraocular Irritation Evaluation of CELLUGEL in Rabbits	Rabbits (NZW)	0.1 mL posterior chamber injection	8	Nonirritating
	IOP and Ocular Irritation Evaluation	Rabbits (NZW)	0.1 mL anterior chamber injection	16 eyes	Nonirritating; No IOP raising potential

X. Summary of Clinical Studies

A. Overview of Clinical Investigations

A total of five clinical studies were conducted using CELLUGEL® during the period of 1989 to 1999. Of these, two studies are presented in key support of safety and efficacy (C-98-22 and C-96-48). The remaining three studies were conducted outside the United States by Vision Biology (VBI) and are supportive of a safe and effective product (C-99-23, C-99-24, C-99-25).

The two key studies presented in support of safety and efficacy are Protocol C-96-48, conducted by Vision Biology, Inc., and Protocol C-98-22, conducted by

Alcon Laboratories, Inc. Clinical Protocol C-96-48 was a controlled, randomized, multicenter study among 396 patients which was designed to demonstrate that CELLUGEL® was: 1) equally effective to the control in its ability to protect corneal endothelial cells and maintain the anterior chamber during surgery and 2) equivalent to the control in its effects on postoperative IOP. Based on the results from this study, Clinical Protocol C-98-22 was designed to specifically address IOP elevation during the expected peak period, 6 hours postsurgery, since this data had not been collected during the previous C-96-48 clinical trial.

In order to obtain a more accurate representation of CELLUGEL®'s effect on IOP during the early postoperative period, no prophylactic medications were administered to patients in the C-98-22 study prior to the 6-hour IOP measurement.

The following table gives an overview of the two clinical studies that were considered key to support safety and effectiveness. These two studies compared CELLUGEL® to a marketed sodium hyaluronate viscosurgical device.

Table 4

CELLUGEL Clinical Studies

Protocol Number	Countries	No. of Sites	Study Duration	Patient Follow-up	CELLUGEL Subjects	Control Subjects
C-96-48	United States	9	1/93-6/95	6 Months	199	197
C-98-22	United States, Canada	9	7/98-2/99	21 Days	149	147
TOTAL Key Studies^a	-	18	1/93-2/99	-	348	344

B. Patient Population and Accountability

1. Demographics

a. Clinical Study C-98-22

No statistically significant differences between CELLUGEL® and the control were found for gender, race, age category and iris color. The treatment groups were similar for mean age for all patients enrolled.

b. Clinical Study C-96-48

No statistically significant differences between CELLUGEL® and the control were found for gender, age category and mean age among all patients enrolled. Information on patient race was not collected in this study. However, all patients who were eligible for the study were included.

Table 5

Key Studies C-98-22/C-96-48 Patient Demographics

Key Studies	Treatment	Enrolled	Male	Female	Mean Age
C-98-22	CELLUGEL	149	53	96	71.7
	Control	147	51	96	73.5
	Total	296	104	192	72.6
C-96-48	CELLUGEL	199	87	112	70.8
	Control	197	80	117	72.1
	Total	396	167	229	71.4
TOTAL		692	271	421	71.9

2. Inclusion Criteria

Clinical Studies C-98-22 and C-96-48

The total study population included 692 patients (male or female), of any race, who were scheduled for the removal of a cataract with the implantation of an intraocular lens. In addition, Clinical Study C-96-48 allowed the inclusion of aphakic patients requiring secondary IOL implantation (1 patient).

3. Exclusion Criteria

a. Clinical Study C-98-22

Patients were excluded from this study if they had other planned surgical procedures or the planned use of an investigational intraocular lens. They were also excluded if they had glaucoma in either eye or ocular hypertension (IOP > 21 mmHg) in the operative eye. Patients with proliferative diabetic retinopathy or uncontrolled diabetes mellitus were excluded from this study, as well as patients with any abnormality

that prevented reliable Goldmann applanation tonometry. In addition, patients with lens pseudoexfoliation syndrome, previous ocular trauma to the operative eye, a history of chronic or recurrent inflammatory eye disease or a congenital ocular abnormality were excluded. Patients were also excluded if they had iris atrophy, significant endothelial guttata or corneal dystrophy.

b. Clinical Study C-96-48

Patients were excluded from participation if they had acute ocular infection or inflammation, chronic uveitis, iritis, iridocyclitis or rubeosis iritis, uncontrolled glaucoma, aniridia, proliferative diabetic retinopathy, iris atrophy, or systemic disease with ocular manifestations.

4. Patient Accountability

All patients who received the randomly assigned study device were evaluable for safety (Intent to Treat data set). A subset of the entire population was also used for some key analyses; this is the Per Protocol data set. The Per Protocol data set included those patients who met inclusion/exclusion criteria and complied with the protocol. In keeping with current standards for the analysis of clinical data, where a patient's fellow eye was enrolled into the study, the patient's second eye was removed from the Per Protocol data set (but remained in the Intent to Treat data set). In addition, a few patients experiencing significant vitreous loss during surgery were excluded from the Per Protocol data set (prior to revealing the treatment codes) as vitreous in the anterior chamber can elevate intraocular pressure and may confound the data. After breaking treatment code, it was observed that an equal number of patients in the CELLUGEL® and the control groups had been excluded.

In general for equivalence hypotheses, the Per Protocol analysis is a more conservative approach. Therefore, primarily Per Protocol analyses have been presented, where equivalence arguments have been made. However, in both studies, the Intent to Treat and Per Protocol data sets support the same conclusions.

Table 6

Key Studies C-98-22/C-96-48 Patient Accountability

Key Studies	Treatment	All Patients (Intent to Treat)	Per Protocol Patients	Patients with missed visits
C-98-22	CELLUGEL	149	140	0
	Control	147	140	0
	Subtotal	296	280	0
C-96-48	CELLUGEL	199	169	20
	Control	197	164	25
	Subtotal	396	333	45
TOTAL		692	613	45

a. Patients with Missed Visits

1) Clinical Study C-98-22

No patient missed a visit during this 3-week study.

2) Clinical Study C-96-48

Forty-five patients missed visits during the course of this 6-month study. The following table lists the reasons for the missed visits. This attrition rate is not unusual for a study of this size and duration. The numbers of patients who missed visits were similar for both CELLUGEL® and the control groups. Although intercurrent illnesses, including those leading to death, occurred in this study, adverse events were not collected for these patients.

Table 7

Reasons for Missed Visits (C-96-48)

Reason for Missed Visit	CELLUGEL		Control	
	No. of Patients	Percent	No. of Patients	Percent
Lost to follow-up	11	5.5%	4	2.0%
Noncompliant with protocol	7	3.5%	10	5.0%
Patient requested withdrawal	0	0%	6	3.0%
Illness	1	0.5%	0	0%
Died	1	0.5%	5	2.5%
Total Patients in Study (N)	199	10.0%	197	12.5%

C. Efficacy Results

1. Endothelial Cell Density

a. Clinical Study C-98-22

No endothelial cell density data were captured in this 21-day IOP study.

b. Clinical Study C-96-48


CELLUGEL® was similar to the control in its ability to protect corneal endothelial cells during cataract/IOL surgery.

Endothelial cell densities were measured by specular microscopy prior to and, again, 6 months following surgery in the C-96-48 study.

Endothelial cell density losses were similar for both CELLUGEL® and the control in both the Per Protocol and Intent to Treat analyses. The mean percent change in cell density from baseline to 6 months was not statistically different between groups.

Figure 1

Change in Endothelial Cell Density (cells/mm²) at 6 Months (Per Protocol)



In the Per Protocol analyses, mean endothelial cell losses, measured at 6 months, were 3.6% and 3.8%, respectively when CELLUGEL (n=138) and control (n=130) were used to maintain anterior and posterior chamber spaces during surgery. At 6 months, CELLUGEL patients lost an average of 119 cells/mm², while control patients had lost an average of 135 cells/mm².

In the Intent to Treat analyses, endothelial cell loss was slightly higher in both groups although, again, the difference between groups was not statistically significantly different (CELLUGEL, 4.1%, n=152; control, 4.4%, n=146). At 6 months, CELLUGEL patients, on average, had lost an average of 127 cells/mm², while control patients had lost an average of 151 cells/mm².

The calculation of mean cell density change from baseline was based upon patient eyes that had a density measurement at both the preoperative baseline visit and the six-month postoperative visit. Patients who discontinued from the study and patients who were missing either a baseline or a 6-month endothelial cell density measurement were therefore excluded from analysis of the mean change from baseline.

2. Anterior Chamber Maintenance

Efficacy was also evaluated by the viscoelastic's ability to maintain a deep anterior chamber depth during surgery. CELLUGEL® was equal to or better than the control in maintaining the anterior chamber depth during surgery as reported by the investigators.

a. Clinical Study C-98-22

The viscoelastic's ability to maintain anterior chamber depth was a subjective evaluation, which was reported by the surgeon. In the Per Protocol data set, CELLUGEL® maintained the anterior chamber depth in a statistically significantly larger proportion of patients than did the control ($p < 0.001$). During anterior capsulotomy, CELLUGEL® maintained the anterior chamber depth in 97.9% of patients compared to 78.6% of the control patients. Thirty (21.4%) shallow anterior chamber depths were reported in the control patients compared to 3 (2.1%) CELLUGEL® patients.

During phacoemulsification, CELLUGEL® maintained the anterior chamber depth in a statistically significantly larger proportion of patients than did the control ($p=0.005$). The anterior chamber was maintained in 99% of the CELLUGEL® patients compared to 92% of the control patients.

During IOL insertion, CELLUGEL® and the control performed similarly (97.1 % Cellugel vs. 92.1 % Control) at maintaining the anterior chamber depth ($p=0.109$).

b. Clinical Study C-96-48

In the Per Protocol data set, the viscoelastic maintained a normal anterior chamber depth in 99.4% of patient eyes in both CELLUGEL® and the control. Only one patient in each group was reported to develop a shallow anterior chamber depth during surgery.

D. Safety Results

1. Intraocular Pressure

a. Mean Intraocular Pressure

1. Clinical Study C-98-22

Clinical Study C-98-22 was designed to evaluate the postoperative IOP profile of CELLUGEL® compared to the control. Patients received either CELLUGEL® or the control during cataract surgery (phacoemulsification) with posterior chamber intraocular lens implantation. No prophylactic IOP-reducing medications were administered at surgery. Immediately following the 6-hour IOP measurement, physicians were allowed to administer IOP-reducing therapy if the IOP was ≥ 30 mmHg. At all subsequent visits, investigators were allowed to prescribe IOP-reducing therapy as needed. IOP-reducing therapies were administered to a similar number of patients in the CELLUGEL® ($n = 24$) and the control ($n = 22$) groups. With the exception of one control patient, all IOP-reducing therapies were discontinued the day following surgery.

CELLUGEL® and the control were statistically equivalent in their effects on postoperative intraocular pressure. This conclusion was based on a statistical test of noninferiority. At each visit, the upper 95% confidence limit for the mean difference in IOP between CELLUGEL® and the control was less than 3.5 mmHg for both the Per Protocol and Intent-to-Treat data sets. Table 8 presents the mean IOPs at each visit for the Per Protocol data set.

Table 8

Mean IOP (mmHg) and Mean IOP Change From Baseline (mmHg) by Visit for
Per Protocol (C-98-22)

Treatment		Visit				
		Baseline	6 Hour	24 Hour	Day 7	Day 21
CELLUGEL	Mean	15.80	22.99	19.29	15.54	15.52
	Std	2.57	8.35	4.91	3.34	3.09
	N	140	139	140	140	140
	Min	9	6	10	6	8
	Max	21	54	36	24	24
	Mean Chg.	-	7.22	3.64	-0.60	-0.98
Control	Mean	15.94	21.63	19.58	15.38	14.96
	Std	2.65	7.91	5.93	3.11	3.14
	N	140	139	140	138	140
	Min	10	2	5	7	8
	Max	21	50	40	24	28
	Mean Chg.	-	5.71	3.49	-0.26	-0.28
Difference (CELLUGEL-Control) ^b		-0.14	1.36	-0.29	0.16	0.56
Upper 95% Confidence Limit ^{ab}		0.47	2.43	0.77	1.23	1.62

^a A one-sided 95% confidence interval was constructed. CELLUGEL® is noninferior to the control if the upper 95% confidence limit is less than 3.5 mmHg.

^b Based upon the difference in Least Squares (LS) Means. The LSMeans may differ slightly from the Arithmetic Means.

2) Clinical Study C-96-48

The postoperative mean IOP results from Study C-98-22 are supported by the mean IOP data from a subpopulation in the C-96-48 study that was similar to the population of C-98-22 (non-glaucoma patients without prophylactic IOP therapy at surgery). There were no statistical differences between the mean IOPs at all visits for these patients (Table 9). [T-test of the largest difference (0.6 mmHg at Day 90) yields $p > 0.05$.]

Table 9

Mean IOP (mmHg) in Nonglaucoma Patients Without Prophylactic IOP Therapy
for Per Protocol (C-96-48)

Treatment		Visit					
		Baseline	24 Hour	Day 7	Day 30	Day 90	Day 180
CELLUGEL	Mean	15.8	19.0	14.7	14.7	13.9	14.4
	Std	2.4	6.5	2.7	2.8	2.5	2.9
	N	63	62	61	62	49	59
	Min	10	7	9	7	8	9
	Max	21	38	19	23	18	24
Control	Mean	16.0	18.7	14.6	15.0	14.5	14.2
	Std	2.9	6.6	2.7	2.7	2.7	2.3
	N	70	70	69	65	56	60
	Min	10	8	9	8	10	9
	Max	28	40	21	22	21	19

b. Frequency of IOPs \geq 30 mmHg

1) Clinical Study C-98-22

With this study design, where prophylactic medications are prohibited, it can be useful to evaluate the frequency of patients presenting in the early postoperative period with IOPs \geq 30 mmHg.

Table 10

Frequency of Patients With IOP \geq 30 mm Hg for Per Protocol (C-98-22)

	Treatment		Visit				
			Baseline	6 Hour	24 Hour	Day 7	Day 21
IOP \geq 30 mmHg	CELLUGEL	%	0.0	15.8%	4.3%	0.0	0.0
		N	0	22	6	0	0
		Total	140	139	140	140	140
	Control	%	0.0	12.2%	8.6%	0.0	0.0
		N	0	17	12	0	0
		Total	140	139	140	138	140

The incidence of IOPs greater than or equal to 30 mmHg were evaluated in the C-98-22 IOP study. At 6 hours following surgery, 15.8% of the CELLUGEL® patients (n = 22) and 12.2% of the control patients

(n = 17) had an IOP \geq 30 mmHg. By 24 hours, a smaller percentage of CELLUGEL® patients had IOP \geq 30 mmHg than the control; 4.3% of the CELLUGEL® (n = 6) and 8.6% of the control (n = 12) patients had an IOP \geq 30 mmHg. These differences are not statistically significant (Fisher's Exact Test yields: p=0.49 at 6 hours and p=0.22 at 24 hours). By the Day 7 examination, there were no IOP elevations \geq 30 mmHg.

2) Clinical Study C-96-48

In the C-96-48 Per Protocol group, the incidences of early IOP elevations \geq 30 mmHg were similar to C-98-22. At 24 hours, 11 CELLUGEL® patients (6.5%) and 11 control patients (6.7%) had an IOP \geq 30 mmHg. In the subgroup of Per Protocol patients without glaucoma who did not receive prophylactic IOP-reducing medication at surgery, the incidences of IOPs \geq 30 mmHg were 11.1% in CELLUGEL® patients (n = 8) and 9.1% in the control patients (n = 7) at 24 hours.

2. Device Failures

There were no device failures or replacements reported during these clinical trials using CELLUGEL®.

XI. Conclusions Drawn from Studies

Results from these clinical studies support the following conclusions:

- CELLUGEL® is clinically equivalent to a marketed control OVD in protecting corneal endothelium cells and maintaining the anterior chamber depth during cataract surgery and IOL insertion.
- CELLUGEL® is clinically equivalent to a marketed control OVD in its effects on postoperative intraocular pressure.
- CELLUGEL® is reasonably safe and effective among patients undergoing cataract surgery and IOL implantation.

XII. Panel Recommendation

In accordance with the provisions of section 515(c)(2) of the act as amended by the Safe Medical Devices Act of 1990, this PMA was not referred to the Ophthalmic Devices Panel, an FDA advisory committee, for review and recommendation because the information in the PMA substantially duplicates information previously reviewed by this panel.

XIII.CDRH Decision

FDA issued an approval order on FEB 24 2000. The applicant's manufacturing facility was inspected on July 2, 1999 and was found to be in compliance with the device Good Manufacturing Practice regulations.

LABELING


CELLUGEL®

OPHTHALMIC VISCOSURGICAL DEVICE

(2% Hydroxypropyl Methylcellulose)

DESCRIPTION:

CELLUGEL® Ophthalmic Viscosurgical Device (OVD) is a sterile, nonpyrogenic, noninflammatory viscoelastic solution of highly purified non-proteinaceous 2% hydroxypropyl methylcellulose (HPMC) with an average molecular weight of 300,000 daltons dissolved in an isotonic, physiological buffer. Each mL of CELLUGEL contains 2% HPMC, 0.525% sodium chloride, 0.075% potassium chloride, 0.048% calcium chloride, 0.03% magnesium chloride, 0.39% sodium acetate, 0.17% sodium citrate, and water for injection. The osmolarity of CELLUGEL is 315 ± 35 mOsM, the pH 7.2 ± 0.4 , and the viscosity $30,000 \pm 10,000$ mPa.s (cps) (at 0.2 sec^{-1} , 25°C).



INDICATIONS

CELLUGEL is indicated for use during surgery in the anterior segment of the eye. It is designed to create and maintain space, to protect the corneal endothelium and other intraocular tissues and to manipulate tissues during surgery. It may also be used to coat intraocular lenses and instruments during cataract extraction and IOL insertion.

CONTRAINDICATIONS

At present, there are no known contraindications to the use of CELLUGEL when used as recommended.

PRECAUTIONS

Precautions are limited to those normally associated with the surgical procedure being performed. As with all ophthalmic viscosurgical devices, a transient rise in IOP in the early postoperative period has been reported in some cases. It is therefore recommended that CELLUGEL be removed from the anterior chamber by thorough irrigation and/or aspiration at the end of surgery to minimize post-operative intraocular pressure increases. Intraocular pressure should be monitored postsurgically and appropriate therapy instituted if significant increases occur. Do not overfill the anterior chamber. In addition to the above, the following precautions should be observed:

- Do not reuse cannulas.
- Use only if material is clear.
- Avoid trapping air bubbles within CELLUGEL before injection.
- **This Product Contains Dry Natural Rubber.**

ADVERSE REACTIONS

In two clinical studies, 348 patients were treated with CELLUGEL and 344 patients were treated with a control Ophthalmic Viscosurgical Device (Healon®*). The incidences of ophthalmic adverse events that were reported in $\geq 1\%$ of the patients are shown in Table 1.

Table 1
Adverse Events

Observation	Treatment	Study 1 ^a		Study 2 ^b	
		N	%	N	%
External Slit-lamp Observations ^c	Cellugel	110	55.3	3	2.0
	Healon	87	44.2	1	0.7
Posterior Capsule Haze	Cellugel	94	47.2	13	8.7
	Healon	87	44.2	13	8.8
Intraocular Slit-lamp Observations ^d	Cellugel	70	35.2	-	-
	Healon	80	40.8	-	-
Macular Degeneration	Cellugel	34	17.1	17	11.4
	Healon	34	17.3	20	13.6
Lid Observations ^e	Cellugel	34	17.1	2	1.3
	Healon	35	17.9	1	0.7
Posterior Segment Observations ^f	Cellugel	26	13.1	2	1.3
	Healon	24	12.2	4	2.7
Nd: YAG posterior Capsulotomy	Cellugel	20	10.1	-	-
	Healon	11	5.6	-	-
Dry Eye ^g	Cellugel	11	5.5	5	3.4
	Healon	8	4.1	2	1.4
Iris Atrophy	Cellugel	10	5.0	-	-
	Healon	6	3.1	-	-
Macular Edema	Cellugel	9	4.5	1	0.7
	Healon	8	4.1	3	2.0

Secondary Glaucoma	Cellugel	Study 1 ^a		Study 2 ^b	
		6	3.0	-	-
Hyphema	Healon	5	2.6	-	-
	Cellugel	5	2.5	-	-
IOL repositioning or exchange	Healon	2	1.0	-	-
	Cellugel	4	2.0	-	-
Intraocular pressure >40 mmHg	Healon	2	1.0	-	-
	Cellugel	3	1.5	6	4.0
Vitreous in the Anterior Chamber	Healon	2	1.0	8	5.4
	Cellugel	2	1.0	1	0.7
Endothelial Damage	Healon	4	2.0	-	-
	Cellugel	2	1.0	3	2.0
Cells (AC cells ≥ grade 3)	Healon	-	-	3	2.0
	Cellugel	-	-	2	1.3
Corneal Edema (≥ grade 3)	Healon	-	-	-	-
	Cellugel	-	-	1	0.7
Nd: YAG anterior synechiolysis	Healon	-	-	2	1.4
	Cellugel	2	1.0	-	-
Retina procedure	Healon	1	0.5	-	-
	Cellugel	2	1.0	-	-
Lid procedure	Healon	4	2.0	-	-
	Cellugel	2	1.0	-	-
Conjunctival cyst or filament removal	Healon	3	1.5	-	-
	Cellugel	2	1.0	-	-
Subjective complaints ^h	Healon	-	-	-	-
	Cellugel	-	-	6	4.0
	Healon	-	-	7	4.8

^a Study 1: CELLUGEL, N=199; HEALON, N=196 (One HEALON patient did not return for follow-up).

^b Study 2: CELLUGEL, N=149; HEALON, N=147.

^c Includes conjunctival injection, conjunctival hemorrhage, superficial punctate keratitis, ecchymosis, arcus senilis, conjunctival chemosis, pinguecula, subconjunctival hemorrhage, hyperemia, conjunctival gape, corneal abrasion.

^d Includes corneal folds, Descemets folds, endothelial folds, striae, guttata, trace endothelial changes, cortical remnants, endothelial pigment, endothelial debris, microcystic corneal edema.

^e Includes blepharitis, dermatochalasis, lid edema, ptosis, collarettes, chalazion.

^f Includes posterior capsular folds/wrinkling, retinal pigment epithelial changes, posterior vitreous detachment.

^g Includes poor tear film.

^h Includes foreign body sensation, ocular pain, diplopia.

Other ophthalmic adverse events considered unrelated to use of the OVD and occurring among patients at a rate of < 1% included: eye discomfort, IOL membrane, puritus, retinal hemorrhage, blurred vision, IOL repositioning with vitrectomy, removal of residual lens cortex, and foreign body removal.

CLINICAL STUDIES

In two controlled, randomized, multicenter clinical studies, 348 patients were treated with CELLUGEL and 344 patients were treated with Healon. A total of 396 patients were enrolled in Study 1 with cell density as the primary endpoint measured at baseline and at the final 6 month visit. Patients who presented with low cell densities were not evaluated. Study 2 was designed to specifically address intraocular pressure elevation during the expected peak period at 6 hours postsurgery with a 21 day follow-up. No prophylactic medications were administered to patients in Study 2 prior to the 6-hour IOP measurement. CELLUGEL and Healon were shown to be clinically equivalent in their effects on postoperative intraocular pressure, based on a statistical test of non-inferiority.

Table 2

Change in Endothelial Cell Density (cells/mm² ± SEM) at 6 months (Study 1)

OVD	Cell Density Change	Percent Change
Cellugel (n=138)	-119 ± 44.3	-3.6%
Healon (n=130)	-135 ± 45.2	-3.8%

Table 3

Frequency of Patients with IOP ≥ 30 mm Hg

OVD	Time Interval		
	Study 1 ^a	Study 2	
	24 Hours	6 Hours	24 Hours
Cellugel	11.1% (n=8/72)	15.8% (n=22/139)	4.3% (n=6/140)
Healon	9.1% (n=7/77)	12.2% (n=17/139)	8.6% (n=12/140)

^a This is a subgroup including only nonglaucoma patients who did not receive a prophylactic IOP reducing medication prior to the 24 hour exam.

Table 4

Mean IOP Change from Baseline (mmHg ± SEM)

OVD	Time Interval		
	Study 1 ^a	Study 2	
	24 Hours	6 Hours	24 Hours
Cellugel	3.0 ± 0.74 (n=72)	7.22 ± 0.69 (n=139)	3.64 ± 0.39 (n=140)
Healon	3.0 ± 0.82 (n=77)	5.71 ± 0.66 (n=139)	3.49 ± 0.51 (n=140)

^a This is a subgroup including only nonglaucoma patients who did not receive a prophylactic IOP reducing medication prior to the 24 hour exam.

HOW SUPPLIED

CELLUGEL is a sterile, nonpyrogenic, single-use, ophthalmic viscosurgical device, supplied in a disposable syringe delivering 1.0 mL, packaged in a sterile peel pouch, and is terminally sterilized by autoclaving. A sterile, disposable, blunt-tipped cannula is provided.

DIRECTIONS FOR USE

FOR INTRAOCULAR USE ONLY. BOTH CELLUGEL AND CANNULA ARE FOR SINGLE-USE ONLY. The syringe assembly is designed only for the injection of the CELLUGEL ophthalmic viscosurgical device it contains. Use of the syringe assembly for aspiration is not advised.

1. Using sterile technique, peel open the pouch containing the sterile syringe or cannula and drop the contents onto a sterile field.
2. Remove cap from syringe tip.
3. It is recommended that the cannula hub be filled with balanced salt solution prior to attaching the cannula to the syringe in order to minimize the introduction of air bubbles into the anterior chamber.
4. Firmly attach the cannula to the tip of the syringe.
5. Remove plastic cartridge from cannula.
6. Purge the remaining air from the system by holding the syringe barrel with one hand and gently depressing the plunger rod with the other hand until CELLUGEL appears at the cannula tip.

CELLUGEL ophthalmic viscosurgical device should be carefully injected into the anterior chamber prior to capsulotomy using standard aseptic techniques. CELLUGEL may be injected into the chamber prior to or following removal of the crystalline lens. Instillation of CELLUGEL prior to lens removal will provide protection to the corneal endothelium from possible damage due to surgical instrumentation during cataract surgery. Additional CELLUGEL may be injected during anterior segment surgery to fully maintain the chamber or replace any volume lost during the surgical procedure.

STORE AT ROOM TEMPERATURE 15° - 30° C (59° - 86° F)

PROTECT FROM FREEZING AND LIGHT.

DO NOT USE THIS PRODUCT AFTER THE EXPIRY DATE WHICH IS PROVIDED ON THE SYRINGE, POUCH, AND CARTON.

STERILE

Caution: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

*Healon® is a registered trademark of Pharmacia & Upjohn Company.