

# Device Errors

## Pulse oximeters

**MONITOR YOUR PATIENT CLOSELY TO AVOID  
SERIOUS INJURIES.**

**BY BEVERLY ALBRECHT GALLAURESI, RN, BS, MPH**  
Nurse-Consultant  
Center for Devices and Radiological Health  
Food and Drug Administration  
Rockville, Md.

A 77-YEAR-OLD WOMAN WAS HOSPITALIZED IN CRITICAL condition for emergency open-heart surgery. The nurse placed a pulse oximeter sensor on the patient's left thumb and taped it tightly. **When** she removed the sensor after the procedure, the thumb was cyanotic and swollen. The patient is now experiencing devascularization, and the circulatory outcome of the digit is unknown.

### What went wrong?

To avoid serious injuries, you need to treat a pulse **oxim-**eter with the same respect you have for more complex medical devices. Assess the patient's peripheral vascular function before you apply the sensor. Place the sensor on a finger that blanches well and is polish-free. Don't place a sensor on the thumb and never tape it tightly. Check the site immediately after you place the sensor and at frequent intervals during a long surgical procedure. If the patient's peripheral circulation is questionable, use a nasal or earlobe sensor instead of a finger probe.

### What precautions can you take?

Follow these recommendations to avoid injuries when using a pulse oximeter:

- Make sure your patient has stable peripheral vascular function before you place a pulse oximeter sensor.
- Clean the sensor site and remove any nail polish.
- If necessary, secure the sensor loosely using paper tape; never tape it tightly.
- Check the site immediately after applying the sensor to ensure that circulation is adequate.
- Include circulation checks with all patient assessments.
- Change the pulse oximeter site every 2 hours or according to your institution's protocol.
- Closely monitor patients with fragile **skin**, such as infants or elderly adults.
- Make sure you're current on how to use any pulse oximeter you're **working** with, including setting and maintaining all alarms. ▀

Although you need to support the adverse event-reporting policy of your health care facility, you may voluntarily report a medical device that doesn't perform as intended by calling MedWatch at 1-800-FDA-1088 (fax: 1-800-FDA-0178). The opinions and statements contained in this report are those of the author and may not reflect the views of the Department of Health and Human Services. *Device Errors* is coordinated by Chris Parmentier, RN.