



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

OCT - 3 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Barbara Pedigo
President
Hunter Research Laboratories, Incorporated
1541 South Vine Street
Denver, Colorado 80210

Re: K002213
Trade Name: Bull's-Eye Disposable Needle Recapping Aid
Regulatory Class: II
Product Code: FMI
Dated: July 18, 2000
Received: July 21, 2000

Dear Ms. Pedigo:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does

not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Statement of Indications for Use

510(k) Number: K002213

Device Name: bull's-eye™

Indications for Use: Bull's-eye™ disposable needle recapping aid helps prevent needle recapping injury by providing a medical grade, cardboard barrier with a target design to promote eye fusion. Bull's-eye™ is octagon shaped to prevent rolling on tray set-ups, it remains on the needle cap through disposal and allows the healthcare worker's hand to remain behind the sharp.

Bull's-eye™ is placed on a sealed needle cap prior to use and remains on the needle cap through disposal. Bull's-eye™ is designed to be used with the dental and medical patient population and on a dental or medical type of needle cap on a dental or medical type syringe.

Rahima Curran
(Division Sign-Off)
Division of Dental, Infection Control,
General Hospital Devices
510(k) Number K002213