

JUL 11 2005

510(k) Summary

Submitter: ClearMedical, Inc.
1776 136th Place NE
Bellevue, WA 98005

Contact: Gene Lim
Ph: (425) 460-2779
Fax: (425) 401-1515

Trade name: ClearMedical Reprocessed Multiple Clip Appliers

Common name: Non-reloadable Multiple Clip Appliers

Classification name: Manual surgical instrument for general use (21 CFR 878.4800)

Product code: GDO – Applier, Surgical, Clip

Predicate device: K771412 – Ethicon Ligaclip MCA clip appliers

Device description: The reprocessed multiple clip applier is an automatic ligating clip applier. The device is preloaded with a minimum of 10 titanium ligating clips that individually advance after each clip application.

Intended use: Reprocessed multiple clip appliers are intended for the ligation of vessels or other tubular structures.

Technological characteristics: Reprocessed multiple clip appliers are used devices that are cleaned, inspected, tested, packaged, and sterilized for an additional single patient use. The technological characteristics of design, material, and functional performance of reprocessed multiple clip appliers are unchanged and remain equivalent to the predicate devices.

Test data: Validation of cleaning, performance, packaging, and sterilization together with biocompatibility testing demonstrate reprocessed clip appliers perform as intended and are safe and effective.

Conclusion: Based on information provided in this submission, ClearMedical Reprocessed Clip Appliers are substantially equivalent to the identified predicate devices and are safe and effective for their intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL 11 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Clear Medical Incorporated
Mr. Mike Kovacs
1776 136th Place Northeast
Bellevue, Washington 98005

Re: K033579

Trade/Device Name: Reprocessed Multiple Clip Appliers Models, MCL20, MCM20,
MCM30, MCS20

Regulation Number: 21 CFR 878.4300

Regulation Name: Implantable staple

Regulatory Class: II

Product Code: NMJ

Dated: May 27, 2005

Received: May 31, 2005

Dear Mr. Kovacs:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2- Mr. Mike Kovacs

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in cursive script that reads "Miriam C. Provost".

Miriam C. Provost, Ph.D.
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K033579

Device Name: Reprocessed Multiple Clip Appliers

Indications for Use:

Reprocessed multiple clip appliers are intended for the ligation of vessels or other tubular structures.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE
IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

Page 1 of 1

510(k) Number K033579