

AUG 27 2004

K040315

MCI Clenching Inhibitor

510(k) Summary of Safety and Effectiveness

Device Generic Name

Occlusal Stent

Device Trade Name

MCI Clenching Inhibitor

Regulation Number and Name

21 CFR 872.5570 Device, Jaw Positioning

Product Code

LQZ

Panel Code

Dental

Applicant Name and Address

MCI-Myohealth Systems

6B Inga Parade, Mt. Martha VIC 3934

AUSTRALIA

Contact Information

Dave Balding

Consultant

26552 Tampico Place

Mission Viejo CA 92691

USA

Phone: (949) 412-0619 or (949) 916-2663

Fax: (949) 916-2663

Device Description

The Myohealth Clenching Inhibitor disrupts the activity of the pericranial musculature by providing contact between the anterior (incisor) teeth and preventing contact of opposing rearward teeth. The modified dental contact provided by this intraoral device modifies the neuromuscular response of the body to reduce contraction of the pericranial musculature and the associated bruxism and temporal mandibular disorders.

The MCI is an anterior device custom fabricated for each patient. The Myohealth Clenching Inhibitor is intended for treating conditions diagnosed by medical professionals.

The MCI Clenching Inhibitor is designed to span 6 teeth, from upper cuspid to cuspid. The larger size of the MCI greatly reduces the risk of inhalation or swallowing this appliance. The MCI Clenching Inhibitor has a thermoplastic lining. This lining ensures excellent long term retention.

K040315

Comparison to Predicate Device

Attribute	Myohealth Clenching Inhibitor	Predicate (NTI Clenching Suppression System, K981546)
Indications for use	Protection of teeth and restorations from the forces of bruxism. Recovery from temporal mandibular disease (TMD) in which bruxism is a contributing factor.	For the prevention of chronic tension and temporal joint syndrome that is caused by chronic clenching of the posterior mandibular and maxillary teeth by the temporalis muscle.
Method of manufacture	Custom fabricated	Custom fabricated
Single patient, intraoral device	yes	yes
Material	Thermoplastic acrylic	Thermoplastic acrylic
Method of action/intended use.	An anterior bite stop that inhibits masticatory musculature.	An anterior bite stop that inhibits masticatory musculature.
Risk of material exposure	None	None



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 27 2004

MCI-Myohealth Systems
C/O Mr. David Balding
Consultant and Authorized Contact
David Balding
26552 Tampico Place
Mission Viejo, California 92691

Re: K040315
Trade/Device Name: MIC-Myohealth Clenching Inhibitor
Regulation Number: 872.5570
Regulation Name: Intraoral Devices for Snoring and Intraoral Devices for Snoring and
Obstructive Sleep Apnea
Regulatory Class: II
Product Code: LQZ
Dated: August 19, 2004
Received: August 23, 2004

Dear Mr. Balding:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K040315

Indications for Use

510(k) Number (if known): K040315

Device Name: MCI-MyoHealth Clenching Inhibitor

Indications For Use:

1. **Protection of teeth and restorations from the forces of bruxism.**
2. **For the prevention of TMJ syndrome through reduction of trigeminally innervated muscular activity.**

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Rano

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K040315

Page 1 of 1