



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 21 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

CR Tech, Inc.
% Mr. Leonard Frier
Official Correspondent
MET Laboratories, Inc.
914 West Patapsco Avenue
BALTIMORE MD 21230-3432

Re: K051494
Trade/Device Name: Medical Image Digitizer, Models:
CDCR 5020s, UCCR 7001 and TCCR 7001
Regulation Number: 21 CFR 892.2030
Regulation Name: Medical image digitizer
Regulation Number: 21 CFR 892.1650
Regulation Name: Image-intensified fluoroscopic
x-ray system
Regulatory Class: II
Product Code: LMA and MQB
Dated: September 8, 2005
Received: September 9, 2005

Dear Mr. Frier:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

21 CFR 876.xxxx	(Gastroenterology/Renal/Urology)	240-276-0115
21 CFR 884.xxxx	(Obstetrics/Gynecology)	240-276-0115
21 CFR 892.xxxx	(Radiology)	240-276-0120
Other		240-276-0100

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510K number: **K051494**

Device Name: Medical Image Digitizer

Indications For Use

The CDCR5020s, UCCR7001 and TCCR7001 are intended to be used in radiological practices as laser digitizers for phosphor storage plate. The images are sent for processing to a PC.

1. Model CDCR 5020s

This device is used on cassette based x-ray systems using the standard Buckey exist on every system.

After exposing the cassette which includes the phosphor plate, the plate is digitized and the digital image is transferred to the PC for processing

2. Models UCCR 7001, TCCR 7001


With these two devices the organ or object for x-ray exposure are positioned in the front of the device (near the phosphor plate which is positioned into the device) and after exposure the digitizing process is taking place, then the Digital image is transferred to the PC for processing.

The only difference in the models are the way of use either Vertical for Upright or horizontal under the Radiographic Table.

ALL ABOVE MODELS ARE NOT TO BE USED FOR MAMMOGRAPHY

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Reproductive, Abdominal &
Radiological Devices

510(k) Number K051494

Prescription Use X

or

Over-The-Counter Use _____