

K051859

SEP 23 2005

ADMINISTRATIVE INFORMATION

Manufacturer Name: Sistema de Implante Nacional, Ltda
Av. Paes de Barros, 485 Mooca
Sao Paulo - SP CEP: 03115-020
Brazil
Telephone +55 11 2169-3000
FAX +55 11 2169-3025

Official Contact: Wladimir Estanquiere

Representative/Consultant: Floyd G. Larson
PaxMed International, LLC
11234 El Camino Real, Suite 200
San Diego, CA 92130
Telephone (858) 792-1235
FAX (858) 792-1236

DEVICE NAME

Classification Name: Implant, Dental, Root Form (DZE) ;
Abutment, Implant, Dental, Endosseous (NHA)
Trade/Proprietary Name: Sistema de Implante Nacional (S. I. N.)
Dental Implant System
Common Name: Endosseous Dental Implant and Abutment

ESTABLISHMENT REGISTRATION NUMBER

The Establishment Registration number for Sistema de Implante Nacional, Ltda is 3004201263. The Owner/Operator number is 9059509.

DEVICE CLASSIFICATION

FDA has classified endosseous dental implants as Class II devices.

CONFORMANCE WITH PERFORMANCE STANDARDS

No performance standards applicable to endosseous dental implants have been established by FDA. However, CP titanium used to manufacture Sistema de Implante Nacional dental implants meet the chemical and mechanical requirements of ASTM F 67 and ISO 5832-2.

PACKAGING/LABELING/PRODUCT INFORMATION

Sistema de Implante Nacional Dental Implants will be packaged in a radiation sterilizable package consisting of a primary container, with implant and auxiliary parts, sealed with a peel-off wrapping and grouped in storage packs.

INTENDED USE

The Sistema de Implante Nacional Dental Implant System is intended to be surgically placed, either immediately or delayed, in the bone of the maxillary and/or mandibular arch to provide support for crowns, bridges or overdentures. Restorations supported by two or more Sistema de Implante Nacional implants may be loaded immediately after implant placement if primary implant stability has been achieved.

DEVICE DESCRIPTION

Sistema de Implante Nacional Dental Implants are threaded, tapered and straight endosseous dental implants made of commercially pure titanium and intended for use with Sistema de Implante Nacional System abutments and instruments. The implants are offered in a multiple of lengths and diameters. They are offered with a machined surface or acid etched.

EQUIVALENCE TO MARKETED PRODUCT

The Sistema de Implante Nacional Dental Implant System is substantially equivalent, for the purposes of FDA's regulation of medical devices, to Class II medical devices that are cleared for marketing in the United States.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 23 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Sin Sistema De Implante Nacional LTDA
C/O Mr. Floyd G. Larson
PaxMed International, LLC
11234 El Camino Real, Suite 200
San Diego, California 92130

Re: K051859

Trade/Device Name: Sistema de Implante Nacional Dental Implant System
Regulation Number: 21 CFR 872.3640
Regulation Name: Endosseous Dental Implant
Regulatory Class: II
Product Code: DZE
Dated: July 6, 2005
Received: July 8, 2005

Dear Mr. Larson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number: K051859

Device Name: Sistema de Implante Nacional Dental Implant System

Indications for Use:

The Sistema de Implante Nacional Dental Implant System is intended to be surgically placed in the bone of the maxillary and/or mandibular arch to provide support for crowns, bridges or overdentures. Implants may be placed immediately after tooth extraction or following bone healing. Restorations supported by two or more Sistema de Implante Nacional implants may be loaded immediately after implant placement if primary implant stability has been achieved.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Ken M. Kelly for HSR

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K051859