

# Session #1: Corrosion


Moderator: Erica Takai, FDA

- Objective 1: Relevance of corrosion testing - Jack Lemons and Brigitta Brott (University of Alabama Birmingham; 50 min)
- Objective 2: Testing commonly performed (HW results; 5 min)
- Objective 3: What was learned from previous testing (HW results & discussion; 30 min)
- Objective 4: Acceptance criteria (HW results & discussion; 60 min)
- Objective 5: Moving forward (HW results & discussion; 90 min)

Objective 1: Describe the relevance of corrosion testing; In vivo corrosion and clinical consequences—Jack Lemons and Brigitta Brott (University of Alabama Birmingham; 50 min)

Objective 2: Identify the types of corrosion testing commonly performed (HW results; 5 min)


# HW: Response demographics regarding types of corrosion testing performed



n = 54


# HW: Response demographics regarding types of corrosion testing performed

**Surface Treatments**



n = 54



**Metal Alloy**



n = 54

# HW: What types of corrosion testing do you perform for cardiovascular implants?


Results from all devices and alloys pooled



\*Note one 'no' response was a non-cardiovascular device

# HW: What types of corrosion testing do you perform for cardiovascular implants?


Results from all devices and alloys pooled



n = 54

n = 10 out of 54 devices reported fretting testing performed separate from fatigue

HW: What types of corrosion testing do you perform for cardiovascular implants?



n = 55

‘Yes’ responses included open circuit potential, explant analysis, and/or immersion tests



Objective 3: Identify what has been learned from previous corrosion testing; (HW results & discussion; 30 min)

HW: Please describe or identify observed corrosion events from in vitro testing (e.g. after fatigue testing) and in vivo experience (if any).


- Total respondents = 15
- In vitro: mostly none observed outside F2129 testing
  - 1 observation of corrosion post-fatigue with overlap
  - 2 observations of nickel release
- In vivo: two 'yes' however, one appears to be due to post-explant device handling
- Few observed in vivo events may be due to low number of returned explants for analysis or under reporting?

# HW: Test parameters for ASTM F2129 testing performed

		N=68
<b>Scan Rate (mV/sec)</b>	<b>1mV/sec</b>	28
	<b>0.167 or 0.2mV/sec</b>	40
		N=64
<b>Solution</b>	<b>PBS</b>	59
	<b>Other (0.9% saline, Hank's)</b>	5
		N=68
<b>Covering/Coating</b>	<b>covered or coated</b>	42
	<b>not covered or coated</b>	26
		N=65
<b>Crevice?</b>	<b>Yes</b>	35
	<b>No</b>	30

- N=11 devices had both pre- and post-fatigue data
- In vivo corrosion observed?  
→ approximately half not analyzed, half not observed


# HW: Breakdown potential for ASTM F2129 testing performed (as manufactured devices):



- Median of values with breakdown


➤ Generally, increase in E<sub>b</sub> observed post-fatigue (median 583mV<sup>12</sup>)

# HW: Stdev reported for individual tests for ASTM F2129 testing performed (as manufactured devices)




➤ Variability within individual tests higher in stainless steel and nitinol

# HW: Rest potential for ASTM F2129 testing performed (as manufactured devices):



Min/max/median

# HW: Eb-Er for ASTM F2129 testing performed (as manufactured devices):



➤ Generally, increase in Eb-Er post-fatigue for NiTi (median 800mV)

# HW: Repassivation potential for ASTM F2129 testing performed (as manufactured vs. post-fatigue devices)

	<b>As-Manufactured</b>				<b>Post-Fatigue</b>
<b>Metal Alloy</b>	<b>All</b>	<b>Stainless Steel</b>	<b>All CoCr</b>	<b>Nitinol</b>	<b>All</b>
<b>Ep (mV)</b>	N=31	N=8	N=5	N=17	N=11
none	N=18	N=6	N=3	N=8	N=5
mean	-48			-164	-73
median	-86			-100	-104
mean stdev	57			23	142
min	-331			-331	-340
max	913			50	350

- <1/2 reported an Ep
- Most reported Ep close to Er (negative #)
- Only 2 reported Ep close to Eb
- No discernable change post-fatigue



HW: If known, please identify and provide values for the in vivo driving forces for corrosion

- Hoar and Mears (1966), Proceedings, Royal Society, Series A, vol 294, p.1439.
  - Potential measured of various alloys (stainless steel, titanium, etc)
  - in goats (femoral plate) and humans (finger pin) over 71 and 90 days respectively
  - Rest potentials in the 100-600 mV range (stainless steel and nickel alloys <300mV)
  - Scratching the surfaces of the metals (in vitro) caused a transient drop in potential 1-30min.

HW: If known, please identify and provide values for the in vivo driving forces for corrosion.

- L. B. Pertile *et al.* (2009) “In vivo human electrochemical properties of a NiTi-based alloy (Nitinol) used for minimally invasive implants.” J Biomed Mater Res 89A, p.1072–1078.
  - OCP of nitinol wires in femoral, iliac and abdominal arteries
  - 6 patients during vascular surgery (measured over ~12min)
  - $-334 \pm 30$  mV/SCE (similar to in vitro measurements)

HW: If known, please identify and provide values for the in vivo driving forces for corrosion


- Shih *et al.* (2000) “Increased corrosion resistance of stent materials by converting current surface film of polycrystalline oxide into amorphous oxide.” Biomed Mater Res, 52, p. 323-332.
  - Nitinol and 316L stainless steel wires with polycrystalline (PO) or amorphous oxide (AO)
  - Implanted in abdominal aorta of dogs, OCP 6hrs
  - -370mV (PO) vs -20mV (AO) for 316L
  - -330mV (PO) vs -30mV (AO) for nitinol
- Nerve cell conduction in  $\sim -10$ -40mV range

## Objective 3. Identify what has been learned from previous corrosion testing (20 min)

- ❖ 3.1. Discuss whether corrosion testing has been predictive of in vivo corrosion (from animal studies or patient explants; 15 min)
- ❖ 3.2. Identify the in vivo driving forces for corrosion (5 min)

Objective 4: Discuss interpretation of  
results and proposed acceptance criteria  
for ASTM F2129  
(HW results & discussion; 60 min)

HW: What parameters do you establish acceptance criteria for / what are those criteria?



- N=3 nitinol indicated no Eb acceptance criterion
- 23 of 26 indicated no Er acceptance criterion

HW: What parameters do you establish acceptance criteria for / what are those criteria?


### Nitinol

<b>Eb-Er</b>	N=13
<b>none</b>	N=7
<b>median</b>	600
<b>min</b>	200
<b>max</b>	718.3

- 16 of 26 indicated no Eb-Er acceptance criterion

- Some responses both with and without set acceptance criteria indicated that they compared results with another marketed device (n = 8)
- Some set acceptance criteria as 'average value' with no more than 1-2 samples below a set minimum value (n = 3)

HW: Is it appropriate to use universal acceptance criteria (across device types and alloys) for ASTM F2129 testing?



n = 17

- Eb-Er might be appropriate but Er varies between alloys
- F2129 results not correlated with in vivo outcomes; need to also consider intended use/location etc.



## Objective 4. Discuss interpretation of results and proposed acceptance criteria for ASTM F2129 (55 min)

- ❖ 4.1. Discuss value of Eb-Er vs Eb as acceptance criteria (15 min)
- ❖ 4.2. Discuss whether specific Eb or Eb-Er values are universally appropriate (20 min)
- ❖ 4.3. Discuss the significance of Ep (10 min)
- ❖ 4.4. Discuss any other corrosion testing concerns (10 min)

## Thought experiment (objective 4):


- ❖ When comparing to a predicate device, when Eb results are the same, but Eb-Er is poorer, is this acceptable?
- ❖ What if Eb is poorer but Eb-Er is better?

## Objective 4. Discuss interpretation of results and proposed acceptance criteria for ASTM F2129 (55 min)

- ❖ 4.1. Discuss value of  $E_b-E_r$  vs  $E_b$  as acceptance criteria (15 min)
- ❖ 4.2. Discuss whether specific  $E_b$  or  $E_b-E_r$  values are universally appropriate (20 min)
- ❖ 4.3. Discuss the significance of  $E_p$  (10 min)
- ❖ 4.4. Discuss any other corrosion testing concerns (10 min)

# HW: Do you have concerns about other limitations of current corrosion test methods?

Responses mainly regarding F2129



n = 17

- Fine as a screening tool but no correlation to *in vivo* outcomes
- Lack of calibration criteria (variability across runs due to test setup)
- Effects of solution (lack of pH fluctuations, composition, de-aeration)



## Objective 4. Discuss interpretation of results and proposed acceptance criteria for ASTM F2129 (55 min)

- ❖ 4.1. Discuss value of  $E_b$ - $E_r$  vs  $E_b$  as acceptance criteria (20 min)
- ❖ 4.2. Discuss whether specific  $E_b$  or  $E_b$ - $E_r$  values are universally appropriate (20 min)
- ❖ 4.3. Discuss the significance of  $E_p$  (10 min)
- ❖ 4.4. Discuss any other corrosion testing concerns (10 min)

Objective 5. Moving Forward:  
Discuss how corrosion should be assessed  
moving forward  
(HW results & discussion; 90 min)

ASTM F2129 as an assessment of clinically  
relevant corrosion (40 min)

# HW: Is ASTM F2129 a sufficient assessment of clinically relevant corrosion?



➤ Only 1 response suggesting to replace F2129 with other tests



## HW: Is there a way to modify?

- Need for correlation to *in vivo* outcomes
- Good for general indication of corrosion resistance but insufficient assessment alone (e.g. add fretting assessment, nickel leach testing)
- Potentiodynamic polarization curves should be assessed more deeply (e.g. Tafel characteristics, repassivation, active pitting, pitting propagation)

# Objective 5. Moving Forward: Discuss how corrosion should be assessed moving forward


- ❖ 5.1. Discuss whether ASTM F2129 is a sufficient assessment of clinically relevant corrosion (10 min)
  - ❖ 5.1.1. If not (30 min):
    - ❖ 5.1.1.1. Discuss whether it is possible to modify F2129 to be more clinically relevant (e.g. acceptance criteria)
    - ❖ 5.1.1.2. Discuss what is a sufficient assessment of clinically relevant corrosion
      - Other standard test methods to evaluate corrosion (e.g. ISO 16429 and 10993-15)?
- ❖ 5.2. Potential modifications to current testing methodologies and assessments (25 min)
  - F2129:
    - Issues surrounding post-fatigue F2129 testing?
    - Others?
- ❖ 5.3. Discuss the utility of galvanic corrosion testing (25 min)
  - Modifications to galvanic corrosion testing?

# HW: What is a sufficient assessment of clinically relevant corrosion?

- F2129 (comparable to predicate devices with acceptable clinical history) + other assessments
- Biocompatibility + 3-6 months animal studies (thrombogenicity and histopathology) + SEM and Auger to check for the passivation layer uniformity and thickness
- Fatigue test (with fretting for a few months) + galvanic corrosion
- Device integrity at expected life

Potential modifications to current testing methodologies and assessments (25 min)

# HW: Best practices for F2129 testing: Do you use a control?






Yes = Predicate device

## How many samples do you use per test?

Metal Alloy	All	Stainless Steel	All CoCr	Nitinol
# Responses	39	8	10	20
Median	6	7	5	10

# HW: Best practices for F2129 testing: Do you induce damage to the covering?



\*Damage induced mainly by simulated use (with and without fatigue), one response intentionally scratched

- All n = 6 stainless steel indicated "Yes"

# Objective 5. Moving Forward: Discuss how corrosion should be assessed moving forward


- ❖ 5.1. Discuss whether ASTM F2129 is a sufficient assessment of clinically relevant corrosion (10 min)
  - ❖ 5.1.1. If not (30 min):
    - ❖ 5.1.1.1. Discuss whether it is possible to modify F2129 to be more clinically relevant (e.g. acceptance criteria)
    - ❖ 5.1.1.2. Discuss what is a sufficient assessment of clinically relevant corrosion
      - Other standard test methods to evaluate corrosion (e.g. ISO 16429 and 10993-15)?
- ❖ 5.2. Potential modifications to current testing methodologies and assessments (25 min)
  - F2129:
    - Issues surrounding post-fatigue F2129 testing?
    - Others?
- ❖ 5.3. Discuss the utility of galvanic corrosion testing (20 min)
  - Modifications to galvanic corrosion testing?

Galvanic corrosion testing (25 min)




HW: Under what conditions do you believe galvanic corrosion testing is needed?

➤ Is galvanic corrosion testing needed if ASTM F2129 results are “good” for a single device containing dissimilar metals?



- F2129 does not evaluate galvanic corrosion (so always should be considered)
- But if materials close in galvanic series (and cathode:anode ratio small) may not be needed
- Confirmatory test

HW: Do you believe there is value in galvanic corrosion testing of overlapped devices of dissimilar metals?



- Yes: if overlap with another device is reasonably likely to occur during clinical use
- Yes: effects of localized transient behavior
- No: Difficulty in determining/ obtaining other device
- No: Most alloys used in devices close together galvanic series (?)

# HW: Best practices for galvanic corrosion testing of a single device containing dissimilar metals

<b>Do you measure the uncoupled potential before and after?</b>	<b>n = 15</b>
Yes	n = 5
No	n = 5
Only before	n = 5
<b>Do you monitor the coupled potential?</b>	<b>n = 14</b>
Yes	n = 11
No	n = 3
<b>What is the endpoint of your test?</b>	<b>n = 14</b>
Time	n = 12
Steady state current	n = 2
<b>Do you perform testing on as-manufactured devices?</b>	<b>n = 14</b>
yes	n = 13
no	n = 1

Time = 12-48 hrs  
(one 4-6 day)

‘No’ response has not  
started testing

## HW: Best practices for galvanic corrosion testing of a single device containing dissimilar metals

What acceptance criteria do you use for steady state current? (n=12)

- $\leq 3 \text{ nA/cm}^2$  (n = 3)
- 2-4 nA/cm<sup>2</sup> (n = 1)
- Mass loss (n = 2)
- None (n = 6)

How many samples do you use per test? (n=13)

- Median = 5

# Objective 5. Moving Forward: Discuss how corrosion should be assessed moving forward

- ❖ 5.1. Discuss whether ASTM F2129 is a sufficient assessment of clinically relevant corrosion (10 min)
  - ❖ 5.1.1. If not (30 min):
    - ❖ 5.1.1.1. Discuss whether it is possible to modify F2129 to be more clinically relevant (e.g. acceptance criteria)
    - ❖ 5.1.1.2. Discuss what is a sufficient assessment of clinically relevant corrosion
      - Other standard test methods to evaluate corrosion (e.g. ISO 16429 and 10993-15)?
- ❖ 5.2. Potential modifications to current testing methodologies and assessments (25 min)
  - F2129:
    - Issues surrounding post-fatigue F2129 testing?
    - Others?
- ❖ 5.3. Discuss the utility of galvanic corrosion testing (20 min)
  - Modifications to galvanic corrosion testing?