

1 this question to you. In your 22-year old child that
2 is being implanted with this intraocular lens because
3 they are minus 13 and too thin for LASIK and contact
4 lens intolerant, what is an acceptable cell loss rate,
5 not coefficient of variation, so that when they are 82
6 they still have a functioning endothelium?

7 What is the number? Is it .9? Is it .6?
8 Is it 1.5? Is it 2?

9 DR. EDELHAUSER: Good question and it's
10 hard to come up with a number because the thing that
11 I think would be ideal to answer that question would
12 be is that if we had a longitudinal study of high
13 myopes and looked and actually measured the cell loss,
14 I think this would be very important. Unfortunately,
15 this is not in the literature. So the -- what you --
16 in order to answer that question, you know, it has
17 jumped around between 1 to 2 percent as to where we
18 stand with it. One can do all kinds of mathematical
19 calculations to see, you know, how many years would
20 the endothelial cells be depleted, half percent, one
21 percent, two percent, and this all assumes a linear
22 decline which I don't think is completely accurate at

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 this particular stage based on the new information.

2 So to answer your question, I don't know
3 what the exact level would be with this without some
4 good longitudinal data to be able to make an absolute
5 judgment on. And I think one of the things that all
6 of the endothelial studies suffer from is that we
7 don't have good epidemiological data on various
8 populations for the corneal endothelium with regard to
9 aging and various types of subsets like high myopia
10 for example.

11 DR. MACSAI: Well, given that lack of
12 security and an absolute number, you know, I'm wary of
13 creating another closed loop anterior chamber IOL
14 disaster that I think most of the cornea surgeons in
15 this room experienced. So what do you think -- I
16 mean, is there a problem with vault and does that
17 correlate with endothelial cell damage? We saw laser
18 flare meter data, not fluoroscopy data, and it looked
19 good but I've posed to the sponsors and I continue to
20 have this concern, vault is good because cataracts are
21 bad but does vault cause posterior chamber -- I mean,
22 posterior iris chafing? Does it release pigment? We

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 don't know. We haven't looked for sample EC lines on
2 gonioscopy.

3 Does that cause some chronic inflammation
4 that over the 60-year life expectancy of this 22-year
5 old, may effect their endothelial cells. Someone
6 needs to, you know, provide some data from the sponsor
7 regarding this concern.

8 DR. VUKICH: Well, there's two things that
9 I think we are putting together. One would be the
10 initial -- that could account for some initial cell
11 loss

12 DR. MACSAI: And then to propose there's
13 an increased rate of loss, there has to be some
14 ongoing irritation to accelerate above baseline. That
15 ongoing accelerated rate, we believe, would be
16 consistent with the morphometric analysis. If we're
17 going to see some sort of insult, whether it be
18 inflammation of which we detected none, whether it
19 would be a mechanical of which again, we would have to
20 postulate some contact with the cornea that we simply
21 have not observed. These chambers have remained well-
22 formed and we have, again, not seen a mechanism by

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 which we can take a posterior chamber lens and equate
2 this into ongoing corneal endothelial trauma. We
3 would really have to propose a new mechanism for a
4 chronic ongoing accelerated loss of endothelial cells
5 that takes into account normal morphology and no other
6 known cause of this accelerated loss. We believe a
7 lot of what we're seeing here is just an extended
8 remodeling period. We have some insult similar to
9 what we'd expect in clear corneal cataract surgery and
10 there is remodeling that stabilizes the population
11 back to its, again, normal redistribution and that, we
12 believe takes as long as three years and we simply
13 can't see an accelerated rate. So I think projecting
14 is difficult but we've certainly accepted the
15 limitations of the data we have and are committed to
16 long-term follow up. It's an important issue and I
17 think it needs monitoring.

18 DR. MACSAI: Have -- you know, when you
19 talk about long-term insult, my concern is not lens
20 corneal touch. My concern is lens iris touch. My
21 concern is that you know, pigment release and has the
22 sponsor in some way separated those with the good high

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 vault, segregate those out, look at their flare meter,
2 look at their angles, look at their transillumination
3 defects, and look at their endothelial cell loss, that
4 particular group, because I think that would help
5 answer the question.

6 DR. SANDERS: Well, we do have data on
7 three lenses were replaced because they were too long,
8 which were the highest vault and if you look at the
9 final endothelial cell densities, 3300, 2400, 2700.
10 They were the highest cell densities at the later time
11 periods so it appears that these cases are not the
12 ones that demonstrate cell loss with time.

13 DR. MACSAI: But they were replaced.

14 DR. SANDERS: Yes, but they were replaced
15 after a fairly long period in the eye.

16 DR. SUGAR: Can I ask a clarification from
17 Marian? Are you implying that pigment release causes
18 endothelial cell loss because I'm not aware of that?

19 DR. MACSAI: I'm not implying pigment
20 release causes endothelial cell loss. I'm implying
21 pigment release implies touch. Touch may insight
22 chronic inflammation and may have some role in this.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 I don't know. I ask the questions of the sponsor
2 because I don't know.

3 DR. WEISS: Dr. McCulley and then Dr. Ho.

4 DR. McCULLEY: I've been around -- Jim
5 McCulley. I've been around since prior to the
6 beginning of clinical specular microscopy. Been
7 through decades of frustrations of trying to listen to
8 people make sense out of and make points based on cell
9 density. And having listened to -- read everything
10 that was provided, having listened to what everyone
11 has said, quite honestly, I'm at a point where it
12 seems to me that what you've presented at least my
13 interpretation of it, would be that we have surgical
14 trauma, endothelial cell loss, and no evidence for
15 anything except continued remodeling. And no evidence
16 for any other mechanism for continued endothelial cell
17 loss or death other than the normal apoptotic death.
18 So I'm not sure where, you know, one could go further
19 with this or what we would ask you to do other than
20 the surveillance that you're doing except to ask is
21 there some other more sensitive way of looking for
22 inflammation which wasn't a part of your PMA. So I'm

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 not even sure how fair that question is.

2 DR. SANDERS: With regard to the
3 inflammation, ocular inflammation was the subject of
4 my PhD thesis so I did quite a bit of work in this
5 area and that's why we included in the PMA five -- I
6 mean, the laser cell flare meter has been basically
7 thought to be too sensitive a measure of inflammation
8 and it's not even allowed for inflammatory studies
9 because it's too easy to show a decrease in
10 inflammation between groups, and five separate studies
11 in the published literature have shown no inflammatory
12 response after the early post-operative period with
13 this implantable contact lens.

14 DR. McCULLEY: Well, then I guess what I
15 would hope is I envision potential hours of discussion
16 about small points relative to endothelial cell loss
17 and cell density in something that is less than an
18 ideal science. So I would hope that the panel would
19 really press Hank, who is the world's expert in my
20 experience on endothelial specular microscopy with any
21 other issues rather than us trying to figure out
22 what's what among ourselves. If we can have -- so I

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 guess my plea is -- to the panel is, please press Hank
2 while he's here to give us the information that will
3 be more expert than we're at to be able to generate
4 amongst ourselves and hopefully have a more efficient
5 discussion of this because to me this is surgical
6 trauma remodeling.

7 DR. WEISS: Dr. Ho.

8 DR. HO: Allen Ho. But is there any
9 evidence that this sub-clinical inflammation has a
10 deleterious effect on the cornea?

11 DR. VUKICH: We have not demonstrated any
12 subclinical information, no.

13 DR. HO: Is there anything in the
14 literature?

15 DR. McCULLEY: I think -- Jim McCulley.
16 They have no evidence for subclinical inflammation and
17 depending on how you define subclinical which
18 presumably would be what we see at the slip lamp,
19 they've gone another step forward, don't have any.
20 What we could do would be again, intuition. My
21 intuition tells me what I've said. It would be
22 intuitively to go back to some of the closed-loop AC

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 IOLs that actually didn't have sub-clinical, they had
2 clinical inflammation that led to loss of endothelium,
3 so I'm not sure that maybe in some of those eyes some
4 of us didn't see the cell and flare that was going
5 along with those AC IOLs but I think if you have
6 chronic inflammation or chronic rise in intraocular
7 pressure, there is data that suggests there is
8 endothelial cell damage over time.

9 But we don't have any of that there and
10 that's one of the things that intuitively leads me to
11 my conclusion, we have no proposed -- we have no
12 support for any mechanism for any continued
13 endothelial cell loss beyond the apoptotic aging.

14 DR. WEISS: Yeah, I would prefer if we
15 could keep the panel discussion in the panel
16 discussion portion and keep the questions while the
17 sponsor is up there because we have limited time. Do
18 you have any other questions specifically for the
19 sponsor?

20 DR. HO: I do. The only patient that had
21 == Allen Ho -- that had severe sustainable loss of
22 vision in this trial was a patient who had a retinal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 detachment and in a group of very high myopes we would
2 expect perhaps without intervention by natural history
3 that you might see retinal detachment.

4 However, one of the predisposing factors
5 to retinal detachment in high myopes is clearly
6 retinal breaks and lattice degeneration. Do you have
7 any data about number one, lattice degeneration
8 retinal breaks pre-operatively and was indirect
9 ophthalmoscopy part of the study procedures pre and
10 post-operatively?

11 DR. VUKICH: A dilated funduscopy
12 examination was required at several intervals
13 throughout the follow-up period and detailed
14 information was collected by the investigators
15 specific to peripheral retinal findings. We don't
16 have that collated specifically but also entry
17 criteria did require a stable retinal exam. Any pre-
18 existing holes or tears or retinal changes that would
19 be considered high risk, of course, were excluded.

20 DR. HO: They were excluded. Stable
21 retinal breaks were included in this or were they
22 treated preoperatively with laser, for example?

1 DR. VUKICH: We do have a patient, I
2 believe, who had -- we had one patient was treated for
3 an acute retinal break.

4 DR. HO: Yeah, I think that's really
5 important to flesh out for a potential consumer of
6 this kind of technology because, you know, that's
7 where you're losing an eye. However, you may not lose
8 that eye based on your intervention. It simply may be
9 natural history. So I think that's -- I would like to
10 see that information. Thank you.

11 DR. WEISS: Do you think the optic size of
12 4.65 had any impact on visual acuity in younger
13 patients who had larger pupils or is this something
14 you didn't look at?

15 DR. VUKICH: Well, visual acuity and
16 quality were two different things. The visual acuity
17 didn't seem to have an impact in terms of the
18 improvement in best spectacle corrected acuity. Those
19 were the patients who actually had the most
20 improvement quality of vision by subject of symptoms.
21 We can stratify that by level and can provide that,
22 yes.

1 DR. WEISS: So you would be able to look
2 at the size of the pupils to see if it had any adverse
3 effect. Dr. --

4 DR. VUKICH: Well, excuse me, let me
5 qualify that by saying, pupil size measurement was not
6 a part of this clinical exam, either preoperatively or
7 during the course of the trial, so we could only
8 stratify it by level of myopia, not by pupil size.

9 DR. WEISS: Okay, so that's an unknown
10 factor.

11 DR. VUKICH: Correct.

12 DR. WEISS: Dr. Grimmer?

13 DR. GRIMMETT: A quick one for Dr.
14 Sanders. Campbell estimated that pigment particles
15 can be as small as one micrometer in size. Does that
16 laser cell meter detect particles that small?

17 DR. SANDERS: Yes, it does. The standards
18 that are used are in the two micron range and those
19 are meant to be certainly large enough. One micron
20 sized particles should be detected by the Kowan
21 machine.

22 DR. GRIMMETT: You had about 20 patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 after the three-month period or so, 25 or something
2 like that up to two years, something like that.

3 DR. SANDERS: Correct, and the cell
4 measurements were essentially below one per area that
5 was seen on average.

6 DR. GRIMMETT: In those 20, okay, thank
7 you.

8 DR. SANDERS: Yes.

9 DR. SLADE: Yeah, Steve Slade, one quick
10 point to address Dr. Macsai's concern about the
11 vaulting, I just want to make it clear that while
12 angle examination, gonioscopy, was not part of the
13 exam, we certainly did slit lamp exams at multiple
14 intervals and at no point did we ever find peripheral
15 touch, so we were looking at grading angles in that
16 fashion and at no point was the vaulting such that it
17 actually caused touch or PAS.

18 DR. WEISS: Thank you. Dr. Mathers?

19 DR. MATHERS: For Dr. Edelhauser, if
20 you're going to postulate that remodeling is the
21 process, it might be helpful to know -- to see these
22 cells and watch them remodel because they're not being

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 created. They've got to be out there. Could you help
2 us by letting us know how many cells you like to see
3 on a cornea to understand the remodeling process.
4 You're looking at 93 here. What would you recommend
5 that we try to look at if we're going to actually
6 understand if remodeling is the issue versus cell loss
7 on a given patient?

8 DR. EDELHAUSER: I think that one, it's
9 important to do more than -- if you want right now the
10 information, more than just central specular
11 microscopy. Obviously, if we have these pooled cells
12 out in the periphery, it would be interesting to see
13 what's happening with those. I mean, and to get a
14 larger cell number, now the -- most of the instruments
15 that we used in specular microscopy you're limited to
16 pretty much about four millimeters in the center,
17 unless you really encourage the patient you can get
18 out to maybe four millimeters off center to look at
19 the periphery. It's not an easy measurement to
20 obtain.

21 DR. MATHERS: But there's a half a million
22 cells in that area, so --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. EDELHAUSER: Yeah, so I mean, one --
2 if one had to say predict the ideal way to really
3 evaluate it, is I think some of the ways that we --
4 that article we published in the AJO is that we did
5 take eight or nine readings across the cornea; one
6 central, four paracentral and for far peripheral and
7 then if you do that, you can -- and then the
8 interesting thing when you do that, Bill, is that you
9 find out that there's a higher percentage of corneal
10 endothelial cells in the superior region. And
11 similarly the German Daus all found the same thing.
12 So you have a 16-percent increase in peripheral
13 endothelial cells in the superior region.

14 DR. MATHERS: Would you recommend that --
15 matching that against controls as a means to obtain
16 this understanding?

17 DR. EDELHAUSER: Well, if we're going to
18 really map out what's happening in the cornea, with
19 any type of surgical situation with remodeling one
20 would have to do that.

21 DR. WEISS: Dr. Bradley has one brief
22 question. I will ask a question and then we're going

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 to have a 10-minute break.

2 DR. BRADLEY: You might need a 10-minute
3 break after my question. I'm bringing -- I'd like to
4 just go back to the issue that Dr. Weiss raised a few
5 minutes ago about pupil size. There seems to be a
6 certain irony here. I mean, one of the motivations
7 for the product is that there are certain people out
8 there whose myopia level is too high although cornea
9 too thin to perform LASIK simply because -- perform
10 LASIK and have the standard 6.5 millimeter diameter
11 optical zone.

12 The replacement product is only having
13 potentially a 4.65 millimeter optical zone. And one
14 of the reasons why we have a large optical zone with
15 LASIK is because we are concerned about pupil size
16 issues. And I'm a bit concerned that we have so
17 little information about pupil sizes of these patients
18 even -- we would anticipate for example, with young
19 adults mesopic light levels that at least half of the
20 light would be passing into the eye outside of the
21 optical zone of the ICL.

22 Under those circumstances, one can only

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 imagine that the image quality would be very poor.
2 Having said all that, the data seems to point that the
3 patients are quite happy with their nighttime driving,
4 your mesopic contrast sensitivity test with a glare
5 source showed perfectly good results and I'm completed
6 confused by that. I wonder if the sponsor could
7 clarify how that could possibly happen with such a
8 small pupil size.

9 DR. VUKICH: Well, we'll start by looking
10 at pupil size. Certainly, when we developed the
11 protocol in 1995, I don't believe that the interest or
12 the understanding of how these pupil sizes could
13 interact with optical quality were fully understood.
14 That said, pupil size we neither an entry criteria nor
15 a parameter that was measured throughout the course of
16 the trial. I think the only way that we can answer
17 that is to go back to the patient's satisfaction
18 surveys and the quality of vision that they report
19 inasmuch as the patients, in fact, didn't seem to be
20 bothered by the theoretical concerns of an optic size
21 smaller than their pupil. Of course, they didn't know
22 this but what they saw they seemed satisfied with.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 I understand and appreciate the concerns
2 even with pupil size. However, there seems to be some
3 variability in the response or the effect of the pupil
4 size that we're understanding now with LASIK where it
5 may not be as much of a correlation as we perhaps,
6 intuitively may expect. So we don't understand the
7 mechanism why a smaller optical size at the level of
8 the lens inside the eye may not have as much influence
9 but yet, we simply have to go back to the results and
10 I believe that they are consistent with patient
11 satisfaction and with the use of this device.

12 To speak to vision quality, there was a
13 subset in a published report looking at vision quality
14 in patients looking at induced aberrations and we
15 found post-LASIK versus ICL, that the ICL patients had
16 one-third as much spherical aberration and half as
17 much coma. And so we certainly believe that it's at
18 least in comparison to LASIK, probably better in that
19 regard at least.

20 DR. WEISS: One last question and this is
21 sort of a bottom line question for Dr. Edelhauser
22 because it seems that the main concern of the panel is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the impact on the endothelium. Would you be surprised
2 if this lens was a contributory factor in causing
3 corneal edema in any of the patients on whom it was
4 implanted?

5 DR. EDELHAUSER: At this stage, no,
6 because the cell density of these patients were well
7 above, you know, 23, 2400.

8 DR. WEISS: I should say eventually. If
9 any of these patients eventually developed corneal
10 edema, in conjunction with having this placed, would
11 that surprise you or do you think that would be
12 totally independent of having this lens placed?

13 DR. EDELHAUSER: Well, when you think
14 about having a lens behind the iris and not rubbing
15 onto the corneal endothelium, it's hard to imagine,
16 you know, the mechanism of what would cause this -- a
17 marked decrease in corneal endothelial cells.

18 DR. WEISS: So you would be -- that as a
19 complication would be surprising to you even 20 years
20 down the line.

21 DR. EDELHAUSER: Yeah.

22 DR. WEISS: Okay.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. SLADE: Just one quick thing, this
2 lens has been implanted outside the U.S., tens of
3 thousands of cases over 10 years and while reporting
4 that experience is not FDA quality, I do believe we
5 would know if this lens ever created a corneal
6 decompensation if the patient had to have a graft and
7 we know of none in that experience.

8 DR. WEISS: Thank you. We're going to
9 take a 10-minute break and I'd ask you to be back here
10 promptly and then we're going to go onto the FDA
11 presentation.

12 (A brief recess was taken.)

13 DR. WEISS: Donna Lochner will be
14 introducing the FDA presentation.

15 MS. LOCHNER: Thank you, Dr. Weiss.
16 Because this is the first phakic IOL to be brought
17 before the panel, I would like to briefly present how
18 FDA's guidance to industry on the design of phakic IOL
19 studies has evolved beginning with the October '98
20 panel meeting. In 1999 ANSI standards and later the
21 ISO meetings began and they currently are held every
22 six months or so. Both the ANSI and ISO standards

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 are expected to be submitted for voting in 2004.

2 Today I'll provide just the highlights of
3 the three panel discussions and then summarize the
4 current ANSI and ISO standards which have incorporated
5 all the major recommendations of the panel with some
6 minor exceptions. FDA issued a draft guidance
7 document in 2000 and expects to issue a final guidance
8 when the ANSI standards are finalized. So this first
9 slide -- I think I went -- this first slide is for the
10 October 23rd, 1998 meeting which, as I said was the
11 first discussion by the panel and at that meeting, the
12 panel recommended that effectiveness criteria
13 generally followed the refractive laser guidance. For
14 example, with respect to the uncorrected VA loss of
15 BSCVA, and also recommended that adverse events in the
16 first year should generally follow the IOL grid for
17 aphakia as a starting point for the study design.

18 The panel recommended a sample size of 500
19 subjects and this was primarily because they felt that
20 as a new indication, new technology, they should take
21 a more conservative approach and the 500 subjects was
22 consistent with what was originally done with IOL

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 aphakia studies. Further, they recommended mesopic
2 contrast and sensitivity testing be done and mesopic
3 pupil size measurements be done. That a questionnaire
4 for visual complaints be administered and that
5 pachymetry, dilated lens and fundus evaluations,
6 topography, keratometry and gonioscopy evaluations be
7 performed.

8 With regards to specular microscopy, the
9 panel recommended a sample size to allow detection of
10 2.5 percent per year and they obtained this figure
11 from the Bourne article that was referred earlier in
12 the discussion this morning. There was a suggestion
13 that all patients be tested but they felt that FDA
14 should try to power the studies to detect the 2.5
15 percent per year. They felt PMA data was needed to
16 three years and if there was a loss or the loss was
17 progressing, a five-year study should be performed.
18 With respect to lens opacities, the panel recommended
19 a clinical grading system and three-year data be
20 collected.

21 The May 12th meeting was held to receive
22 the panel's input prior to publication of FDA's draft

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 guide and at that meeting, the panel generally
2 endorsed our proposals to power the studies to be able
3 to detect a 1.5 percent loss in the specular
4 microscopy study per year and the 1.5 percent figure
5 came after iterating several hypothetical annual
6 losses from a phakic IOL taking an average endothelial
7 cell densities at different age ranges from the
8 literature and determining the age at which the
9 hypothetical annual loss would result in corneal
10 decompensation for the various age groupings. From
11 there we assigned a standard deviation of five percent
12 and sort of arrived at -- which was sort of arrived at
13 as being a reasonable loss so that even young adults
14 would be in their 70s prior to decompensation and that
15 the sample size would still remain reasonable for
16 these studies.

17 The panel endorsed this approach and also
18 asked for data analysis to include a stratification by
19 age. And they further recommended that the analysis
20 look at the mean rate of loss and a frequency analysis
21 to show the percent of patients losing greater than 10
22 percent over the course of the study. With respect to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 lens opacities, the panel again recommended a
2 preoperative and post-operative clinical grading
3 system and at this meeting they also -- there was
4 quite a bit of discussion about control group and felt
5 that that was recommended. The panel also again
6 emphasized gonioscopy and dilated fundus exam.

7 After another two years of meetings with
8 ANSI and ISO we brought a composite of the standards
9 to the panel but with a focused review of endothelial
10 cell density, lens opacity and the contra-sensitivity
11 study. We assigned primary reviewers for each of
12 these three topics and also invited speakers to
13 address endothelial cell design and lens opacity
14 clinical study design issues. The panel recommended
15 that the cell density studies be able to detect the
16 1.5 percent annual loss and this, again, was based
17 upon entry criteria on cell density and acceptable
18 density for the life of the patient. Depending upon
19 the standard deviation, they commented that this will
20 equate to about 200 to 300 eyes. They recommended use
21 of a central reading center or other methods with
22 similar precision and validity. They recommended the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 three-year data was needed for the PMA and also that
2 an intermediate measure between the two and three-year
3 point might be needed to help to establish linearity.

4 Depending upon the three-year data, the
5 panel recommended that additional two years post-
6 marketing study may be needed. And finally, again,
7 the frequency analysis was requested. With respect to
8 lens opacities, again, the panel recommended a
9 clinical rating system and the three-year data also
10 was needed to address the issue of lens opacity and
11 that consideration will be given to longer term at
12 least a five-year post-marketing study. Once a PMA
13 has been reviewed, the panel felt it was useful to
14 look at laser flare and high resolution ultrasound for
15 source of any opacities. And they felt that two or
16 more lines loss with glare or one line without glare
17 would be the level that would be considered clinically
18 significant for any opacity.

19 They further recommended that contrast
20 sensitivity testing be done on all patients to
21 document the severity of any future opacity. With
22 respect to the contrast sensitivity discussion, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 major recommendation that came out of that was that
2 the panel felt a clinically significant decrease in
3 contrast sensitivity should be set at .3 log units and
4 again, the panel emphasized gonioscopy and further
5 stated at this meeting that consideration should be
6 given to collection of data post-market depending upon
7 how the PMA data looked.

8 Again, as I said, all of this culminated
9 in the current draft ANSI and ISO standards with
10 recommendations for a three-year, 300-subject
11 preoperative control study. Safety end points from
12 the FDA's aphakic IOL grid are also used as control
13 data in these standards and now I'll just briefly go
14 through the current recommendations and the most
15 current versions of these standards and that is that
16 the following evaluations be performed; in corrected
17 CVA, distance and near, BSCVA distance and near,
18 manifest and cycloplegic refractions, a subject
19 questionnaire, a slit lamp exam including aqueous cell
20 and flare, gonioscopic exam, corneal edema, pupillary
21 irregularities, iris atrophy and pigment dispersion.

22 These standards recommended a dilated

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 fundus exam, that IOP testing be performed, mesopic
2 pupil size be measured and that pachymetry,
3 preoperative axial length, anterior chamber duct
4 measurement and kerotometry be performed. With
5 respect to specular microscopy, the standards assume
6 a 10 percent surgical loss and recommend that the
7 studies be able to deduct a two-percent loss per year.
8 The standards recommend that all 300 subjects be
9 tested so that at least 200 good images would be
10 obtained.

11 They recommend use of a central reading
12 center and they recommend that 100 to 150 cells be
13 counted. With respect to lens opacities, again the
14 standards recommend a clinical grading system and they
15 recommend that a change in contrast sensitivity
16 performance from preop to each post-op visit at which
17 an opacity is observed be performed to document any
18 significance to the opacity. The standards recommend
19 contrast sensitivity be performed under mesop and
20 mesopic with glare and the sample size recommended is
21 61 subjects.

22 Now, I'd like to thank and acknowledge the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 PMA review team for this application. Dr. Alexander,
2 who is the lead reviewer for the PMA, Dr. Eydelman,
3 the clinical reviewer, Dr. Gray who performed the
4 statistical review, Don Calogero, our jack of all
5 trades who performed engineering, contrast sensitivity
6 and specular microscopy reviews. Susanna Jones
7 reviewed the toxicology. Susan Gouge, microbiology,
8 Charles Sawyer, patient labeling, Pam Reynolds
9 performed the bio-research monitoring review and
10 Vertleen Covington on the quality systems or good
11 manufacturing practices review. And last but not
12 least, I have to give a special thanks to Sally
13 Thornton, who due to the expedited nature of this PMA
14 really had to do above and beyond the amount of normal
15 running around and we couldn't have gotten here today
16 without her excellent support.

17 Now, Dr. Eydelman will present the
18 clinical questions.

19 DR. EYDELMAN: Good morning. This PMA is
20 truly precedent setting and I wanted you to be aware
21 of it for several reasons. First of all, there are
22 currently no phakic intraocular lenses approved in the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 U.S. There are also no currently approved devices
2 requiring intraocular surgery for correction of
3 refractive error. Thirdly, there are no current FDA
4 approved devices for the correction of myopia greater
5 than 15 diopters. In addition, FDA approved IOLs for
6 use only in adults 60 years of age and older until
7 this year.

8 Currently, responses may require lowering
9 age for indication to all adults by reference to our
10 recent publication. This is the first time,
11 therefore, that you're going to be considering a PMA
12 for an IOL intended solely for implantation in young
13 adults. As you heard, this PMA received an expedited
14 review status. That truly meant much shorter
15 turnaround time for both the sponsor and us. To make
16 a point of it, I want you to be aware that the last
17 major clinical amendment wasn't received by FDA till
18 September 3rd.

19 As a result of all this, I haven't been
20 able to receive the sponsor's final panel presentation
21 until today, so please forgive any redundancies that
22 I might have in my presentation. As you have all

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 seen, this was a very large PMA with numerous analysis
2 and I will not try to summarize all of it. I'm merely
3 trying to bring your attention to some information
4 which is relevant to the questions that we ask for
5 your consideration.

6 Regarding lens opacification, there were
7 five eyes in the whole PMA that developed nuclear
8 opacities of two plus at the LOCS scale at two to
9 three years. There were 14 cases of ASC opacities of
10 trace or more. Eleven of them occurred at or before
11 the six months and three cases at one year to 26
12 months post-op. In view of these, do you believe that
13 the three-year follow up is sufficient to establish a
14 lens opacification profile associated with this
15 device? If not, what is your recommendation?

16 Eleven out of the 14 cases of ASC appeared
17 at or before the six-month visit suggesting surgical
18 trauma. Combining surgical experience with V3 and V4
19 models, 50 percent of 87.5 percent if you exclude the
20 problematic site number 15, of early ASC cases
21 occurred within the first eight surgical cases. In
22 the Canadian trial performed by three inexperienced

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 surgeons, 22.5 percent of cases developed ASC
2 opacification.

3 The Dominican Republic study which was
4 performed under supervision of a surgical proctor,
5 demonstrated a rate of 4.8 percent. In light of these
6 findings, do you believe surgeon experience to be an
7 important factor in ASC development, secondary to
8 surgical trauma? If yes, do you believe that future
9 users of this lens should be required to undergo
10 special training?

11 Vault measurements in the study were
12 clinical estimates comparing the slit lamp appearance
13 of the corneal thickness to the interval centrally
14 between the crystalline lens and the ICL. Five
15 hundred micron corneal thickness was assumed for
16 conversion from a percentage of corneal thickness to
17 microns. All measurements in an individual case at
18 every visit were averaged to derive at a vault
19 measurement. So as you can see, it was not a very
20 precise measurement estimate. However, it was done.

21 Patients were graded as having poor vault
22 if investigators consistently graded the space

1 between ICL and crystalline lens as less than 10
2 percent of the central corneal thickness and that
3 equated to about 50 microns. Twenty-four cases of the
4 V4 cohort with this technique were determined to have
5 poor vault, 16.7 percent of them or four out of 24 V4
6 cases with poor vault, subsequently developed ASC
7 opacification in contrast only two percent of cases
8 with good vault had ASC.

9 All three cases of significant ASC
10 opacification of late onset defined as greater than
11 six months in V4 cohort were in the eyes with poor
12 vault. In V3 cohort, 41 percent of cases with poor
13 vault developed ASC versus nine percent of cases with
14 good vault. Gonvers, et al, in his recent publication
15 further supported the relationship of vaulting to
16 cataract information. In the PMA the sponsor
17 recommended replacement of the ICL only in cases of
18 poor vault that exhibited early ASC in areas of ICL
19 touch in subjects with UCVA worse than 20/50. Do you
20 agree with this recommendation? If not, what would
21 you recommend?

22 In the clinical trial, sizing was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 determined by the horizontal white-to-white and ACD,
2 anterior chamber depth measurements. Inherent
3 measurement error associated with caliper measurements
4 was judged by the sponsor to be plus or minus .1
5 millimeter. Anterior chamber depths in the study was
6 measured by ultrasound, Orbscan and IOL master. From
7 the literature review, the sponsor concluded that
8 results may differ by as much as .3 millimeters
9 between different measurement methods.

10 Our own literature review revealed lack of
11 correlation of white-to-white measurements and the
12 sulcus-to-sulcus dimension. We also believed that the
13 literature shows that none of the external
14 measurements, including anterior chamber depth and
15 axial length, have been able to accurately predict
16 internal ocular dimensions. The sponsor believes that
17 this literature evidence currently available is
18 anecdotal and they further point out that all the
19 safety and efficacy data available were obtained with
20 a current sizing algorithm based on white-to-white and
21 ACD measurements.

22 It's interesting to note that looking at

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the distribution of the ICL implanted, 50 percent were
2 performed with 12.5 millimeters, versus 7.6 percent
3 was 11.5 millimeter lens. In the overall PMA cohort,
4 1.5 percent of the lenses were replaced due to
5 inappropriate sizing. Do you believe that the method
6 currently recommended by the sponsor for determination
7 of the overall diameter of the ICL to be inserted is
8 appropriate? If not, what do you recommend?

9 As you heard previously, we asked the
10 sponsor to break up their cohort into four refractive
11 groups. Fifteen to 20 diopter group contained 31 eyes
12 at three years. I want to make sure that you're aware
13 that while preliminary discussion for refractive laser
14 guidance for myopia greater than seven diopters was
15 held at the '97 panel meeting. There was no consensus
16 reached on several issues and therefore, there is no
17 currently available guidance for acceptable safety and
18 efficacy outcomes for high myopes after refractive
19 surgery. For eyes with MRSE greater than 15 diopters,
20 in the ICL cohort, there were 3.8 percent or two eyes
21 that lost greater than two lines, 3.8 percent that
22 lost 2 lines and 17.3 percent that lost 1 line. If

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 you calculate it out, it turns out that at 15 diopters
2 of myopia, magnification factor account for a one-line
3 loss being equivalent to a two-line loss. Therefore,
4 we ask the sponsor to include that in the analysis of
5 their high myopia group.

6 Thus, if you add it up, total loss of one
7 line or greater was 25 percent for the small cohort.
8 Some additional safety outcomes for these eyes were
9 retinal detachment at 3.8 percent, ASC opacification
10 of 5.8 percent and as of 9/15, only -- the sponsor
11 informed us that only one eye of these was clinically
12 -- had clinically significant ASC and that is 1.9
13 percent. Clinically significant nuclear cataract in
14 7.7 percent, ICL removal/cataract extraction performed
15 in 3.8 percent and again, 3.8 percent had an increase
16 of greater than two diopter cylinder.

17 As you heard, currently limitation of ICL
18 power is minus 20 diopters. Inadvertently a lot of
19 eyes with MRSE greater than 15 diopters were targeted
20 for under-correction. Eleven point five percent of
21 them were targeted for greater than three diopters,
22 28.8 for greater than two and 65.4 for greater than

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 one. Looking at predictability, 23.3 percent had
2 accuracy within half diopter, 53.3 was within one
3 diopter. Combining the targeted under-correction was
4 a predictability that you saw resulted in rather large
5 range for resultant MRSE for this group at three
6 years. As you can see, it ranged from minus .85
7 diopters to plus .5 with 10 percent of the eyes ending
8 up greater than four diopter myopia, 26.6 greater than
9 three diopters.

10 Looking at all eyes with preop MRSE
11 greater than 15 diopters 38.7 percent of them were
12 able to achieve 20/40 or better. There were no eyes
13 available that were targeted for emmetropia and had
14 preop of 20/20 or better. While all eyes in this sub-
15 group were -- while there were no eyes that were --
16 there were no patients that were unsatisfied, looking
17 at very extremely satisfied patients, you see that for
18 the group of greater than 15 diopters, the
19 satisfaction percentage drops somewhat to 75 percent.

20 Does the safety and efficacy data for eyes
21 with preoperative myopia of greater than 15 to 20
22 diopters support approval of this refractive range?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 If approval for eyes with preoperative MRSE greater
2 than 15 to 20 is recommended, is the term "correction
3 of" as it relates to this refractive range,
4 appropriate in the indication statement? If not, what
5 alternative term do you recommend?

6 Any time we at FDA consider risk benefit
7 analysis for each of the refractive groups, we have to
8 consider two factors. First, is a safety and efficacy
9 profile for each refractive group with the device in
10 question. In addition, we look at safety and efficacy
11 profile for the currently approved or alternate
12 devices available; in this case, glasses, contacts,
13 LASIK, for each of the refractive groups? With this
14 in mind, does the safety and effectiveness outcomes
15 support approval of STAAR ICL for the eyes with the
16 following preoperative MRSE, minus 3 to minus 7,
17 greater than 7 to 10, and from greater from 10 to 15
18 diopters? Twenty patients in overall PMA cohort
19 required treatment other than IOP-lowering meds in the
20 early post-op period. Seventeen of them requiring
21 additional iridotomies, and three requiring additional
22 irrigation/aspiration procedure. In these 20 eyes,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 IOP ranged as high as 65, with IOP spikes observed
2 between one and 21 days post-op. Most of them,
3 however, were seen in one to two days post-op.

4 Incidents of early post-op spikes was
5 stratified by study site and was shown to range
6 between zero to 20 percent. The differences were not
7 found to be statistically significant. Do you believe
8 that specific recommendations regarding early post-op
9 follow up are needed in the labeling? I want to bring
10 your attention to the fact that the labeling you
11 currently have is not -- did not undergo final FDA's
12 review. We always correct all the inconsistencies.
13 Patient symptoms and quality of vision assessment
14 stratified by refractive groups would automatically be
15 included. Demographics is always included.

16 What we are asking your input on is issues
17 unique to ICL that need to be communicated in
18 physician and patient labeling, possibly as a warning
19 or precaution. In addition, we're asking you to
20 consider issues that will be common to all phakic
21 IOLs, such as possible requirement for exclusion of
22 subjects with low endothelial cell density as a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 function of age. This would be consistent with ANSI
2 PIOL draft standards recommendation for clinical
3 studies. It would, however, imply access to specular
4 microscope for all implanting surgeons.

5 In addition, recommendations for
6 gonioscopy and mesopic pupil size assessment preop and
7 post-op in all patients. This is consistent, once
8 again, with our standards recommendation for all
9 clinical studies. Overall, we want to know what
10 additional labeling recommendations do you have. Now,
11 I would like to introduce Dr. Gerry Gray who will
12 review all of the endothelial cell data analysis and
13 when the Chair is ready, I'll be happy to project all
14 questions as they appear in your handout.

15 DR. GRAY: Good morning. My name is Gerry
16 Gray. I'm the team leader for cardiovascular and
17 ophthalmic statistics and I was the statistical
18 reviewer for this PMA. My comments are going to be
19 restricted to the specular microscopy sub-study. This
20 is an overview of the design. We've heard it several
21 times. We're talking about endothelial cell counts
22 and measurements on endothelial cells based on

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 photographs from a specular microscope and all the
2 images were read at a core center with one reader.

3 The study was originally designed to have
4 a preoperative and then three-month, one-year, two-
5 year follow up. During the course of the study it was
6 modified to add three and four-year visits. And the
7 purpose was to investigate the effects on endothelial
8 cells. There were a total of 306 eyes that were
9 enrolled in this sub-study and it had at least one
10 count. I'm just going to go through a little bit
11 about the accountability of the eyes because it gets
12 a little confusing here.

13 The pattern of missing is not quite
14 standard where everyone has a preop visit and then
15 people start to drop off after that. It's a fair
16 amount different. In fact, there were -- 94 of the 306
17 patients had no preoperative visit. Six people had
18 preop and one subsequent. Thirty-four had preop, two
19 subsequent and 172 had preop and then three of them
20 were after that, and the small numbers after that tell
21 you where the person's last visit was.

22 And all this accountability information is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 based on a data set that was submitted to me by the
2 sponsor for analysis. So actually, I think it was a
3 SAS formatted data set. A couple of more
4 accountability combinations; 154 patients had preop
5 and three-year visits, 57 comes up a couple of times.
6 It's not the same 57 patients but 57 had three and
7 four-year visits, 57 had preop and four-year visits.
8 A total of 67 people had all the visits up to three
9 years and a total of 37 had all visits up to four
10 years. So there's 37 patients out of these 306 that
11 had all the visits.

12 So here's a plot that we've seen before.
13 It's the raw results from the data -- from the study,
14 excuse me. The year or the time has been jittered a
15 little bit to show the distribution there. There are
16 preop measurements and then three months, one year,
17 two years, three years and four years. The dashed
18 blue line here just simply connects the means at those
19 time points. And when we look at this, there's really
20 two questions that are key here. The first one is how
21 -- at what point in time can we say that any effect of
22 the actual surgical procedure, whether it would be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 just lost due to surgical trauma and/or some amount of
2 remodeling, at what point in time would we say that is
3 negligible and we can ignore it and use the data after
4 that to get some estimate of what long-term loss might
5 be? So that's the first key question that we need to
6 think about.

7 And then the second thing is what happens
8 off to the right-hand side of this graph, what happens
9 after five, 10, 20 years down the road? Just to set
10 the stage a little bit, this is -- these numbers here
11 are the mean cell counts for various cohorts of
12 patients that you might think about using in this
13 study. The first cohort is all eyes. That's just all
14 306 eyes that were measured whenever, the baseline
15 preoperative measurement, the mean was 2657 and it
16 steadily declined after that to 2355 at the four-year
17 point.

18 The next cohort, I couldn't fit it in very
19 well, so I call it pre and two plus. Those are all
20 the patients that have a preoperative measurement and
21 then they had at least two measurements after that.
22 So that's 206 of the 306 patients and you can see

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 there, it's fairly similar actually to what we get
2 with all eyes. The next two cohorts are somewhat
3 different. The cohort that only has three and four-
4 year measurements and this is a cohort that in the
5 analysis presented to us by the sponsor for the three
6 to four-year loss they used. You'll note that the
7 main difference here is at the three-year point that
8 measurement of 2355 is somewhat lower and it's
9 actually in fact, lower than the average measurement
10 they got at four years for those 57 patients.

11 And then finally, an even smaller subset
12 was everyone who had all the visits and that shows a
13 similar pattern to the three and four-year one, and I
14 presume these are the numbers that were used to make
15 that plot that came up in the sponsor's presentation.
16 So over the duration of the study, over the three and
17 four years we're talking about here, the estimates of
18 cell loss are fairly stable regardless of how you
19 calculate them. At three years, the range of
20 estimates is 8.5 to 8.9 percent. If you use the 154
21 patients who had preoperative and then three-year --
22 a three-year visit, the estimate is 8.7 percent. And

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the competence interval for that is anywhere between
2 a 10.3 and 7.1 percent loss. In raw numbers that's
3 220, 235 cells per millimeter square and that
4 calculation includes anything that happened to the
5 patients between preop and the three-year point which
6 would be any initial operational loss, any kind of
7 remodeling, any normal loss due to aging over that
8 period.

9 And at four years, we've added on a little
10 bit here and it's anywhere from 8.4 to 9.7 percent
11 loss. Okay, now the big question, of course, is
12 what's the steady-state long-term loss that we can
13 expect to see. What's the long-term rate of change in
14 the endothelial cell density we might think we would
15 see? And it turns out that this estimate depends
16 mostly on those -- on the question of how long we
17 believe the effects of the implantation persist, at
18 what time point can we say whatever remodeling or
19 operative loss we have seen is negligible at this
20 point. And that translates into which of the cohorts
21 we actually used to do that estimation. As you saw in
22 the previous slide, the table of cell densities, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 two cohorts on the bottom that only had -- that had
2 three and four-year measurements had a markedly lower
3 three-year cell count than the others and that's the
4 main difference in terms of what you get out in the
5 estimates.

6 The analysis that was presented to us by
7 the sponsor in this PMA was basically using the
8 percent change between the three and four-year time
9 points, using only those patients who had both three
10 and four-year measurements. That's the 57-patient
11 cohort and it properly did some statistics to account
12 for a correlation within a patient between eyes. And
13 the net result there is an estimated percent change of
14 .07 percent, that is a slight gain. In fact, it was
15 one cell per millimeter squared with a confidence
16 interval between minus 1.4 and positive 1.6 percent.

17 Now other cohorts you'll recall, have
18 relatively higher three-year counts and you can do a
19 lot of different kinds of analyses but the bottom line
20 is that the various analyses using those other cohorts
21 and using all time points or time points other than
22 just the three and four-year, produce a change of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 around minus two percent per year. If you go the
2 fancy statistics route and do random coefficients
3 regression, you get a loss of minus 1.9 percent per
4 year. If you believe that whatever -- that the time
5 cutoff for the operational and/or remodeling change is
6 three months and just use the data after three months,
7 and go through an analysis, it's exactly like the one
8 done by the sponsor, in other words, just use the
9 changes from time point T to T plus one, you get an
10 estimate of minus two percent per year.

11 If you believe that any trauma or
12 remodeling is done after two years and you use the two
13 to three-year difference plus the three to four-year
14 differences, you get an estimate of minus 1.8 percent.
15 And the confidence intervals change a little bit. The
16 one for the -- using the regression is probably the
17 smallest because it has a model to help it make the
18 balance smaller, but those are fairly consistent
19 estimates compared to the difference between them and
20 the one that only uses the three and four-year data.

21 So the key question, of course, is where
22 is that cutoff between operative and/or remodeling

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 loss and whatever you might call steady-state, long-
2 term loss. And all I can do is statistics, right? I
3 don't have the clinical knowledge but I have the data,
4 so using the data that we do have, the question here
5 from the statistical point of view is we see that
6 there's some amount -- in many of the cohorts, there's
7 some amount of leveling off after the three-year point
8 between three and four years and the question is, is
9 that statistically significantly different than
10 whatever the slope we saw between three months and one
11 year, one year and two years, two years and three
12 years. And the answer to that is no. If you'll
13 recall the previous plot, it showed the dotted line
14 that connected the means, it looked pretty much like
15 a straight line and the statistics confirm that.
16 There's no strong evidence that the rate of
17 endothelial cell loss between three and four years is
18 any different than the rate -- the annual rate before
19 that. So in the data we have, there's not strong
20 evidence that it's different. Of course, we only have
21 57 people at four years and that could be do to just
22 random fluctuation or we just don't have a big enough

1 sample at four years to have much statistical power
2 but that's what we have.

3 And just in case you care about the
4 details, this was all based on a piecewise linear
5 model that assumes there's a preoperative loss between
6 zero and three months and then after that, it's steady
7 decline either to three years and then a change to
8 four years or it's straight from three months on. But
9 the implication of all this from the data we have is
10 that is that the steady state loss should be
11 estimated using all the data after three months. And
12 if you'll recall from a previous slide, even if we
13 want to go to two years, it doesn't make that much
14 difference here.

15 And so my best guess is due to long-term
16 loss would be that we have -- first of all, there's a
17 mean preoperative measure of 2651 and with the first
18 three months, the absolute loss is about 1.9 percent,
19 so about a two-percent loss over the first three
20 months, and then after that, the rate of loss per year
21 is about 1.9 percent.

22 If we extend this model a little bit to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 include a three and four-year slope, which again was
2 not warranted by the statistics probably, you do get
3 a pretty similar estimate to what the sponsor had
4 between three and four years of an actual slight gain.
5 So here's the results from the two different fits, the
6 two main different kinds of fit that I'm talking
7 about. First of all, there's a blue line here that's
8 just like the one you saw in the previous plot that's
9 pretty much overlaid by the black line. The black
10 line is the fit that I was describing where we had a
11 linear drop at the three months and then a straight
12 line after that. And the green line out at the end is
13 the analysis that was presented to us by the sponsor
14 which is just using the patients who have three and
15 four-year data. And you look at this plot and you
16 say, well, that's not that much different because you
17 know, the only thing different is maybe the difference
18 between the mean at three years there, but the problem
19 is that we don't really care that much at the four-
20 year point. What we care is what happens after 10,
21 20, 30 years and when you make the plot -- when you
22 show the time span we're talking about those are quite

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 a bit different results.

2 And if you believe the three to four --
3 using the three to four-year data, we're basically a
4 flat line, slight increase over time on the
5 endothelial cell density. If you believe that the
6 loss is going to continue linearly at 1.9 percent per
7 year forever, then after about 20 years you're at the
8 1500 cells per millimeter squared and somewhere around
9 35 years you're down to 800. I don't have any -- I
10 don't show any errors around these lines, in the error
11 bars. If you know much -- if you know about errors
12 for regression the errors go, they move outward the
13 further away you get from the center of the data and
14 if I put them on here, they would -- these estimates
15 are pretty much meaningless I think after 15 to 20
16 years. You don't have very much confidence at all in
17 them.

18 And that brings me to, of course, the
19 caveats that the statisticians always give about
20 extrapolation. It's always a questionable exercise to
21 extrapolate beyond the range of the data we have and
22 especially when we're talking about the range we have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 here. It's highly -- any extrapolation you would make
2 would be highly dependent on the model we use and the
3 assumptions we want to make and both those lines that
4 you saw previously assume that whatever linear trend
5 you saw between three and four years is going to
6 continue forever beyond that.

7 And it's probably in this case a lot more
8 important to think about if it's necessary to obtain
9 good long-term data and if so, how to go about doing
10 that. Okay, now, I'm going to switch gears a little
11 bit and talk about individual patients because maybe
12 more important than the average cell loss through time
13 which is described by the linear fits are questions
14 like what proportions of patients will show a cell
15 loss greater than some critical amount. In other
16 words, what proportion of patients will have cell
17 densities less than 1500 or 800 cells per millimeter
18 squared in 10, 20 or 30 years.

19 And from my point of view, the problem is
20 you can't really answer this with much confidence
21 using the data we have here. But let me just
22 summarize what we do have here. If you'll recall one

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of the previous -- the fancy statistical model I used
2 previously actually gives me an estimate for each eye
3 of what the post-operative ECD change for that eye is
4 and then after that, what's the annual change through
5 time, and so you have a distribution of those
6 estimates for each eye.

7 And using that, you can get -- you can
8 create tables like this that tell you something like
9 in this case four and a half -- excuse me, 10 percent
10 of the patients will have an initial loss of four and
11 a half percent or more and 10 percent of the patients
12 will have an annual loss of 2.9 percent or more.

13 Now, that's based on again, I'm making some assumption
14 that whatever we've seen in the first three or four
15 years is going to continue however far in the future
16 you want to go. Okay, and finally, there were some
17 co-variants that seemed to be significant predictors
18 of endothelial cell loss, notably is the anterior
19 chamber depth which was a statistically significant
20 predictor of cell loss regardless of how you analyze
21 it really. The sponsor presented analysis in the PMA
22 that showed the used binned data, in other words, they

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 broke the ACD into four different groups based on
2 three, three and a half, four millimeter cuts and then
3 presented the cell loss for each of those groups.

4 A bunch of other co-periods didn't appear
5 to be significant predictors of cell loss. Just to
6 help put the ACD effect into context, I created this
7 graph here that takes -- for each eye, you take all
8 the possible annual differences that you got for that
9 eye and calculate from those the percentage loss for
10 that eye and then average those for that one eye. So
11 on the Y axis is for each eye now an annual percentage
12 ECD change that we see in the four years -- after
13 three months. I threw out the first few months
14 because that seemed to be somewhat different. And
15 then platted on the X axis is the ACD measurement for
16 that eye.

17 And the point is that, remember the
18 average ACD is around 3.5 and the average cell --
19 annual cell loss was right here, it's around two
20 percent and down here it says estimated slope is 1.6,
21 so you know that the difference between -- if this
22 right here is about two percent loss, and someone

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 that's a half a unit to the left is going to have a
2 loss that's about 0.8 percent more, 2.8 percent, and
3 someone who is a half a unit to the right is going to
4 have about 0.8 percent, less cell loss. They're loss
5 is going to be about 1.2 percent per year. This is
6 just an attempt to kind of put the -- take the
7 statistical significance of the ACD effects and try to
8 put it in some terms that might be hopefully relevant.

9 So after all that, there's two main
10 questions here for the panel. The first one is that
11 the mean change between three and four years in that
12 57-patient cohort that had both of those was an actual
13 gain of .1 percent in endothelial cell density, so is
14 there sufficient data to support the conclusion that
15 the losses in the first three years are reflective of
16 surgical trauma with some prolonged remodeling period
17 that culminates in a stabilization after three years
18 and if not, what minimum eyes in follow up would you
19 try to make a recommendation that we might need to
20 make that assessment?

21 The second question relating to the
22 anterior chamber depth eyes with the smaller anterior

1 depth of 2.8 to 3 had a greater loss of endothelial
2 cells than the eyes with a greater than 3 millimeter
3 ACD. So the question is, do the outcomes of the ACD
4 analysis provide some assurance of safety in this
5 device for eyes in the lower end and then the upper
6 end of the ACD spectrum? Thank you very much for your
7 attention.

8 DR. WEISS: Thank you. We will now have
9 questions for the FDA from the panel. I'm just going
10 to start off, just to clarify for myself about the
11 endothelial cell loss in terms of determination
12 whether it levels off or increases between three to
13 four years versus whether it continues dropping. From
14 what I understood you to say, if you look at the
15 cohort of 57 which is what the sponsor was looking at
16 between three to four years, you could possibly say
17 that it was going to level off, but if you look at the
18 other cohorts, it does not show that. Am I
19 misinterpreting it or is that basically --

20 DR. GRAY: That's correct. Your estimated
21 amount of endothelial cell loss depends primarily on
22 which cohort you use and the one cohort -- the cohort

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that has either the three and four-year measurement
2 that has three and four-year measurements has a lower
3 three-year count and therefore, you get basically a
4 flat line after that.

5 DR. WEISS: So we have a choice of
6 basically looking at the cohort of 306 and if we look
7 at the cohort of 306, it does not support leveling off
8 between three to four years. If we look at the cohort
9 of 206, it does not support leveling off at three to
10 four years. And if we look at the cohort of 37?

11 DR. GRAY: Well, when you say "support" it
12 might mean a different thing to you than to me. When
13 you get down to the 57 or 37 patients, there is more
14 of a leveling off but on the other hand, there's more
15 air because you have fewer patients. So I didn't
16 actually do the test with the 37 patient cohort, but
17 my guess is that you couldn't say statistically that
18 there was a difference, but I didn't actually do that.

19 DR. WEISS: But certainly for the larger
20 groups, which would have more statistical strength, it
21 shows no leveling off.

22 DR. GRAY: That's correct. I personally

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 used -- concentrated on the group that had a
2 preoperative measurement and then two or more
3 measurements after that because that was the one that
4 I -- in order to do these tests you have to be able to
5 fit a model of some sort.

6 DR. WEISS: So we're talking about if you
7 look at the group of 206, which had the preoperative
8 measurement and measurements at each of these time
9 points, or at some of these time points, at least on
10 two.

11 DR. GRAY: Two or more, yes.

12 DR. WEISS: At two or more of those time
13 points. If you looked at that group, this did not
14 support leveling off between three to four years.

15 DR. GRAY: From a statistical point of
16 view doing the test for leveling, that's correct, it
17 did not support it.

18 DR. WEISS: Okay, thank you. Dr.
19 Grimmett?

20 DR. GRIMMETT: Michael Grimmett. Dr.
21 Gray, I appreciate your comments. On the group of 37,
22 you may not have run the analysis at the end but did

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 you calculate the rates of or the confidence intervals
2 for the endothelial cell loss, what it ranges between
3 for the 37 eyes at year four? Did you show that? I
4 mean, I know for the 57 it was a 90 percent confidence
5 interval was 1.4 something. Did you do the same thing
6 for the 37 eyes? It's probably wider, right?

7 DR. GRAY: No, I didn't do that. It would
8 most likely be wider because of the sample size is
9 three-quarters. So that would increase it by some
10 amount, yes.

11 DR. GRIMMETT: Okay, thank you.

12 DR. WEISS: Dr. Bradley?

13 DR. BRADLEY: Dr. Gray, on one of your
14 last slides there you showed us the relationship
15 between anterior chamber depth and cell loss and you
16 did a linear regression that 1.6 percent per
17 millimeter.

18 DR. GRAY: Yes.

19 DR. BRADLEY: Did you do the analysis to
20 find out how much of the variance was explained by the
21 linear model? That becomes quite an important number
22 for us.

1 DR. GRAY: Well, that was part of the
2 analysis but I don't have that number here on me. The
3 reason I -- I guess their point is that there is a
4 statistically -- when you ask how much of the
5 variation is explained, there is a statistically
6 significant -- that slope is significantly different
7 than zero, okay, so from a statistical point of view
8 there is -- that's a significant slope. And what I
9 was trying to get at was that what's the clinical
10 relevance of that and that's where -- why I made the
11 plot that calculated the 1.6 percent per year. But I
12 don't have that number on me.

13 DR. BRADLEY: Yeah, but it's the clinical
14 significance that's driving my question here in a
15 sense that the linear regression might be highly
16 significant but it may explain a very tiny amount of
17 the variance and thus making policy based upon a
18 parameter which explains only a tiny amount of the
19 variance is really meaningless. So if we had that
20 number or after the meeting somehow that number could
21 be available, that might help policy.

22 DR. WEISS: Dr. Bandeen-Roche, Dr.

1 McCulley and then Dr. Mathers.

2 DR. BANDEEN-ROCHE: Thank you for your
3 presentation. I just have a brief clarification
4 question which is that the numbers that you cited for
5 the four-year mean cell counts differ from the
6 calculations that I cited earlier. And so for
7 instance the three, four-year mean that you cited
8 three years and four years is 2355 and 2356 and
9 reading from Volume 4 of 4, page MD19, those numbers
10 are cited as 2455 and 2456. Now, this in a way sounds
11 like a little point but it goes to the
12 representativeness, the relative representativeness of
13 the various cohorts. So I don't know whether it's
14 clear which one of those is right.

15 DR. GRAY: I'm not, those all differ by
16 exactly 100?

17 DR. BANDEEN-ROCHE: Yes, yes.

18 DR. GRAY: So my first guess is somebody
19 has a typo because that's probably not just a
20 coincidence that they're both exactly 100 off. These
21 calculations that you see here, the mean cell, the
22 sponsor sent me a data set at the end of July, July

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 25th, that has the endothelial cell counts that I
2 later discovered they were rounded -- these are mean
3 so they were rounded off to the nearest cell, the one
4 I got. And that -- the numbers you see here are what
5 I calculated using the data set that I was sent.

6 Now, if the three and four years -- if the
7 two-year number is correct of 2428, then I would say
8 2455 and 2456 are probably not correct, because that
9 would mean that there was an increase between two
10 years and three years as well.

11 DR. BANDEEN-ROCHE: Okay, thank you.

12 DR. WEISS: Dr. McCulley?

13 DR. MCCULLEY: Yeah, I've already
14 expressed a little bit of skepticism about the
15 emphasis being put on cell density but I know those
16 are the numbers you had when you did your analysis,
17 but from a clinical standpoint just over the years,
18 I'm a little skeptical about putting too terribly much
19 weight on something that can vary depending on where
20 you take the count and the variability over time, the
21 reproducibility, so I remain a little skeptical in
22 that regard based on my clinical experience and what

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 I've seen in reviewing papers and hearing
2 presentations over many years.

3 So I guess then my question is, did you do
4 any statistical analysis assessing the size and shape
5 variation over time of the cells?

6 DR. GRAY: No, I did not do that. I used
7 the results that we were submitted to us by the
8 sponsor which seemed to indicate there was really not
9 an issue. So I didn't --

10 DR. McCULLEY: Not, an issue, I'm sorry,
11 meaning what, that there wasn't a change over time?

12 DR. GRAY: There did not seem to be a
13 change through time for either the percent hexagonal
14 or the CV and I didn't dig into that further. I used
15 the same thing that you got in the submission, which
16 is the analysis that the sponsor did.

17 DR. McCULLEY: Yeah, I mean, in the
18 absence of data, I don't really know for sure what's
19 right here and your extrapolation caveats, I think,
20 are good and it would be nice to have the very long-
21 term data, but at least from a cell density
22 standpoint, my impression is that the critical cell

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 density for corneal edema is 800 plus/minus 400
2 roughly tremendous range and tremendous variability.
3 And that these other factors seem to play a very
4 critical role and it would be more comforting for me
5 to know that we had more data to support the size and
6 shape didn't change over time. The numbers just
7 aren't -- or the density isn't the only thing and
8 there's tremendous variability in the measurement
9 methodologies.

10 DR. WEISS: One thing, and I hope that we
11 can pull this perhaps on the lunch break is one
12 difficult item is for the August 2002 panel meeting
13 when we had some of the people who were working with
14 sponsor actually consult and guide the panel as far as
15 what the requirements should be for such a study, I do
16 not recall any such emphasis on hexagonality and
17 coefficient of variation. The number -- the cell
18 density is what was emphasized. Dr. Grimmett can
19 comment in terms if your recollection is any
20 different.

21 DR. GRIMMETT: Yeah, Mike Grimmett. I was
22 the assigned primary reviewer for endothelial analysis

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 at that meeting in August of '02 and in the
2 presentation I made and included in the outline were
3 the references that Dr. Edelhauser was citing
4 regarding the sensitivity of pleomorphism and
5 polymegathism so it was covered. I don't think the
6 sponsor emphasized it or the presenters emphasized it
7 but I did cover it in my presentation, making very
8 similar comments to what Dr. Edelhauser said.

9 DR. WEISS: Dr. Mathers?

10 DR. MATHERS: Thank you for the clarity of
11 your presentation. I thought it was very helpful. In
12 the written work that we were given beforehand, you
13 note that the -- by your model one you had an
14 endothelial cell density loss in absolute numbers of
15 about 49 cells per year and 20 percent of the
16 population actually had a cell loss of 100 cells per
17 year. That's what you're saying. Am I correct in
18 assuming then that that 20 percent of the population
19 in this population would then have an endothelial cell
20 loss rate of about 3.8 percent per year by that
21 calculation? If the 4.9 is average and the average is
22 1.9 by your model 1, it seems to me that would give a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 20 percent of this group that were having a loss of
2 3.8 percent per year. I mean, that's the logical
3 conclusion.

4 DR. GRAY: That is a conclusion that I
5 didn't actually calculate. It's very difficult -- the
6 problem is it's hard enough to estimate the mean
7 function here and now we're trying to estimate the
8 line below which only 10 percent of the people are
9 going to be. And that actually is not -- is even
10 harder statistically.

11 DR. MATHERS: Right, okay.

12 DR. GRAY: The best estimate I can do
13 right now, based on the data we have are what I gave
14 in the presentation, which is that 25 percent of the
15 people will have 2.3 percent or more. Now, if I
16 understand your confidence limits on that, it would be
17 pretty wide.

18 DR. MATHERS: Right.

19 DR. GRAY: I'm not sure exactly what they
20 are. I haven't -- I don't have them on me.

21 DR. WEISS: Seeing no other -- Dr. Macsai?

22 DR. MACSAI: I have three brief comments

1 and I thought all your presentations were great, thank
2 you. The first, they all revolve around endothelium
3 but the first is to Donna. In all your presentations
4 about ANSI and the guidance documents, nowhere did you
5 mention a history of contact lens work and in light of
6 all this discussion about endothelial cell remodeling,
7 I would ask the agency to consider adding that so that
8 that -- I think it's a critical piece of information
9 to help us in the future on any intraocular device.

10 So I didn't see it. Maybe it's there.

11 MS. LOCHNER: It was discussed at some of
12 the earlier panel meetings and the end result was that
13 I think given considerations to the population you're
14 treating and that there is going to be contact lens
15 wear and what's the practical thing to impose on a
16 clinical study, in the end the panel didn't give that
17 emphasis, but I do hear what you're saying and I
18 appreciate the comment.

19 DR. MACSAI: I'm simply asking for history
20 so that you could segregate out --

21 MS. LOCHNER: Oh, yes, yes.

22 DR. MACSAI: -- who wore lenses and who

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 didn't preoperatively. It helps analyze this
2 endothelial cell data.

3 MS. LOCHNER: Yes, and I think many
4 studies will be able to do that.

5 DR. MACSAI: Okay, the second two
6 questions are for Dr. Gray. In this data set you
7 received from the sponsor, do you know if patients who
8 had exchanges at the time of implantation or
9 subsequent to the time of implantation were excluded
10 because that would skew this data, I think
11 significantly?

12 DR. EYDELMAN: I think I actually touched
13 on this in my review. I believe there were two
14 different analysis. In the overall analysis by the
15 sponsor, the data for the eyes that underwent
16 secondary procedure were included, but they were
17 excluded in the analysis where they were determining
18 ACD significance.

19 DR. MACSAI: But were they excluded in
20 measuring endothelial cell density long term?

21 DR. EYDELMAN: They weren't excluded from
22 continuation of collection of data if that's what

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 you're asking. We don't have the analysis for those
2 eyes separated out.

3 DR. MACSAI: Well, do we have an analysis
4 of the eyes that had the lens put in once and only
5 once and never touched again and what happened to the
6 endothelial cells?

7 DR. EYDELMAN: I believe that would be the
8 analysis where the tables for the ACD depth
9 significance were performed.

10 DR. MACSAI: And then I would ask Dr.
11 Gray, looking at those tables, does your slope still
12 hold to the green versus the black slide number 15 or
13 whatever it was, 13, sorry?

14 DR. GRAY: I guess I'm -- first of all,
15 I'm not entirely sure because I don't recall the exact
16 -- I didn't actually do that analysis both ways to
17 compare but the key difference between the estimates
18 that we saw was the fact that the 37 or the 57
19 patients had a lower, a much lower count at the three-
20 year time point than the other group and that's what
21 is driving most of the difference. All the other
22 methods of analysis and different groups of patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 that you include, if you get beyond just the three and
2 four-year data, you have a switch and so all of a
3 sudden, it's about two percent, 1.8, 1.9, 2 and so it
4 really comes down to a question of what time point you
5 think the remodeling is over or whatever happens
6 during the surgery is done with and beyond that, we
7 can consider steady state. And then you get into the
8 whole issue of what does that even mean and how can we
9 extrapolate 20 years down the road which is sort of
10 unanswerable, I think, with the data we have.

11 DR. MACSAI: Maybe I'm not getting
12 something here.

13 DR. EYDELMAN: Let me just try to add, we
14 don't have exactly what you're asking for, Dr. Macsai.
15 We don't have the analysis of just the eyes that had
16 secondary intervention, the endothelial cell separated
17 out. What I do want to point out were that there were
18 few eyes to start out with and chances are some of
19 them did not have the analysis all together. As far
20 as I'm aware, PMA did not contain breakdown for the --
21 on this issue. Certainly your recommendation can look
22 upon it after the panel.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. WEISS: Since we're running 40 minutes
2 behind and we haven't gotten into a discussion, I'm
3 going to have one brief comment by Dr. McCulley, and
4 then we're going to go to five minutes of questions
5 for the sponsor and then break for a 45-minute lunch.

6 DR. McCULLEY: Okay, a critical question
7 seems to be in humans, how long does it take for the
8 endothelium to remodel after an injury and is it
9 degree of injury dependent, is it age dependent? I
10 don't know the answers to those questions but that
11 seems to be absolutely -- the answer to that seems to
12 be absolutely critical in knowing how to interpret the
13 cell density and the cell shape and size change. Do
14 we know that? Do we know how long it takes to -- and
15 maybe when the sponsor comes back, Hank will know.
16 But that's a key question to all of this.

17 DR. WEISS: I want to thank FDA for an
18 excellent analysis and presentation. Sponsor, would
19 you be able to answer or address some of these issues?
20 So you have five minutes to answer all our questions.
21 While the sponsor is setting up, when we break for
22 lunch, I'll just point out, this will be abbreviated.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 It will be 45 minutes, not an hour as listed in
2 deference to the fact that we are running over
3 significantly at this early point in time.

4 MS. THORNTON: Are you ready, Dr. Vukich?

5 DR. VUKICH: Pardon me?

6 MS. THORNTON: Are you ready?

7 DR. VUKICH: I believe so. For some
8 reason, I believe the projector was changed out from
9 underneath us. Okay. We would like to just take a
10 moment to respond to a couple of the questions that
11 were requested of the sponsor. For the number of
12 sites that were contributing to the four-year
13 analysis, this data was collected at eight of the nine
14 sites that were collecting specular micrographs. We
15 were able to calculate the confidence interval for the
16 37-eye consistent cohort of eyes at all of the
17 intervals and that will be the graph that follows.

18 There was clarification that we will need
19 from Dr. Bandeen-Roche on her request for information
20 on an overlies of one of our cohorts, but it may take
21 a little more time than we have available and a little
22 more clarification on exactly what we would like to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 provide. This is the 90-percent confidence interval
2 of the mean for the 37-eye cohort and at four years,
3 which I think is the point of interest. It was 2244
4 to 2509. I see we're taking notes here. Okay, good.
5 This is the entire cohort then for the endothelial
6 cell density. For a point of clarification, this
7 cohort did include all patients and these were also --
8 who were examined and did include patients who had
9 secondary procedures so in some essence it does look
10 at a worse case scenario.

11 A separate analysis of the data,
12 subtracting those patients out has been done. We can
13 tell you that it shows no difference in our
14 estimation. We were hoping it would, but it didn't.

15 There was one final question that we'd
16 like to address and that was from Dr. Bradley. There
17 was a question concerning pupil size and quality of
18 vision. We wanted to point out that our contrast
19 sensitivities were all done under mesopic illumination
20 at 3 candelas per meter squared. Although we did not
21 have pupil size to correlate with that, there would be
22 some assumption that the pupils would be at least

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 smaller than under photopic conditions and that with
2 and without glare there was no demonstrable difference
3 at post-operative contrast sensitivity and in fact, at
4 four of the five measured intervals, there was
5 actually an improvement in contrast sensitivity so we
6 hope that speaks to the quality of vision at least
7 under mesopic conditions.

8 Finally, we'd like to thank the members of
9 the FDA panel for their thoughtful and thorough review
10 of all of this information. Thank you.

11 DR. WEISS: Thank you for making it brief.

12 DR. McCULLEY: Does Hank have an answer to
13 my question?

14 DR. WEISS: We'll find out. Can you make
15 it -- can you give a brief answer and if the answer
16 is, we don't have the information, then that is the
17 answer.

18 DR. EDELHAUSER: I think that is the
19 answer. We don't really have the information. The
20 only really data that we can rely on is probably the
21 keratoplasty data from Bill Bourne which showed a
22 market drop-off, you know, and that's not really the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 data we're after. So we don't have the data.

2 DR. WEISS: No data. Forty-five minutes
3 for lunch and then we'll be starting promptly.

4 (Whereupon, the proceedings in the
5 above-entitled matter went off the record at 12:22
6 p.m. and went back on the record at 1:14 p.m.)

7 DR. WEISS: Can everyone from the panel
8 take their seat, please. We're going to continue the
9 Committee deliberations on this PMA with presentations
10 from Primary Panel Reviewers, beginning with Dr.
11 Marian Macsai-Kaplan. I will remind Panel Members and
12 Sponsor, and FDA, et cetera, that we are now about an
13 hour behind, so I would suggest or request that all
14 comments be short, to the point, and have the purpose
15 of moving this PMA ahead.

16 DR. MACSAI: I'm done.

17 DR. WEISS: With that non-intimidating
18 introduction, I have Dr. Macsai.

19 DR. MACSAI: I would like to first
20 acknowledge a few things. One is, that the Sponsors
21 did an amazing job on a really fast track PMA, and
22 that the FDA did an outstanding job in getting us this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 information as fast as it could be gotten. And I want
2 to really thank Sally for being in such close
3 communication. This was a difficult PMA to review I
4 think for all of the reviewers.

5 The Sponsor has gone through a lot, and so
6 has the FDA, so I'm going to try and limit my
7 comments, but I have a few things I just feel obliged
8 to say.

9 First of all, you saw in the distribution
10 of the patients enrolled in the study, that the vast
11 majority were Caucasian. And from previous devices we
12 looked at, we realized that we do need to look at the
13 affects in non-Caucasian patients. The Sponsor did
14 supply data from the Dominican Republic data set, and
15 I think it would be important for that to be included
16 in anything made available to the public, segregated
17 by refractive error, to help the non-Caucasian
18 population with their expectations.

19 Second of all, exclusion criteria were
20 included, and 65 eyes with pre-existing conditions
21 were included in the study. The results of what
22 happened to those 65 eyes should also be made

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 available by the Sponsor to the Agency, because from
2 those 65 eyes, we may glean information that would
3 help patients who might be treated in an off-label
4 manner.

5 In addition, in the exclusion criteria,
6 limbal pathology was not included, and must be
7 included if a white-to-white measurement is required
8 to size this IOL.

9 Another additional criteria that must be
10 included for exclusion is what the lower limit of
11 endothelial cell counts are per age group. And I
12 would reference Dr. Grimmert's excellent review for
13 that.

14 I'm going to now address efficacy, and
15 then the questions put forward by the Agency.
16 Efficacy of this device is really good, very good.
17 And I'm going to just limit by comments by saying that
18 I was happy to see the efficacy of this data in the 3
19 to 7, 7 to 10, and 10 to 15 diopter groups, and leave
20 the over 15 diopter group for later in my discussion.

21 I would have some questions why a
22 refractive surgeon might use this in a minus 3 diopter

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 group, and until I personally see data that this is
2 superior to refractive surgery already out there, I
3 would personally wonder about that issue.

4 Regarding the specular microscopy data,
5 which was my question 1 in the original questions
6 provided by Dr. Eydelman to us, I feel uncomfortable,
7 plain and simple. I feel uncomfortable because we
8 haven't set a limit of what is the minimal number of
9 endothelial cells that a patient needs to have. We're
10 talking about implanting a device in a 22 year old
11 patient, taking worse case scenario, as the Sponsor
12 said earlier.

13 We've segregated out the patients that had
14 complications, replacements, removal, and if you take
15 a 22 year old and assume that they don't become in
16 need of a cataract until they're 62, assuming they're
17 myopic, they have a higher prevalence of nuclear
18 sclerotic cataracts, you're talking about the device
19 remaining in place for 40 years. And at 40 years,
20 according to Dr. Gray's chart, they're going to drop
21 to a dangerous limit. And so my discomfort comes from
22 the fact that the surgeons who participated in this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 trial are the best of the best. They have the best
2 hands, they have the best experience. I've had the
3 privilege of being taught by some, and observing them,
4 and they are really the best there is, so we're taking
5 a device and releasing it to Joe Q. Average surgeon,
6 and this device will be seen as sort of a drive-
7 through procedure, I'm afraid, where you drive in, you
8 get your IOL, you drive out, you move to Outer
9 Mongolia, and we don't know what happens to you. And
10 we don't know what's going to happen in 10, 20, 30, 40
11 years to the endothelium. So I, of course, having
12 experienced the closed-loop AC IOL induced pseudo
13 phakic bullous keratopathy, am concerned about this
14 device and its effect on the endothelium. And that,
15 to me, is the biggest issue with this PMA. Everything
16 else is really pretty small in comparison to that.

17 Along those lines, we were asked to look
18 at the anterior chamber depths. And I think the
19 Sponsor has shown, Dr. Gray has shown, everyone has
20 shown that in the hands of the best, with an anterior
21 chamber depth less than 3, this device induces a 50
22 percent higher endothelial cell loss. So at this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 time, my recommendation would be that this device not
2 be labeled to be used in an eye under 3 millimeters
3 anterior chamber depth. And that if the Sponsor has
4 further data, that can, of course, be looked at in the
5 future.

6 Question 2 is the nuclear opacities.
7 Nuclear opacities in this population developed at two
8 time courses, early-on, probably surgically-related.
9 Later on, probably nuclear sclerosis developing in
10 these high myopes.

11 I didn't have a big problem with this, but
12 it brings very much to the surface the training of
13 surgeons who are going to use this device. If you
14 look at the Canadian data in those three inexperienced
15 surgeons, there was a 22.5 incidence of anterior
16 subcapsular opacities, while the surgeons that were
17 proctored in the Dominican Republic only had a 4.8
18 percent incidence of anterior subcapsular opacity
19 development. So clearly, that technique used in the
20 Dominican Republic has some effect, so the Sponsors
21 are now left with a huge challenge; how do you take
22 Joe Q. Average surgeon and make him good enough to use

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 this device?

2 And some of my suggestions would be that
3 this device, this Collamer ICL is very similar to the
4 Collamer posterior chamber intraocular lens and Toric
5 intraocular lens that is currently available, and has
6 been for years, for cataract surgery. And that any
7 surgeon who wants to implant this device must first
8 become proficient using that intraocular lens and
9 loading it in the shooter, which is the exact same,
10 and implanting it in the eye. And only after they're
11 proficient with that device, should they then be able
12 to use this device. And they should be proctored one-
13 on-one in the use of this device.

14 But it brings to mind another concern,
15 which is, if you look at the analysis of the
16 investigational sites, one surgeon at one site had a
17 significantly higher number of complications, and a
18 significantly higher number of IOL removals and
19 exchanges. And remember, we're dealing with the best
20 of the best, so I raise this question to the Sponsor,
21 pending release to the general public, how is the
22 Sponsor going to monitor this? If the Sponsor has to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 supply these IOLs to someone who's exchanging them too
2 often, or repositioning them too often, the Sponsor
3 seemingly should have some kind of tracking method for
4 this, and further training required prior to the
5 release of this device. And it's a big, onerous task,
6 but we're talking about putting this in young people
7 with clear lenses, so I think that there's a degree of
8 responsibility the Sponsor will have on this
9 post-approval.

10 Regarding the Agency's question about
11 removal, and if there's areas of touch, and if the
12 uncorrected vision is worse than 20/50, I thought
13 these were fine caveats, but I would also raise the
14 question to both the Agency and the Sponsor, if there
15 is an anterior subcapsular cataract in the visual
16 pathway, should that also be added as a reason for
17 removal?

18 Question 3 regarded the use of the
19 horizontal white-to-white in the anterior chamber
20 depth measurements to determine the sizing of the ICL.
21 I too, like Dr. Grimmett, went back to my operating
22 room and looked at what I had available to measure

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 white-to-white, and it's just a little, I think,
2 Castroviejo caliper, and mine goes by 1 millimeter
3 increments.

4 I, like Dr. Vukich and Dr. Slade, was
5 trained in a time that we did extra caps, we measured
6 white-to-white. I think my residents have done five
7 extra caps in their entire training. I don't think
8 they know how to measure white-to-white. I think the
9 Sponsor is going to either have a huge task of
10 teaching them how to do it, or find a better
11 technique. And for that, I would recommend
12 consideration of the Orbscan, which we now know has
13 been shown in the Wang article from the Development of
14 Ophthalmology Journal to be reproducible. It also
15 supplies your anterior chamber depth.

16 I'm not endorsing that product. I hold no
17 interest in that product, but it's out there, and it
18 would give a reliable reproducible measurement for the
19 beginning surgeon. Regardless though, if the patient
20 has limbal pathology, you cannot ascertain a
21 white-to-white measurement; therefore, that is an
22 exclusion criteria in my mind for this device.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Question 4. There are currently no
2 devices approved in the U.S. for correction of myopia
3 greater than 15 diopters. True. So I feel once again
4 very uncomfortable here.

5 First of all, clearly this device in that
6 population does not correct myopia, it only reduces
7 it. So in light of Dr. Eydelman's question, we have
8 to change "correction of" to "reduction of". But I
9 worry that we, as a panel, are going to arbitrarily
10 set a standard by approving this in this age range.

11 I look to the Agency, and ANSI in their
12 wisdom for guidance, and my feeling is once a guidance
13 document is developed in this population, minus 15 to
14 minus 20, and the Sponsor has this engineering thing
15 worked out, that at that time, once the guidance
16 document is set, if the device meets the guidance
17 document criteria, approval is a no-brainer. But at
18 this time, we have no guidance, and I'm uncomfortable
19 with arbitrary approval, which would set a standard,
20 because I am certain there will be more phakic IOLs to
21 come in the future.

22 Question 5, does safety and effectiveness

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 data support approval of the STAAR ICL for the eyes
2 with the following pre-operative MRSE, minus 3 to
3 minus 7, minus 7 to minus 10, minus 10 to minus 15.
4 And in general, my response to this question is yes.
5 However, there remains this outstanding issue
6 regarding endothelial cell loss, sizing of the IOL,
7 cataract information. I'm not uncomfortable with the
8 cataract formation, sizing of the IOL is fixable. And
9 I guess I feel if Dr. Edelhauser doesn't have the
10 answer for endothelial cell loss, I don't know who
11 will. And so, we're functioning in a big old gray
12 zone. And maybe a warning that might be appropriate
13 is that endothelial cell count must be done on these
14 patients pre-operatively, and should be done on these
15 patients post-operatively for a very long time. And
16 if there is a decrease long-term in endothelial cell
17 count, not from an otherwise obvious condition, such
18 as a high fema, trauma, iritis, that perhaps this
19 device should be explanted to protect these patients
20 from pseudo phakic bullous -- from bullous keratopathy
21 at some time in the future.

22 The Sponsor Question 6, management of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 acute intraocular pressure rises in post-operative
2 period. Well, I'm disappointed that gonioscopy was
3 not performed post-operatively in these patients, and
4 I think that Dr. Lochner's presentation has addressed
5 this issue. A mistake was made in the development of
6 this PMA protocol, and it will have to be rectified in
7 the future. But perhaps if the PIs were made farther
8 in advance - I don't know, one week seems awfully
9 early to me - the PI would have healed, and not of
10 them might have been included. And there wouldn't be
11 a need for reopening in the future.

12 In addition, I think the Sponsor must
13 mandate that the surgeon check the pressure within 4
14 to 6 hours after placement of the device, and again in
15 24 hours, so that if it's the viscoelastic, this can
16 be addressed.

17 Question 7, Sponsors have reported that a
18 number of patients noted glare and/or halos
19 post-operatively. Again, I'm disappointed because
20 though Dr. Schallhorn might feel pupil does not make
21 a difference, and I know this lens is much farther
22 inside the eye, I think we could have learned a great

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 deal from that information. And I would ask the
2 Agency to mandate pupil measurements in the future, so
3 that our patients can have a better idea of what to
4 expect from a device. Without it, we can't answer the
5 question, so we're kind of left -- we need to include
6 the data about glare and halos, what patients
7 experienced. We need to include the data about the
8 quality of vision pre-operatively. It was poor at
9 11.6 percent of patients pre-operatively, but at 36
10 months, it was still poor in 5.8 percent of patients.
11 And that's a little disconcerting, because if you read
12 the recently published paper where they compared an
13 eye with an ICL and an eye with LASIK, those patients
14 were doing great. And I have no doubt that the
15 refractive quality with this device for patients will
16 be better than a minus 10 LASIK. And that the higher
17 order aberrations will be less with this device than
18 a minus 10 LASIK. But I'm still wondering why 5.8
19 percent of the patients rated their vision poor. Who
20 were they, and why was it poor? So that concludes my
21 presentation. Thank you.

22 DR. WEISS: Thank you very much, Dr.

1 Macsai. We're going to have Dr. Joel Sugar, who's the
2 second primary reviewer.

3 DR. SUGAR: Thank you. I'm going to just
4 skip through various parts of my review. Of course,
5 I want to thank and compliment the Sponsor and the FDA
6 reviewers for the excellent job they did in both
7 putting the data together, and then analyzing the
8 information.

9 The accountability was good in the study.
10 The efficacy was good up to the minus 15 diopter
11 range, and beyond that range, certainly reduction of
12 myopia should be the indication, or the labeling
13 should be for reduction of myopia, not for correction
14 of myopia. The stability was good.

15 In terms of safety, the loss of lines of
16 best corrected visual acuity, I thought was very
17 acceptable. I think that you can play games about the
18 fact that the minification has changed and, therefore,
19 you should lose less lines, but what matters to the
20 patient is how well they see. And if they don't lose
21 lines of vision, even though they should have
22 theoretically gained a line of vision, I think they're

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 still benefitted.

2 I was concerned about the patients who
3 required enlargement of their laser iridotomies
4 post-operatively because of elevated interocular
5 pressures. In my review, I had the wrong time periods
6 because I measured from the baseline examination, not
7 from the day of treatment. I'm concerned about the
8 Sponsor developing a better means of assessing the
9 iridotomies, both their spacing and their size, so
10 that these patients won't have the pressure elevations
11 as high as 65, as were noted in the presentations.

12 The retinal detachments, I think were
13 acceptable given the population that was being
14 assessed. The cataracts, I think, were acceptable
15 given the population that was being assessed.
16 Although I have concern about the recommendation for
17 removing the lens when anterior subcapsular cataract
18 is seen at an acuity of 20/50 or greater, I would be
19 more concerned about removing it when there's
20 progression of cataract. If, however, I had the data
21 that I don't have, which is, is going in and taking
22 the IOL out, putting a new one cause more progression

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 of the ASC or not, and I don't think we've been
2 presented with any information to tell us whether that
3 does or does not happen.

4 I'm also concerned in terms of the issue
5 of cataracts, since these are patients who will
6 develop cataracts in the long run, like all of us. Is
7 axial length measurable through the IOL easily or not?
8 Does a new algorithm have to be developed for
9 ultrasonic, or whatever technique is used for
10 measuring axial length?

11 People who have their axial length
12 measured, their anterior chamber depth measured
13 ultrasonically could presumably have that data, their
14 axial length captured concurrently and presented to
15 the patient. And it would, I think, make sense, since
16 this is an implant, that the patients be given a card
17 with the data on the lens implanted. But also, if
18 there's data on their axial length, that that be
19 captured, unless it's easy to measure their axial
20 length with the IOL in place, and it would be nice to
21 know that from the Sponsor. It would also be nice to
22 know whether exchanging the lens in and of itself

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 induces another order of complications.

2 Endothelial cell loss has, I think, been
3 very well discussed, and I guess I do feel that,
4 contrary to what I wrote in my review, that anterior
5 chamber depth less than 3 should probably be
6 contraindicated for this lens.

7 There are a few other minor issues.
8 There's some in the labeling that I mentioned in my
9 review. For example, in the brochure it says that
10 surgeons should never touch the center of the optic
11 with instruments when it's in the eye. I don't know
12 if that's because of concern about leaving imprints on
13 the lens, or it's because pushing the lens, pushing
14 the IOL into the crystalline lens could induce
15 cataract. It would be worth having a statement in the
16 brochure saying why that's an issue.

17 The statements made, again, in the
18 labeling, that this device has "been shown to improve
19 the overall quality of vision", I think that's too
20 broad a brush to paint this with. I think you need
21 specific data saying that some patients have overall
22 vision improvement, some don't, and give data.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The brochure should also, I think, have a
2 picture of the device, and a picture of the
3 positioning so that even if someone's taking a course,
4 they will have some hard copy information, should a
5 question arise about lens positioning; although it
6 seems pretty obvious.

7 In terms of the specific questions, is
8 there sufficient data to suggest that there's
9 remodeling? I think that there is. I'm concerned
10 that we capture more data in four years, and
11 definitely capture data at five years on endothelial
12 cell loss. I don't think that we should wait for that
13 information to approve the product.

14 I already talked about the anterior
15 chamber depth. Do I believe surgeon experience is an
16 issue? Absolutely, and that's been addressed by the
17 Sponsor, saying that there will be mandated training.
18 I also talked about the anterior subcapsular cataract,
19 that we need more information on what secondary
20 interventions do.

21 Do I believe the method for determination
22 of overall diameter is appropriate? I think that it

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 is. I think that white-to-white is not as difficult
2 to measure as has been implied. While Orbscan gives
3 it a .1 millimeter on a standard printout, it gives
4 you the white-to-white up to .1 millimeter, I don't
5 think that that -- and that's been shown to be
6 reproducible, I don't know that it's been shown to be
7 any better than manual white-to-white measurements.
8 And certainly, hasn't been shown with this device to
9 provide any advantages. And it's a substantial
10 expense for the average practitioner, who may not have
11 the Orbscan.

12 We talked about the greater than 15. I
13 think that the device should be approved for
14 correction of myopia up to minus 15 diopters, and for
15 reduction of myopia beyond that level. And I think
16 that ends my review. Thank you.

17 DR. WEISS: Thank you very much, Dr.
18 Sugar. The last reviewer, Dr. Grimmett.

19 DR. GRIMMETT: I'm pleased to have the
20 privilege to make a few comments about the
21 application. I apologize for any redundancy. I
22 didn't have any of the talks before during the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 preparation of my talk. Additionally, part of my
2 purpose and mission is to get the information in the
3 public record, so that interested patients in the
4 future can search relevant issues regarding this
5 device.

6 You've obviously all read my review, the
7 cure for insomnia, and I will try to highlight just a
8 few of those issues, but will not go over the data in
9 excruciating detail. You can be happy about that.

10 Before I dig into the PMA, I'd like to go
11 over a few background issues regarding the application
12 to help us in our overall analysis. First, I want to
13 review a few issues related to phakic IOL lens vault.
14 Proper lens vaulting is clearly critical to the
15 success of this phakic IOL. Excessive vault over the
16 crystalline lens will push the iris forward, and has
17 the following potential complications; angle closure,
18 angle synechiae, iris chafing with potential
19 complications of pigment dispersion and pigmentary
20 glaucoma, iris sphincter erosion, iris translumination
21 defects, and alteration of the normal aqueous dynamics
22 that is pupillary block.

1 On the flip side, a poor vault in the --
2 over the crystalline lens has the potential to induce
3 cataracts due to IOL crystalline lens contact.
4 Moreover, if the IOL is too short, it's theoretically
5 possible for it to be mobile, with possible rotation
6 or anteroposterior movement. Clearly, the vault has
7 to be just right to minimize complications, and the
8 tolerances are expected to be low.

9 With an older version of the ICL, Version
10 3, the Sponsor believes that poor lens vault led to a
11 higher right of anterior subcapsular opacities,
12 quoting results from Sanders, in the Journal of
13 Refractive Surgery in 2002. The current application
14 states that Version 4 has an additional .13 to .21
15 millimeters of anterior vault, as compared to Version
16 3. And while I didn't find data in the PMA to
17 substantiate that, the Sponsor clarified today that's
18 a design issue.

19 In the literature, Gonvers & Associates
20 examined central vaulting with digitized slit lamp
21 photographs in 75 eyes. They had 24 V3s and 51 V4s.
22 At three months, the central vaulting of the 24 V3s

1 was slightly less than the central vault of the 51
2 V4s, but the difference in their study was not
3 statistically different. And they concluded, "The
4 change in design between models V3 and V4 did not
5 achieve its goal, which was an increase in vaulting."
6 I just bring that up because I didn't see any data in
7 this application to substantiate the assertion in
8 vivo. Certainly, it's important to keep in mind that
9 increased vaulting may reduce cataractogenesis at the
10 expense of iris and angle complications.

11 In the application, when looking at
12 vaulting, one gets the impression that the phakic IOL
13 vault is a static situation, but I don't -- this
14 couldn't be further from the truth. Stable phakic IOL
15 vaulting on a day-to-day basis is probably not
16 achievable for numerous reasons. Number one,
17 accommodation has been shown to decrease anterior
18 chamber depth by about a quarter of a millimeter,
19 increase the lens thickness by .28 millimeters, and it
20 decreases the radius or curvature of the anterior
21 surface of the crystalline lens.

22 Number two, lens vaulting may differ,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 depending on whether the patient is supine or prone;
2 that is, gravitational effects. And number three, the
3 light reflex has been shown to cause a reduction in
4 the phakic IOL anterior capsular distance. Therefore,
5 on a day- to-day basis, the actual lens vault is
6 probably a dynamic variable.

7 Here's an ultrasonic image from Kim and
8 colleagues in AJO in 1998. The third image on the top
9 shows accommodation on a 30 centimeter target, and
10 displays a decreased distance between the IOL and the
11 crystalline lens right there, due to changes in lens
12 thickness and radius of curvature.

13 The fourth image shows a relationship of
14 the phakic IOL to the crystalline lens in total
15 darkness right here. And then the relevant change
16 when shining a penlight on this eye. In this
17 particular case, there's IOL lens contact with simply
18 a light reflex. Based upon these data, perfectly
19 static phakic IOL crystalline lens relationships on a
20 day-to-day basis are improbable.

21 Moreover, stable IOL vault over the
22 lifetime of the eye is probably not achievable either

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 for numerous reasons. One, the soft IOL material may
2 flatten with time. Dr. Vukich, I believe, mentioned
3 European or outside the United States data over 10
4 years, that it may not. There is an article in the
5 literature that indicates that it may. I believe it's
6 from Arne.

7 Number two, aging has been shown to
8 increase the lens thickness by 1.24 millimeters from
9 age 40 to age 65. Number three, plate phakic IOLs may
10 rotate or have mobility. And number four, the ciliary
11 sulcus diameter has been shown to decrease by
12 approximately 1 millimeter in diameter from age 40 to
13 80.

14 All of these day-to-day and lifetime
15 issues may lead to intermittent or permanent IOL
16 crystalline lens contact, and may lead to
17 cataractogenesis, pigment dispersion, subclinical
18 inflammation, and/or disruption of the normal aqueous
19 humor dynamics. Given these factors, I can't imagine
20 that ICL positioning will be stable and problem-free
21 for the lifetime of a given patient, especially since
22 this device is intended for young recipients.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 Let's talk about issues related to the
2 sizing of these IOLs. The sizing of the ICL myopic
3 lenses was determined by the horizontal white-to-white
4 and the anterior chamber depth measurements in the
5 following fashion. For anterior chamber depths 2.8 to
6 3-1/2, they added half millimeter to the
7 white-to-white, and for anterior chambers greater than
8 3-1/2, they added 1 millimeter to the white-to-white.
9 For in-between sizes, there was a rounding down and
10 rounding up protocol.

11 Hence, STAAR's sizing methodology is based
12 upon white-to-white measurements. However, valid
13 scientific evidence exists saying that white-to-white
14 measurement do not correlate to the sulcus dimension.
15 So white-to-white measurement does not -- is not a
16 good surrogate marker of the variable of interest, the
17 sulcus diameter.

18 Here is just one piece of information from
19 Reinstein's study, in which he examined white-to-white
20 values with calibrated photographs and
21 sulcus-to-sulcus dimensions with high frequency
22 ultrasound. All this information is in the public

1 domain. It's right off the Internet.

2 The top value shows that of myopic eyes,
3 plotting white-to-white on the X axis, and
4 sulcus-to-sulcus on the Y axis, that there's no
5 correlation for myopic eyes. The same was true for
6 hyperopic eyes.

7 These data imply that a one-size fits all
8 phakic IOL would seemingly have just a good chance of
9 success or failure as basing the ICL upon the
10 horizontal corneal diameter.

11 Let's go ahead and look at a few examples
12 of basing the ICL on white-to-white measurements to
13 display this fact. Here's a case where white-to-white
14 is 11-1/2 OU. Put the ICL based on that, bravo, it
15 looks pretty good - adequate lens vault in both eyes,
16 left and right, so we're pleased with ourselves on
17 this case.

18 The next one we have an asymmetric
19 white-to-white, 11-1/2 on the right, and 12 on the
20 left. However, despite differing white-to-white
21 measurements, the lenses were over-sized in both by
22 about the same amount, rather than an asymmetric

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 amount, and the vault is excessive, causing angle
2 closure, as you can probably see.

3 Here's a case where the same
4 white-to-white existed on both sides, but the vault
5 was excessive on the right, and non-existent on the
6 left, with lens IOL touch. I simply would say that
7 because there's valid scientific evidence indicating
8 there's a lack of correlation between white-to-white
9 and sulcus-to-sulcus, that physician labeling should
10 include relevant material facts indicating the lack of
11 the correlation. In fact, in knowing this data now,
12 it's amazing to me that the vault data within the
13 application is as good as it looks.

14 We'll review a few issues related to
15 glaucoma. And please pardon me, Dr. Coleman. I will
16 defer to your judgment on these issues. I'm just the
17 cornea guy. Projected glaucoma risks for this device
18 include pigment dispersion syndrome, angle narrowing,
19 and angle closure.

20 Regarding pigment dispersion syndrome,
21 it's important to realize that STAAR's study cohort
22 fit squarely within the known risk factors for pigment

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com