

1 [Slide]

2 As a control group--I would like to speak
3 a little bit about our MIDCAB experience in
4 Hanover. We have now enrolled more than 700
5 patients. Out of the first 500 patients we did
6 angiographic follow-up in 6-7 percent of these
7 patients. The first group, which was the big one
8 with 297 patients, had a pre-discharge angiogram.
9 What was pretty interesting was that in about 6
10 percent of these patients we had a highly
11 significant problem at the site of the anastomosis,
12 as you can see here. As a Swedish colleague
13 presented his data with the same problem four years
14 ago at the ASCTS and recommended just to wait
15 because this is part of the healing response, we
16 changed our politics, which you will see on the
17 next slide, and just let the whole situation be as
18 it was; waited for 3-6 months, reevaluated these
19 patients and saw that the degree of stenosis or the
20 number of intimal hyperplasia went down without any
21 intervention from 6 percent to 1 percent.

22 [Slide]

23 There is another example here and, as you
24 can see again, there is a highly significant
25 stenosis here at the pre-discharge angiogram;

1 perfect anastomosis 3-6 months later when we
2 reevaluated the patients. What we learned here is
3 that the healing response is still evolving in the
4 earlier time frame. We changed our angiographic
5 follow-up from pre-discharge to a 6-month follow-up
6 so the remaining 203 patients were evaluated 6
7 months after surgery instead of having a
8 pre-discharge angiographic follow-up.

9 [Slide]

10 So, what is the Hanover experience now
11 with anastomotic devices? Just to give you a quick
12 overview, Hanover does approximately 2,000 open
13 heart procedures per year. It is a large teaching
14 institution. We are affiliated with several
15 research centers so we are exposed to new
16 technologies and clinical trials. The studies I
17 have performed were with Ventrica, St. Jude and
18 Converge. In addition, I have a little experience
19 also with Cardica and Coalescent, however, I just
20 want to present you the data where I have
21 angiographic follow-up.

22 [Slide]

23 St. Jude--we had a prospective, randomized
24 trial with 11 patients where every patient received
25 two proximal anastomoses. One was hand-sewn and

1 the other one was an automatic anastomosis. This
2 is the strongest study design that you can create.
3 The data that we saw was that in 11 patients who
4 were enrolled in the study and came back after 6
5 months and there were 10 postoperative angiograms
6 showing that only 3 grafts were patent. We had 6
7 occlusions and 1 highly significant stenosis at the
8 site of the anastomosis, with a consequent PTCA and
9 stent after graft. Even though the patients were
10 asymptomatic, the study was stopped because the
11 data did not look the way we wanted to have it.

12 [Slide]

13 Even though all patients were
14 asymptomatic, due to several reasons that we can
15 discuss, I think independent of the cause a
16 prospective six-month angiographic evaluation was
17 sufficient in our study to detect performance
18 issues of the device.

19 [Slide]

20 Ventrica was part of a multicenter trial
21 that we did with two other centers. We enrolled
22 100 patients, 48 came from Hanover--

23 DR. TRACY: Can I ask you to start
24 wrapping up?

25 PROF. KLIMA: Yes. The most important

1 information I want to give you here is that in the
2 first 48 patients we had a pre-discharge angiogram
3 and a 6-month angiogram to study efficacy of the
4 device and performance of the anastomosis after 6
5 months.

6 [Slide]

7 The last study we did was with Converge.
8 It was also a multicenter trial. We had 8 weeks of
9 follow-up with good data. However, I think these
10 data just showed us something about the feasibility
11 of the device. It does not give sufficient
12 information about how good the anastomosis will be
13 at 6 months.

14 [Slide]

15 In conclusion, I would say that
16 prospective, multicenter trials should be used to
17 evaluate the performance of anastomotic devices.
18 Retrospective clinical endpoints are not sufficient
19 to give you any information about how good the
20 device really is. A comparison to historical
21 controls should be acceptable. Angiographic
22 follow-up is the gold standard and should be used
23 to evaluate an anastomotic device. I think, as we
24 saw from our experience, a six-month angiographic
25 follow-up is sufficient to address the performance

1 of an anastomotic device. Thank you.

2 DR. TRACY: Thank you. Any questions?

3 Dr. Krucoff?

4 DR. KRUCOFF: Doctor, do you feel like
5 your conclusions are equally applicable to a device
6 used for proximal anastomosis as opposed to a
7 device used for distal anastomosis?

8 PROF. KLIMA: Yes, I think so because we
9 used this kind of protocol in proximal and distal
10 anastomotic devices. I know there are different
11 mechanisms causing graft failure, anastomotic
12 failure, but the six-month follow-up is the period
13 where I think wound healing has finished and the
14 problems which can come up are really device
15 related.

16 DR. FERGUSON: I like very much the idea
17 that you used the patient as his own control. The
18 question is did you have some randomized way in
19 which it went to one vessel distally and another?

20 PROF. KLIMA: Yes. Either the patient
21 needed two distal anastomoses or four distal
22 anastomoses so that you had either a single graft
23 with a connector or not a connector or a sequential
24 graft with a connector or not a connector. We
25 preoperatively randomized which graft would be

1 connected with the automatic connector or would
2 have hand-sewn anastomosis.

3 DR. TRACY: Dr. Yancy?

4 DR. YANCY: Just to be clear, the limited
5 data you showed us suggested 60 percent occlusion
6 for the Symmetry device in the randomized effort
7 you did. Is that correct?

8 PROF. KLIMA: Yes, with six occlusions out
9 of ten patients after six months and patients were
10 asymptomatic.

11 DR. YANCY: So, how has that affected your
12 clinical use of the Symmetry device?

13 PROF. KLIMA: Well, after 11 patients,
14 that was the last implant of that device.

15 DR. TRACY: No other questions?

16 [No response]

17 Thank you. Mr. Foley?

18 MR. FOLEY: Dr. Klima gave our report.

19 DR. TRACY: Thank you. At this point,
20 those are all the people who had specifically
21 requested presentations but I would like to ask the
22 audience if there is anyone else who wishes to
23 address the panel on today's topic or any other
24 topic.

25 [No response]

1 We will then close the open public hearing
2 and we will take a five-minute break.

3 [Brief recess]

4 **Open Public Discussion**

5 DR. TRACY: We will try to start the open
6 committee discussion at this point. We would
7 appreciate it if people who were here this morning
8 speaking would remain for as much of the discussion
9 as possible. Once again, thank you for being here,
10 speakers, and please remain if it is at all
11 possible so that the committee has a chance to ask
12 you any additional questions we didn't get to.

13 At this point we will begin the open
14 committee discussion and I would like the panel
15 members to keep in mind the series of questions
16 that were presented to us by the FDA earlier. So,
17 we are trying to discuss things that will be
18 relative to ultimately coming up with answers to
19 these particular questions.

20 At this point, are there any opening
21 questions or comments from the committee? Dr.
22 Aziz?

23 DR. AZIZ: Obviously, when we talk about
24 patency we want to get some idea of what the
25 anastomosis is like in the OR. A number of people

1 are using variable Doppler flow related probes or
2 different types to give us some idea. Then,
3 obviously, doing angiograms a week later sometimes
4 with those ultrasonic systems you obviously can't
5 get an idea of the patency though you may get some
6 idea of the flow. I mean, ideally, I think
7 obviously one would want to know what the
8 anastomosis looks like at the time you complete it.
9 I believe there are some techniques afoot now that
10 allow you to angiographically evaluate what the
11 anastomosis looks like. Do you any of you know the
12 device I am talking about? Are you guys aware of
13 that? Because I think that is what should be used.

14 DR. TRACY: If you are going to make any
15 comments, just come forward to the podium and
16 please identify yourselves.

17 DR. WOLFE: There is an editorial--

18 DR. TRACY: Sir, please state your name.

19 DR. WOLFE: My name is Randall Wolfe. I
20 am a cardiothoracic surgeon at the University of
21 Cincinnati. There is published, in The Journal of
22 Thoracic and Cardiovascular Surgery, about two
23 years ago, an article on acute assessment for
24 intraoperative assessment of coronary grafts. I
25 was the author of that editorial. It outlines the

1 different techniques.

2 The newest technique, which may be the
3 best, is a 13 MHz probe that can be placed directly
4 over the distal anastomosis and gives data both on
5 the flow and the anatomic construction of the
6 graft.

7 DR. AZIZ: Because, I mean if there is any
8 problem that is intraoperative, then that could be
9 fixed there.

10 DR. WOLFE: That is the idea of the
11 editorial, to ask surgeons to please do
12 intraoperative assessment so we can remove all the
13 technical errors.

14 DR. AZIZ: I mean doing an angiogram a
15 week postop to identify the problem makes it really
16 in a sense difficult to fix at that point. So, in
17 our discussions I think we should really be aware
18 of some of the newer things that are coming down
19 the pike that would allow us to detect and correct
20 it, if at all possible, rather than waiting for a
21 week later.

22 DR. WOLFE: I would be happy to get that
23 to you if you wish.

24 DR. EDMUNDS: Randy, if you are going to
25 use ultrasound won't the clips or the metal

1 confound the signal?

2 DR. WOLFE: Actually, the answer is no.
3 The qualification is there is some degradation of
4 the signal if you have a complete ring but in
5 general you can see quite well. Now, with
6 angiography sometimes there can be a little bit
7 with a ring right at the anastomosis but with the
8 ultrasound you can see the anastomosis.

9 DR. EDMUNDS: Well, I know that when it is
10 hand sutured. We have a paper in the Annals for
11 that. But I thought that any time you have some
12 metal, you know, like looking at a valve, you get
13 reflections.

14 DR. WHITE: We look at stents all the time
15 with IVIS and I think the amount of degradation can
16 be handled. I mean, you can see the lumen pretty
17 well.

18 DR. AZIZ: I think the other thing, you
19 know, some of these, let's say, graft failures
20 clearly are related to the site of the anastomosis
21 but I think rheology also obviously plays a role in
22 terms of your competitive flow in the vessel. If
23 you are doing an LAD and the LAD is not that
24 stenotic or appears stenotic, I mean that could be
25 a contributing factor. Do some of these ultrasound

1 devices help you detect if you have competitive
2 flow?

3 DR. WOLFE: Yes, the 13 MHz probe
4 information was published. I can also provide the
5 panel with that information if you would like
6 because it gives you physiologic data as well as
7 anatomic data.

8 DR. AZIZ: And if you detect that there is
9 competitive flow, what then?

10 DR. WOLFE: In fact, there have been some
11 studies that have shown--actually, there was a
12 single center German study that showed that if you
13 remove the patients that have low flow at the time
14 of their bypass graft you can actually get a better
15 handle on patency because those low flow grafts do
16 have a higher incidence of occlusion in the first
17 week.

18 DR. TRACY: I think there are several
19 things if one were trying to design a trial to
20 figure out how these things work. The immediate
21 issue is can you put the thing on in the first
22 place. It seems like the issue that you are
23 addressing is assessing the acute patency issue.
24 Maybe I could ask Dr. Emery a question. It seems
25 as though there were a variety of technical things

1 that you were describing to overcome an inherent
2 problem with the device, tacking things on, etc.
3 Is that part of the initial assessment of acute
4 patency? Do you look at something, do you know
5 something by looking, and how do you pass on that
6 type of information to another surgeon to deal with
7 this acute issue?

8 DR. EMERY: We address that in several
9 points. First of all, I think the training for the
10 use of these devices was not adequate. It went
11 from deploying three or four of these devices in a
12 pig aorta with artificial pulsations to taking it
13 to off-pump beating heart surgery which involves
14 many other considerations. Just the difference
15 between off-pump beating heart surgery and on-pump
16 beating heart surgery is a whole different mental
17 and physical attitude for the surgeon. I think my
18 other colleagues here would agree. Then you apply
19 a device that changes your operative protocols from
20 distal to proximal, for instance; different ways to
21 measure grafts and the different quality of vein
22 grafts. So, the training I think is important.

23 Then, as you discover technical issues you
24 need to carefully modify what you do to make these
25 work, and tacking of the grafts was one thing. As

1 I mentioned, I over-measure the grafts specifically
2 because I have more fear of a short graft, which is
3 more dangerous to occlude, than a longer graft that
4 you can lay out in various patterns. The length
5 has very little relation to the total blood flow
6 through the graft. It is a very minor portion of
7 Poiseuille's equation for flow. So, you can have
8 your graft a centimeter too long and it can be a
9 very appropriate graft. You just have to be sure
10 it doesn't move, flop or kink itself because the
11 place it will kink is at the nearest fixed point,
12 which is either the distal anastomosis or, more
13 commonly, the proximal anastomosis. So, these
14 little technical issues arose over using these
15 things over time and trying to evaluate what I was
16 doing.

17 DR. TRACY: Dr. Krucoff?

18 DR. KRUCOFF: Just scanning the questions
19 the agency has posed to the committee, it seems to
20 me that we are in a technology where it might be
21 worth using stent study design from the time that
22 plain balloon angioplasty first addressed stenting
23 to current drug-eluting stent platforms and think
24 about the array of technical--there are obviously
25 some structural elements; there are technical

1 implantation techniques. There is a sort of acute
2 procedural outcome; there is near-term and there is
3 long-term outcome however we have gone about
4 defining that over the years. It seems to me to be
5 quite relevant to the purpose of these devices, to
6 make the surgery faster, easier, off-pump.

7 You know, I think you could design a
8 pretty clear series of targets of what are you
9 after in using these devices and then what would
10 demonstrate safety and efficacy along the way. The
11 one thing that, for instance, to me seems very
12 clear is that I do think the proximal anastomosis
13 issues are very different than the distal
14 anastomosis issues and to separate that would seem
15 to me to be a very obvious place to start.

16 Then, getting descriptors together from
17 the literature--you know, what is it about age,
18 diabetes, the number of grafts, the diameter of the
19 graft targets the surgical community has found to
20 be predictive, published as predictive, would
21 create a propensity population. Then, depending on
22 what the objective of a give device is I think you
23 could begin to draw down on when do you want to
24 assess it; do you want to assess it; why you are
25 doing the procedure in a way that may help you do a

1 better procedure. Do you want to assess it in an
2 initial cohort where you do an angiogram before
3 they are out of the hospital? Do you want a larger
4 trial and get out to six months, ten months, or
5 have a discussion about the location for the
6 specific device what is the timing that makes the
7 sort of best primary endpoint.

8 I think maybe we could find that the range
9 of approaches here sort out into something very
10 similar to what we have done over the past twelve
11 years with stents.

12 DR. TRACY: Dr. White?

13 DR. WHITE: Mitch, I completely agree with
14 you. Particularly for the proximal anastomosis
15 issues we are talking about stent-like designs. As
16 an angiographer, as we have just been talking
17 about, I would also make a plea that we consider
18 non-invasive imaging for these devices and we
19 explore the limits of non-invasive imaging because
20 there is a risk of angiography. There is a finite
21 risk of angiography and if there was a way that we
22 could satisfy ourselves about the patency and about
23 the anatomy, morphology, then I would prefer any
24 non-invasive tool to do this than to actually ask
25 for routine angiography. There may not be a way to

1 do that. I don't know enough about the
2 non-invasive imaging to be conclusive but that is
3 the reason why I asked about the MR. But I think
4 we shouldn't just sit back and say everybody has to
5 have an angiogram at six months if we are trying to
6 do no harm. I think if there was a non-invasive
7 way to get that information I would like to promote
8 that.

9 DR. TRACY: Of course, I agree with some
10 non-invasive way to look at these things in the
11 long run. I am concerned that there are different
12 time periods. There is the acute issue where you
13 are plugging the thing in and then how you assess
14 that intraoperatively, then what are the time
15 frames and what are the correct follow-ups. I
16 think your only tool during angiography with the
17 stents is with IVIS. Is that correct?

18 DR. WHITE: There are tools. There is the
19 intracardiac echo machine that has low frequencies,
20 9 Hz and so forth, that is analogous to TEE in many
21 respects, that can be used to look nicely at the
22 ascending aorta but that is a venous invasive exam
23 which might be preferable to an arterial invasive
24 exam.

25 DR. YANCY: By the same token, you are

1 particularly using like a 4 tesla magnet MR and can
2 see the proximal stenosis of grafts pretty well so
3 you can get some structural data.

4 DR. TRACY: Dr. Hirshfeld?

5 DR. HIRSHFELD: I think I would sort of
6 like to echo what Dr. Krucoff and Dr. White said.
7 I think there are two core questions. The first
8 question is do these devices confer an important
9 advantage over traditional hand-sewn anastomoses?
10 The second is do they have any downsides either
11 early or late in terms of how patients are in the
12 long term?

13 So, it seems to me that the way to assess
14 these devices is a combination of documenting early
15 patency, which I think can be done non-invasively
16 with CT angio or MR, and then at some point
17 documenting the morphology of these anastomoses
18 with a technique that has high enough resolution
19 such as selective angiography.

20 DR. TRACY: Dr. Mack, did you have a
21 comment?

22 DR. MACK: Michael Mack. I have three
23 comments regarding the issues that have been raised
24 here. First of all regarding intraoperative
25 patency, I think that is going to be extremely

1 problematic in terms of designing any trial. First
2 of all, you could do intraoperative angiography.
3 The problem with that is the image intensifiers
4 that are available in most operating rooms, as well
5 as the experience of most surgeons performing
6 angiography, isn't going to make that practical.

7 You could do IVIS. You could do
8 intracardiac echo or you could do epiaortic
9 scanning. Again, the problem is you are using an
10 unproven technique with significant user
11 variability to prove an unproven technology and I
12 just don't see that as being a realistic way of
13 evaluating a new device in this situation.

14 Secondly, alternatives to angiography
15 after surgery for follow-up, we have experience
16 with both MRI and EBCT. EBCT is great to show
17 whether the graft is occluded or not occluded but
18 it is not accurate in terms of degree of stenosis,
19 which I think is a harbinger in these grafts of
20 potential occlusion later on so it will not pick
21 that up.

22 MRI, at least with stents in our
23 experience and I assume you could make the jump to
24 anastomotic connectors, creates a flow void in the
25 area where the stent is and you cannot diagnose

1 instant restenosis with cardiac MRI.

2 Thirdly, I think what a real problem of
3 this whole trial design is, is that there are a
4 number of devices out there that have been given
5 previous guidance as to this was a 510(k) pathway
6 and that six-month angiographic follow-up would be
7 an appropriate endpoint to determine safety and
8 efficacy. I realize that the function of this
9 panel, first and foremost, is patient safety and
10 secondly efficacy, but there are a whole number of
11 studies that are out there done that have six-month
12 angiographic follow-up and a very real, practical
13 question is what happens to those studies and what
14 happens to those devices? Do they just get thrown
15 away and we start all over again?

16 I would say that if that is the case, I
17 think that a number of these devices will never
18 make it to market and the companies will go under.
19 I realize that that is not your purview but I think
20 it is a very practical consideration of the problem
21 that exists right now. Thank you.

22 DR. TRACY: Dr. Zuckerman?

23 DR. ZUCKERMAN: Dr. Mack has raised some
24 important points. I would indicate to both the
25 industry and the investigators here that the

1 ongoing studies are developing important data sets
2 for FDA. The question that may need to be raised
3 for some of these companies that are in the process
4 of doing these studies is, depending on what advice
5 this panel gives us today in addition to what is
6 ongoing now, what might be additionally required.
7 But the agency is by no means saying that what
8 companies and investigators have performed to date
9 needs to be thrown out.

10 DR. TRACY: Dr. Kato?

11 DR. KATO: I have a question for Dr.
12 Martin. Can you comment a little bit about the
13 time frame of when you performed your
14 catheterizations in terms of whether 6 months or 12
15 months, in your experience, was of any value or do
16 you have to take this out further in your
17 experience?

18 DR. MARTIN: Well, I will be honest with
19 you, most of my presentation was obviously
20 anecdotal based on the one patient I saw. My
21 problem is the fact that you are subjecting
22 patients to a metallic device without the benefit
23 of platelet inhibition. We know with recent FDA
24 warnings about subacute thrombosis in drug-eluting
25 stents how important that is, and these are

1 patients that don't get Plavix usually
2 intraoperatively because of the bleeding
3 consequences, or may not get it for 48 hours
4 because of the absorption issues.

5 So, the damage is already done when you
6 install these devices. Basically, the platelet
7 adhesion and the basis or the nidus for fibril
8 intimal hyperproliferation starts immediately and
9 that is exactly what I postulated when I saw this
10 first patient, first case. Basically, that is what
11 I think the problem is. I mean, you can't put
12 metallic devices in an aorta. Like Dr. Weinberger
13 said, it is a different pathology virtually.

14 But the point is you are trying to auger a
15 hole, start some sort of clotting cascade with your
16 thumb while you are getting this device ready to
17 implant so you are hoping it clots, but then you
18 don't have anything on board to keep it from
19 over-clotting for instance. I think some of the
20 early cases that had thromboses of the grafts were
21 related to that and I think the nidus begins at the
22 time of implantation with these devices, period,
23 because they are metallic.

24 Now, as far as my experience, I have a
25 very low threshold for cath'g these patients. I

1 mean, basically anginal syndromes, thallium--yes,
2 we are trying to get all these patients back. The
3 problem is that a lot of these patients are taken
4 care by a multitude of cardiologists and, you know,
5 they each have their own individual thresholds for
6 doing non-invasive testing, for instance, in
7 bringing the patients in and, of course, economic
8 considerations as well.

9 But, I mean, my opinion is the same as it
10 was two years ago and now--as I was talking to Dr.
11 White--what do we do with these patients? Can we
12 use drug-eluting stents? Well, not really because
13 a lot of these patients, as was shown in the
14 article two months ago in JAC by Cavendish's group,
15 have a high propensity for restenosis because of
16 the propensity that we see in a lot of stent
17 patients for instant restenosis.

18 So, I hope I have answered your question
19 but the bottom line is I think the injury starts at
20 the time of surgery whether technical or not
21 technical. Every film I looked at on the screen
22 showed evidence of proximal hyperproliferation or
23 intimal endothelial proliferation, every one of
24 them, every graft that was shown up there, and I
25 think that is part and parcel of this whole

1 problem.

2 DR. KATO: Well, my question was really in
3 your experience is it going to be 6 or 12 months
4 for your study.

5 DR. MARTIN: That is a great question but
6 basically we normally treadmill people within six
7 months, use thallium treadmills. Normally most of
8 them are cath'd within six months. That is when we
9 use sort of a definitive marker. But, as I said
10 before, there are lots of patients that can't get
11 back for various reasons and there are patients
12 that have to be done over and over again. Some of
13 these patients have been done six, seven and eight
14 times in various circumstances.

15 DR. TRACY: Dr. Bridges?

16 DR. BRIDGES: I have a question for anyone
17 either on the panel or in the audience. Do we have
18 six-month angiographic data on the Symmetry device?
19 Because one of the key questions here is what is
20 the length of time that it takes to discover a
21 problem with these devices. Is that data available
22 because I haven't seen it published and I was
23 wondering if anyone can comment on that.

24 DR. SLAUGHTER: Mark Slaughter. I can't
25 specifically answer for Symmetry; I think I know

1 the answer. But I would like to respectfully
2 disagree with Dr. Martin so it doesn't go unspoken.
3 To say, you know, a stent is a stent and metal is
4 metal and these things are all the same is not
5 true. And, there is lots of data for the distal
6 anastomosis and Converge. They do get Plavix
7 starting within 24-48 hours. And, the German
8 experience which he put up is 60 days and they have
9 100 percent angiographic follow-up and it is 97
10 percent patency, and they have morphology. That
11 is, there is no stenosis; there is no narrowing.

12 So, to say that because you have something
13 in the lumen and you didn't have Plavix
14 automatically means it has to fail is not true.
15 Within our own experience with now half of our
16 patients back and angiograms 6-7 months, we have no
17 occlusions and we have essentially no stenosis
18 whatsoever. So, to just sort of blanketly say that
19 without, you know, preoperative medical therapy if
20 you have this foreign body it automatically means
21 failures and it will never work is not true.

22 By the same token, there are issues
23 related to that as to material design, where it is
24 placed, and what type of injury occurs are
25 important issues and need to be answered along the

1 way. I think the issue of six months, I mean, I
2 think someone needs to step up and answer but the
3 idea is in any other study so far with these
4 devices at six months the majority of all failures
5 are identified, or you have changes in morphology
6 that would subsequently predict later failure.

7 DR. EDMUNDS: To answer your question,
8 Charles, blood never sees the metal. It sees a
9 granular layer of proteins that are adsorbed onto
10 the metal or any foreign body and that is what the
11 reaction is. The mosaic of those proteins is
12 empirically derived but it never sees the metal.

13 DR. TRACY: Prof. Klima?

14 PROF. KLIMA: Uwe Klima, Hanover. As I
15 presented in my talk, yes, we have postop
16 follow-up, angiographic follow-up with the Symmetry
17 device six months after surgery and this was
18 sufficient enough to detect a real, real problem at
19 the site of the anastomosis even those these
20 patients were asymptomatic. I think this is a very
21 clear message to say that even though you have
22 asymptomatic patients you might have a significant
23 problem at the site of the anastomosis and I think
24 the six-month interval from surgery until you put
25 the patient back into a cath is sufficient to

1 detect any significant problem at the site of the
2 anastomosis.

3 DR. TRACY: Dr. Zuckerman?

4 DR. ZUCKERMAN: Although you have reported
5 very important six-month angiographic follow-up
6 results, what percent of your cohort had six-month
7 angiography?

8 PROF. KLIMA: My total cohort?

9 DR. ZUCKERMAN: Yes.

10 PROF. KLIMA: Well, 67 percent of our
11 MIDCAB patients, 100 percent of the Ventrica
12 patients, 100 percent of the Converge patients and
13 100 percent of the St. Jude patients.

14 DR. ZUCKERMAN: Right, and so in terms of
15 this advisory panel appreciating some of the
16 problems that we have and whether six months can be
17 the gold standard, our paradigm has usually been at
18 least about 80 percent angiographic follow-up,
19 assessed in an independent core lab, to really make
20 statements about the totality of the data. So, I
21 would ask people to consider Dr. Bridges' question.
22 Even if certain data have been spoken of here, has
23 the angiographic follow-up been sufficient or is
24 there still a paucity of data, reviewed
25 independently, to make any decision regarding Dr.

1 Bridges' question about what is the length of
2 follow-up.

3 DR. TRACY: Could I ask you, did you have
4 other assessments on those patients, other
5 non-invasive assessments performed before the
6 angiogram or was this just part of the routine
7 follow-up?

8 PROF. KLIMA: Well, intraoperatively we
9 did an ECHO measurement of flow. However, I am not
10 a big believer about the accuracy of how good the
11 ECHO is. It gives information about is the graft
12 patent or is it not patent but it does not give you
13 very good information about is there 50 percent
14 stenosis, yes or no. If we have a situation where
15 we see that there is no flow in the graft by
16 ultrasound, yes, we go back and we do anastomosis
17 again. All the other information you get is yes or
18 no; it is patent or it is not patent, no more than
19 this.

20 DR. WHITE: She was asking about thallium
21 stress tests. Did they detect your asymptomatic
22 occlusions? Did any other non-invasive functional
23 tests detect asymptomatic occlusions?

24 PROF. KLIMA: Well, there was a clinical
25 interview, so to say, when we brought the patients

1 back about how they would perform. We put them on
2 an exercise stress test and about half of the
3 patients showed ischemia when they were on the
4 bicycle, so to say. But in the clinical situation
5 of daily living--you know, how patients move,
6 whatever they do, they were asymptomatic.

7 DR. WHITE: So, the important data then is
8 that the non-invasive test for ischemia did not
9 completely identify the set of patients who had
10 occlusions.

11 PROF. KLIMA: Exactly. That is why I
12 would really recommend to you to have a six-month
13 angiographic follow-up and if your data look very
14 good I think you can go on with the clinical tests
15 a year later, bring back the patients and see how
16 your patients are doing. In case they are having
17 troubles you certainly have to reevaluate those
18 patients. Either you do it with an angiogram,
19 which I would not do because the angiogram itself
20 has some morbidity, so to say, or some problems
21 coming up sometimes with an angiogram. I think
22 valuable information would also be with an MRI or
23 CT scan.

24 DR. TRACY: Dr. Hirshfeld?

25 DR. HIRSHFELD: I wanted to ask Dr. Emery,

1 because Dr. Emery spoke quite a bit about the
2 technical aspects of using these devices, we have
3 heard reports from a couple of groups of very
4 experienced surgeons and they believe that they
5 have an excessive device failure rate. You have
6 obviously great experience using this device and I
7 wondered if you could offer an opinion as to
8 whether or not the excessive device failure rates
9 may be attributable to technical considerations in
10 terms of how the device is used, and how you would
11 distinguish between that an inherent
12 characteristics of the device.

13 DR. EMERY: I was just thinking that our
14 experience is very much different from Dr.
15 Schoettle's and I would like to get together with
16 him actually and see our differences. On the
17 angiograms, I just reviewed them and there are two
18 out of all the ones I showed. Dr. Martin had had a
19 proximal stenosis; the others were widely open
20 proximally.

21 One of the issues that Ms. Marders dealt
22 with was the dehiscence which was a disastrous
23 event. I think that is a technical problem because
24 the graft is too short and when the patient's heart
25 fills or they cough or retch after surgery, they

1 can pull on that device and pull it off. It is
2 clearly not as strong a suture or as strong as the
3 device Dr. Hausen described. So, I think you need
4 to establish some in vitro criteria for these
5 valves before implantation, not just animal data
6 but the parameters of their strength, just like you
7 proposed for cardiac valve prostheses before they
8 go into clinical trials. I do think the technical
9 issues need to be dealt with because there is no
10 data as to what surgical techniques contribute to
11 failure and that certainly could be one reason,
12 plus regional differences in patient
13 characteristics. It is clear that there are
14 differences in our regional populations--incidence
15 of diabetes, obesity, male versus female, quality
16 of the veins, these issues all come up and these
17 are not addressed in any of the trials that were
18 presented, except for Dr. Mack's trial which came
19 out with diabetes as a very important issue in
20 these device failures.

21 DR. WHITE: Are you still using this
22 device?

23 DR. EMERY: I do for specific indications,
24 people over 80 and in people that have evidence of
25 calcified aorta. Dr. Jim Harte's series is a big

1 proponent for epiaortic scanning before coronary
2 surgery. If people have aortic disease, then I
3 will use it. I do not use the smaller size unless
4 it is absolutely necessary for very specific and
5 patient-oriented reasons. It has its indications
6 and contraindications I think, like we try and do
7 with all of our procedures, processes and
8 medication in the field of medicine.

9 DR. TRACY: Dr. Krucoff?

10 DR. KRUCOFF: I just wanted to step back
11 to the notion that the range of devices we are
12 talking about, in addition to their individual
13 differences, extend from some that are approved and
14 on the market, some that are already in the course
15 of clinical trials that were discussed and
16 designed, and then others that I guess are thinking
17 about the future. It seems to me there are some
18 common pieces that we could focus on them.

19 One of them is called informed consent. I
20 think whether a trial is being planned or is
21 already in motion or whether a surgeon is planning
22 to use a device that is on the market in a human
23 being, we could look at the informed consent
24 document and process across the board and make sure
25 that it contains the appropriate information

1 reflecting current knowledge about the potential
2 risks or benefits.

3 The other is medication. I think there is
4 no question that while maybe loading people with
5 Plavix preoperatively has some downsides, once a
6 gadget like these goes in and hemorrhage and other
7 circumstances is not an issue, anti-platelet
8 therapy across the board, whether you are already
9 doing a trial, is worth probably making sure that
10 you can do as well as possible.

11 Lastly, in the safety assessment is the
12 use of a data and safety monitoring board. What we
13 again very directly from the stenting world
14 recognize is that angiographically, whatever your
15 time window is, you are going to generate clinical
16 events because if at six months people aren't
17 symptomatic but you see 90 percent stenosis, the
18 tendency to do something about that is going to
19 impact the population.

20 So, at least insofar as, for instance,
21 with a good data and safety monitoring board
22 supervision program you could move the angiographic
23 endpoint to a later point and if you have a trial
24 that has done all the paperwork but hasn't enrolled
25 a patient it may be a little different than a trial

1 that has enrolled 210 out of your 215 planned
2 patients. But I think for anything that is in
3 motion, recognizing that an angiographic
4 endpoint--and is pretty clear from what we have
5 heard that detailed morphology of how these things
6 heal is probably going to be an angiographic
7 endpoint somewhere along the line. But if we push
8 that time window later clinical events will occur
9 more consistently with the natural history of the
10 intervention and be less driven by a response to an
11 angiogram.

12 So, I think as long as there is a safety
13 board or there is a safety surveillance mechanism
14 that can make sure patients aren't being harmed,
15 there is a way to envision using angiographic
16 endpoint but pushing out later to allow clinical
17 events to also tell you how patients really respond
18 to these devices.

19 I think from consent to data safety and
20 where you put your time window, frankly, whether it
21 is for a device that is approved and on the market
22 or for a trial that is already enrolling or for a
23 trial that is being planned, that will at least
24 give a backbone that we have learned a lot with
25 stent-like devices that you could think about

1 structuring trials.

2 DR. TRACY: Dr. Edmunds?

3 DR. EDMUNDS: Bob Emery, is it your
4 feeling that the technical limitations of the
5 anastomotic device are much more constrained than
6 they are with the hand-sutured anastomosis
7 proximal, one question; distal for the second? And
8 we are talking just vein grafts.

9 DR. EMERY: Not necessarily, Dr. Edmunds.
10 I think, like anything, we have to learn how to do
11 it and I am sure my first few proximal anastomoses
12 when I was a cardiac surgical fellow were not as
13 good as they are right now after twenty years. I
14 think the same is true with the devices and there
15 may be a need for more mentoring on the first
16 several cases so someone can suggest what is good
17 and what is bad, much as is done with other more
18 complex devices. These devices seem intuitively
19 simple but obviously we are coming to the
20 conclusion here that they are not intuitively
21 simple. They have a lot of subtleties in their use
22 both in terms of healing, deployment and events
23 that occur, and those all need to be addressed, not
24 just the fact that you push a button and you have
25 an anastomosis. That is clearly very simple.

1 These things, like any other techniques in our
2 medical profession, we need to learn how to do
3 them.

4 DR. EDMUNDS: But you showed us angiograms
5 in which there were problems, that you needed to
6 put it on the side of the aorta towards the
7 pulmonary artery, which you don't need to do for
8 hand-sewn anastomoses. It has to be totally at
9 right angles. It can't be at a little bit of a
10 hood. Those are the constraints that seem to be
11 much more confining than would be in a hand-done
12 anastomosis. Is that true or not true?

13 DR. EMERY: Yes, and they have to be
14 defined. That is what I was talking about for some
15 of the in vitro tests and even the in vivo tests
16 before the application to our clinical patients.
17 It is not just patency that is important; it is how
18 you use it and define the limitations of the device
19 as part of the preclinical issues that go into
20 preparation for a clinical trial.

21 DR. EDMUNDS: Therefore, the device
22 anastomosis is less robust than is a hand-sutured
23 anastomosis. It is less tolerant of small error.

24 DR. EMERY: That is possible. You know,
25 it is possible but, on the other hand, sewing an

1 anastomosis in a diseased aorta is less tolerated
2 by the patient and there is value in these devices
3 even though they may be more complex. Also, on one
4 angiogram I showed side-by-side anastomoses where
5 one was hand-sewn because I was over-aggressive in
6 pulling the vein graft and pulled the connector
7 right off the aorta, which, you know, caused a
8 little bit of a stir in the operating room, of
9 course. The other connector was just fine so I
10 just hand-sewed the other anastomosis. The patient
11 returned at three months with both occluded. I
12 can't explain why there was a difference in
13 something like that occurring.

14 DR. TRACY: Dr. Yancy, then Dr.
15 Weinberger, Dr. Kato and then lunch.

16 DR. YANCY: Putting some thought into a
17 clinical trial design and thinking of endpoints
18 that would be evaluated in that design, I have not
19 heard a discussion this morning about the presumed
20 advantages for this device, specifically the time
21 of the proximal anastomosis, decreasing that
22 variable to the extent that that is clinically
23 relevant. The second would be the CNS event
24 circumstances because that was purportedly one of
25 the major reasons for developing the technology.

1 Are we to assume that that has been proven or are
2 those issues? Because in my judgment, if we are
3 addressing questions of efficacy with regards to
4 the integrity of coronary perfusion and we have
5 missed the purported benefit, then I think there
6 are some serious questions that have to be
7 addressed that are more global than trial design.
8 Are there people in the audience that can
9 specifically comment on the CNS issues and the
10 timing issues, especially if we follow the most
11 recent discussion that more care and consideration
12 needs to be made in generating these proximal
13 anastomoses?

14 DR. BRIDGES: I don't think there is any
15 data other than speculation that I am aware of,
16 unless somebody else knows otherwise, that prove
17 those points that you made. I mean, I think the
18 supposition is that there would be less incidence
19 of stroke and that may very well be true but I
20 don't think that that data is available.

21 DR. YANCY: Well, I respect that but there
22 is a comment in our packet that 30,000 of these
23 devices have been implanted and we have any number
24 of events reported to the FDA, and I think if we
25 are dealing with supposition and speculation the

1 trial design needs to be a bit more rigorous than
2 we may have thought earlier today.

3 DR. FERGUSON: Can I respond to that?

4 DR. TRACY: Yes.

5 DR. FERGUSON: I think you are mixing a
6 little bit apples and oranges because the genesis
7 for many of these mechanical devices has been
8 MIDCAB and off-pump where putting on an anastomosis
9 by hand is much more difficult, sometimes
10 impossible. It obviates the use of an aortic clamp
11 which is not possible if you are doing it
12 hand-sewn. So, I think that responds to the issue
13 of CNS issues because those issues are due to
14 clamping the aorta so you can do the anastomosis of
15 the hand-sewn.

16 DR. AZIZ: But then they should be able to
17 look at the cases they have done and see what the
18 incidence of neurological problems is in the 30,000
19 cases or so that have been mentioned.
20 Theoretically, it should be less than what was
21 expected. I mean, somebody should have that data.

22 DR. TRACY: Dr. Weinberger?

23 DR. WEINBERGER: I would like to sort of
24 echo and expand a little bit on what Dr. Krucoff
25 said. I think that in the past couple of years the

1 FDA, together with the cardiology community, has
2 developed a paradigm for analyzing new endovascular
3 stents. If you go back to 2000 or 2001 in the
4 early days of drug-eluted stents, the studies that
5 were designed included both clinical endpoints and
6 angiographic substudies, and the angiographic
7 substudies were done to power angiographic
8 endpoints which are completely different than
9 powering clinical endpoints.

10 I think that both are important. I think
11 the surgeons and a lot of us are focusing on
12 morphological endpoints right now, and I think the
13 study done to get a morphological endpoint when you
14 have a continuous variable could be done with a
15 much smaller number of patients.

16 On the other hand, the FDA standard in the
17 past has always been clinical benefit to patients
18 or at least clinical equivalence to previous
19 devices and I think we need to gather clinical
20 endpoints for those particular pieces of
21 information. So, what we really need is to define
22 what the rates of major adverse cardiac events are
23 in hand-sewn saphenous vein surgery at one year or
24 two years, some time point, and use a clinical
25 endpoint for that piece of the information.

1 We are all sort of focusing on an
2 angiographic endpoint which, at best, is going to
3 be a surrogate to suggest what is going to be
4 happening clinically and probably should represent
5 only a substudy of an ultimate study that is done
6 to approve this.

7 Now, I am sure that the manufacturers are
8 going to howl and say that we are really subjecting
9 this to PMA standards, but I think that given where
10 we are and the fact that we have introduced what
11 appears to be problems at unexpectedly high
12 frequency the data necessary to at least inform
13 patients that we are not exceeding the previous
14 event rates are necessary.

15 DR. TRACY: Dr. Kato and then we should
16 really break for lunch.

17 DR. KATO: Why don't I let Dr. Frater talk
18 and I can reserve my comments until after lunch.

19 DR. TRACY: That is fine.

20 DR. FRATER: I just want to try to respond
21 to a couple of the issues that were raised. You
22 did raise the question of angiographic evidence at
23 6 months on the connector in a miscellaneous group
24 of patients in which angiograms were available.
25 There were 221 of them. The occlusion rate between

1 6-12 months was 20 percent. The range was 2.3
2 percent to 58.3 percent. That is one of our
3 problems, dealing with a range that wide. How can
4 you possibly account for it?

5 It is perhaps fair just to mention that
6 there are studies already published in abstract
7 form in which the outcomes judged by MACE followed
8 by angiography are nowhere near what you heard this
9 morning. I will just quote one from Brady,
10 University of East Carolina, with 400 patients with
11 650 veins and they found, based on clinical events,
12 3 veins that had problems of occlusion or stenosis.
13 The difference between these various studies is
14 enormous.

15 In terms of neurological episodes, those
16 400 patients studied from 2001 to September 2003
17 had a 1.7 incidence of neurological adverse events
18 in the postoperative period. Dr. Schoettle could
19 tell you about his connector patients and their CNS
20 adverse events but I will leave him to do that if
21 he wishes.

22 So, I won't spend more time; we clearly
23 don't have time, but there are other similar good
24 series to report that we are faced with
25 extraordinary diversity in results from competent

1 people.

2 DR. TRACY: Thank you. Why don't we break
3 for lunch and let's be back here at 1:20.

4 [Whereupon, at 12:20 p.m., the proceedings
5 were recessed for lunch, to reconvene at 1:20 p.m.]

1 A F T E R N O O N P R O C E E D I N G S

2 DR. TRACY: If everybody would please take
3 their seats, I would like to resume the open
4 committee discussion because I know there are a lot
5 of additional issues that need to be addressed this
6 afternoon.

7 Again, I would like to thank presenters
8 for hanging around to help us clarify some issues.
9 It is very helpful to get the input of everybody
10 who has made the time to come here. This is a
11 difficult issue that we are wrestling with because
12 there are an awful lot of variables that are
13 involved in the discussion that we are here to
14 have. At the end of the day we are expected to be
15 able to answer some questions or begin to answer
16 questions that have been posed to us from the FDA.

17 I would just like to at this point sort of
18 summarize what I see as some of the complexities
19 that we are looking at right now. I think there is
20 a variety of variable here that we are discussing.
21 One is what type of devices are we talking about.
22 There is a variety of devices that are out there
23 and there is a variety of devices that are coming
24 down the pike.

25 I think an issue that we have to address

1 is, is the analysis of each of these different
2 types of devices going to be the same. My
3 suspicion is that that will not be true but is
4 there some sort of paradigm we can come up with
5 that would help to analyze what type of process
6 should be gone through to look at different types
7 of anastomotic devices?

8 I think one thing that is clear that we
9 have heard so far today is that the target vessel
10 is very important--whether something is being
11 connected to an LAD, whether it is an aortic
12 anastomosis or artery or vein anastomosis. these
13 are variables that are inherent in this type of
14 device and the analysis or process to look at these
15 things may be different if we are talking about
16 distal or proximal anastomosis, and how do we
17 handle that in a given patient population?

18 I think we have heard an awful lot about
19 operative variables, some of which may be
20 technically overcomable, some of which may not be.
21 Maybe a MIDCAB is very different from a coronary
22 bypass graft placed through a thoracotomy. How do
23 we handle that?

24 So, let's focus on some of those things as
25 we think about our trial designs and what

1 parameters we can look at both acutely, what is the
2 appropriate mid-term analysis on these patients,
3 and what does long-term analysis and long-term
4 follow-up really mean. Are we looking at an
5 anatomic outcome or are we looking at a functional
6 outcome? And what other outcomes are appropriate
7 to consider?

8 One of the original reasons perhaps for
9 developing any of these anastomotic devices was to
10 avoid some of the neurologic outcomes that come
11 from cross-clamping the aorta. So, what other
12 outcomes do we need? We sort of focused on
13 angiographic but are there other neurologic or
14 other types of outcomes we need to think about?

15 The final area that I would like people to
16 sort of focus in on is that it does seem to be
17 almost a separate cardiovascular surgeon training
18 program to learn how to use these things, and what
19 is it that needs to be done to train appropriately
20 for the use of these devices?

21 So, if you can try to keep in the back of
22 your mind the questions that have been posed to us
23 and let's try and focus our way through some of
24 these issues.

25 DR. KRUCOFF: As a point of order, I don't

1 know if this is inappropriate or appropriate but
2 would it be possible to actually go through the
3 afternoon discussion by starting with the questions
4 instead of doing the discussion and then ending
5 with the questions? We are not voting.

6 DR. TRACY: We can do that but I think
7 there may be some additional issues. You have the
8 questions listed here so if there are additional
9 things that you think are relevant before we answer
10 the questions, I want to give plenty of time for
11 discussion and input from the audience on this.
12 Yes?

13 DR. MACK: Michael Mack. Do you want me
14 to start?

DR. TRACY: Sure.

15 DR. MACK: I have a couple of comments
16 based upon what you just said. The first alludes
17 to what you said, you need a whole separate
18 training program to implant these devices, and a
19 lot of attention was focused on Bob Emery's
20 presentation about the permutations that are
21 necessary for this. I think this addresses the
22 issue that implantation of a medical device by
23 surgery or by a catheter is not a pill and you
24 don't just give a pill and a placebo and that is
25 the variable. There are intricacies associated

1 with the implantation of a device. It is not that
2 the use of an anastomotic device is necessarily
3 more constraining or that the margin of error is
4 narrower, but it is just different and just because
5 you have been sewing coronaries for twenty years
6 doesn't mean you can automatically put in a device
7 the first time.

8 I will go to the analogy of drug-eluting
9 stent devices. There are a whole bunch of adverse
10 events with subacute thromboses that happen
11 immediately upon approval of this. It probably had
12 to do with the fact that it is a device, a
13 particular type of stent that most cardiologists
14 weren't implanting. It was just different. It is
15 stiffer. The delivery platform wasn't as usable.
16 How high do you inflate the balloon? Do you
17 overlap stents? What is the amount of coverage
18 that you have? All those little intricacies just
19 had to be brought up to speed and were different
20 than cardiologists had been doing before.

21 It is the same thing with anastomotic
22 devices. you find out that you have to put a
23 suture to tack a vein graft to keep it straight
24 coming off, whereas you don't with a sutured one,
25 and there are little permutations like that that

1 you find when there is broader experience with it.

2 The second has to do--it was either Chris
3 or Clyde that mentioned this, we talked a lot about
4 adverse events about the procedure. Well, what
5 about the benefit? Well, I think there is benefit
6 potentially both with the device itself but what it
7 is ultimately going to lead to.

8 One is the device itself. Maybe you
9 actually have a more reliable anastomosis. That is
10 not proven. But this is a brick in a wall of less
11 invasive surgery in general and both minimal access
12 surgery and off-pump surgery. The reasons for that
13 are the following: The reason that minimal access
14 surgery is so limited is the most difficult thing
15 to do through a little incision is to sew. If
16 minimal access surgery is every going to get any
17 place, this is a necessary brick in the wall to
18 move that along.

19 Similarly, you can argue whether there is
20 a benefit to off-pump surgery or not. The
21 preponderance of evidence would seem to indicate
22 that there is. But it is stuck at 25 percent of
23 CABG in the United States. Why? Because most
24 surgeons aren't comfortable putting the heart on
25 end for 20 minutes to sew on a beating heart. The

1 potential benefit of an anastomotic device, if
2 proved to be safe and efficacious, is immediately
3 catalyzing off-pump surgery. If it takes a minute
4 to do the anastomosis most surgeons are going to
5 feel comfortable with that.

6 So, I think the ultimate benefit is that
7 it is a potentially more reliable anastomosis, not
8 clamping the aorta, not having neurological events
9 but ultimately what it allows this whole field to
10 progress to is that it is not going to be the same
11 operation as for the last fifty years.

12 DR. TRACY: Thank you. With that as a
13 jumping point and trying to go along a little bit
14 with the lines Mitch indicated, we have heard the
15 historic comparison to CABG and anastomosis as sort
16 of the gold standard. Is that the right gold
17 standard or is this technique different enough that
18 we don't use sutured CABG anastomoses as the gold
19 standard? Dr. Aziz?

20 DR. AZIZ: Just before answering that
21 question, you know, we have a technique which is
22 very safe. It is not perfect; the vein grafts may
23 not last long but I think that is usually related
24 to vein biology. So, I mean, we have a certain
25 standard and if we are going to adapt a new

1 technology I think at a minimum it can't be worse
2 than what we have, and I think we have to be
3 careful on that. I think we need to sort of keep
4 that in the back of our mind. For example, this
5 particular device that was presented, the St. Jude
6 device--clearly, the results are worse than what we
7 have with hand-sewn. So, I think, you know, we
8 have to make sure that whatever we do is not worse
9 than what we have. It is not like we are in a
10 void.

11 Regarding the sort of controls, as someone
12 else had also mentioned, using the same patient as
13 his own control would obviate the need to have
14 randomized studies in the sense of, you know,
15 having two different groups of patients. You could
16 do one anastomosis with one of these techniques
17 either for the distal or proximal and one doing it
18 the old-fashioned way or the regular way with a
19 suture technique. When you evaluate by angiogram,
20 that in itself could be its own control.

21 DR. FERGUSON: I have kind of done a 180
22 on this because I was impressed with--I think it
23 was Bob Emery but others too. I came in with the
24 idea that using the traditional statistics for a
25 control, matched control would be a good thing to

1 do. But the more I have listened today, the more I
2 think that the patient population has changed so
3 much in the last twenty years that those historical
4 controls are truly historical and probably not
5 relevant to what is going on today given the kinds
6 of patients that are being operated on today, and
7 so forth. I think, Bob, you mentioned something
8 about that too. So, I would only comment that if
9 we truly need control for experimental devices,
10 then I think that we need to make those controls
11 modern-day standard heart/lung operations.

12 DR. TRACY: Dr. Blumenstein?

13 DR. BLUMENSTEIN: I am really seconding
14 what you are saying and maybe stating a little more
15 explicitly that I don't see we can do anything by
16 randomize in some sense. I think that you have to
17 have the power of randomization to give you the
18 stochastic equivalence between those treated with
19 the experimental intervention and some kind of a
20 control group. Whether it is matched or whether it
21 is a two-group study, and so forth, will have to be
22 discussed but there are just too many extraneous
23 factors that cannot be controlled. I am just blown
24 away by the number of them. So, you have to have
25 randomization to control the experimental

1 intervention and let those other factors be
2 stochastically controlled with randomization.

3 DR. EDMUNDS: Just to maybe start a little
4 bit of discussion on this issue because this really
5 is the heart of the issue, let's just look at the
6 problems of a non-randomized control group and what
7 are the alternatives. Well, propensity matching
8 and univariate and multivariate logistic
9 regression. But a lot of the technical problems
10 are not things that we collect data on, as Bob
11 Emery has brought out. Yet, they are very relevant
12 to the success or failure of patency, let's say,
13 and that will confound any kind of propensity
14 matched control group or taking it out of a
15 logistic regression equation because you don't have
16 the data in there in the first place. So, I can't
17 think of a control group other than a concurrent
18 prospective, randomized control. Maybe somebody
19 else can but I can't.

20 DR. TRACY: Mitch?

21 DR. KRUCOFF: Because of the breadth of
22 range we are looking at, we really do have to be
23 careful about babies and bath water. I think if we
24 take this as a new device, then a characteristic
25 first question would be is it safe? The first

1 in-man type of experience would be one way to go.
2 I think we have seen at least one elegant
3 illustration of using a patient as their own
4 control where you have even numbers of grafts that
5 can be characterized and randomized and carefully
6 followed angiographically.

7 That is living in the morphologic side of
8 are we hurting people anatomically by putting these
9 things in. That would also be a very good
10 opportunity for a bunch of smart surgeons to look
11 at what are the technical features that we would
12 want to capture in a larger, more definitive trial
13 and explore that a little bit in a small but very
14 intensively designed study. I think the one
15 feature there that certainly stood out for me out
16 of this morning's presentations is the potential to
17 use a patient as their own control over a series of
18 non-LAD grafts.

19 I think as you move into a more definitive
20 "okay, this thing is safe," now is it effective,
21 and is it effective for what? Is it effective as a
22 stepping stone toward fully robotic surgery, or is
23 it effective to not have to cross-clamp the aorta
24 in 80 year-old people with crunchy aortas? I mean,
25 it depends on the question which the FDA is very

1 particular about; it depends on the label; it
2 depends on the indication.

3 But I think if we do this in steps there
4 are points where patients can serve as their own
5 control, where historical descriptors may be
6 useful, and then I think ultimately there is a
7 point where you are going to have to do real
8 meaningful and probably randomized trials to show
9 that it is effective, building along the way.

10 DR. HAUSEN: Bernard Hausen, Cardica.
11 What is the worry, Dr. Ferguson? That we are over-
12 or underestimating with our historic controls the
13 true incidence of occlusions?

14 DR. FERGUSON: It is a different patient
15 population. That is my point. We had this problem
16 yesterday a little bit about using a group of
17 patients that are 20, 25 years out as our control
18 group. I don't think we can do that. Somehow or
19 other, we have to have concurrent. Now, how that
20 is done is up for grabs here I think. Go ahead.

21 DR. HAUSEN: Let's say Dr. Frater is right
22 and patients are getting sicker; they are
23 definitely getting older; the veins are getting
24 poorer and all this will result into poorer outcome
25 of vein grafts. You could hypothesize that

1 nowadays, if I did a brand-new control group with
2 vein grafts hand-sewn, my results would be worse
3 than what Dr. Mack presented to us.

4 DR. TRACY: I think that you are right but
5 the problem is that we don't have a proximal
6 control here. We don't have something that we can
7 look at that we feel is comparable.

8 DR. HAUSEN: I understand that but if you
9 are worried that when I am comparing my average
10 patency of a device to a control, the worry would
11 be that it gets approved because it is
12 statistically not worse or similar to control where
13 my control would underestimate the true prevalence
14 of occlusions. I don't think that is the real
15 world.

16 DR. FERGUSON: That is not a definition of
17 a control in my world.

18 DR. HAUSEN: I know. I know, but from a
19 regulatory point of view you are trying to prevent
20 products, or the FDA is trying to prevent products
21 from getting onto the market that look better than
22 they really are and that pose a patient risk. So,
23 if you have a control group that is worse than what
24 the real world is showing us right now--

25 DR. BRIDGES: Can I interrupt? I think we

1 got the point. I think part of the problem here
2 though is that one of the safeguards of a
3 randomized trial, or at least a prospective
4 trial--I mean, yes, perhaps the patients today are
5 "worse" than the patients were then but what is to
6 prevent someone from doing a study where they are
7 selecting patients from today who aren't worse and
8 applying the technology selectively to groups of
9 patients who are not worse than some group of
10 historical controls? That is the whole reason why
11 prospective trials--it is just to mitigate against
12 that kind of deficiency.

13 So, to simply stand up and say, well, you
14 know, the patients that we are operating on today
15 are worse than they were before, therefore, we
16 don't need prospective trials I think is overly
17 simplistic.

18 DR. TRACY: Dr. Sapirstein?

19 DR. SAPIRSTEIN: Given the template or the
20 sample I provided where we use a 95 percent point
21 estimate of patency and a lower confidence limit of
22 5 percent, 95 percent lower confidence limit of 90
23 percent, how much better can you do on a randomized
24 control? The only reason I bring this up is
25 because of the sample size that you would require

1 and the difficulty of subjecting patients to
2 angiography.

3 DR. KRUCOFF: Come on, Ralph, how are you
4 not going to submit patients to angiography? In
5 this status of a new device how are you going to
6 avoid doing angiograms in these patients? You are
7 going to have to do angiograms in these patients.

8 DR. SAPIRSTEIN: Yes, I think so. I think
9 that you have to do angiography on a new device but
10 on standard procedure, a LIMA to LAD, are you
11 justified in doing an angiographic evaluation of
12 the control? Maybe you are. It is a point of
13 discussion.

14 DR. TRACY: Dr. Edmunds?

15 DR. EDMUNDS: Before you stand down,
16 because it is very expensive to do randomized,
17 prospective, controlled trials, you have to decide
18 in a power analysis how much of a difference is it
19 that you want to see, and presumably a composite
20 primary outcome, to be meaningful, and do your
21 power analysis on that basis. Now, if you really
22 want to say that a one percent difference
23 multiplied by 260 million people, not all of whom
24 have coronary-artery disease, is what is needed,
25 therefore, you are going to have to have a very,

1 very small difference, like half a percent, nobody
2 can afford that trial.

3 DR. SAPIRSTEIN: Exactly.

4 DR. EDMUNDS: So, that power analysis and
5 what you define as a meaningful difference is
6 critical I think to any kind of regulation that you
7 set up for these companies.

8 DR. SAPIRSTEIN: Absolutely. That is why
9 we put up that template. With 80 percent power, a
10 5 percent alpha, can we do--

11 DR. EDMUNDS: Well, what difference?

12 DR. ZUCKERMAN: [Not at microphone;
13 inaudible]... the expected observed rate for a
14 venous study would be 95 percent with 80 percent
15 power, alpha 5, and setting a delta of 0.05,
16 meaning that you are projecting a lower rate--

17 DR. EDMUNDS: I am not sure.

18 DR. ZUCKERMAN: --you need 150 patients.
19 If we expand that delta to 0.07 it could be a lower
20 rate of 88 percent and you would only need 59
21 patients but perhaps the delta should be tighter
22 [not at microphone; inaudible]...consider this
23 approach given that for having just normal CABG
24 patients come back for follow-up angiography, as
25 Dr. White was alluding to, perhaps can be very

1 difficult.

2 DR. EDMUNDS: Well, the five percent
3 difference is probably in the realm of
4 reasonableness.

5 DR. TRACY: Let me be clear what we are
6 talking about. Have we moved away from a
7 randomized trial to a trial where we are using
8 historic 85, 95 percent?

9 DR. WHITE: Can I support that?

10 DR. TRACY: You support moving away and
11 using the historic CABG data?

12 DR. WHITE: Right. The reason I say that
13 is that obviously I am an angioplaster so I am
14 anti-intellectual, but the point is that I do like
15 randomized trials and I won't argue with you that
16 that is the best of all worlds, but I think we can
17 gain confidence that we are doing as well or
18 better, and we have such extensive historical data,
19 as Dr. Mack presented to us today and other data
20 that is available, that if the manufacturer chooses
21 a relatively high bar for patency, it could be done
22 as a single arm because we would be sure that they
23 would be 90 percent or better and that would
24 certainly meet historical controls and it would be
25 hard to argue that they were worse than historical

1 controls.

2 If the manufacturer, on the other hand,
3 decided that that bar was too high and they wanted
4 a randomized, controlled trial because the control
5 group really was going to perform at a lower rate,
6 then I think that is a decision we ought to leave
7 to the investigators and to the company.

8 DR. TRACY: Dr. Blumenstein?

9 DR. BLUMENSTEIN: I am not sure exactly
10 what this represents. I assume this means a
11 single-arm trial with a criterion of success of,
12 say, 95 percent and precluding the possibility--it
13 is a sort of non-inferiority situation--precluding
14 the possibility that the outcome is truly 90
15 percent or worse with a 5 percent delta, or
16 whatever.

17 The danger here for a company undertaking
18 such a trial with such a high bar is that if you
19 take that sample size of, say, 150 and you have
20 guessed wrong about what the success rate is, then
21 you are in trouble.

22 DR. WHITE: And we have been here once
23 before when companies come back to us, missing
24 their bar and asking us to make that exception--

25 DR. BLUMENSTEIN: Right.

1 DR. WHITE: --and it is up to us not to
2 make that exception if they miss that bar.

3 DR. BLUMENSTEIN: But there are two
4 reasons why this is a real problem. Number one is
5 that if you had started with a lower success rate,
6 supposing it was reasonable to do so and we all
7 accepted that it was reasonable to do so, the
8 sample size goes up as the base reference success
9 rate moves closer to 0.5. So, by putting the
10 success rate close to 1.0, then you are getting a
11 smaller sample size than you would had you put the
12 base rate closer to 0.5. So, by setting the bar
13 high, which is a good thing, you are in fact
14 anti-conservative with respect to the sample size.

15 DR. TRACY: Dr. Zuckerman, you look like
16 you have something to add here.

17 DR. ZUCKERMAN: I just want to state that
18 Dr. Blumenstein has correctly summarized what this
19 trial design is trying to do, say, for a LIMA
20 trial. We are trying to show that the new device
21 has an observed performance rate or patency rate
22 above 0.90. Therefore, as he was--

23 DR. TRACY: Dr. Zuckerman, stay near the
24 microphone.

25 DR. ZUCKERMAN: --if the observed rate is

1 0.93 instead of 0.95 they can do a sample size
2 calculation and it will go up. But the question
3 is, is 0.90 the right bar for this performance
4 goal, Dr. White.

5 DR. WHITE: I would be very happy. Given
6 what we know about the historical controls for
7 hand-sewn grafts, the meta-analysis and data that
8 we have looked at, that if I was confident that the
9 device could do 90 percent or better at one-year
10 patency or six-month patency, I would feel pretty
11 comfortable about that.

12 DR. KRUCOFF: How about 89? How about 88?
13 How about 87? I mean, come on, Chris, the
14 likelihood of creating an ambiguous data set that
15 ultimately then would occupy all of us for a
16 day--you know, I think we have to recognize that
17 the likelihood for ambiguity if we court certain
18 marginal design structures is so intense--

19 DR. WHITE: I don't have a problem with
20 89.

21 DR. KRUCOFF: What if they don't make
22 their endpoint? What about that they find out that
23 in 50 year-old women, as a retrospective subgroup,
24 it was 100 percent? You know, the vulnerability to
25 starting in the wrong place to then, in this venue,

1 having a long discussion about completely ambiguous
2 data is huge and I think the key is don't start in
3 the wrong place.

4 DR. TRACY: I think the problem that will
5 come up is that there may be times where it would
6 be inappropriate to have a historic control, but it
7 may be just as inappropriate to have the patient
8 serve as their own control--a single vessel study.
9 There are all kinds of reasons why either one of
10 these is going to be not optimal fit for any given
11 device that is being tested. So, I think there has
12 to be a little bit of breadth. I think that gets
13 back to the idea that you cannot possibly have an
14 exact study that has to be done by each one of
15 these different devices because they are different
16 devices. So, the patients that are going to be
17 enrolled are going to be different from one study
18 to another. I don't know how to get around that
19 except to accept different controls in different
20 studies.

21 DR. EDMUNDS: But a Type 2 error is just
22 as bad as a Type 1 error.

23 DR. TRACY: Dr. Emery?

24 DR. EMERY: I think two issues come to
25 mind. One is that patients serving as their own

1 control will likely only be utile in proximal
2 anastomotic studies and this would separate the
3 distals from the proximals.

4 The second thing is that, on the side of
5 patients as their one control, you have to be sure
6 you are asking the proper question from the study.
7 For instance, there may be biology of aging. We
8 know from prosthetic valves that a bioprosthesis in
9 an elderly patient will last a lot longer than a
10 bioprosthetic valve in a younger patient. The same
11 may be true for connectors because of the biology
12 of wound healing and inflammation. We don't know
13 that. So patients serving as their own control
14 would be a good means of separating out the
15 patients at risk and the patients not at risk based
16 on their individual biology.

17 In our paper published in Circulation, the
18 younger patients that had the aggressive restenosis
19 were at risk for the connectors. On the other
20 hand, if you are looking at stroke risk, which has
21 been brought up, that is not a good study for
22 patients as their own control because, perforce,
23 you have to put a clamp on the aorta for all
24 patients in the study to randomize them and that is
25 one of the highest risk things that we are trying

1 to avoid with connectors, particularly in the aged.

2 So, it is important to address the
3 question before the study is designed to examine
4 what answer you are looking for, the value to the
5 patient population or the patency and effectiveness
6 of the connector itself.

7 DR. BRIDGES: I have a question related to
8 what Wolf Sapirstein said, which is that one
9 question the panel should consider is, if we are
10 advocating a randomized trial, how do we feel about
11 the risk to control patients being
12 subjected--assuming that the patient can't serve as
13 his own control for example, how would we feel
14 about control patients being cath'd at six months
15 or one year?

16 DR. TRACY: I think that was Dr. White's
17 point, that it is hard to get patients back. We
18 have historic data on those patients regarding
19 long-term patency.

20 DR. FERGUSON: Some of the hardest data we
21 have seen today is from Dr. Klima, at Hanover, and
22 they don't seem to have a problem over there when
23 they set up a design trial for any of these issues.
24 I am not quite sure why we do.

25 DR. YANCY: Well, just to take some

1 exception to that, any study design that comes up
2 has to clear local IRBs that are under increasing
3 scrutiny with regards to what they will subject
4 particularly control individuals to. So, I think
5 that Chris' point is not an inconsequential
6 consideration. Certainly from the stent era we do
7 use MACE as an endpoint and I wonder how that
8 compromises this study design if we use MACE.

9 DR. TRACY: So, we have the scenario of a
10 historic control. We know what the patency rates
11 are. We have the scenario of the patients serving
12 as their own control, which will limit the type of
13 device that can be studied. There will be limits
14 there. And, we have a scenario of a randomized,
15 controlled study group that would be brought back
16 for intervention. I think each of these has some
17 problems. Dr. White?

18 DR. WHITE: I think that we need very hard
19 endpoints. I think hard endpoints make for
20 smaller, more firm studies. So, I really think we
21 want patency. I think we know already that
22 asymptomatic patients show up this way. We have
23 heard that perhaps non-invasive testing isn't an
24 adequate screen so we want the hard endpoint. If
25 we ask for patency at whatever interval, six

1 months, one year, then I think we take care of Dr.
2 Edmunds' concern about the procedural issue because
3 if you are having closure at procedure you will see
4 the closure at six months as well, and you can then
5 backtrack that.

6 I think if you do single-arm studies with
7 objective criteria, which has its own problems with
8 a slippery slope, at least you can subject those
9 patients to a hard endpoint of patency, and I think
10 that minimizes the number of patients who are
11 studied. It gives you the hardest endpoint for the
12 follow-up. And, I think it really puts the
13 patients who are at risk for this problem as the
14 subjects of the study as opposed to the control
15 population.

16 DR. TRACY: Dr. Hirshfeld?

17 DR. HIRSCHFELD: I think there are really
18 two questions. The first question is do the
19 devices perform equivalently to standard hand
20 suturing? That, it seems to me, can be answered in
21 a study with a relatively small population of
22 patients and a study design analogous to what the
23 Hanover people presented. That has a very hard
24 endpoint and has a great deal of statistical power,
25 and that could settle this question with a

1 relatively modest number of patients, I would
2 think.

3 The second question is the question of
4 what global benefit does the patient derive from
5 this, and those are all surrogate endpoints. They
6 are neurological endpoints and they are MACE
7 endpoints. That could be addressed with a larger
8 randomized trial in which all the endpoints that
9 were collected were the surrogate endpoints in
10 terms of perioperative neurologic events and
11 long-term cardiovascular MACE and wouldn't really
12 require angiographic follow-up.

13 So, it would seem to me that one possible
14 approach to this issue and I think Dr. Mack's
15 comment about the importance of continuing
16 development in this field, independently of how the
17 current generation devices actually sort out, is a
18 very important point. But the issue of whether or
19 not we are hurting patients in terms of graft
20 patency can be answered with a relatively small
21 trial, and then the issue of whether or not
22 patients are deriving benefit in terms of surrogate
23 endpoints could be answered with a larger but less
24 complicated trial.

25 DR. TRACY: Dr. Blumenstein, do you have

1 any comments on that type of a design?

2 DR. BLUMENSTEIN: I am sorry, I wasn't
3 paying attention. I am doing some trial size
4 computations.

5 DR. HIRSHFELD: I am sure that if you try
6 to do a trial size computation there is no way you
7 can concentrate on anything else. But what I was
8 thinking was that if you did a study design
9 analogous to the study that was described from
10 Hanover in which each patient serves as their own
11 control for the performance of the device in terms
12 of graft occlusion, that is a study that has all of
13 the variables controlled through the randomization
14 process and has the highest possible degree of
15 precision in terms of the outcome event, and should
16 be able to be completed with a relatively modest
17 sample size.

18 DR. BLUMENSTEIN: That is exactly what I
19 am trying to compute.

20 DR. HIRSHFELD: But that doesn't answer
21 the other question which is are patients better off
22 in terms of events because of the use of the
23 device, or are they equivalently off, or are they
24 worse off. That would require a larger trial in
25 which patients are randomized between the device or

1 not the device and the endpoints are all the
2 clinical endpoint events, neurological events and
3 cardiac MACE.

4 DR. BLUMENSTEIN: Well, we can discuss
5 this later I suppose but when designing the trial,
6 in order to size the trial you have to pick a
7 primary endpoint. Basically, I think what we are
8 stuck with here in this discussion, without getting
9 into a great amount of detail, is the dichotomous
10 outcome or either success or failure assessed at,
11 say, seven, eight months down the road, giving a
12 window for six-month angiography or something like
13 that.

14 The secondary endpoints or the other
15 benefits that might accrue then have to be assessed
16 in the context of the sample size computed for that
17 structure of the study and within the framework of
18 that study.

19 I think what you are asking is, if you
20 were to do a two group study, then you would have a
21 much better idea of the overall impact of the
22 benefit of this and a two group study might be
23 something like a treatment failure-free survival
24 study where you are using a time-to-event endpoint
25 and you randomize into two groups. In that case,

1 what you gain is that the patients in each group
2 are assessed with respect to all the outcomes and
3 are not muddled by having one vessel randomized to
4 the experimental procedure and one vessel
5 randomized to the control procedure. Maybe that is
6 a good follow-on study, a postmarketing type study.

7 But it seems to me that if you can achieve
8 the matched vessel type of study where each patient
9 is serving as their own control. If that study is
10 doable administratively, if you feel that you can
11 do the quality control, implement the intervention,
12 the randomization, etc., then I think that gets at
13 the answer about whether there is success with
14 respect to the vessels and a comparison to a
15 control which would take away the host factors.

16 DR. SAPIRSTEIN: [Not at microphone;
17 inaudible]...safety together with effectiveness,
18 and much of the effectiveness of these devices has
19 yet to be determined by virtue of facilitating
20 other changes to the performance of CABG which
21 haven't in themselves been determined. For
22 instance, off-pump, beating heart, MIDCAB and all
23 those sort of things. One of our most important
24 considerations is patency right now. We are
25 interested in safety and MACE events but right

1 now--well, not right now but at this stage of the
2 questions we are really interested in
3 effectiveness.

4 DR. TRACY: I think that is really very
5 important because it doesn't matter how safely you
6 can do something that doesn't work.

7 DR. SAPIRSTEIN: That is right.

8 DR. TRACY: So, I think that a definition
9 or determination of patency is critical to this.

10 DR. HIRSHFELD: To put this in context, we
11 have a currently approved, marketed device about
12 which some data has surfaced that raises the
13 question that its performance may actually be
14 inferior to conventional techniques. The data to
15 date are inconclusive but there is enough of a
16 cloud that that question has arisen. So, I think
17 our first obligation is to answer the question as
18 to whether or not it is equivalent or not to
19 conventional hand-suturing techniques. That, I
20 think, is something that can be answered with a
21 highly focused trial. Because if it turns out that
22 it is inferior, then all the other questions, all
23 of a sudden, are moved way to the background. If
24 the trial can demonstrate that it is equivalent or
25 superior to conventional techniques, then the next

1 question arises which is, well, does that benefit
2 in performance translate into clinical outcomes.

3 DR. TRACY: Dr. Maisel?

4 DR. MAISEL: In my view, I don't think
5 that one size fits all. I think there is more than
6 one way that we could feel comfortable with the
7 efficacy, and certainly a randomized clinical trial
8 would meet that standard. But I certainly can
9 imagine an observational trial with 10,000 patients
10 and no events or, you know, 1,000 patients and 5
11 events. There certainly could be an observational
12 trial, a one-arm trial that I think would meet that
13 standard. So, you know, the challenge of drawing a
14 line in the sand is a big challenge but I don't
15 think it is impossible.

16 DR. TRACY: Dr. Yancy?

17 DR. YANCY: This is tangential and I am
18 happy to table this item at the Chair's
19 prerogative, but Dr. Hirshfeld's reference
20 resonates clearly with the concern that has been
21 building in my mind over the last several hours,
22 and that is that there is a device on the market
23 about which questions have been raised. Part of
24 the questions perhaps have to do with the umbrella
25 under which it received approval. Part of the

1 questions may have to do with the very issue we are
2 discussing right now in terms of trial design. It
3 is hard to know the denominator. We know there are
4 a number of events that have come to the FDA.

5 It is also incumbent upon us to support
6 the development of the minimally invasive surgical
7 procedures. So, I think the analogy of babies and
8 bath water is correct. But I am beginning to
9 wonder if it is anywhere within our purview to
10 address this one specific device, or maybe it has
11 to be set aside to another entity, but it seems as
12 if it is not unprecedented to make statements about
13 concerns regarding the safety of a device that is
14 already marketed, not necessarily to withdraw it
15 but indicate that we have seen a signal that raises
16 questions. I am not saying it is a bad device
17 because, again, the statistics of looking at this
18 are difficult. If this is inappropriate, then I
19 will withdraw this and beg apology, but it is a
20 question that I would like to pursue.

21 DR. TRACY: If we could let the FDA speak
22 to that.

23 MS. FLEISHER: Dina Fleisher, from the
24 FDA. This issue has been addressed by the FDA.
25 There was a whole team here from our Office of

1 Surveillance that has been taking care of this and
2 addressing it in the purview of the limits that we
3 have at FDA. For the interest of this panel
4 meeting, that really isn't probably relevant to
5 this particular discussion of what future designs
6 we would like to have.

7 DR. YANCY: One follow-up question?

8 DR. TRACY: Yes?

9 DR. YANCY: Are there activities you plan
10 to do, without compromising your deliberations,
11 that you can tell us? I mean, is there a time
12 sequence or is there a statement to be prepared?
13 Is there anything that you can say in that regard
14 or is that an inappropriate question?

15 MS. FLEISHER: We could probably do a
16 one-minute summary if that is helpful to you.

17 DR. YANCY: I am happy to do that off-line
18 since I may be the only person who has this
19 concern.

20 [Several members reply, "no, you're not."]
21 Then I think it would be important to us.

22 MS. FLEISHER: I will let the lead
23 reviewer from that office actually address that.

24 MS. HOANG: I am Quynh Hoang, from the
25 Office of Surveillance and Biometrics. As Dina

1 Fleisher has indicated, the FDA is aware of a
2 number of MDR reports regarding the Symmetry device
3 and we have convened a committee, a cross-center
4 committee, to look at the issue. Your question is
5 specific to what FDA plans to do about the adverse
6 events that we have seen with the Symmetry, and the
7 only thing I can say is that we are working with
8 the company.

9 DR. TRACY: Thank you. Dr. Slaughter?

10 DR. SLAUGHTER: Dr. Slaughter. Just two
11 issues, using a patient as their own control and
12 trying to identify two arteries has been done
13 before, and that was the ARIA trial using a
14 synthetic conduit. Although it seems appealing, it
15 is extremely difficult to administer. The issue
16 that comes up is that even preoperatively you can
17 identify two vessels that you think are similar.
18 The problem is you are doing that on a preoperative
19 angiogram for which the stenosis may be different
20 in each, and the question is were they under-filled
21 because the dye is not getting there or are they
22 really that small? So, you are really using a very
23 crude comparison of external diameter. Then, when
24 you get there you sort of flip a coin and pick one.
25 You open it; you can size it; you open the other

1 one and it turns out it is different.

2 The other issue is still that although
3 they may look the same size, you know, before you
4 start there are other issues as to the size of the
5 bed they supply and things like that. It becomes
6 very difficult technically to administer and it can
7 prolong the operation significantly when you are
8 trying to randomize them intraoperatively. So,
9 although it seems appealing, and certainly to use
10 the patient as their own control just in general,
11 and if you do complete angiography as opposed to
12 selective angiography it is more doable, but to
13 truly try and randomize two what you think are
14 matched vessels in an individual patient is
15 technically very difficult.

16 DR. BRIDGES: I don't think you need to
17 necessarily match the vessels in each patient. I
18 think that the randomization will take care of
19 that. I mean, you don't match patients that are
20 randomized. You don't say, well, I am going to
21 randomize these two patients and I think they are
22 equivalent. The randomization itself--it could
23 very well be that the first patient will have one
24 OM that has a tight stenosis and one that doesn't
25 and then in the next patient the converse will be

1 true. So, as long as you are randomizing the
2 technique that you use within the patients as you
3 go along, that should be a valid study design,
4 unless I am missing something.

5 DR. SLAUGHTER: Then you are back to where
6 you have a control that is not getting the device
7 and you have to have follow-up with invasive
8 monitoring and you are subjecting a control to an
9 invasive test.

10 DR. BRIDGES: No, what I am suggesting is
11 the same design that the group in Hanover used. In
12 that case, unless I misunderstood it, each patient
13 either had two or four anastomoses performed. If
14 they had two, then a random decision was made as to
15 which of those two would be done with one
16 technique, hand-sewn versus the device and each
17 patient was randomized in that fashion so that
18 there was no control group that was angiogrammed.
19 These are all patients who had an anastomotic
20 device used that served as their own control. So,
21 within those patients it is not relevant whether
22 the two arteries that were randomized were
23 equivalent. That is not the question.

24 DR. SLAUGHTER: I understand. I think
25 they are close but I think it does make a

1 difference because there is a difference in
2 distals. Unless you quantify the information when
3 you start to review it, it is very difficult to
4 assess patency.

5 DR. BRIDGES: Sure.

6 DR. EDMUNDS: Well, to show off the
7 culture at the University of Pennsylvania, I am
8 going to disagree with my colleague, Dr. Bridges,
9 and agree with Dr. Slaughter. There is plenty of
10 data to show that the degree of proximal stenosis
11 in a native artery influences patency. Moreover,
12 there is a difference in patency rate between
13 target arteries--circumflex, right, LAD--well, LAD
14 is not in this equation. And, it will be a
15 logistic nightmare, Charles, to try to get patients
16 in which you can randomize the arteries. You just
17 won't be able to get enough patients without
18 getting into a whole lot of institutions, I don't
19 think.

20 DR. BRIDGES: Sure. I am not trying to
21 say it would be easy or simple to interpret, but I
22 am simply making the point that I think what has
23 been suggested is following the design of the
24 Hanover study and, as I understand it, that is the
25 way they conducted the study and that is my point.

1 My point is to say that their study did not involve
2 trying to match arteries. They simply randomized
3 technique selection within patients who had an even
4 number of arteries. I certainly acknowledge that
5 that is not the same question as to whether the
6 degree of stenosis or the circumflex versus the
7 LAD, or whatever, will have different patency
8 rates. Hopefully, that would be taken care of by
9 enrolling enough patients but there still would be
10 challenges with that kind of a study. I am not
11 advocating it; I am simply trying to confirm what
12 the design of that study was.

13 DR. TRACY: I am going to try and move us
14 on to some of the next points, but just to briefly
15 summarize, the clear message is that not having
16 some type of control or some comparative basis is
17 not acceptable. Just having a study where the
18 device alone is being analyzed and then hoping that
19 something will be picked up at a later point is not
20 an adequate endpoint. There needs to be some type
21 of control built into the study but we are having a
22 hard time grappling with what that control should
23 be. Should it be intra-patient; should it be
24 patient-to-patient; or should be looking at a hard
25 outcome compared to historic control, patency

1 outcome compared to historic control plus
2 randomization on surrogate endpoints such as
3 neurologic events.

4 I think that we would need, and I don't
5 think we can do that sitting here, a statistical
6 analysis to figure out what size would be needed
7 for each of these different models--unless he is
8 very good and very fast. I would like to move on
9 to some of the other issues and maybe, when Dr.
10 Blumenstein is read, to come back to that if we
11 can.

12 DR. KRUCOFF: Cindy, just before you leave
13 that, let me make one comment about size. The
14 other side of what we have heard a lot of wisdom on
15 today is that you can drive the density of
16 endpoints into a higher direction to do a smaller
17 trial. So, if you know that actually arteries that
18 have greater plaque burden are more likely to fail,
19 or if you have diabetics they are more likely to
20 have problems, you can use propensity scores to
21 actually target a higher risk population and do a
22 smaller trial. So, you get into the conundrum of
23 you have to find those patients--

24 DR. TRACY: I think we have had some
25 disagreement that a propensity score would be

1 effective in comparing at least to historic
2 controls so--

3 DR. KRUCOFF: I am not talking about
4 historic controls. I am talking about how
5 frequently you can anticipate a failure endpoint in
6 a population to do a prospective study and do 100
7 patients and have an answer rather than 500 or
8 1,000 and have an answer.

9 DR. TRACY: Right. Well, we will leave it
10 there for the moment. Let us know when you have
11 finished your analysis, Dr. Blumenstein.

12 The second set of questions, with regard
13 to device placement and device design, please
14 address the following: Given the considerable
15 differences between the proximal and distal CABG
16 anastomoses, what, if any, differences in study
17 criteria should be required?

18 My observation would be you can't do a
19 proximal without doing a distal. I am not clear,
20 Dr. Sapirstein, on exactly what the question is
21 here.

22 DR. SAPIRSTEIN: There has been a lot of
23 discussion about the connector device which is
24 strictly for the proximal anastomosis of a venous
25 conduit to the aorta and doesn't involve anything

1 about suturing to the distal coronary. Now, there
2 are a lot of devices which specifically address the
3 distal anastomosis to the coronary artery and not
4 to the aorta, and other devices which are just
5 going to anastomose an artery, the LIMA to the LAD
6 or something of that nature. So, there is a lot of
7 variation in the indications for a participant
8 device.

9 DR. WHITE: Well, I would say that the
10 first criteria, the primary criteria has to be
11 patency, and patency whether it is proximal or
12 distal doesn't matter to me. I want to know that
13 the graft and the device is patent and I am not
14 confident that we have non-invasive ways to screen
15 these patients for patency. So, that brings us
16 back to a hard endpoint, patency. So, I am in
17 favor of applying the same criteria to the distal
18 as I would to the proximal anastomosis.

19 DR. TRACY: John?

20 DR. HIRSHFELD: I am saying this with some
21 fear that I will subsequently be exposed as being
22 simplistic, but I think that the issues are the
23 same and I think that the trial design, whether we
24 are examining proximal or distal anastomoses,
25 should be able to be the same, with the same

1 departments.

2 DR. WHITE: No, they are not the same.

3 DR. HIRSHFELD: Okay, I knew somebody
4 would bring that up.

5 DR. TRACY: Dr. Edmunds, do you want to
6 make some comments?

7 DR. EDMUNDS: I would suggest for starters
8 that for proximal saphenous vein, aortal saphenous
9 vein anastomosis the following for endpoints: An
10 aortic complication, changes in neurocognitive
11 state or stroke, patency whether it is an occlusion
12 or stenosis, hemorrhage, and acute revision, you do
13 it at the table or just after you closed, and any
14 device-related death as endpoints. I have left out
15 myocardial infarction.

16 DR. TRACY: I think those were the reasons
17 for wanting the MIDCAB in the first place, to
18 reduce some of those neurologic types of outcomes.
19 The one thing that is in there though that is
20 similar to distal would be patency. But I agree,
21 those are additional concerns for the proximal
22 anastomosis. Dr. Ferguson?

23 DR. FERGUSON: I think we can have the
24 same endpoints but the devices themselves are so
25 dissimilar and used for such different reasons that

1 we have to be cognizant of that.

2 DR. KRUCOFF: I agree with Tom. I think
3 it has been made very clear that even the technique
4 of putting something into a beating heart where you
5 are doing a distal anastomosis may be totally
6 different than the technique of connecting
7 something to an aorta. There may be some overlap
8 in the endpoints but these are separate trials.

9 DR. SAPIRSTEIN: A question that comes up,
10 Hank, is, for instance, do you require a six-month
11 patency evaluation for a proximal and a six-month
12 for a distal, or would you require different rigor
13 of evaluation for a proximal, which is a large
14 vessel, a large anastomosis, relatively crude?
15 That is one of the questions that sort of
16 distinguishes the two areas.

17 DR. EDMUNDS: Wolf, I am going to defer to
18 the collective wisdom in the room.

19 DR. AZIZ: From the data that we saw with
20 the other device, by six months we were seeing
21 problems so I think even though the time period may
22 be different I think six months should be at least
23 the minimum time frame.

24 DR. TRACY: Dr. Kato?

25 DR. KATO: I would agree with Dr. Aziz

1 that six months should determine patency both for
2 proximal and distal anastomoses. The only
3 difference with the distal is that the criteria of
4 any aortic complication or stroke won't be there.

5 DR. TRACY: Dr. Yancy?

6 DR. YANCY: Going back to the question of
7 being simple, I think simple is good, John, and the
8 study design is quite similar, in my judgment. It
9 is just the endpoints and surrogates and the
10 variables that we follow are different, some
11 referable to the aortic anastomosis, others
12 referable to the technical concerns. But
13 ultimately the study designs are not terribly
14 different. It is just the length of follow-up and
15 what we follow.

16 DR. TRACY: Yes?

17 PROF. KLIMA: Uwe Klima, from Hanover. I
18 think the myocardium does not care whether the
19 blood flows through a proximal anastomotic device
20 or through a distal device anastomotic device so
21 the endpoint really should be patency rate after
22 six months and does the anastomosis look good
23 whether it is running through a proximal or through
24 a distal anastomosis.

25 DR. TRACY: Dr. Weinberger?

1 DR. WEINBERGER: Yes, I am sitting here
2 and I am wondering whether people are beginning to
3 blur the difference between graft patency and graft
4 stenosis. I think that the numbers that people are
5 talking about, looking at six months, makes sense
6 in the context of stenosis. If you are talking
7 about patency, the numbers sound to me like looking
8 just at six months you are not going to have enough
9 events. Maybe we can get the surgeons to comment
10 on that. You care ultimately about how the graft
11 feels when it has done healing. The failure rates
12 that we are talking about really are complete
13 occlusions. If you look at six months, is the
14 acute response to the surgical process, given these
15 new connectors, completely finished by that point?

16 DR. AZIZ: You may not get complete
17 adhesion but you may even get severe stenosis so
18 even that is bad.

19 DR. WEINBERGER: I agree that it is bad,
20 but we are talking about endpoints being graft
21 failure. Is that right?

22 DR. EDMUNDS: Well, let's just stay with
23 the aortic saphenous vein proximal anastomosis and
24 then deal separately with distal. Okay? The
25 reason is that I think that you can't say when

1 healing ends, except when you change your address
2 permanently. So, I think you just pick a period in
3 time and you say that is our endpoint. If we see
4 25 percent stenosis in 50 percent of the device
5 applications and we see 3 percent stenosis in the
6 control group, that probably would work out, if we
7 got enough power, to significance. Now, how
8 important that is, that is not the question. That
9 is not the outcome either. It is the stenosis.
10 That is what I am advocating.

11 DR. WEINBERGER: So, just to concretize
12 this, we are talking about a continuous variable at
13 six months and we are talking about not measuring
14 graft patency rate, although that might be a
15 secondary endpoint. We are talking about the
16 difference in the distribution of patients with
17 respect how tight the stenosis is proximally. Is
18 that right?

19 DR. EDMUNDS: The question you are asking
20 is the same as what is the patency rate of the
21 rapamycin stent at five years? No one knows. We
22 haven't been five years. It could be 50 percent.

23 DR. TRACY: So, I guess the answer here is
24 that you find a time, a period of time where you
25 look and that becomes the data point that you have.

1 DR. KRUCOFF: Well, there is an issue with
2 a biomarker. In fact, what you are describing is
3 minimal luminal diameter so as a continuous
4 variable you take your point, six months for the
5 proximal anastomosis in the device group. The
6 minimal luminal diameter is smaller than for the
7 non-device group. That may be a statistically
8 significant difference. Then a separate question
9 becomes is it so small that they are having angina,
10 infarctions or deaths? And, that is a trial design
11 that is, again, pretty well established but it is a
12 biomarker where the statistical significant
13 difference is not one and the same as the
14 clinically meaningful endpoint.

15 DR. TRACY: Let me skip ahead to part b),
16 are there certain aspects of the clinical study
17 design, example, length of follow-up and endpoints,
18 that should be required for all devices
19 irrespective of device form and function? for
20 example, the U-clip performance closely duplicates
21 that of a suture, whereas the Symmetry has greater
22 similarity to a stent.

23 So, where we are I believe is that there
24 would be specific aspects of a clinical trial
25 design that would be different whether you are

1 talking about proximal or distal, but there has to
2 be some concrete endpoint, and probably you could
3 come up with a generic time frame--six months, nine
4 months, whatever you choose--whereby the patency or
5 the percent stenosis would be analyzed. The other
6 endpoints might be more appropriately analyzed, of
7 course, at different times--neurologic events
8 acutely or catastrophic events, obviously, would be
9 analyzed acutely.

10 DR. BRIDGES: I think to answer this
11 question, 2.b, my opinion would be that these
12 devices are so different and so protean in their
13 designs that I think it would be difficult for us
14 to prospectively or a priori try to figure out what
15 we expect the differences to be. So, I would think
16 that the design of a study for all distal devices
17 ought to be the same, at least as a starting point.
18 I mean, you might look at a device and say, well,
19 that looks like a suture and, therefore, we don't
20 need maybe as large a study group. But I think the
21 most straightforward approach would be to have the
22 same endpoint irrespective of what our opinion
23 about how we think the device will function is.

24 DR. KATO: But the implication really is,
25 is one device similar to another one, which is a

1 discussion, as Dr. Frater who I guess is not here
2 anymore brought up before about, you know, is this
3 device definitely going to be a PMA or is it going
4 to be a 510(k)?

5 DR. TRACY: It is not in the purview of
6 this group to change whether something is a 510(k)
7 or PMA. That is completely off the table even as a
8 question to this group. The question really is
9 study design. I think that you have to design a
10 study that is going to capture important endpoints
11 with it is a 510(k), PMA, or whatever you call it.
12 A concrete point in time to say let's look and see
13 what the anatomy looks like, with it is a
14 continuous variable like stenosis or absolute
15 patency/occlusion, that is a time you choose that
16 makes biologically some sense--six months, nine
17 months, and the other variables appropriate to
18 those types of outcomes. Does that seem to make
19 sense?

20 DR. YANCY: So, in that context I would
21 like to support what Dr. Bridges said, that there
22 is sufficient variability and sufficient degree of
23 unknown that the answer to the proposed question is
24 that there should be some consistency in the time
25 of follow-up--

1 DR. TRACY: Right.

2 DR. YANCY: --irrespective of the
3 mechanism of action.

4 DR. TRACY: Right. I think that is a good
5 point. Part c) to that question, it is rarely
6 possible to determine the cause of conduit failure.
7 Can you suggest criteria to determine whether a
8 failure is device related? Dr. Yancy?

9 DR. YANCY: That is the essence of a
10 clinical trial design. You have to have a clinical
11 trial that rules out the confounders, or at least
12 equilibrates them between the intervention and
13 reference group so that you can have only one
14 remaining variable and presume that it is, in fact,
15 the device. I mean, that is the whole purpose for
16 a clinical trial.

17 DR. ZUCKERMAN: That is perhaps the first
18 part of the equation, but FDA was also wondering
19 whether the panel agreed with some of Dr. Krucoff's
20 earlier comments. For example, should some of the
21 methodologies from the stent trials be adapted to
22 these CABG trials, meaning should there be an
23 independent clinical events committee looking at
24 these events to try to determine cause of graft
25 failure. Two, is it worthwhile to, for example,

1 have an independent core angiographic laboratory
2 look at the angiograms, at least in a percentage of
3 patients, to confirm what the sites are calling
4 device versus non-device related?

5 DR. EDMUNDS: Bram, do you have the
6 criteria for the stent trial? I am not sure that I
7 know what those were.

8 DR. ZUCKERMAN: A set of criteria has been
9 developed over the last ten years in order to
10 minimize potential sources of bias and to make sure
11 that we are comparing apples with apples, and they
12 include randomization, use of an independent
13 clinical events committee and, because of the
14 importance of correct angiographic analysis and the
15 known ability of site investigators to look at
16 angiograms and develop measurements that can't be
17 confirmed, the use of independent core labs. Those
18 are the main ones.

19 DR. MAISEL: I think an independent core
20 lab is obviously going to be a critical part of any
21 assessment of graft stenosis. I do think that an
22 independent data and safety monitoring board or
23 clinical events committee is also a critical
24 component. I think getting at the issue of cause
25 of device failure or cause of failure of graft

1 patency, whether it is due to the device or not, is
2 going to be extremely difficult based on
3 angiographic or clinical criteria, and I think, you
4 know, the trial designs we talked about are really
5 the only way to get at that.

6 DR. AZIZ: I think the cause of graft
7 patency or failure sometimes may be clear-cut. For
8 example, if you are doing an angiogram and you see
9 90 percent stenosis at the proximal end, that may
10 be okay but if you do the angiogram at six months
11 and the graft was stenosing at three months and the
12 flow was slowing by six months, it may be occluded.
13 So, you know, I don't think it will always be
14 clear-cut but sometimes I think it will be but
15 sometimes it may not be. All you would say is that
16 the graft is occluded and you may not be able to
17 say if it was, let's say, a rheology type problem
18 or whether it was because of stenosis. So,
19 hopefully, by doing a randomized study that might,
20 you know, become clearer.

21 DR. KATO: But one question for you, from
22 what I hear from the angiographers, when you get
23 stenosis or an occlusion from one of these devices,
24 the clot or the platelets go right up to the edge
25 of the aortal-vein junction. Correct? Because

1 that is relatively atypical for a hand-sewn graft
2 with a good cobra head which is usually open and
3 then there is clot distal to that, which suggests
4 that it is not a proximal anastomosis problem.

5 DR. WHITE: I think that when you close a
6 graft--I was interested in some of the images that
7 were shown today, and I think John was too, about
8 some of the stenoses in the middle of the graft.
9 What I think you are saying is that you don't know
10 if it is intimal hyperplasia from within the stent
11 part of the device or whether it is the middle part
12 of the graft that actually goes down because as
13 soon as the graft goes, then the thrombus
14 propagates. In the hand-sewn stuff you have to
15 make that hood and we always have a little nipple.
16 These devices seem to go right up flush. But I
17 don't know if you can tell where the problem
18 started angiographically after the occlusion has
19 occurred. It could be with the device or it could
20 be in the middle of the graft.

21 DR. KATO: But don't you think if the clot
22 or the platelet plug goes right up to the end of
23 the graft, you know, to occlude the device that
24 most likely that is where it started? In the
25 standard hand-sewn anastomosis--you are right, you

1 still have that little nipple there.

2 DR. TRACY: Are we struggling with the
3 fact, exactly what the FDA is saying, that it is
4 difficult to determine the cause of conduit
5 failure? I am not sure that anything we have
6 discussed so far has really gotten us any closer,
7 except to say that having an independent DSMB and
8 core lab to look at it. I am not sure what they
9 are looking at though. I mean, we don't understand
10 the pathology, as far as I am hearing from this
11 discussion. Unless, Dr. Emery, you have something
12 that would clarify this.

13 DR. EMERY: I don't think so. I agree
14 that when you get either a proximal or distal
15 stenosis, thrombosis propagates to the next major
16 branching which in vein grafts is either the
17 proximal or distal anastomosis.

18 What I was going to address here was the
19 issue of distal anastomoses because you need to
20 remember that the conduits are different. Of all
21 the data we have seen today that is variable, the
22 data that is least and almost invariable is the
23 early and one-year patency of internal mammary
24 artery grafts. As Michael Mack has said, we need
25 to think like cardiologists where we solve problems

1 rather than fight about them. The problem of the
2 LIMA to the LAD is fairly straightforward.

3 Historical data is new. It is done on
4 MIDCAB and done on off-pump surgeries. There are
5 multiple papers and very recent literature and it
6 is very solid data, and you are evaluating an
7 artery to an artery anastomosis with a
8 stent-like--I don't want to say a stent, a
9 stent-like process there so interventional
10 angiography could be applied to a mammary to the
11 LAD. In saphenous vein grafts that is not true
12 because the conduit is variable for various
13 individuals. So the criteria for establishing
14 patency leans more towards randomization and, as
15 Dr. Edmunds has addressed, there are different
16 patencies for different systems, different run-offs
17 and different degrees of stenosis with the vein
18 graft in particular. That does not appear to be
19 true with mammary artery grafts. So, you may have
20 to vary your study criteria for the conduit that is
21 being utilized in a particular study.

22 DR. TRACY: John?

23 DR. HIRSHFELD: Thinking ahead to what the
24 data will be, the first endpoint is going to be
25 patent versus not patent and that is going to tell

1 us whether there is a performance difference
2 between the device and the traditional technique.
3 Patent versus not patent is not going to say
4 anything about the etiology of the loss of patency.
5 There likely will be a group of grafts that are
6 patent but stenotic and in those grafts the
7 location of the stenosis will probably give us some
8 insight as to what the etiology of the impending
9 failure of that graft would be, and we would
10 probably be able to tease out the answer to the
11 etiology of increased overall graft failure from
12 that data.

13 DR. TRACY: So, it really becomes
14 observational, once again emphasizing the need to
15 have angiographic follow-up so we can over time
16 figure out what these failures are related to. Dr.
17 Yancy?

18 DR. YANCY: I think there is just one
19 other variable referable to the question that Dr.
20 Zuckerman has on the table that captures everything
21 we have been dealing with. The one concept we
22 haven't addressed is that of bias. There is no way
23 that this can be a blinded protocol because the
24 surgeon knows what he or she has done. And one of
25 the strengths of having adjudication committees,

1 DSMBs and core labs is that it really mitigates
2 that whole variable. So, I think the answer to
3 what you have raised, Bram, is an emphatic yes, we
4 do need to populate that study with the usual
5 complement of oversight committees, maybe even more
6 so because of concerns we might have.

7 DR. BRIDGES: I agree with that and the
8 point that Dr. Hirshfeld made. I think it might be
9 important to make sure as we go through this
10 evaluation process that some attention is paid not
11 just to patency versus non-patency, stenosis versus
12 non-stenosis, but that a deliberate attempt be made
13 to indicate the distribution of stenoses within
14 grafts and that might provide additional useful
15 information. I mean, typically that information is
16 not available when we read studies of this nature
17 but, you know, if you found a different
18 distribution of stenoses that might provide the
19 clues as to whether it was device related.

20 DR. WHITE: I would like to say that I
21 think we can fulfill two out of three of Bram's
22 criteria, that is the core objective criteria, the
23 core lab and DSMB. I would still like to leave the
24 door open for a trial that is a single-arm
25 objective performance criteria trial. That would

1 still lend itself to DSMB oversight, high rate of
2 angiographic follow-up and independent assessment
3 of the angiographic endpoints, and worry about the
4 slippery slope issues that are real that will come
5 up if you come close but don't make the endpoint.

6 DR. SAPIRSTEIN: We have used the
7 generally applied criteria for conduit patency as
8 greater than 50 percent, not more than 50 percent
9 stenosis. If there is stenosis greater than 50
10 percent, we have considered that an obstructive
11 lesion, a failure. Is this an acceptable endpoint?

12 DR. FERGUSON: I am no angiographer by a
13 long shot, but I think we know from the rheology
14 that 50 percent is a flow-limiting lesion--

15 DR. SAPIRSTEIN: Yes.

16 DR. FERGUSON: --so I think that is a good
17 place to start.

18 DR. WHITE: I would like to start at the
19 place where intervention is contemplated, and I
20 think that failure would be defined as the level at
21 which you would be willing to reintervene either
22 with reoperation, angioplasty or stent. I think
23 that either 50 or 70 percent is going to be that
24 threshold.

25 DR. KRUCOFF: I think as a true continuous