

FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
OPHTHALMIC DEVICES PANEL

Ninety-fourth Meeting

*This transcript has not
been edited and FDA
makes no representation
regarding its accuracy*

Tuesday,
January 12, 1999

2146 '99 JAN 21 P2:53

Main Conference Room
Office of Device Evaluation
9200 Corporate Boulevard
Rockville, Maryland

IN ATTENDANCE:

Voting Members

JAMES P. McCULLEY, M.D., Chair
MARK A. BULLIMORE, MCOptom, Ph.D.
EVE J. HIGGINBOTHAM, M.D.
JANICE M. JURKUS, O.D.
MARIAN S. MACSAI, M.D.
JOSE S. PULIDO, M.D.
JOEL SUGAR, M.D.

Consultants, Deputized to Vote

KAREN BANDEEN-ROCHE, Ph.D.
MICHAEL R. GRIMMETT, M.D.
ALICE Y. MATOBA, M.D.
WOODFORD S. VAN METER, M.D.
MING X. WANG, M.D., Ph.D.

Non-Voting Members

RENEE A. MIDDLETON, Ph.D., Consumer Representative
MARCIA S. YAROSS, Ph.D., Industry Representative

Panel Executive Secretary

SARA M. THORNTON

C O N T E N T S

PAGE

Call to Order

James P. McCulley, M.D., Chair

5

Introductory Remarks

Sara M. Thornton, Executive Secretary

5

OPEN COMMITTEE DISCUSSION

Conflict of Interest Statement and
Appointment to Temporary Voting Status

Sara M. Thornton

9

Branch Updates

James F. Saviola, O.D., Chief
Vitreoretinal and Extraocular Devices Branch

11

Morris Waxler, Ph.D., Chief
Diagnostic and Surgical Devices Branch

14

Donna R. Lochner, Chief
Intraocular and Corneal Implants Branch

15

Introduction of PMA P980031

Donna R. Lochner

18

Ashley A. Boulware
Team Leader and Biomedical Engineer
Intraocular and Corneal Implants Branch

18

C O N T E N T S

PAGE

Presentation by the Sponsor,
KeraVision, Inc., of PMA P980031,
Intacts Intrastromal Corneal Ring Segments

Thomas M. Loarie	20
Darlene E. Crockett-Billig	21
David J. Schanzlin, M.D.	26
Michael A. Lemp, M.D.	37

Clinical Review

Malvina B. Eydelman, M.D. Division of Ophthalmic Devices	58
---	----

Primary Panel Reviewers

Joel Sugar, M.D.	73
Woodford S. Van Meter, M.D.	80
Michael R. Grimmett, M.D.	87

Questions to the Sponsor	95
--------------------------	----

Panel Discussion	153
------------------	-----

Closing Comments by the Sponsor	223
---------------------------------	-----

Vote on Recommendation	232
------------------------	-----

P R O C E E D I N G S (8:09 a.m.)

DR. McCULLEY: I'd like to call to order the 94th meeting of the Ophthalmic Devices Panel.

I would like to turn the floor now to Sara Thornton for introductory remarks.

MS. THORNTON: Good morning. I'd like to welcome you all to the meeting. I'm Sara Thornton, the executive secretary for the Ophthalmic Devices Panel. Before we proceed with today's agenda, there are a few announcements I'd like to make.

I'd like to remind everyone that you're requested to sign in on the attendance sheets in the registration area just outside this meeting room. You may pick up an agenda and information about today's meeting and how to obtain summary minutes or panel transcripts there as well.

You should make a note that the panel meeting tentatively scheduled for March 11th and 12th of this year has been canceled. The status of the May 3 and 4 meeting will be available early in March. Check our Web site or you can call me. It should be on the hotline, though, and we'll let you know as soon as we know.

Messages for the panel members and FDA participants, information and special needs should be directed through Ms. Ann-Marie Williams or Ms. Shirley

1 Meeks who are available at the registration table or in the
2 registration area. The phone number for contacting you
3 here is (301) 443-8011. For those of you with cell phones
4 and pagers, we ask that you turn them off or put them on
5 vibration mode while you are in this room.

6 Your lunch sheets, which I hope all of you have
7 had an opportunity to fill out and submit the required
8 monies to Ms. Williams, I hope you've done that now, but if
9 you haven't, would you finish that business up at the
10 break, please, so we can get those orders over and you'll
11 have your lunch when you have the time to eat it.

12 Lastly, will all meeting participants please
13 speak into the microphone and give your name clearly so
14 that the transcriber will have an accurate recording of
15 your comments.

16 At t his time, I'd like to introduce our panel
17 to you. We are missing at the moment a couple of folks and
18 I will introduce them when they come in. But I would like
19 to begin the introductions. I think it would be good to
20 start with Dr. McCulley here to my right.

21 DR. McCULLEY: Jim McCulley, professor and
22 chairman, Department of Ophthalmology, University of Texas
23 Southwestern Medical School in Dallas.

24 DR. HIGGINBOTHAM: Eve Higginbotham, professor
25 and chair, Department of Ophthalmology, University of

1 Maryland, Baltimore.

2 DR. SUGAR: Joel Sugar, professor of
3 ophthalmology, University of Illinois at Chicago.

4 DR. BULLIMORE: Mark Bullimore, the Ohio State
5 University College of Optometry.

6 DR. JURKUS: Jan Jurkus, professor, Illinois
7 College of Optometry in Chicago.

8 DR. ROSENTHAL: Ralph Rosenthal, director,
9 Division of Ophthalmic Devices.

10 DR. YAROSS: Marcia Yaross, director of
11 regulatory affairs, Allergan, and industry representative
12 to the panel.

13 DR. BANDEEN-ROCHE: Karen Bandeen-Roche,
14 associate professor of biostatistics, Johns Hopkins
15 University, Baltimore.

16 DR. GRIMMETT: Michael Grimmett, assistant
17 professor of clinical ophthalmology, University of Miami
18 School of Medicine, Bascom Palmer Eye Institute.

19 DR. WANG: Ming Wang, assistant professor of
20 ophthalmology, Vanderbilt University.

21 DR. VAN METER: Woodford Van Meter, private
22 practice, cornea and external disease, Lexington, Kentucky.

23 DR. MATOBA: Alice Matoba, associate professor
24 of ophthalmology, Baylor College of Medicine, Houston,
25 Texas.

1 DR. MACSAI: Marian Macsai, professor of
2 ophthalmology, West Virginia University School of Medicine.

3 MS. THORNTON: Thank you, and welcome to you
4 all.

5 I would like to point out a typo in the first
6 page of the agenda. My sincerest regrets to Dr. McCulley.
7 He is the Chair. He is not the Interim Chair.

8 DR. MCCULLEY: As you pointed out, I'm still in
9 training.

10 MS. THORNTON: Yes.

11 DR. MCCULLEY: As I pointed out, I'll be in
12 training the day I leave.

13 MS. THORNTON: Now I think I can turn it back
14 over to Dr. McCulley, and we will proceed with the first
15 portion of our meeting.

16 DR. MCCULLEY: I'd like to open the public
17 hearing session of this meeting. We had only one person
18 who indicated he wished to speak before who has canceled.
19 If there is anyone in the audience who would like to come
20 to the podium to make comments, you're welcome to do so.

21 (No response.)

22 DR. MCCULLEY: Seeing none, that closes the
23 open public hearing.

24 There will be an opportunity for public
25 comments toward the end of the deliberations of the PMA and

1 prior to the panel vote. That comment period will be for
2 comments relative specifically to the PMA and deliberations
3 that have taken place.

4 We'll now begin the open committee discussion.
5 I'll turn the floor back to Ms. Thornton.

6 MS. THORNTON: First, I'd like to read the
7 conflict of interest statement for the Ophthalmic Devices
8 Panel meeting of January 12, 1999. The following
9 announcement addresses conflict of interest issues
10 associated with this meeting and is made part of the record
11 to preclude even the appearance of an impropriety.

12 "To determine if any conflict existed, the
13 Agency reviewed the submitted agenda and all financial
14 interests reported by the committee participants. The
15 conflict of interest statutes prohibit special government
16 employees from participating in matters that could affect
17 their or their employer's financial interest. However, the
18 Agency has determined that participation of certain members
19 and consultants, the need for whose services outweigh the
20 potential conflict of interest involved, is in the best
21 interest of the government.

22 "We would like to note for the record that the
23 Agency took into consideration certain matters regarding
24 Drs. Janice Jurkus, James McCulley, and Ming Wang. These
25 panelists reported past and/or current involvement in firms

1 at issue but in matters not related to today's agenda.
2 Since their interests are unrelated, the Agency has
3 determined that they may participate in the committee's
4 deliberations.

5 "In the event that the discussions involve any
6 other products or firms not already on the agenda for which
7 the FDA participant has a financial interest, the
8 participant should excuse himself or herself from such
9 involvement and the exclusion will be noted for the record.
10 With respect to all other participants, we ask in the
11 interest of fairness that all persons making statements or
12 presentations disclose any current or previous financial
13 involvement with any firm whose products they may wish to
14 comment upon."

15 I would like now to read the appointment to
16 temporary voting status. "Pursuant to the authority
17 granted under the Medical Devices Advisory Committee
18 Charter dated October 27, 1990, as amended April 20, 1995,
19 I appoint the following individuals as voting members of
20 the Ophthalmic Devices Panel for the duration of this
21 meeting on January 12, 1999: Dr. Karen Bandeen-Roche, Dr.
22 Michael R. Grimmatt, Dr. Woodford S. Van Meter, Dr. Alice
23 Y. Matoba, Dr. Ming X. Wang.

24 "For the record, these persons are special
25 government employees and are consultants to this panel or

1 consultants or voting members of another panel under the
2 Medical Devices Advisory Committee. They have undergone
3 the customary conflict of interest review and have reviewed
4 the material to be considered at this meeting." Signed,
5 Dr. D. Bruce Burlington, Director of the Center for Devices
6 and Radiological Health, dated December 15, 1998.

7 Thank you.

8 DR. McCULLEY: Okay. We will shortly begin the
9 division updates. But prior to that in just a moment, Ms.
10 Thornton has another announcement to make. When we begin
11 that, we will begin with Dr. Rosenthal. But prior to that,
12 Sally?

13 MS. THORNTON: Yes, I would just like to
14 announce for the record that Dr. Pulido and our consumer
15 representative, Dr. Renee Middleton, are not here and will
16 not be here for this meeting, and their absences were not
17 anticipated.

18 Thank you.

19 DR. McCULLEY: Dr. Rosenthal?

20 DR. ROSENTHAL: I have no comments, Mr.
21 Chairman.

22 DR. McCULLEY: Branch updates. Dr. Saviola?

23 DR. SAVIOLA: Good morning. I have one update
24 to report to the panel. Our branch has completed a review
25 of the Wesley-Jessen Precision UV (Vasurfilcon A) soft

1 contact lens that requested a modification to the labeling
2 to include an additional indication statement, and a
3 modified UV lens labeling note and warning statement.

4 Since the lens is marketed for both daily wear
5 and extended wear, a 510k, K982988, and a PMA supplement,
6 P940013/S006, were both reviewed and cleared for marketing
7 as of last week.

8 The additional indication statement reads as
9 follows. It's a single sentence. "Precision UV lenses
10 help protect against transmission of harmful UV radiation
11 to the cornea and into the eye." The remainder of the
12 indication statement was unchanged.

13 The revised warning is: "UV-absorbing contact
14 lenses are not substitutes for protective UV-absorbing
15 eyewear such as UV-absorbing goggles or sunglasses because
16 they do not completely cover the eye and surrounding area.
17 You should continue to use UV-absorbing eyewear as
18 directed."

19 The new phrase in the warning, "because they do
20 not completely cover the eye and surrounding area," was
21 added to explain the limitations of the lens.

22 The new note and the old note are on the
23 overhead and I will read them for the record. The new note
24 reads: "Long-term exposure to UV radiation is one of the
25 risk factors associated with cataracts. Exposure is based

1 on a number of factors such as environmental conditions --
2 altitude, geography, cloud cover -- and personal factors --
3 extent and nature of outdoor activities. UV-absorbing
4 contact lenses help provide protection against harmful UV
5 radiation. However, clinical studies have not been done to
6 demonstrate that wearing UV-absorbing contact lenses
7 reduces the risk of developing cataracts or other eye
8 disorders. Consult your eyecare practitioner for more
9 information."

10 Previously, the note had stated: "The
11 effectiveness of wearing UV-absorbing contact lenses in
12 preventing or reducing the incidence of ocular disorders
13 associated with exposure to UV-light has not been
14 established at this time."

15 During the review of these documents, the
16 recommendation of the Promotion and Advertising Staff in
17 the Office of Compliance were sought and considered by our
18 branch. Also, since the note, warning, and indications
19 appear in patient labeling, a review of the proposed
20 changes was conducted by human factors reviewers in the
21 Office of Surveillance and Biometrics, Division of Device
22 User Programs.

23 This was a literature-based application. There
24 have been a number of studies published since the original
25 policy letter was issued in 1987. However, since there are

1 no clinical data to support the claim that a UV-absorbing
2 lens will reduce the incidence of cataracts, that statement
3 is not included in the revised labeling.

4 We are currently working on a policy to advise
5 other firms how they may proceed to change their product
6 labeling for other UV-absorbing contact lenses and we will
7 be issuing guidance to the industry in the near future.
8 Our goal is to take the least burdensome regulatory path.
9 However, any modification to a currently approved
10 indication statement would require submission of a new 510k
11 and/or PMA supplement, depending on the lens' indication of
12 daily or extended wear.

13 That concludes my update.

14 DR. McCULLEY: Thank you. I should have
15 pointed out that Dr. Saviola is chief of the Vitreoretinal
16 and Extraocular Devices Branch.

17 Now I would like to recognize Dr. Morris
18 Waxler, chief, Diagnostic and Surgical Devices Branch.

19 DR. WAXLER: Good morning. I want to thank the
20 staff of the Diagnostic and Surgical Devices Branch, the
21 medical officers, and others in the Division of Ophthalmic
22 Devices without whose diligent work device applications
23 would not be approved. Thanks also to panel members who
24 have been helpful on many occasions.

25 On December 17, 1998 the Food and Drug

1 Administration approved PMA P970053 for the Nidek EC-5000
2 Excimer Laser System for PRK treatment for the reduction or
3 elimination of mild to moderate myopia, -0.75 to -13.00
4 diopter spherical equivalent, less than or equal to 0.75
5 diopter astigmatism.

6 The FDA has approved a number of sponsor-
7 investigator clinical trials to investigate a wide range of
8 important variables in the PRK and LASIK treatment of
9 myopia and hyperopia, with and without astigmatism. These
10 studies are, or have been, conducted using refractive
11 lasers from many of the manufacturers. A large amount of
12 information has been learned about PRK and LASIK. Many of
13 the recent applications from sponsor-investigators are
14 redundant and therefore have been disapproved. The Agency
15 will approve additional sponsor-investigator clinical
16 trials for excimer lasers only if the sponsor has submitted
17 a study that has a unique and scientifically sound
18 investigational plan.

19 Thanks.

20 DR. McCULLEY: Does that conclude your report?

21 DR. WAXLER: Yes.

22 DR. McCULLEY: Thank you.

23 Next, Donna Lochner, chief, Intraocular and
24 Corneal Implants Branch.

25 MS. LOCHNER: Thank you. First, I'm pleased to

1 announce that Staar Surgical Company's PMA
2 P880091/Supplement 14 for their Toric Intraocular Lens was
3 approved by the FDA on November 4th of 1998. This PMA was
4 reviewed by the panel in July of 1998.

5 Second, I'd like to make the panel aware that
6 the Division plans to release in the very near future a
7 guidance document entitled "Accountability Analysis for
8 Clinical Studies for Ophthalmic Devices." This document
9 will be available shortly on FDA's Web site and will be
10 mailed to the panel when it is released. The guidance is
11 intended to provide general information about the analysis
12 of accountability in ophthalmic device investigational and
13 marketing applications and notifications. The guidance
14 defines common terms used in reporting accountability,
15 describes in a fair amount of detail the key factors to
16 consider in presenting loss to follow-up analyses, and
17 provides suggested formats for reporting accountability.

18 The Division hopes terminology and methods of
19 presentation can be somewhat standardized so that the
20 Division, sponsors, and the panel can more effectively
21 analyze these data and so that a common understanding of
22 accountability may result. We will appreciate any comments
23 you may have when the document is released.

24 Third, the Intraocular and Corneal Implants
25 Branch also plans to release very shortly the next draft of

1 the intraocular lens guidance document. Of particular note
2 in the release of this draft version is that it contains an
3 update of the historical control, or FDA grid, clinical
4 data. Again, we plan to mail a copy to the panel when it
5 becomes available, and we will welcome any comments you may
6 have on any part of the guidance.

7 And fourth and last, I'd like the panel and
8 public to be aware that FDA has issued a requirement that
9 all medical devices or packaging containing natural rubber
10 must include a statement in the product labeling stating
11 the presence of natural rubber. FDA is requiring this
12 statement because medical devices with natural rubber may
13 pose a significant health risk to some consumers who are
14 sensitized to natural latex proteins. The Branch expects
15 to ensure that this requirement is being met as we review
16 new labeling, especially for viscoelastic products that use
17 latex rubber stoppers in their packaging.

18 Thank you.

19 DR. McCULLEY: Thank you.

20 We will now move on to the introduction of PMA
21 980031. Donna Lochner, Ashley Boulware -- oh, pronounced
22 "Bouler"? It's not spelled that way.

23 (Laughter.)

24 DR. McCULLEY: Anyway, sorry, Ashley. I'll
25 learn. Remember, I'm still in training.

1 Would you please introduce the PMA?

2 MS. LOCHNER: Yes, I'd like to just acknowledge
3 the hard work of the review team on this PMA, in
4 particular, Ashley Boulware, the team leader for the PMA,
5 and I'd like to read into the record the review team who
6 reviewed this document. First, Ashley Boulware was the
7 team leader and the engineering reviewer; Dr. Malvina
8 Eydelman, the clinical reviewer; Dr. Gene Pennello, the
9 statistician; Susanna Jones, toxicology; Susan Gouge,
10 microbiology; Jean Toth Allen, bioresearch monitoring;
11 Carol Clayton, labeling; Ronald Swann, from the Office of
12 Compliance, performed the good manufacturing practices
13 review; and Mervin Parker from our PMA policy staff.

14 That concludes my comments. I will turn the
15 floor over to Ashley Boulware.

16 MS. BOULWARE: Mr. Chairman, members of the
17 panel, ladies and gentlemen, PMA P980031 has been submitted
18 for the KeraVision Intacts Intrastromal Corneal Ring
19 Segments. The Intacts consist of two half circles referred
20 to as "ring segments." The Intacts are machined from
21 Perspect CQ polymethylmethacrylate using techniques similar
22 to those employed in intraocular lens manufacturing. The
23 Intacts are available in three thicknesses -- 0.25
24 millimeter, 0.30 millimeter, and 0.35 millimeter. The
25 product has a fixed outer diameter of 8.10 millimeter and a

1 width of 0.8 millimeter. Each segment has an arc length of
2 150 degrees and a cross-section that is hexagonally shaped.

3 A single positioning hole is drilled into the
4 superior end of each segment to aid in the surgical
5 manipulation. The Intacts are inserted between the layers
6 of the corneal stroma through a small incision made in the
7 periphery of the cornea. When surgically placed at
8 approximately two-thirds depth into the corneal stroma, the
9 Intacts reshape the corneal curvature by increasing the
10 thickness of the cornea in the periphery. This peripheral
11 thickening causes the interior curvature of the cornea to
12 flatten, thereby correcting for myopia by lowering the
13 optical power of the eye. The degree of corneal flattening
14 achieved using the device is directly related to the
15 thickness of the Intacts implanted.

16 The proposed indication is for the reduction or
17 elimination of myopia of minus 1.00 diopter to minus 3.00
18 diopters at the spectacle plane in patients who are 21
19 years of age with documented stability of refraction as
20 demonstrated by a change of less than or equal to 1.00
21 diopters for at least six months prior to the preoperative
22 examination, and with preoperative myopia error ranging
23 from minus 1.00 diopters to minus 3.50 diopters with 1.00
24 diopter or less of astigmatism.

25 The Intacts are implanted using a set of

1 surgical instruments designed by KeraVision specifically
2 for this procedure. Hand-held surgical instruments are
3 generally regulated as Class I devices exempt from 510k
4 submission. FDA is currently considering the appropriate
5 classifications for the KeraVision instruments.

6 The primary panel reviewers for P980031 are
7 Drs. Grimmett, Sugar, and Van Meter. The sponsor has been
8 advised of the questions and concerns raised by the primary
9 panel reviewers and FDA's clinician, Dr. Malvina Eydelman.
10 Following the sponsor's presentation, Dr. Eydelman will
11 summarize issues from her clinical review.

12 Thank you, Mr. Chairman.

13 DR. McCULLEY: Thank you.

14 We'll move on to the sponsor presentation. A
15 point, one, is that you have one hour for your
16 presentation, and then just a point of clarification for
17 me. I'd been advised that you had a movie of the surgical
18 procedure that you were going to show and we were going to
19 allow time for that. You're not going to show that? Okay.

20 I will now turn the floor to the sponsor.

21 MR. LOARIE: Good morning. I am Tom Loarie,
22 chairman and chief executive officer of KeraVision.
23 KeraVision was founded in 1986 to develop new solutions to
24 correct common vision problems. We appreciate the
25 opportunity to appear before you today with clinical

1 results for our first product, KeraVision Intacts for the
2 correction of myopia.

3 KeraVision Intacts, shown here, are designed to
4 reshape the interior surface of the cornea for the
5 correction of myopia without compromising the central
6 optical zone.

7 The outline for today's presentation is as
8 follows. Darlene Crockett-Billig, KeraVision's vice
9 president of regulatory affairs and clinical research, will
10 present a regulatory overview. Dr. David Schanzlin,
11 professor of ophthalmology, University of California at San
12 Diego and the chief clinical investigator for KeraVision
13 trials, will provide an overview of Intacts technology and
14 an efficacy assessment. The safety assessment will be
15 delivered by Dr. Michael Lemp, who is chairman of
16 KeraVision's data and safety monitoring board. Dr. Lemp
17 will also present an assessment of reversibility and
18 adjustability before concluding today's presentation with a
19 risk-benefit assessment of the Intacts product.

20 Now I would like to introduce Darlene Crockett-
21 Billig.

22 MS. CROCKETT-BILLIG: Thank you, Tom.

23 I'd like to begin our regulatory overview with
24 a description of the Intacts product. This non-laser,
25 vision correction alternative reshapes the cornea

1 mechanically by adding material rather than by removing
2 tissue from the cornea. Thickness of the device is a
3 parameter which determines the refractive effect. Three
4 thicknesses were evaluated in this PMA cohort. The
5 proposed indication for this technology is for patients
6 having -1.00 to -3.50 diopters of myopia with 1.00 diopter
7 or less of astigmatism.

8 Intacts consist of two clear segments, each
9 having an arc length of 150 degrees. They are manufactured
10 from PMMA, which has a proven history as an ocular implant
11 material. Intacts are placed in the corneal periphery
12 outside the central optical zone. Because Intacts are an
13 additive technology, they can be removed, if desired, which
14 results in a reversal of the refractive effect. Through
15 exchange procedures, they are also potentially adjustable.

16 I will now review the device's regulatory
17 history. In April of 1991, IDE No. G910034 was approved
18 for the initial product design. A 360 degree ring. Ten
19 subjects were enrolled in this Phase I blind eye trial.
20 This coincided with treatment of our first ten sighted eye
21 subjects in Brazil. Our Phase II sighted eye trial with a
22 360 degree ring was approved in February of 1993. May of
23 1995 initiated the Phase IIa trial for our new 150 degree
24 segment design with five thicknesses and a predicted range
25 of -1.00 to -6.00 diopters.

1 As a result of these clinical data, a new
2 nomogram was devised and the thicknesses were bifurcated
3 into two separate trials. The 0.40 and 0.45 sizes were
4 evaluated in an expanded Phase II trial using the revised
5 nomogram. This trial was initiated in August of 1996.

6 The Phase IIIa trial for the 0.25, 0.30, and
7 0.35 millimeter Intacts was initiated in October of 1996
8 for -1.00 to -3.50 diopters.

9 In November of 1997, the Phase IIIb trial for
10 the 0.40 and 0.45 Intacts was initiated along with a new
11 0.21 millimeter size. Data collection is currently ongoing
12 for this trial. A modular PMA with twelve month data on
13 the combined Phase II and Phase III cohort of 410 eyes for
14 three thicknesses was submitted in July of 1998 and
15 accepted for filing in August of 1998.

16 As I mentioned, the PMA cohort combines Phase
17 IIa and Phase IIIa data for the 0.25, 0.30, and 0.35
18 thicknesses. Ninety eyes from Phase II were implanted and
19 359 eyes from Phase III, for a total of 449 eyes. The
20 available post-operative follow-up ranges from 12 to 24
21 months.

22 Eleven clinical sites, listed here,
23 participated in the Phase II and Phase III trials.

24 The Phase II and Phase III trials have been
25 conducted with the oversight of a data and safety

1 monitoring board. Dr. Michael Lemp is the chairman of this
2 board, and other members include Dr. Gary Foulks and Dr.
3 Thomas Clinch.

4 The inclusion criteria for the PMA cohort was
5 typical of most refractive surgery protocols, except that
6 we required all subjects to have a best spectacle-corrected
7 visual acuity of 20/20 or better. Likewise, the exclusion
8 criteria was typical of most refractive surgery protocols,
9 with consideration given to conditions affecting wound
10 healing and overall corneal health.

11 A comprehensive list of ophthalmic tests were
12 performed to evaluate the safety and effectiveness of the
13 Intacts procedure. The tests included visual acuity, slit
14 lamp, tonometry, refractions, and mesopic contrast
15 sensitivity. Tests were conducted according to the
16 schedule shown here. To help ensure standardization,
17 extensive site training and certification were conducted
18 prior to initiation of these trials.

19 The following tests were performed on subgroups
20 for the PMA cohort. Spectral microscopy, central corneal
21 sensation, A-scan, automated visual field, and slit lamp
22 photography. All patients enrolled at a site were required
23 to participate in the designated subgroup tests.

24 As we look at the demographics for the 449
25 evaluable implant eyes, we see that 51 percent were females

1 and 49 percent were males, the mean age was 39 years, the
2 racial demographics are typical of contemporary refractive
3 surgery studies in the United States.

4 Excellent accountability was maintained
5 throughout the 12 month follow-up period. At the month 12
6 exam, 97.6 percent of subjects were evaluated.

7 The distribution of the PMA cohort is presented
8 here. Of the 449 eyes who received Intacts, 410 eyes were
9 evaluated at month 12. The remaining 39 subjects who were
10 not evaluated at month 12 was for the following reasons --
11 20 had the Intacts removed prior to month 12, 5 subjects
12 were lost to follow-up, 5 subjects missed the month 12
13 exam, 4 subjects completed the month 12 exam but missed the
14 analysis cutoff date, 4 subjects underwent an exchange
15 procedure, and 1 subject had a single segment.

16 Subject enrollment was balanced by thickness
17 within each of the 11 sites to achieve an overall balanced
18 distribution.

19 DR. McCULLEY: While you're changing, let me
20 point out to the panel that there is hard copy of her
21 presentation in your folder, if anyone wants it and hasn't
22 found it.

23 MS. CROCKETT-BILLIG: The initial Phase II
24 nomogram, shown here, was based on Eye Bank eye research.
25 Based on the preliminary data from Phase II, the nomogram

1 was revised for Phase III. The actual range of
2 preoperative refractive error used for each sickness is
3 specified here. As you can see, the actual range was from
4 -0.75 to -3.875 for subjects in the PMA cohort.

5 I'd like to conclude by defining the study
6 population evaluated in our PMA safety and efficacy
7 analyses. Four hundred and fifty-four eyes were evaluated
8 for safety. These data will be presented by Dr. Michael
9 Lemp. For efficacy, I've already described the subject
10 disposition which provides us with 410 evaluable eyes at
11 month 12.

12 I would now like to introduce Dr. David
13 Schanzlin, who is the director of keratorefractive surgery
14 at the University of California in San Diego and chief
15 clinical investigator for KeraVision's clinical trials. He
16 will present an overview of Intacts technology and provide
17 an efficacy assessment.

18 DR. SCHANZLIN: Thank you, Darlene.

19 Before I begin, I would like to state that I
20 have served as a consultant to KeraVision for the last 12
21 years and I am the chief clinical investigator for the
22 studies reported in this PMA cohort. My consulting fees
23 are paid directly to the University. I recently acquired a
24 small equity interest in the company, and KeraVision paid
25 my way to this meeting. In my position on the full-time

1 faculty of the University of California at San Diego, I
2 head up the refractive surgery service and I routinely
3 perform various refractive procedures including PRK and
4 LASIK.

5 Now let's review the basics of this new
6 refractive technology for the correction of myopia. As
7 shown here, the KeraVision Intacts are placed in the
8 peripheral corneal stroma outside of the optical entrance
9 pupil. The fact that the central optical zone is preserved
10 is one of the primary advantages of the Intacts approach.
11 Preservation of the central visual axis coupled with the
12 fact that the Intacts can be removed allow for a refractive
13 effect that is reversible and potentially adjustable.

14 The KeraVision Intacts consist of two clear
15 micro thin PMMA segments, each having an arc length of 150
16 degrees. Each segment is precision lathe-cut to plus or
17 minus 0.1 millimeter, and has a hexagonal cross-section
18 that lies along a conic section. With a fixed outer
19 diameter of 8.1 millimeters and an inner diameter of 6.8
20 millimeters, the Intacts have a large, clear central
21 optical zone. Here, at the superior edge of each segment,
22 you can see a small positioning hole which aids in surgical
23 manipulation of the segments. As mentioned, three
24 thicknesses, or sizes, are presented in this PMA for the
25 correction of myopia.

1 Let's look at the role thickness plays in
2 determining the device's refractive effect. Intacts act as
3 passive spacing elements that change the arc length of the
4 anterior corneal curvature. As shown here, placement of
5 the Intacts in the periphery of the cornea causes local
6 separation of the corneal lamellae. This results in
7 shortening of the corneal arc length which has the net
8 effect of flattening of the corneal curvature. When the
9 Intacts thickness is increased, greater amounts of local
10 separation occur, resulting in increased corneal flattening
11 and greater refractive effect. Hence, the refractive
12 effect achieved by the device is directly related to
13 thickness.

14 To demonstrate how the Intacts create central
15 corneal flattening, we have prepared this animation. Here
16 you can see with the animation thickening of the ring
17 segments. Watch the center of the cornea, how it flattens.
18 Let's look at it again, however, and really this time,
19 rather than looking at the center of the cornea, view the
20 area in the mid-periphery and this time watch how the mid-
21 periphery also flattens within the area inside of the ring
22 itself. The importance of this observation relates to the
23 maintenance of the cornea's natural prolate shape. Intacts
24 are, to our knowledge, the first procedure for the
25 correction of myopia that maintains normal corneal

1 asphericity.

2 The concept of arc shortening was verified by
3 our early Eye Bank eye research. A nearly linear
4 relationship between the device thickness and the change in
5 corneal curvature was established, as demonstrated in this
6 slide. This linear relationship has been confirmed and
7 further refined by our clinical studies. With Intacts,
8 each successive increase of 0.05 millimeters in thickness
9 imparts an additional 0.70 diopters of corneal flattening.

10 Let's review quickly the surgical procedure for
11 the Intacts procedure. First of all, a comment about the
12 ease. This surgery is really very easy to master. Most
13 surgeons are able to complete this surgery in between ten
14 to fifteen minutes. In my hands, this is similar to the
15 time that it takes me to do a LASIK procedure.

16 The procedure begins first by prepping and
17 draping the eye in order to fully isolate the eyelashes
18 from the surgical field. Next, using a geometric center as
19 the reference, the corneal surface is inked with markings
20 that indicate the incision placement site and the final
21 positioning of the Intacts segments. Following pachymetry
22 over the incision, a 15 degree diamond knife is set at two-
23 thirds depth and a 1.8 millimeter peripheral incision is
24 made along the incision mark. The vacuum centering guide
25 is aligned on the central corneal mark and suction is

1 applied.

2 Next, we make the peripheral lamellar channels
3 using a clockwise dissector. The dissector is a blunt
4 instrument that is designed to preferentially dissect up
5 rather than down if it ever gets out of corneal plane. A
6 glide is inserted into a previously created lamellae pocket
7 at the incision site, and the dissector is advanced under
8 the glide to ensure that the lamellar dissection begins at
9 the proper depth. The dissector is then rotated to create
10 the peripheral lamellar channel.

11 DR. McCULLEY: Suction is maintained throughout
12 that?

13 DR. SCHANZLIN: Suction is maintained
14 throughout, and suction times run about one minute, one
15 minute and fifteen seconds for both dissections.

16 DR. McCULLEY: Intraocular pressure during that
17 time?

18 DR. SCHANZLIN: It's lower than LASIK. It's
19 running around 80. The pupil occasionally will dilate and
20 occasionally vision will blur.

21 The counter-clockwise dissection is then made
22 in a similar fashion.

23 The Intacts segment is grasped with special
24 forceps and inserted into the peripheral lamellar channel.
25 There is no resistance encountered to this motion since the

1 channel has already been created. The segment is placed in
2 its final position with a blunt Sinsky hook. It is
3 positioned under the peripheral corneal markings, thus
4 assuring its proper position not only in the outer extant,
5 getting to the 8 millimeter outer diameter, but also
6 keeping it free of the incision site.

7 The second segment is then inserted in a
8 similar fashion. The wound is gently reapproximated with a
9 single 11.0 nylon suture, the knot is buried, topical
10 antibiotic and corticosteroid drops are instilled, and a
11 clear shield is applied.

12 Let's now review the efficacy data from our PMA
13 cohort of 410 subject eyes. To demonstrate the efficacy of
14 KeraVision's Intacts for the correction of myopia, we
15 assessed the variables summarized here -- uncorrected
16 visual acuity, predictability based on cycloplegic
17 refraction spherical equivalent, and stability of
18 refractive effect. And, indeed, all efficacy endpoints
19 were exceeded at month 12, both those specified in the two
20 Intacts clinical protocols as well as those in the FDA
21 guidance document for refractive surgery lasers which
22 included predictability based on manifest refraction
23 spherical equivalent. We recognize, however, that the
24 laser guidance document was written for a range of -1.00 to
25 -7.00 diopters and our data range was from -1.00 to -3.50

1 diopters. However, all of the data presented here today
2 substantial exceed all of these endpoints.

3 Let's look at each of these endpoints in
4 detail. The first protocol endpoint was uncorrected visual
5 acuity. At one year, 97 percent of the patients had
6 unaided visual acuity of 20/40 or better, 74 percent were
7 20/20 or better, and a remarkable 53 percent of the
8 patients had unaided visual acuity of 20/16 or better.

9 The visual recovery at the Intacts procedure is
10 quite rapid. On day 1, 57 percent of the patients were
11 20/25 or better uncorrected vision, 34 percent were 20/20
12 or better, and 13 percent were already 20/16. At month 1,
13 at month 6, and month 12 we see continued improvement, and
14 by month 12 over half of the subjects had achieved an
15 impressive uncorrected vision of 20/16 or better. To our
16 knowledge, no other refractive surgical procedure has such
17 a high percentage of patients achieving better than 20/20
18 visual acuity.

19 Our next protocol endpoint was predictability
20 of refractive effect based on cycloplegic refraction
21 spherical equivalent. As shown here, 68 percent of the
22 subjects were within half a diopter of intended correction,
23 and 90 percent were within a diopter of intended
24 correction. The protocol endpoints for predictability of
25 refractive effect were exceeded by these Intacts

1 performance.

2 As I just noted, 90 percent of the subjects
3 were within a diopter of the intended outcome at month 12.
4 Of the 17 undercorrected subjects, only 1 was more than 2
5 diopters from intended. Similarly, only 1 of the 23
6 overcorrected subjects was more than 2 diopters from
7 intended.

8 When we look at the deviation from plano, 92
9 percent of subjects were within a diopter of plano, only 1
10 of 23 undercorrected subjects had deviation from plano of
11 more than 2 diopters, and, similarly, only 1 of 8
12 overcorrected subjects had a deviation from plano of more
13 than 2 diopters.

14 As previously noted, the protocol endpoint for
15 predictability based on cycloplegic refraction spherical
16 equivalent was exceeded. Likewise, the Intacts
17 predictability based on manifest refraction spherical
18 equivalent is similar and also exceeded FDA guidance
19 document for the refractive surgical lasers.

20 In any longitudinal analysis, it is important
21 to look at a consistent population. Therefore, the same
22 408 subjects were analyzed at each exam point to document
23 the percentage of eyes within plus or minus a half diopter
24 and within plus or minus of 1.00 diopter of intended at 3,
25 6, and 12 months after surgery. We see that the

1 percentages achieved at month 3 remain constant over time.

2 The stability of manifest refraction over time
3 was assessed using a constant population of the 384
4 subjects who had manifest refraction spherical equivalent
5 results available at all exams from month 1 through month
6 12. We can clearly see that stability was achieved three
7 months after surgery and was maintained through month 12.

8 Both the protocol and the FDA guidance
9 endpoints for stability of refractive effect specify that
10 95 percent of subjects should have a change of less than or
11 equal to 1.00 diopter of manifest refraction spherical
12 equivalent between two refractions performed at least three
13 months apart beginning at month 3. As we can see from this
14 table, stability was reached at three months after surgery
15 and confirmed at subsequent test intervals. So this
16 endpoint was met for the cohort overall and for each of the
17 three Intacts sizes at all time intervals.

18 Throughout the remainder of this presentation I
19 will discuss the results within the recommended prescribing
20 range for each of the Intacts thicknesses. To understand
21 how the recommended prescribing range was derived, it is
22 useful to examine the distribution of refractive errors for
23 the cohort subjects by Intacts thicknesses.

24 As you can see, there was some overlap in the
25 preoperative cycloplegic refractions between each of the

1 thicknesses. The preoperative range for the 0.25
2 millimeter size was -0.75 to -1.875 diopters. For the 0.30
3 millimeter size the range was -1.50 to -3.00 diopters. And
4 it was -2.00 to -3.875 for the 0.35 millimeter Intacts.

5 To better assist ophthalmologists in selecting
6 the appropriate Intact size for their individual patients,
7 we have defined a recommended prescribing range, RPR, for
8 each thickness. As indicated by the wide bracketing lines
9 shown here, the ranges are contiguous but do not overlap.
10 The ranges are based on the nominal predicted correction
11 for each thickness plus or minus 0.35 diopters. For the
12 0.25 millimeter Intacts, the recommended prescribing range
13 is a spherical equivalent of between -1.00 and -1.625. For
14 the 0.30 millimeter, we have an RPR of -1.75 to -2.25. And
15 for the 0.35 millimeter Intacts, the RPR is -2.375 to -3.00
16 diopters.

17 Of the 410 subjects in the overall PMA cohort,
18 317 were within this prescribed recommended prescribing
19 range. The number of evaluable eyes within the recommended
20 prescribing range for each thickness is also shown here and
21 the distribution was approximately equal.

22 This table lists the uncorrected visual acuity
23 performance by thickness within the RPR. We had 56 percent
24 of subjects achieving excellent visual acuity at the level
25 of 20/16 or better for all Intacts thicknesses. As we look

1 at the data by thickness, 63 percent of the 0.25 millimeter
2 subjects, 54 percent of the 0.30 millimeter subjects, and
3 49 percent of the 0.35 millimeter subjects achieved this
4 high level of uncorrected visual acuity. Both the protocol
5 and FDA guidance endpoint called for 85 percent of the eyes
6 achieving uncorrected visual acuity of 20/40 or better.
7 Clearly, this criterion was easily surpassed by all three
8 Intacts thicknesses.

9 More importantly, performance at the level of
10 20/20 or better was strong for each of the Intacts
11 thicknesses, with 84 percent, 82 percent, and 66 percent
12 achieving this level. And 78 percent of the RPR cohort
13 overall achieved 20/20 or better visual acuity.

14 The Intacts demonstrated excellent
15 predictability based on cycloplegic refraction spherical
16 equivalent correction achieved at month 12. A tabulation
17 of the predicted versus the achieved correction for the PMA
18 cohort and for the RPR groups is shown here. As you can
19 see, there's excellent correlation between these numbers
20 for both the cohort and the RPR groups.

21 Predictability within the prescribing range was
22 also good for each individual Intacts size. Both the
23 protocol and FDA guidance endpoints specified that an
24 outcome would be considered predictable if 75 percent of
25 subjects were within a diopter of their intended

1 correction. This criterion was clearly exceeded for all
2 three thicknesses.

3 Similarly, all three Intacts thicknesses
4 surpassed the endpoint specified that 50 percent of
5 patients should be within a half a diopter of intended
6 correction. Having satisfied these endpoints for plus or
7 minus 0.50 and plus or minus 1.00 diopter, we conclude that
8 the Intacts procedure is, indeed, predictable.

9 In summary, we believe that the data from this
10 PMA cohort have established the efficacy of the KeraVision
11 Intacts for the correction of myopia with spherical
12 equivalence from -1.00 to -3.00 diopters. We propose a
13 recommended prescribing range for each Intacts thickness
14 that is contiguous and non-overlapping. These ranges,
15 shown here, are based on the nominal predicted correction
16 plus or minus 0.35 diopters.

17 Thank you for your attention. I would now like
18 to introduce Dr. Michael Lemp, president of University
19 Ophthalmic Consultants of Washington and chairman of the
20 data and safety monitoring board for the Intacts clinical
21 trials.

22 Michael?

23 DR. LEMP: Good morning. Thank you, David.

24 Before I begin our assessment of safety, I
25 would like to state that I have served as a consultant for

1 KeraVision for five years and that also I have no ownership
2 interest in the company. I would also like to point out
3 that as chairman of the DSMB board it is my opinion and
4 that of the other members of the DSMB board that this
5 clinical study was exceptionally well designed, conducted,
6 and managed.

7 Now I'll turn to a consideration of the safety
8 issues. To demonstrate the safety of the KeraVision
9 Intacts for the correction of myopia, we assessed the
10 safety variables which are summarized here. Indeed, all
11 the safety endpoints were met at month 12, both those
12 specified in the Intacts clinical profile as well as those
13 in the FDA guidance document for refractive surgery lasers.

14 Now let's look at each of the endpoints in
15 detail. The first criterion was best spectacle-corrected
16 visual acuity. No subject had a best corrected visual
17 acuity of 20/40 or worse at month 12, thus satisfying the
18 endpoint. Best corrected acuity was maintained as 98
19 percent of the subjects were within nine letters of their
20 preoperative value. Four subjects lost ten or more
21 letters, and six gained ten or more letters at 12 months.

22 Again, the protocol endpoint for best corrected
23 maintenance was satisfied for the cohort overall and for
24 each of the individual Intacts thicknesses. There was no
25 statistically significant difference among the Intacts

1 sizes.

2 This slide details the four subjects who had a
3 loss of ten or more letters or two or more lines. It is
4 important to note that each of these subjects was 20/20 or
5 better at the last reported exam.

6 This slide compares best corrected acuity at
7 month 12 to preop by level of visual acuity. All subjects
8 had a best corrected acuity equal to or better than 20/32
9 at the month 12. The one subject who was 20/32 had lost
10 only nine letters from the preop and so was not listed on
11 the previous slide. By the next exam, this subject had
12 returned to a preop baseline of 20/20.

13 But let's look here at the increase in the
14 number of subjects having 20/10, 20/12.5, and 20/16 or
15 better visual acuity. Of the subjects, 90 percent were
16 20/16 or better at month 12, and 99 percent were 20/20 or
17 better. The improvement in best corrected acuity was
18 statistically significant.

19 As I just noted, this improvement in best
20 corrected acuity from preop was statistically significant.
21 An analysis of the month 12 data compared to preop
22 indicated that 19.5 percent of the subjects had an increase
23 of five or more letters or one or more lines at best
24 corrected acuity. Additionally, 34.6 percent had an
25 uncorrected visual acuity equal to or better than their

1 preop best corrected visual acuity.

2 Now let's look at induced cylinder. Per the
3 protocol endpoint, less than 5 percent of subjects were to
4 have an induced cylinder greater than 2.00 diopters. No
5 subject had a cylinder greater than 2.00 diopters. Here we
6 see that 92 percent of subjects were within 0.75 diopter of
7 their preoperative cylinder. This is considered to be
8 within the range of measurement error for repeatability for
9 manifest refractions per Zadnik, et al. Of those subjects,
10 15, who had a cylinder increase of greater than 1.00
11 diopter, 100 percent had a best corrected acuity of 20/20
12 or better, with 93 percent 20/16 or better. The
13 uncorrected acuity for all patients was 20/32 or better,
14 with 80 percent 20/25 or better. It is important to note
15 that no subject had greater than 2.00 diopters of induced
16 cylinder.

17 Mesopic contrast sensitivity both with and
18 without glare was performed at all Phase III sites. A
19 functional acuity contrast test chart was used with a view-
20 in tester. The spatial frequencies tested were 1.5, 3, 6,
21 12, and 18 cpd.

22 The endpoint was met for all spatial
23 frequencies tested both with and without glare. There was
24 no mean decrease greater than 0.1 log unit. While some of
25 the changes in spatial frequencies at the 1.5 and 6 cpd

1 without glare were statistically significant, they were not
2 considered functionally significant. No clinically
3 relevant changes were seen for any of the spatial
4 frequencies evaluated.

5 Specular microscopy was a Phase III subgroup
6 test with four sites participating. The Konan Robo non-
7 contact specular microscopy unit was used. Three regions
8 were assessed -- a central region, a 6 o'clock peripheral
9 region, and a 10 o'clock peripheral region. Cell density,
10 coefficient of variation, and percent hexagonality data
11 were analyzed separately at a reading center at Emory
12 University under the direction of Dr. Edelhauser.

13 In this photograph of an Intacts subject, we've
14 illustrated three regions which we evaluated, they're in
15 yellow, the central, 6, and 10 o'clock peripheral regions.
16 It's important to note that the instrument used was
17 designed primarily for central measurements. Measurements
18 in the corneal periphery are normally difficult to obtain
19 due to the angle of approach. The change in contour of the
20 corneal surface and the peripheral area adjacent to the
21 Intacts presented even more than the usual difficulties in
22 obtaining the 10 o'clock images. These measurements were
23 still taken, however, primarily to rule out changes in the
24 cell morphology near the Intacts.

25 The protocol endpoint for endothelial cell

1 density specifies that the mean density should not decrease
2 by more than 10 percent of the preoperative value. As this
3 slide shows, all three regions clearly met the endpoint
4 criteria, and all observed losses were substantially less
5 than the specified maximum of 10 percent.

6 This slide compares endothelial cell density
7 for the initial treated eyes and the untreated fellow eyes
8 at month 6. Since untreated fellow eyes were eligible for
9 surgery after month 6, this comparison was not available
10 for later exam points. If we look at those eyes with a
11 decrease greater than 10 percent, the central region shows
12 that three of the initial or treated eyes and two of the
13 fellow eyes had these readings. In the peripheral 6
14 o'clock region, seven of the treated eyes and six of the
15 fellow eyes also had this. And in the peripheral 10
16 o'clock region, six treated eyes and five fellow eyes.

17 There were no statistically significant
18 differences between the initial treated eye and the
19 untreated fellow eye for any region evaluated. These data
20 indicate that variability in the readings appears to be a
21 reflection of the reproducibility of the measurement
22 technique and not an indication of an ongoing safety issue.

23 The only data we were able to find in the
24 literature on longitudinal peripheral endothelial cell
25 densities assessed in a serial fashion was a paper by

1 Trocme, et al. in which peripheral cell densities were
2 assessed preoperatively and following PRK. As you can see,
3 the loss from preop to month 12 for the PRK was reported at
4 6.9 percent compared to 1.9 percent for the Intacts in the
5 peripheral 10 o'clock position.

6 The coefficient of variation was actually
7 improved at both month 6 and month 12 exams with the
8 improvement at month 12 being statistically significant at
9 the 0.02 level. The fact that the coefficient of variation
10 was actually improved speaks against any morphological
11 change associated with the Intacts procedure. No
12 statistically significant change was seen in the percent
13 hexagonality analysis. Also, there were no cellular
14 morphologic changes that would indicate that the
15 endothelium was compromised due to the Intacts procedure.

16 We now turn to the intraoperative clinical
17 findings. The Intacts procedure was successfully completed
18 98.9 percent of the time. The intraoperative adverse
19 events rate was 0.2 percent. One subject experienced an
20 anterior chamber perforation. This was related to either
21 an incorrect diamond knife setting and/or a deviation in
22 the pocketing technique that resulted in an overly deep
23 pocket. The event was not related to the product, and the
24 subject has fully recovered.

25 There were six subjects who experienced

1 intraoperative ocular complications as defined in the
2 protocols. three of these six complications were related
3 to a deviation in the surgical protocol, two were related
4 to subject movement, and the remaining intraoperative
5 complication was related to an allergic reaction to the
6 antiseptic material used. None of these intraoperative
7 complications was related to the device.

8 There were five safety related adverse events
9 during the twelve month reporting period. Two anterior
10 chamber perforations, one which was just discussed as an
11 intraoperative event during the initial implant procedure,
12 and a second one which occurred during an exchange
13 procedure. Both incidents were related to deviations from
14 the surgical procedure; that is, incorrect knife settings
15 and pocketing technique. One subject had an infection
16 during the early post-operative course which completely
17 resolved. One subject had a decrease of two or more lines
18 of best corrected acuity over two consecutive exams which
19 subsequently resolved. And one subject had a shallow
20 placement of the temporal segment which was subsequently
21 removed. The nasal segment remains in place and the
22 subject has had a good visual outcome.

23 Looking at a breakdown of the adverse events by
24 protocol, we see that three of these adverse events
25 occurred in Phase II, for an adverse event rate of 3.3

1 percent. The two adverse events in Phase III resulted in
2 an adverse event rate of 0.6 percent, and a combined
3 adverse event rate for the PMA cohort is 1.1 percent.

4 This slide provides the visual outcomes for the
5 five adverse event subjects. The case of infection was
6 resolved and the subject is stable with 20/16 best
7 corrected vision. As a result of this early Phase IIa
8 adverse event, an infection prophylaxis guideline was
9 implemented in March of 1997. No infections occurred in
10 the Phase IIIa subsequent to the guidelines'

11 implementation. As previously noted, the subject with a
12 shallow placement of the temporal segment still has the
13 nasal segment in place. Her uncorrected acuity is 20/20,
14 and best corrected acuity is 20/12.5. The subject with a
15 loss of two or more lines of best corrected acuity for two
16 consecutive exams was stable at month 12 interim exam with
17 no loss of best corrected acuity. The subject with an
18 anterior chamber perforation after the initial surgery has
19 a best corrected acuity of 20/12.5 and is scheduled for a
20 second surgery. The subject with an anterior chamber
21 perforation during the exchange procedure subsequently had
22 the Intacts removed. The subject recently underwent a
23 successful exchange procedure and the best corrected acuity
24 is now 20/20.

25 It is important to note that the best corrected

1 acuity for all of these adverse events patients was 20/20
2 or better at their most subsequent and most recent exam.
3 No additional adverse events have occurred for the
4 contralateral eyes in the PMA cohort. And the total
5 enrollment now stands at 735 eyes, and the combined adverse
6 event rate is 0.7 percent.

7 Central corneal sensation was a subgroup test
8 involving nine sites. A Cochet-Bonnet anesthesiometer was
9 used. A subject was considered to have a significant loss
10 if the decrease in central corneal sensation was 20
11 millimeters or more from the preoperative exam. To our
12 knowledge, there are no studies available which validate
13 the reproducibility of this methodology.

14 With that in mind, let's review the data. At
15 month 12, 5.5 percent of the subjects had equal to or
16 greater than 20 millimeters reduction in central corneal
17 sensation. It's interesting to note that six of the
18 thirteen subjects here were from one site, perhaps
19 indicating some subjective variability in the method. It
20 is also important to mention that no subject had a complete
21 loss of corneal sensation, or had an epithelial defect, or
22 had any clinical sign of neurotrophic keratitis. At their
23 most recent exam, all of the 13 subjects have returned to
24 their preoperative baseline.

25 We evaluated the reduction in central corneal

1 sensation by thickness and found no statistically
2 significant relationship.

3 The incidents of clinically significant
4 complications noted between month 9 and month 12 was low.
5 We saw the following. Four subjects with a loss of equal
6 to or greater than ten letters or two lines of best
7 corrected acuity; three have returned to within ten letters
8 or two lines, the remaining subject was 20/20 with an
9 eleven letter loss. The one subject who had a best
10 corrected acuity worse than 20/25 has subsequently returned
11 to baseline. Fifteen subjects had an induced cylinder
12 greater than 1 diopter. The range was between 1.25 and
13 1.50 diopters for twelve of the fifteen, and 1.75 diopters
14 for three of the fifteen.

15 As noted earlier, all of these subjects had a
16 best corrected acuity of 20/20 or better with 93 percent
17 20/16 or better. All cylinder subjects had an uncorrected
18 acuity of 20/32 or better with 80 percent 20/25 or better.
19 No subject had more than 2 diopters of induced cylinder.

20 As we look at neovascularization, I'd like to
21 note that all of the cases were in the region of the
22 incision site. What we saw between the month 9 and month
23 12 interval were six subjects with pannus, four of whom had
24 had preexisting pannus. Five subjects had a single deep
25 vessel. None of these vessels were progressive. All five

1 subjects had a history of contact lenswear, and one subject
2 had a deep vessel preoperatively, two had preexisting
3 pannus, and one of the five had the vessel resolved by the
4 month 12 exam. One subject was reported to have a
5 persistent epithelial defect located temporally which was
6 attributed to dry eye during the month 12 test interval.
7 This was not considered to be related to the device. The
8 onset was at month 9 and the incident had been resolved by
9 month 12. One subject had an onset of Uveitis at the month
10 9 interim exam, it had completely resolved by the month 12
11 exam. Thirteen of the 237 subjects experienced a reduction
12 in central corneal sensation of equal to or greater than 20
13 millimeters at the month 12, and all 13 have returned to
14 within 20 millimeters of their preoperative central corneal
15 sensation reading.

16 Twenty-three, or 7 percent of subjects, have
17 reported having "always" and "severe" visual symptoms at
18 month 12. These symptoms included difficulty with night
19 vision, blurry vision, double vision, glare, halos, and
20 fluctuating distance vision. Although exact comparisons
21 are difficult to obtain due to differences in reporting
22 methodologies, these numbers are similar to those reported
23 for other approved refractive technologies.

24 We evaluated the data by thickness for the 23
25 subjects who reported visual symptoms, and you can see the

1 data here. It's important to note that 11 of the 13
2 subjects with 0.35 millimeter Intacts had a 1 diopter or
3 more deviation from plano. As you recall from Dr.
4 Schanzlin's presentation, many, 29 percent, of our enrolled
5 0.35 subjects were outside the proposed recommended
6 prescribing range for this thickness.

7 The question then becomes what is the
8 significance of these visual symptoms. We looked at the
9 following data to answer this question. How many of the 23
10 subjects with "always" and "severe" symptoms at month 12
11 elected to have their second eye treated? Eleven, or 48
12 percent, had a contralateral eye procedure. We also saw
13 that three of the twenty-three subjects went on to request
14 a removal procedure, and four had exchange procedures.

15 We analyzed the data to determine probable
16 reasons for visual symptoms. Post-operative deviation from
17 plano had a statistically significant relationship with
18 frequency and severity of difficulty with night vision,
19 diplopia, and blurry vision. A statistically significant
20 relationship was seen between the month 12 manifest
21 refractive cylinder and the frequency and severity of
22 halos, diplopia, fluctuating vision, and blurry vision.
23 Other variables associated with visual symptoms included
24 mesopic pupil diameters of 7 millimeters or greater, and
25 preoperative RGP contact lens wear most likely associated

1 with some degree of corneal warpage.

2 Finally, the results of a self-administered
3 patient survey indicated that 90 percent of the subjects
4 were satisfied with their initial Intacts procedure, and 95
5 percent of those who had had bilateral procedures.

6 Let's move on to a consideration of
7 reversibility and adjustability. As previously mentioned,
8 one of the unique features of this additive technology is
9 that it can be removed, if desired. And because there's no
10 surgical invasion of the central visual axis of the cornea,
11 a reversal of refractive effect is possible. Intacts can
12 be easily removed in a brief outpatient procedure. It
13 takes approximately about five minutes. No clinically
14 significant complications or sequelae have been associated
15 with the removal procedure. And this feature provides a
16 unique option for restoring a patient's eye to its previous
17 optical performance.

18 The proposed claim of reversibility is based
19 primarily on two criteria, preservation of a subject's best
20 corrected visual acuity and the ability of the subject to
21 return to within 1.00 diopter of the preoperative
22 refraction. These criteria were derived from the FDA
23 guidance document.

24 Of the 34 removal procedures, one, and
25 infection case, was removed for safety reasons, fourteen

1 were removed due to dissatisfaction with the correction
2 achieved, sixteen for dissatisfaction due to visual
3 symptoms, and three for other removals. One subject in
4 this latter category had the Intacts removed from two eyes
5 due to FAA restrictions for pilots having "experimental"
6 procedures. And this subject had an uncorrected acuity of
7 20/16 in one eye and 20/25 in the other. One subject was
8 explanted due to an anterior chamber perforation during the
9 exchange procedure and has subsequently undergone a
10 successful exchange procedure.

11 Of the 34 removals, comparison to the
12 preoperative status is provided for 28 subjects with three
13 months of postremoval data available. Per the protocol,
14 patients were exited after month 3, hence our selection of
15 this time period for postremoval analysis. All subjects
16 were within ten letters or two lines or better of their
17 preoperative best corrected acuity, and all subjects were
18 20/20 or better.

19 Postremoval predictability data indicates that
20 81 percent of the subject eyes returned to within a half a
21 diopter of their preoperative manifest spherical
22 equivalent, and 96 percent returned to within 1.00 diopter.
23 The one subject with a manifest spherical equivalent
24 greater than 1.00 diopter compared to preop actually had an
25 improvement in the manifest spherical equivalent.

1 Here we see the 26 subjects, or 96 percent,
2 that returned to within 1.00 diopter of their preop
3 manifest spherical equivalent. The one subject who was not
4 within the 1.00 diopter actually had a decrease in his
5 refractive error.

6 In looking at the manifest refraction cylinder
7 at month 3 postremoval exam, we see that 100 percent of the
8 subject eyes returned to within 1.00 diopter, and 93
9 percent returned to within a half a diopter based on
10 manifest refraction cylinder by the month 3 postremoval
11 exam.

12 This slide shows the stability of the manifest
13 refraction within 1.00 diopter over time. All subjects
14 were stable within 1.00 diopter manifest spherical
15 equivalent by the month 3 postremoval.

16 The uncorrected acuity reversibility data
17 indicate that 93 percent of subject eyes returned to within
18 two lines or better of their preoperative uncorrected
19 vision, and 82 percent returned to within one line or
20 better of their preoperative uncorrected acuity.

21 Reversibility of visual symptoms is defined as
22 returning to the same or better level of severity and
23 frequency as compared to the preoperative level of the
24 referenced symptoms. By symptom, 81 to 96 percent of
25 subjects had a reversal in their visual symptoms upon

1 removal. To assess reversibility, we looked at several
2 refractive criteria.

3 Accounting for subjects who were within or
4 improved from preop, 89 percent were within a half a
5 diopter of manifest spherical equivalent, a half a diopter
6 of cylinder, and one line of best corrected acuity; 100
7 percent were within three-quarters of a diopter of manifest
8 refraction, three-quarters of a diopter of cylinder, and
9 one line of best corrected acuity. Obviously, 100 percent
10 were well within the criteria derived from the FDA guidance
11 document listed here.

12 In review, 100 percent of subjects maintain a
13 best corrected acuity of 20/20 or better, all were within
14 1.00 diopter or better of their preoperative manifest
15 spherical equivalent, 100 percent had a stable postremoval
16 refraction from month 1 to month 3, and 93 percent were
17 within a half a diopter of their preoperative manifest
18 refraction cylinder, and a majority of visual symptoms were
19 reported at levels equivalent or better than their
20 preoperative levels.

21 Adjustability of the Intacts refractive effect
22 is achieved through an exchange procedure in which the
23 Intacts of one thickness are removed and substituted for
24 new ones of a different thickness. Twelve, or 2.7 percent,
25 of the subject eyes had an exchange procedure as of the

1 November 12, 1998 clinical update. All 12 exchanges were
2 for undercorrection. The subjects were eligible for an
3 exchange procedure after six months and provided they met
4 the specified protocol criteria.

5 This slide provides a comparison of pre-
6 exchange and post-exchange deviation from plano. Post-
7 exchange, over half of the subjects shifted to within 1.00
8 diopter of plano. We saw fewer subjects with a deviation
9 greater than plano after their exchange procedure, and no
10 subject greater than 2.00 diopters post-exchange. Also, 73
11 percent of the subjects were 20/40 or better, and 55
12 percent were 20/25 or better after their exchange
13 procedure.

14 All subjects had increased refractive effect
15 following an exchange procedure, 73 percent had an improved
16 uncorrected acuity, and three exchange subjects went on to
17 have their Intacts removed due to continuing
18 undercorrection.

19 In conclusion, let's consider the overall risks
20 and benefits of the procedure. The performance data for
21 Intacts were excellent, allowing us to conclude that the
22 Intacts effectively reduce myopia between 1.00 and 3.50
23 diopters for subjects with 1.00 diopter or less of
24 astigmatism. The uncorrected visual outcomes were
25 outstanding, and best corrected visual acuity was

1 maintained. It is important to note that no clinically
2 significant harm to any subject occurred. The rapid visual
3 recovery was demonstrated and that effect appears to be
4 reversible upon removal.

5 Reviewing the procedure itself, we see that it
6 does not compromise the central optical zone, that the
7 cornea's natural physiological shape is maintained; i.e.,
8 an aspheric or prolate surface. Since Intacts correct
9 myopia by mechanical rather than surgical remodeling of the
10 cornea, no tissue removal is required and no long-term
11 corticosteroid therapy is required. And finally, the
12 procedure is relatively easy to learn and yields consistent
13 refractive results.

14 As with any procedure, there are some risks.
15 Patients may need another procedure if the results are not
16 satisfactory. Visual symptoms may occur following the
17 procedure, visual symptoms primarily related to the post-
18 operative deviation from plano and cylinder. Patients may
19 experience induced astigmatism which typically decreases
20 over time, however, the visual acuity does not seem to be
21 significantly affected. Infection is a risk with any
22 surgical procedure. It is important to note that the
23 adverse event rate however was quite low.

24 All safety endpoints have been met. It has
25 been demonstrated that Intacts are well-tolerated in the

1 cornea. And no patient has had a clinically significant
2 loss of best corrected acuity.

3 All the efficacy endpoints have been met or
4 exceeded. Excellent visual acuities have been achieved.
5 The performance is predictable. And the refractive effect
6 is stable.

7 The proposed indication for use of the
8 KeraVision Intacts is for the reduction or elimination of
9 myopia of -1.00 to -3.00 diopters at the spectacle plane in
10 patients who are 21 years of age or older, with a
11 documented stability of refraction as demonstrated by a
12 change of less than or equal to 1.00 diopter for at least
13 six months prior to the preoperative exam, and with
14 preoperative myopic error ranging from -1.00 to -3.50
15 diopters with 1.00 diopter or less of astigmatism.

16 And so in summary, the KeraVision Intacts we
17 believe represents a safe and effective alternative for the
18 correction of myopic errors. The benefit-risk ratio is
19 very favorable with excellent visual outcomes, the ability
20 to return a patient's eye to its preoperative optical
21 performance, and a very low adverse event rate. The
22 clinically relevant point for any refractive surgery
23 procedure is how well do patients see. And these patients
24 see extremely well. Over half of the subjects have
25 corrected visual acuities of 20/16 or better. As chairman

1 of the Data and Safety Monitoring Board, I believe this
2 study has demonstrated the safety and effectiveness of the
3 Intacts for the indicated use.

4 This concludes our presentation. Thank you for
5 your attention.

6 DR. McCULLEY: Thank you.

7 We will recess and take a 15-minute break.
8 Please note the clock. It is 9:37.

9 (Recess.)

10 DR. McCULLEY: We're now going to begin with
11 the clinical review by Dr. Malvina Eydelman. But prior to
12 that, Ms. Thornton has a couple of introductions.

13 MS. THORNTON: Before we proceed, I would like
14 to introduce to you and to our panel and our staff here the
15 two people that we've waited for and that came in for the
16 properdin of the sponsor's presentation. They have had a
17 great deal of difficulty getting here and we really
18 appreciate their perseverance.

19 The first person I would like to introduce to
20 you is Dr. Renee Middleton. She is our interim consumer
21 representative who has graciously agreed to attend in the
22 place of our panel consumer representative, Ms. Lynn
23 Morris. Dr. Middleton comes to us from the Ear, Nose, and
24 Throat Devices Advisory Panel. She is an assistant
25 professor in the Department of Counselling and Counselling

1 Psychology, and director of Human Resources and Outreach in
2 the College of Education at Auburn University in Auburn,
3 Alabama. Thank you very much for your perseverance. We
4 really appreciate your coming here and taking up our
5 consumer representative banner for today.

6 The other person I'd like to introduce to you
7 is a voting member of our panel, Dr. Jose Pulido, professor
8 and chair of the Department of Ophthalmology at the
9 University of Illinois Eye and Ear Infirmary in Chicago,
10 Illinois. Welcome to you also, Dr. Pulido.

11 Thank you, Dr. McCulley.

12 DR. MCCULLEY: Malvina?

13 DR. EYDELMAN: Good morning. I would to thank
14 the sponsor for providing me with a copy of their
15 presentation prior to this meeting, allowing me to avoid
16 redundancy in my presentation. Today I will, therefore,
17 only highlight some points for panel consideration and will
18 not present a comprehensive review of the clinical study in
19 this PMA.

20 The sponsor has eloquently summarized their
21 specular microscopy outcomes in the presentation you just
22 heard. I would like to bring to your attention some
23 additional outcomes. For all implants, the mean central
24 cell loss at 12 months was 0.41 percent. Since contact
25 lens users are known to experience a transient rise in

1 endothelial cell density after discontinuation of contact
2 lens wear that might offset the initial decrease in
3 endothelial cell density, the sponsor was asked to stratify
4 the cell density outcomes by preoperative use of contact
5 lenses. The mean central endothelial cell loss over one
6 year for contact lens wearers was 0.15 percent, and 0.98
7 percent for those who wore glasses preoperatively.

8 No operative eye in this study had a decrease
9 in cell density of 10 percent or more in the central
10 region. Stratifying the central loss by implant thickness
11 revealed no statistically significant differences.

12 Peripheral cell density outcomes, however,
13 differed significantly from central. For all implants, the
14 mean loss was 1.8 percent at 6 o'clock, and 1.9 percent at
15 10 o'clock position. Thirteen eyes had a decrease in
16 peripheral cell density of 10 percent or more. Six
17 subjects had a decrease at the 6 o'clock position, and
18 seven at ten. Out of the thirteen eyes with a decrease in
19 the peripheral endothelial cell density of 10 percent or
20 more, four had 0.25 millimeter implant, three had 0.30, and
21 six had 0.35 millimeter implant.

22 Cell density analysis stratified by thickness
23 revealed statistically significant differences between the
24 thicknesses for change between preop and month 6 as well as
25 preop and month 12 for 10 o'clock region only. At 10

1 o'clock location, between preop and month 6, 12 percent of
2 eyes implanted with 0.35 millimeter Intacts had cell
3 density loss of greater than 10 percent as compared to 2
4 percent for 0.25 and 3 percent for 0.30 Intacts. At the
5 same location, between preop and month 12, 15 percent of
6 eyes implanted with 0.35 ring had cell density loss of
7 greater than 10 percent as compared to 2 percent for 0.25
8 Intacts, and 3 percent for 0.30.

9 Furthermore, the mean peripheral cell change at
10 twelve months for the 10 o'clock position also showed
11 statistically significant differences between the
12 thicknesses. There was 4.95 percent loss associated with
13 0.35 millimeter implant as compared to 0.24 percent
14 increase with 0.30, and a 1.5 percent loss associated with
15 0.25 ring.

16 Given that the device in question is a corneal
17 implant, the sponsor was requested to perform several
18 analyses to try to establish lack of progression of
19 endothelial cell loss. The mean endothelial cell density
20 at the preoperative month 6 and month 12 exams was compared
21 for each of the corneal regions. There was a statistically
22 significant decrease in endothelial cell density between
23 month 6 and month 12 for all three regions. Sponsor
24 compared findings at the central region to the expected
25 decrease in normal eyes reported by Bourne, et al., of 0.6

1 percent per year and found the differences from month 6 to
2 month 12 not to be statistically significant.

3 Peripherally, at 6 and 10 o'clock, mean change
4 from month 6 to month 12 was 1.4 percent and 1.1 percent
5 respectively. The expected rate of decrease in the
6 peripheral corneal region of normal subject eyes could not
7 be found in the literature, so a comparison of the mean
8 peripheral change to the rate of loss in normal eyes was
9 not possible.

10 The panel is being asked to consider the
11 outcomes of all the endothelial cell density analyses and
12 comment on whether sufficient assurance of safety is
13 currently available or whether additional analysis and/or
14 additional longitudinal follow up is needed prior to
15 conclusion.

16 Visual symptoms in the study were assessed for
17 frequency of occurrence as well as magnitude. Analysis of
18 visual symptoms reported as "often" or "always" between the
19 preoperative and month 12 exams for all thicknesses
20 combined reveals some interesting findings. Glare
21 frequency increased from 1.7 to 9.7 percent, halos from 0.3
22 percent to 11.9 percent, difficulty with night vision from
23 5.9 percent to 17.3 percent, blurry vision from 1.5 to 12.2
24 percent, sensitivity to light increased to 7 percent from
25 3.4, double images from none preoperatively to 6.7 percent

1 at month 12. All of the visual symptom increases was
2 statistically significant with the exception of fluctuating
3 distance vision.

4 There was statistically significant changes
5 between the frequency of double images, fluctuating near
6 vision, and fluctuating distance vision between
7 preoperative and month 12 exams among different
8 thicknesses.

9 Double images were reported as occurring at a
10 frequency of "often" or "always" at the month 12 exam by
11 3.7 percent of subjects with 0.25 Intacts, 5.5 percent with
12 0.30 Intacts, and 10.9 percent with 0.35 Intacts.
13 Fluctuating near vision was reported as occurring at a
14 frequency of "often" or "always" by 1.8, 1.8, and 7.3
15 percent for the 0.25, 0.30, and 0.35 millimeter Intacts
16 respectively. Fluctuating distance vision was reported as
17 "often" or "always" by 0.9, 1.8, and 7.3 percent of
18 subjects with the three rings, as seen on this slide.

19 Magnitude of visual symptoms in this study,
20 combining moderate and severe symptoms, is presented here.
21 Glare occurred at 16.1 percent, halos at 17.2 percent,
22 difficulty with night vision 19.1 percent, and blurry
23 vision at 15.9 percent. Results for light sensitivity,
24 double images, and fluctuating vision of moderate or severe
25 magnitude are seen here.

1 Among Intacts thicknesses, a statistically
2 significant difference in symptom magnitude was seen as
3 well. Double images were reported at a magnitude of
4 moderate or severe at month 12 by 6.1, 6.9, and 15.3
5 percent of subjects with the 0.25, 0.30, and 0.35
6 millimeter Intacts respectively. Fluctuating distance
7 vision was reported at a magnitude of moderate or severe by
8 4.7, 8.8, and 14.3 percent of subjects with the three
9 different sizes.

10 The panel is being asked to consider visual
11 symptoms data in their recommendations for safety outcomes
12 of each of the three thicknesses of the Intacts. Aspects
13 of this data the panel feels are important to be present in
14 the labelling need to be identified.

15 Reductions in central corneal sensation that
16 were 20 millimeters or more from the preoperative baseline
17 was seen for 9 percent of subject eyes at month 6, 5
18 percent at month 12, 3 percent at month 18, and 2 percent
19 of the eyes at month 24. Of the subject eyes with
20 reduction of central corneal sensation that were 20
21 millimeters or more, 13 out of 24 subject eyes examined at
22 month 6, and 6 out of 13 subject eyes examined at month 12
23 were observed at a single clinical site.

24 Review of the implant videotapes for these
25 subject eyes suggest that the reduction in central corneal

1 sensation may have been related to the surgical technique.
2 It appears that in each of these cases there may have been
3 some excessive manipulation of the incision site during the
4 pocketing procedure which may have contributed to the loss
5 of central corneal sensation. The panel's recommendation
6 for appropriate labelling addressing corneal sensation loss
7 is requested.

8 The sponsor has already presented the
9 predictability data for the overall device as well as for
10 each of the thicknesses. This slide is meant just to draw
11 your attention to the fact that the predictability of the
12 0.35 millimeter implant to achieve both plus or minus half
13 a diopter and plus or minus 1.00 diopter of intended is
14 statistically lower than the thinner implants.

15 The average achieved corrections for each
16 implant thickness in the PMA cohort were as follows: 0.25
17 millimeter, 1.48 plus or minus 0.52 diopters; for 0.30
18 millimeter, 2.07 plus/minus 0.56; and for 0.35 millimeter,
19 2.76 plus/minus 0.69 diopters. As you have already heard,
20 the sponsor has revised the recommended prescribing range
21 for each implant size. While the sponsor's nominally
22 predicted corrections of 2.00 and 2.70 diopters closely
23 mimic the mean achieved correction for both the 0.30 and
24 0.35 millimeter implants, for the 0.25 millimeter Intacts
25 the sponsor is recommending a nominally predicted

1 correction of 1.30 diopters even though the mean achieved
2 correction was 1.50 diopters. Panel members are being
3 asked to consider this in their recommendations for
4 labelling.

5 The combined recommended prescribing range for
6 all three thicknesses currently spans from -1.00 diopter to
7 -3.00 diopters of myopia. The proposed indication for use
8 for this device reads as follows: "KeraVision Intacts are
9 intended for the reduction or elimination of myopia of
10 -1.00 to -3.00 diopters at the spectacle plane in patients
11 who are 21 years of age or older, with documented stability
12 of refraction as demonstrated by change of less than or
13 equal to 1.00 diopter for at least six months prior to
14 preoperative examination, and with preoperative myopic
15 error ranging from -1.00 to -3.50 diopters with 1.00 or
16 less of astigmatism."

17 There were 40 eyes treated with 0.35 millimeter
18 Intacts that had preoperative CRSE greater than 3.00
19 diopters. Outcomes for these eyes are presented on this
20 slide. As you can see, while only 43.6 percent of these
21 attained CRSE within half a diopter of preop, 87.5 percent
22 achieved 20/40 or better uncorrected visual acuity. Please
23 consider these outcomes in your recommendations regarding
24 the appropriate upper myopic range for the indication for
25 use statement.

1 As you have already heard, out of 34 eyes in
2 the removal cohort of this PMA there has been one safety
3 related removal due to infectious keratitis, 16 removals
4 occurred due to dissatisfaction with visual symptoms
5 achieved. Visual symptoms that prompted removal were night
6 vision, glare, halos, and double images, a finding that's
7 consistent with the outcomes of subjective assessment of
8 visual symptoms. Fourteen removals occurred due to
9 dissatisfaction with the correction, seven of these
10 occurred due to dissatisfaction with correction secondary
11 to induced astigmatism. Among the three removals
12 classified as "Other," one was associated with a deferred
13 exchange procedure secondary to anterior chamber corneal
14 perforation with a diamond knife during an attempted
15 exchange procedure.

16 The mean duration of ICRS implant time prior to
17 removal was 10.3 months. The removals due to
18 dissatisfaction related to visual symptoms tended to occur
19 at an earlier time point than removals due to
20 dissatisfaction with outcomes related to correction. Even
21 though refractive stability for this device was established
22 to occur by six months, most removals occurred after month
23 6. This was one of the factors that prompted FDA to
24 request a minimum of 12 months follow-up prior to PMA
25 submission.

1 This submission provides distribution of
2 Intacts removal by thickness. Five removals occurred in
3 eyes with 0.25 millimeter Intacts, 10 in 0.30, and 19
4 removals in eyes with 0.35 millimeter Intacts. Thus, 56
5 percent of removals occurred in the 0.35 millimeter Intacts
6 eyes. This is consistent with the earlier discussions of
7 greater visual symptoms and decreased accuracy of
8 correction achieved with the 0.35 millimeter device.

9 Sponsor's reversibility claim is based upon
10 data from 28 eyes that have completed a month 3 postremoval
11 exam. Six of these exams were not conducted within the
12 designated exam window and ranged in time from two to nine
13 months postremoval. Thus, three months postremoval data
14 within the designated exam window is available for 22 eyes.
15 No data analysis beyond three months postremoval is
16 currently available.

17 Twenty-seven out of 28 eyes analyzed by the
18 sponsor had manifest refractive cylinder data. All of
19 these returned to within plus or minus 1.00 diopter
20 manifest refraction cylinder by the month 3 postremoval
21 exam. For the one eye with no MR cylinder value available
22 at month 3 there was an increase in cycloplegic refractive
23 cylinder of 1.25 diopters. Analysis of the seven subjects
24 who had an increase of greater than 1.00 diopter of
25 manifest cylinder with the implant shows that upon removal

1 the change in the resultant cylinder as compared to preop
2 varied from minus 0.25 diopters to plus 0.75 diopters.

3 Removal subjects from Phase IIa reported in the
4 severity of visual symptoms. At the preoperative exam, the
5 severity of all key symptoms were rated as none or mild.
6 At the month 3 postremoval exam, two subjects reported a
7 total of three key symptoms as severe. All remaining key
8 symptoms were reported as none or mild. Out of the 22 eyes
9 with three months postremoval data, 19 had both month 1 and
10 month 3 exams. All of these 19 eyes were within plus or
11 minus 1.00 diopter of MRSE, and 12 of these eyes were
12 within plus or minus of half a diopter MRSE.

13 Panel members are being asked for their
14 assessment of the reversibility data. Specifically, is
15 MRSE plus or minus 1.00 diopter in 19 eyes between one and
16 three months postremoval exams sufficient to establish
17 stability at three months postremoval? Is the data
18 currently available for the removal cohort sufficient to
19 make the reversibility claim for this device? If not, how
20 much data and what length of follow-up is needed before
21 such a claim can be established?

22 There have been a total of 12 exchange
23 procedures as of 11-12-98. As you have already heard, all
24 subjects to date were exchanged for reasons related to
25 undercorrection. It is important to point out that 11 of

1 the 12 exchanges performed involved switching to 0.40 and
2 0.45 millimeter Intacts, thicknesses which are not
3 currently available and are not subject for evaluation in
4 this PMA.

5 In the analysis of the post exchange outcomes,
6 sponsor has combined data from all the most recent post
7 exchange exams. These ranged from seven days to eighteen
8 months post exchange. The outcomes of this most recent
9 post exchange data analysis shows a range of the change in
10 MRSE achieved as a result of exchange to span from 0.12
11 diopters to 1.50 diopters. Unfortunately, the sponsor has
12 not yet been able to identify a precise nomogram for
13 predicting the ideal size of the Intacts needed for
14 exchange. As a result, residual post exchange MRSE in the
15 12 eyes range from minus 2.00 diopters to plus 0.25
16 diopters. Three subjects had their Intacts subsequently
17 removed due to continuing undercorrection.

18 The panel is being asked to consider data from
19 these exchange procedures in consideration of sponsor's
20 claim for adjustability of refractive effect.

21 Current proposed labelling for this device
22 reads as follows: "An analysis was performed to assess the
23 percentage of subjects who achieved enhanced visual
24 performance at month 12. The results indicate that 21.5
25 percent of patients had a post-operative UCVA better than

1 the preoperative BSCVA, 32.7 percent of patients achieved
2 their post-operative UCVA better than 20/20 and equal to or
3 better than the preoperative BSCVA, 19.5 percent of
4 patients achieved their post-operative BSCVA that was
5 better than the preoperative BSCVA. These results indicate
6 that the KeraVision Intacts may be capable of increasing
7 the resulting power of the cornea to provide for enhanced
8 visual capabilities for some patients. Patients with the
9 KeraVision Intacts may be able to achieve better visual
10 acuity than they could achieve with their prior methods of
11 vision correction. This finding may occur because the
12 natural prolate shape of the cornea has been maintained
13 inside the central optical zone for these patients."

14 There were 80 eyes with a month 12 BSCVA that
15 had better than their preoperative BSCVA. Phase II data is
16 analyzed in terms of lines change, and Phase III in terms
17 of letter change. For the 18 eyes in this category in
18 Phase II, the mean BSCVA line change was 1.1 plus or minus
19 0.3, and the mean BSCVA letter change for 62 eyes in Phase
20 III subjects was 6.2 plus or minus 1.7 letters.

21 This enhanced visual performance statement is
22 based upon a change of less than two lines of Snellan VA.
23 In another analysis in this PMA, in calculating the amount
24 of BSCVA loss following the procedure, the sponsor did not
25 count the one line loss as clinically meaningful. The

1 sponsor believes that this inconsistency can be explained
2 by skewed preoperative distribution and the fact that the
3 ETDRS chart only measures to a visual acuity of 20/10, thus
4 limits accordingly a gain of 10 or more letters for those
5 subjects with preoperative BSCVA better than 20/16.

6 This slide contains a summary of the month 12
7 UCVA compared to the preoperative BSCVA for the other two
8 categories of the enhanced visual performance claim. The
9 mean line change is 1.1 and the mean letter change is 3.5
10 for subjects with post-op UCVA better than preop BSCVA.
11 The mean changes for the last group are 0.3 lines and 2.7
12 letters.

13 In light of this discussion, panel members are
14 being asked to comment on the enhanced visual performance
15 section in the labelling.

16 This concludes my comments and I will now
17 restate the questions for panel consideration.

18 Question 1. "Do the outcomes of the
19 endothelial cell density analysis presented in this PMA
20 provide reasonable assurance of safety for all three
21 thicknesses of the Intacts? What, if any, additional data
22 are needed to make this decision?"

23 Question 2. "Do the assessments of visual
24 symptoms provide reasonable assurance of safety for all
25 three thicknesses of the ICRS?"

1 Question 3. "Do the reports of corneal
2 sensation losses provide reasonable assurance of safety for
3 all three thicknesses? What, if any, additional data are
4 needed to make this decision?"

5 Question 4. "The range for the average
6 correction achieved with 0.25, 0.30, and 0.35 millimeter
7 Intacts is from minus 1.48 diopters to minus 2.76 diopters.
8 Does the achieved correction data support requested
9 indication for patient population with preoperative myopic
10 error ranging from minus 1.00 to minus 3.50 diopters with
11 1.00 or less of astigmatism?"

12 Question 5. "Do the safety and effectiveness
13 outcomes support approval of 0.25, 0.30, and 0.35
14 millimeter Intacts? Is distinct labeling warranted for any
15 one of the three proposed thicknesses?"

16 Question 6. "Is the current data in the
17 removal cohort sufficient to support reversibility claim?
18 If not, what is the minimum number of eyes and the minimum
19 length of follow-up that you recommend for this
20 assessment?"

21 Question 7. "Is the current data on exchange
22 procedures sufficient to support the claim for
23 adjustability of refractive effect? If not, what is the
24 minimum number of eyes and the minimum length of follow-up
25 that you recommend for this assessment?"

1 Question 8. "The sponsor would like to make a
2 claim of 'enhanced visual performance' in their labeling.
3 Do you feel that the data in this PMA support this claim?"

4 Thank you very much. This concludes my
5 presentation.

6 DR. McCULLEY: Thank you.

7 Do any of the three primary reviewers have
8 visual aids? Okay, then can we turn the lights back up
9 again.

10 If we could go to the primary reviewers, we're
11 going to take them in the order that they appear on the
12 agenda. For those of you who don't have an agenda in front
13 of you, that will be Dr. Sugar, then Dr. Van Meter, then
14 Dr. Grimmett.

15 Dr. Sugar?

16 DR. SUGAR: Thank you.

17 The sponsors are to be complimented for their
18 organized and complete presentation and their responses to
19 the FDA requests, and, likewise, Dr. Eydelman is to be
20 complimented for her complete and very insightful review.
21 I have to apologize for the redundancies in this process,
22 so I'm going to review some information that has already
23 been reviewed.

24 Four hundred and fifty-two patients were
25 enrolled in the initial portion of the study, first eye

1 portion, 449 of which underwent implant procedures.

2 Ultimately, many second eyes were operated upon, leaving a
3 total of 735 eyes undergoing implants. Most of the data
4 was reviewed, however, for just the first 449 eyes.

5 Demographics were consistent with other refractive studies,
6 and the accountability was very high with 410 of 420, or
7 97.6 percent, available for study at twelve months. The
8 other 29 not eligible included 20 who had the implants
9 removed, 4 exchanged, 1 who had a single segment removed,
10 and 4 who had analyses out of window.

11 The efficacy exceeded guidelines for all
12 measures, although as Dr. Schanzlin mentioned, it should be
13 kept in mind that the guidelines are for correction of
14 myopia up to minus 7.00 diopters and this study attempted
15 to correct only up to 3.00 diopters of myopia. At month
16 12, for the 0.25 millimeter segments 83.7 percent had 20/20
17 or better acuity, and 99.3 percent 20/40 or better
18 uncorrected visual acuity; for the 0.30 millimeter segments
19 this was 77.5 percent and 97.1 percent; and for the 0.35
20 millimeter segments 60.6 percent and 93.4 percent. And
21 with the ultimate recommended prescribing ranges, the
22 figures are even a little bit higher.

23 An attempted versus achieved scattergram was
24 never presented and data on specific patients with greater
25 than 3.00 diopters of myopia preoperatively were not

1 presented to me for analysis. This would be useful.
2 Predictability, however, as defined by plus or minus 1.00
3 diopter of manifest refractive spherical equivalent, was 92
4 percent at twelve months and 69 percent plus or minus 0.5
5 diopters and this decreased with increasing segment
6 thicknesses.

7 Stability appeared to be achieved between three
8 and six months, and this was acceptable.

9 The sponsors make an efficacy claim for
10 enhanced visual performance, which Malvina just reviewed,
11 or hyperacuity. The three definitions Malvina went
12 through. They conclude that "these results indicate that
13 the keraVision Intacts may be capable of increasing the
14 resolving power of the cornea to provide for enhanced
15 visual capabilities for some patients. Patients may be
16 able to achieve better visual acuity than they could
17 achieve with their prior methods of vision correction."
18 However, it should be noted that 13.9 percent lost one or
19 more lines of BSCVA while 19.5 percent gained one or more
20 at twelve months. Admittedly the ability to gain lines is
21 limited since many patients were 20/12.5 or 20/10
22 preoperatively. The percentage of contralateral,
23 unoperated eyes that also gained best corrected visual
24 acuity was not presented.

25 Given that examiners and subjects were not

1 masked as to treatment and there was no specific wording
2 that I'm aware of for testing so that examiners and
3 subjects would likely to have been more eager to have good
4 outcomes in the treated eye, this phenomenon may not be
5 real. The fact that a number of examiners apparently
6 commented on how well patients saw may reflect a difference
7 between this procedure and other refractive procedures
8 rather than truly enhanced performance since the visual
9 axis was not incised or ablated.

10 In terms of safety, this data has been reviewed
11 and safety appears to not be an issue in terms of contrast
12 sensitivity or loss of best spectacle-corrected visual
13 acuity.

14 Induced cylinder likewise decreased as time
15 went on, although the increased likelihood of induced
16 cylinder with increased ring thicknesses should be
17 mentioned in the package inserts.

18 Endothelial cell density and morphology were
19 assessed in a subgroup. Cell density overall decreased by
20 a mean of less than 1.5 percent which was not statistically
21 significantly different from the 0.6 percent cell loss per
22 year found by Bourne in normal patients. But 13 of 110
23 patients however did show a greater than 10 percent cell
24 loss at one year. While the mean cell loss is not
25 clinically significant, patients with a 10 percent cell

1 loss are of concern and follow up endothelial analyses on
2 these patients may be worthwhile. It is conceivable that
3 especially with the thicker implants that there is
4 continuing flection of the cornea and flection of the
5 posterior surface of the cornea that may induce continuing
6 endothelial cell loss. Perhaps a post-PMA analysis of
7 these patients and of future patients would be worthwhile.

8 Changes in pachymetry were not significant.
9 Channel deposits and lamellar channel haze and peripheral
10 vascularization do not appear to be significant. There was
11 central clouding noted at 1 and 2-plus level in two
12 patients and this rapidly resolved and appears not to be
13 significant.

14 The corneal sensation data is difficult to be
15 certain about. There were 25 of 259 patients, or 9.3
16 percent, who lost 20 millimeters or more of sensitivity at
17 six months while 13 of 237, or 5.5 percent, showed 20
18 millimeters or more loss at twelve months. Interestingly,
19 seven of these thirteen did not show the loss at six months
20 but did show it at twelve months, which is certainly hard
21 to explain. The changes however do not appear to be
22 practically important.

23 Visual fields were unchanged. Intraocular
24 pressures were not an issue.

25 Subjective symptoms are an issue. The severity

1 of symptoms was assessed in Phase II subjects while
2 severity and frequency were assessed in Phase III. At
3 twelve months in Phase III important changes were noted in
4 difficulty with night vision, blurry vision, halos, double
5 image, and glare. Double vision and fluctuating vision
6 increased with ring thicknesses and frequency of them
7 increased with ring thicknesses. This information must be
8 reflected in the physician and patient labeling.

9 And 90 percent of patients were somewhat or
10 strongly satisfied with outcomes at twelve months. Levels
11 were higher, 95 percent for patients undergoing bilateral
12 implants, but this increased percentage presumably reflects
13 the fact that patients who are satisfied with their first
14 eye are more likely to go on and have surgery on their
15 second eye.

16 Adverse events included one infection with
17 staph epidermidis which resolved with removal of the
18 segments. It is uncertain whether removal of the segments
19 was necessary or not. All patients with adverse events
20 recovered 20/20 or better best spectacle-corrected visual
21 acuity in the study.

22 Reversibility of effect, 36 of 449 patients, or
23 8 percent, underwent implant removals; 12.2 percent in
24 Phase II, and 7 percent in Phase III. Twenty-one removals
25 were because of symptoms of glare, fluctuating or blurred

1 vision, halos, or poor night vision; seven were due to
2 under or overcorrection; five due to induced astigmatism,
3 one due to infection; one due to perforation, and one
4 because of this mentioned patient who had his removed
5 because the FAA did not accept them.

6 No patient lost best spectacle-corrected visual
7 acuity or was worse than 20/20, and all patients with
8 available data were plus or minus 1.00 diopter of their
9 preoperative refraction, and 82 percent were plus or minus
10 a half diopter. Two patients however had severe trouble
11 with night vision, double images, or fluctuating vision.
12 All others had either no or mild symptoms. So the
13 procedure does appear to be reversible.

14 In terms of lens exchanges, in my review I
15 found that 12 eyes underwent exchanges, all for
16 undercorrection. Three of these were later removed because
17 of persistent undercorrection. All exchanges but one were
18 with segments of 0.40 millimeters or 0.45 millimeters,
19 sizes which were not part of this PMA.

20 Improvement in manifest refractive spherical
21 equivalent was noted in all patients but the magnitude
22 varied from only 0.12 diopters with a 0.35 to a 0.45
23 exchange to 1.50 diopters with a 0.35 to 0.40 exchange.
24 Patient numbers, data on the larger ring sizes, and
25 predictability of outcomes are insufficient to allow

1 approval of wording suggesting that undercorrection can be
2 improved through ICRS exchanges.

3 With these findings mentioned, I recommend
4 approval with conditions, and we will discuss this more
5 later. This ends my review.

6 DR. McCULLEY: Thank you.

7 Dr. Van Meter? Feel free to be as detailed or
8 as summarizing as you feel comfortable with.

9 DR. VAN METER: This is Woodford Van Meter.
10 Mr. Chairman, do you want us to go through the answers of
11 the questions proposed now or will we do that later?

12 DR. McCULLEY: No, we'll do that later. But if
13 you have something incorporated in your review that
14 addresses it, then do that. But we're not going to go
15 through those specific questions now.

16 DR. VAN METER: Woodford Van Meter. I also
17 would like to thank the sponsors for a thoroughly organized
18 presentation, and thank Dr. Eydelman for a very detailed
19 clinical review that made this reviewer's job much easier
20 than it otherwise might have been.

21 The accountability of this study I thought was
22 very reasonable with 452 patients initially enrolled, 449
23 eyes comprising the patient cohort, as previously
24 mentioned. Thirty-nine eyes did not have a month 12
25 examination, leaving 410 eyes which completed the month 12

1 exam as the basis for the data presented. Of the 39 eyes
2 missing the month 12 exam, there were 20 removals, 4
3 exchanges, 1 single segment removal, and only 5 patients
4 out of this cohort lost to follow-up. The data were
5 presented for 97.6 percent of the 410 implanted eyes. This
6 is very reasonable accountability.

7 In the effectiveness category, there were three
8 efficacy endpoints reviewed, which were uncorrected acuity,
9 predictability, and stability. At twelve months,
10 uncorrected visual acuity showed 74 percent of subjects
11 20/20 or better, and as might be expected, the findings
12 were slightly better with the lower thickness implants than
13 the higher thickness implants. The patient numbers were
14 comparable between the three sizes. Overall, 15.2 percent
15 of patients saw 20/40 or better preoperatively, but only
16 13.8 percent were in the lower myopia group that received
17 the 0.25 millimeter rings.

18 I think it is significant that FDA guidance
19 calls for 85 percent uncorrected visual acuity of 20/40 or
20 better, but we're comparing apples and oranges because the
21 guidance document was really issued for lasers which
22 correct up to 7.00 diopters and in this particular device
23 correction only up to 3.50 diopters is reasonable.
24 However, the data did show a significant improvement in
25 uncorrected visual acuity compared to preop in all three

1 groups.

2 The predictability analyzed for manifest and
3 cycloplegic spherical equivalent again showed better
4 predictability at twelve months with a 0.25 millimeter
5 height than the other two. The revised prescribing
6 indications when stratified by ring height showed that the
7 manifest refraction plus or minus 1.00 diopter was achieved
8 in 95.5 percent with a 0.25 height, 85 percent with a 0.35
9 height. Again, this reaches FDA guidance document
10 parameters with the proviso already mentioned.

11 I would like to make a comment however that the
12 reversibility of this implant device probably makes
13 additional risk acceptable and these devices do appear to
14 be reversible.

15 As might be expected, the predictability is
16 slightly better in the lower myopia range. And as was
17 mentioned in the revised document that we received in late
18 December, the fact that the greater scleral elasticity in
19 younger patients probably accounts for the improved effect
20 in younger patients rather than older patients as was
21 initially suspected.

22 Certain subjects appeared to achieve an
23 enhanced level of vision following the implants -- 21
24 percent achieved a post-operative uncorrected acuity that
25 was better than preoperative best spectacle-corrected

1 acuity, 19.5 percent of patients achieved a post-operative
2 best corrected acuity better than preop. There was
3 approximately 60 percent overlap between patients in these
4 two groups.

5 I believe that labeling is reasonable to
6 describe the enhanced visual capabilities observed in some
7 patients, but the lack of predictability for which patients
8 are likely to see better and the obvious problems from
9 fitting these patients that are undercorrected with any
10 contact lens to correct residual refractive error, which
11 might be your best option if they don't want to wear
12 glasses, make labeling a very critical issue that I think
13 we should discuss. I think it's fair for the company to
14 represent the improved effect, but I don't think any claims
15 can be made regarding this effect.

16 The sponsor states stability by manifest
17 refraction was obtained by month 3 because 97 percent of
18 eyes had MRSE at month 6 plus or minus 1.00 diopter of the
19 month 3 examination. I don't think stability is a major
20 concern. In the inclusion criteria patients had to have a
21 manifest refraction plus or minus a diopter within six
22 months prior to surgery. Therefore, stability after the
23 procedure is no worse than stability before the procedure.

24 Safety issues are discussed in terms of best
25 spectacle-corrected acuity, induced cylinder, contrast

1 sensitivity, and endothelial cell counts. With best
2 spectacle-corrected acuity, there was very minimal risk of
3 losing best spectacle-corrected acuity. At month 12, three
4 patients had a BSCVA worse than 20/20. All patients had
5 BSCVA 20/20 or better at 18 and 24 months. The number of
6 subjects that lost BSCVA was comparable to the number of
7 patients that gained BSCVA, suggesting a normal statistical
8 distribution in effect. At six months, seven subjects lost
9 two or more lines of BSCVA, none had BSCVA worse than
10 20/32, and nine had improved BSCVA. At twelve months, four
11 patients had worse BSCVA but none worse than 20/25, and six
12 actually had improved BSCVA.

13 For induced cylinder, the percentage of
14 patients that had an induced cylinder greater than 1.00
15 diopter at twelve months was 3 percent with a 0.25 implant,
16 7.3 percent with a 0.30, and 11.7 percent with a 0.35.
17 Again a slant that would not be unexpected. Since no
18 subject had greater than 2.00 diopters of cylinder and
19 since the procedure is potentially reversible, I believe
20 this complication can probably be accepted as a learning
21 phenomenon of implantation surgery. And it should be
22 covered in labeling that the possibility of induced
23 cylinder with the higher ring thicknesses, i.e., 0.35
24 millimeters, probably is worth including in the
25 preoperative consent.

1 Endothelial specular cell counts. The slightly
2 increased endothelial cell loss at twelve months in the
3 0.35 millimeter group is of some concern. The percentage
4 of patients with greater than 10 percent loss at month 12
5 was 2 percent with a 0.25 implant, 3 percent with a 0.30,
6 and 15 percent with a 0.35 implant. However, given the
7 limitations of endothelial cell counts, the variations in
8 cell counts in even normal individuals, and the lack of
9 concrete data on what peripheral cell loss means in normal
10 subjects -- which to my mind makes me wonder why we collect
11 peripheral cell loss because we're not certain what we want
12 to do with this data -- in addition to the known
13 fluctuations in cell counts in contact lens wearers, the
14 seriousness of the problem, or whether it is a problem at
15 all, is not known.

16 I believe that additional follow up should be
17 collected on the 0.35 millimeter subjects to gather more
18 data to determine what the potential problem might be down
19 the road. But I note that even at the 10 o'clock cell
20 count, which is where most of the problem seems to have
21 appeared, is taken through the edge of the implant and it
22 is reasonable to suspect that there is some distortion in
23 the measurement here. So I think postmarket follow-up for
24 the 0.35 implants with endothelial cell counts is probably
25 reasonable.

1 Contrast sensitivity testing showed minimal
2 adverse effect of the implants on glare and low contrast
3 acuity. I thought it was acceptable.

4 There were five adverse events reported in the
5 four hundred and fifty-four patients. One patient with
6 infectious keratitis, one patient had shallow placement of
7 the stints, two patients with anterior chamber perforations
8 -- I will not repeat these because they've been discussed.

9 The corneal sensation loss probably means that
10 labeling should include some notation that partial loss of
11 corneal sensation that may be temporary is possible. The
12 clinical significance however of this loss of corneal
13 sensitivity is not known.

14 Twenty-nine of the 31 removals in the cohort
15 were due to patient dissatisfaction with the device. There
16 have been 34 removals in 725 subject eyes implanted through
17 May of 1998 for a 4.7 percent removal rate. The sponsor
18 has presented sufficient data on reversibility, with an 87
19 percent return to manifest refraction of plus or minus a
20 half a diopter, to make reversibility a reasonable claim.

21 That concludes my presentation except for the
22 proposed questions and I will hold that until later.

23 DR. McCULLEY: Dr. Grimmett, you have a very
24 scholarly and outstanding review. If you could stress
25 areas where in more detail where you bring up new points

1 and where you're in agreement, if you could state that
2 you're in agreement, if it needs fine tuning, then state
3 the fine tuning.

4 DR. GRIMMETT: Yes, I refer everyone to my
5 detailed comments of 13 December. I will not belabor the
6 point, and to avoid redundancy I will just state some
7 highlights of issues that I felt important. But I will
8 refer everyone to the 13 December comments that I made. I
9 will try to summarize my comments into four main issues --
10 regarding the data in general, safety issues, efficacy
11 issues, and regarding labeling for the thicker implants.

12 Regarding the data in general, I echo the
13 comments of Drs. Sugar and Van Meter regarding
14 complimenting the sponsor for an extraordinarily
15 comprehensive data presentation. I appreciated it. I have
16 concerns that the follow-up percentage in month 18 and 24
17 are low overall and have the potential for sample bias. So
18 my comments will be restricted to month 12 data. I don't
19 believe that the data from those longer intervals may
20 represent the entire cohort for sample bias reasons.

21 Regarding safety issues, I will comment on five
22 highlights. First, regarding distance best spectacle-
23 corrected visual acuity loss without glare. It is
24 important to note by way of background that a loss of 2 or
25 more Snellan lines, that is ten letters, has been

1 recognized as being clinically meaningful because a gain or
2 loss of one line can occur from one exam to another in
3 unoperated eyes, particularly in individuals who can see
4 20/10, 20/12, or 20/16. Using this definition of greater
5 than two lines of loss as clinically meaningful, the
6 patients in this study had an acceptably low rate of loss;
7 that is, 4 out of 410, or 1 percent, loss greater than two
8 lines at 12 months.

9 The second safety issue regarding the corneal
10 endothelium, my major concern concerns the 0.35 millimeter
11 implant. There was a statistically higher rate of mean
12 cell loss, that is, approximately 5 percent, at the 10
13 o'clock position. I calculated instantaneous annual cell
14 loss rates and it is a footnote in the longer document,
15 footnote 9, and determined that at a 5 percent cell loss
16 rate, if we started with a cell count of 2,700 cells per
17 square millimeter, if that were continuous annually, in
18 order to reach 800 cells, which is the threshold for which
19 we may start seeing clinically significant corneal edema,
20 it would take approximately 23 and three-quarters years to
21 reach that. That would be concerning for younger patients
22 because that may occur during their lifetime. Hence, due
23 to that reason, I believe that longer term follow-up is
24 warranted for the 0.35 millimeter implant peripherally.
25 Perhaps this is best achieved by postmarket surveillance.

1 The third issue on safety is regarding corneal
2 sensation. It has been previously discussed that 13 out of
3 237, or 5 percent, had decreases greater than 20
4 millimeters in central corneal sensation. But, very
5 importantly, no patient with a central corneal sensation
6 loss had corneal staining of any kind. Dr. Lemp in his
7 slides indicated that all 13 of these patients returned to
8 normal, which is quite reassuring. However, I would state
9 that in my opinion appropriate labeling should include a
10 statement regarding the potential for altered corneal
11 sensation in select patients, if it's not already done.

12 Regarding the subjective assessment by
13 subjects, it has been pointed out that the frequency of all
14 visual symptoms increased at month 12 to a maximum of 17
15 percent for difficulty with night vision. And regarding
16 the magnitude of visual symptoms, up to approximately one
17 in five patients complained of moderate or severe visual
18 symptoms at month 12. These findings suggest that optical
19 quality has been altered in select patients.

20 I agree that labeling should reflect tables
21 that show the frequency data with "often" and "always" in
22 the categories and in the magnitude data "moderate" and
23 "severe." Simply reporting the "always" and "severe" in my
24 mind downplays important visual symptoms that occur in a
25 fair number of subjects within this study. Assuming that

1 labeling reflects these data, I think that the visual
2 symptoms are acceptable.

3 Regarding the reversibility of this data, I
4 appreciated the updated data from the manufacturer that
5 increased the follow-up percentage to 82.4 percent, 28 out
6 of 34 eyes, at three months following removal. That's up
7 from 62 percent follow-up in the original submission that I
8 reviewed. This essentially removes the specter of sample
9 bias. I believe that greater than 80 percent follow-up is
10 reasonable.

11 I believe that appropriate labeling regarding
12 reversibility should reflect the fact that the data is
13 limited to three months. We're only talking about 28 out
14 of 34 eyes that had the procedure removed. The long-term
15 stability of the refraction is not known. I have concerns
16 that three months is not long enough to fully establish the
17 stability of the refraction postremoval. And, number
18 three, unless data is available to the contrary, the
19 suitability of these eyes for further refractive procedures
20 is unknown. I did not locate data that indicates that
21 they're suitable for further work.

22 The third point regarding efficacy issues,
23 first regarding uncorrected visual acuity. Uncorrected
24 visual acuity has been previously shown to be an important
25 predictor of patient satisfaction, and, indeed, 97 percent

1 of these patients had 20/40 or better uncorrected visual
2 acuity at month 12. I thought that was excellent.

3 Regarding the stability of the manifest
4 refraction spherical equivalent, the majority of the
5 patients achieved, that is greater than 80 percent achieved
6 plus or minus a half diopter from visit to visit. The
7 procedure does appear to be stable.

8 I believe that the labeling should specifically
9 state the percentage of patients changing plus or minus a
10 half diopter from visit to visit since the level of
11 preoperative myopia is low. Hitting the mark 95 percent of
12 the time plus or minus 1.00 sounds very impressive until we
13 consider the fact that some of these patients went into
14 this procedure with myopia as low as minus three-quarters
15 of a diopter. A two diopter spread for those patients I
16 don't believe is acceptable. As long as the labeling
17 reflects what is plus or minus a half, I think that would
18 be sufficient.

19 Regarding the predictability, for the thinner
20 implants, 0.25 and 0.30, the procedure is capable of
21 achieving plus or minus a half of intended in approximately
22 three-quarters of patients at 12 months. Similar to my
23 comments on the stability, I believe the labeling should
24 indicate predictability within a half diopter of intended.
25 And I had a comment regarding that the mean achieved

1 correction should equal the nominally predicted correction
2 but I believe one of Dr. Schanzlin's slides altered the
3 prescribing range for the 0.25 millimeter implant. The
4 data reflect a mean achieved correction of approximately
5 1.5. In the original submission, the nominally predicted
6 correction was 1.3. I believe that should be reconciled,
7 but it may have already been done.

8 Regarding the hyperacuity phenomenon, I believe
9 in order to be consistent throughout the report, or any
10 report, that the definition of a clinically significant
11 gain or loss of Snellan vision should remain the same
12 within the report. I don't think it's fair to report loss
13 greater than two lines and report gains greater than one
14 letter, because in doing so you're minimizing the downside
15 and emphasizing the upside. I would go with the generally
16 accepted clinical definition of a two line cutoff.

17 Moreover, Dr. Sugar's point regarding the potential for
18 bias in the post-operative vision measurements is valid and
19 may, indeed, skew the data toward better post-operative
20 vision measurements. Hence, I endorse only reporting
21 changes in the measured visual acuity greater or equal to
22 two lines.

23 The sponsor has an argument that the data may
24 be skewed because it is difficult and/or impossible to
25 measure patients with a two line gain for those patients

1 starting with 20/10 or 20/12.5. That is a valid point.
2 However, I think the issue is moot. The number of eyes by
3 my review that started in that circumstance is only
4 approximately 20 eyes total. That's only 5 percent of the
5 eyes in the total study started with that circumstance. So
6 I don't believe that issue matters.

7 Regarding adjustability, I believe the number
8 of exchanges is too low for claims of adjustability.
9 Additionally, I believe the implant sizes utilized for many
10 of the exchanges are outside the range of those reported in
11 this study. Therefore, the data do not support claims of
12 adjustability.

13 My final issue is with regard to the labeling
14 of the 0.35 millimeter implant. I believe that the thicker
15 implant will require different labeling. As detailed in my
16 16-page document of 13 December, there are eight issues
17 that I believe should be specific to the 0.35 millimeter
18 implant.

19 One, the endothelium, the 0.35 millimeter
20 implant, as previously mentioned, had a 5 percent mean loss
21 over 12 months. That was a statistically significant
22 difference from the others. I believe that should be
23 mentioned.

24 Two, regarding induced cylinder greater than
25 1.00 diopter, there is an increasing trend with higher ring

1 thicknesses and that was a statistically significant
2 difference.

3 Three, regarding explants, there is a trend for
4 a higher rate of explantation due to refractive or optical
5 aberrations with the thicker implant. That is 19 of 34
6 eyes with 0.35 millimeter were explanted.

7 Four, with regard to the frequency of visual
8 symptoms, the thicker implant has an increasing trend for
9 double images and fluctuating distance and near vision.
10 That difference was statistically significant.

11 Five, regarding the magnitude of visual
12 symptoms, the thicker implant has increasing trend for
13 double images and fluctuating distance vision. That also
14 was statistically significant.

15 Six, regarding uncorrected visual acuity, the
16 thicker implant has lower success rates for achieving 20/20
17 or 20/40 uncorrected vision. That was statistically
18 significant.

19 Seven, regarding the stability of the manifest
20 refraction spherical equivalent, the thicker implant has
21 lower proportion of subjects plus or minus a half diopter
22 from visit to visit. That finding was statistically
23 significant.

24 Eight, finally, regarding predictability, the
25 thicker implant has a lower predictability for both plus or

1 minus a half of intended and plus or minus one diopter of
2 intended. That also was statistically significant.

3 I'll conclude my comments there and wait for
4 the questions. Thank you for your attention.

5 DR. McCULLEY: Thank you. That was excellent.
6 I thank all three of the reviewers for excellent and
7 thoughtful reviews. It made the job much easier for all
8 the rest of us.

9 At this point, what we'd like to do is recall
10 the sponsor to the table, and the purpose of this portion
11 of the proceedings is for panel to ask sponsor specific
12 questions relative to the PMA and for sponsor to respond
13 directly to the questions posed. It is not an opportunity
14 for the sponsor to introduce additional issues, unless
15 they're brought up by panel.

16 So if the sponsor would please return to the
17 table. Is sponsor ready?

18 MS. CROCKETT-BILLIG: Almost.

19 DR. McCULLEY: Let me know when you're ready.
20 You have your water bottle. Are you ready?

21 MS. CROCKETT-BILLIG: Sponsor's now ready.

22 DR. McCULLEY: We'll open the floor now for
23 panel to query the sponsor.

24 Dr. Van Meter?

25 DR. VAN METER: I'd like to ask two questions.

1 One, it has to do with the surgical procedure. When you
2 are doing an exchange, is it necessary to use the
3 instrument that separates the corneal stromal lamella in
4 the same way that you do a primary procedure or can you
5 pull one implant out and just place larger ones in place
6 without having to recarve the channel?

7 MS. CROCKETT-BILLIG: This is Darlene Crockett-
8 Billig from KeraVision. It is not necessary to redissect
9 the channel. The segments are simply removed and new
10 segments of a different size are reinserted into the
11 existing channel.

12 DR. VAN METER: Okay. Well, then the one
13 perforation that occurred on an exchange procedure, did
14 that involve the device actually being driven through
15 Descemet's membrane into the anterior chamber?

16 MS. CROCKETT-BILLIG: No. That adverse event
17 was related to an incorrect diamond knife setting. When
18 the initial incision was reopened, it was related to a
19 diamond knife setting.

20 DR. VAN METER: So the perforation occurred on
21 the initial incision reopening.

22 MS. CROCKETT-BILLIG: Yes, that's correct.

23 DR. VAN METER: Then my second question is
24 there was one perforation that I presume was understood to
25 actually dissecting the lamella, the one perforation that

1 occurred in the primary procedure. Was that a diamond
2 knife problem or a dissection problem?

3 MS. CROCKETT-BILLIG: I'll have our director of
4 clinical, Dan Beck, respond to that question.

5 MR. BECK: This is Dan Beck, director of
6 clinical for KeraVision. In reference to your question of
7 one of the two operative posterior corneal perms, yes, upon
8 review of the videotape in that particular case, we
9 observed that there was additional dissecting of the
10 incision that really deviated from the protocol, and we
11 attribute the dissection ultimately to that underlying
12 cause.

13 DR. VAN METER: So the incision was too deep?

14 MR. BECK: It may have been too deep. It may
15 have been incised too deep. That was not obtainable from
16 the video.

17 DR. VAN METER: I guess I'm trying to
18 conceptually understand how the perforation occurred, and
19 that was when you actually put in the device that creates
20 the channel.

21 MR. BECK: Now, it did not -- it's hard to tell
22 from the video when exactly the perforation did occur. It
23 may have been ultimately due to a too deep diamond knife
24 setting or too excessive manipulation of the incision prior
25 to insertion.

1 DR. VAN METER: Because one of the concerns I
2 think a surgeon would have is that the effect of this is
3 clearly dependent on how deep the device is inserted, and
4 it probably makes a difference that approximately two-
5 thirds thickness is a nice approximation, but whether you
6 want 200 microns or 300 microns or 400 microns between the
7 device and Bowman's probably makes a difference in the
8 effect. I presume this is something that you all have
9 looked at.

10 MS. CROCKETT-BILLIG: I'd like to call Dr.
11 David Schanzlin to the podium to address that.

12 DR. SCHANZLIN: We have looked at thickness
13 based on --

14 DR. McCULLEY: Identify yourself.

15 DR. SCHANZLIN: I'm sorry. Dave Schanzlin,
16 representing KeraVision.

17 We have looked over the last several years at
18 all of our results, looking at the effect of depth of the
19 implant versus effect, and there doesn't seem to be an
20 effect as far as efficacy. There is a concern, of course,
21 that you be deeper than 50 percent in these cases, and so
22 that is a concern for the surgeon who's doing the surgery.

23 The incident discussed here, where the surgeon
24 set the diamond knife at two revolutions -- in other words,
25 he was off by more than 300 microns -- is something that

1 certainly should be addressed in training surgeons how to
2 do that.

3 Then the other question you had was on the
4 exchange procedure. If an exchange is done early, within I
5 would say the first six to nine months, a Sinsky hook can
6 easily open the incision, and you can find the channel and
7 then grab it with the positioning hole with the Sinsky
8 blunt section up, and draw out the ring with ease. If it
9 is longer than that, you sometimes will have to reincise
10 the incision if it's scarred, and that's what happened in
11 this case where a diamond knife was used.

12 DR. VAN METER: Okay.

13 DR. McCULLEY: So your perforations were with
14 your incisions, not with the insertion of the segments.

15 DR. SCHANZLIN: True. The other case that was
16 presented here, the question was was it the incision or was
17 it the initial pocketing with the Saurez spreader? Was it
18 a deep incision to begin with and then the Saurez then
19 caused the dissection? It's hard to tell from the video
20 which of those two steps, but it was not from the device
21 itself.

22 DR. VAN METER: Thank you.

23 DR. McCULLEY: Dr. Wang?

24 DR. WANG: Ming Wang. I've been formulating
25 this question all morning after hearing the presentation

1 and all the reviewers' comments. Under this question, I
2 have three subquestions under it. The question is does the
3 0.35 millimeter group present as a distinct, correlatively
4 different group compared with the other two?

5 My question is supported by three comments.
6 Number one, there appears to be grossly doubling of visual
7 symptoms and frequencies in the .35 group compared with the
8 other two. For example, diplopia in the .3 group is 6.9
9 percent and the .35 group is 15.3 percent and endothelial
10 loss is higher. So therefore, there appears to be a
11 grossly doubling of problems with the .35 group.

12 The second subquestion, everybody mentioned
13 about possibly the necessity of squeeze, so to speak, of
14 the guidelines of the guidance document for a lower range
15 of correction. The reference has been for the 1 to 7
16 diopter correction, and the issuing question is 1 to 3
17 diopter. My question is does a squeeze in the guideline in
18 the guidance document for this lower range result in a
19 questioning of this .35 group quality?

20 For example, the guidance document asks for 75
21 percent within 1 diopter, though grossly overall, all three
22 groups included, you have 91 percent, which is far above
23 the 75 percent within the 1 diopter, but the .35 group, 83
24 percent, which is higher than 75 percent, but not that
25 high. Then 50 percent within .5 diopter final

1 predictability, though all three groups combined is 71
2 percent, which is much higher than the 50 percent, but the
3 .35 group again is 62 percent, which is marginally higher.

4 So the question again is if we need to squeeze
5 the guidance document for the low range, would that result
6 in the questioning of this .35 group?

7 Third, and the last subquestion in the .35
8 group, is is there any data -- I understand it's ongoing --
9 about the .4 and .45 groups? I understand it's not
10 included in the present review. Is there any data or
11 indication that would suggest there might be problems with
12 the thicker group, and because it resulted in more stress,
13 therefore, since it is a continuum from .35 to .4 to .45,
14 which would again result in a questioning of this .35
15 group, because it is the closest to the .4 and .45 groups.

16 Therefore, in summary, the question is does the
17 .35 group present as a distinct, correlatively different
18 group compared with the other two?

19 DR. McCULLEY: Sponsor?

20 MS. CROCKETT-BILLIG: I'm just taking notes on
21 the question. All right. To make sure I understood what
22 you said --

23 DR. McCULLEY: Can you speak more into the
24 mike, please?

25 And everyone, please, each time you speak,

1 please identify yourself.

2 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
3 KeraVision. To repeat the question as I understand it, you
4 asked basically does the 0.35 group present a distinct and
5 correlative difference as compared to the other two
6 thicknesses that we evaluated?

7 Then underneath that you had three
8 subcategories. The first category that you discussed was
9 having to do with the performance based on visual symptoms.
10 The second category had to do with basically the
11 performance of the .35 as compared to the other two
12 thicknesses. The third, you requested information on the
13 ongoing trials with the .40 and .45 thicknesses to see if
14 that was suggesting a safety problem.

15 One moment. We're bringing up the data.

16 As you recall from Dr. Schanzlin's
17 presentation, 29 percent of the subjects in the 0.35 group
18 were undercorrected. They were outside of the RPR range.
19 Here's a slide that was presented by Dr. Lemp providing the
20 visual symptoms at month 12. Now, as we go to look at by
21 ring thickness --

22 DR. McCULLEY: Dr. Rosenthal?

23 DR. ROSENTHAL: I frankly think that the issue
24 was well discussed by both the sponsor and the FDA and the
25 primary reviewers, and I'm not sure we're going to get much

1 more out of reviewing their data.

2 DR. McCULLEY: I think one issue that was --

3 DR. ROSENTHAL: This is more a comment about a
4 philosophical issue relating to the larger segment, and I
5 appreciate the comments made, but I'm not sure a response
6 is going to elucidate any greater clarity for the panel.

7 DR. McCULLEY: I think that there are some
8 philosophical points that have been brought up by the
9 reviewers that I agree I'm not sure that we'll resolve.
10 There was one point that was brought up that I think might
11 be of use, certainly, and that is do you have any insights
12 from the .40 and .45 that might help us assess the .35, and
13 whether there are trends that then continue in the negative
14 direction past the .35 that would cause concern?

15 DR. ROSENTHAL: May I? This is Dr. Rosenthal.
16 I think the PMA does not really address the two greater
17 segments and --

18 DR. McCULLEY: But Ralph, I think the PMA does
19 mention the .40 and the .45. Sponsor brought it up. I
20 think that information from that could be of value to the
21 panel in assessing the safety and effectiveness of the .35.

22 DR. PULIDO: I agree.

23 DR. McCULLEY: They may not have data. If they
24 don't have data, fine, but I think it's reasonable for us
25 to ask if they have information that can help us in our

1 deliberations. If they don't have the data, they don't
2 have it, but I think it's perfectly appropriate for us to
3 ask this, since they introduced it into the discussions and
4 it is reality.

5 DR. ROSENTHAL: This is Dr. Rosenthal. I
6 appreciate your need to want to know, but it's really not
7 part of the PMA and a formal presentation is not being
8 given to you. They can give you a general statement, but I
9 think it would be inappropriate to start bringing
10 formalized evaluation to us when no one has had a chance to
11 even look at it, including our own reviewer.

12 Certainly, I would think if you want to make a
13 general statement about it, please do.

14 DR. MCCULLEY: Could I clarify the point then
15 that it would be inappropriate for sponsor to bring or
16 present data, but if there are general comments that can be
17 made, that we would appreciate those general comments.

18 MS. CROCKETT-BILLIG: All right. Darlene
19 Crockett-Billig. On the data that we have available in our
20 ongoing Phase IIb and Phase IIIb trials with the 0.4 and
21 0.45 millimeter thicknesses, the endpoints are met, as
22 specified in our protocol for the .40 and .45 thicknesses.
23 When we look at within the recommended prescribing range,
24 we have 98 percent of patients 20/40 or better and 67
25 percent 20/20 or better.

1 DR. McCULLEY: I guess the points that Dr. Wang
2 brought up were more targeted toward symptoms, endothelial
3 loss, and the like.

4 MS. CROCKETT-BILLIG: We're collecting that
5 data. We don't have that available at this point in time
6 on the thicker ring sizes.

7 DR. McCULLEY: Fair enough.

8 Dr. Macsai?

9 DR. MACSAI: This is Dr. Macsai. I have a
10 series of questions for the sponsor if the chair will
11 allow.

12 The sponsor has talked about the comparison of
13 post-op vision with preop best spectacle-corrected vision,
14 and I am aware that Dr. Eydelman brought up in her review
15 the question of preop contact lens vision. Does the
16 sponsor have comparison of post-op uncorrected visual
17 acuity with preop contact lens visual acuity, either soft
18 or rigid?

19 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
20 Your question is does the sponsor have available
21 uncorrected visual acuity data comparing outcomes with
22 preop contact lens wearers versus post-op contact lens
23 wearers?

24 DR. MACSAI: Yes, that is my question. In
25 other words, did you measure best contact lens-corrected

1 visual acuity preop?

2 MS. CROCKETT-BILLIG: No, we did not.

3 DR. MACSAI: So then you don't have that
4 comparison?

5 MS. CROCKETT-BILLIG: The protocol specified
6 for best spectacle-corrected visual acuity.

7 DR. MACSAI: My next question is in regard to
8 the incidence or rate of fluctuating distance vision,
9 double images, and halos. With halos, if you look at the
10 often and always percent of patients with these symptoms
11 post-operatively, there is 11.8 percent patients with post-
12 op complaint of halos, compared to .3 percent preop, a 37
13 times increase. Yet if you look at the same symptom under
14 mesopic conditions with pupil size less than 7 millimeters
15 or greater than 7 millimeters, there appears to be no
16 statistically significant difference in the patients'
17 perception of halos. Therefore, I ask you if these halos
18 are directly due to this device.

19 The same thing can be seen for double images,
20 where the percent noted by patients is 6.7 percent
21 often/always post-operatively, zero percent preoperatively,
22 and approximately equal under mesopic conditions, and for
23 fluctuating distance vision the increase is from 1.1
24 percent preop to 3.3 percent post-op, or three times, where
25 it is approximately equal with the pupil being 7

1 millimeters, greater or less than 7 millimeters. So is it
2 the device that's causing these visual symptoms?

3 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
4 KeraVision. We did our analysis on visual symptoms --

5 DR. MACSAI: Well, I saw your data when you
6 presented it, but it doesn't answer the question.

7 MS. CROCKETT-BILLIG: We found that the
8 statistical significance was primarily related to post-
9 operative deviation from plano.

10 DR. MACSAI: So it's the post-operative
11 deviation from plano that's causing the halos, double
12 images, and fluctuating vision? Because these things don't
13 seem to change dependent upon pupil size, which you would
14 expect them to.

15 MS. CROCKETT-BILLIG: I'd like Dr. Michael Lemp
16 to address that question, please.

17 DR. LEMP: Thank you. Michael Lemp,
18 representing KeraVision.

19 Dr. Macsai, in response to your question, what
20 we found was the thing that was statistically associated
21 with these symptoms was deviation from plano and induced
22 cylinder, and you correctly point out that in some of these
23 symptoms there did not appear to be a statistical
24 association between pupil size and this.

25 If you have residual cylinder, you would expect

1 to have some of these symptoms irrespective of pupil size,
2 and so the induced cylinder was statistically associated
3 with these symptoms, in addition to deviation from plano,
4 so that the pupil size probably may play a role under
5 certain conditions, but the farther away you get from the
6 refractive effects, you may have these symptoms induced
7 without the influence of the pupil size.

8 DR. MACSAI: So Dr. Lemp, are these symptoms
9 due to induced cylinder or uncorrected cylinder?

10 DR. LEMP: Well, we think that the statistical
11 association that the analysis of the data showed was an
12 association between undercorrection and induced cylinder,
13 so both of those show an association.

14 DR. MACSAI: Dr. Lemp, could the induced
15 cylinder be corrected with spectacles? Because it seems
16 that contact lens fitting when this device is in your
17 cornea might be a challenge.

18 DR. LEMP: We don't have data on the contact
19 lens fitting aspect of that, but it's the fact that even
20 the patients with induced cylinder had very good visual
21 acuities with this and indicated that their acuity as such,
22 uncorrected, was not significantly affected, even though
23 it's probable that some of the visual symptoms that they
24 complained of were related to this. I don't have specific
25 data on the rest of your question.

1 DR. MACSAI: Thank you.

2 DR. McCULLEY: Dr. Macsai, how many questions
3 do you have?

4 DR. MACSAI: I have three more.

5 DR. McCULLEY: Okay.

6 DR. MACSAI: Is that okay?

7 DR. McCULLEY: For now.

8 DR. MACSAI: The next question is how does the
9 Intacts affect the physician's ability to examine the
10 peripheral retina of these patients, what sort of
11 distortion is there, and can these patients have a
12 gonioscopic examination? Can you put a three-mirror on
13 these corneas?

14 MS. CROCKETT-BILLIG: I'd like to ask Dr. David
15 Schanzlin to address that question.

16 DR. SCHANZLIN: Yes, the gonioscopy is not
17 difficult to do. You do notice with the thicker rings, .35
18 and the other ones that we aren't discussing today, that
19 there is a slight indentation at the 7 to 8 millimeter
20 optical zone. You'll see a little indentation on the
21 endothelium there, but it does not preclude visualization
22 of the angle in these cases. Examination of the retinal
23 periphery is also not obscured by this.

24 This is at a 7 millimeter optical zone. When
25 the pupil is very widely dilated, you get an interesting

1 shadow when you're doing indirect ophthalmoscopy, but not
2 something that you can't work around, so I don't think this
3 is any significant problem for a retinal surgeon, if that's
4 your question.

5 DR. MACSAI: Is the shadow eliminated when you
6 use a contact lens?

7 DR. SCHANZLIN: Yes.

8 DR. MACSAI: There were three patients or three
9 eyes that had the Intacts removed due to continuing
10 undercorrection in the exchange data. Did those three eyes
11 go on to have other refractive procedures?

12 MS. CROCKETT-BILLIG: This is Darlene Crockett-
13 Billig from KeraVision. We don't have that information
14 immediately available in terms of what the three explants
15 or the three subjects who had exchange procedures and went
16 on to have subsequent removals, if they had a second
17 refractive procedure.

18 Dr. David Schanzlin would like to comment.

19 DR. SCHANZLIN: Your question was specific, but
20 I think before someone else brought up the issue about
21 other refractive procedures, and patients who have had
22 removals have gone on to have PRK, patients have gone on to
23 have LASIK procedures without complication, and the cases
24 that at least we know about have good visual outcomes.

25 DR. MACSAI: The last question I have is simply

1 one that I found a bit confusing. Why did the sponsor
2 change the definition of the adverse event from Phase II to
3 Phase III regarding visual results? It made it very
4 difficult for review of that.

5 MS. CROCKETT-BILLIG: This is Darlene Crockett-
6 Billig, KeraVision. The sponsor did not change adverse
7 event definition between Phase II and Phase III for visual
8 symptoms or visual results. What we did is in Phase II the
9 protocol defined a removal or an exchange procedure as an
10 adverse event.

11 DR. MACSAI: Right.

12 MS. CROCKETT-BILLIG: And when we prepared the
13 Phase III protocol, those were additional procedures that
14 you would have, but they were not considered a safety-
15 related adverse event, and that was the reason for the
16 difference. That is our approved Phase III protocol that
17 the trial was conducted under.

18 DR. MACSAI: Thank you.

19 DR. McCULLEY: Are you through?

20 DR. MACSAI: Yes.

21 DR. McCULLEY: That was very good, Dr. Macsai.
22 Dr. Middleton?

23 DR. MIDDLETON: Renee Middleton. I'm
24 representing the consumer component.

25 Just really one follow-up, and then a specific

1 question. You are conceding then to the fact that you
2 really don't have firm data with respect to your claim
3 about suitability for additional refractive surgical
4 procedures? That you don't have any firm data on that?

5 MS. CROCKETT-BILLIG: This is Darlene Crockett-
6 Billig. I would like to call Dr. David Schanzlin to
7 address that question.

8 DR. MIDDLETON: He mentioned, just in his
9 comment there --

10 DR. McCULLEY: Did you have data submitted to
11 the FDA?

12 MS. CROCKETT-BILLIG: No. We do not have data.
13 It was not part of our protocol.

14 DR. McCULLEY: Then, as I understand it, you
15 may not introduce new data at this point.

16 MS. CROCKETT-BILLIG: We do know that our
17 patients have gone on to have other refractive surgeries.
18 We do have that information. We can tell you what
19 procedures they had, but we did not collect this data.
20 This is not our information.

21 DR. MIDDLETON: And your comment was just for
22 those that you're aware of, so there may be others that
23 you're not aware of, but for the ones that you're aware of
24 the procedures were --

25 MS. CROCKETT-BILLIG: Yes. Dr. Schanzlin will

1 elaborate on that point.

2 DR. McCULLEY: Based on data submitted to the
3 FDA.

4 DR. SCHANZLIN: Well, we don't have data
5 submitted. We have data here, but we don't have data
6 submitted.

7 DR. ROSENTHAL: This is Dr. Rosenthal. May I
8 say that I've written down the issue relating to the
9 labeling regarding refractive suitability for future
10 procedures, and we will certainly address that with the
11 sponsor.

12 DR. McCULLEY: But you've received no data.

13 DR. ROSENTHAL: We've received no data.

14 DR. McCULLEY: So as I understand it, then it
15 cannot be introduced at this time.

16 DR. ROSENTHAL: Right, but we have flagged it
17 and it will be discussed with the sponsor.

18 DR. McCULLEY: And I think the panel has
19 expressed its concern, and we'll probably come back to that
20 again when we make our final recommendations.

21 DR. ROSENTHAL: All right. Thank you.

22 DR. VAN METER: Woody Van Meter. Mr. Chairman,
23 would it be fair to ask Dr. Schanzlin to answer that
24 question, though, that he was prepared to do?

25 DR. McCULLEY: If it involves introducing data

1 not submitted to the FDA, as I understand it, it cannot be
2 submitted.

3 DR. ROSENTHAL: This is Dr. Rosenthal. I
4 think, if I heard Dr. Schanzlin correctly, he did say that
5 anecdotally he had reports that some refractive procedures
6 went well, so I'm not sure you really --

7 DR. McCULLEY: You know, we'd all like to know,
8 but I think this is not appropriate for these proceedings.

9 Dr. Middleton, did you have another question?

10 DR. MIDDLETON: Just a second question for my
11 understanding. You had six eyes, I guess, that underwent
12 exchange procedures because of the undercorrection, and I
13 was just interested in knowing what factors play a role in
14 undercorrection. Why might that occur?

15 MS. CROCKETT-BILLIG: Dr. Schanzlin, would you
16 please address that question?

17 DR. SCHANZLIN: Well, recall in this study
18 there were two phases, Phase II and Phase III studies, and
19 if you remember the graph that I showed of the refractions
20 of patients going into the study, clearly some of the
21 patients with the thicker rings were beyond the performance
22 of the product. So that automatically would bias them
23 toward an undercorrection. The risk factors for
24 undercorrection are primarily --

25 DR. MIDDLETON: In terms of the patient?

1 DR. SCHANZLIN: Yes. That, of course, should
2 be better following the RPR, the recommended prescribing
3 guidelines that we've presented today.

4 DR. MIDDLETON: How many of those were eye
5 contacts or contact lenses?

6 DR. SCHANZLIN: How many of the patients --

7 DR. MIDDLETON: Did any of those wear contact
8 lenses?

9 DR. SCHANZLIN: Of the ones who were
10 undercorrected?

11 DR. MIDDLETON: Yes.

12 DR. SCHANZLIN: You mean preoperatively, before
13 surgeries?

14 DR. MIDDLETON: Yes.

15 DR. SCHANZLIN: Can you help me with the actual
16 number of patients that we submitted?

17 MS. CROCKETT-BILLIG: We have that information.
18 We'd have to look it up. We'll get back to you on that
19 question.

20 DR. McCULLEY: Dr. Pulido?

21 DR. PULIDO: Jose Pulido. Just a quick
22 question. Apparently, there was a difference in the number
23 of the percentage of explants that were done at the
24 different centers. Emory and Northwest Corneal Services
25 had about 15 percent explantation rate versus almost half,

1 7 percent, explantation rates elsewhere.

2 Is there a big learning curve on this? Is that
3 something that needs to be put in the labeling?

4 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
5 KeraVision. When we analyzed the data, we looked at that
6 as a factor, and basically I think it comes down to patient
7 selection and making sure that you're selecting patients
8 that are within the correct range. Our data did not
9 indicate that there was an effect with learning curve per
10 se.

11 DR. PULIDO: Jose Pulido again. I don't
12 understand what you mean by patient selection.
13 Conceivably, you had strict criteria that patients would be
14 selected similarly at all the different sites.

15 MS. CROCKETT-BILLIG: Similarly has to do with
16 selection of the patients in terms of are they going to be
17 a suitable refractive surgery patient, period, or not.

18 Dr. Schanzlin, would you comment on that
19 further?

20 DR. McCULLEY: No one-upmanship, Dr. Schanzlin.
21 (Laughter.)

22 DR. SCHANZLIN: Can you say the question again
23 so I reformat the answer for you?

24 DR. PULIDO: Jose Pulido. My concern was the
25 basically doubling of explantation rates at two centers

1 versus the other centers, and why should there be such a
2 marked difference in explantation rates. Is there a
3 problem with learning curves that needs to be put in the
4 labeling? Why is there such a major difference between all
5 the different centers?

6 DR. SCHANZLIN: I don't see it as much a
7 problem with the surgeon having a worse result, but rather
8 a lower threshold for having a patient who is not pleased
9 with the result and just saying, well, we can take it out.

10 One of the unique things of the ring is that we
11 can remove it, and so when patients have some of the
12 symptoms that you've heard about today, maybe double vision
13 at night or whatever, one of the first questions we ask
14 clinically is, well, would you like us to take the ring
15 out? It's really a nice decision point question for the
16 patients, because very few of them, as you know, went on to
17 have rings removed. The data we have today is really the
18 data from those who selected to keep the rings in, with the
19 exception of the removal rate cases there.

20 So in the cases that you're talking about
21 there, and I don't want to go into individual surgeons,
22 because that's really not what we should do in a study like
23 this, but I think the surgeon involved in at least one of
24 those two centers has a very high rate, and yet other
25 surgeons at the same center involved in the study do not

1 have such a high rate. So it isn't that he had worst
2 results. I think he had a lower threshold for removals.

3 But addressing the bigger concern is what is
4 the learning curve on this? It doesn't take long to teach
5 a surgeon to do this procedure. In fact, in the study, and
6 I won't go into names again, one of the surgeons who came
7 in in the latest series in the Phase III studies has data
8 that far excels my cases, and that's perhaps more than
9 statistical variation, or maybe it just is, but I don't
10 think that that's a major concern based on the way you
11 addressed the question.

12 DR. McCULLEY: Dr. Schanzlin, you used "I
13 think" several times. Have you analyzed your data to
14 demonstrate that in fact what you said you thought is in
15 fact the case, that it's surgeon threshold, rather than any
16 other factor? If you have not, I think that would be
17 useful for the FDA to have.

18 DR. SCHANZLIN: Yes, we have. Do we have that
19 to present or --

20 DR. McCULLEY: Well, you don't have to present
21 it. You can just say that --

22 MS. CROCKETT-BILLIG: Yes, we have looked at
23 it.

24 DR. McCULLEY: And that is indeed the case,
25 that it's surgeon threshold that seems to be the issue in

1 terms of the excessive removal rate at two centers, not
2 anything else?

3 MS. CROCKETT-BILLIG: That's correct.

4 DR. McCULLEY: Dr. Higginbotham?

5 DR. HIGGINBOTHAM: Dr. Eve Higginbotham,
6 University of Maryland. Dr. Schanzlin, please stand.
7 Thank you.

8 (Laughter.)

9 DR. HIGGINBOTHAM: Simon says.

10 I would like clarification of two points that
11 are clinically related. On page 166 of your document,
12 Table 64, it's apparent that 54 of your procedures of 449
13 procedures, or 12 percent, took longer than 20 minutes to
14 perform, and I think the range went up to 44, as indicated
15 on Table 63.

16 What part of that surgical time, would you say,
17 was related to an increase in pressure when the suction cup
18 was applied? Is it the majority, most of that surgical
19 time, or is it still about a minute and 15 seconds, as you
20 stated earlier?

21 DR. SCHANZLIN: Well, perhaps the team here can
22 draw out the suction time data or give me a range on that.

23 Certainly, the first surgeries you do have
24 longer suction times. In training surgeons elsewhere in
25 the world, I can tell you that early surgeons tend to be

1 one and a half to three minutes maximum on the suction
2 time, and I'm usually down to around a minute to a minute
3 and a few seconds, and perhaps there are surgeons in the
4 room that are faster than I am.

5 But that certainly is one of the things that is
6 addressed in training, is the need to be cautious of that,
7 and also to be mindful that we do raise the intraocular
8 pressure during that process. The intraocular pressure
9 rise, though, is not as high as with LASIK and other
10 procedures. We don't need to raise as much vacuum on the
11 eye, because we only need to resolve the torque of the
12 dissection. We don't need to create a firm eye like the
13 microkeratomes do.

14 DR. HIGGINBOTHAM: Thank you.

15 I have one more question, just a very quick
16 question, and forgive me if this is precisely indicated in
17 your document, but for your post-operative visual field
18 assessment, was that a 120 screening program, as it was in
19 the preoperative assessment?

20 DR. SCHANZLIN: I believe it was.

21 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
22 That is correct.

23 DR. HIGGINBOTHAM: So none of these patients
24 actually underwent any full threshold or SITA testing? Is
25 that correct?

1 DR. SCHANZLIN: That's correct.

2 DR. HIGGINBOTHAM: Okay. Thank you.

3 DR. McCULLEY: Just one point of clarification.
4 You're saying that you raise the pressure with your suction
5 range to approximately 80 millimeters of mercury?

6 DR. SCHANZLIN: It's approximately 80, yes.

7 DR. McCULLEY: That's approximately what it is
8 in LASIK as well.

9 DR. SCHANZLIN: Well, I don't want to argue the
10 LASIK data, but I believe it's much higher than that in
11 LASIK.

12 DR. McCULLEY: Well, measuring with
13 pneumotonometry with capability up to 100, I would argue
14 with you, but we won't argue that. I just think we need to
15 leave that maybe then out of the discussion of the
16 comparison without solid data on either side.

17 Dr. Bandeen-Roche?

18 DR. BANDEEN-ROCHE: Karen Bandeen-Roche. The
19 first part of my question I think also will go to Dr.
20 Schanzlin, and it's a follow-up to Dr. Pulido's question, a
21 two-part question.

22 So the first part is that I also noticed
23 variation by site, and not only with explants, but in a
24 fair variety of the outcomes of visual acuity and
25 predictability, et cetera, so given the range of outcomes

1 in which there was variability by site, does your answer
2 still stand that there is good data to support that that's
3 selection, rather than learning curve?

4 DR. SCHANZLIN: I think so. Perhaps this is
5 really a statistical question. Are the sites poolable?
6 And I think that there is no data that suggests that we
7 could exclude any one site, because they really are not
8 that different. It's only when you look at the results.
9 That's why that statistically -- you know, I think you have
10 to consider the group as a whole, and not worry whether one
11 surgeon that has done 200 of these has a little bit better
12 data than somebody who's done just a few. I don't think
13 that there's a difference to that degree.

14 DR. BANDEEN-ROCHE: I would comment that
15 poolability analyses do have fairly low power often, and so
16 really large differences would be needed to fail that
17 analysis.

18 The second part of my question I think is
19 related, and it's just an analysis question. Were any of
20 the analyses done to determine P values or confidence
21 intervals that took into account the clustering within
22 sites or within positions, the effect of which would be to
23 decrease precision?

24 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
25 KeraVision. I'd like one of our statisticians to answer

1 that question, please. Don Young will answer that question
2 for us.

3 MR. YOUNG: Don Young, San Jose State
4 University. I am consultant to KeraVision. I have no
5 equity interest.

6 With regard to your question on whether or not
7 there was any clustering correction or adjustment, no,
8 there wasn't. We did not do any of those techniques.

9 DR. McCULLEY: I think Dr. Sugar was the person
10 I was aware of next.

11 DR. SUGAR: Two unrelated questions. First,
12 concerning specular microscopy, in the analysis of the
13 patients with the greater cell losses, was there a sense in
14 looking at the photographs, videotapes, or whatever you had
15 that there was cell loss that was emanating from the point
16 of flexure of the posterior cornea, Descemet's membrane,
17 and would this suggest the possibility of ongoing cell loss
18 in the deeper implants or the thicker implants?

19 I guess Dr. Edelhauser, or I'm not sure who
20 would be the person to answer.

21 MS. CROCKETT-BILLIG: I'd like to call Dr.
22 Edelhauser, please, to address that question.

23 DR. EDELHAUSER: Dr. Edelhauser, Emory
24 University, consultant and reading center for specular
25 microscopy for KeraVision.

1 The specular that we performed was with the
2 non-contact Robo, and indeed, getting the peripheral 10
3 o'clock reading was probably the most difficult reading to
4 get, because the image sometimes would be adjacent to the
5 ring, and most of the speculars that we analyzed had
6 between 70 and 150 cells, but indeed you could see in the
7 image that came up on the Robo that in that peripheral
8 region there could be a dark area where some of the cells
9 could be missing that we could not read. So the chances of
10 variation in the 10 o'clock region would be much greater
11 than, say, central.

12 DR. SUGAR: Variation in accuracy of your
13 capturing data or variation in cell density?

14 DR. EDELHAUSER: Well, we were masked when we
15 read this. I would say that our accuracy, when we analyzed
16 it, we looked at it in terms of -- well, I didn't do the
17 analysis. We sent the data there, but when we analyzed on
18 the Robo, you touched the center of each one of these
19 cells, and to get a maximum of 70 to 150 cells, so I would
20 say what cells showed up were as accurate as the instrument
21 could provide.

22 DR. SUGAR: And did you have a sense, even with
23 the later analyses, that there were cells that were
24 different in that area?

25 DR. EDELHAUSER: I would say no. The cells

1 looked the same in terms of what we saw on the Robo. The
2 percent hexagons were the same and the coefficient of
3 variation were the same.

4 DR. SUGAR: Thank you.

5 DR. McCULLEY: Were there more dark spots that
6 you couldn't interpret?

7 DR. EDELHAUSER: Well, the dark spots were
8 particularly off in the periphery, whereas in the central
9 reading, because this instrument will seek the apex of the
10 cornea, you were always right on the money at all times.

11 DR. McCULLEY: Any other questions for Dr.
12 Edelhauser while he's at the podium?

13 (No response.)

14 DR. McCULLEY: Thank you.

15 Dr. Sugar, you had a second question?

16 MS. CROCKETT-BILLIG: Dr. Holladay wanted to
17 make a comment.

18 DR. HOLLADAY: Jack Holladay, consultant to
19 KeraVision. The one comment I wanted to make about that,
20 we used the Robo also, and what you have to realize is when
21 it takes a measurement, the endothelial surface has to be
22 almost flat over the area that the reflection comes from in
23 order for you to be able to see all the cells in focus,
24 because it's in one focal plane, but what happens is when
25 the plastic is in the cornea, that indentation of the back

1 surface where there are endothelial cells makes it
2 impossible for that particular endothelial camera to get
3 that entire plane in focus.

4 So what he's talking about is that you can see
5 one portion of that plane and those cells look the same,
6 and then where it begins to bend off it gets dark and you
7 can't see anything. So the only indication we have is that
8 cells in the area that you can see are normal.

9 The second aspect is that the 10 o'clock that
10 you measure preop, post-op has to be right inside the ring,
11 so the chances of hitting the same point pre and post-op
12 are very difficult to do, and that's why I think you're
13 seeing this variation. It doesn't have anything to do with
14 the cells.

15 DR. SUGAR: Thank you.

16 The other question I have is about surgeon
17 education. You presented in your package what the program
18 would be, and I'm just trying to verify that surgeons would
19 be required before receiving the product to take the
20 course.

21 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
22 Yes, that is correct.

23 DR. SUGAR: Thank you.

24 DR. McCULLEY: Dr. Bullimore?

25 DR. BULLIMORE: Yes. This is Mark Bullimore.

1 I've heard a lot of data presented on visual acuity and a
2 lot of different sort of slices and dices done. My
3 question is this. When you compare the best spectacle-
4 corrected visual acuity pre and post, and you had a table
5 in your presentation, and you can probably find it quicker
6 than me, but what I want to know is what is the mean change
7 in best spectacle-corrected visual acuity in terms of
8 number of letters? You alluded that there's a change in
9 the distribution, but when you've quantized the data in
10 terms of lines of visual acuity, it's not immediately
11 obvious to me what that might be. So that's one thing I'd
12 like to know.

13 The other thing is, in terms of the induced
14 cylinder, I can't remember what type of astigmatism is
15 being induced by the procedure. Are we introducing a with
16 the rule or against the rule astigmatism or have we got a
17 shotgun effect? Is the sort of astigmatism usually this
18 way or this way? I feel like the Pope now.

19 (Laughter.)

20 DR. BULLIMORE: But given the construction of
21 the device, I would have thought you'd get some trends
22 occurring.

23 So take those questions in either order you
24 would like them.

25 DR. McCULLEY: And I had one, too, that really

1 fits right in with that. What was the association of the
2 induction of cylinder relative to the 30 degree gap and
3 relative to the segment? Which is what you were saying,
4 too.

5 DR. BULLIMORE: You said it more eloquently.

6 DR. McCULLEY: I didn't have to cross anybody,
7 though.

8 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
9 KeraVision. In regards to Dr. Bullimore's first question
10 comparing what was the mean change for BSCVA, for both
11 protocols for the PMA cohort the main change was 0.3 lines,
12 plus or minus 0.83 lines.

13 DR. BULLIMORE: Sorry. Say that again.

14 MS. CROCKETT-BILLIG: A 0.13 line change.

15 DR. BULLIMORE: It's 0.13 line. So that's just
16 under one letter, basically, translating it.

17 MS. CROCKETT-BILLIG: Yes.

18 DR. BULLIMORE: And that would be consistent
19 with a magnification effect that you might anticipate with
20 the change of the refractive correction moving from the
21 spectacle plane to the corneal plane, so one could explain
22 the change in best spectacle-corrected visual acuity almost
23 entirely on the basis of a magnification effect without
24 even considering what the optical effects of the procedure
25 might be.

1 What about the axis of cylinder? Have we got
2 that?

3 MS. CROCKETT-BILLIG: I'd like Dr. Holladay to
4 address that question.

5 DR. HOLLADAY: Jack Holladay, Houston, Texas,
6 consultant. You're absolutely right. We would expect for
7 the average correction to effect a one letter
8 magnification, so for the overall cohort, with everyone
9 mixed together, you didn't see any significant change from
10 pre to post-op when you took all of the patients. The
11 point was there was a subgroup of those, though, in which
12 there were 20 percent that had more than a one line change,
13 which was statistically significant at a very small P
14 value.

15 The point was that that number we think is a
16 good number. The ones that decreased we also knew could go
17 down as many lines as you wanted to because they weren't up
18 at the upper end. As you saw, three of the four patients
19 that dropped lines were 20/12, so we believe that there is
20 a claim for 20 percent of the patients, which we can't
21 predict, that actually can have a statistically significant
22 improvement in one line in vision.

23 In terms of the whole cohort of the 494, there
24 is no difference between the pre and the post-operative
25 value. So all we'd like to do, if it's possible, is

1 mention that there may be a group of patients who actually
2 do have a one line improvement in their vision.

3 DR. BULLIMORE: I guess I have the same problem
4 with that statement that the reviewers seemed to have, and
5 that is, you know, you can't say that when you have an
6 average variable change that some people got worse and some
7 people got better, and then really just choose to ignore
8 the people who got worse. I mean, that's just sort of
9 statistics that we expect to come from inside the Beltway
10 and not from this part of the country.

11 (Laughter.)

12 DR. BULLIMORE: I accept your statement that
13 some people did get better and I accept the statement that
14 some people get worse, but I'd like to sort of look at the
15 distribution of the whole, and the .13 lines seem to me to
16 be a bit much better indication in terms of what's
17 happening to the sample as a whole and what is the most
18 likely thing that's going to happen to a patient that I
19 would refer or, were that patient me, what would be my
20 expectations of my own vision. So that's the only point I
21 want to make.

22 What about the axis of cylinder? Have you got
23 some data on that?

24 DR. HOLLADAY: Can I make one other comment
25 about that in response?

1 DR. McCULLEY: Well, I think that this is an
2 important point, yes, and I think let's stay on it until
3 maybe anyone else that has comments about it as well,
4 because I think one thing that we should have learned over
5 the years, if we don't make our advisory panel opinion
6 clear, then we don't have the effectiveness that we're
7 being asked to provide.

8 DR. HOLLADAY: I don't disagree with any of
9 those comments at all. As a matter of fact, to pick a
10 subgroup and to say that they are statistically
11 significant, although you can't predict who those are, is
12 certainly a good point.

13 The only reason I raise that is I've been
14 involved with several other PMAs that come from other
15 refractive procedures, but that we've never seen any
16 numbers in which anyone has ever had an improvement in
17 visual acuities in a small percentage, and that's the only
18 reason we brought that up.

19 DR. BULLIMORE: The claim's not been made, I
20 believe, in terms of the things that were presented to me
21 in the literature, the claim's not been made with reference
22 to other refractive surgery procedures. The claim has been
23 made with respect to other corrective technologies, such as
24 spectacles and contact lenses, and the data doesn't support
25 that claim, based on the fact that you've only got a .13

1 line or less than a letter change in visual acuity. So
2 that's really what I'm speaking to. I'm not wanting to
3 compare it to anything else. I just want to get the data
4 straight.

5 DR. McCULLEY: And their study was not done to
6 compare, so that's an inappropriate line of reasoning for
7 these discussions.

8 Karen?

9 DR. HOLLADAY: The cylinder issue --

10 DR. McCULLEY: We'll come back to it. We're
11 still on this enhanced acuity issue.

12 DR. BANDEEN-ROCHE: Yes, Karen Bandeen-Roche.
13 My follow-up question would be that as far as I understand,
14 a one line change in acuity is fairly constant with the
15 measurement error associated with visual acuity, and from
16 Dr. Sugar's report there was approximately -- I mean,
17 slightly less, but not a very different percentage who
18 declined by one letter as who improved by one letter.

19 So my question would be what is the evidence
20 that the hyperacuity phenomenon represents something other
21 than what would be expected by the measurement error of the
22 procedure?

23 DR. McCULLEY: Would sponsor like to respond to
24 that or let the comment stand?

25 MS. CROCKETT-BILLIG: One moment, please.

1 We're getting information together.

2 DR. McCULLEY: I'm sorry, sponsor. We really
3 need for you to come forward.

4 MS. CROCKETT-BILLIG: All right. What we found
5 was that we found 19.5 percent of our subjects did have an
6 increase in their BSCVA of five or more letters or one or
7 more lines. That's what our analysis showed. We have the
8 P value. We found that the increase in BSCVA was
9 statistically significant and with a P value of 0.002.

10 DR. McCULLEY: Do you have comments relative to
11 those that were just made in criticism of your approach to
12 this?

13 MS. CROCKETT-BILLIG: I'd like Dr.
14 Chiacchierini to address that comment, our consulting
15 statistician.

16 DR. CHIACCHIERINI: Yes, I am Dr. Richard
17 Chiacchierini. I was the former director of the Division
18 of Statistics in the Center for Devices and Radiological
19 Health. I'm now the vice president of statistical services
20 at C.L. MacIntosh & Associates. I have no financial
21 interest in the company other than my consulting fee.

22 The whole issue of determination of whether or
23 not something is within sampling variation is the precise
24 nature of statistical testing. The fact that we have a
25 statistically significant improvement is beyond doubt.

1 There is a statistically significant improvement.

2 The issue before this body is whether that is a
3 clinically important improvement, so I think what we need
4 to do is distinguish between statistical significance,
5 because had this been in the sampling error of the
6 distribution, it would not have been statistically
7 significant. It was in fact statistically significant.

8 Now, the fact is, is this a clinically
9 important improvement? That is, a one line increase? That
10 I think is what we're discussing.

11 DR. BANDEEN-ROCHE: This is Karen Bandeen-
12 Roche. If I could just ask for a clarification, could you
13 please just remind us exactly what was the test for which
14 you're quoting statistical significance? I don't mean the
15 method, but what was the denominator, what was being
16 compared, et cetera? Just refresh our memory on that,
17 please.

18 DR. McCULLEY: While we're waiting for sponsor
19 to pull a response together, just a point of clarification
20 for us. In the past, we have taken two lines as being
21 clinically significant in our guidance documents and in our
22 discussions. Whether that's appropriate or not, that has
23 been our tendency and we're talking about one line, not two
24 lines, here. So if the issue is clinical significance and
25 that's your argument, then we have taken two lines in the

1 past as being clinically significant.

2 DR. CHIACCHIERINI: The line change denominator
3 was based on 410 patients.

4 DR. McCULLEY: Do you have a response to that,
5 Karen?

6 DR. BANDEEN-ROCHE: No, although, as I
7 understand it, that's mean change? Mean change. So I
8 think combining Dr. Bullimore's and my comments, I wouldn't
9 change my comment.

10 DR. CHIACCHIERINI: It would be the median
11 change. It was a non-parametric test.

12 DR. McCULLEY: Your comments still stand? And
13 would you state clearly for us what your view on this is,
14 so that we're real clear?

15 DR. BANDEEN-ROCHE: Well, my view is that in
16 terms of the percentage with improved visual acuity, not
17 the median change, but the percentage with a line or more
18 visual acuity improvement, I am not convinced that it's
19 beyond what would be expected by measurement error of the
20 procedure.

21 DR. McCULLEY: Dr. Macsai?

22 DR. MACSAI: Since we're discussing this issue,
23 Dr. McCulley, I would like to bring again to the panel's
24 attention and to the sponsor's attention that we not
25 compare apples with oranges. We do not have the best

1 corrected visual acuity on these patients preop. Their
2 best corrected visual acuity would be in a rigid, gas-
3 permeable contact lens, not a spectacle.

4 So you don't know for sure, statistically
5 significant or not, whether this is of any significance
6 clinically to the patient, and I think it would be
7 inappropriate to make a statement as such, because you need
8 to know what their best corrected visual acuity is preop if
9 you are going to say that post-op they have had an
10 enhancement or improvement of their best corrected visual
11 acuity.

12 DR. McCULLEY: Okay. We're still on this
13 point, prior to coming back to the astigmatism question
14 that Dr. Bullimore raised. Are there other questions on
15 this issue? Dr. Wang?

16 DR. WANG: Ming Wang. I'd like to second Dr.
17 Joel Sugar's comment about lack of masking of the
18 possibility of enhanced visual acuity of the unoperated
19 eye. It's relevant to the hyperacuity issue.

20 DR. McCULLEY: I think that the panel seems to
21 have a fairly consistent opinion relative to that issue.
22 Can we now go to Dr. Bullimore's question about astigmatism
23 gap segment, et cetera?

24 MS. CROCKETT-BILLIG: Yes. Dr. Bullimore,
25 would you repeat that other question to make sure I have it

1 right?

2 DR. BULLIMORE: Yes. I seem to remember
3 somewhere in the materials I was sent the fact that the
4 main or maybe the median surgically induced change in
5 refractive -- or surgically induced change in cylinder,
6 whatever the appropriate term is -- was on the order of
7 half a diopter. My question was is that randomly
8 distributed or was there a tendency for it to be, say, with
9 the rule or against the rule, clustered around a certain
10 meridian?

11 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
12 KeraVision. The induced cylinder that we did have was
13 primarily with the rule cylinder.

14 DR. BULLIMORE: I guess, as a follow-up to
15 that, and this is a question as much for the panel as the
16 sponsor, is that information that would be useful to have
17 in the labeling? For example, were a patient to present
18 with astigmatism at axis 180, one might counsel them
19 differently than if they had astigmatism at axis 90. I'll
20 leave that for discussion later.

21 DR. McCULLEY: Other comments in that regard?
22 Dr. Macsai?

23 DR. MACSAI: Does the sponsor have any data as
24 to whether or not induced astigmatism was orthogonal or
25 non-orthogonal -- i.e., regular or irregular -- in light of

1 the visual symptoms represented by some of the recipients
2 of the device?

3 MS. CROCKETT-BILLIG: I'd like Dr. Schanzlin to
4 answer that question, please.

5 DR. SCHANZLIN: Dr. Schanzlin. It's primarily
6 orthogonal. That is, when you refract these patients,
7 they're easily correctable with spectacles.

8 DR. McCULLEY: Other comments? Dr. Jurkus?

9 DR. JURKUS: This is Dr. Jurkus. I had a
10 question regarding a number of patients who reported
11 difficulties with night driving, halos, and glare, and they
12 were then treated with pilocarpine, and then if I read the
13 table correctly in Volume 3, that treatment was then
14 discontinued. I was wondering if you could give any
15 information regarding what their subjective symptoms were
16 after the discontinuation of pilocarpine.

17 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
18 The patients primarily went back to where they were prior
19 to treatment with pilocarpine; i.e., the treatment really
20 wasn't effective in reducing the visual symptoms for most
21 subjects.

22 DR. McCULLEY: Dr. Matoba?

23 DR. MATOBA: Alice Matoba. I had a question
24 about the 34 patients in whom you have removed the
25 implants. Fourteen were because they were dissatisfied

1 with their vision and 16 because they had visual symptoms
2 with which they were dissatisfied. Were those numbers
3 excluded from your final analysis, the 449 patients for
4 whom you've given us the data? And if so, how would you
5 justify excluding them?

6 MS. CROCKETT-BILLIG: The 449 subjects were
7 included in the safety analysis. We had 449 implants for
8 safety.

9 DR. MATOBA: But did you exclude the 34 in whom
10 you removed the implants? They're not included in the 449,
11 is that correct?

12 MS. CROCKETT-BILLIG: At month 12, we had 410
13 subjects, and they were not included in that information,
14 but they were part of the overall safety of the 449.

15 DR. MATOBA: Overall safety, but not patient
16 satisfaction, for example.

17 MS. CROCKETT-BILLIG: They didn't have an
18 implant at month 12 to have the patient to take the survey.

19 DR. MATOBA: Right, but somehow it doesn't seem
20 valid to just drop them out of analysis, since they weren't
21 satisfied.

22 MS. CROCKETT-BILLIG: We did an alternate
23 analysis, a carry forward analysis, where we looked at all
24 subjects, carrying them through the whole time period as if
25 they would have been in the trial for the whole 12 months,

1 and I'd like Don Young to address that question.

2 MR. YOUNG: Don Young, San Jose State. As
3 Darlene has mentioned earlier, we did a carry forward
4 analysis on all patients if they dropped out of the cohort
5 because of removals and the various other things that
6 you've cited, and those results are presented in the
7 appendix to the cohort report. Looking at the differences
8 in the results, there are no substantial differences in the
9 results.

10 DR. MATOBA: Did you look at the different
11 sizes? For example, the majority of those patients, 56
12 percent, had the .35 millimeter implant. If you looked at
13 different sizes of implants, was there still no significant
14 difference if you subdivided and categorized the data?

15 MR. YOUNG: As far as I can recall, there was
16 no difference.

17 DR. McCULLEY: So your efficacy data did not
18 include the 34.

19 MS. CROCKETT-BILLIG: That's correct.

20 DR. McCULLEY: Your safety data did.

21 MS. CROCKETT-BILLIG: That's correct.

22 DR. McCULLEY: So there has to be clear
23 labeling that in the efficacy data that includes an
24 additional 4.7 percent explant rate.

25 Dr. Macsai?

1 MS. CROCKETT-BILLIG: Darlene Crockett-Billig,
2 KeraVision. Only 20 eyes were explanted at the month 12
3 time period, though.

4 DR. MACSAI: This is Dr. Macsai. I'm a little
5 bit confused. In comparing visual symptoms that are often
6 or always mentioned by these patients, in my previous
7 comment to you talking about pupil size, mesopic, and in my
8 discussion with Dr. Lemp regarding halos, double images,
9 and fluctuating distance vision, those percentages or those
10 increases -- let's just take halos, preop, one out of 357
11 patients, .3 percent. Post-op, 39 out of 328, 11.8
12 percent, a 37 times increase in often/always complaints of
13 halos, not a variation with pupil size.

14 Is 39 out of 328 accurate if you have excluded
15 the dissatisfied patients, the 16 with visual symptoms that
16 caused explantation? Should that numerator in fact be 55?

17 MS. CROCKETT-BILLIG: We didn't do the analysis
18 in that fashion in the data that we submitted. Now, we
19 could certainly go back and look at that.

20 DR. McCULLEY: I'm a little confused, too.
21 This is McCulley. Your symptom percentage does not include
22 those patients who were sufficiently symptomatic to have
23 explant.

24 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
25 Our symptom percentage was based on the subjects that had

1 the implants at month 12. There were 20 subjects who had
2 removals prior to month 12, and those subjects would not be
3 in that specific N for visual symptoms at month 12.

4 DR. McCULLEY: So we just have to be very
5 careful about making sure that when the final statistics
6 are done that the percentages that are presented by one
7 mechanism or another make it clear that there are those
8 that had the symptoms that still had the implant and those
9 that had the symptoms that were so severe that they had the
10 implant removed, so that those numbers don't misrepresent.

11 Dr. Wang?

12 DR. WANG: Ming Wang. I just have a quick
13 comment and a question for the sponsor. The comment's
14 about the .35 grouping possibly needing to be treated
15 separately in response to Dr. Rosenthal. I think probably
16 it will be more than philosophical, in that Dr. Sugar, Dr.
17 Grimmatt, and Dr. Van Meter have all mentioned incidences
18 of possible separate labeling needing for that size.

19 I'm concerned about the doubling of problems in
20 that size compared with other sizes. For example,
21 diplopia, 15.3 versus 6.9 percent. So are we hitting a
22 point of nonlinearity in which that size may present a
23 problem physiological to the cornea?

24 My question to the sponsor is is there any
25 stratification -- this is sort of related to Dr. Matoba's

1 question -- of explant rate and complications according to
2 size?

3 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
4 We have looked at the explant rate according to size, and
5 we're pulling the information up.

6 DR. McCULLEY: I think that was presented,
7 wasn't it? That's already been presented. I think we'll
8 deal with that. Your points are very good, but I think in
9 the normal proceedings of the panel we'll deal with those
10 issues when we come to labeling recommendations. Does the
11 FDA need anything? I think we can deal with it. The
12 issue's real and valid, and it will be dealt with when we
13 deal with labeling issues. It really doesn't require a
14 response from sponsor at this time.

15 Dr. Grimmett?

16 DR. GRIMMETT: Dr. Michael Grimmett. I wasn't
17 here on the panel when your protocol was designed and
18 developed, so this is a background question. When the
19 testing was selected at a 7 millimeter pupil for the
20 mesopic testing, how did you all choose 7 millimeters as
21 your target pupil size and did you analyze the data for 6
22 millimeters, for example, or that just wasn't the way it
23 was done, and did you take into account the Stiles-Crawford
24 effect and other issues like that?

25 MS. CROCKETT-BILLIG: I'd like Dr. Schanzlin to

1 answer that question, please.

2 DR. SCHANZLIN: We just did the cutoff at 7
3 because that's the optical zone where the ring comes into
4 effect, so when we were designing this, we were looking
5 mainly at what is the visualization of the ring. If you
6 have someone whose pupil would dilate to 9, they might
7 actually see the shadow of the ring, so that was why we
8 looked at that cutoff.

9 DR. GRIMMETT: As I recall in George Waring's
10 lengthy book on RK, due to the Stiles-Crawford effect,
11 under scotopic conditions the effective pupil size was
12 approximately 5.5 or something like that.

13 I was just interested in how 7 was chosen. It
14 was just due to the physical dimension.

15 DR. SCHANZLIN: It was physical dimension of
16 the implant.

17 DR. HOLLADAY: Jack Holladay. I wrote that
18 chapter for George.

19 (Laughter.)

20 DR. HOLLADAY: And it's true that when
21 everything's open -- and that goes all the way back. I
22 don't want George to get credit for that either. It goes
23 all the way back to Stiles and Crawford, and it just shows
24 that when you actually do take into the Stiles-Crawford
25 effect that anything beyond a 5.5 millimeter size pupil has

1 no input or any light intensity that forms the main image
2 on the macula.

3 It doesn't, however, relate to peripheral
4 images, shadows, and other things that you see. So the 7
5 millimeter choice was just specifically to see whether it
6 was inside or out of the ring, but there's nothing magic
7 about that 5.5 millimeter, other than when you get bigger
8 than that the contribution to the macular image is small,
9 but it still forms a lot of peripheral images, and that's
10 why we do those, to see if there's peripheral glare or
11 unwanted images. That's the reason.

12 DR. GRIMMETT: Thank you.

13 DR. McCULLEY: Are there other questions from
14 the panel? I have two. One is, you presented no
15 information on topography.

16 DR. SUGAR: Yes, they did.

17 DR. McCULLEY: They did?

18 DR. SUGAR: It wasn't presented today, but it's
19 in there.

20 MS. CROCKETT-BILLIG: Yes, the information is
21 in the panel package.

22 DR. McCULLEY: Okay. Well, my problem with
23 topography has been interpretation of it, so I guess my
24 question then is, not having reviewed that specific hard
25 copy, can you give us a summary of what your topographical

1 analysis showed?

2 MS. CROCKETT-BILLIG: Yes. I would like Dr.
3 Holladay to address that. He's independently reviewed our
4 topography.

5 DR. BULLIMORE: Can we have Dr. Waring?

6 (Laughter.)

7 DR. McCULLEY: Jack, I know you love this
8 stuff. Just bottom line.

9 DR. HOLLADAY: Can we show a slide?

10 DR. McCULLEY: A slide? A slide you can show
11 me and everyone else.

12 DR. HOLLADAY: Well, let me just tell you,
13 rather than showing you, because we're all getting late in
14 the day anyway, but I'll tell you this. The topography
15 data pretty much helped us correlate exactly with what the
16 vision data was, and that was this. We saw the same
17 predicted corneal acuity changes on the topography as we
18 saw in the regular data, and we saw a very high correlation
19 between the changes in those two, meaning this, that about
20 97 percent of the corneas preoperatively were at the
21 optical limit at 20/10. Post-operatively, that dropped
22 down to about 94 or 95 percent at 20/12 and 20/16, the same
23 sorts of things that we saw in the visual acuity change.

24 Again, the with the rule astigmatism that was
25 seen showed up on the topography also, and the aspheric

1 changes that we saw, both in the quantitative maps and on
2 the surfaces, showed up, too. So what it did is it helped
3 -- and that's why I was kind of interested in that enhanced
4 acuity. When you get more prolate, sometimes you can get
5 better optical images and we were unable to correlate those
6 in terms of proving that those were the patients that got
7 better, so that fell through, but the fact is every
8 topographic change in terms of quality, shape, and
9 astigmatism and aspheric change paralleled exactly what
10 their data showed in terms of visual acuity.

11 Any other specific questions about that?

12 DR. McCULLEY: Are there any others? Karen?

13 DR. BANDEEN-ROCHE: Karen Bandeen-Roche. I had
14 a few more specific questions. Hopefully, they shouldn't
15 take very long.

16 The first one goes to patient satisfaction. I
17 was curious. Of those patients who had a good
18 contralateral candidate eye, what was the percentage that
19 actually opted to have the other one implanted?

20 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
21 We'll pull that information for you.

22 DR. BANDEEN-ROCHE: Okay. The second one is
23 very quick. It goes to the rate of having a removal. Are
24 you fairly satisfied that that has stabilized or are there
25 yet patients who might yet opt to have their ring removed

1 in a reasonable expectation?

2 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
3 Yes, we feel the removal rate has stabilized. Actually, we
4 encouraged the patients, even though they were unhappy, to
5 stay in the trial for at least six months to see if their
6 situation would stabilize over time. So we feel that the
7 information is accurate and, again, that is stabilized.

8 DR. BANDEEN-ROCHE: Thank you.

9 The third question goes to data for corneal
10 sensation. I believe in six months there were 259 patients
11 evaluated. At 12 months, that had dropped to 237. So my
12 question is is there association between corneal sensation
13 and those who dropped out? In other words, are we possibly
14 masking a higher loss of corneal sensation from those who
15 weren't evaluated?

16 MS. CROCKETT-BILLIG: Darlene Crockett-Billig.
17 No, that is not the case.

18 DR. McCULLEY: Dr. Wang?

19 DR. WANG: Ming Wang. You have reported
20 adverse events, such as anterior chamber perforation,
21 function decrease of two lines, and also shallow implant
22 segment removed. I was wondering, is there any association
23 of these complications with the experience of the surgeon?
24 It would be important for me as a surgeon to know that
25 perhaps these complications are less likely once one is

1 well trained, if it's indeed offered to the general public
2 and to all surgeons.

3 MS. CROCKETT-BILLIG: I would like to ask Dr.
4 Schanzlin to address that comment.

5 DR. SCHANZLIN: I think properly trained these
6 are easily addressed. Certainly, the perforations we
7 talked about earlier, that's certainly an issue of just how
8 to set a diamond knife blade.

9 At least one surgeon, one where the shallow
10 implant was there, it was a shallow incision to start with,
11 and if you look at what happens when you do a shallow
12 incision and try and make the dissection, you get a lot of
13 resistance. You can easily learn that when you see a wave
14 in front of the advancing dissector and you're struggling
15 to advance it, that means you're too shallow, so you should
16 back off and you can rechannel deeper and get a deeper
17 incision.

18 The other thing is that if you see post-
19 operatively that the ring on day 1 is at less than 50
20 percent depth, that's the time to intervene, not wait six
21 months like we do in trials and see what happens.

22 So these are things that certainly can be
23 addressed in training that can only improve the results, I
24 believe.

25 DR. McCULLEY: Dr. Pulido?

1 DR. PULIDO: Jose Pulido. We know that with
2 radial keratotomy there is inherent instability of the
3 cornea with altitude and there is a decrease in the tensile
4 strength of the cornea. Do those things exist with this
5 implant?

6 DR. McCULLEY: Dr. Schanzlin?

7 DR. SCHANZLIN: Well, there a couple of things
8 in RK that we're concerned about. One is what happens with
9 altitude, and our patients have traveled in airplanes. I
10 don't have anyone that's been to Mount Everest yet, but
11 certainly going up in planes to that equivalent, 8,000
12 feet, and also on the West Coast, people doing a lot of
13 skin diving, it doesn't seem to cause fluctuations in
14 vision based on that.

15 The other large concern, of course, is what
16 happens with time with this, and we've closely looked at
17 the data to look at what happens. Is there any drift,
18 hyperopic drift, over time? It seems like it's very
19 stable, three months, six months, nine months, 12 months.
20 So I don't anticipate we have a danger there.

21 DR. McCULLEY: Dr. Higginbotham?

22 DR. HIGGINBOTHAM: Dr. Eve Higginbotham. A
23 similar question, Dr. Schanzlin. We have seen changes in
24 refraction in women that are perimenopausal and post-
25 menopausal. Did you analyze the data for those factors in

1 terms of the women?

2 DR. SCHANZLIN: Yes, we did, and am I right in
3 saying that there was no difference between males or
4 females, no difference between females premenopausal or
5 post-menopausal? That's true.

6 DR. McCULLEY: I think this last question that
7 I have, Dr. Schanzlin, would be most likely for you as
8 well. You mentioned that there was localized staining in 4
9 to 7 percent of patients at various time points. Where was
10 that located on the cornea? Was there a pattern or was it
11 random?

12 DR. SCHANZLIN: Well, the 4 percent, of course,
13 you're talking about the pooled data and it's pretty much
14 random, usually at the limbus, 7 o'clock and 5 o'clock,
15 close to the limbus.

16 DR. McCULLEY: Outside the ring.

17 DR. SCHANZLIN: Outside the ring.

18 DR. McCULLEY: So none in the center that you
19 specifically commented on? Again, getting back to the
20 question of any adverse effect from the decrease in corneal
21 sensation, have you looked at the central corneal
22 epithelium with specular confocal to try to determine
23 whether it's healthy or not?

24 DR. SCHANZLIN: We have not and certainly have
25 not submitted that. I believe Penny Aspell presented, I

1 think it was at CLEO, some time ago looking at ultrasonic
2 microscopy of epithelial thickness close to the ring
3 centrally, and I believe, if I recall, the conclusion was
4 that there was a mild thickening of the epithelium just
5 inside of the ring, but again, she's not here to present
6 that data and we did not submit it. It's a one-center
7 study that has not been verified anywhere else.

8 DR. McCULLEY: Okay. Thank you.

9 Any other questions from the panel before we
10 break for lunch?

11 (No response.)

12 DR. McCULLEY: Seeing none, Ms. Thornton has a
13 couple of comments, and then we're going to take a one-hour
14 lunch break.

15 MS. THORNTON: I'd like to announce that the
16 lunch for the panel will be in Room 20H just down the hall.
17 That room has been reserved for panel and for FDA. The
18 sponsor has Room 20G, which has been reserved for their
19 lunch.

20 Please take with you your cups, your cans, your
21 paper products, and deposit them in the rubbish bins just
22 outside the door at the request of the management.

23 Thank you.

24 (Whereupon, at 12:26 p.m., the meeting was
25 recessed for lunch, to reconvene at 1:26 p.m.)

AFTERNOON SESSION

(1:40 p.m.)

1
2 DR. McCULLEY: Can I ask everyone please to
3 take your seats? Let me ask the panel just for a general
4 feel. Do we want to have open discussion and then go to
5 the questions, or do we want to incorporate our discussion
6 in the questions? In the questions. There seems to be
7 strong sentiment for that.

8 Do you guys want to project your questions,
9 FDA? Malvina? I'll start reading the first question.

10 "Do the outcomes of the endothelial cell
11 density analysis presented in this PMA provide reasonable
12 assurance of safety for all three thicknesses of the ICRS?
13 What, if any, additional data are needed to make this
14 decision?"

15 I think the primary reviewers seemed to have a
16 pretty consistent feel about this. Joel, can I ask you to
17 open discussion?

18 DR. SUGAR: I think all three of our reviews
19 suggested that some further analysis of the patients that
20 have had difficulties, especially those who had 0.35 rings,
21 ought to be obtained, further longitudinal data.

22 DR. McCULLEY: So postmarket study. Is that
23 the appropriate word? We have to be careful in terms of
24 the words. So it's postmarket study.

25 Dr. Grimmett, do you have anything further to

1 add to that? I think there was consensus.

2 DR. GRIMMETT: Dr. Grimmett. I echo those
3 comments.

4 DR. McCULLEY: Does anyone disagree with that?
5 (No response.)

6 DR. McCULLEY: Okay. My impression is that if
7 we are going to request that, that it is advantageous for
8 us to be as specific and narrow in our recommendation as
9 possible, not make an open-ended recommendation. Does
10 anyone have a specific recommendation as to how one would
11 like to see the study carried out?

12 Dr. Sugar?

13 DR. SUGAR: Sugar. I would like to see on the
14 patients who had -- I don't want to use the term
15 "significant" -- on patients who had cell losses in the 10
16 percent range, peripheral or central, that further data be
17 obtained on those patients, if possible. Also, that
18 prospectively -- and I don't know how to get specific with
19 this, but that a subset of future patients receiving the
20 0.35 ICRS have endothelial data obtained on them.

21 DR. McCULLEY: Okay. I gave some thought to
22 the specificness of it, and it's in the same vein as what
23 you've stated. It's that we request two-year follow-up on
24 the patients that have had 10 percent or more cell loss.
25 What we have as our controls in that are going to be the

1 other two size rings, that a comparable number of patients,
2 a matched number of patients with the 0.25 and 0.30 ring be
3 done as well as the control, because that's our comparison,
4 and that we request an additional 25 patients with the 0.35
5 ring for future study, starting from preop through two
6 years.

7 Dr. Higginbotham?

8 DR. HIGGINBOTHAM: Eve Higginbotham. Would you
9 also specify from the same site, considering our previous
10 discussion about site differences in outcome?

11 DR. McCULLEY: I wouldn't have a problem with
12 that. Well, I'm not sure that you really want same site,
13 because we want real world. I think let's leave that -- I
14 would suggest we leave that detail out.

15 DR. HIGGINBOTHAM: I'm suggesting a set of
16 case-controls from the same site as the patients --

17 DR. McCULLEY: For the prospective patients.

18 DR. HIGGINBOTHAM: Exactly.

19 DR. McCULLEY: Because the others are already
20 scattered, presumably, among more than one site. Okay, I
21 think that's reasonable.

22 Marian?

23 DR. MACSAI: Dr. McCulley, I'm not certain what
24 number of patients are needed in each group to ascertain
25 statistical significance in endothelial cell loss, so I

1 question why we would just randomly ask for 25 more.
2 Perhaps someone from the agency, such as Dr. Eydelman or
3 Dr. Boulware, could help me out here?

4 DR. McCULLEY: Recommend a number that would
5 give us some degree of confidence? Okay. So instead of
6 saying a specific 25, a number of patients to be studied
7 prospectively that will give us a statistical degree of
8 comfort.

9 DR. MACSAI: To detect what?

10 DR. McCULLEY: I know you two were talking at
11 the time. Did you hear what I said?

12 DR. EYDELMAN: Part of it.

13 DR. MACSAI: I would think what we would be
14 interested in is a greater than 6 percent cell loss rate in
15 the periphery, greater than that found by Bourne in the
16 normal --

17 DR. McCULLEY: Bourne was 0.6 centrally.
18 Again, what I suggested was a 0.25, a 0.30, and a 0.35 as
19 the control, the first two for control relative to the
20 0.35, which is the one that is carrying our significant
21 concern, and we will assume that there's -- I mean, that's
22 making the assumption that there is not a significant
23 change in the 0.25 and 0.30. But we can't ask them to
24 define what happens in the normal population in --

25 DR. MACSAI: Right, I understand.

1 DR. McCULLEY: So I think that's putting too
2 much on them. But I think if the prospective study is
3 statistically -- the sample size is statistically
4 determined with the three ring diameters --

5 DR. MACSAI: The sample size may already be
6 there. Excuse me, Dr. Macsai. The sample size may already
7 be enrolled. They may not need to enroll more patients.

8 DR. McCULLEY: Marian, what I said --

9 DR. MACSAI: That's what I'm trying to say.

10 DR. McCULLEY: There were two parts to the
11 suggestion. There are those that have demonstrated
12 significant cell loss to date, that they be followed to two
13 years. That data is submitted.

14 DR. MACSAI: Yes.

15 DR. McCULLEY: That sample size is
16 predetermined. That a to-be-determined sample size for
17 each of the three ring sizes be prospectively enrolled from
18 preop through a to-be-determined length of time
19 additionally to support whether there is endothelial cell
20 loss in the periphery or not.

21 DR. MACSAI: I guess I don't understand the
22 second part of that recommendation. If they have three
23 groups, they've done endothelial cell counts that have been
24 monitored by a central unit for 0.25, 0.30, 0.35
25 millimeters. Why not just follow those patients out

1 longitudinally?

2 DR. McCULLEY: If that sample size would give
3 statistical significance, then I think that would be
4 reasonable. If it won't, then I would request the
5 prospective study. I was anticipating that that sample
6 size would be very difficult to get statistical
7 significance on.

8 DR. EYDELMAN: Dr. Eydelman. Panel would need
9 to help us specify the difference that you're trying to
10 detect as to the percent loss between 0.35 and 0.30 and
11 0.25 rings, because we have done statistics previously in
12 order to establish a 10 percent loss, and you need at least
13 80 eyes. So just comparison of the groups, we need another
14 variable to detect the sample size -- i.e., different
15 within what percentage? What would be your equivalency
16 threshold between the two endothelial cell populations?

17 DR. McCULLEY: You mean a percentage greater
18 than?

19 DR. EYDELMAN: Yes.

20 DR. McCULLEY: Well, we've talked about the
21 variability of specular before, and the 10 percent has been
22 the number.

23 DR. EYDELMAN: So 10 percent over what length
24 of follow-up for the second group? I understand for the
25 first group you want two years. What length for the second

1 group?

2 DR. McCULLEY: Is the sample size large enough
3 as it exists if followed to two years to give statistical
4 confidence as to whether there is a progressive loss or
5 not? That's a question.

6 DR. EYDELMAN: Progressive loss of what
7 magnitude?

8 DR. McCULLEY: Greater than 10 percent. Ten
9 percent or greater.

10 DR. EYDELMAN: Progressive 10 percent over what
11 period of time?

12 DR. McCULLEY: Two years.

13 DR. EYDELMAN: Yes. If you follow the same
14 sample as they have currently for two years as opposed to
15 one year, they were able to determine that there was no
16 loss of 10 percent at one year.

17 DR. McCULLEY: Okay. So the recommendation
18 could be as specific as following the current population
19 that have lost more than 10 percent for two years, for a
20 total of two years data gathering.

21 DR. EYDELMAN: That's the first part of your
22 recommendation. And the second part?

23 DR. McCULLEY: Well, the second part, again,
24 may need to drop. If you can reach a reasonable degree of
25 comfort one way or the other with the current sample size,

1 then I don't think the second part is necessary. If you
2 can't, then I think that the second part will help.

3 Dr. Pulido?

4 DR. PULIDO: How about looking at all the
5 enrolled patients in the 0.35 millimeter group over the two
6 years and comparing that to a randomly selected group from
7 the 0.30 and 0.25 millimeters over the same two-year length
8 of time?

9 DR. McCULLEY: I thought I had a better
10 suggestion than you. I clearly didn't. So I think that
11 sounds reasonable as well.

12 DR. MACSAI: Dr. McCulley?

13 DR. McCULLEY: Yes, Dr. Macsai.

14 DR. MACSAI: We've spoken many times that these
15 patients who have refractive surgery frequently have
16 preoperatively worn contact lenses, and upon cessation of
17 their contact lenses they realize an increase in their
18 endothelial cell count. So perhaps, in fact, you would
19 have to stratify that data to clarify pre-Intact
20 implantation contact lens wear and not.

21 DR. McCULLEY: They did that before, and I
22 would assume that they would continue to do that.

23 DR. MACSAI: But to make sure that those groups
24 were large enough for statistical validity.

25 DR. McCULLEY: Matched appropriately.

1 Is that specific enough for you? I mean, I'm
2 not sure we got it nailed.

3 DR. EYDELMAN: Neither am I.

4 DR. McCULLEY: So you need more help?

5 DR. EYDELMAN: If you can.

6 DR. McCULLEY: Yes, Marcia.

7 DR. YAROSS: Marcia Yaross. I would suggest
8 that the sponsor be given an opportunity to demonstrate
9 that by postapproval follow-up of the existing enrolled
10 population, they can make the case for determining, with
11 adequate statistical significance, that there is no safety
12 issue before they be asked to enroll additional prospective
13 subjects.

14 DR. McCULLEY: Okay, that's well stated. I
15 mean, that was, in effect, the intent. That's better
16 stated. That's the general comment, but then that leaves
17 it for FDA and sponsor to determine exactly who is going to
18 be looked at. I think for how long, we're saying two years
19 at a minimum. And whether you're going to take all of the
20 0.35 group or whether you're going to take the 0.35 group
21 that's had 10 percent or greater loss, get a matched group
22 in the 0.25 and 0.30, I guess we really didn't hit that,
23 although there's no reason you can't do both. It could be
24 a two-pronged study. You could take all of the 0.35, those
25 that have had 10 percent or more, match with 0.25 and 0.30

1 population, and follow those to two years as well.

2 Dr. Matoba?

3 DR. MATOBA: Alice Matoba. I think we still
4 have to specify what percentage cell loss we're looking
5 for. Otherwise you could end up with data that doesn't
6 show significant loss, but you haven't proven that there is
7 no -- the converse hasn't been proven. So I think we still
8 have to specify what percent loss per year we are looking
9 for.

10 DR. McCULLEY: Talking about greater than 10
11 percent? Or from baseline at two years? One of the
12 problems with specular is it's not as precise as we'd like
13 it to be.

14 Dr. Macsai, you were --

15 DR. MACSAI: Well, I was going to stick my neck
16 out for a change.

17 (Laughter.)

18 DR. MACSAI: Dr. Grimmett's calculations, which
19 I double-checked somewhere here, he calculated 5 percent
20 loss per year. At 23 and three-quarters years, you're
21 toast, or your cornea is toast. So that means if you're 25
22 when this is implanted, when you're 48.75, your endothelial
23 cell count has dropped to 800 and you're on a transplant
24 list.

25 DR. McCULLEY: Dr. Matoba?

1 DR. MATOBA: Alice Matoba. I had a question
2 about those calculations, and that is, did you assume that
3 it was 5 percent cell loss for the whole cornea over that
4 period of time? But really, it was like a focal area that
5 showed 5 percent loss, not the central or 6 o'clock. So I
6 don't think you can assume 5 percent loss per year over
7 all. You were probably over-estimating the effect of that
8 finding.

9 DR. MACSAI: Right.

10 DR. McCULLEY: Dr. Van Meter?

11 DR. VAN METER: I would like to reiterate that
12 point. Given the fact that this was a focal measurement
13 and the fact that our ability to measure cells in this
14 particular area of concern is more questionable than our
15 ability to measure central corneal epithelial cells, I
16 would favor just continuing to follow these patients and
17 see where they go. It might be possible to follow all of
18 the 0.35 cohort, but I'm not sure that -- again, I think
19 we're collecting data and we're not sure how to manage the
20 data we get.

21 DR. McCULLEY: I think that's true, that before
22 the fact, we're not necessarily sure -- Dr. Macsai?

23 DR. MACSAI: Dr. Macsai. In fact, you're
24 correct, Dr. Van Meter. This is perhaps data where we're
25 not sure, but we're talking safety. We're talking long-

1 term safety of a new device. If the cohort is already
2 enrolled, they've already been measured centrally at 6 and
3 10, continue to follow that cohort and see what happens.

4 DR. VAN METER: I thought that's what I said.

5 DR. McCULLEY: Okay. We had agreement before
6 that we were going to recommend a postapproval study, and
7 what we're trying to do is fine-tune what that specific
8 recommendation is. What I seem to be hearing is that the
9 recommendation is to follow for up to a maximum of two
10 years, or for two years as a minimum for now, the 0.35
11 enrolled cohort. Does that accurately reflect the
12 sentiment?

13 DR. ROSENTHAL: This is Dr. Rosenthal. And
14 what about the 0.25 and the 0.30? Do you want us to sample
15 those and continue to follow some of those?

16 DR. MACSAI: Yes.

17 DR. ROSENTHAL: I think that's a reasonable
18 request. I think the real issue is that you don't know
19 what's really going to happen.

20 DR. McCULLEY: That's right.

21 DR. ROSENTHAL: It's fun to follow them, but
22 what are we going to do? I think that rather than start on
23 another path of another prospective study, to look at these
24 patients and look at them carefully I think is a very wise
25 recommendation and one which I think the company could

1 probably live with quite happily, and would want to live
2 with.

3 DR. McCULLEY: They'd probably rather not do
4 it.

5 Are there any other comments about the
6 specular?

7 (No response.)

8 DR. McCULLEY: So Question 1 is dealt with.

9 Question 2. "Do the assessments of visual
10 symptoms provide reasonable assurance of safety for all
11 three thicknesses of the ICRS?"

12 DR. SUGAR: That's also the same as Question 5,
13 I think.

14 DR. MATOBA: No, it's not.

15 DR. McCULLEY: Question 2 is limited to
16 symptoms. Question 5 brings up other issues, and that then
17 raises a question about how best to approach this issue
18 relative to the labeling and the different ring sizes. Dr.
19 Grimmett had eight specific points that he brought up about
20 the 0.35. Others, Dr. Wang and others brought up issues
21 relative to the 0.35 and labeling. So do we want to
22 combine 2 and 5?

23 DR. BULLIMORE: This is Dr. Bullimore. I think
24 that there's a consensus that we do need to include this in
25 the labeling, but basically it's okay and we should move on

1 and cover it under Question 5.

2 DR. McCULLEY: We can do that.

3 DR. MACSAI: I beg to differ.

4 DR. McCULLEY: Dr. Macsai?

5 DR. MACSAI: Excuse me. I respectfully
6 disagree with Dr. Bullimore. I think that the visual
7 symptoms that these patients are discussing with the
8 sponsor represent a quality of vision, which is maybe
9 different from visual acuity. I understand we have
10 limitations in our ability to measure them, but it's
11 important to include the patients who had the segments
12 removed due to their visual symptoms, and that's not so in
13 the table that we have. So that the statistics are lower
14 in these tables than, in fact, if you include those
15 patients.

16 DR. ROSENTHAL: This is Dr. Rosenthal. I think
17 what you're saying is that there should be some reference
18 in the labeling to the issue of symptoms.

19 DR. MACSAI: Yes.

20 DR. ROSENTHAL: Both for those who have
21 maintained the ring in the cornea and for those who have
22 had it taken out, so that patients and doctors know what
23 they have to face when they are contemplating putting in
24 that ring. You would like that reflected in the labeling.

25 DR. MACSAI: Yes.

1 DR. BULLIMORE: This is Dr. Bullimore. I agree
2 with Dr. Macsai.

3 DR. McCULLEY: Okay. There are multiple
4 issues. Let me do this. Let me ask Dr. Grimmett to read
5 -- we're not going to downplay symptoms, we're not going to
6 downplay signs or anything else. But let me ask Dr.
7 Grimmett at this point to read his eight concerns about the
8 0.35 ring that would presumably be the substance of the
9 label warning or the labeling.

10 DR. PULIDO: Point of order.

11 DR. McCULLEY: Okay.

12 DR. PULIDO: Are we proceeding, then, to
13 Question 5?

14 DR. McCULLEY: No. We're moving 5 to 2. We're
15 combining 2 and 5 and taking it at 2. The two are tied, so
16 we're moving them all to one basket and we're going to deal
17 with them now.

18 Dr. Grimmett?

19 DR. GRIMMETT: Dr. Grimmett. Dr. McCulley, did
20 you just want me to list the headers for the eight
21 sections? Like the first was endothelium. Do you want me
22 to just list them like that?

23 DR. McCULLEY: I think so.

24 DR. GRIMMETT: Okay. The eight issues that I
25 believed were different for the 0.35 millimeter implant

1 based on the data: number 1 was endothelium; number 2 was
2 induced cylinder greater than 1 diopter; explantation;
3 frequency of certain visual symptoms; magnitude of certain
4 visual symptoms; uncorrected visual acuity; stability of
5 manifest refraction spherical equivalent; and
6 predictability.

7 DR. McCULLEY: So if we recommended putting in
8 the labeling that with the 0.35 ring there were increased
9 concerns about those eight issues, would that deal
10 effectively with what we want to have addressed?

11 DR. SUGAR: I think the specific data should be
12 included. This is Sugar. So that there should be a table
13 listing for each ring size, the frequency of whatever it is
14 of those eight points that we're looking at. This would be
15 for the physician's package.

16 For the patient's package, I think it would
17 have to be put into some kind of more user-friendly
18 wording.

19 DR. McCULLEY: So recommend including raw data
20 or summary data.

21 DR. JURKUS: This is Jan Jurkus. The other
22 thing I would suggest is not only including just the ring
23 size but also the refractive data, that the 0.35 ring is
24 for the 3 diopter myope, and if you are this moderate type
25 of myope, these might be problems that you have,

1 particularly for the patients, since there are a number of
2 other methods that can correct that patient that might have
3 fewer symptoms. The patient might not know what a 0.35
4 ring size is.

5 DR. McCULLEY: Oh, okay. But in the labeling
6 it would have to be clear to them what refractive error was
7 being addressed.

8 DR. ROSENTHAL: This is Dr. Rosenthal. I agree
9 that one should put in a table. I think the issue is do
10 you emphasize the areas -- I mean, I would feel it
11 appropriate to emphasize the eight areas that Dr. Grimmett
12 mentioned, to present the data and then to say that these
13 are the following problems that can exist. That's how I
14 would proceed.

15 DR. McCULLEY: Both a statement and the
16 numbers.

17 DR. ROSENTHAL: As well as the numbers.

18 DR. McCULLEY: Yes, yes. I agree.

19 Dr. Pulido?

20 DR. PULIDO: Just a question for the panel. I
21 think we all feel very comfortable in approving -- or at
22 least I do. I feel comfortable in approving the 0.25
23 millimeter and the 0.30 millimeter, but I have questions
24 still about approving the 0.35 millimeters. We're moving
25 along as if we are going to approve the 0.35 millimeter as

1 well. Do we want to first find out if people want to
2 approve all three of them?

3 DR. McCULLEY: Okay. Point raised. Good
4 point. Is there sentiment -- I mean, this gets the cart
5 before the horse, and that's one of the problems with going
6 down the questions. But since we're bringing up issues
7 that are different with the three different sizes, then I
8 guess we're working on labeling for 0.35 as though it's
9 going to be part of the total recommendation. Is there a
10 sentiment to continue to include 0.35 as part of our
11 overall discussion, or should it be taken out into a
12 separate discussion? That's trying to say politically what
13 was already alluded to. Is there sentiment to continue
14 discussion of all three ring sizes, with some possible
15 differential labeling language?

16 DR. SUGAR: Yes.

17 DR. BULLIMORE: Yes.

18 DR. VAN METER: Yes.

19 DR. McCULLEY: Any dissent?

20 DR. MACSAI: I'd like to raise the question to
21 our medical officer. This is Dr. Macsai. Will longer-term
22 follow-up of the 0.35 segments that are implanted provide
23 more data on stability, more data on the rate of explants
24 and on the endothelium loss? Simply, do we need to look at
25 that segment for more than 12 months? If we looked at that

1 segment a year from now, would we feel more comfortable
2 approving it?

3 DR. ROSENTHAL: This is Dr. Rosenthal. You've
4 already asked that it be looked at for an additional year
5 for the endothelial, all the patients with the 0.35.

6 DR. MACSAI: That have endothelial cell counts.
7 That's a subgroup.

8 DR. ROSENTHAL: Oh, a subgroup. I see. So now
9 you're suggesting possibly looking at all of them?

10 DR. MACSAI: I'm asking the panel.

11 DR. McCULLEY: We're getting two questions on
12 the floor at the same time. Let's try to get resolved the
13 first question before we go on to whether we want
14 additional postmarket.

15 Dr. Van Meter?

16 DR. VAN METER: I'd like to directly answer the
17 question. I think it's reasonable to follow those patients
18 that have already had specular -- those 0.35 millimeter
19 implants that have already had specular started and follow
20 them longer and see where this goes, because I don't think
21 the additional data on additional patients -- the
22 recruitment of more patients is not going to give us data
23 that's helpful. I would also like to point out that the
24 procedure is potentially reversible and the visual symptoms
25 can be a labeling issue. But I think it's very reasonable

1 to include this in the discussion.

2 DR. McCULLEY: Okay. So is the agreement that
3 we will continue to discuss all three ring sizes, keeping
4 the option of differential labeling recommendations?

5 DR. MACSAI: Yes.

6 DR. McCULLEY: Is there dissent to that?

7 (No response.)

8 DR. McCULLEY: Unanimous -- oops. Karen?

9 DR. BANDEEN-ROCHE: This is just a question.
10 Karen Bandeen-Roche. In terms of a bad scenario, we follow
11 for two years and something is alarming, what happens at
12 that point? Is approval rescinded? What happens?

13 DR. McCULLEY: That's an FDA policy-driven
14 issue, I believe, and not under our purview. Not to try to
15 skirt your issue, but it's not our issue.

16 Dr. Wang?

17 DR. WANG: Ming Wang. I just want to echo Dr.
18 Pulido's comment just now and clarify for myself. Are we
19 going to treat 0.35 separately or still discussing as a
20 whole?

21 DR. McCULLEY: Just a moment ago we had
22 unanimous consent that we were going to discuss all three.

23 DR. WANG: Okay.

24 DR. McCULLEY: We just decided that.

25 DR. WANG: As one premarket approval.

1 DR. McCULLEY: Yes. But we reserve the right
2 and the ability to have labeling that is specific for
3 individual ring sizes, if we choose.

4 Renee?

5 DR. MIDDLETON: Renee Middleton. My comment is
6 with respect to all three ring sizes. While we're talking
7 about the labeling under the patient booklet, when the
8 question is asked what are the risks and we have the chart
9 there with night vision and blurring vision, one of the
10 things that stands out for me is the fact that it says you
11 may have discomfort and/or pain. From a person who really
12 dislikes pain, that stands out to me. I think that it
13 would be helpful to indicate to the patient did everybody
14 experience pain, how much pain, and you'd have the
15 information. They'd report the information on pain and
16 discomfort, 15 percent of the individuals.

17 Somewhere that should be reported to give me a
18 sense that I may not have to worry about pain. But if I
19 do, it's only like 15 percent of the individuals that they
20 sampled.

21 DR. McCULLEY: Good point. So our
22 recommendation would be, along with the -- in the labeling,
23 to indicate the severity and frequency of pain that a
24 patient is apt to encounter.

25 DR. MIDDLETON: Because two days of pain can be

1 significant.

2 DR. McCULLEY: Oh, you bet.

3 So is the sentiment, then, that we would have
4 the general statement relative to the eight points that
5 were brought up? Did we have fluctuating vision? We had
6 25 percent of patients with the 0.35 ring had fluctuating
7 vision of greater than 0.5. We didn't have a stability in
8 your eight. So there would be nine.

9 DR. GRIMMETT: Yes, there was.

10 DR. McCULLEY: Seven? Less stable. You got
11 it. Sorry.

12 Does that cover all of the additional issues we
13 had with the 0.35 ring? There was increased -- no, there
14 wasn't.

15 DR. GRIMMETT: This is Dr. Grimmatt. I do have
16 some typed summary comments that I could simply hand over
17 and you could copy them and it enumerates all this.

18 DR. McCULLEY: Okay. Is there agreement on the
19 eight points, the general statement in the table with
20 percentages? Is there disagreement?

21 (No response.)

22 DR. McCULLEY: Okay. Question 3. "Do the
23 reports of corneal sensation losses provide reasonable
24 assurance of safety of all three thicknesses of the ICRS?
25 What, if any, additional data are needed to make this

1 decision?"

2 Dr. Grimmer?

3 DR. GRIMMETT: Dr. Grimmer. I believe the
4 answer is yes, reasonable assurance of safety. As either
5 Dr. Schanzlin or Lemp mentioned, of the 13 that had corneal
6 sensation losses, the 5 percent at 12 months, subsequently
7 they all regained less than 20 millimeters of difference on
8 the esthesiometer. Additionally, those patients with
9 corneal sensation loss at 12 months had no evidence of
10 corneal staining. So I think the answer is yes.

11 DR. McCULLEY: Would it be reasonable to put in
12 the labeling that a decrease in corneal sensation has been
13 noted and it's of unknown clinical significance at this
14 point? Just to put something in the labeling.

15 DR. GRIMMETT: I agree. It's been noted in
16 select patients and we have reason to believe it may be
17 temporary in nature.

18 DR. McCULLEY: And of no yet-determined
19 clinical significance.

20 Any other comments on that point?

21 (No response.)

22 DR. McCULLEY: Okay, Question 4. "The range
23 for the average correction achieved with 0.25 millimeters,
24 0.30, and 0.35 ICRS is from -1.48 to -2.76 diopters. Does
25 the achieved correction data support requested indications

1 for patient population with preoperative myopic error
2 ranging from -1 to -3.50 with 1 diopter or less of
3 astigmatism?"

4 Dr. Matoba?

5 DR. MATOBA: It might be less confusing to say
6 up to -3 spherical equivalent with up to 1 diopter of
7 astigmatism, because I think you can only implant up to -3
8 spherical equivalent. I think they go to -3.5 because you
9 might have 1 diopter of astigmatism, but the way it's
10 stated is a little vague.

11 DR. McCULLEY: Yes, I was confused by that.
12 But I see sponsor shaking their head, which they're not at
13 liberty to un-confuse us. There is the issue about the
14 3.50, the 3, the spherical equivalent, the sphere, and so
15 forth.

16 DR. MACSAI: Spectacle plane, corneal plane.

17 DR. McCULLEY: Is that all it is?

18 DR. MACSAI: No, I think it's not.

19 DR. McCULLEY: What's the -3.50 at the
20 spectacle plane or the corneal plane, or vice-versa?

21 DR. MACSAI: It's not the spectacle plane
22 versus the corneal plane. My understanding from reading it
23 is it's the spherical equivalent, because -3 plus 1
24 equals -3.

25 DR. McCULLEY: That was your understanding of

1 it, and that's what Dr. Matoba said, and that makes sense,
2 but in reading the audience, which I'm not supposed to do,
3 I'm not sure that's what was intended. We can make our
4 recommendation.

5 Dr. Eydelman?

6 DR. EYDELMAN: I believe the indication
7 statement is referring to spherical equivalent because the
8 tables in my review, which came from the sponsor's
9 submission, were in terms of spherical equivalence,
10 preoperative CRSE cycloplegic spherical equivalent.

11 DR. McCULLEY: Okay. So it is a -3 upper limit
12 spherical equivalent.

13 DR. BULLIMORE: No, 3.50 spherical equivalent.

14 DR. McCULLEY: Dr. Eydelman?

15 DR. EYDELMAN: Back in my slides, what I tried
16 to break out was -- if you give me a second, I'll try to
17 get to that slide. This is referring to greater than 3
18 diopters of spherical equivalent.

19 DR. BULLIMORE: This is Dr. Bullimore. You
20 have a slide that's headed "Appropriate Refractive Range"
21 and it gives the RPR for all three thicknesses and then the
22 proposed indication. I think that clarifies it. This is
23 the one.

24 DR. EYDELMAN: Does that help?

25 DR. BULLIMORE: Yes.

1 DR. McCULLEY: I think we're understanding
2 amongst ourselves. I'm going to exercise the prerogative
3 of asking a person from the sponsor to approach the podium,
4 if you wish, to clarify this beyond what we have stated,
5 which is, in effect, it would be -- the upper limit of
6 approval would be a -3 spherical equivalent.

7 DR. ROSENTHAL: This is Dr. Rosenthal. Mr.
8 Chairman, if you look at Dr. Eydelman's next slide, she has
9 shown you the patients greater than -3 diopters spherical
10 equivalent as well.

11 DR. McCULLEY: I know.

12 DR. ROSENTHAL: All right. So I'd like you to
13 consider that.

14 DR. McCULLEY: Okay. Is it 3 or 3.50 spherical
15 equivalence? A simple question, I think.

16 DR. BULLIMORE: No, that's not a simple
17 question.

18 DR. McCULLEY: I knew I should never say that.

19 DR. BULLIMORE: If you go back a slide, the
20 proposed indication says "for reduction or elimination of
21 myopia of -1 to -3" blah blah blah, "in patients" quack
22 quack quack, "with preoperative myopic error ranging from
23 -1 to -3.50." The way I read that is that the range of
24 likely correction is -1 to -3, but it would be labeled such
25 that it could be done on patients up to -3.50, on the

1 understanding that they would be slightly, on average,
2 undercorrected. The first clause is the correction. The
3 second is the patient range that would be eligible.

4 DR. McCULLEY: And we had data that was shown
5 on those greater than -3 that actually came in under
6 guideline.

7 DR. BULLIMORE: Yes.

8 DR. McCULLEY: Okay. So let's just be
9 absolutely sure that we're all on the same page.

10 Dr. Holladay?

11 DR. HOLLADAY: Jack Holladay. One comment.
12 Could we go back to that prior slide? It's exactly what
13 we're saying. We want to get from 1 to 3 diopters of
14 effect in the reduction of myopia. The patients in the
15 study went up to 3.50, and actually a little higher, that
16 got that effect. Our question may be that we might even
17 consider dropping that last line because there are
18 patients, for example, that may be -4.50 or -4 who you want
19 to be -1 and still get a 3 diopter effect.

20 So what we'd like to say is "for the reduction
21 or elimination of myopia from -1 to -3," period. And the
22 range of corrections doesn't matter.

23 DR. McCULLEY: That becomes a practice of
24 medicine issue. So your request, then, is -1 to -3
25 spherical equivalent.

1 DR. HOLLADAY: Correction.

2 DR. McCULLEY: Correction, and dropping the
3 second portion of this.

4 Dr. Eydelman?

5 DR. EYDELMAN: If we follow that suit, we have
6 to somehow assure ourselves that, since we don't have any
7 data, that if we implant the 0.35 millimeter ring in a 5
8 diopter myope, we're still going to get 3 diopter
9 correction. I don't know that.

10 DR. MACSAI: We have no clue about that. That
11 data hasn't been presented, and as you said earlier, Jim, I
12 believe that's a practice of medicine issue regarding
13 monovision, et cetera, or presbyopia. But really upon
14 explanation, reduction of -1 to -3 diopters of myopia -- do
15 you understand?

16 DR. McCULLEY: Yes. It is in this patient
17 population. The second part is qualified.

18 DR. MACSAI: You need to have it, because
19 otherwise you could have -4 plus 2.

20 DR. McCULLEY: Good point. Then the question
21 comes back, are we comfortable with it being -3.50?

22 DR. MACSAI: Yes.

23 DR. McCULLEY: Dr. Pulido?

24 DR. PULIDO: Yes. From Table 4 of page 40 of
25 Dr. Eydelman's review, 99.8 percent were less than 4

1 diopters, and 99.1 percent were from 0.5 up to 3.5
2 diopters. So I think by putting it as presently labeled
3 would be very reasonable.

4 DR. McCULLEY: General sentiment to that, now
5 that we're unconfused and all understand what we're talking
6 about? Is there any disagreement?

7 (No response.)

8 DR. McCULLEY: Are you satisfied, Dr. Eydelman?

9 DR. EYDELMAN: Yes.

10 DR. BULLIMORE: This is Dr. Bullimore. I'm in
11 violent agreement with myself.

12 (Laughter.)

13 DR. BULLIMORE: The only issue I want to raise,
14 and I raise it tentatively, is the issue of overcorrection.
15 Given the fact that we've set the lower limit of the range
16 at -1, and the average effect for the 0.25 millimeter
17 device was about 1.5, if a patient comes in as a -1 and
18 gets the 0.25, their most likely expectation is that
19 they're going to be overcorrected by half a diopter. Of
20 course, there's quite a reasonable chance that they will be
21 overcorrected by quite a bit more. Do we want to include
22 something in the labeling about the possibility of
23 overcorrection in low levels of myopia? I just offer that.
24 If it dies, I'll die with it.

25 DR. McCULLEY: Other comments?

1 DR. SUGAR: I think that's very reasonable.

2 DR. McCULLEY: Okay. Any other discussion one
3 way or the other?

4 (No response.)

5 DR. McCULLEY: So basically what you're saying
6 is that there should be a warning in the labeling that a -1
7 or thereabouts with the 0.25 ring has a risk of being 0.5
8 diopters or more overcorrected.

9 DR. BULLIMORE: Yes.

10 DR. McCULLEY: That would be in the labeling.
11 Is there agreement?

12 DR. SUGAR: Yes.

13 DR. McCULLEY: Disagreement? Dr. Macsai.

14 DR. MACSAI: Do we have data to suggest that
15 "or more"?

16 DR. BULLIMORE: The main change with the 0.25
17 ring, which presumably would be the one you'd put in a -1
18 myope, was 1.5 diopters. So the average change is going to
19 take the patient to +0.50. So there's a 50 percent chance
20 they're going to be that side of the line, and a 50 percent
21 chance they'll be on the other side.

22 DR. SUGAR: The standard deviation is 0.52.

23 DR. MACSAI: Okay. Thank you.

24 DR. McCULLEY: Other discussion on that point?

25 Karen.

1 DR. BANDEEN-ROCHE: Dr. Bandeen-Roche. It
2 might make it clearer just to state that the mean
3 correction was -1.48. Therefore, there is some reasonable
4 chance of overcorrection in people with myopia much less
5 than that.

6 DR. McCULLEY: Okay. Any further comments?

7 DR. BULLIMORE: I live.

8 DR. McCULLEY: Despite your violent agreement
9 with yourself? If you don't get into a violent
10 disagreement with yourself. Actually, you used to have a
11 full head of hair, right?

12 (Laughter.)

13 DR. McCULLEY: Okay. We've done 5. Wait.
14 We're not through with 4? Oh, 5. Yes, you're right.

15 "Do the safety and effectiveness outcomes
16 support approval of 0.25, 0.30, and 0.35 ICRS? Is distinct
17 labeling warranted for any one of the three proposed ring
18 thicknesses?"

19 We covered much of that before, but there are
20 some things we didn't.

21 Dr. Grimmett?

22 DR. GRIMMETT: Dr. Grimmett. Since we lumped 2
23 and 5 together, I held my comment on 2 until we got to 5
24 here.

25 (Laughter.)

1 DR. GRIMMETT: I want to reiterate what I
2 already commented on, but just to reemphasize it. It's my
3 belief that both in the patient booklet and in the
4 physician information, if it's not already done, which it
5 may be, that the frequency and magnitude data reflect often
6 and always categories combined, as well as moderate and
7 severe categories combined. I think it downplays visual
8 symptoms to only report the severe and always categories.
9 I think that's very important.

10 DR. McCULLEY: My interpretation of the data
11 was that there was significant increase in symptoms that
12 were of increasing magnitude going from 0.25 to 0.35, but
13 they were real at 0.25.

14 DR. MACSAI: Absolutely.

15 DR. GRIMMETT: Because if you combine those
16 figures, up to approximately 1 in 5 can have visual
17 symptoms that certainly reflect a new set of higher order
18 visual aberrations occurring in those patients. I believe
19 that's important to know.

20 DR. MACSAI: Higher, 37 percent.

21 DR. ROSENTHAL: I don't think it's fair to
22 combine them. I think it's fair to put them both in, but
23 I'm not sure it's fair to combine them, just like you
24 wouldn't combine zero and 1 to make a half. You put in
25 moderate and severe. You don't leave them out. But to say

1 moderate or severe, I think it's fairer to give the patient
2 and the doctor some idea of these very difficult terms.
3 What is moderate and what is severe for Dr. Macsai may not
4 be moderate and severe for me. But at least it's better
5 than combining moderate and severe.

6 DR. McCULLEY: So it just means more data.

7 DR. GRIMMETT: I completely agree. More
8 information is better than just listing the two columns, so
9 people can have a balanced appraisal of what's going on.

10 DR. McCULLEY: Karen?

11 DR. BANDEEN-ROCHE: Dr. Bandeen-Roche. This is
12 a global issue not connected to one of the three sizes. It
13 might be redundant. It's the issue of those patients who I
14 would say fail altogether, to make sure that there is some
15 labeling that is very clear about that, which means
16 patients whose surgeries had to be aborted, just didn't
17 succeed in the first place, and patients who end up
18 explanting ultimately being reported as a percentage for
19 which that happens and not being mixed up about do visual
20 symptoms include these people or not. We just have to be
21 very clear, and I'd like to see an up front statement of
22 this may not work for you at all, the percentages, thus and
23 such, being clear about what that means.

24 DR. McCULLEY: So we need a strong statement in
25 labeling about the 4.7 percent incidence of explantation,

1 and that those patients represent additional problems above
2 and beyond those reported in the symptom percentages.

3 DR. BANDEEN-ROCHE: Right, and a very small
4 percentage for whom the surgery never succeeded in the
5 first place, 5 out of 454, something like that.

6 DR. McCULLEY: I don't remember the exact
7 number, but okay.

8 Dr. Wang?

9 DR. WANG: Ming Wang. The question is, is
10 distinct labeling needed just for the 0.35 or not, or is it
11 sufficient just providing that table including all three?
12 My feeling is that 0.35 is somewhat different animal, and
13 we may need to put a sentence for the patient's sake,
14 rather than just present a table neutrally, all three
15 sizes, but have one sentence or two warning about this,
16 because it is doubling in the visual symptoms of 0.35
17 compared with the other two. It's non-linear.

18 DR. McCULLEY: The sentiment has been there,
19 and that was well stated, that the additional problems with
20 0.35 should not be just buried in the table for a person to
21 sort out and realize that there is significant increase at
22 the 0.35 ring.

23 DR. SUGAR: That was the Grimmatt proposal.

24 DR. McCULLEY: Right, Grimmatt-Wang.

25 Dr. Macsai?

1 DR. MACSAI: I would also like to suggest that,
2 as far as the effectiveness outcomes, that not only the
3 plus or minus 1 diopter of intended correction be included,
4 but also the plus or minus 0.5 diopter, because we're
5 talking about such a low range of myopia in these patients.
6 If you're -1, plus or minus 1 doesn't mean much.

7 DR. McCULLEY: You're talking about in the
8 predictability and stability label warning?

9 DR. MACSAI: Not warning. Labeling.

10 DR. McCULLEY: Okay, labeling.

11 DR. GRIMMETT: This is Dr. Grimmett. That may
12 already be done in the labeling that they presented, both
13 the data for plus or minus 0.5 and plus or minus 1.

14 DR. MACSAI: Right.

15 DR. McCULLEY: But I think the point is, rather
16 than a naive patient who is a -1 not taking that into
17 account, that there needs to be added wording for that
18 person to be sufficiently aware.

19 Any other comments on Question 5?

20 (No response.)

21 DR. McCULLEY: Question 6. "Is the current
22 data in the Removal Cohort sufficient to support
23 reversibility claim? If not, what is the minimum number of
24 eyes and the minimum length of follow-up that you recommend
25 for this assessment?"

1 Dr. Van Meter, do you want to start that?

2 DR. VAN METER: Woodford Van Meter. I think
3 the data support the reversibility claim. I don't think
4 that the adjustability claim can be supported, but I think
5 the reversibility claim is adequately supported by the
6 data.

7 DR. McCULLEY: Okay, we're just talking about
8 reversibility now. We'll come to adjustability the next
9 point. Reversibility is on the table.

10 Karen?

11 DR. BANDEEN-ROCHE: Karen Bandeen-Roche.
12 Whether reversibility is supported depends on what we
13 consider to be adequate reversibility. I computed some
14 rough lower confidence bounds based on the sample of 21,
15 and then assuming the same percentages in 34. Just to very
16 briefly summarize, there is very weak evidence, even at 21
17 over 21, that reversibility for BSCVA is higher than 94
18 percent, very strong evidence that it's higher than 70
19 percent, and then a whole range in-between. So, yes, it
20 really depends on what we think is an acceptable percentage
21 re-achieving their preop status.

22 DR. McCULLEY: Right. I'm not sure about that.
23 Is there some way that you could phrase this to take the
24 reality as presented under consideration, rather than an
25 all or nothing?

1 DR. BANDEEN-ROCHE: Well, the 95 percent
2 confidence bounds that I -- these are exact lower
3 confidence bounds that I roughly computed, back of the
4 envelope sort of a calculation, were 0.87 for best
5 corrected visual acuity, 0.79 for post-removal
6 predictability, 0.73 for manifest refraction cylinder.
7 Again, that was at 21. I realize there's been some
8 additional data since then, but if you extrapolated that
9 out to 34, those quantities would be 0.91, 0.82, and 0.75,
10 respectively.

11 DR. McCULLEY: Could you put a recommendation
12 into words that would be appropriate for labeling?

13 DR. BANDEEN-ROCHE: I think that it's beyond my
14 purview to suggest clinically what is an acceptable target
15 for reversibility. So I think that has to precede a
16 recommendation on my part.

17 DR. MACSAI: Dr. McCulley?

18 DR. McCULLEY: Dr. Macsai.

19 DR. MACSAI: I just want to raise the question
20 that there appears to be reversibility of visual symptoms
21 and vision, but I didn't see any data about reversibility
22 of changes in sensation, endothelial cells, contrast
23 sensitivity in those patients.

24 DR. McCULLEY: Dr. Grimmett brought that up in
25 his review. The symptoms were not completely reversible,

1 nor was the refractive error completely reversible. But,
2 then again, the point that Karen is bringing up is what is
3 acceptable? What kind of guideline would we have for
4 acceptability stating that we agree it's reversible? It is
5 reversible to a degree. Is it reversible to a sufficient
6 degree that that statement should not be qualified in some
7 way?

8 Dr. Van Meter, then Dr. Sugar.

9 DR. VAN METER: The recent data that was
10 submitted at the end of December with the final packet --
11 and perhaps we can do this just by reiterating some
12 statements from that final report. No subject lost best
13 corrected acuity. Ninety-six percent of patients returned
14 to plus or minus 1 diopter of their manifest refraction.
15 The one that was outside that particular bracket was
16 better. Stability appears to have been shown at three
17 months, and all of the patients that had cylinder appeared
18 to return when the implants were taken out. I think you
19 might say that the loss of glare symptoms or endothelial
20 cell counts, we really don't have any data on that, and I
21 think you can say that is not known. But I believe we have
22 reasonably good data.

23 DR. McCULLEY: I thought we had the data on
24 symptoms but we didn't on endothelium or corneal sensation,
25 slit lamp appearance and so forth.

1 DR. SUGAR: This is Sugar. For the symptoms,
2 there were two patients. All patients had zero to 1 for
3 all of the subjective symptoms, except for two patients,
4 who had three symptoms -- severe difficulty with night
5 driving, double images, and fluctuating vision. I think
6 that it's appropriate to have a statement that it is
7 reversible and list the parameters under which it is so for
8 all patients 20/20 or better best spectacle-corrected
9 visual acuity. Two of 21, or whatever the denominator
10 becomes, had symptoms that were moderate or severe. These
11 are two patients with those three symptoms.

12 DR. McCULLEY: So just put data in the
13 labeling.

14 DR. SUGAR: And just put the specifics there.

15 DR. McCULLEY: Karen?

16 DR. BANDEEN-ROCHE: As long as it's very clear
17 the amount of data that contributed to this, I think I'd
18 like to make it even stronger to get people to stop and
19 think about what strength of evidence that means. Now,
20 maybe for the physician's booklet a lower confidence bound
21 would do it. I still think it's worth at some point FDA
22 and all of us thinking about what is acceptable. But in
23 the meantime, there needs to be some measure whether 21 out
24 of -- well, I'm sorry, it's now 28 out of 28, whether
25 that's very powerful evidence or very weak evidence, or

1 something in between. I would suggest put the confidence
2 bounds in there for the physicians.

3 DR. ROSENTHAL: This is Dr. Rosenthal. You
4 know, it becomes a philosophical issue what is reversible.
5 If patients have unrealistic expectations that this can be
6 put in and just removed with complete assurance that there
7 will be no problems, that would be -- I would consider that
8 a reversible situation. But I don't think that's
9 absolutely the case. Hence, one could, instead of using
10 the word -- I'm just throwing this up for panel's
11 consideration -- instead of using the word "reversible,"
12 use some other word saying that it can be removed, and
13 therefore there's a 96 percent chance that you will get
14 back to normal, but there's a chance that there's going to
15 be some serious problems. That's a realistic evaluation,
16 rather than calling something reversible.

17 DR. McCULLEY: I think that makes sense.

18 Renee?

19 DR. MIDDLETON: As a consumer, if you told me
20 that something was reversible, I would assume that you
21 meant that I would go back to where I was, your former
22 statement. That would be my assumption. So if that's not
23 the picture that you want to leave with the consumer, then
24 perhaps some other wording is appropriate, or more language
25 to clarify what you mean by reversibility.

1 DR. McCULLEY: Or a 96 percent probability of
2 reversibility. Something.

3 DR. MACSAI: Removable.

4 DR. McCULLEY: Well, it is removable, but it's
5 largely reversible. We need a better word because, as
6 indicated, in our human mind, when we hear reversible, we
7 assume return to baseline. It has a probability but not an
8 absolute certainty of doing that. How do we want to have
9 that -- how can that be honestly, effectively conveyed?
10 Marcia?

11 DR. YAROSS: This is Marcia Yaross. What I
12 would suggest, instead of struggling over the semantics, is
13 just have a factual description of what were the results
14 observed and basically state that in a group of X patients
15 who had this device removed, X percent recovered their
16 vision to the previous levels. I'd also point out as a
17 benchmark in terms of effectiveness numbers in the ranges
18 of 84 to 88 percent achieving an effect have been
19 considered effective by this panel in the past. So 100
20 percent is not necessary to determine that something is
21 effective.

22 DR. McCULLEY: As long as it's done so that
23 reversibility is not what stands out and there's a footnote
24 that might be ignored. So if we take the word out, you're
25 suggesting to --

1 DR. YAROSS: I'm suggesting a factual
2 description of what happened to those subjects that had the
3 device removed.

4 DR. ROSENTHAL: This is Dr. Rosenthal. There's
5 no question that that's what would be done. The question
6 is that the company is going to want to say it is
7 reversible. That is a very powerful claim that they will
8 want to make in their marketing of the product, as compared
9 to the other modalities of refractive surgical
10 intervention. So I think it's a very important issue what
11 we allow them to say and what we allow them not to say.
12 Certainly in the labeling we would spell it all out, but
13 are you going to let them say it's reversible?

14 DR. YAROSS: Marcia Yaross again. That's why I
15 was trying to propose that something in the 80 percent is
16 what's been our criterion. So if we agree on what
17 reversibility should mean, then you don't need to be 100
18 percent to say you're effective in some regard, based on
19 the history of this panel.

20 DR. McCULLEY: Just one point here. We're
21 getting into the world of lay people and their definitions
22 of words, not the scientific community, with an arbitrary
23 softening of a definition. I think Renee's point is very
24 well taken, and I think that the population, if they hear
25 reversible, they're going to assume it is absolutely

1 reversible. This is not only caveat emptor but caveat
2 venditor. It's who's buying and who's selling. So
3 informed consent and the risk of someone selling something
4 that is reversible that ends up not being truly reversible,
5 the patient is going to be unhappy because they use the
6 word reversible as they would use it in their normal daily
7 lives.

8 Eve, I think you had something to say.

9 DR. HIGGINBOTHAM: I was going to offer
10 "somewhat reversible."

11 (Laughter.)

12 DR. HIGGINBOTHAM: The suggestion is that this
13 is only up to three months of data, and I just question
14 whether or not that's long enough. I know in glaucoma we
15 like to see at least six months. Because the refractive
16 stability seems to get better up to 12 months, I don't know
17 if three months is really enough.

18 DR. McCULLEY: Again, we need a thesaurus to
19 come up with a better word.

20 Renee?

21 DR. MIDDLETON: I think we can see that there
22 is a difference between misleading and outright lying, but
23 they're not doing either. I think if they simply said -- I
24 would be comfortable if you said reversibility. I don't
25 have a problem with that as a consumer, as long as you tell

1 me in what framework you're referring to reversibility as,
2 and that's up front to the patient. Then I wouldn't feel
3 like they were misleading me. With the data that's there,
4 I don't think they're intending to mislead, unless you just
5 said reversibility and leave all that information out.

6 DR. McCULLEY: Okay. Dr. Pulido?

7 DR. PULIDO: Dr. Middleton, I agree on a lot of
8 things with you, but on this one I disagree because people
9 are going to fixate on that word and forget the rest. I
10 think we still should struggle with this and make it a more
11 equitable term than reversible.

12 DR. McCULLEY: Dr. Wang?

13 DR. WANG: Ming Wang. I would like to suggest
14 maybe language something like this: "This device can be
15 removed and the majority of symptoms are reversible," and
16 provide a table, or "is largely reversible." Some
17 language, immediately followed with a table.

18 DR. McCULLEY: Dr. Eydelman?

19 DR. EYDELMAN: If the panel chooses to come up
20 with different language, we're still going to need the
21 recommendation from the panel as to what would constitute a
22 reversible claim that the sponsor can then come back with
23 enough data, because that claim I have a feeling will come
24 back. So even if we choose not to call it reversible at
25 present, can the panel then make a recommendation as to how

1 much data would be needed for the sponsor to be able to
2 make that claim?

3 DR. McCULLEY: Reversible as has been suggested
4 is an all or nothing in the mind of the public, and it's
5 going to have to be clarified if it's anything less than
6 absolutely reversible. Just on that point, is there
7 agreement to that point, to use the word reversible as it
8 is used? We can't make this a scientific term that we put
9 under a guideline. It's used in common, everyday language,
10 and I think we understand how it's understood.

11 Dr. Sugar?

12 DR. SUGAR: I think this affects how the
13 company markets their product, and I think people market
14 refractive surgery as in most circumstances you will no
15 longer need glasses for distance vision, and I think it's
16 reasonable to say that this is reversible in most
17 circumstances, because I believe that it is.

18 DR. McCULLEY: Okay, as long as it's qualified,
19 and then how much qualification, how specific should the
20 qualification be. It just has to be qualified and not lost
21 as a footnote.

22 Karen?

23 DR. BANDEEN-ROCHE: Dr. Bandeen-Roche. Might
24 you also want to put a comment in there for the lay person
25 that, however, only thus and such many patients have been

1 evaluated for only three months of follow-up and more data
2 will solidify -- something like that, you know?

3 DR. McCULLEY: That makes sense as well.

4 Dr. Matoba?

5 DR. MATOBA: Alice Matoba. I agree with Dr.
6 Sugar regarding his comments about the reversibility.
7 Also, I wanted to suggest that maybe we consider adding a
8 statement that the average residence of the implants was
9 about 10 months and that we don't have information
10 regarding potential reversibility following long-term
11 residence in the cornea.

12 DR. McCULLEY: Good point.

13 Dr. Pulido?

14 DR. PULIDO: What about the endothelial cell
15 loss? That is not reversible.

16 DR. McCULLEY: We don't know. We don't know
17 whether it's significant or not yet, so we really can't
18 address that. We don't know that it's there yet.

19 Dr. Macsai?

20 DR. MACSAI: It seems to me that you and Dr.
21 Pulido and Dr. Higginbotham and Dr. Sugar have just
22 answered this question. We do not know if the endothelial
23 cell change is reversible.

24 DR. McCULLEY: Well, we don't know if it's real
25 yet.

1 DR. MACSAI: Well, that's right. So we don't
2 know, and we only have three-month data. I think
3 "removable" is a great word. There's no eraser at the end
4 of a laser, but this is removable. Take it out. I mean,
5 that's a big selling point.

6 DR. McCULLEY: Not nearly as big a selling
7 point as reversible.

8 (Laughter.)

9 DR. MACSAI: But we don't know.

10 DR. ROSENTHAL: This is Dr. Rosenthal. The
11 company is not going to put "this is largely reversible,"
12 or "under circumstances mostly reversible." The company is
13 going to want to say it's reversible. I must say, I have a
14 little problem with it, because I do genuinely feel that it
15 means it's reversible, you can erase it completely and go
16 back to where you were, and I'm not sure that's totally
17 true.

18 DR. McCULLEY: Well, we know it's not totally
19 true, and the panel is in agreement with you.

20 DR. MACSAI: We don't know, so we should not
21 say.

22 DR. YAROSS: Mr. Chairman, can I offer a
23 possible clarification? This is Marcia Yaross. Or an
24 alternative? And that might be to say that the refractive
25 effect is reversible. It has a reversible refractive

1 effect, and that leaves open the fact that some of these
2 other safety issues, it is not yet known.

3 DR. McCULLEY: To me, that obfuscates and it
4 implies other things, so I don't think that's a good
5 solution to it. I'd like to hear one. It's out there.

6 Renee?

7 DR. MIDDLETON: I think I liked what Dr. Sugar
8 proposed, but I hear Dr. Rosenthal saying that the sponsors
9 are not going to accept that.

10 DR. ROSENTHAL: This is Rosenthal. I didn't
11 say they weren't going to accept it. What I said was, as
12 you know, if you're a sponsor and you've got a product, you
13 want to try to sell it using the best buzz words you can
14 use. That's my only point. I don't want the company to be
15 in any doubt about what the panel feels they should be
16 allowed to say. That's what my point is.

17 DR. MIDDLETON: Can I ask Dr. Sugar to repeat
18 the language that he would recommend? I believe it was
19 him.

20 DR. SUGAR: I don't think that Dr. Rosenthal is
21 going to leave the sponsor in any doubt.

22 (Laughter.)

23 DR. MIDDLETON: But what was the language that
24 you proposed?

25 DR. SUGAR: Reversible in most circumstances.

1 DR. MIDDLETON: Reversible in most
2 circumstances. I would be comfortable with that as a
3 consumer. But if you have -- I'm speaking from the role of
4 a consumer. If you have serious concerns about that
5 language and the sponsors would refuse to use that
6 language, I would say that they can't then say it's
7 reversible if they're going to be uncomfortable with
8 proposed language similar to that.

9 DR. McCULLEY: How about "reversible in most
10 but not all circumstances"?

11 DR. MIDDLETON: That's similar.

12 DR. ROSENTHAL: I get the gist, and I think the
13 big issue will be with our compliance section if there are
14 issues relating to complaints about ultimate claims.

15 DR. McCULLEY: Dr. Eydelman?

16 DR. ROSENTHAL: But I understand, and I think
17 the company understands as well.

18 DR. EYDELMAN: I know we're all trying to
19 discuss the best term, but if we can just step back, I
20 wanted to make sure that I'm clear on the panel's feeling
21 about the appropriateness of stability between -- only two
22 measures were taken, one and three months. This is
23 different than all other stabilities. It's different than
24 our guidance, et cetera. So before we conclude with the
25 best term, I wanted to make sure I understand what the

1 feeling is about that.

2 DR. McCULLEY: Good point. I mean, we're
3 cutting a corner on our guidelines. The definition of
4 stability is within a diopter three months apart.

5 DR. EYDELMAN: At least two measurements at
6 three months.

7 DR. MACSAI: We need six-month data.

8 DR. EYDELMAN: Granted that's not the same
9 device, but --

10 DR. SUGAR: Preoperatively at three months
11 post-removal, 100 percent of patients were within -- or 96
12 percent of patients were within plus or minus 1 diopter.
13 That meets the guideline. And the one patient who wasn't
14 was actually less myopic than they started.

15 DR. EYDELMAN: No, but you had a big dip in the
16 curve. That definition came into assuming some kind of a
17 linear curve. So you need to look at the points post-
18 removal.

19 DR. McCULLEY: Post-removal. One month and
20 three month were the time points that were measured.

21 Dr. Bullimore?

22 DR. BULLIMORE: I think there's a danger here
23 that the sponsor is put in an untenable position, and I
24 want to speak out on their behalf. They produced a device
25 that seems to be reasonably safe, reasonably effective, and

1 in 5 percent of people it needs to be explanted. We're
2 arguing over those 5 percent as to whether it should be
3 reversible. What do we want them to do? Do we want them
4 to produce something that fails in 10 percent so we have
5 more data to go on?

6 (Laughter.)

7 DR. BULLIMORE: We're trying to make the best
8 of the data that we have, and to sort of say we need to go
9 out further and argue about the stability I think is really
10 disingenuous on the part of this panel, and I think we
11 should just put closure on this item and move on.

12 DR. McCULLEY: I think it's very important with
13 what we do to be certain that patients are adequately
14 informed. If we tell a patient that they are getting a
15 procedure done to them that is reversible, then we do them
16 a disservice if it is not 100 percent reversible. That's
17 not to say it's not an acceptable or good procedure. It's
18 that the patient needs to be effectively informed of what
19 they can expect. So I don't think that's being
20 disingenuous.

21 DR. BULLIMORE: No, but I think we're getting
22 outside of what is reasonably expected of this panel and
23 what is reasonably expected on the part of the agency of
24 this panel. The agency has compliance people. This seems
25 to be in their court. Whether this group of people can

1 agree on what's reversible and what's not seems to be a
2 moot point. The advertising, the promotion of these
3 devices, we're basing it on a small number of people. We
4 have the data here. Let the FDA do what they want with it.

5 DR. McCULLEY: I think the FDA is asking us for
6 our opinion.

7 DR. BULLIMORE: I think they've heard enough
8 opinion. We can sit here until 7:00 discussing this and
9 not get any closer.

10 DR. McCULLEY: I think we'll get closer. We're
11 not done yet.

12 Alice?

13 DR. MATOBA: Alice Matoba. I wanted to say
14 that I don't think anything in medicine is 100 percent. If
15 you're going to hold people to 100 percent reversibility as
16 a standard, no one will be able to achieve that. I think
17 they've come close enough to showing reversibility that
18 they can use that word in some way.

19 DR. McCULLEY: I don't have a problem with
20 that. It just has to be qualified, and it has to be
21 qualified so the lay person will understand the
22 qualification and not skirt over it. I think Joel's
23 reversible in -- what was it?

24 DR. SUGAR: In most circumstances.

25 DR. McCULLEY: In most but not all, and then

1 that underlines that it's not all for the human mind, not
2 to skip over and not pay attention. I would be comfortable
3 with that statement.

4 Dr. Grimmer?

5 DR. GRIMMETT: We said this before, but I'll
6 throw it in because it has to do with reversibility.
7 Regarding claims that patients may elect to have an
8 alternative refractive surgical procedure performed,
9 there's no data at this point to support that claim. That
10 was an Amendment A attachment to page 25. So I think those
11 claims need to be removed.

12 DR. McCULLEY: Right. Okay. Another point
13 that I would agree with as well. We don't have data to
14 make this statement one way or the other.

15 Dr. Sugar?

16 DR. SUGAR: Just to get back to Malvina's
17 point, in Volume 8, Section I, page 11, they have from day
18 14 to month 1, 94 percent changed less than 1 diopter, and
19 from 1 month to 3, 24 out of 24 changed less than 1
20 diopter. So I think that there is evidence that there's a
21 plateau.

22 DR. EYDELMAN: Right, but from day 14, I
23 believe there was something about 14 eyes.

24 DR. SUGAR: Eighteen eyes.

25 DR. McCULLEY: Dr. Pulido?

1 DR. PULIDO: I just would like to ask if maybe
2 you, Mr. Chairman, would want to ask the panel if we feel
3 comfortable with Dr. Sugar's recommendation so we can move
4 on.

5 DR. McCULLEY: Okay. Dr. Sugar, would you
6 restate your recommendation? And we'll take a panel poll.

7 DR. SUGAR: Sure. This procedure in the
8 labeling could be stated as reversible in most but not all
9 circumstances.

10 DR. McCULLEY: Is there comfort with that
11 statement?

12 DR. PULIDO: How about removable and reversible
13 in most circumstances?

14 DR. SUGAR: It's always removable.

15 DR. MACSAI: Out to 10 months.

16 DR. McCULLEY: We have to have the
17 qualification that this is based on data that is on a small
18 number of patients with short-term follow-up. So it has to
19 be qualified. We go along with the statement as though it
20 appears that it's okay, but we don't have long enough data
21 or enough patients yet.

22 MS. LOCHNER: Excuse me, Dr. McCulley. There
23 was one question that was raised earlier that I really
24 haven't heard the panel discuss, and it speaks to this
25 question of reversibility. I think it was Dr. Bandeen-

1 Roche who said that the number of eyes that were evaluated,
2 whether it's the 21 or the 34, gives you a lower confidence
3 interval on this reversibility claim that is at some
4 points, I think you said, as low as 70 percent. I guess I
5 would like to hear what the panel -- that's sort of setting
6 a threshold of an absolute amount of equivalence to support
7 a claim such as this, and I hadn't really heard whether
8 that was felt to be adequate, the fact that the sample size
9 you have only gives you assurance down to about 70 percent.

10 Now, I don't suggest that the sponsor have more
11 failures, but there certainly could be more eyes enrolled
12 and the sample size itself increased so that you have
13 better confidence in the data. And I'm sure I didn't get
14 the absolute percentages correct.

15 DR. BANDEEN-ROCHE: This is Karen Bandeen-
16 Roche. I just wanted to clarify that wasn't for all
17 outcomes. That was for, I think, manifest stability at
18 plus or minus 0.5. It was whatever one had the lowest
19 achieved percentage.

20 DR. BULLIMORE: This is Dr. Bullimore. Just a
21 point of clarification. What is the 0.71? Is that a
22 proportion?

23 DR. BANDEEN-ROCHE: It's a lower confidence
24 bound on the proportion reversible.

25 DR. BULLIMORE: Okay. So the sample or the --

1 DR. BANDEEN-ROCHE: The actual proportion in
2 that case was 0.88.

3 DR. BULLIMORE: And the lower 95 percent --

4 DR. BANDEEN-ROCHE: The 95 percent confidence
5 bound is 0.75, given --

6 DR. BULLIMORE: So it's a proportion.

7 DR. BANDEEN-ROCHE: That's the unit.

8 DR. BULLIMORE: That's the unit. That's what I
9 want to know. Okay.

10 DR. PULIDO: And therefore it's still
11 reversible in most cases.

12 Can we move on?

13 DR. McCULLEY: I don't think we want to set any
14 mark at a 70 percent confidence that would come back to
15 haunt us.

16 MS. LOCHNER: No, I just think hearing your
17 opinion on that was all we really wanted.

18 DR. McCULLEY: We're trying to hedge, and I
19 think we're trying to accept that sponsor wants the word
20 somewhere in there reversible, and it seems to be a
21 reasonable request, but it has to be qualified.

22 DR. BANDEEN-ROCHE: If I may, Karen Bandeen-
23 Roche. I just want to make sure I'm being absolutely clear
24 because Dr. Grimmett thinks that I have not been.

25 (Laughter.)

1 DR. BANDEEN-ROCHE: I thank you for this.

2 We're talking about a 95 percent confidence
3 bound. The value of that bound is 0.75.

4 DR. McCULLEY: I knew I didn't understand.

5 Any other comments on that?

6 (No response.)

7 DR. McCULLEY: Thank you very much. So are we
8 okay on this point? Is the FDA okay on this point? Are we
9 tired of this point?

10 Question 7. "Is the current data on exchange
11 procedures sufficient to support claim for adjustability of
12 refractive effect? If not, what is the minimum number of
13 eyes and the minimum length of follow-up that you recommend
14 for this assessment?"

15 Let's take this question in two parts. I think
16 there was fairly uniform agreement among the primary
17 reviewers. Is the current data on exchange procedures
18 sufficient to support the claim of adjustability?

19 DR. SUGAR: No.

20 DR. BULLIMORE: No.

21 DR. MACSAI: No.

22 DR. McCULLEY: Does anyone have any countering
23 argument to that stated opinion?

24 (No response.)

25 DR. McCULLEY: Okay, now to the more difficult

1 part. If not, which we're not, what is the minimum number
2 of eyes -- and the record should show that that was a
3 unanimous no -- what is the minimum number of eyes and the
4 minimum length of follow-up that you'd recommend for this
5 assessment?

6 Dr. Sugar.

7 DR. SUGAR: I don't think we can come up with a
8 number because the predictability was so low in the few
9 patients that they did that we don't have any idea of what
10 the efficacy is going to be. We don't know where to grab a
11 number.

12 DR. McCULLEY: We respectfully wish not to get
13 put in that corner.

14 Any other comments?

15 (No response.)

16 DR. McCULLEY: Question 8. "The sponsor would
17 like to make a claim of 'enhanced visual performance' in
18 their labeling. Do you feel that the data in this PMA
19 support this claim?"

20 DR. SUGAR: No.

21 DR. MACSAI: No.

22 DR. McCULLEY: The sentiment being stated is
23 no. Is there dissent to that sentiment?

24 (No response.)

25 DR. McCULLEY: Unanimous no, that it does not

1 support the claim.

2 Does the FDA need further comment from panel?

3 There's been a lot of comment that has all been in one
4 direction. I realize just as you may have fights in the
5 future relative to reversibility, which is going to be a
6 major marketing issue, that you may have future fights
7 about this issue. Do you have sufficient information
8 relative to why the panel feels the way it does, or do you
9 need more?

10 DR. EYDELMAN: I believe we have enough.

11 DR. ROSENTHAL: This is Dr. Rosenthal. May I
12 add that if there really is a conflict between the sponsor
13 and the agency, the sponsor can make a proposal to do a
14 study to which the claim that they wish to promote will be
15 acceptable and the study can be done. So there's no
16 impasse as long as the company cooperates.

17 DR. McCULLEY: Those of you who looked at the
18 physician training, we didn't ask the sponsor about what
19 they were putting as a recommended upper limit for the
20 suction ring being applied with a pressure in the range of
21 80. In the physician instruction portion, was a limit
22 stated and recommended? Because I can imagine if a person
23 is not adequately warned and not watching the clock, that
24 an eye could be highly -- five minutes?

25 DR. MACSAI: Five minutes?

1 DR. McCULLEY: Five minutes?

2 Dr. Higginbotham, would you like to comment?

3 DR. MACSAI: What's that based on?

4 DR. HIGGINBOTHAM: This is Dr. Higginbotham.

5 Five minutes at 80 millimeters of mercury is certainly a
6 long time for any eye, particularly when you've only
7 assessed these eyes by super threshold testing. So I would
8 have some particular concerns about that duration.
9 Certainly the literature suggests four minutes is the
10 absolute maximum, but something less than that even would
11 be more acceptable. So if the ideal time and what your
12 experience has dictated in the procedures you've done to
13 date mostly being a minute and a half, I would limit it to
14 a minute and a half as a goal.

15 DR. McCULLEY: Based on experience with LASIK,
16 I think one can go past a minute and a half before
17 panicking. But, boy, I would never go to five minutes. I
18 think this is going to be important for labeling and
19 instruction of physicians, that they know the clock has to
20 run. I'm not even sure that docs doing some of these other
21 procedures that pressure eyes are effectively watching, but
22 we haven't had the opportunity to make the comments to
23 them.

24 What would you consider -- a minute and a half
25 is too short.

1 DR. HIGGINBOTHAM: Well, I guess I'll yield to
2 those clinicians that are currently doing LASIK in the same
3 population. So I would yield to your wisdom, Dr. McCulley.

4 DR. McCULLEY: Well, it's seat-of-the-pants
5 rather than wisdom. But I start to worry big time at two
6 and a half, and I wouldn't go past three. But that's not
7 based on any good, solid data, unless the person from whom
8 I adopted that, and I couldn't tell you who it was, had
9 good, solid data.

10 Dr. Wang?

11 DR. WANG: Ming Wang. Accompanying all LASIK
12 surgery practice, a two-minute threshold was based on case
13 reporting in one of the international conferences, LASIK
14 beyond two minutes. Two minutes is the common upper limit
15 in LASIK.

16 DR. McCULLEY: So you use two minutes.

17 DR. WANG: Yes.

18 DR. McCULLEY: Based on one case?

19 DR. WANG: Yes, based on one case, I believe in
20 South America.

21 DR. McCULLEY: Okay.

22 DR. HIGGINBOTHAM: I missed that case report.
23 Sorry.

24 DR. McCULLEY: Does anyone else have any
25 comments about this? I think we would be uncomfortable

1 with five minutes. I think that this is going to be
2 extremely important for the company to pass on to the
3 surgeon.

4 Dr. Macsai?

5 DR. MACSAI: This is a technical question. You
6 put the suction on and cut clockwise and counter-clockwise.
7 If you extend beyond one and a half or two minutes, can't
8 you take it off, put it back on, then do the other half?
9 Why not?

10 DR. HIGGINBOTHAM: Dr. Higginbotham. I guess
11 my comment interpreted the five minutes as being continuous
12 application of vacuum for five minutes. But certainly if
13 you allow the nerve to perfuse intermittently, as we do
14 when we're massaging eyes post-retrobulbar, I think that
15 would be sufficient. But again, we're going on more
16 anecdotal information at this point.

17 DR. McCULLEY: But it would be better to have
18 one and a half minutes with a break, and then one and a
19 half minutes again, rather than three.

20 DR. HIGGINBOTHAM: Agreed.

21 DR. McCULLEY: So I think direction to FDA to
22 pay very close attention to working out something with the
23 company for guidelines relative to the length of time for
24 continuous suction, and if that is exceeded, how long the
25 eye should be allowed to perfuse prior to reestablishing

1 suction.

2 Dr. Pulido?

3 DR. PULIDO: Can I ask a separate question at
4 this time?

5 DR. McCULLEY: Yes.

6 DR. PULIDO: I'm still very concerned about the
7 marked difference in results between different sites, and I
8 would like to see something in labeling, a statement in the
9 labeling to the effect that there can be marked variability
10 in results in explantation rates between different
11 surgeons.

12 DR. McCULLEY: Boy, I'll tell you, that opens
13 up, I think, from some different vantage points that we
14 won't talk about in public forums, some real potential cans
15 of worms.

16 Dr. Higginbotham?

17 DR. HIGGINBOTHAM: I would have some concerns
18 about putting that in, just because the dialogue that
19 occurs between a patient and physician and the decision
20 that follows to take a device out can be independent of the
21 technical outcome of that procedure but maybe some other
22 subjective reasons that are more patient based. So I would
23 not tie the sponsor's hands with that kind of statement.

24 DR. McCULLEY: Dr. Bullimore?

25 DR. BULLIMORE: I have another suggestion to

1 make.

2 DR. McCULLEY: Along these lines? Let's be
3 sure this point is resolved. Is it still on this point?

4 DR. BULLIMORE: No, it's not on this point.

5 DR. McCULLEY: Okay. Any other comment on this
6 point? Dr. Macsai?

7 DR. MACSAI: I agree with Dr. Pulido's
8 observation. However, it seems to me inherent in any
9 surgical procedure that there's going to be variability
10 between surgeon to surgeon, whether it's an intracorneal
11 ring, an intraocular lens, or a laser. There is going to
12 be variability, and I would ask that the agency determine
13 if this was a -- I assume from their bio, as we just
14 learned about, that this was not an issue.

15 DR. McCULLEY: Other comments on this issue?
16 Dr. Bullimore?

17 DR. BULLIMORE: It's not on this issue.

18 DR. McCULLEY: No, no. We're going to a new
19 issue now.

20 DR. BULLIMORE: Oh, I get a new issue. Oh,
21 goody.

22 In the proposed indication, the sponsor says
23 that refractive stability preoperatively is defined as a
24 manifest refraction change of 1 diopter or less for at
25 least six months prior to the preoperative examination.

1 That seems a very generous definition of stability compared
2 to what currently resides in the guidance document, and I
3 would propose that we adopt something more consistent with
4 other such standards.

5 DR. McCULLEY: Dr. Rosenthal?

6 DR. ROSENTHAL: Unfortunately, the definition
7 of preoperative stability has varied, and this sponsor
8 worked on the definition which was presented to you, and I
9 think the study was done with that definition and it would
10 have to be carried through.

11 DR. McCULLEY: I understand what you're saying,
12 and I understand what Dr. Bullimore is saying, and I had
13 the same concern.

14 DR. ROSENTHAL: I know.

15 DR. McCULLEY: It's back to the point of we
16 have now fine-tuned, we have moved forward, we now have
17 something we think is better, so now do we carry forward
18 something that wasn't? Yet that's how they did the study.
19 There must be some way to deal with both points
20 effectively. For instance, it could be stated that, in
21 effect, that's how the study was done that led to approval.
22 However, current FDA guidelines are 0.5 millimeter change
23 no greater than in the past year. So we bring into this
24 previous, now not-so-good decision where we are now.

25 DR. BULLIMORE: The sponsor can collect data

1 on, say, -1 to -20 diopters, and this panel might vote for
2 approval on -1 to -2. So just because a study has been
3 done under a different set of entry criteria which is
4 distinct from proposed indications for the public release
5 of the device, I would argue that we should adopt a more
6 conservative definition of preoperative stability.

7 DR. McCULLEY: I think in current practice a
8 diopter variation or close to it in a six-month period
9 would not represent a patient that was an appropriate
10 person for keratorefractive surgery. That's not adequate
11 stability in the real world now. This doesn't disadvantage
12 sponsor, can put both pieces of information in, hopefully.

13 Dr. Eydelman?

14 DR. EYDELMAN: I just had a response to Dr.
15 Bullimore's comment. When sponsor does a study on a large
16 refractive range and then we choose to limit it, it's
17 usually based on data. Unfortunately, we don't have any
18 data on a tighter preop stability.

19 DR. McCULLEY: So can we use my argument rather
20 than his?

21 DR. BULLIMORE: Oh, please.

22 (Laughter.)

23 DR. EYDELMAN: So is the panel recommending the
24 indication for use statement to be changed, or the panel's
25 recommendation somehow will get reflected somewhere else in

1 the labeling? I just want to be clear on --

2 DR. McCULLEY: You have a regulatory issue that
3 I can't tell you how to deal with in how the study was
4 done. But wherever that was stated, I think it needs to be
5 not asterisked and put at the bottom of the page but
6 imprint right after it that even though the study was
7 initially done with this parameter, current practice is to
8 have refractive stability defined as no greater than a half
9 diopter change over the preceding year.

10 Is there agreement with that? Disagreement?
11 Dr. Yaross?

12 DR. YAROSS: Marcia Yaross. I would just
13 question is it the sponsor's job to say what is current
14 medical practice if it deviates from how the study was
15 done. I think that can cause some confusion.

16 DR. McCULLEY: Well, but I think it can also be
17 harmful for the FDA to put its mark of approval on
18 something that we currently don't consider adequate
19 stability.

20 Dr. Jurkus?

21 DR. JURKUS: Jan Jurkus. I would have a real
22 problem with the 1 diopter. Say, for example, if we had
23 someone that might have been a quarter diopter myope and
24 now was a 1.25 diopter myope six months later because of
25 functional changes. They had the surgery done and they

1 might continue to become more myopic. So I think the half
2 diopter is certainly more realistic.

3 DR. McCULLEY: So again, our recommendation is
4 not to disadvantage sponsor, and I think both viewpoints or
5 circumstances can be dealt with effectively.

6 Are there other comments from panel? Karen?

7 DR. BANDEEN-ROCHE: Karen Bandeen-Roche. This
8 is purely a clarification. I'm just concerned that my
9 comments in the record are very unclear as they stand. So
10 just one sentence. The parameter for which I was computing
11 a confidence bound is the proportion of subjects who can be
12 expected to achieve some reversibility criterion, such as
13 the proportion of subjects who can be expected to return to
14 within 0.5 manifest refractive cylinder post-removal.

15 DR. McCULLEY: Okay. Other comments?

16 DR. MIDDLETON: I have two comments with
17 respect to the patient booklets. Is it appropriate to make
18 those now?

19 DR. McCULLEY: I think it's appropriate now.

20 DR. MIDDLETON: In the patient booklet under
21 "Warnings," they list for near-sightedness that you should
22 discuss with your doctor if your near-sightedness is
23 changing, diabetic, et cetera. I think that they should
24 indicate why, and they do so in the physician's booklet.
25 They simply need to say the safety and effectiveness of the

1 ring has not been established in patients with these
2 conditions. Rather than just listing those, don't leave a
3 question in the consumer's mind. You should say why.

4 Then the second point, are you a good candidate
5 for the KeraVision ring, I think that in the last bullet
6 you're asking the patient to be informed as to the risks
7 and the benefits as compared to all other available
8 treatments. You should indicate what those treatments are,
9 such as spectacles, contact lenses. Again, you do so in
10 the physician's booklet. I don't think it's too much to
11 ask or too much information to put there for the
12 individual. In the very first bullet, the requirement is
13 to be at least 21 years of age. You don't say to be the
14 appropriate age and leave it to the individual to find out
15 what that age is. You say 21. So by the same token, other
16 available treatments, just list them. It's not a lot:
17 spectacles, contact lenses, other surgical procedures. The
18 same way you do in the physician's booklet I think you
19 should do here.

20 Those are just two points that I would address.

21 DR. McCULLEY: Thank you.

22 Are there other comments from panel? Dr.
23 Higginbotham.

24 DR. HIGGINBOTHAM: I'd like to follow Dr.
25 Middleton's lead and suggest that in the patient booklet --

1 this is to humor me as well, I suppose -- on page 6, under
2 "Warnings," I think we should explicitly state that you
3 should discuss with your doctor if you have elevated
4 pressure or a history of elevated pressure or a suspicion
5 of glaucoma. I think as I read the booklet my
6 understanding was that glaucoma was to be rolled into this
7 overriding term called "disease." It may not come to the
8 patient's mind if they don't actually see that word. So
9 just a suggestion.

10 DR. McCULLEY: Other comments?

11 (No response.)

12 DR. McCULLEY: At this point I want to open the
13 floor to another public hearing session. Anyone in the
14 audience who would like to approach the podium other than
15 the sponsor to make comments about the deliberations
16 relative to this PMA specifically are invited to -- and it
17 also can't be panel, Dr. Bullimore.

18 (Laughter.)

19 DR. McCULLEY: You'll lose more hair.

20 (Laughter.)

21 DR. McCULLEY: Please approach the podium at
22 this time.

23 (No response.)

24 DR. McCULLEY: Seeing none, the open public
25 hearing session is closed.

1 We now have the opportunity for five minutes of
2 closing comments from sponsor.

3 DR. LEMP: Thank you. Dr. Lemp representing
4 the sponsor.

5 First of all, the sponsor would like to thank
6 the panel and the agency for a very careful evaluation of
7 this device. We realize that we provided you with a lot of
8 data, and it's difficult to analyze a lot of data, and we
9 appreciate all the time and effort that you put in on this.

10 What I'd like to do is just make a couple of
11 comments to a couple of the areas concerning comments that
12 the panel had to try to clear things up.

13 The first comment relates to the vacuum time
14 that Dr. Higginbotham brought up. It's important to
15 recognize that the mean vacuum time on these procedures was
16 1.4 minutes, with a standard deviation of 0.43. So most of
17 these procedures were about a minute and a half, just a
18 little over a minute and a half.

19 A comment on your other question, this can be
20 removed and reapplied, so it doesn't have to be continuous.
21 So that's one point.

22 Another point that I would make has to do with
23 the reversibility claim on this. I would just reiterate
24 one thing that went by on a slide perhaps quickly and was
25 not easy to assimilate all of these things.

1 If you use the most rigid criteria that we used
2 for reversibility, which was return to within a half a
3 diopter of the preop manifest spherical equivalent, a half
4 a diopter of the preop cylinder, and one line plus or minus
5 of visual acuity, 24 out of 27, or 89 percent met that
6 criteria. That relates to some of the data points that Dr.
7 Yaross brought up. If you use three-quarters of a diopter
8 for the spherical equivalent and the cylinder, and one
9 line, it's 100 percent of that N of 27.

10 The other points that I would bring up relate
11 to endothelial cell density, and I think it's important to
12 just clarify one or two points about that.

13 If you look at the eyes that had greater than
14 10 percent cell loss and you look at the six-month data,
15 which we had for both the treated eye and the fellow eye,
16 and that's as far out as we could get from the fellow eye
17 because the fellow eye was eligible for implantation beyond
18 that point, the number of patients who had greater than a
19 10 percent cell loss was 12 percent for the treated eyes
20 and 9 percent for the fellow eyes at the six-month point.
21 If you go out to 12 months, it was 14 percent for the
22 treated eye -- that's one more eye than we found at the
23 six-month data -- and four of those six eyes were the same.

24 Now, if you look at the data which we also
25 presented, which was the paper by Trocme et al., which

1 looked at the longitudinal study of peripheral cell data in
2 patients with PRK, which was the only study we could find
3 that had any data on peripheral cell endothelial cell
4 densities in a longitudinal study, they found a loss in the
5 periphery of -6.9 percent, or 221 cells. If you look at
6 our study, at the 10 o'clock position, which was the area
7 of concern, we had at the same time a loss of -4.8 percent
8 at the 10 o'clock position, which was 138 cells.

9 Now, in terms of the variability of this test
10 that's already been alluded to, Dr. Edelhauser was a
11 principal author on a paper that's just been published
12 looking at endothelial cell density after LASIK. They were
13 looking at central endothelial cell density, and even in
14 the central endothelial cell density studies, there's a 10
15 percent variation, which relates to what Dr. McCulley was
16 relating to earlier. That's 300 cells variation. So even
17 in the central area, that seems to be the limit of the
18 technology that we can do. That becomes considerably more
19 problematic when you get to the periphery of the cornea,
20 and particularly when you get to the periphery in a cornea
21 that has its curvature changed in the area that you're
22 looking at.

23 If you look at some of the specular
24 photomicrographs, you can actually see the area where you
25 can't see the cells as well, presumably when you start to

1 get that change, and it really adds to the question of the
2 reproducibility. As Dr. Edelhauser has pointed out, you
3 can't be sure you're coming back to the same area.

4 Now, in terms of your concern about
5 longitudinal studies, I would say that the sponsor also
6 agrees with some of your concerns. The real question is,
7 is there any real cell loss here? We don't think there is.
8 We think this is noise in the system, but it doesn't hurt
9 to do some longitudinal studies. I'm not sure you're going
10 to find anything, however, by doing longitudinal studies of
11 the peripheral cornea, because we don't know that we're
12 measuring anything there, that we have any reproducibility
13 in this.

14 There have been a number of studies dating back
15 quite a few years in terms of what happens when you do have
16 real change in the periphery of the cornea secondary to
17 injury, such as after surgery or when you create that
18 injury, and you see it reflected in central endothelial
19 cell counts because you get an equilibration of cells
20 coming in, and we would think that it would make more sense
21 if you're going to follow this that you follow the central
22 cell counts because they would be much more reliable, more
23 likely to have meaningful data than the peripheral cell
24 counts because we have nothing to go on in terms of
25 longitudinal studies for normals, and the methodology in

1 Cell appears to be inherently flawed.

2 So those are some of the comments we would make
3 about those studies.

4 DR. McCULLEY: Mike, you have five minutes
5 total. You're honing in on that.

6 DR. LEMP: Okay. Well, actually, that's about
7 what I really wanted to say, except that as a sponsor I
8 think all of us feel that we have provided reasonable
9 assurance for the safety and effectiveness in all the ring
10 sizes for the study. Thank you for your careful
11 consideration.

12 DR. McCULLEY: Thank you, sponsor.

13 We now have five minutes for the FDA to make
14 closing comments before we move towards voting, a motion.

15 (No response.)

16 DR. McCULLEY: Okay. You need to read the
17 voting options.

18 MS. THORNTON: The panel voting options. The
19 panel has a copy of these in their folder if they want to
20 follow along.

21 The Medical Device Amendments to the Federal
22 Food, Drug and Cosmetic Act require that the Food and Drug
23 Administration obtain a recommendation from an outside
24 expert advisory panel on designated medical device
25 premarket approval applications that are filed with the

1 agency. The PMA must stand on its own merits, and the
2 recommendations must be supported by safety and
3 effectiveness data in the application, or by applicable
4 publicly available information.

5 Safety is defined in the Act as reasonable
6 assurance based on valid scientific evidence that the
7 probable benefits to health under conditions of use
8 outweigh any probable risks. Effectiveness is defined as
9 reasonable assurance that in a significant portion of the
10 population, the use of the device for its intended uses,
11 and the conditions of use, when labeled, will provide
12 clinically significant results.

13 Your recommendation options for the vote are as
14 follows:

15 Approval. If you recommend approval, there are
16 no conditions attached.

17 You can recommend approvable with conditions.
18 You may recommend that the PMA be found approvable subject
19 to specified conditions, such as resolution of clearly
20 identified deficiencies which have been cited by you or FDA
21 staff. Prior to voting, all the conditions are discussed
22 by the panel and listed by the panel chair. You may
23 specify what type of follow-up to the applicant's response
24 to the conditions of your approval recommendation you want;
25 for example, FDA or panel. Panel follow-up is usually done

1 through homework assignments to the primary reviewers of
2 the application or to other specified members of the panel.
3 A formal discussion of the application at a future panel
4 meeting is not usually held.

5 If you recommend post-approval requirements --
6 for example, post-approval follow-up studies or postmarket
7 surveillance to be imposed as a condition of approval --
8 then your recommendation should address the following
9 points: the purpose of the requirement, the number of
10 subjects to be evaluated, and the reports required to be
11 submitted.

12 The third option is not approvable. Of the
13 five reasons the Act specifies for denial of approval, the
14 following three reasons are applicable to panel
15 deliberations: the data do not provide reasonable
16 assurance that the device is safe under the conditions of
17 use prescribed, recommended, or suggested in the proposed
18 labeling; reasonable assurance has not been given that the
19 device is effective under the conditions of use prescribed,
20 recommended, or suggested in the labeling; and based on a
21 fair evaluation of all material facts in your discussions,
22 you believe the proposed labeling to be false and
23 misleading.

24 If you recommend that the application is not
25 approvable for any of these stated reasons, then we ask you

1 to identify the measures you think are necessary for the
2 application to be placed in an approvable form. If FDA
3 agrees with the panel's not approvable recommendation, we
4 will send a "Not Approvable" letter. Please note: This is
5 not a final agency action on the PMA. The applicant has
6 the opportunity to amend the PMA to supply the requested
7 information. The panel at a future meeting may review the
8 amended application.

9 Please note that following the voting, the
10 chair will ask each panel member to present a brief
11 statement outlining the reasons for their vote.
12 Traditionally, the consumer representative and the industry
13 representative do not vote, and Dr. McCulley as chairperson
14 votes only in the case of a tie.

15 Thank you, Dr. McCulley.

16 DR. McCULLEY: Dr. Sugar, would you like to
17 make a motion?

18 DR. SUGAR: I'd like to move for approval with
19 conditions, the conditions being the changes in the
20 labeling that were discussed as we went through the eight
21 questions. In addition, after approval, I would like to
22 suggest that data continue to be accrued on endothelial
23 cell counts as we discussed earlier also.

24 DR. McCULLEY: Is there a second?

25 DR. HIGGINBOTHAM: Second.

1 DR. McCULLEY: Is there further discussion?

2 DR. VAN METER: To specifically elaborate --

3 DR. McCULLEY: Dr. Van Meter.

4 DR. VAN METER: Dr. Van Meter. The endothelial
5 cell counts you would do for the 0.35 millimeter implants
6 only, and continue to follow those patients?

7 DR. SUGAR: No. What we discussed was all of
8 the 0.35's and a selected group of the two others as a
9 control.

10 DR. VAN METER: Okay, but that select group of
11 others are the ones who have already had --

12 DR. SUGAR: All of this is patients who have
13 already had surgery but will be followed subsequent to now.

14 DR. VAN METER: Thank you.

15 DR. McCULLEY: Dr. Grimmett?

16 DR. GRIMMETT: Just as a point of
17 clarification, as long as we're thinking about a study on
18 peripheral cell counts, as long as we're there, central
19 cell counts too? Is that correct?

20 DR. SUGAR: The statement was endothelial cell
21 counts presumably derived the same way they were previously
22 in the study, which was three points -- 6, 10, and central.

23 DR. GRIMMETT: Thank you.

24 DR. McCULLEY: Other discussion? Dr.
25 Bullimore.

1 DR. BULLIMORE: In addition to the Grimmett
2 list of eight, I'm assuming that there are --

3 DR. McCULLEY: That was the eight questions,
4 not the Grimmett eight. We're talking about the panel's
5 response. I'm not going to try to reiterate all of those
6 because I will leave something out.

7 DR. SUGAR: They also include the comments from
8 Dr. Middleton on changes in the patient brochure.

9 DR. McCULLEY: Yes, the additional one. But we
10 reached consensus on each point as we went through the
11 eight questions. The various consensuses -- I guess that's
12 a word -- that we agreed upon are part of the conditions.

13 Seeing no further requests for discussion,
14 there's a call for the question.

15 All in favor of the motion, please signify by
16 raising your right hand.

17 (Show of hands.)

18 DR. McCULLEY: I count 11 ayes.

19 All opposed?

20 (No response.)

21 DR. McCULLEY: Eleven is the total number of
22 voting members. Therefore, the motion is unanimously
23 recommended for approval, with conditions as stated.

24 Now we'll go around and each voting member will
25 be asked to state the reasons for voting the way in which

1 you voted. We'll start with Karen, Dr. Bandeen-Roche.

2 DR. BANDEEN-ROCHE: Karen Bandeen-Roche. I
3 voted approvable with conditions because I believe that the
4 sponsor has demonstrated the device to be safe and
5 effective, subject to the labeling and other conditions we
6 discussed.

7 DR. McCULLEY: Dr. Grimmett?

8 DR. GRIMMETT: Dr. Grimmett. Based on my
9 review, I believe the sponsor has shown reasonable safety
10 and effectiveness with the labeling specifications that we
11 discussed and the additional data that's needed on the
12 endothelium.

13 DR. McCULLEY: Dr. Wang?

14 DR. WANG: Ming Wang. I think this is a
15 reasonable and relatively safe option for refractive
16 correction in the range that we discussed, and I recommend
17 approval with conditions, with endothelial counts and other
18 considerations included.

19 DR. VAN METER: Woody Van Meter. I voted
20 approvable with conditions for the same reasons.

21 DR. McCULLEY: Dr. Matoba?

22 DR. MATOBA: Alice Matoba. I voted for
23 approval for the same reasons.

24 DR. McCULLEY: Dr. Macsai?

25 DR. MACSAI: Dr. Macsai. I voted for

1 approvable with conditions basically because the perfect
2 refractive surgical procedure does not exist. However, the
3 KeraVision Intacts bring us closer to it. The Intacts
4 appear to be reasonably effective and safe. There are
5 concerns of continued endothelial cell loss. Some patients
6 noted a decrease in their visual quality, with halos,
7 double vision, and fluctuations in distance vision
8 regardless of pupil size. However, as opposed to other
9 refractive surgical procedures, the KeraVision Intacts are
10 removable.

11 (Laughter.)

12 DR. McCULLEY: Dr. Pulido?

13 DR. PULIDO: I voted approvable with conditions
14 for similar reasons mentioned previously. Also, I'd like
15 to state that I do thank the sponsors for giving us the
16 data even though sometimes it was difficult to flush out.
17 But they did give us good data to evaluate, which made our
18 evaluations easier.

19 DR. HIGGINBOTHAM: Dr. Higginbotham. Based on
20 the data presented both by the sponsor and FDA, and I must
21 say that the data was well presented by both, and also the
22 prospect of the long-term study on the cornea, I voted for
23 approvable with conditions. Thank you.

24 DR. McCULLEY: Dr. Sugar?

25 DR. SUGAR: I voted for approvable with

1 conditions, and I pretty much stated my reasons. This
2 appears to be a safe, effective, and reversible procedure.

3 (Laughter.)

4 DR. BULLIMORE: This is Dr. Bullimore. I voted
5 for approvable with conditions for the same reasons.

6 DR. McCULLEY: Dr. Jurkus.

7 DR. JURKUS: Dr. Jurkus. I voted approvable
8 with conditions for the same reasons.

9 DR. McCULLEY: Does FDA have any concluding
10 remarks?

11 MS. THORNTON: I would just like to take the
12 time to thank all of you for a very, very thorough
13 preparation for this meeting. I know it was very
14 difficult. There were a lot of volumes. It was quite
15 heavy. It's in the interest of physical fitness that we
16 send you these things.

17 (Laughter.)

18 MS. THORNTON: And to thank the sponsor for
19 their time, and I appreciate your cooperation with this.

20 DR. McCULLEY: I'd like to thank the sponsor,
21 the FDA, and the panel members for I think a well done
22 session that did some good. Thank you.

23 We stand adjourned.

24 (Whereupon, at 3:35 p.m., the meeting was
25 adjourned.)