

sgg

AT

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
FOOD AND DRUG ADMINISTRATION  
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

GENERAL AND PLASTIC SURGERY DEVICES PANEL  
MEDICAL DEVICES ADVISORY COMMITTEE

OPEN SESSION

Wednesday, June 16, 1999

10:00 a.m.

Gaithersburg Hilton  
Salons D and E  
620 Perry Parkway  
Gaithersburg, Maryland

MILLER REPORTING COMPANY, INC.  
507 C Street, N.E.  
Washington, D.C. 20002  
(202) 546-6666

## PARTICIPANTS

Thomas V. Whalen, M.D., Acting Chairperson  
David Krause, Ph.D., Executive Secretary

## VOTING MEMBERS

Phyllis Chang, M.D.  
Benjamin O. Anderson, M.D.  
Susan Galandiuk, M.D.  
David L. DeMets, Ph.D.

## DEPUTIZED VOTING MEMBERS

Thomas B. Ferguson, M.D.  
Blake Hannaford, Ph.D.  
Michael D. Crittenden, M.D.  
Mark A. Talamini, M.D.  
Nancy N. Dubler, LLB.  
Cedric F. Walker, Ph.D., PE

## NONVOTING MEMBERS

Maxine F. Brinkman, RN, Consumer Representative  
James W. Burns, Ph.D., Industry Representative

## C O N T E N T S

Conflict of Interest and Opening Remarks, David Krause, Ph.D.	4
Panel Introduction	7
Update Since Last Meeting, Stephen Rhodes	10
"Postmarket Evaluation at FDA's Center for Devices and Radiological Health, Larry G. Kessler	12
Presentation of Award, Mr. Jim Dillard	30
Applicant Presentation, Intuitive Surgical Incorporated:	
Introduction, Michael Daniel	32
Technology, Fred Moll, M.D.	34
Statistical Analysis, Dan Bloch, Ph.D.	40
Clinical Study and Results, Barry Gardiner, M.D.	49
FDA Presentations	
Preclinical and Technical Aspects, Dwight Yen	98
Clinical Aspects, Roxi Horbowyj, M.D.	101
Statistical Aspects, Harry Bushar, Ph.D.	110
FDA Questions, Dwight Yen	119
Panel Deliberations	
Preclinical Overview, Blake Hannaford, M.D.	121
Clinical Overview, Mark A. Talamini, M.D.	131
Statistical Overview, David L. DeMets, M.D.	138
Discussion of FDA Questions by the Panel	159
Concluding Panel Deliberations and Vote	184

P R O C E E D I N G S

**Conflict of Interest and Opening Remarks**

DR. KRAUSE: I would like to get the meeting started, please. Good morning, everyone. We are ready to begin this meeting of the General and Plastic Surgery Devices Panel. My name is David Krause. I am the executive secretary of the panel. I am also a reviewer in the Plastic and Reconstructive Surgery Branch in the DGRD.

I would like to remind everyone that you are requested to sign in on the attendance sheets, which are outside the door, on the tables. You may also pick up an agenda, a panel meeting roster and information about today's meeting there. The information includes how to find out about future meeting dates through the advisory panel phone line and how to obtain meeting minutes or transcripts.

Before I turn the meeting over to Dr. Whalen I am required to read a number of statements into the record. I have three statements to read. One is the deputization of temporary voting members. The second is deputization of acting chair, and the third is the conflict of interest statement.

This is the temporary voting status: Pursuant to the authority granted under the Medical Devices Advisory Committee Charter, dated October 27th, 1990, as amended on April 20th, 1995 and October 10th, 1997, I appoint the



1 following people as voting members of the General and  
2 Plastic Surgery Devices Panel for this meeting on June 16,  
3 1999: Blake N. Hannaford, Nancy A. Dubler, Mark A.  
4 Talamini, Cedric F. Walker, Thomas B. Ferguson and Michael  
5 D. Crittenden.

6 For the record, these people are special  
7 government employees and are consultants to the Center for  
8 Devices and Radiological Health under the Medical Devices  
9 Advisory Committee. They have undergone customary conflict  
10 of interest review. They have reviewed the material to be  
11 considered to this meeting. And, it is signed by Dr.  
12 Elizabeth Jacobson, Acting Director, Center for Devices.

13 I would now like to read the appointment of the  
14 temporary panel chair person. I appoint Thomas V. Whalen,  
15 M.D. to act as temporary chairperson for the duration of the  
16 General and Plastic Surgery Devices Panel meeting on June  
17 16th, 1999. For the record, Dr. Whalen is a special  
18 government employee and is a voting member of the General  
19 and Plastic Surgery Devices Panel. Dr. Whalen has undergone  
20 customary conflict of interest review. He has reviewed the  
21 issues to be considered at this meeting. It is signed by  
22 Dr. Feigel who is the present Director of the Center for  
23 Devices.

24 Finally, I would like to read the conflict of  
25 interest statement. The following announcement addresses

1 conflict of interest issues associated with this meeting,  
2 and is made part of the record to preclude even the  
3 appearance of an impropriety. To determine if any conflict  
4 existed, the agency reviewed the submitted agenda and all  
5 financial interests reported by the panel participants. The  
6 conflict of interest statutes prohibit special government  
7 employees from participating in matters that could affect  
8 their or their employers' financial interests. However, the  
9 agency has determined that participation of certain members  
10 and consultants, the need for whose services outweighs the  
11 potential conflict of interest involved, is in the best  
12 interest of the government.

13           The agency took into consideration matters  
14 regarding Drs. Cedric Walker, Blake Hannaford, Susan  
15 Galandiuk and Benjamin Anderson. These individuals reported  
16 financial interests in firms at issue, but in matters not  
17 related to the topic to be discussed by the panel. The  
18 agency has determined, therefore, that they may participate  
19 fully in today's deliberations.

20           In the event that the discussions involve any  
21 other products or firms not already on the agenda for which  
22 an FDA participant has a financial interest, the participant  
23 should excuse him or herself from such involvement, and  
24 their exclusion will be noted for the record.

25           With respect to all other participants, we ask in

1 the interest of fairness that all persons making statements  
2 or presentations disclose any current or previous financial  
3 involvement with any firm whose products they may wish to  
4 comment upon.

5 At this time, I would like to turn the meeting  
6 over to Dr. Whalen.

7 **Panel Introduction**

8 DR. WHALEN: Thank you, Dr. Krause. My name is  
9 Thomas V. Whalen, and I am head of the Division of  
10 Pediatrics Surgery at Robert Wood Johnson Medical School at  
11 Camden, where I also hold the role of Associate Professor of  
12 Surgery and Pediatrics.

13 Today, our panel will be making recommendations to  
14 the Food and Drug Administration on a premarket approval  
15 application. Our next item of business is to introduce our  
16 panel members, who are giving of their time to help the FDA  
17 in these matters, and the FDA staff who are here at this  
18 table. I would like to ask that each person introduce him  
19 or herself, stating his or her specialty, position title,  
20 institution and his or her status on the panel as voting  
21 member, industry or consumer representative or deputized  
22 voting member. I would like to ask Dr. Burns to begin and  
23 we will go around the table.

24 DR. BURNS: I am Jim Burns. I am Vice President  
25 of Biomaterial and Surgical Products Research at Genzyme

1 Corporation, and I have a Ph.D. in bioengineering, and my  
2 expertise is in biomaterial science.

3 MS. BRINKMAN: I am Maxine Brinkman, Director of  
4 Women's Services, from Mercy Medical Center in Mason City,  
5 Iowa. I have an RN and a masters in public health.

6 DR. WHALEN: Could you each mention your status on  
7 the panel?

8 DR. BURNS: I am the industry rep. for this panel.

9 MS. BRINKMAN: I am the consumer rep.

10 DR. DEMETS: My name is David DeMets. I am a  
11 biostatistician, University of Wisconsin, and I am a voting  
12 member of the panel.

13 DR. FERGUSON: I am Tom Ferguson, Cardiothoracic  
14 Surgery, Washington University School of Medicine. I am a  
15 deputized voting member.

16 DR. HANNAFORD: My name is Blake Hannaford. I  
17 have no middle name, which is what the "N" stands for I  
18 think --

19 [Laughter]

20 I am a Professor of Electrical Engineering at the  
21 University of Washington, in Seattle. I am also Adjunct  
22 Professor of Bioengineering and Adjunct Professor of  
23 Surgery, and I have a Ph.D. in electrical engineering and  
24 computer science. I am a deputized voting member of the  
25 panel.

1 DR. GALANDIUK: Susan Galandiuk. I am a  
2 colorectal surgeon and I am an Associate Professor of  
3 Surgery at the University of Louisville, and I am a voting  
4 member of the panel.

5 DR. CRITTENDEN: I am Mike Crittenden. I am a  
6 cardiothoracic surgeon at the West Roxbury VA, and  
7 affiliated with the Harvard Medical School. I am a  
8 deputized voting member.

9 DR. ANDERSON: I am Ben Anderson. I am an  
10 Associate Professor of Surgical Oncology at the University  
11 of Washington, in Seattle. I am the Medical Director of the  
12 Breast Care Program there, and a member of the Fred  
13 Hutchinson Cancer Research Center.

14 DR. CHANG: I am Phyllis Chang. I am an Associate  
15 Professor in the Department of Surgery, Section of Plastic  
16 Surgery, with a joint appointment in the Department of  
17 Orthopedics Surgery for Hand and Microsurgery, in the  
18 University of Iowa, College of Medicine. I am a voting  
19 member.

20 DR. TALAMINI: Mark Talamini. I am an Associate  
21 Professor at Johns Hopkins University School of Medicine. I  
22 am a general surgeon, gastrointestinal surgeon primarily,  
23 and Director of Minimally Invasive Surgery at Hopkins, and I  
24 am a deputized voting member.

25 MS. DUBLER: I am Nancy Dubler. I am Director of

1 the Division of Bioethics at Montefiore Medical Center, and  
2 a Professor of Bioethics at the Albert Einstein College of  
3 Medicine, and I am a deputized voting member.

4 DR. WALKER: Cedric Walker, Professor of  
5 Biomedical Engineering at Tulane University, in New Orleans.  
6 I am a deputized voting member.

7 DR. WITTEN: Celia Witten, Division Director of  
8 DGRD in ODE in CDRH and FDA. I am a representative of FDA,  
9 not a panel member.

10 DR. WHALEN: Thank you. I would like to note for  
11 the record that the voting members present constitute a  
12 quorum, as required by 21 CFR, Part 14.

13 Next, we are going to hear from Mr. Stephen Rhodes  
14 who will give the panel an update since our last meeting in  
15 November. Mr. Rhodes?

16 **Update Since the Last Meeting**

17 MR. RHODES: Thank you. I am Stephen Rhodes. I  
18 am the Branch Chief of the Plastic and Reconstructive  
19 Surgery Devices Branch. In November the panel made  
20 recommendations on the classification of five wound  
21 dressings, gauze, hydrophilic wound dressings, hydrogel  
22 wound dressings, occlusive wound dressings and porcine wound  
23 dressings.

24 The FDA is working on a final rule classifying  
25 four of those wound dressings, the gauze, the hydrogel wound

1 dressings, the occlusive wound dressings and the hydrophilic  
2 wound dressings, as Class I exempt devices. These are  
3 products that do not have any drugs, any biologics or any  
4 animal source materials in them. We are still working on  
5 classifying the porcine wound dressings. We appreciate your  
6 assistance in helping us to classify these products, and we  
7 are pleased that we are reaching closure on classifying all  
8 these important classes of products.

9 Other activities in the branch -- we have  
10 completed a guidance document on surgical meshes. We are  
11 working on updating several guidance documents, a guidance  
12 document for breast implants which was developed in 1995,  
13 and we are also updating the guidance documents for  
14 interactive wound dressings and noninteractive wound  
15 dressings. We are also developing a new guidance document  
16 for sutures. Lastly, we are working on a final rule,  
17 calling for PMAs for safety and efficacy data for saline-  
18 filled breast implants for mammary prostheses. We are  
19 expecting that that final rule will be published this year,  
20 and we are looking forward to the panel's participation in  
21 the discussion of those PMAs in the future.

22 That concludes my update. Thank you.

23 DR. WHALEN: Thank you, Mr. Rhodes. We will next  
24 hear from Mr. Larry Kessler regarding postmarket evaluation  
25 at FDA's Center for Devices and Radiological Health. Mr.

1 Kessler?

2 "Postmarket Evaluation at FDA's Center  
3 for Devices and Radiological Health"

4 MR. KESSLER: Good morning. My name is Larry  
5 Kessler. I am the Director of the Division of Surveillance  
6 and Biometrics at FDA. I am not a voting member..

7 [Slide]

8 Around two years ago, maybe a little less, Dr.  
9 Susan Alpert, who is the Director of the Office of Device  
10 Evaluation, asked me to give a brief presentation to a group  
11 of the panel chairs who collected together to review issues  
12 and common to panels. And, I talked a little bit about  
13 postmarket evaluations at FDA and a number of the panel  
14 chairs asked me to come, give presentations to the members  
15 of the panels, not just the chairs.

16 So, Dr. Tom Gross, who is the head of my Division  
17 of Postmarket Surveillance, and myself have been going  
18 around this year to give presentations to all the panels,  
19 trying to give you a little flavor for how we work on the  
20 postmarket side. So, in 15 minutes I am going to describe  
21 everything we think you need to know, in a brief primer, on  
22 postmarket evaluation.

23 As I am going to show you in a minute, there are a  
24 variety of different mechanisms which use similar names and  
25 may sometimes get confusing to you, as panel members, and



1 are even sometimes confusing to FDA staff who work with them  
2 on a week to week basis.

3 [Slide]

4 I am going to describe a few of the methods of  
5 device postmarket evaluation at the Center. I am going to  
6 present the challenges we face in accomplishing postmarket  
7 evaluation roles, and just try and describe the pivotal role  
8 that you, as the advisory panel, can plan in helping us  
9 accomplish the goals for making sure that devices are safety  
10 and efficacy, not only as they reach market but as they  
11 spend time on the market.

12 As you well are aware, most of the time a device  
13 is in use is postmarket, not premarket, and that is the bulk  
14 of our concern. I believe you have handouts on this in case  
15 you want to take notes, and I hope you have questions.

16 [Slide]

17 This is a brief schematic of how we view the  
18 world. This is the world of the Center for Devices viewed  
19 from the context of the Office of Surveillance and  
20 Biometrics.

21 For the most part, from design through clinical  
22 testing the bulk of this work, the best majority of this is  
23 done by manufacturers, clinicians, their patients and not by  
24 the FDA. This is where the bulk of the work is done in the  
25 premarket phase. FDA gets increasingly involved as any

1 device or modification of a device goes through lab and  
2 bench testing and clinical testing. Once it reaches FDA  
3 review, we have the clinical community that we hope we are  
4 interacting with a lot. The advisory panels are one  
5 mechanism for that. Once a device is cleared for market, we  
6 have a variety of tools that we can use in the postmarket  
7 period to evaluate continued safety and effectiveness issues  
8 during the postmarket period.

9 I am going to talk quite a bit about the medical  
10 device reporting program, and then about postmarket  
11 surveillance authorities and postapproval studies. You are  
12 probably most familiar with postapproval studies but I will  
13 try and talk about all three and weave them together a  
14 little bit. I am not going to talk very much about our  
15 epidemiology program or field inspection program, which are  
16 also important components of postmarket evaluation but I  
17 don't have that much time today and I can always come back  
18 if you are interested.

19 [Slide]

20 What do we care about in the postmarket period?  
21 Clearly, long-term safety is one issue that matters to us,  
22 matters to you, matters to everybody involved. We are also  
23 very interested in the postmarket period about looking at  
24 performance of a device in community practice. You will see  
25 premarket applications that frequently come from well-

1 controlled studies, clinical settings where there may be  
2 device experts; there is a great deal of attention to  
3 protocol. Once devices get out into the community they  
4 don't achieve the same level of performance often, and  
5 sometimes you see problems you would not have seen.

6 In a minute I will tell you a little bit about the  
7 MDR program, but if you think you can imagine everything  
8 that can happen to devices you clear, you can't. You must  
9 see these reports. It is amazing what people would think  
10 about doing when devices have been approved and cleared for  
11 market.

12 Effects of change in user setting -- as you well  
13 know, an enormous amount of product has left the hospital  
14 walls and moved to the bedside, and it has been particularly  
15 true over the last 16 years since the imposition of DRGs by  
16 HCFA and the insurance world. People are trying to move  
17 patients out of hospitals faster and faster.

18 Some of this is very positive from a therapeutic  
19 standpoint, but it means increasingly sophisticated  
20 technologies making it to the bedside. We, on the  
21 postmarket end, tend to see an increase in the different  
22 kinds of problems that we see. This includes serious  
23 injuries and deaths from a variety of products that you  
24 would not tend to see in a hospital setting but you will see  
25 in a home at bedside because you have people who are less

1 expert and who don't have the resources to deal with  
2 problems once they arise. So, whether in fact some piece of  
3 technology should move from the hospital to the bedside is  
4 something that we will often ask in the postmarket of a  
5 product.

6 [Slide]

7 One of our most quoted and used mechanisms for  
8 looking at product problems is the Medical Device Reporting  
9 Program. Manufacturers must by law report deaths and  
10 serious injuries, if a medical device may have caused or  
11 contributed to the event, as well as malfunctions. Now, in  
12 the European Union these are called "near incidents." A  
13 product that fails that did not happen to cause a death or  
14 serious injury but could if it occurred without such  
15 fortuitous circumstances is a mandated reportable event by a  
16 manufacturer.

17 Since the Safe Medical Devices Act of 1990, all  
18 user facility in this country -- all hospitals, all nursing  
19 homes, ambulances, surgical care centers, they are supposed  
20 to also report all deaths to the FDA, and all deaths and  
21 serious injuries to manufacturers of devices if they know  
22 the manufacturer of a product.

23 We get roughly 95 percent of our 80,000 reports  
24 from manufacturers and less than 5 percent come from  
25 hospitals, nursing homes, etc., and most of the events we

1 get from manufacturers are happening in hospitals, nursing  
2 homes, etc. So, we get the reports from manufacturers.  
3 Many of the manufacturers you see here stand up before you  
4 are critical components of our ability to monitor postmarket  
5 data. Hospitals and nursing homes -- some of them do well  
6 and many do not.

7 [Slide]

8 Beginning in about 1992, we were receiving at FDA  
9 over 100,000 medical device adverse events a year. We were  
10 really swamped with work. It was really a challenge because  
11 we had roughly 15 analysts to look over 100,000 reports and  
12 it was really a daunting task.

13 Information in these reports is supposed to  
14 include device specifics, event description, event date,  
15 patient characteristics, etc., but these reports often have  
16 limited information. Sometimes a manufacturer will come to  
17 us and they will say, "I'm going to send you a report but  
18 it's not complete." By law, they are supposed to complete  
19 all the information they can get. We will ask them why not.  
20 They will call the hospital who will tell them, "my lawyer  
21 tells me not to tell you what really happened." I am sorry  
22 for the legal guys out there; it's just something that  
23 happens in the walls of the hospital, and it sometimes is a  
24 challenge for us.

25 But often the reports are very useful and we use

1 them to prompt a lot of actions. What I have done is select  
2 a couple of products that are relevant to this panel where  
3 we have taken some actions based on the data from the MDR  
4 program. For example, we get reports of operating room  
5 fires with ESUs frequently. It is a common problem. You  
6 may have seen a piece that I think "20/20" did about a year  
7 ago, and it is something that we see repeatedly, and we  
8 think it is a problem that can be reduced by more Public  
9 Health action and we proposed to take some actions this year  
10 to make clinicians and the public aware of how we can reduce  
11 fires related to ESUs.

12 More recently, in 1997, a product was approved by  
13 FDA, I believe through the 510(k) program, for a collagen  
14 impregnated synthetic sling. We noticed some increased  
15 reports in 1998 and conducted several inspections with the  
16 manufacturer who eventually, upon more careful examination  
17 of their complaint data and our MDR reports, decided to  
18 recall the product in January of this year, and they will no  
19 longer market or manufacture the device.

20 DR. TALAMINI: Excuse me, sir --

21 MR. KESSLER: Yes?

22 DR. TALAMINI: -- as a surgeon, I have to know  
23 what an ESU is.

24 MR. KESSLER: Electrosurgical cautery unit.

25 DR. TALAMINI: Sorry.

1 MR. KESSLER: We are talking about things that are  
2 used in a very rich oxygen environment around drapes.

3 Based on the MDR system, we will also conduct  
4 directed inspection of facilities, and we will also conduct  
5 postmarket studies. We did one on polyurethane foam-coated  
6 breast implants a few years ago -- successful.

7 [Slide]

8 I am going to spend a couple of minutes on these  
9 two authorities because these two will probably interest you  
10 more than anything else. The Safe Medical Devices Act  
11 passed Section 522 postmarket surveillance. We have had  
12 this authority for about nine years now. We also are  
13 allowed to call for postapproval studies under the PMA.  
14 Both of these are mandated by law. Postapproval refers to  
15 PMA products only, whereas, Section 522 covers, as modified  
16 by FDA in May, 1997, Class II or III products whose failure  
17 may present a public health problem. I have stated that  
18 quite generally; the statutory language is a little more  
19 specific.

20 We have two mechanisms at FDA for calling for or  
21 requiring manufacturers to conduct studies in the postmarket  
22 period. These are principally done for safety concerns.  
23 They don't have to be done only for safety concerns; we can  
24 ask effectiveness questions but the principal spirit of the  
25 postapproval and postmarket and studies is for safety issues

1 but there are effectiveness and performance questions which  
2 are sometimes asked as well, and sometimes they are not  
3 separable for a certain kinds of device; its safety is its  
4 effectiveness issue. We see both authorities as  
5 complementary to postmarket.

6 [Slide]

7 To call for a study, we need to figure out what  
8 are the criteria to ask for such a study because studies in  
9 the postmarket period are quite challenging, and I will  
10 explain that in a minute. First of all, we need to  
11 understand the critical public health question, and these  
12 postmarket study requests can either occur from "for cause"  
13 situations or new expanded conditions of use or an evolution  
14 of technology.

15 We have to consider before we actually require a  
16 postapproval or postmarket study other mechanisms we might  
17 use to collect the same relevant data. This could include  
18 the Medical Device Reporting Program, secondary database  
19 analysis of, say, the HCFA data that we have in house, etc.

20 I want to talk for a bit about practicality and  
21 feasibility on my next slide and, most importantly, how will  
22 the data be used? Before we ask for a postmarket study, we  
23 need to figure out what we are going to do with it. Are we  
24 going to modify the label? Are we going to ask for a  
25 recall? Change the way we approve a device? What about



1 communication to the public, health professionals, etc?

2 [Slide]

3 There is a wide variety of study approaches that  
4 we should use in the postmarket period. With the  
5 institution of required or discretionary postmarket  
6 surveillance after SMDA, FDA tended to use a fairly heavy  
7 approach to postmarket studies and asked for things at the  
8 bottom end of this list. We were asking for randomized,  
9 case control or well-documented trials, and we had a lot of  
10 problems getting those accomplished.

11 With the recent FDAMA, the clear signal from  
12 Congress was to expand the kinds of things we might allow  
13 manufacturers to do to provide valuable postmarket data.  
14 So, in our recent Federal Register notices about our  
15 postmarket program we have expanded the kinds of designs  
16 that we might ask for not only from detailed randomized  
17 studies, but also possibly something as simple as detailed  
18 review of complaint history in files of manufacturers or  
19 literature, as well as non-clinical testing of devices.

20 [Slide]

21 But postmarket challenges have been very  
22 frustrating. They are frustrating for at least these four  
23 reasons: First, rapid evolution of technology make the  
24 studies that we are doing sometimes obsolete. By the time  
25 we get a protocol in house, the manufacturer goes and asks

1 clinicians for data, they begin to collect data, they are  
2 already on the second or third generation of that product  
3 and it is questionable of what utility the postmarket study  
4 might be.

5           It is not true of all medical device products but  
6 it is true of some. I still remember -- I spent, I guess,  
7 ten years at the National Cancer Institute working in cancer  
8 prevention and control, and when I arrived in 1984 one of  
9 the first things I was told there is, "the Pap smear is  
10 about to become obsolete, guys. It's been used for thirty  
11 years and we have much greater technology that's coming down  
12 the road, and it's going to be all gone by the time you hit  
13 1990." As you all know, it is still one of our major tools  
14 for cancer prevention and detection. It is true for  
15 devices. We hear things are going to be obsolete but you  
16 would be surprised at how long they hang on.

17           There is a lack of incentives for industry. The  
18 plain fact of this is when we ask industry for postmarket  
19 studies, it is generally not good news and their interest in  
20 helping us with those data is relatively minimal. It is  
21 unlikely that the end of the postmarket market study is  
22 going to be "what a wonderful product." It is more likely  
23 to be "we've discovered some problem and we want to tell you  
24 about them." Although this may be good marketing in the  
25 long-run, for industry it is not exciting news generally.

1 So, they are generally reluctant. Some of them have dragged  
2 their heels on some of the postmarket studies we have asked  
3 them for.

4 There is a general lack of interest in the  
5 clinical community. Those of you at Hopkins and the Fred  
6 Hutch and Washington U don't make a lot of publications off  
7 products that are already marketed, but it is nice to have  
8 sexy new technology for which you can help get information  
9 out in the literature. So, getting the clinical community  
10 interested in postmarket studies has been a frustrating  
11 challenge not only for us but also for industry. They find  
12 it hard.

13 Finally, sometimes on both ends, industry and the  
14 FDA end and on the advisory panel end, we have not clearly  
15 specified the public health question of importance for which  
16 to conduct a postmarket study.

17 [Slide]

18 That is the challenge I will leave you with. When  
19 considering postmarket studies, whether postapproval or 522  
20 -- and FDA can figure out the right mechanism to use; they  
21 each have different strengths and weaknesses. Postmarket  
22 studies under 522 can be done for 510(k) products.  
23 Condition of approval only apply to PMA products. We need  
24 to ensure that this is of primary importance, that it is not  
25 a secondary question. I reviewed, in 1997/8, almost a dozen

1 panel recommendations for postapproval studies or for PMA  
2 products. For half of them a postapproval study was asked  
3 for, and in most of those I could not detect clearly  
4 specified the question the panel was asking for us to ask a  
5 manufacturer, and that is a real challenge for us. So, if  
6 you are going to ask for a study, you really need to specify  
7 what you are trying to answer. That is critical.

8           Try and figure out what the clinical or regulatory  
9 relevance is for the question. What are you going to do  
10 with the data? Are you going to suggest we change the  
11 label? Are you going to suggest we put out a public health  
12 advisory to the clinical community, the lay community,  
13 hospitals? Who is going to get this information? Will we  
14 conduct a recall? If you discover, say, a long-term safety  
15 problem that was above some threshold, tell us the issue  
16 that you are interested in; clearly specify the question  
17 that will help us specify the question and work with the  
18 company to develop a design that makes sense. If you can do  
19 that, you can help us a lot.

20           [Slide]

21           Finally, I hope the future of the Medical Device  
22 Reporting program and postmarket surveillance -- I haven't  
23 talked a lot about the future of MRD, it is less relevant  
24 here -- we are trying to use a wider variety of design  
25 approaches to give the industry more flexibility and us more

1 flexibility to answer the right question correctly. We work  
2 in collaboration with the industry and the clinical  
3 community, and we are trying to obtain expanded access to  
4 different data sources, including registries, and we are  
5 beginning to hold a variety of workshops and conferences  
6 with industry and with the academic community to get better  
7 data in the postmarket period. Our data are somewhat  
8 limited, but with better specified questions, better  
9 specified studies and increased access to wider data sources  
10 which exist in the public domain we think we can answer the  
11 questions that we need to fill our mission in monitoring  
12 postmarket safety and effectiveness of devices.

13 Thank you for your time. I will take questions if  
14 you have any.

15 DR. CRITTENDEN: I have two questions. Could you  
16 talk about off-label use, and the second question is, there  
17 seem to be 100,000 adverse events per year and I just want  
18 to know how close that is to what the estimates are in terms  
19 of the actual number.

20 MR. KESSLER: I will answer the second question  
21 first because it is easy. We have lousy estimates -- let me  
22 go to the first question and then I will come back because  
23 the second one I have a longer answer for. Off-label use --  
24 we don't explicitly track off-label use in terms of studies  
25 unless you or we suggest what we should look for. However,

1 that is one of the great things about the Medical Device  
2 Reporting program.

3           We get routinely evidence of off-label use, and I  
4 will give you a great example -- it is not so much relevant  
5 for this panel, although I am sure we have dozens of them  
6 for you. It is an example I presented to the panel chairs.  
7 Soon after coronary stents were developed and started to be  
8 used, we would get reports of deaths or serious injuries  
9 related, and one of the egregious uses of stents is one  
10 clinician, wanting to keep a large part of an artery patent  
11 used 27 stents stringed in a row. They were labeled single  
12 stent only. Stringing a few of them together was something  
13 some clinicians were starting to do but one clinician got  
14 carried away, let's say, and the patient did pass away. The  
15 patient may have passed away anyway; it may not have been  
16 directly related to the use. But the MDR system is rife  
17 with those examples, and those are some of the ones on which  
18 we act. We will see a use that was not intended by the  
19 manufacturer which will provoke death or serious injury, and  
20 we will use public health advisories and safety alerts to  
21 the clinical community to tell them what we think is going  
22 on. So, that is the answer to the first question.

23           The second answer, the General Accounting Office  
24 made an estimate about a decade or so ago that we get about  
25 one percent of the real events that are happening. So, if

1 we get 100,000 events there are a million, and I think that  
2 estimate could be plus or minus, another factor of five or  
3 ten, quite frankly. It could be way up or way down. It is  
4 not very good. We have lousy denominator data, and even if  
5 we ask manufacturers to supply it, which we were going to do  
6 in 1996 but have walked away from it for a variety of  
7 reasons, it is not clear we would get good numerator data.  
8 You need both.

9 [Slide]

10 So, Congress, in its wisdom, supported the FDA by  
11 providing the legislative mandate for FDA to begin  
12 conducting a sentinel or device surveillance network  
13 program. We did a pilot study two years ago, which just  
14 ended last year, where we used 24 facilities and  
15 concentrated on those facilities in getting them to report  
16 faithfully, and we can get much better tracking of what is  
17 going on. The new sentinel system is supposed to begin in  
18 the next year or two, depending on the appropriations. We  
19 have been directed in statute to develop a sample user  
20 facility system so we can have a lab where we can look at  
21 good estimates of adverse events and denominator uses.  
22 Right now we don't have that. That is where we are going,  
23 and you will be hearing more and more about this system over  
24 the years provided Congress can find the wherewithal to  
25 appropriate the funds. But that is the way we need to go.

1 It is a really exciting program for us, and we think it is  
2 going to make a big difference in being able to find out how  
3 devices are really used.

4 DR. CRITTENDEN: Thank you. You have partially  
5 answered my question, but if you had all the data you wanted  
6 do you have the resources to really evaluate it properly at  
7 this stage?

8 MR. KESSLER: If we got 100 times the reports?  
9 No, no. If we got all the reports I wouldn't need to. I  
10 will just tell you about something we have done in the past  
11 two years, and give you a flavor for how we can work around  
12 this. We get 100,000 reports, but many of them are  
13 absolutely repetitive. It is the same thing over and over.  
14 The best example that relates to you is capsular contracture  
15 of breast implants. Once we get the first few hundred of  
16 those, the next 20,000 are not very interesting.

17 Now, I have to say this, and I am glad there is a  
18 consumer advocate here, we have a problem convincing women  
19 with breast implants and some of their associated legal  
20 counsel that more MDR reports of capsular contracture aren't  
21 a scientific way to help us evaluate a problem. So, what we  
22 have done is to move to a program that we call summary  
23 reporting. Once we see a problem over and over, we work  
24 with the industry, we work with the advisory panels, we work  
25 with the clinical community and try to figure out is there



1 something to be done. When there isn't we have the  
2 wherewithal to either have reports sent in, in a summary  
3 fashion so we get, say, one table four times a year of the  
4 same kind of complaint and, in that way, you can get a  
5 thousand complaints reviewed in a minute. So, that is one  
6 avenue we have taken. We can also give an exemption to  
7 companies -- don't report that. Keep it in your complaint  
8 file and in your quality system reporting and we can handle  
9 it there.

10 But the magic behind the sentinel system is that  
11 if we can get a good sample of between 5-10 percent of the  
12 hospitals in the country, the way the National Nosocomial  
13 Infection Surveillance System does in CDC, you can track the  
14 real problems fairly fast and you don't need a 100 percent  
15 sample. So, if we can get a good 5 or 10 percent of the  
16 hospitals reporting, and the rate in the sentinel pilot was  
17 20 times current report rates, we would get a lot more  
18 reports but from a much smaller sample. So, then we could  
19 handle it. I have 15 outstanding analysts. Almost all of  
20 them are nurses; one is not, but they are fabulous, and they  
21 pick up really unusual problems. They find needles in  
22 haystacks -- really unusual.

23 Other questions?

24 DR. WHALEN: Thank you, Mr. Kessler. Before  
25 proceeding to the next listed item on the agenda, for an

1 important unlisted piece of business I turn the floor back  
2 to Dr. Krause.

3 DR. KRAUSE: Sometimes we have the fun activity of  
4 handing out plaques and things. So, Jim?

5 **Presentation of Award**

6 MR. DILLARD: I am Jim Dillard. I am the Deputy  
7 Director in the Division of General and Restorative Devices.  
8 It is my distinct honor to get to stand up here and  
9 recognize one of you, although we would like to recognize  
10 all of you, but one of you for the distinct accomplishment  
11 of being on this panel now for your requisite period of  
12 time, and perhaps it may be this individual's last panel  
13 meeting, although this particular individual has been so  
14 good for us that usually when we have a particularly good  
15 member we keep them around and call them back from time to  
16 time, and have them perform their services on a semi-regular  
17 basis.

18 So without further ado, I would like to recognize  
19 Dr. James Burns, and your distinct recognition says, "in  
20 recognition of distinguished service from the General and  
21 Plastic Surgery Panel and Medical Devices Advisory  
22 Committee," and it is signed Elizabeth Jacobson, Acting  
23 Center Director, and also Jane Haney, Commissioner of the  
24 Food and Drug Administration. So with that, Dr. Burns,  
25 thank you for all of your service and, hopefully, you will

1 hang this plaque next to some of those approval orders you  
2 get from your own company's products.

3 [Applause]

4 I will turn the meeting back to you, Dr. Whalen.

5 DR. WHALEN: Thank you, Mr. Dillard, and  
6 congratulations, Dr. Burns.

7 We will now proceed with the open public hearing  
8 session of the meeting. All persons who wish to address the  
9 panel should speak clearly into the microphone as the  
10 transcriptionist is dependent on this means of providing an  
11 accurate record of this meeting.

12 We are requesting that all persons making  
13 statements during this open public hearing disclose whether  
14 or not they have any financial interests in any medical  
15 device company. Thus, before making a presentation to the  
16 panel, in addition to stating your name and affiliation,  
17 please state the nature of your financial interest, if any,  
18 or state that there is none. At the time that the agenda  
19 was drafted for this meeting there had been no formal  
20 requests made to FDA to present at this particular juncture  
21 of the meeting. Since there are no formal requests, I ask  
22 now that if there is anyone who wishes to address the panel  
23 that they please raise their hand.

24 Since there are no requests to speak in the open  
25 public hearing, we will now proceed to the open committee

1 discussion. I would like to remind the public observers at  
2 this meeting that while this portion of the meeting is open  
3 to your observation, public attendees may not participate  
4 except at the specific request of the panel. There will be,  
5 however, a further opportunity for the public to comment  
6 near the end of the meeting. We will now proceed to the  
7 sponsor's presentation.

8       **Applicant Presentation, Intuitive Surgical Incorporated**  
9                               **Introduction**

10               MR. DANIEL: Good morning, ladies and gentlemen,  
11 members of the panel and FDA. On behalf of Intuitive  
12 Surgical, we are pleased to have the opportunity to present  
13 today data and information in support of safety and  
14 effectiveness for the Intuitive surgical endoscopic  
15 instrument control system. My name is Mike Daniel. I am  
16 the vice president of regulatory and clinical affairs for  
17 the company.

18               [Slide]

19               We have with us today our president and CEO,  
20 Lonnie Smith; our founder and medical director, Dr. Fred  
21 Moll who will provide background for the system and discuss  
22 the technology. We have Dr. Guthart here with us today to  
23 answer any technical questions. He is the director of  
24 systems engineering. We have two of our four investigators.  
25 We have our principal investigator, Dr. Barry Gardiner from

1 Oakland, California, and we have our second U.S. clinical  
2 investigator, and that is Dr. Alan White from Tacoma,  
3 Washington. Dr. Gardiner will present the bulk of the  
4 clinical data today. We also have our statistician, Dr. Dan  
5 Bloch from Stanford, and Dr. Bloch will provide a brief  
6 statistical overview.

7 [Slide]

8 Intuitive Surgical was founded in 1995, with the  
9 purpose of developing computer-assisted technology. We have  
10 technology licenses from Stanford Research Institute, IBM  
11 and MIT. We have approximately 100 full-time employees,  
12 half of whom are engineers. We are currently marketing the  
13 system in Europe.

14 [Slide]

15 A brief summary of the regulatory status -- we  
16 obtained clearance from FDA via the 510(k) process in July  
17 of '97 for the following indication: Assistance in the  
18 accurate control of blunt dissectors, retractors,  
19 stabilizers and endoscopes in endoscopic surgical  
20 procedures. That would include laparoscopic as well  
21 thoracoscopic.

22 [Slide]

23 What we are here today to do is for a proposed  
24 additional indication, and that is, assistance in the  
25 accurate control of graspers, sharp dissectors, needle

1 holders, and electrocautery and accessories, and we are  
2 going to limit this to laparoscopic surgical procedures  
3 because that is the data we are presenting.

4 [Slide]

5 A quick overview of the regulatory process to date  
6 -- we submitted the 510(k) in December of 1996, and we  
7 obtained clearance in July of 1997. At about the same time  
8 the agency determined that clinical data would be necessary  
9 for sharp dissection, electrocautery and, in their minds and  
10 our minds conjunctively, surgical procedures as opposed to  
11 assistance.

12 We submitted an IDE and obtained conditional  
13 approval in July, 1998, and then proceeded with our clinical  
14 study and completed the 30-day follow-up required in  
15 December of 1998. We submitted a 510(k) in January of this  
16 year and last month FDA made the decision to convert that  
17 510(k) to a PMA.

18 [Slide]

19 Clearly, today's objective is to present data in  
20 support of that PMA approval.

21 [Slide]

22 I would like to introduce now our medical  
23 director, Dr. Fred Moll, who will describe the technology.

24 **Technology**

25 DR. MOLL: Good morning. I would like to first

1 thank the panel for the opportunity to present here today.  
2 By way of background, I am an M.D. and received residency  
3 training in general surgery. About five years ago I was  
4 introduced to the pioneering work of the Stanford Research  
5 Institute in the area of computer-assisted surgery. I  
6 became convinced that this approach could have a positive  
7 impact in minimally invasive technique, and it is this  
8 interest that led to the founding of Intuitive Surgical.

9 Before we present the study results today, I want  
10 to take you through a few minutes to introduce Intuitive's  
11 technology and the reasons for its development.

12 [Slide]

13 First a clarification of terms, computer-assisted  
14 surgery can be defined as a technique in which the surgeon's  
15 motion is assisted by a computer. How is this technique  
16 fundamentally different than conventional surgery? It  
17 interrupts the direct mechanical connection between the hand  
18 and the instrument, and inserts an electrical interface.

19 How can it improve existing technique? To explain  
20 this I want to discuss the advantages and disadvantages of  
21 both open and conventional laparoscopic technique, and  
22 suggest how computer-assisted surgery can provide clinical  
23 value.

24 [Slide]

25 In open surgery a large incision permits full

1 range of motion of the surgeon's hands and wrists inside the  
2 body cavity and close to the target anatomy. The  
3 disadvantage of this technique is simply that it requires a  
4 large incision.

5 [Slide]

6 In contrast, the advantages of conventional  
7 laparoscopy over open surgery relate directly to the  
8 dramatic reduction in incision size afforded by this  
9 technique. Advantages include reduced postoperative pain,  
10 shorter recovery, reduced hospital stay and, in many cases,  
11 lower health care costs.

12 [Slide]

13 However, the disadvantages of conventional  
14 laparoscopy relate to the fact that dexterous tasks become  
15 more difficult. The surgeon experiences reduced control  
16 because his or her hands are outside the body cavity and at  
17 the end of long instruments. Also, the body wall constrains  
18 movement and freedom of motion is reduced. Thus, control is  
19 necessarily transferred from the fingers and wrists to the  
20 surgeon's shoulders and elbows. Counter-intuitive motion is  
21 also a problem due to the body creating a fulcrum effect.  
22 The surgeon is forced to move his hands left to move the  
23 instrument tip right. Finally, the loss of eye-hand  
24 alignment and depth perception create challenge to good  
25 surgical technique.



1 [Slide]

2 The great attribute of computer-assisted surgery  
3 is that in very specific ways it allows enhanced instrument  
4 control. This enhancement is accomplished by, first,  
5 increasing range of motion by allowing Intuitive control of  
6 an articulating wrist inside the body cavity. Second,  
7 computer-assisted technique increases control by  
8 electronically shortening instruments and eliminating the  
9 fulcrum effect. Further precision is enhanced by the  
10 addition of motion scaling and tremor reduction. Finally,  
11 the Intuitive surgeon console provides 3D vision and eye-  
12 hand alignment.

13 [Slide]

14 I would like to show a video to illustrate some of  
15 these points. First, open surgery.

16 [Video presentation]

17 As you see in this video, the surgeon enjoys full  
18 range of motion for his or her hands and wrists because of  
19 the large incision used in open technique.

20 In laparoscopy, however, the surgeon's hands are  
21 at a significant distance from the tissue holding the end of  
22 long instruments. Here he or she is forced to control the  
23 instruments with shoulder and elbow movements rather than  
24 using the more dexterous movements of the fingers and  
25 wrists. In addition, the surgeon operates from a two-

1 dimensional video that is not in the same visual field as  
2 his or her hand movements.

3 In computer-assisted surgery, as shown here,  
4 Intuitive manipulators are able to provide precise motion at  
5 the instrument tip. The surgeon, positioned at a console,  
6 can control three manipulators with attached articulating  
7 tools such that the surgeon directs the tip of the  
8 instrument inside the body cavity. This increases range of  
9 motion, control and dexterity. In addition, the system  
10 console provides 3D vision and eye-hand alignment.

11 With the addition of motion scaling and tremor  
12 reduction, the surgeon's intended movements, at the bottom  
13 of the screen, are transmitted precisely to the system's  
14 instrument tips, at the top of the screen. The result is  
15 enhanced dexterity and precision. In addition, the system  
16 provides a means for efficient instrument change and the  
17 manipulators are able to remove the system from the  
18 operative field very quickly if necessary to do. With  
19 computer-assisted technique the instrument tips are afforded  
20 a much larger range of motion, provided in a variety of tips  
21 for different clinical purposes, and the system has a master  
22 that feeds electronically to instrument manipulators. These  
23 manipulators articulate at the instrument tip to provide  
24 control.

25 This, again, is the three-dimensional viewer that

1 provides eye-hand alignment and, as you can see, the motion  
2 of my fingers at the bottom of the screen are transmitted  
3 very precisely to motions of the instrument tip at the top  
4 of the screen.

5 This is an example of how instruments are changed  
6 in the system. If the manipulator needs to be moved quickly  
7 out of the surgical field, that is also quite easy.

8 The Intuitive System is being used clinically in  
9 general, gynecologic and cardiac surgery outside the U.S.,  
10 where it is approved for these indications. As seen in this  
11 footage, the enhanced dexterity offered by the system has  
12 helped accomplish precise suture placement in procedures  
13 such as mitral valve repair, as seen here.

14 [Slide]

15 In summary, the Intuitive computer-assisted  
16 surgery system enhances instrument control, first, by  
17 augmenting freedom of movement via articulating instruments  
18 which provide seven degrees of freedom inside the body  
19 cavity; second, by transferring hand control to instrument  
20 tips, thereby electronically shortening the surgeon's  
21 instrument; third, by eliminating counter-intuitive  
22 movement; and, fourth, by providing motion scaling, tremor  
23 reduction, coaxial eye-hand alignment and three-dimensional  
24 imaging.

25 [Slide]

1           The net effect is that the system enhances  
2 laparoscopic technique by recapturing many of the advantages  
3 of open surgery, thus, allowing more dexterous surgery in a  
4 minimally invasive format.

5           Having described the dexterous attributes of the  
6 system, it is important to understand that the study that we  
7 will present was not designed to quantify dexterity. We do  
8 believe that this study gives ample evidence that the  
9 Intuitive device is safety and efficacious, and that the  
10 enhanced dexterity can be appreciated from the video  
11 footage.

12           [Slide]

13           I would like now to ask Dr. Dan Bloch to present  
14 the statistical analysis of our clinical trial. Dr. Bloch?

15                           **Statistical Analysis**

16           DR. BLOCH: I have no financial interest in  
17 Intuitive Surgical. I have no stock options. I am being  
18 reimbursed for my expenses to attend this meeting, and I do  
19 get paid for the hours that I do consult for Intuitive  
20 Surgical.

21           [Slide]

22           What you are going to hear today are the results  
23 of two different studies, two separate studies, each one  
24 controlled. The design of the studies was as randomized,  
25 controlled clinical trials. The objective for both studies

1 was to show equivalence.

2           So, I am going to start by just giving a brief  
3 overview of tests of equivalence and how we defined the  
4 methods that we used for these two studies. The objective  
5 of the equivalence study is to show that the study group is  
6 clinically not significantly worse than the control group  
7 with respect to predetermined endpoints. Another way of  
8 saying this is that the aim is to demonstrate that outcomes  
9 of the treatment group and the control group are close  
10 enough so that the treatment and control group do not differ  
11 in a clinically important way.

12           What is different with clinical equivalence tests  
13 is that one has to define what differing in a clinically  
14 important way means. The methods that we employed to  
15 determine these parameters for the different endpoints were,  
16 first, prior to study we estimate the outcome average for  
17 the control group, denoted by the letter C in this slide,  
18 and we obtained that from the literature. Then in relation  
19 to that average value that we obtain as an estimate from the  
20 literature, then in consultation with our expert clinicians  
21 we determined the maximum amount that the test group average  
22 can attain and still be clinically acceptable, as denoted by  
23 T. Then simply the amount that is not considered to be  
24 clinically important is defined as either the difference or  
25 the ratio of these averages. I have used the Greek symbol

1 for delta to denote this difference, which statisticians  
2 love to use for parameters. Finally, the test of  
3 equivalence simply compares the observed data to this  
4 parameter delta.

5 [Slide]

6 We employed two different types of statistical  
7 tests for equivalence. The first is appropriate for  
8 success/failure data, yes/no data. For example, in our  
9 study the primary outcome for safety is complication. From  
10 that data we obtain a proportion for the group, either the  
11 treatment group or the control group, and using the  
12 published methods Blackwelder we simply take the difference  
13 between the two, and that is the delta.

14 Specifically, as an example, one of the studies  
15 has to do with what I will call the lap missing procedure.  
16 Dr. Barry Gardiner will explain these in more detail. I  
17 will defer questions of clinical meaning to him. I don't  
18 feel I am competent to answer those kinds of questions. But  
19 as an example, for lap missing study, there, from the  
20 literature and our best interpretation of what is right, the  
21 complication rate in the control group is estimated to be  
22 two percent. The acceptable complication rate for the test  
23 group is nine percent, so simply the difference of nine  
24 minus two is delta.

25 What this means is that if the difference between

1 the treatment and the control group is over seven percent,  
2 then we would say that is not clinically acceptable. That  
3 is too big to call the two treatment groups equivalent.

4 [Slide]

5 For continuous outcomes we use the ratio of the  
6 expected group averages rather than the difference. This is  
7 based on a method that is published by Dr. Fieller. As an  
8 example from the second study, having to do with lap choles,  
9 for procedure duration the literature estimate for the  
10 amount of time for the surgical procedure was 89 minutes,  
11 and the clinically acceptable procedure duration deemed by  
12 the clinicians was at least double that, and I have written  
13 down 178 minutes here, which is double 89. Here the delta  
14 is simply the ratio, 178 divided by 89.

15 This methodology is in the protocol for continuous  
16 outcomes, and it is especially useful when the actual  
17 control observations differ from the prior assumptions.  
18 Recall that our assumptions having to do with the control  
19 means come from the literature, and when we do the  
20 experiment it might be quite different from that. We hope  
21 we got a good estimate because our estimate of delta depends  
22 on it.

23 But, as an example, procedure duration for the lap  
24 chole control arm was 69 minutes. That actually was the  
25 average, quite different from 89 minutes. The ratio is

1 especially useful when such a situation occurs, that is, the  
2 ratio of two still makes sense in this case.

3 [Slide]

4 Some of the delta values in the tables were  
5 incorrectly specified. In what follows I will try to  
6 present a fair and balanced presentation of the assumptions  
7 and try to clearly indicate that the deltas were  
8 predetermined. There was a disconnect between the text and  
9 the tables. Unfortunately, I did not have an opportunity to  
10 review the tables and text before the protocol was submitted  
11 to the FDA, and did not discover the discrepancies until  
12 later, however, upon careful reading of the protocol as a  
13 whole, I think it is quite clear what our intent is and this  
14 is what I wish to review now in the next slides.

15 Regarding the next two bullets, as I have already  
16 indicated, sometimes the control averages were quite  
17 different that we saw in this study results and that does  
18 have a bearing on the actual test of statistical  
19 equivalence. We will see that later. Finally, the testing  
20 of equivalence methods that I have described are not valid  
21 if we have low outcomes, low average outcomes and we will  
22 see that for conversions and complications, in fact, we have  
23 no outcomes and you can't then estimate the standard error.

24 [Slide]

25 Again, I just want to defer to Dr. Gardiner any



1 clinical interpretations regarding the outcomes that I am  
2 going to present.

3 [Slide]

4 The primary outcome variable in both groups was  
5 conversions, that is the need to stop using the method that  
6 the patient was randomized to and use a different method to  
7 complete the procedure.

8 I am going to present the conversion results in  
9 two ways, one which includes data before the system was even  
10 brought to the operating table, and what is involved with  
11 that will also be presented by Dr. Gardiner. He will make  
12 that clear as to what part of the surgical procedure does  
13 take place before the Intuitive System is even in place.

14 In that case, the protocol did not specify test of  
15 equivalence. Only note here that in this case, in this  
16 period of the intervention there were two conversions in the  
17 lap chole arm and one in the lap Nissen arm, both using the  
18 Intuitive procedures, two and one occurrence respectively in  
19 these proportions, 3.5 percent and 1.7 percent.

20 [Slide]

21 Our intent was to test for equivalence for  
22 conversions when the Intuitive instruments were actually  
23 used, and that is what I mean by "during the ISI segment."  
24 Clearly, you can't have an Intuitive System actually in  
25 place in the control group but there is a corresponding

1 surgical interval and that is what is meant here.

2 In this particular case there were no conversions  
3 either in the control group or in the Intuitive group for  
4 either study. The N/A for the p value simply means that one  
5 can't calculate a test of substantial equivalent; you can't  
6 divide by zero; the standard error is zero.

7 [Slide]

8 The other primary outcome for efficacy is only  
9 applicable in the lap Nissen study, called the DeMeester  
10 score, and this will be described by Dr. Gardiner. The  
11 delta, which is in the middle of the table -- now, this is  
12 continuous outcome -- is 14.8 divided by 6. Notice at the  
13 top the estimated control value was equal to 6; the 14.8 is  
14 what DeMeester published as the thresholds. If a person has  
15 a DeMeester score over 14.8, that is considered to be  
16 abnormal and a score lower than 14.8 is considered to be  
17 normal.

18 Notice that at the bottom of this transparency I  
19 have something called the minimum delta. This is the delta  
20 so that if you take the evidence from the trial into  
21 consideration the p value would be exactly equal to 0.05.  
22 That is, if a delta is prespecified and equals 1.71 or  
23 larger, then the difference between the observed ratio,  
24 which is directly above that bottom line, 1.02, and 1.07,  
25 the data would support equivalence at the 0.05 level.

1 Again, I want to remind you that low p values means that we  
2 have established equivalence.

3 Notice in the middle of the table again that the  
4 actual ratio, 2.5, is considerably larger than 1.71 and, of  
5 course, that is why the p value is so small. It is less  
6 than 0.001.

7 [Slide]

8 The primary safety outcome is complications.  
9 Again, we have the period prior to the instrumentation being  
10 used, even being brought to the table. This was not  
11 intended for a test of equivalence and that is why the delta  
12 is not specified. But for information, we did observe 2/57  
13 cases in the lap chole group versus 3/60 in the lap Nissen  
14 groups in the Intuitive treated arms versus zero in the  
15 control arm for the lap chole and 5 percent in the lap  
16 Nissen arm.

17 [Slide]

18 What was meant is to compare complications as a  
19 primary outcome again during the segment in which the  
20 Intuitive instrumentation was actually in place. As with  
21 conversions, there were no occurrences in either arm for  
22 either study.

23 [Slide]

24 One of our secondary outcomes is procedure  
25 duration. I would like to read this quote that is in the

1 protocol: "Our clinical consultants maintain that for  
2 laparoscopic procedures the operative time would need to  
3 double -- that is where I have 89 times 2 or 178 in the  
4 previous example -- before we would begin to be clinically  
5 significant. For purposes of final analysis, a conservative  
6 delta -- and I would like to underline "conservative" -- of  
7 45 for lap chole and 50 for lap Nissen minutes will be used.

8 [Slide]

9 You will notice in the middle of this table that  
10 there are two deltas. There is a ratio delta, and I said  
11 conservative, and a ratio delta equal to acceptable. In the  
12 protocol, and this is confusing, we have 1.5 as the ratio  
13 delta to be considered for the lap chole arm, but we also  
14 have in this quote that we could use the delta equal to 2.

15 Now I want to draw your attention to the bottom of  
16 this transparency where, again, the minimum deltas are  
17 listed, 1.76 and 1.71. Again, the idea is very simple. If  
18 the delta that we prespecified was greater than 1.76 then  
19 the evidence of the data would support equivalence. If we  
20 had chosen a value that was lower than that -- and,  
21 remember, this is a parameter you choose up front -- then it  
22 was too conservative. But we admit that we were choosing a  
23 conservative one so I hope you will bear with this in this  
24 regard. Notice that the p value is less than 0.001. We  
25 established equivalence used in the deltas of 2.

1 [Slide]

2 The final outcome that I want to refer to is  
3 postoperative stay, again a secondary outcome. The  
4 estimated values from the controls in the literature were  
5 two days and three days on average for the lap chole and lap  
6 Nissen studies. This is a continuous outcome so we are  
7 using a ratio delta, and the clinicians deemed that an  
8 average of three days for the lap chole group was  
9 acceptable, hence the ratio of 1.5 -- 3 divided by 2, and  
10 also the ratio of 1.5 for lap Nissen, which means that up to  
11 4.5 days would be clinically acceptable for postop stay for  
12 lap Nissen cases, again with a delta of 1.5. In both cases  
13 equivalence is established, noting again that the p values  
14 were less than 0.05.

15 [Slide]

16 I think now I would like to stop this presentation  
17 and Dr. Gardiner will present the study findings and speak  
18 to the clinical interpretation of the data.

19 **Clinical Study and Results**

20 DR. GARDINER: Good morning. I would like to  
21 thank the FDA and all the panel members for the opportunity  
22 that they have given us to make this presentation today.

23 I am a general surgeon, and my focus today is  
24 going to be on the clinical importance of the data rather  
25 than the statistical one. I am a surgeon in northern

1 California with an active clinical practice in advanced  
2 laparoscopic surgery. I have been a consultant to Intuitive  
3 Surgical since April of 1996 and at the outset, and in the  
4 interest of full disclosure I should mention that the  
5 company does compensate me for my time through a consulting  
6 agreement, and I do have a small equity position in the  
7 company.

8 In July of last year the company began a clinical  
9 trial, clinical study to demonstrate that the Intuitive  
10 System can be used safely and effectively in performing  
11 laparoscopic surgery.

12 [Slide]

13 I was the principal investigator for that trial  
14 and over a three-month period we completed 129 laparoscopic  
15 cholecystectomies and laparoscopic funduplications using  
16 this system. We did 12 training cases and 117 randomized  
17 cases as part of the clinical trial.

18 As we begin to examine some of this data today, we  
19 are going to see that there are some differences in some of  
20 the data and some of the values between the control groups  
21 and the study groups. When these differences are put into  
22 clinical context, and the appropriate clinical context, I  
23 think you will believe and see that none of these  
24 differences are actually clinically important, and that the  
25 data does adequately support the conclusion that this system

1 can be used safely and effectively in both basic and  
2 advanced laparoscopic procedures.

3 [Slide]

4 This was a prospectively randomized clinical trial  
5 concurrently run that included 12 training cases and 233  
6 randomized patients that were operated on between July and  
7 October of last year. Now, there were two procedures that  
8 we evaluated in this trial. One was the cholecystectomy, a  
9 procedure that is probably the most commonly done and one of  
10 the most basic of laparoscopic procedures. The Nissen  
11 fundoplication is a procedure requiring advanced  
12 laparoscopic skills and suturing techniques. In each arm of  
13 the study, both the cholecystectomy arm and the Nissen arm,  
14 the patients were randomized into one of two groups, a  
15 control group in which the procedure was done with  
16 conventional laparoscopic technique and traditional  
17 instruments and a study group in which the procedures were  
18 done using the Intuitive surgical device.

19 In prospective collaboration with the FDA, the  
20 company chose to do this trial at Hospital Torre Medica, in  
21 Mexico City. This is a three-month trial and the company  
22 was unable to locate any center in the United States in the  
23 available time that had a sufficient number of patients with  
24 sufficiently severe, untreated gastroesophageal reflux  
25 disease to meet the inclusion criteria. We had a full-time

1 clinical monitor present on site throughout the study to  
2 ensure accurate and complete data collection and patient  
3 follow-up.

4 [Slide]

5 Four completely separate surgical teams  
6 participated in this study, each with their own assistant  
7 and their own scrub nurse. These four surgeons were all  
8 very highly experience laparoscopic surgeons. On average,  
9 each of them had done over 1000 laparoscopic  
10 cholecystectomies and 400 Nissen funduplications using  
11 conventional laparoscopic technique and traditional  
12 laparoscopic instruments. I think you will see as we get  
13 into the data that this will need to be considered as you  
14 interpret some of the data that we present today.

15 [Slide]

16 This is a photograph of the operating room at  
17 Hospital Torre Medica. The surgeon's console is over on the  
18 left. The patient is on the right, undergoing a  
19 cholecystectomy. The Intuitive device is in place here and  
20 the assistant is positioned between the patient's legs with  
21 the scrub nurse to his direct right.

22 [Slide]

23 In the cholecystectomy arm of this trial only  
24 patients with symptomatic gall stones, documented on  
25 ultrasound, were included. Patients were randomized



1 preoperatively, not intraoperatively. As a result, there  
2 were no intraoperative exclusions specified in the protocol.  
3 So, in essence, this data is being analyzed on an intent-to-  
4 treat basis. Morbidly obese patients were excluded, as were  
5 those requiring emergent surgery because of acute  
6 cholecystitis. In an attempt to try and reduce variability  
7 between treatment groups and control groups, patients with  
8 the suspicion of common duct stones that would need a common  
9 duct operation in addition to a cholecystectomy were  
10 excluded preoperatively. Finally, any patient that had a  
11 relative or absolute contraindication to having their  
12 disease treated laparoscopically were excluded from the  
13 trial preoperatively.

14 [Slide]

15 In the fundoplication arm only patients with  
16 symptomatic gastroesophageal reflux disease were included.  
17 This required that they both had endoscopic-proven, biopsy-  
18 proven esophagitis, and that the DeMeester score was over  
19 14.8. For those of you not familiar with this test, it is a  
20 measure of the severity of acid reflux from the stomach up  
21 into the esophagus.

22 As in the cholecystectomy arm, in an attempt to  
23 try and reduce patient variability between control and study  
24 groups, we excluded patients with morbid obesity, patients  
25 with intrinsic esophageal disease, and patients with

1 periesophageal hernias since that is fundamentally a  
2 different condition and, in addition, like the laparoscopic  
3 chole series, any patient with a relative or absolute  
4 contraindication for laparoscopic was excluded. As in the  
5 cholecystectomy arm, there were no intraoperative exclusions  
6 defined.

7 [Slide]

8 Remember, this study was done from a clinical  
9 point of view to demonstrate that the Intuitive System is  
10 safety and effective in performing laparoscopic surgery.  
11 That is what the object was. The endpoint in both the  
12 cholecystectomy and fundoplication arms to establish safety  
13 was equivalent complications, equivalent complications  
14 occurring in the study group compared to those occurring in  
15 the control group. I believe you will see that we met this  
16 endpoint from a clinical standpoint and, in fact, in many  
17 ways patients in the study group fared considerably better  
18 than those in the control.

19 We also had a series of endpoints related to  
20 effectiveness. In the cholecystectomy arm the primary goal  
21 was to successfully remove the gallbladder without  
22 conversion to either conventional laparoscopic technique or  
23 open surgery.

24 There was a series of secondary endpoints --  
25 equivalent procedure time, equivalent postoperative hospital

1 stay, and comparable quality of life scores. This is a  
2 psychological well being test that is a measure of quality  
3 of life.

4 [Slide]

5 The endpoints for safety and effectiveness for the  
6 fundoplication arm were identical to those in the  
7 cholecystectomy arm, with one exception. There was an added  
8 primary effectiveness endpoint that was the equivalent  
9 reduction in the DeMeester scores following surgery between  
10 the control and the study groups. Otherwise, the Nissen arm  
11 and the cholecystectomy arm were the same.

12 [Slide]

13 You are going to see some differences between the  
14 control group and the study group in several of these  
15 effectiveness endpoints, primarily as it relates to  
16 procedure duration. I believe we understand why some of  
17 those differences occurred, if not all of them, and we will  
18 go into that in a little bit more detail in just a little  
19 bit. But you will see that none of these differences are  
20 clinically important, and certainly not sufficient to  
21 conclude that this device isn't clinically effective.

22 The evaluation of safety of this surgical device  
23 is based on a comparison of complications that occurred  
24 between the study group and those that occurred in the  
25 control group. It, therefore, becomes very important that

1 we understand why each and every one of these complications  
2 occurred in this clinical trial. To do so, we need to know  
3 precisely what part of each procedure the study device was  
4 used for.

5 [Slide]

6 In each cholecystectomy the abdomen was initially  
7 insufflated and the trocars were placed in a traditional and  
8 usual fashion. In the control group all the remaining steps  
9 of the cholecystectomy were done with conventional  
10 laparoscopic instruments, using traditional laparoscopic  
11 techniques. In the study group it was only after the  
12 trocars were in place that the Intuitive System was  
13 introduced into the patient. Once that had been done, the  
14 remaining steps of the operation were done with the  
15 Intuitive System. This is the reason that trocar injuries  
16 have nothing to do with either the use of or the safety of  
17 this device because the trocars were placed before the  
18 Intuitive System was ever introduced into the patient.  
19 After the gallbladder has been dissected from the liver bed  
20 the system is removed, the gallbladder is removed with  
21 traditional instruments, and the trocars are then removed  
22 and the incision is closed in the usual fashion.

23 [Video]

24 The study device in this video is being used on  
25 the left and standard instruments are being used on the

1 right. The gallbladder is initially freed from the  
2 surrounding structures here. Clips sufficiently large to  
3 occlude the cystic duct and the artery that goes to the  
4 gallbladder were unavailable in the Intuitive System at the  
5 time this study was done.

6 For this reason, once the duct and artery were  
7 completely dissected they were clipped on the patient side,  
8 and you may actually see a clip here and over here, using  
9 the conventional clip applier in both the study group and  
10 the control group, and then both were doubly tied using  
11 intracorporeal knotting and suturing techniques. The  
12 gallbladder was then separated from the liver bed, after the  
13 cystic duct and artery are divided, in the study case being  
14 done with blunt dissection and cautery with control. You  
15 can see the advantage that the articulated end of this  
16 instrument has in terms of directing the tip of the  
17 instrument to the tissue. In the study group the Intuitive  
18 System is removed after the gallbladder has been extracted  
19 using conventional instruments.

20 [Slide]

21 In assessing the comparative complications in the  
22 fundoplication arm of this study, it is just as important as  
23 it was in the cholecystectomy arm that we understand what  
24 part of each procedure the device was used for, and you will  
25 see why that is important as we get into the data. In each

1 fundoplication procedure the abdomen, just as in the  
2 cholecystectomy, was insufflated in the usual fashion;  
3 trocars were placed, as they would be in traditional  
4 laparoscopic surgery.

5           Now, harmonic scalpel was not available on the  
6 study device at the time this trial was done and, therefore,  
7 in both study groups and control group the short gastric  
8 vessels that tether the stomach to the spleen were  
9 cauterized, divided and taken down as the next step in the  
10 operation. So, we put the trocars in, divide the short  
11 gastric vessels with the harmonic scalpel, and this is all  
12 done using traditional laparoscopic instruments and is  
13 identical in both groups, whether it is the control or the  
14 study group.

15           In the control group the remaining steps of the  
16 operation are then completed using traditional instruments  
17 and conventional technique. In the study group it was only  
18 after the division of the short gastric vessels that the  
19 system was introduced into the patient. This is the reason  
20 that neither trocar injuries nor harmonic scalpel injuries  
21 produced during the operation have anything to do with  
22 either the use of or the safety of this device because the  
23 device hasn't been brought to the table prior to the trocars  
24 and the short gastric vessels being divided with the  
25 harmonic scalpel.

1           Once in place, the system is then used to do the  
2 rest of the operation and the fundoplication is completed,  
3 anchored to the crura, and then once the operation is  
4 finished the system is removed and then the trocars are  
5 removed and the incision is closed in the traditional and  
6 conventional fashion.

7           [Video]

8           A study case is being done on the left, a control  
9 case being done on the right. The initial step of the  
10 operation in both study group and the control group is to  
11 divide the short gastric vessels with the harmonic scalpel,  
12 and that is being done here. So, you can see that the  
13 system at this point has not been brought into the field.  
14 It hasn't even been brought up to the operating table at  
15 this point, and the control and the study groups are being  
16 done with the same instruments.

17           Once the Intuitive System is brought into the  
18 field, the right crus is dissected and, once the hiatus has  
19 been freed up, the fundus is brought underneath the  
20 esophagus and that is happening here. This is the  
21 esophagus. Then the fundal wrap is sutured together in  
22 front of the esophagus here using intracorporeal suturing  
23 techniques and intracorporeal knotting. You can see once  
24 again the advantage that this articulated instrument gives  
25 you here in terms of having the tip of the instrument

1 directed toward the tissue in the proper orientation. The  
2 crura are then closed and the wrap is then anchored to the  
3 crura, and then in the study group the Intuitive System is  
4 removed and the trocars, and wounds are closed in the  
5 traditional fashion.

6 [Slide]

7 There was a total of 245 patients enrolled in this  
8 trial, 129 patients in the study group and 116 controls. Of  
9 the patients that were enrolled on the Intuitive side, there  
10 were 12 training cases done with the Intuitive System,  
11 leaving a total of 233 patients that were actually  
12 randomized, essentially equally distributed between the  
13 control and the study groups.

14 [Slide]

15 There were no meaningful differences between the  
16 study group and the control group with regard to age, body  
17 mass index or PGWB scores, and on the fundoplication arm,  
18 which is the only arm that has relevance to the DeMeester  
19 score, the preoperative DeMeester scores were essentially  
20 identical.

21 [Slide]

22 We completed 245 procedures in three months,  
23 roughly divided equally between the two arms of the study  
24 and between the control and study groups within each arm.  
25 There is some distribution and variation between the number



1 of cases surgeons did and the distribution of those cases,  
2 but those differences were based on patient availability at  
3 the time the surgeon was present and these differences are,  
4 as you can see, relatively minor.

5           There were only 12 training cases, 5 on the Nissen  
6 side and 7 on the cholecystectomy side, that were done by  
7 these 4 surgeons prior to randomization of the patients into  
8 the trial. This meant that at the outset of the study none  
9 of us had any appreciable clinical experience with this  
10 device. The average number of study cases that each surgeon  
11 did during the trial was about 29 so that even by the end of  
12 the trial our experience was still quite limited with this  
13 device, and this may well have contributed significantly to  
14 the procedure duration, as you will see later on.

15           [Slide]

16           Let's take a look at the results from this trial.  
17 Remember that safety was defined in both arms of this study  
18 as equivalent complications between the control group and  
19 the study group. I will only talk about complications that  
20 we observed in this study.

21           [Slide]

22           This is an overview of every adverse event that  
23 was experienced in these 245-patient study, ranked in  
24 decreasing order of severity from top to bottom. The  
25 patient who sustained the gastric perforation had by far and

1 away the most serious complication in this series, and this  
2 complication ultimately resulted in her death. This was a  
3 patient in the control group, and the complication was due  
4 to an unrecognized injury to the stomach, felt to be related  
5 to the use of the harmonic scalpel. This was recognized 24  
6 hours following surgery. Multiple reoperations were  
7 required for abscess drainage, and this patient ultimately  
8 developed adult respiratory distress syndrome and died on  
9 the 83rd postoperative day. That is the patient in the  
10 control group, the most serious complication.

11           There were two treatment failures, also both in  
12 the control group. Both were due to migration of the wrap  
13 around the esophagus, having the wrap migrate up into the  
14 mediastinum. The first of these, this perioperative wrap  
15 migration, occurred 8 hours following surgery. It was  
16 recognized and the patient was immediately taken back to the  
17 operating room and re-laparoscoped. The wrap was taken down  
18 but the fundus was found at the time of the second operation  
19 to have already developed an area of necrosis. This was  
20 resected laparoscopically, and the patient was discharged on  
21 the 21st postoperative day and has done well. So, that is  
22 the first treatment failure, the perioperative wrap  
23 migration.

24           The second treatment failure, also in the control  
25 group, was a late wrap migration. This patient is an opera

1 singer and the patient's reflux symptoms recurred about two  
2 months following the original operation. The evaluation  
3 revealed that the wrap had migrated into the chest, and this  
4 patient required reoperation three months later. This  
5 patient is now reported to be asymptomatic and doing well.  
6 This wrap migration problem is something that we are  
7 becoming increasingly aware of as a recognized complication  
8 of the Nissen procedure.

9           There were two trocar injuries to the bowel in  
10 this study, one in a training case and one in a randomized  
11 patient in the trial. Neither one was related to the study  
12 device since the trocars were put into these patients before  
13 the system was even brought up to the operating room table,  
14 and that is why I spent so much time emphasizing what part  
15 of this operation was done by the study device.

16           The injury to the small bowel was recognized 24  
17 hours after the original procedure and, again, this is a  
18 training case. This patient required reoperation and repair  
19 of the enterotomy but she had an uneventful postoperative  
20 recovery and was discharged on the 8th postoperative day,  
21 and she is doing well.

22           The injury to the stomach -- this one -- occurred  
23 in one of the randomized patients, randomized to the study  
24 group. This was recognized intraoperatively, repaired  
25 laparoscopically, and the patient experienced no untoward

1 event.

2           There was one patient in the study group that  
3 developed bleeding from the trocar site, is this one, and  
4 was returned to the operating room. This was actually my  
5 patient. The wound was re-explored locally and a vessel in  
6 the abdominal wall was suture ligated. On review of the  
7 tape following this operation, it was clear that the  
8 bleeding was related to the initial insertion of the trocar  
9 site through the abdominal wall and not related to the study  
10 device itself. This patient also had an uneventful  
11 postoperative recovery.

12           There were three minor complications that occurred  
13 on the study side of this trial. There were two minor  
14 serosal tears to the stomach related to use of the harmonic  
15 scalpel. These were reinforced laparoscopically with  
16 sutures and were of no consequence. Once again, these  
17 occurred before the system was brought up to the patient.

18           There was a superficial umbilical wound infection.  
19 That was also my patient. The gallbladder ruptured in this  
20 patient as we were bringing it up through the umbilical  
21 incision and stones and some bile were spilled in the  
22 incision and I think, most likely, that is what led to this  
23 infection. It was a trivial infection and was managed with  
24 just warm compresses.

25           [Slide]

1           In summary, there were nine complications in these  
2 245 patients. By far and away, the three most serious ones  
3 occurred in the control group, including one death and two  
4 treatment failures on the fundoplication arm of this trial.  
5 There were five complications that occurred in the study  
6 group, and none was related to use of the study device.  
7 Four occurred before the system had been brought up to the  
8 table, and two of these were related to trocar placement and  
9 two were related to use of the harmonic scalpel, and there  
10 was one minor wound infection with a known antecedent cause.  
11 There was one trocar injury to the small bowel that occurred  
12 in a training case prior, again, to introduction of the  
13 system into the patient.

14           [Slide]

15           Concluding this, from a clinical standpoint, it  
16 doesn't appear that there were any device-related  
17 complications associated with the use of this system. All  
18 of the complications were related to the use of conventional  
19 laparoscopic instruments and, in fact, the control group  
20 actually fared appreciably worse in this regard than did the  
21 study group. The data establishes that the safety endpoints  
22 for this device have been met.

23           [Slide]

24           Well, we have addressed the issues of safety of  
25 this device, now let's deal with its effectiveness. We have

1 defined the primary effectiveness of this device by its  
2 ability to successfully remove the gallbladder and complete  
3 the fundoplication without conversion to either conventional  
4 laparoscopic or open surgical technique. Then as a  
5 secondary endpoint there was another endpoint for  
6 equivalency which was a reduction in the DeMeester score on  
7 the Nissen side.

8 [Slide]

9 Let's look and see how we did with regard to these  
10 effectiveness endpoints. There were 12 training cases done  
11 with the study device; 117 patients were randomized into two  
12 groups to have their operations performed with the Intuitive  
13 System. All 12 training cases and 114 of the 117 randomized  
14 patients had their operation successfully completed using  
15 the Intuitive device.

16 There were three remaining cases that had been  
17 randomized into the two study groups but whose procedures  
18 were not completed according to this random assignment.  
19 These three cases were converted either to an open  
20 laparotomy or conventional laparoscopic before the system  
21 was ever brought to the operating table. So, there were no  
22 conversions after the system had actually been placed into  
23 use.

24 [Slide]

25 So we can truly understand what happened in these

1 patients, because I think it does go to the issue of  
2 effectiveness and these primary effectiveness criteria,  
3 let's go into a little more detail about what happened in  
4 these three cases. There were three patients that had been  
5 randomized into the study arm but in whom the Intuitive  
6 System actually was never even brought to the operating room  
7 table.

8           The first case was actually my patient. We  
9 inserted the laparoscope and immediately found severe  
10 macronodular cirrhosis, and this patient had extensive  
11 portal hypertension. I just felt there was an absolute  
12 contraindication to proceeding with any form of laparoscopy  
13 and I basically removed the laparoscope and immediately  
14 converted this patient to an open laparotomy and completed  
15 the patient in the traditional open fashion.

16           There were no intraoperative exclusions defined in  
17 the protocol so technically, for the purposes of this  
18 analysis, this patient has been considered a conversion.  
19 The study device, however, was never actually inserted into  
20 this patient, and had the pathology of this patient been  
21 detected preoperatively he wouldn't have been allowed to  
22 participate in the clinical trial in the first place.

23           The other two cases that were not completed  
24 according to random assignment were completed using  
25 conventional laparoscopy. Let's talk about those two. One

1 of these patients had severe acute cholecystitis, severe  
2 acute inflammation and scarring in the porta hepatis, and  
3 the other had a very significant nodular cirrhosis. These  
4 two cases were Dr. White's, and at this point in the trial  
5 he had done one training case with the system. These had  
6 been his first and his third randomized patients in this  
7 trial using this device. Given the severity of the  
8 pathology he was dealing with, and his inexperience with the  
9 system, he made a clinical judgment that it was in the best  
10 interest of his two patients that he use instruments and  
11 techniques with which he was already very familiar.

12 He converted both of these cases without ever  
13 bringing the system to the patient's table because of his  
14 position on the learning curve and the severity of this  
15 pathology. Dr. White is here today with us, if any of you  
16 have additional questions about either of these two cases.

17 [Slide]

18 With these three cases in mind, let's review what  
19 the protocol actually says about conversions. The endpoint  
20 for primary effectiveness of this device was defined in the  
21 protocol as successful completion of surgery without  
22 conversion to either conventional laparoscopy or open  
23 technique. That is what the protocol says.

24 Now, this definition was well intentioned, but in  
25 retrospect it fails to take into account that conversions



1 can occur and may occur for reasons other than those related  
2 to the device itself, and that is exactly what happened in  
3 these three cases, and it wasn't anticipated. This has led  
4 to a curious situation of having us rely on cases, in which  
5 the system was never actually used, to draw conclusions  
6 about how well the device works.

7           Since what we are trying to do here is to evaluate  
8 how effectively a surgical device can be used in surgery, it  
9 would seem that a better definition for conversion would be  
10 terminating the use of the device once the surgeon has  
11 started operating with it. Interpretation of the protocol  
12 in this fashion would lead to a more meaningful evaluation  
13 of the effectiveness of this device since it would rely on  
14 cases in which the device was actually used instead of those  
15 in which it was never used.

16           If the protocol is interpreted in this fashion  
17 there were no treatment-related conversions in this trial.  
18 This seems, to me, to be a more clinically appropriate  
19 definition and, given this interpretation, this device  
20 clearly met its primary effectiveness endpoint with respect  
21 to conversion.

22           On the fundoplication arm of the study there was  
23 an additional primary endpoint for effectiveness, and that  
24 was equivalent reduction of DeMeester scores following  
25 surgery in the control and the treatment groups.

1 [Slide]

2 There was no difference between the study group  
3 and the control group in either the percent reduction of the  
4 DeMeester score following surgery, 65 percent and 69 percent  
5 respectively, or in the number of patients whose  
6 postoperative DeMeester score returned to normal, in this  
7 case 39 patients in both study and control groups.

8 [Slide]

9 This slide shows you a graphic representation of  
10 the DeMeester score data from the preoperative data to the  
11 postoperative data, with the control group on the left side,  
12 the study group in the middle and DeMeester's classic data  
13 that was published in June of 1995 on the right. It is  
14 clear that there is no real substantial difference in any of  
15 these three groups.

16 [Slide]

17 So, I think the data does show that this device is  
18 effectiveness in carrying out basic and advanced  
19 laparoscopic surgery. Every procedure that was begun with  
20 the device was successfully completed with the device. In  
21 the fundoplication arm of the study the reduction of  
22 DeMeester scores was equivalent. There were no conversions  
23 that occurred that in any way we could attribute to a system  
24 failure, that the system wouldn't do or couldn't do what it  
25 needed to do or what it was designed to do, or that the

1 surgeon started with the system and couldn't get the  
2 operation finished with the system. Those were the issues  
3 that we felt we were addressing in conversions.

4           So, in my opinion, from a clinical point of view,  
5 had any of those issues occurred it would have gone straight  
6 to the question of efficacy of this device, but none of  
7 those issues did occur. What did happen were conversions  
8 related to discovering a contraindication to laparoscopy or  
9 the insecurity of a surgeon being randomized to operate with  
10 a brand-new device in a very difficult situation. The data  
11 shows the device has met its primary effectiveness endpoint.

12           [Slide]

13           There were three secondary effectiveness endpoints  
14 that were defined in the protocol, procedure duration,  
15 postoperative hospital stay and quality of life scores. In  
16 looking at the procedure duration, on average it took us  
17 about 40 minutes longer to do the cholecystectomies in the  
18 study group with the study device than with control  
19 instruments, and it took about 50 minutes longer to do the  
20 Nissens with the study device compared to standard  
21 traditional laparoscopy.

22           Now, I don't believe those differences are  
23 sufficient to have a negative clinical impact on patient  
24 care or treatment outcomes. They are just not sufficient to  
25 be clinically important. But as a surgeon, and having spent

1 a significant amount of time operating on this system, I  
2 have to ask what in the world was it that caused this  
3 difference. Why did it take us longer to complete the same  
4 procedures using the Intuitive System compared to the  
5 control traditional laparoscopic instruments? It is  
6 actually a very interesting question.

7 [Slide]

8 I have looked at a great number of the tapes of  
9 all four of us operating, and in my opinion there are a  
10 number of factors that have contributed to these longer  
11 procedure times. Probably the most significant factor was  
12 that this was a completely new device being used clinically  
13 for the first time by all four of these surgeons. As you  
14 watch these tapes, it is pretty clear that all four of us  
15 were gaining confidence in the use of the system as the  
16 study progressed. We were obviously more comfortable, all  
17 four of us, by the end of the study than we were at the  
18 beginning, but it is also obvious that all four of us were  
19 still learning how to use this device. Historically,  
20 operating times for every laparoscopic procedure, including  
21 the lap chole, have fallen as experience has been gained,  
22 and it will be no different for this device.

23 Secondly, it is clear when you actually watch  
24 these videos, and I have watched a lot of them, that the  
25 movements made by the surgeons with the device are

1 definitely more deliberate; they appear to be more  
2 meticulous, more precise, but they take longer to make.

3           The third factor I believe is that having the  
4 surgeon sit next to the patient at a console changes the  
5 dynamic between the surgeon and the assistant. Surgeons  
6 using this system need not only to learn how to use the  
7 system but they will also need to learn how to use their  
8 assistant in a more effective and the most efficient manner.  
9 This is going to take some time and some experience, I  
10 think, with the system to do that, and I think it is going  
11 to be a unique and different process for each individual  
12 surgeon as they work through the experience curve.

13           Finally, time is required to set up the device, 7  
14 minutes for the cholecystectomy and an average of 14 minutes  
15 for the fundoplication, and this needs to be added to the  
16 procedure duration.

17           So, as you might expect, procedure duration was,  
18 in fact, affected by the use of this device and it is a  
19 multifactorial issue. But in the final analysis, even if  
20 the procedure duration stays the same, which it almost  
21 surely will not, it just isn't sufficient to be clinically  
22 important.

23           [Slide]

24           Now let's look at length of stay. According to  
25 the protocol this outcome has been analyzed on the basis of

1 mean length of stay. On this basis there is no significant  
2 difference between the control group and the study group in  
3 the cholecystectomy arm of this study, being 1.3 days for  
4 each group. In the fundoplication arm, this arm here, the  
5 length of stay is artificially high in the control group,  
6 1.4 days for the study group and 3 days for the control  
7 group. This longer length of stay is because there were 2  
8 patients in the control group who suffered major  
9 postoperative complications, and when you take these 2 cases  
10 into account and you analyze the data that way there just  
11 aren't any clinically important differences between control  
12 and study group in either arm of this study with regard to  
13 length of stay.

14 [Slide]

15 The improvement in quality of life scores  
16 basically show no significant differences between study and  
17 control groups but there were some additional observations  
18 that were noted during this study that, although they were  
19 not defined in the protocol as either safety or  
20 effectiveness endpoints, deserve some comment.

21 [Slide]

22 Dysphagia is a common and well-recognized  
23 complaint in the early postoperative period in patients  
24 undergoing fundoplication. There were eight patients in  
25 this trial that underwent dilatation following surgery,

1 three in the control group and five in the study group, and  
2 all of those patients are now asymptomatic. These  
3 differences are not statistical significant.

4 Blood loss in this trial was trivial in both  
5 control and study groups and in both arms of the study.  
6 That having been said, there was an observed difference in  
7 the average blood loss between control and study group in  
8 both arms of the study. There was a difference of about 13  
9 cc in the cholecystectomy arm and not quite 12 in the  
10 fundoplication arm. Now, this amount of blood, just to put  
11 it into perspective, is a little over 2 teaspoons, and it is  
12 actually less blood than the amount of blood that was drawn  
13 from these patients to do their preoperative lab testing.  
14 This amount of blood loss is just clinically unimportant.

15 [Slide]

16 The system performed, with regard to reliability,  
17 very well and it performed as it had been designed  
18 throughout the course of this clinical trial. The system-up  
19 time was 99.7 percent. There were 3 system faults that  
20 occurred in 3 separate cases for a total down-time of 13, 20  
21 and 12 minutes respectively. The system behaved as it had  
22 been designed to and went into an immediate "safe" state.  
23 The device was successfully rebooted in all instances and  
24 the procedures were completed without further incident.  
25 There were no adverse event outcomes as a result of these 3

1 faults and the software has been modified to eliminate them.

2 [Slide]

3 In conclusion, there were no device-related  
4 complications that occurred in this study. So, the safety  
5 endpoint has been met. All of the procedures that were  
6 begun with the Intuitive System were successfully complete  
7 with the Intuitive System, and the reduction in the  
8 DeMeester scores in the study group was equivalent to that  
9 in the control. So, the primary effectiveness endpoints  
10 have been met. Based on the clinical data, there is a  
11 reasonable assurance that the Intuitive System is both safe  
12 and effectiveness when used in accordance with its intended  
13 use.

14 [Slide]

15 Now, from a regulatory point of view, this device  
16 has already been cleared through the 510(k) pathway for  
17 blunt dissection, retraction, stabilization and manipulation  
18 and control of endoscopes. With this limited array of  
19 tools, the system is basically an assisting device.

20 [Slide]

21 What the company is seeking is to have this  
22 additional indication for use added to the labeling for this  
23 device that would allow tissue grasping, sharp dissection,  
24 suturing and use of the electrocautery.

25 [Video]



1           This additional indication and the tools it covers  
2 will enable surgeons to actually use this system clinically  
3 to safely and effectively perform laparoscopic surgery.

4           Thank you very much.

5           DR. WHALEN: Thank you, Dr. Gardiner. I would ask  
6 that each of the sponsor's presenters be near or at the main  
7 podium to answer any questions the panel may have, and while  
8 we may be asking those questions I would like to ask the  
9 others of the sponsor's group at the table if they would be  
10 kind enough and begin clearing from there in anticipation of  
11 FDA assuming that position very shortly.

12           Are there any of the panel members who wish to ask  
13 any questions or make any specific comments to the sponsor  
14 at this juncture? Dr. Galandiuk?

15           DR. GALANDIUK: I have two questions. Looking  
16 through the materials that were sent to us ahead of time, it  
17 listed in the device descriptions that it was for use by a  
18 "professional" and I was just concerned that this might be  
19 of use by nurses, and I think it should be intended for use  
20 by surgeons and I think that should be stated in the label.

21           The second question is for Dr. Gardiner.  
22 Considering that he, as an expert who used this device, took  
23 longer to do his operations with this, whether is the extent  
24 of the learning curve? And, if this new indication is  
25 approved, how will that be addressed by the company in terms

1 of teaching surgeons how to use this or giving them, you  
2 know, the 12 training cases that the investigators had here  
3 or some significant experience so that they know how to use  
4 this?

5 DR. GARDINER: I think the learning curve is  
6 always a difficult issue. Having been involved in the  
7 development of cholesection, for example,  
8 laparoscopically, it took me 65 or 70 cases before I was  
9 really comfortable that I could do those operations with  
10 equal facility to open surgery. Now, I don't think it is  
11 going to take anywhere near that degree of training on this  
12 device. It is pretty clear to me that you are safe and  
13 effectiveness sitting down and operating with it initially,  
14 but there is no question you are going to get better as you  
15 go -- I don't know, 10 or 15 cases maybe. But I think that  
16 the learning curve is going to continue much beyond that.  
17 In terms of training, what the company is going to do about  
18 that I don't know. We probably ought to have Dr. Moll  
19 handle that.

20 DR. MOLL: Just to remind you, there were 12  
21 training cases in this study distributed between 4 doctors.  
22 I think in one sense surgeons never have enough training  
23 but, clearly, training is a very important part of this  
24 story and will be a very important part of how this system  
25 is introduced. There is no surgical device that is

1 introduced and is immediately picked up by the surgeon and  
2 used properly without training.

3 I won't go into specific plans about how the  
4 system, if sold in the United States, will be trained. I am  
5 probably not the right person to do that, but it is at the  
6 top of our mind and we will have very clear plans for  
7 introducing a training protocol together with the sale of  
8 this device.

9 DR. GALANDIUK: How is the training conducted in  
10 Europe?

11 DR MOLL: In Europe it begins with bench-top  
12 training. In other words, you remember from the video that  
13 the system is used in the lab and the surgeons and the  
14 surgical assistants are brought in and have a very thorough,  
15 didactic overview of the system, its capabilities and sort  
16 of a general philosophy of how it works and how the system  
17 architecture is designed. Next, the surgeon and the  
18 surgical team have an opportunity to spend a lot of time on  
19 the bench-top suturing and performing dexterous tasks, and  
20 just using the system, moving their arms around and  
21 understanding how it is placed in the operating field, how  
22 it is draped, and all the things that go along with proper  
23 setup. Once the proper didactic and bench training is  
24 performed, we then sponsor a number of cadaver labs and  
25 animal labs where the system is used both in animal tissue

1 and in cadaver tissue to understand, again, proper setup and  
2 the limitations of the system, and how it is ideally set up  
3 and how it works most effectively.

4           So, it is bench-top training, didactic training,  
5 animal training, cadaver training extensively before it is  
6 introduced clinically. Obviously, it doesn't stop there.  
7 When we introduce it clinically we have both a proctor that  
8 is required to be on site to walk the new investigators and  
9 the new users of the system through the first cases. We  
10 also script very carefully what they do and what they do not  
11 do with the system, such that people walk before they run  
12 with the system. Most of the time they begin with open  
13 procedures. In other words, they use the system in an open  
14 case before they are asked to do minimally invasive  
15 technique. So, it is an extensive process of making sure  
16 that before there are procedures that challenge the  
17 surgeon's familiarity with the device he has an extensive  
18 period of getting comfortable and he understands the basics  
19 of the system.

20           Having said that, as Dr. Gardiner mentioned,  
21 surgery is all about experience, and there obviously is  
22 going to be a learning curve with this system as with any  
23 other surgical device. But we take very seriously the need  
24 for training as a part of the story.

25           DR. ANDERSON: I have just a few questions. Dr.

1 Moll, in terms of the prior experience with these tools in  
2 practice -- I can see that because of the articulation the  
3 system has a lot more nooks and crannies than your standard  
4 laparoscopic tools. Are there any problems with cleaning  
5 these devices in comparison to those other tools, and things  
6 getting caught inside those delicate articulations?

7 DR. MOLL: Yes, they are based on a cable and  
8 pulley system so there are areas of the wrist where that  
9 mechanism is exposed. Having said that, they have been  
10 designed to be cleaned and sterilized such that they have  
11 flush ports where fluid can be run through the inner parts  
12 of the system so they can be thoroughly cleaned and  
13 sterilized. Our experience in now 400 or more cases is that  
14 we have not seen any problem with cleaning or sterilization  
15 of these instruments.

16 DR. ANDERSON: Does that mean that there do need  
17 to be special instructions about how to clean these devices  
18 in comparison to other devices?

19 DR. MOLL: I think absolutely there need to be  
20 special instructions in the labeling as to how the  
21 instruments are best cleaned and sterilized.

22 DR. ANDERSON: I also have a question for Dr.  
23 Gardiner, if I could. I just want to make a comment. You  
24 talked about the protocol and when you randomize. I think  
25 you reported your data correctly by doing it from the time

1 of randomization because otherwise patient selection gets  
2 introduced that screws up the randomization principle. But  
3 I do have a question from a technical, practical standpoint.  
4 The two procedures that you are talking about are not  
5 procedures commonly associated with major bleeding but we  
6 can imagine that with time surgeons are going to try  
7 different techniques, and thinking of a laparoscopic  
8 splenectomy maybe you could get into a major vessel injury  
9 during this and you, as a surgeon, are sitting across the  
10 room from this. You are not scrubbed in at that point. Your  
11 hands are on the controls. What if you got into major  
12 bleeding and you needed to get into that abdomen quickly?  
13 How can you imagine that happening in a way that is safe and  
14 not life-threatening for the patient?

15 DR. GARDINER: Well, we are really not sitting  
16 across the room, we really are, certainly in vision, sitting  
17 next to the patient; not in the sterile field but next to  
18 the patient. You know, bleeding is a problem with any  
19 laparoscopic procedure and, certainly, that was a lot of our  
20 concern when we first started doing laparoscopies -- what  
21 happens if you have to convert? Well, you make an incision  
22 and convert. As I think you can see from the video, it  
23 takes just a second or two to get that system out of the way  
24 if you need to convert, and it is a matter of gowning and  
25 just stepping up to the table.

1 DR. ANDERSON: But the difference is you are not  
2 scrubbed in; you are not sterile, unlike a standard  
3 laparoscopic procedure. So how would that happen? Do you  
4 need to have somebody else standing by? Your assistant?

5 DR. GARDINER: The assistant is standing by. This  
6 is not intended to be used, you know, as a stand-alone  
7 instrument. So, you do have an assistant at the table.

8 DR. CRITTENDEN: The assistant needs to be a  
9 surgeon?

10 DR. GARDINER: I think that is a good question.

11 DR. CRITTENDEN: Given the scenario you are  
12 talking about, it sounds like it.

13 DR. GARDINER: Yes, you know, certainly whether  
14 the assistant is a surgeon or a PA, that certainly varies  
15 from community to community. Certainly, at this point in  
16 time we would envision that a surgeon would be an assistant.

17 DR. TALAMINI: In this case, wasn't the surgeon at  
18 the console scrubbed, with gloves and already sterile, or  
19 did I misread that in the protocol?

20 DR. GARDINER: That is true.

21 DR. TALAMINI: So, the surgeon would just need to  
22 take the gloves off or gown off and they would be sterile  
23 already, or not?

24 DR. GARDINER: Right. No, that is correct. You  
25 probably would want to strip the gown off or put a new gown

1 on.

2 DR. TALAMINI: I have two questions, and as a  
3 clinical investigator these are a little hard for me to ask  
4 but I think I need to ask them. First of all, perhaps for  
5 Dr. White or Dr. Gardiner, the two cases that you described  
6 in detail where it was elected, once the laparoscope was  
7 placed, to not use the system, the fact that those two cases  
8 were withdrawn early in the study, did that skew the  
9 results? In other words, were there equally tough cases  
10 later on that were completed using the device? That is one  
11 question.

12 The second question is it seems to me, with a  
13 general knowledge of the literature, that the gallbladder  
14 rupture rates and two trocar injuries in 200 patients are on  
15 the high side for those incidents. I guess I feel like I  
16 need you to help me with why that doesn't contaminate the  
17 rest of the study. I understand that many of those things  
18 occurred when the system was not in place, yet they are  
19 there and they are part of the results, and I need you to  
20 help me with why that shouldn't contaminate our  
21 interpretation of the remainder of the results when the  
22 system was in place.

23 DR. GARDINER: Well, none of those complications  
24 occurred while the system was in place. So, these are  
25 complications that occurred during the traditional standard



1 conventional laparoscopic part of the operation. So, it has  
2 really nothing to do with the system. Why they happened? I  
3 don't know. I mean, I can tell you that these surgeons were  
4 very experienced. Actually, two of the injuries occurred to  
5 Dr. Cadriere. I have watched him operate. He has an  
6 extraordinarily fine reputation, and I don't know why they  
7 happened to him in this study. I can't answer that. But I  
8 don't think that they had any relationship at all to the use  
9 of the system so it is really two kind of separate  
10 questions, it seems to me.

11           With regard to all of these complications, I think  
12 it goes to the serosal injuries to the stomach, as well as  
13 gallbladder perforation. We have been very rigorous in what  
14 we looked at and if there was a drop of bile that leaked out  
15 during the case, it was counted as a rupture. So, I think  
16 that we have perhaps, if anything, skewed the data to a more  
17 critical side and I think that is probably the explanation  
18 for the gallbladder issue. But these injuries that occurred  
19 in the traditional part of the operation -- I struggled with  
20 why they happened and I don't know.

21           I mean, I can tell you that with that case that I  
22 had with the abdominal wall port site bleeding -- I can't  
23 remember when that has happened to me but, here we are,  
24 dealing with it. So, I don't know.

25           DR. FERGUSON: Could I ask a question of Dr.

1 Gardiner? It gets back to the gowning and gloving aspect.  
2 The four surgeons who did these cases, did they gown and  
3 glove as we read in the protocol, or did they not do that?

4 DR. GARDINER: No, they are gowned and gloved but,  
5 you know, you are not sterile obviously --

6 DR. FERGUSON: Not scrubbed?

7 DR. GARDINER: Right.

8 DR. FERGUSON: So, it is not a matter of putting  
9 on a gown and gloves, like you might do in some cases,  
10 scrubbing first and then gowning and gloving and being able  
11 to change quickly if you needed to.

12 DR. GARDINER: That is correct. I mean, you could  
13 certainly do that.

14 DR. FERGUSON: The reason I am pursuing this is  
15 that I think the issue is important as the device is  
16 disseminated because this is a true difference from a  
17 surgical situation, and I think it needs to be emphasized.

18 DR. GARDINER: We would certainly agree with that.  
19 We did not address the issue of those two cases of Dr.  
20 White's.

21 DR. WHITE: Mr. Chairman, panel members, Alan  
22 White, Tacoma, Washington. I am a general surgeon by  
23 training. I practice in a multi-surgical group, six of us,  
24 and our chief emphasis is in laparoscopic breast and  
25 colorectal surgery. In addition, I am the medical director

1 of the Multi-Care Health System Endosurgical Institute, in  
2 Tacoma, Washington, where we give advanced laparoscopic  
3 training and proctoring and precepting to surgeons in  
4 private practice.

5 Specifically, would you repeat the question in  
6 regard to the two patients that were converted in my series?

7 DR. TALAMINI: My only concern, sir, was whether  
8 the removal of those two cases early in the study had the  
9 potential to skew the data. The other side of the coin  
10 would be were there equally difficult cases later on in the  
11 study that were accomplished with this device that would  
12 eliminate that possibility of the data being skewed by early  
13 removal of those cases?

14 DR. WHITE: First off, I am not a statistician so  
15 I can't address the skewing of the data. I would defer that  
16 to our statistician. Specifically though, in the conduct of  
17 the procedure and what we saw and did as this, my first and  
18 third randomized events, let me take you there.

19 Specifically, I had done the day before a training  
20 procedure using the ISI system for laparoscopic  
21 cholecystectomy and a laparoscopic Nissen. Both were  
22 standard cases with no obvious or significant and severe  
23 pathology. On my very first case on the first randomized  
24 day, my second day in Mexico City, I put the scope into the  
25 first case and encountered a cholecystitis that, having done

1 probably 2000 laparoscopic cholecystectomies, I would count  
2 in the top five as far as severity, a case that had both  
3 severe acute and chronic inflammation.

4           At that time, I held the procedure in abeyance. I  
5 called the on-site coordinator and we discussed it at that  
6 time and, based on my understanding of the protocol at that  
7 time, I was under the impression that we could make an  
8 intraoperative exclusion. That is, we found something that  
9 should have been found ahead of time and we were not aware  
10 of. Discussing it at that time, I opted to say clinically  
11 this is a case that I don't believe I can do  
12 laparoscopically. Therefore, my option is made to convert.

13           Conversion is a term that we have had to struggle  
14 with in laparoscopy since its inception, that is, conversion  
15 to an open event is not bad. At that point in time, with my  
16 level of training with the system being rudimentary or at  
17 least not sophisticated, I opted not to use the system on  
18 that case but, rather, proceed to an open event, and I opted  
19 to use standard laparoscopic technique, and went ahead and  
20 proceeded with the case in that way.

21           The third case was much akin to the one Dr.  
22 Gardiner reported. Again, that same day, a case where, on  
23 placement of the scope, both macronodular cirrhosis and  
24 portal hypertension were found and, again, with that same  
25 reasoning the same series of events occurred.

1 Later in the trial I probably did  
2 cholecystectomies and Nissen funduplications as severe as  
3 anything that was seen, obviously, on that day. But on that  
4 day, clinically I was not willing to proceed with the  
5 device.

6 DR. TALAMINI: Thanks. Obviously, that was the  
7 wise thing to do.

8 DR. WHITE: I don't know about wise but it was  
9 what happened and it is reported.

10 MS. DUBLER: I have two questions and I think they  
11 are probably for Dr. Moll. The first has to do with the  
12 structure of the study itself. You said that in the three  
13 months of the study there was no site outside of the one you  
14 chose that had sufficient numbers of cases. Was that  
15 correct?

16 DR. MOLL: That is correct.

17 MS. DUBLER: Why did you limit it to three months?  
18 What was so important about that that you could not collect  
19 other sites and see if there were any possible differences  
20 among the sites in how the study proceeded?

21 DR. MOLL: There were financial considerations and  
22 there were considerations about the availability of the  
23 system. We had one system that was available and able to be  
24 used for this trial. So, rather than accomplish a  
25 multicenter study, we believed that in the interest of

1 efficiency we could have surgeons come to one site and use  
2 the system rather than -- which would have been impossible  
3 at the time -- move the system from site to site, because we  
4 only had one system. So that was sort of the overriding  
5 feature.

6 I don't have the data but we very carefully looked  
7 at the number of procedures performed in this country,  
8 laparoscopic Nissen funduplications, and although I am not  
9 sure of the total number it is spread out among a number of  
10 centers, and to accomplish this number of procedures under  
11 one roof anywhere else seemed virtually impossible in a  
12 reasonable period of time that we could financially afford  
13 to do. I think that was the major point.

14 MS. DUBLER: I have a second question, and that  
15 has to do with the informed consent process that you used  
16 and the result of that or the reflection of that in the  
17 informed consent form that was presented to your patients.  
18 In the form that I reviewed just today there are no benefits  
19 suggested for this procedure.

20 DR. MOLL: Correct.

21 MS. DUBLER: And there has been a discussion today  
22 of equivalence but not of benefit and, therefore, my  
23 question to you is what is the benefit of this procedure?

24 DR. MOLL: Yes, I tried to address a little bit of  
25 that up front, but it is a very good question and a question

1 that clearly comes out of the study the way it was designed.  
2 This was an equivalence study. We had no intention of  
3 trying to prove that this surgery was somehow better than  
4 conventional technique. So, I think it gets back to the  
5 capabilities of this system in surgery, in minimally  
6 invasive surgery, that although the benefits in a routine  
7 gallbladder surgery may be very hard to describe and, in  
8 fact, there probably aren't any, I think you would agree  
9 that if a system is safe and efficacious from the standpoint  
10 of being able to do the procedure and deliver more dexterity  
11 intraoperatively than the surgeon is accustomed to with  
12 conventional technique, then we believe that efficacy and  
13 the benefit will come from this technique.

14 Now, there is anecdotal evidence that has to  
15 remain anecdotal at this time about what the system can  
16 offer. I think Dr. White can describe for you a situation  
17 in one of his cases where he needed, because of the dilation  
18 of the cystic duct, to suture ligate the duct. I don't want  
19 to speak for him but I think is very difficult, if not  
20 impossible to do with conventional laparoscopic technique.  
21 It is those sorts of situations where articulation and  
22 control can really add true benefit to a routine  
23 laparoscopic procedure, and I think as procedures get more  
24 complex the benefits of the system become more apparent.

25 MS. DUBLER: Let me just have one brief follow-up

1 question. Was this protocol reviewed by an IRB in the  
2 States, even though it was instituted in Mexico?

3 MR. DANIEL: Yes, a very similar -- I can't say  
4 identical; I would have to look it up, but a very similar  
5 protocol was, in fact, reviewed by Summit Medical's IRB.

6 MS. DUBLER: By whose IRB?

7 MR. DANIEL: Summit Medical, the hospital in  
8 Oakland, California where Dr. Gardiner operates. We had  
9 anticipated a possibility of proceeding at that location  
10 early on, before we got up to the numbers that we ended up  
11 getting. My memory is that there is virtually no difference  
12 there but, of course, we did have an ethics committee in  
13 Mexico.

14 MS. DUBLER: I just want to point out what is  
15 puzzling to me, given the discussion this morning, given the  
16 learning curve, given the teaching cases, etc., that there  
17 is no mention in that discussion of the fact that this is a  
18 procedure that surgeons are learning to use, which brings  
19 with it its own level of risk; nor are there any potential  
20 benefits suggested, which puts an IRB in a peculiar  
21 situation of weighing articulated risks against no possible  
22 benefit that is articulated, in which case it is  
23 questionable whether they should permit it to go forward.  
24 So, I think that some attention, if this is approved, toward  
25 an informed consent process that reflects experience and



1 possible risk and benefit in a more comprehensive way would  
2 probably be advisable.

3 MR. DANIEL: You may help me because you have it  
4 in front of you and I do not, but I remember wording that  
5 indicated clearly that the patient's surgeon had a great  
6 deal more experience with the conventional tools than he or  
7 she would have with the Intuitive System. So, we tried to  
8 impart the idea that the surgeon did not have anywhere near  
9 close to the experience with the Intuitive System as we the  
10 conventional.

11 MS. DUBLER: Just one very brief question on the  
12 experience and how people will be trained. It would seem to  
13 me that that experience would need to be quantified in some  
14 way, and I wondered whether you have considered a  
15 certification process for the use of this system.

16 DR. MOLL: Yes, we have not closed on exactly what  
17 sort of formal training will be required for the system, and  
18 I think that is something that deserves a lot of discussion  
19 but it is not something that we have addressed to date.

20 I want to go back and just make sure that there is  
21 no misunderstanding on the sterile field issue. The  
22 surgeons that performed this trial and the surgeons in our  
23 experience in Europe is that the surgeon is gowned and  
24 gloved as he sits at the console. So, if there is a need --  
25 and in my recollection I don't remember this occurring in

1 either the clinical study or the European cases, but if  
2 there is a need to step to the sterile field, the surgeon is  
3 scrubbed and he needs to merely don gloves and a gown to  
4 step to the operative field. So, that does take some time  
5 but it can be very rapid if, in fact, there is a need to do  
6 that.

7 DR. CHANG: Just to clarify that, it was my  
8 impression, and I would like Dr. Gardiner or Dr. White to  
9 clarify that but it was my impression that the clinicians  
10 were scrubbed but I would like Dr. Gardiner or Dr. White to  
11 say, yes, they were scrubbed during this clinical trial, or  
12 perhaps this should be a strong recommendation as a safety  
13 feature that a surgeon be available should intra-abdominal  
14 bleeding occur.

15 DR. GARDINER: I think I would agree with you.  
16 No, we were scrubbed because the surgeons put the trocars  
17 in. So, what would happen is we would scrub and gown just  
18 as we would in any traditional laparoscopic procedure, step  
19 to the table, place the trocars, then step away and put  
20 gloves on so that we didn't contaminate the instrument, and  
21 then do the procedure.

22 DR. MOLL: Just as a sort of backup to that, this  
23 was obviously an issue that we struggled with in the design  
24 phase of this system and it obviously would be possible to  
25 drape and design a console so that it is sterile so that a

1 surgeon can sit down and remain in a sterile field. We  
2 chose not to do that for some very good reasons, that if the  
3 surgeon is away from the sterile field it greatly increases  
4 the chances of contamination. So, once the surgeon, in the  
5 system that we have currently, sits down it is clear that he  
6 is not sterile and if he chooses to go back to the operative  
7 field he needs to become sterile. So, that was a very clear  
8 decision on our part and I think it is the right decision.

9 DR. CHANG: Dr. Moll, my actual question is that  
10 in the protocol conventional clip was used for the cystic  
11 artery and duct. So, we truly don't have a sense of the  
12 effectiveness if the ligatures by using the Intuitive  
13 System. So, in fact, it is belt and suspenders. Are there  
14 plans for adding the clip as an instrument? Would this be  
15 an addendum to your PMA, or would surgeons rely on the  
16 traditional laparoscopic instruments to add that clip in  
17 addition to the suture ligature?

18 DR. MOLL: We will at some point add the  
19 capability to deliver large clips to the system. We  
20 actually have a system that now can deliver small clips but  
21 not large clips. It was sort of a story as to two competing  
22 philosophies on what is the best method of ligation of the  
23 cystic duct and artery that led to this protocol and,  
24 arguably, you know, if we had to do it again we might have  
25 done it differently to be more clear about what was doing

1 the occluding, but maybe I can ask Barry to comment on that.

2 DR. GARDINER: You are absolutely correct about  
3 not making any conclusion about the integrity of the  
4 occlusion of the cystic duct or artery because that was done  
5 by a clip. In the investigators meeting in sorting this  
6 out, three of them used traditional clips and that is the  
7 way they did this operation. I happened to tie, that is the  
8 way I did the operation and I was agitating for the  
9 ligatures and the rest of the investigators felt we ought to  
10 do this the way we do traditional surgery, with clips, and  
11 that is the way the lap chole is done and that is the way it  
12 ought to be done. So, that is the process that we went  
13 through to get there. As Dr. Moll said, I think in  
14 retrospect we probably would have preferred to do it with  
15 just ligatures.

16 With regard to the integrity of the suturing, of  
17 the knot tying, we do have documentation of that with the  
18 DeMeester score. If you look at the DeMeester data in the  
19 control and the study group, there is no difference between  
20 the DeMeester score reduction, and if the knots that this  
21 system tied still weren't intact those DeMeester scores  
22 would have been different. So, I think we have established  
23 the suturing capability by the DeMeester score. That was  
24 the reason we chose that particular endpoint.

25 DR. WHALEN: I am going to ask that other

1 questions be held in abeyance because we need to get the  
2 lunch time. We will have more time for those questions. We  
3 will reconvene in this room at 1:30 for a closed session, at  
4 which time Intuitive Surgical will have the opportunity to  
5 present or discuss any trade secret data that they may or  
6 may not have. At two o'clock the public session will reopen  
7 and at that time FDA presenters will bring to the microphone  
8 the information that we need.

9 [Whereupon, at 12:35 p.m., the proceedings were  
10 recessed, to be resumed in closed session at 1:30  
11 p.m., followed by an open session at 2:00 p.m.]

## 1 AFTERNOON SESSION

## 2 FDA Presentations

3 MR. YEN: Good afternoon, Dr. Whalen and members  
4 of the panel. FDA will now present a summary of the reviews  
5 of this submission.

6 [Slide]

7 For your convenience, a copy of all of the  
8 presenters' slides are at your disposal. I am the lead  
9 reviewer, Dwight Yen, and will give a brief summary of the  
10 regulatory and engineering review. Dr. Horbowyj will  
11 present the clinical review and Dr. Bushar will present the  
12 statistical analysis.

13 [Slide]

14 As the sponsor has indicated, the initial system,  
15 the initial control system was cleared in July of 1997 for a  
16 limited set of instruments, including blunt dissectors and  
17 retractors, based on substantial equivalence to other  
18 instruments and laparoscope holders already on the market.

19 At the beginning of this year Intuitive Surgical  
20 submitted a new 510(k) for the same control system but  
21 adding new instruments, including forceps, scissors,  
22 scalpels, clip appliers, needle holders and electrocautery,  
23 so that the system can also be used to perform grasping,  
24 cutting, electrocautery and suturing.

25 FDA found the system is no longer substantial

1 equivalent, or NSE, to an instrument holder because the  
2 change in intended use afforded with the new instruments  
3 raised different types of questions of safety and  
4 effectiveness. Therefore, this device is classified by  
5 statute into Class III and is considered a PMA.

6           The sponsor has already presented a description of  
7 their device and principles of operation. FDA considers the  
8 device description information factual and adequate.  
9 Therefore, very briefly, the system consists of three  
10 components -- the surgeon's console, the patient table side  
11 surgical cart and the system electronics.

12           [Slide]

13           I have reviewed the hardware design, and to the  
14 credit of the designer, they have made safety one of the  
15 highest priorities, starting with the design, and have  
16 implemented numerous safety-related functions in the  
17 hardware.

18           Our software group at the Office of Science  
19 Technology provided consultation and review of the software  
20 and has determined that the software design and development  
21 meets FDA guidelines for software in medical devices.

22           I will elaborate on the next slide more on this  
23 performance testing. Like other surgical tools, material  
24 use and the manufacturing of the instruments are considered  
25 limited duration patient contact materials. The sponsor has

1 met FDA guidelines for biocompatibility of each patient  
2 contact material that was used.

3 The instruments are single use and reusable  
4 devices that are required to be sterile. Information for  
5 sterilization and instruction to clean and resterilize the  
6 instruments is adequate and has been validated.

7 Finally, the system meets the FDA guidelines for  
8 electrical medical safety standards and electromagnetic  
9 compatibility standards.

10 In the area of performance testing, two hysteresis  
11 studies were performed to evaluate each degree of freedom of  
12 the control arm for reproducibility and precision. Six  
13 animal studies were conducted to evaluate system setup, port  
14 placement, vision and system performance. Four cadaver  
15 studies were conducted to evaluate port placement, access  
16 issues, tissue manipulation and suturing. Surgeons were  
17 invited to use the system to perform a variety of surgical  
18 procedures using a porcine heart. Feedback was received in  
19 terms of the hand grip design and vision to the system.

20 A prototype system was used in Belgium, in 1997,  
21 for clinical a feasibility study on five patients to  
22 initially test the system and, of course, the clinical trial  
23 was conducted in Mexico City last year.

24 [Slide]

25 The clinical trial was completed over a three-



1 month period, between July and October of 1998. As the  
2 sponsor has earlier indicated, three system faults were  
3 experienced during surgery resulting in a 12-13 delay in two  
4 cases and a 20-minute delay in the third case.

5 The sponsor described that in all three cases the  
6 system responded in a fail-safe manner. The instruments  
7 stopped responding to surgeon input. The definition of  
8 fail-safe here is that the system entered the appropriate  
9 error-landing state without uncontrolled motion. This  
10 ensured patient safety and allowed the instrument to be  
11 removed safely. While the exact cause of the failures has  
12 not been determined, the sponsor was able to pinpoint the  
13 failures to an electronic board. A modification has been  
14 made to reduce the likelihood of future faults of this kind.  
15 In each case the system had to be restarted to complete the  
16 procedure. A second modification has been made to add a  
17 reset button to the interface panel that will reduce the  
18 time to restart the system and continue the procedure.

19 This concludes my presentation. I will return to  
20 review the panel questions after Dr. Horbowyj and Dr. Bushar  
21 have completed their presentations. I would now like to  
22 introduce Dr. Horbowyj.

23 **Clinical Aspects**

24 DR. HORBOWYJ: Good afternoon. My name is Roxi  
25 Horbowyj. I am a general critical care surgeon and the

1 clinical reviewer for this application.

2 [Slide]

3 So, I will be presenting the FDA clinical  
4 perspective on the Intuitive System's laparoscopic clinical  
5 study. A lot of the aspects have already been gone over in  
6 good detail by the sponsor and so I will really just go over  
7 highlights, including the objective, the design, procedures,  
8 endpoints, sample size determination, the target population,  
9 outcomes and end with a brief summary.

10 [Slide]

11 As you know, the Intuitive Surgical System allows  
12 surgical tasks to be performed with software-assisted three  
13 ports. Two are for the surgeon's hands and one is for the  
14 laparoscope. In conjunction also conventional laparoscopic  
15 instruments may be used for any additional ports or for any  
16 instruments that are not adapted to the system.

17 [Slide]

18 The objective of this study, as you have heard,  
19 was to demonstrate that the Intuitive surgical endoscopic  
20 instruments would be equivalent in safety and to standard  
21 laparoscopic equipment, which was in the control, when used  
22 to perform general laparoscopic tasks such as grasping,  
23 cutting, blunt and sharp dissection, approximation,  
24 ligation, electrocautery and suturing.

25 [Slide]

1           Taking these tasks into account, the design of the  
2 study was intended to provide valid scientific data that  
3 would allow reasonable clinical assessment of device safety  
4 and effectiveness, independent of the regulatory path to  
5 market.

6           So, the study was prospective, concurrently  
7 controlled, carried out by multiple investigators with  
8 single mask, and randomized, in this case preoperatively,  
9 after inclusion and exclusion criteria were met and informed  
10 consents were signed. Follow-up was for 30 days.

11           [Slide]

12           The procedures that were chosen had these surgical  
13 tasks in mind. Laparoscopic cholecystectomy was chosen  
14 because it is a well-established, widely practiced  
15 procedure, usually straightforward and it is excisional.  
16 So, any trauma caused to tissue would most likely not be  
17 seen in the patient as the tissue that is manipulated --  
18 most in this procedure is removed from the patient.

19           Laparoscopic Nissen fundoplication was undertaken  
20 because it is technically more challenging, and it is a  
21 reconstructive procedure, therefore, any effects on tissue  
22 would most likely be retained within the patient and could  
23 possibly present sequelae.

24           [Slide]

25           Sample size was determined upon consideration of

1 multiple variables. We looked at the literature-reported  
2 complication rates with these procedures; the literature-  
3 reported cohort sizes; the sample size that might be needed  
4 for learning curve assessment, the sample size needed for  
5 clinically reasonable assessment of safety and  
6 effectiveness; as well as sample size that would be  
7 determined by statistical calculations.

8 [Slide]

9 The endpoints that were chosen were as much as  
10 possible endpoints that have been previously validated and  
11 were objective, and represented safety and effectiveness.

12 [Slide]

13 The endpoints, therefore, chosen by the sponsor  
14 were conversion rate were defined in the protocol as  
15 conversion of ISI to conventional instruments or conversion  
16 of control to open technique. It was recognized, and would  
17 be recognized, that patient anatomy and pathology, software  
18 or hardware failure, and surgeon or surgical team position  
19 on the learning curve could contribute to conversion.

20 Procedure duration was defined as being from skin  
21 incision to skin closure; postoperative hospital stay in  
22 days; DeMeester score at 30 days, specifically for  
23 laparoscopic Nissen fundoplication.

24 Quality of life was evaluated using the  
25 psychological well being score at 30 days, as well as

1 preoperatively. This score particularly has been used in  
2 gastrointestinal procedures in the past. It has been used  
3 evaluating Nissen fundoplication done both in the open  
4 technique and with laparoscopic technique.

5 Other measures of safety that we usually look at,  
6 such estimated blood loss, were considered, and then  
7 specific to these procedures would be bioleak and dysphagia.

8 [Slide]

9 The target population, as you heard described,  
10 were otherwise healthy adult patients with gallbladder  
11 disease or gastroesophageal reflux disease confirmed by this  
12 protocol, who were expected to benefit from non-emergent  
13 laparoscopic cholecystectomy or laparoscopic Nissen  
14 fundoplication, and who were willing to participate in a  
15 clinical study.

16 [Slide]

17 The outcomes can be looked at preoperatively,  
18 intraoperatively and postoperatively. Preoperative data  
19 shows, as the sponsor has shown, that the control and  
20 investigational device study populations were clinically  
21 comparable for demographics and inclusion and exclusion  
22 criteria.

23 [Slide]

24 Intraoperative data, such as review of video  
25 tapes, shows that laparoscopic cholecystectomy demonstrated

1 grasping, blunt dissection, cautery dissection, as well as  
2 suture placement around cystic duct and arteries.

3 Evaluation of knot integrity, as was mentioned  
4 earlier by one of the panel members, was precluded by  
5 conventional clip placement on the patient side of the  
6 cystic artery and duct.

7 As you have heard, two investigational device  
8 randomized patients were converted to and completed with  
9 control due to patient pathology and possibly because of  
10 surgical team position on the learning curve.

11 [Slide]

12 Similarly, for laparoscopic Nissen fundoplication  
13 video review demonstrated grasping, blunt dissection,  
14 cautery dissection, needle suture placement and suture tie  
15 for tissue approximation. As you have heard, one  
16 investigational device randomized patient was converted to  
17 and completed with open technique.

18 [Slide]

19 Looking at procedure duration and estimated blood  
20 loss has brought up some discussion. Looking at the means,  
21 for both procedures differences are seen, as well as in the  
22 standard deviation and the range. Similarly, there are some  
23 differences that are seen in estimated blood loss.

24 When looking at the data per investigator and  
25 comparing investigator to investigator, generally what I

1 have found was that there was variability per investigator  
2 as well as from investigator to investigator. There seemed  
3 to be more variability for the investigational device  
4 compared to control, and there seemed to be approximation of  
5 investigational device procedure duration to control with  
6 time or number of cases. But this pattern was variable and  
7 the degree of approximation to control varied from  
8 investigator to investigator.

9 [Slide]

10 This is an example which I hope will help to  
11 illustrate this. This is a single investigator who  
12 completed 57 cases. This includes his training cases. This  
13 slide shows both laparoscopic cholecystectomy as well as  
14 Nissen fundoplication for both device types used by this  
15 investigator. This investigator shows some of the points  
16 that I wanted to demonstrate but it is not necessarily the  
17 person who had the largest number of cases.

18 In the blue and the pink are the investigational  
19 device procedures. The triangles represent Nissen  
20 fundoplication and the circles being laparoscopic  
21 cholecystectomies. The yellow and the orange are control  
22 devices.

23 I think what this demonstrates most is the trends.  
24 There are variabilities in all but the controls for  
25 procedure duration are low and some are more constant

1 throughout the study. These cases are plotted in time as  
2 completed by the investigator. So this would have been his  
3 training case on laparoscopic cholecystectomy and this would  
4 have been his training case on Nissen fundoplication. These  
5 would have been his final two cases performed during the  
6 study. What you also see is that, for example with Nissen  
7 fundoplication performed with the investigational device, as  
8 time goes by, as he does more cases, the procedure duration  
9 time approximates control. You don't quite see that as much  
10 with laparoscopic cholecystectomy performed with the device,  
11 perhaps because there just were fewer cases performed by  
12 this investigator, and so we didn't see that approximation.

13 [Slide]

14 As you have heard, there were several unexpected  
15 events, three episodes of intraoperative software and  
16 hardware shutdown into safe mode during the middle of the  
17 clinical study. These did extend operative time 12-20  
18 minutes, with no known patient sequelae. They were  
19 recovered, however, with active engineering intervention,  
20 and they required system and system use modification.

21 [Slide]

22 Further, outcomes in the postoperative time show  
23 that for control and investigational device study  
24 populations adverse events rates, quality of life at 30  
25 days, DeMeester scores at 30 days and postoperative length



1 of hospital stay were clinically comparable.

2 [Slide]

3 In summary, the ability to perform surgical tasks  
4 with the investigational device and the control device in  
5 laparoscopic cholecystectomy and laparoscopic Nissen  
6 fundoplication has been demonstrated in the study  
7 population, specifically grasping, blunt dissection, cautery  
8 dissection, suture tie placement around tubular structures,  
9 needle suture placement and suture tie for tissue  
10 approximation. Notice, sharp dissection is not listed here.

11 [Slide]

12 Unexpected system shutdown into safe mode occurred  
13 requiring active engineering intervention and system as well  
14 as system use modification.

15 There was some increase in procedure duration and  
16 variability in estimated blood loss compared to control, and  
17 there were non-device failures associated with conversion of  
18 two investigational device randomized laparoscopic  
19 cholecystectomies to control device which were completed  
20 with control device, and this may be attributable to surgeon  
21 or surgical team position on the learning curve for the  
22 device use. There were no conversions due to device  
23 software or hardware failure during the study.

24 Thank you. Dr. Bushar will now present the  
25 statistical aspects of this study.

## Statistical Aspects

1  
2 DR. BUSHAR: Thank you, Dr. Horbowyj. My name is  
3 Harry Bushar. I did the statistical review of this  
4 submission. What I did is I actually entered the sponsor's  
5 data and did all of the statistical analyses that I thought  
6 were necessary. I also checked the sponsor's work and I  
7 deferred to the sponsor's results over mine. In no way do I  
8 differ with the sponsor on the technical analyses. The  
9 differences lie in the way I approached the data and the way  
10 I looked at the protocol. The main thing is that I took an  
11 intent-to-treat approach, which means I looked at all the  
12 patients. I didn't exclude anybody. They were randomized;  
13 they were in the study. Also, I didn't make any decision  
14 about what was or what was not to be included. If an event  
15 occurred and was reported, I used it in the analysis.

16 [Slide]

17 I did look at the clinical equivalence studies  
18 that the sponsor completed. These were the laparoscopic  
19 cholecystectomy, and I will refer to that as LC from now on,  
20 and the laparoscopic Nissen fundoplication, or LNF study.

21 The clinical trial design, as Dr. Horbowyj  
22 explained, was preoperatively randomized, patient masked,  
23 concurrently controlled by conventional laparoscopic  
24 instruments, and it was multi-investigator, that is, four  
25 surgical teams within one site in a hospital in Mexico City.

1 [Slide]

2 What I will be talking about here is mainly the  
3 statistical equivalence testing, which was promised in the  
4 protocol, and the objective of statistical equivalence  
5 testing is to show that the Intuitive group in this case is  
6 not worse than, that is better than or equal to the control  
7 group within delta, where delta is some positive value. The  
8 reason for that is that everything I am going to be  
9 measuring -- larger is worse. So, I am trying to show that  
10 Intuitive is not larger or not worse than control. To quote  
11 the sponsor, delta is equal to some predefined and clinically  
12 meaningful difference above which the two different  
13 methodologies are no longer considered to be substantially  
14 equivalent. This is sort of the essence of what I am  
15 talking about. I get into a lot of details but the main  
16 thing that the panel has to focus on is what is clinically  
17 relevant.

18 [Slide]

19 To put this in terms of equations, the statistical  
20 equivalence test hypotheses are -- well, the null hypothesis  
21 is what you are trying to reject because this means not  
22 equivalent. It means that the mean of the outcome for  
23 Intuitive is greater than or equal to the mean of the  
24 outcome for the control plus some delta.

25 The alternative hypothesis, which you would

1 achieve by rejecting the null hypothesis, is one of  
2 equivalence and that is just the opposite of what the first  
3 equation says, that the mean outcome of the Intuitive is now  
4 less than the mean outcome of the control plus some delta.

5 This is not just a matter of subtracting the  
6 means, which these hypotheses might seem to indicate. The  
7 means are from a sample and we are trying to project to the  
8 population. The delta applies to the population. That is  
9 why it is not a simple difference. You have to take the  
10 mean into consideration; you have to take the standard  
11 deviation and the sample size into consideration to actually  
12 perform this statistical test.

13 [Slide]

14 To begin, I will show you the clinical endpoints  
15 or deltas from the LC protocol table on page 387. What I  
16 have done, I have gone down the left-hand side of that table  
17 which says what the effectiveness outcomes and safety  
18 outcomes are, and down the right-hand column which says what  
19 the delta should be. This is not what the sponsor said was  
20 intended. There are different values in the text; there are  
21 different values in the footnotes, but I am going strictly  
22 with this table.

23 To begin with, with the primary effectiveness  
24 outcome the conversion rate had a delta of 3.5 percent. The  
25 secondary effectiveness outcomes consisted of two, the

1 procedure duration with a delta of 45 minutes, and hospital  
2 stay with a delta of 0.4 days. The safety outcome referred  
3 to overall complication rate. I am not leaving any  
4 complications out. I am not saying this is the best way to  
5 look at the data. This is one way to look at the data.  
6 This is what was stated in the protocol and that is what I a  
7 going to be covering. Here the delta is 3.5 percent, the  
8 same as it is for conversion rate.

9 [Slide]

10 If I apply that to the LC clinical findings and do  
11 the statistical test for equivalence, for the conversion  
12 rate there were two conversions in the Intuitive. So that  
13 gives me 3.5 percent mean for the Intuitive and a zero  
14 percent mean for the control. The predetermined delta is  
15 3.5 percent. You don't even need a statistical test there.  
16 This is not equivalent. There is no way you can show that  
17 it is less than 3.5 percent.

18 However, the minimum delta that would just allow  
19 equivalence to be established, in other words, if you used  
20 the equations that the sponsor provided, you can back-  
21 calculate a delta of 7.5 percent, which means that if that  
22 were the delta you used, which was not, then you could  
23 conclude from this data, from this sample, that the  
24 difference was less than 7.5 percent.

25 This is a secondary analysis and something like

1 this really generates a new hypothesis which should be  
2 tested with the new studies. So, this is just a way of  
3 looking at how we interpret the data when it doesn't show  
4 what it was intended to show.

5 [Slide]

6 Moving on, for procedure duration we got 109  
7 minutes in the Intuitive group and 67 minutes in the control  
8 group, with a predetermined delta of 45 minutes. This leads  
9 to non-equivalence. We can't reject the null hypothesis.  
10 However, the minimum delta that would just allow equivalence  
11 to be established is not much greater than that. You don't  
12 have to double it. You just add 6 minutes; you are up to 51  
13 minutes. You can actually show from this data, if that were  
14 the delta, that the Intuitive-control difference would be  
15 expected to be less than 51 minutes in the population.

16 [Slide]

17 Continuing, the LC procedure duration learning  
18 curve, which has been mentioned and is presented in an  
19 analysis by the sponsor, Dr. Bloch did this using general  
20 estimating equations and he showed that the linear slope of  
21 the procedure duration curves is 0.69 minutes per procedure  
22 for Intuitive and 0.83 minutes per procedure for the  
23 control.

24 Now, the above reduction rates are not statistical  
25 different from zero so that the conjecture that learning

1 here would be expected to eventually catch up, in other  
2 words, Intuitive would eventually catch up with the control,  
3 does not appear to be plausible from this data. But it is  
4 going down. But I think the point here is that surgeons are  
5 getting better as they do more and more Intuitive  
6 operations; they are also getting better and better as they  
7 do more control operations within this clinical study.

8 [Slide]

9 The LC clinical findings with statistical test for  
10 equivalence for hospital stay -- here the hospital stay for  
11 Intuitive was 1.3 days and the hospital stay for control is  
12 1.2 days. Even a predetermined delta as low as 0.4 leads to  
13 equivalence. There is equivalence there.

14 [Slide]

15 Going on now to safety, looking at the  
16 complication rate -- again, I am looking at all  
17 complications -- there is 3.6 percent in the Intuitive group  
18 and zero percent in the control group. The predetermined  
19 delta is 3.5 percent. We see that as not equivalent.  
20 However, if we boost that delta up to 7.7 percent we could  
21 just get equivalence.

22 [Slide]

23 Here is another variable that I looked at and the  
24 sponsor looked at also. This is the LC blood loss. This  
25 was not mentioned in the protocol but it was collected by

1 the sponsor and they did analyze it. The means were 19.2 ml  
2 for the Intuitive and 5.6 ml for the control, and this is  
3 statistical significantly different, with a very low p  
4 value. So, there is more blood loss on average in the  
5 Intuitive than in the control but whether or not 20 ml means  
6 anything is, of course, a matter of clinical interpretation.

7 [Slide]

8 Continuing now with the next study, the LNF, again  
9 I am beginning with the clinical endpoints, deltas, from the  
10 LNF protocol table on page 388. The primary effectiveness  
11 outcome is stated to be conversion rate, but then the  
12 sponsor goes on to discuss the DeMeester score, and I think  
13 the delta given there, 2.22, was meant to apply to the  
14 DeMeester score but I just used the same delta as for the  
15 previous study, 3.5 percent. I am not sure what was meant  
16 by what was in that table. Anyway, I am doing the same  
17 thing for this study that I did for the previous study when  
18 it comes to conversion rate.

19 For secondary effectiveness outcomes, the  
20 procedure duration is now slightly increased, from 45 to 50  
21 minutes. The hospital stay, comparing LC to LNF, is now up  
22 from 0.4 to 0.5 days. For safety outcome, the overall  
23 complication rate, the delta here was doubled, from 3.5 to 7  
24 percent. So, this is what the sponsor wrote in the table.

25 [Slide]



1           When I begin with procedure duration I see I get  
2 137 minutes for the mean for the Intuitive and 89 minutes  
3 mean for the control. With a predetermined delta of 50  
4 minutes this is not equivalent. However, I have to increase  
5 that by 13 minutes to just get equivalence.

6           [Slide]

7           Looking at the learning curve now, here things  
8 look a little bit different. The linear slopes, again using  
9 the sponsor's analysis for GEE, general estimating  
10 equations, you get the Intuitive reduction at 4.5 minutes  
11 per procedure and a control reduction at 2.1 minutes per  
12 procedure. These reduction rates are statistical  
13 significant so they really are coming down, but there is no  
14 statistically significant difference between Intuitive and  
15 control. Even though one is double the other, statistically  
16 they don't show up as being different. So, you may or may  
17 not conclude how easy it is going to be for Intuitive to  
18 catch up to control. Obviously, they are both coming down.  
19 They are eventually going to catch up to each other but we  
20 haven't seen beyond 19 procedures at most per surgeon. So,  
21 we really don't know where this is going but this looks good  
22 because the Intuitive is coming down apparently harder than  
23 the control.

24           [Slide]

25           To look at the hospital stay, here we get 1.4 days

1 for Intuitive and 3 days for the control. The predetermined  
2 delta of 0.5 days is not equivalent, and the reason for that  
3 anomaly, as the sponsor mentioned, is that there were two  
4 outliers in the control group which greatly inflated the  
5 variance. When I looked at this using non-parametrics I  
6 also got pretty much the same answer, and it would be  
7 possible to achieve equivalence by raising the minimum delta  
8 up to 0.93 days.

9 [Slide]

10 Going on to complication rate, here it was very  
11 high. This is all complications. It was 19 percent for  
12 Intuitive and 15 percent for control. The predetermined  
13 delta is 7 percent. This leads to the conclusion that this  
14 is not equivalent. The minimum delta that would just allow  
15 equivalence to be established is 15 percent.

16 Here I did something a little different. Another  
17 way to approach "what if" when you don't get to reject your  
18 null hypothesis is, well, suppose we want to get 19 percent  
19 Intuitive and 15 percent control, suppose we just raise the  
20 sample size and keep the delta. Well, the sample size would  
21 have to go through the roof before this would become  
22 equivalent. It would have to be 847.

23 So, the point of this is that if you try to do  
24 equivalence testing when the Intuitive is coming out  
25 slightly worse than the control it becomes a very difficult

1 task to do. In other words, it may be better to look at  
2 complication rates and conversion rates solely in terms of  
3 qualitatively looking at the rates and making sense out of  
4 it rather than trying to do a formal statistical analysis  
5 using equivalence testing.

6 [Slide]

7 Again, with the LNF blood loss, here it was very  
8 similar in numbers to what we saw with the LC, 18 ml for  
9 Intuitive and 9 ml for control. However, in this case it is  
10 not statistically significantly different.

11 [Slide]

12 As far as my conclusions go from the clinical  
13 study of statistical equivalence, with the one exception,  
14 that is the LC hospital stay, statistical equivalence was  
15 not demonstrated, using the clinical endpoints or deltas  
16 from the LC and LNF protocol tables.

17 Therefore, deltas may be increased, if clinically  
18 feasible to do so -- that is your decision -- or the sample  
19 sizes per group may be increased in an attempt to establish  
20 statistical equivalence.

21 That ends my presentation. I would like to turn  
22 the podium back over to Dwight Yen to provide the questions.

23 **FDA Questions**

24 MR. YEN: At this time, I would like to read the  
25 panel questions that were provided to the panel members back

1 in May. You also have a one-sheet summary of these  
2 questions in front of you.

3 Although not limited to these questions, these are  
4 the particular questions that we are asking for input from  
5 the panel during the panel deliberation following the  
6 presentation from the primary panel members. At that time  
7 we would be happy to clarify any additional issues that you  
8 may have.

9 The first question, please discuss what you  
10 consider are the benefits and the risks of this device for  
11 the intended use based on the preclinical and clinical  
12 performance data presented for laparoscopic cholecystectomy  
13 and laparoscopic Nissen fundoplication.

14 Question number two, please discuss whether or not  
15 you believe the net risk-benefit ratio adequately supports  
16 the use of this device for general laparoscopic surgeries.

17 Question number three, clinical use of this device  
18 during this study was limited to relatively healthy, adult  
19 patients, expected to benefit from elective LC and LNF  
20 procedures. Please discuss concerns, if any, regarding the  
21 use of this device in general surgical procedures in  
22 populations that may be vulnerable to increased blood loss  
23 and/or procedure duration associated with device use, for  
24 example, patients requiring emergent intervention, pediatric  
25 patients, elderly or small adults.

1           Question number four, the sponsor of this device  
2 has claimed that the device is fail-safe. That is, the  
3 device has adequate safeguards built into both the hardware  
4 and software such that any failure of the device during  
5 surgery will not introduce unacceptable risk to the patient  
6 or surgeon. Please discuss the adequacy of device fail-safe  
7 design for this intended use.

8           The last question, limited clinical experiences  
9 indicate that surgeons and surgical teams need to be  
10 properly trained to use this device. Please discuss the  
11 types of training that will be warranted to assure that the  
12 device will be safe and effective.

13           Thank you.

14           DR. WHALEN: Thank you, Mr. Yen. We will now have  
15 the panel deliberations and comments. As stated in the  
16 agenda, we will begin with three scheduled panel member  
17 presentations, which will include Dr. Hannaford on a  
18 technical overview of the PMA, followed by Dr. Talamini with  
19 a clinical overview, and ending with Dr. DeMets on the  
20 statistics of the submission. Dr. Hannaford?

21                           **Panel Deliberations**

22                           **Preclinical Overview**

23           DR. HANNAFORD: I really like using Power Point  
24 but it is still an immature technology and I am afraid to  
25 count on it for an important meeting. So, I went to Kinko's

1 last night and printed it out. Once you take a floppy disc  
2 to Kinko's I always destroy it because you get so many  
3 viruses at Kinko's. So, now I don't have a floppy to stick  
4 in your machine. So, I will just go with the old-fashioned  
5 system.

6 [Slide]

7 This is my first panel session. I am just a  
8 consulting member of the panel, and what I have done is  
9 divided this talk into three segments. The first is a very  
10 generic overview of this technical area. In the second  
11 segment I am talking about the ISI system as I understand it  
12 from the material that I got, and not having the benefit of  
13 anything I have heard today, of course. The third section  
14 is sort of what I recommend as action. So, if you want you  
15 can stop me if this isn't the appropriate time for that  
16 since I don't quite know the protocol yet.

17 [Slide]

18 This is an area which I think is new to medical  
19 devices but not new to technology, and it started in the  
20 Manhattan Project era when they had to handle very dangerous  
21 nuclear materials for building nuclear weapons, and they had  
22 to do dexterous things to them in order to shape them into  
23 bomb components.

24 So, they invented these mechanical teleoperators  
25 or through-the-wall devices, called waldos, as early as the

1 1940s. But these devices were characterized by a fully  
2 mechanical implementation. So, there was a mechanical link  
3 between the operator's hands and the manipulators on the  
4 other side.

5 But then, for a couple of reasons, it was  
6 advantageous to make that an electronic link instead of a  
7 mechanical link, for example, if you had to exceed a 15 ft  
8 distance or you had to have a moving relationship between  
9 the two sides. So, they developed them in an electronic  
10 form.

11 What you are seeing, I think, today is sort of a  
12 1990s version of these systems which are still in wide use.  
13 Even the mechanical systems are still in wide use and work  
14 very well. Then, in research we saw the introduction of  
15 computers into these systems and now we are seeing  
16 telerobots controlled on the Internet in very beginning  
17 demonstrations.

18 [Slide]

19 You heard the term bilateral, some of you did.  
20 This topic was discussed in the closed session in response  
21 to my questions. This is a well-known term; it is not a  
22 secret. This term falls into the class of bilateral  
23 teleoperators because it creates this virtual physical  
24 connection between, in this case, the surgeon and the  
25 patient's tissue.

1           Here are some technical issues in that kind of  
2 system that have to be addressed, and certainly have been in  
3 this product for the most part. Mechanical design has to be  
4 done very carefully to achieve those attributes so there is  
5 any chance of this force information passing back and forth  
6 in both directions.

7           The visual registration between the operator's  
8 hands and the motion of the robotics system is very  
9 important for performance.

10          Force feedback is very important as well, at least  
11 in some cases. When they initially remotely controlled  
12 these mechanical devices in the '50s or late '40s, it was a  
13 position only control. That was totally rejected by the  
14 operators because they didn't have any feel.

15          Now we are in a similar situation as open surgery  
16 converted, say, ten years ago to laparoscopic surgery. The  
17 surgeons have lost the ability to, say, palpate tissues and  
18 make certain kinds of discriminations.

19          Finally, a technical issue is control properties  
20 of this system. Because it is computer controlled or  
21 electronically controlled, what is the performance and the  
22 quality and stability of that system?

23                 [Slide]

24          To get more into control, bilateral control means  
25 a control system where there are two points where the system



1 interacts with the physical world, or ports is a term that  
2 can be used. Here, we are talking about the surgeon  
3 touching a handle at one end and the tool touching the  
4 patient's tissue at the other end. So, there is some kind  
5 of control equation buried in there, somewhere, that has two  
6 inputs, one from each of these ports, and two outputs to  
7 each of these ports. If that is done properly it creates  
8 this virtual physical link. So, if I move the system until  
9 it comes into contact I should feel some force that stops  
10 contact.

11 I think it is interesting because I think we will  
12 see more -- as far as I know, this is the first commercial  
13 device employing this technology and I expect to see more of  
14 them in the future. But it is not something where there are  
15 established standards or protocols for validating  
16 performance and safety yet.

17 [Slide]

18 What are some aspects of performance? Can the  
19 operator distinctly feel some contact with a soft or hard  
20 object? Does the operator feel impeded or damped or some  
21 kind of resistance when they are moving in free motion?  
22 When I move in real free motion I don't feel anything  
23 significant. Is the system free of vibration or  
24 oscillation? Is the force accurately controlled at the end  
25 effector when it comes in contact with something? And, what

1 is the force scale? In other words, what is the ratio of  
2 the surgeon's force to the tool force?

3 [Slide]

4 Stability -- this is really the most challenging  
5 engineering area in these systems, and it is an emergent  
6 property of the entire system. You can't predict stability  
7 or analyze it without thinking about everything, including  
8 the surgeon's and the patient's biomechanical properties.  
9 So, it depends on all of these different things, and other  
10 ones.

11 [Slide]

12 Generally speaking, unstable systems are unsafe  
13 because they are not responding to the surgeon's inputs;  
14 they are doing their own thing. But there are safety  
15 measures that do work. Instability is really a cause of the  
16 system, say, applying too much force or moving too fast.  
17 So, if you have a safety system which detects those  
18 conditions, that is a valid way to address those risks.

19 [Slide]

20 The ISS system itself -- its name, as we have  
21 heard a couple of times, is an instrument control system,  
22 intended to assist in the accurate control of endoscopic  
23 instruments. But just to clarify, it is not an assist in  
24 the sense of power steering on your car where you have a  
25 mechanical link to the steering gear and a power assist.

1 This is something different from that because there is no  
2 physical link in between.

3 It certainly is intended to follow the surgeon's  
4 motions and commands and, to some extent, to couple the  
5 response of the body back to the surgeon, but it controls  
6 all the energy flow into and out of the surgical tool. I am  
7 thinking of mechanical energy here, force, displacement and  
8 their product.

9 [Slide]

10 This is sort of a simplified block diagram  
11 abstracted out of the block diagrams that I have seen. The  
12 point is that this is the surgeon's side. There is a  
13 mechanism that contains motors and sensors. The motors are  
14 capable of pushing back on the surgeon, and the sensors can  
15 detect the surgeon's motion. Then, there are similar things  
16 over here, on the patient side, actuators, meaning motors  
17 that can drive the instrument in all its different  
18 directions and sensors which measure what the instrument is  
19 doing, and there is redundancy in these for safety. But the  
20 key thing is these are connected not by any physical link  
21 but by the computer's hardware and software.

22 [Slide]

23 I would classify that as a bilateral impedance-  
24 controlled, force-reflecting teleoperator. Again, this is  
25 based not on complete information but what I have in the

1 distribution. The worst case for stability -- as I pointed  
2 out, the stability of this system depends on what it is in  
3 contact with and sometimes, for some architectures, the  
4 worst case will be when it is in contact with something  
5 rigid. For other architectures it is when the system is in  
6 free motion. So, this is not a pejorative term about this  
7 system; it is just how you have to analyze it for stability.  
8 That case is when both sides are free. The best case is  
9 when both sides are touching something or are constrained by  
10 something. Each possible architecture has a worst case and  
11 best case. Again, this is not either pro or against the  
12 system we are talking about today.

13 [Slide]

14 I guess I will continue then with what I think  
15 about the system and recommend for this panel. First of  
16 all, there is the issue of level of concern -- again, I am  
17 new to this area but from what I read, it applies primarily  
18 to the software, or exclusively to the software part of the  
19 system, and some of the documents, I guess, are hand-me-  
20 downs from the previous filing so they mention the lowest  
21 level of concern but, clearly, just to "disambiguate" that,  
22 and I think that is the case here, as we change from blunt  
23 tools to sharp tools the same amount of physical energy can  
24 do a lot more damage. So, I agree that a concern of  
25 moderate is warranted.

1           Again, I think in principle the system could with  
2 those sharp tools do a lot of damage but I think that the  
3 level of concern applies to software only. The hardware in  
4 the system is a good check on the software, as I understand  
5 it. So, I think describing the level of concern of software  
6 as moderate makes sense here.

7           [Slide]

8           Issues that I have -- first of all, in looking at  
9 the videos I thought I saw instances of instability, and I  
10 certainly have not done any careful, exhaustive study of a  
11 large sample of material but I did see cases where I saw the  
12 end-effector -- not on all the tapes -- vibrating and that  
13 gives me a concern.

14           Kind of along with that, there is really very  
15 little documentation at all about the control system in what  
16 I have seen. So that gives me a concern. That includes the  
17 gain-stability margins, instances of instability, and even  
18 specifications analysis or testing is completely absent.  
19 Now, this may have to do with proprietary information but I  
20 think that is a concern.

21           [Slide]

22           These are some things I would like to see  
23 documented on the system, namely the control equations, or  
24 at least the nature of the inputs and outputs used in the  
25 control; the stability analysis; analysis of the worst case

1 loads for stability; the maximum stable gains. These are  
2 technical things but I really think what you are seeing here  
3 is a brand-new introduction of technology into surgery. Of  
4 course, I am very excited about it. I agree with the  
5 statements that the sponsor has made about the future  
6 potential of this system. So, I think it is in everybody's  
7 interest to document this first foothold of a very important  
8 new technology.

9           Then, the gain values that are used in production  
10 compared to, say, this maximum stable gain. The other two  
11 important related features are tremor filtering and scaling,  
12 which are mentioned with literally about three words in the  
13 filing but not documented in any detail. For example, what  
14 frequencies are passed through and what frequencies are  
15 attenuated when tremor is filtered by the system.

16           [Slide]

17           Finally, the questions -- I really don't have any  
18 expertise pertaining to the first three questions so I will  
19 just confine it to questions four and five. I think the  
20 system does show --

21           DR. WHALEN: Dr. Hannaford, forgive me for  
22 interrupting but we are going to ask for each panel member  
23 to make comments upon the questions as they come up.

24           DR. HANNAFORD: Fine. In that case I am finished.

25           DR. WHALEN: Thank you very much. Dr. Talamini?

## Clinical Overview

DR. TALAMINI: My name is Mark Talamini, Associate Professor at Johns Hopkins, and this is a non-statistician, non-engineer viewpoint of what I had before me to look at before this meeting. I am an active clinical surgeon. I do about 300 cases a year. About half are laparoscopic, advanced laparoscopic, and about half are open. So, I am not a laparoscopic evangelist for these kinds of technologies. On the other hand, my primary research interest is in the physiology of this type of surgery.

[Slide]

A lot of what I put together for this is already brief but I will make it briefer because I think it now falls into the category of beating a dead horse. We have already talked about advantages of laparoscopic surgery. I think they are clear and well-documented, both in surgeons' experience, patients' experience and in the literature at this stage.

[Slide]

We have talked today already also about the disadvantages of laparoscopic surgery, the things that are taken away when surgeons are operating in this fashion. As has been alluded to, there are substantial things that surgeons give up when they do this type of surgery.

[Slide]

1 In terms of the two procedures that were studied  
2 and have been talked about a lot today, again I will give  
3 you the general surgeon's perspective. A laparoscopic  
4 cholecystectomy, quite frankly, is a pretty easy operation  
5 to do. That is the reason that this operation was  
6 introduced in this country in late 1989 or 1990, but by 1992  
7 three-quarters of the cholecystectomies done in the State of  
8 Maryland were being done this way. Think about that. The  
9 entire general surgery work force trained to do this  
10 operation in two years. So that tells you two things: It  
11 tells you that it is a pretty easy operation and that  
12 surgeons are pretty adaptable.

13 [Slide]

14 Laparoscopic Nissen is different. The  
15 laparoscopic cholecystectomy is an exterpative operation.  
16 You are taking something out so there are things to ligate,  
17 things to divide. The Nissen is a different operation. It  
18 is more difficult technically. It takes a longer period of  
19 time. It is a reconstructive operation where the surgeon is  
20 required to rearrange tissues so it takes a different set of  
21 skills, a different level of skills. Because of that, the  
22 potential complications are more substantial. I think it is  
23 good that we have information for both of these types of  
24 procedures to look at today because they are different, and  
25 they look at different aspects of surgery.



1 [Slide]

2 The whole idea of computer-assisted or robotics in  
3 surgery is an interesting one. On this slide I have put  
4 down some ideas that do not pertain specifically to the  
5 application for today but to this whole field in general.  
6 Those are that laparoscopic surgery really is ideal for the  
7 application of this sort of thing, high technology  
8 applications, for a number of reasons. It depends on the  
9 image stream and that image stream can be digitized,  
10 manipulated and altered. Information can be added to it.  
11 But, really, the promise for us, surgeons, is that high  
12 technology would not replace us but enhance what we do.

13 Now, the HMOs might like to replace us. There  
14 might be a cheaper alternative, I don't know. But we  
15 believe really that the promise is that as we move forward  
16 and apply technology, that is out there in other industries,  
17 to what we do, we will be able to do things better than we  
18 do them now. In the liver we will be able to see where  
19 blood vessels are beneath the surface of the liver, even  
20 though they can't be felt or observed in any other way. The  
21 idea of being able to tell what type of tissue might be  
22 tumor and what type is not by using different ways of  
23 evaluating the surfaces of those tissues, force feedback.

24 Telemedicine I think is a controversial area  
25 because I don't know of any patient of mine or one here, in

1 Gaithersburg, that would like to have their surgeon other  
2 than at the operating room table. But it has been an area  
3 of growth, and certainly one published in the literature,  
4 and it is one that the military is very interested in.

5 [Slide]

6 In terms of the current application and, again, I  
7 don't want to beat a dead horse because all of these issues  
8 have been alluded to and I think we are about to discuss  
9 them, the learning curve effect is important and it sounds  
10 like, from our presentation so far, there are at least two  
11 different ways to look at the learning curve effect with  
12 this application. The experience of laparoscopic  
13 cholecystectomy is an important one in that surgeons did  
14 learn that technology fairly quickly.

15 The blood loss aspect, again, I think it is  
16 significant but we are talking about very small levels of  
17 blood loss. I was surprised that in the study the  
18 investigators were actually able to measure blood losses  
19 that small given these types of procedures. It is very  
20 difficult for me at the end of an operation to know if I  
21 have less than 50 cc of blood loss, let alone discriminate  
22 between 10 and 20 or 30. So, it does appear that there is a  
23 significant difference between the control groups and the  
24 study groups but it is at a level that has me scratching my  
25 head a little bit as a clinical surgeon.

1           The complication rates and, again, we have already  
2 talked about this -- are certain complication rates more  
3 pertinent in evaluating this than others? The investigators  
4 have certainly given their opinions on which are more  
5 important and which are not.

6           [Slide]

7           This part is not beating a dead horse. This is  
8 what I think is new for me as a surgeon. For me, in  
9 evaluating the mass of material that was sent to my office,  
10 the video data was by far the most important. That is  
11 probably because as a surgeon that is what I do every day.  
12 So, in looking at those videos I could think about and  
13 experience what I do every day.

14           I think the statistical analysis and the study  
15 obviously is important and it is meant to do a certain  
16 thing, but for me to see exactly what surgeons were able to  
17 accomplish during operations was extremely important.

18           I have summarized here what I thought I saw this  
19 system able to accomplish, and I looked at many hours of  
20 videotapes, more than I think the rest of the panel was able  
21 to. The system did have the ability to make quite precise  
22 motions, and I saw in a number of instances that the  
23 surgeons were able to pick a target, very, very small, and  
24 go in and exactly pick it up or precisely manipulate it, to  
25 my eye, more effectively than I usually can do with a

1 chopstick instrument at a 1.5 ft distance.

2 I also was impressed in a number of instances,  
3 when bleeding did occur during these operations with this  
4 system, that the surgeon was able to take this instrument  
5 and very precisely grab right where that blood vessel was  
6 and control that bleeding, which again left me scratching my  
7 head a little bit about the differences in the statistical  
8 analysis for blood loss because on the videos, it seemed to  
9 me, most of the time the surgeons were able to very  
10 precisely control bleeding.

11 Also, despite the issue of force feedback, I think  
12 the strategy of giving very good visualization and good  
13 magnification and three-dimensional vision allow the  
14 surgeons to very gently dissect tissues, in my view, looking  
15 at these videos.

16 I thought they were able to tie the knots very  
17 effectively. There were some videos where I saw the sutures  
18 slip just a little bit, others where I didn't see that and I  
19 imagine there may have been some technical changes from one  
20 to the other in the design of the instruments or the  
21 settings in terms of the amount of pressure they were  
22 putting. But, certainly, in some the knots were nice and  
23 tight and there was no slippage; in others there was a  
24 little bit of slippage -- not slippage of the knot itself  
25 but slippage of the grasper on the suture. But I could see

1 that the surgeon recognized that and worked to make sure  
2 that the knot was still tight despite that.

3 Finally, the other thing that I saw that I think  
4 is perhaps the power of the system from the surgeon's point  
5 of view is the similarity of this system to actual hand  
6 motions. Of course, that is because of the multiple degrees  
7 of freedom.

8 That is all I have from a formal point of view.  
9 From a surgeon's perspective, the other point that I would  
10 make is about this whole issue of training. I want to just  
11 make two points in terms of training. One is the  
12 laparoscopic cholecystectomy experience. You know, lap  
13 chole was never exposed to a randomized, prospective trial.  
14 I think general surgeons were actually fairly lucky in that  
15 it turned out to be a safer, better operation because that  
16 whole thing happened so fast that nobody really ever had the  
17 time to do an important randomized, prospective trial which,  
18 of course, is the gold standard for all of us, clinicians.

19 Despite that, surgeons learned the procedure and,  
20 with the exception of some increase in complications, that  
21 has really been a success story for general surgery, yet, it  
22 was not legislated. There was no committee saying, yes, you  
23 can do this or, no, you can't do that.

24 So, I have that on the one hand of my experience  
25 history. On the other hand of my experience history I have

1 the whole laser thing. Remember, when we started  
2 laparoscopic cholecystectomy it was called laser  
3 laparoscopic cholecystectomy, oftentimes primarily for  
4 marketing. The laser turned out not to be very effective,  
5 in fact, perhaps dangerous in some instances, and that was  
6 device based. So, that was a different situation. I don't  
7 understand what was going on with the FDA at that point. I  
8 think these were all already approved devices being used for  
9 cholecystectomies.

10 So, I am not quite sure in my own mind yet where  
11 to put this application in that experience of mine, with lap  
12 chole on the one hand, which worked out very well, and the  
13 laser application, which was the use of a device with that  
14 operation that I think did not work out so well and has now  
15 fallen by the wayside. So, I offer you those as my  
16 thoughts. Thanks.

17 DR. WHALEN: Thank you, Dr. Talamini. The final  
18 scheduled presentation, Dr. DeMets on the statistical  
19 overview.

#### 20 Statistical Overview

21 DR. DEMETS: Thank you. I am going to have a very  
22 low tech presentation, and a very brief presentation, partly  
23 because I was heading out the door on a week-long trip when  
24 Dr. Krause tipped me off that I might have to make some  
25 remarks. So, I have resorted to the low tech approach.

1           What I thought I would do is spend a little time  
2 trying to set the stage for the panel about this active  
3 control equivalence business because I think it is something  
4 that either gives you a headache or makes you glassy eyed a  
5 little bit.

6           [Slide]

7           First of all, let me say it is a very challenging  
8 problem. Statisticians have spent most of their time  
9 thinking about how researchers detect differences, find  
10 differences, and the process that we are talking about is  
11 somewhat different than that. There is not a lot of  
12 literature on it and there is not a single agreed upon way  
13 to do this, although I will have a suggestion.

14           One of the things that happens to you is that  
15 things get reversed. If you think about superiority trials,  
16 which is what we typically are thinking about, versus non-  
17 inferiority or maybe equivalence, the process gets reversed.  
18 Let me try to explain it very simply. Dr. Bushar got to it  
19 just a little bit. In the typical situation you have a null  
20 hypothesis of no difference in event rates, or mean values,  
21 or something, and you have an alternative hypothesis where  
22 you are trying to detect some size difference and you power  
23 your trial to detect that. That is the usual process we go  
24 through.

25           In the active control versus experimental

1 situation or the equivalence you are flipping things around.  
2 The null hypothesis is that, in fact, there is a difference  
3 of a certain size; that the two treatments differ by some  
4 amount. I call that delta minimum, and I want to come back  
5 to that. What you are trying to do is establish that, in  
6 fact, the difference that exists is less than that delta.  
7 That is why the delta business is floating throughout the  
8 discussion of today.

9 In either case of superiority or equivalence you  
10 have to have a large enough study to detect the differences  
11 that you think are important. The size of that difference  
12 is not a statistical question; it is a clinical question, a  
13 clinical decision question. It has nothing to do with  
14 statistics. How you determine it and how you detect it and  
15 rule it out is a statistical problem. So, I guess, that is  
16 point number one.

17 [Slide]

18 Two is in the material and even, in some ways, in  
19 the discussion today. There has been some back and forth  
20 about what scale of reference do we use. In fact,  
21 throughout the FDA these days there are several pockets of  
22 discussion going on, whether it is devices or whether it is  
23 drugs or biologics -- I don't know about devices so much,  
24 but there is no agreed upon scale of reference.

25 What do I mean by that? Well, regardless of



1 whether it is binomial, or yes or no, or success or failure  
2 outcome or some continuous variable, one can express those  
3 differences as either an absolute difference, absolute  
4 difference in rates, absolute difference in events or  
5 measurements, or on a relative scale, such as a relative  
6 risk, the ratio of the success rates or the failure rates,  
7 or percent of change relative to the control arm. It  
8 doesn't really matter which scale we choose. It doesn't  
9 matter whether it is continuous or binomial. But we  
10 probably some day, within this panel or the FDA in general,  
11 should kind of try and get some consensus whether the  
12 language is conducted in English or in French because that  
13 is really all that is involved here -- what language, which  
14 scale do we want to talk about so we are not confusing each  
15 other every time we pick up a different variable.

16           The fourth point is that to think about this kind  
17 of a problem you have to specify what I call the minimum  
18 delta. What size delta for a particular variable -- and you  
19 have to remember which variable and what kind of scale,  
20 absolute or relative, and that is how you express the delta  
21 -- do you think is minimally of clinical interest? In other  
22 words, where would you walk away from the process thinking  
23 it doesn't matter to me? That is a clinical decision that  
24 you have to figure out. It doesn't depend so much on  
25 statistical issues at this point, and that is a key issue.

1 So some of the discussion that was presented by the FDA  
2 review had to do with what difference was stated versus what  
3 we could have detected, and the presentation by the sponsor  
4 also alluded to that a little bit.

5 [Slide]

6 The only way that I can keep my thoughts straight  
7 about this whole process is reflected best in a paper that  
8 Tom Fleming wrote almost ten years ago, when we were  
9 debating and discussing AIDS trials, but it is very generic  
10 methodology. I was disappointed in the presentation this  
11 morning because I thought, based on what you handed out, you  
12 came very close to this method.

13 [Slide]

14 Let me try and explain this picture, which is the  
15 only thing I had time to xerox before I walked out the door.  
16 I think it is a very simple, clear way to think about what  
17 is on the table for our discussion. In the classical  
18 setting -- and I am doing this on a ratio, a relative scale.  
19 This is a figure from Tom Fleming's article. So, 100  
20 percent on the scale says that the two treatments have  
21 absolutely the same effect, whether it is event rates,  
22 failure rates, success rates or measurements. What you want  
23 to do is to show that the difference that you observe, which  
24 in his little example was 125 percent, plus or minus two  
25 standard errors actually exceed delta. That would give you

1 a significant result, with a p value less than 0.05.

2 So, this is the confidence interval. It excludes  
3 100 percent. You say, well, we are better than that. So,  
4 we have an improvement. If we don't have a confidence  
5 interval that excludes 100 percent or a ratio of 1, you say  
6 we can't rule out that they are the same. So, we can't  
7 claim superiority.

8 In our setting, remember, we are flipped around.  
9 We are down here. So, here 100 percent says that the two  
10 therapies are exactly equal and we don't want to be worse  
11 than that by very much. That is this delta minimum. So,  
12 there is some delta below the 100 percent that we will pick.  
13 And, he has it sitting right here in his diagram. Of  
14 course, there is an opposite side delta which is not so  
15 interesting for the discussion here.

16 So, if you do a study, do the confidence interval,  
17 and in this case the lower limit of the confidence interval  
18 is above that delta minimum, you have ruled out that you  
19 could be worse off by this much or below. You have ruled  
20 out you are not worse than the standard you set and, in  
21 fact, you are even a little bit better than that. On the  
22 other hand, this confidence limit includes that delta. So,  
23 you haven't ruled out that you could be that much worse.

24 So, this is a very simple, comfortable approach  
25 that expresses the data you got, and by looking at that

1 confidence interval you can see what is in and what is out.  
2 What have you ruled out and what can't you rule out?

3           Then you say, okay, what delta did I really care  
4 about? Where do I clinically walk away from this? Thus, I  
5 would say if there is anything that is converging as a point  
6 of view among biostatisticians that I hang out with, this  
7 approach is seemingly where we are heading. The reason I  
8 was disappointed in the presentation this morning was that  
9 in the materials handed out to us, not presented but handed  
10 out, there are confidence intervals for the experimental  
11 group and the control arm and you can see the variability  
12 and the overlap and what you can rule out. It is not  
13 expressed as differences, which is what I prefer, but at  
14 least it is expressed in terms of confidence intervals.

15           That is kind of how I sort of think about this.  
16 So, the presentations about what delta are we arguing about  
17 really has to do with what you think is clinically  
18 important. The way this process should go, in my opinion,  
19 is that you should say what is the question, and what  
20 measurement are you going to use for that question, and what  
21 scale. What is clinically relevant? Select the scale in  
22 which you are going to present things. Select the delta  
23 minimum that, again, is based on clinical considerations,  
24 not statistical. Size your trial accordingly so you can  
25 detect that minimum. Present your results in this

1 confidence interval approach, and then say do I rule out the  
2 delta minimum that I a priori specified, which presumably  
3 has some clinical relevance, and you can make your decision  
4 whether you make your case.

5 I was somewhat curious, and perhaps that could be  
6 clarified, that the sample size for this study seems to have  
7 been said in advance at 50 per arm, without a lot of  
8 consideration of what is the question. And, the deltas are  
9 sometimes specified and sometimes not. So, I am a little  
10 unclear about that.

11 I think that the issue about counting all patients  
12 in studies is one that we settled long ago. That is, you  
13 have to account for all patients that you entered, and to  
14 try, as tempting as it might seem, to talk your way out that  
15 some patients didn't sort of fit the criteria has a problem  
16 which was alluded to earlier. That is, if you get rid of  
17 all the tough cases then, of course, things get better.  
18 This debate in the surgical community took place intensely  
19 in the CABG debates about the patients who didn't get CABG  
20 shouldn't be counted. Well, okay, if you get rid of all the  
21 sick patients then, of course, I am going to look better.  
22 We have been through that debate and I think we have  
23 convinced ourselves that it is very difficult -- you don't  
24 have comparability anymore and you cannot call it a  
25 randomized trial.

1           My final comment is that the learning curve, based  
2 on the data I observed, seems to be longer than being in  
3 Mexico City for one day and then the real stuff the next  
4 day, even though I know surgeons are very adaptable. If you  
5 look at the data, it takes ten to twelve, to twenty patients  
6 before the data starts settling down. So, perhaps another  
7 comment on this design is that the learning curve for the  
8 experiment that we have seen wasn't long enough.

9           But my main point is to go back to the way you  
10 think about the data that you have before you, which is to  
11 look at what delta is of clinical value and importance to  
12 you. So, I would actually draw your attention in the  
13 discussion to what was not presented today, and maybe the  
14 sponsor could pull those slides up. But that is how I tried  
15 to sort this out. It is a tough problem. If you are having  
16 trouble with this as a clinician, don't worry; the  
17 statisticians have also been having trouble addressing this  
18 problem.

19           That is all I have for formal remarks.

20           DR. WHALEN: Thank you, Dr. DeMets. Before the  
21 panel tackles the first question, as promised this morning,  
22 this is the opportunity that the panelists have to ask  
23 questions either of the sponsor or of the FDA if, indeed,  
24 any of those questions linger. Dr. Burns?

25           DR. BURNS: I have a question about exactly what

1 it is that we are reviewing for approval. Are we looking at  
2 the instruments as well as the system for approval? Is that  
3 all part of the system?

4 DR. WITTEN: I am not sure what your question is.

5 DR. BURNS: Well, there are instruments that are  
6 used with the system, and those are not standard  
7 laparoscopic instruments, as I understand it.

8 DR. WITTEN: It is the instruments plus the  
9 indication change from the already cleared part of the  
10 product.

11 DR. BURNS: So, under approval it is the  
12 instruments that would go with the system.

13 DR. WITTEN: Well, it is the system plus the  
14 indication that they now want. They presented what the  
15 indication was they have and what indication they want. So,  
16 it is the system plus the indication.

17 DR. WALKER: I guess we still haven't heard the  
18 answer to the question we are both asking, which is no tests  
19 were done on harmonic scalpels. Are we being asked to  
20 approve the harmonic scalpel as a factor of this device  
21 today, or will you come back to the FDA later for harmonic  
22 scalpels? Can the sponsor enumerate for us the factors that  
23 we are being asked to approve today?

24 DR. MOLL: I wish we had that list back up. There  
25 were some very specific tools to accomplish electrocautery,

1 graspers, needle drivers, scissors, and that it.

2 DR. WALKER: No scalpels of any kind?

3 DR. TALAMINI: It says sharp endoscopic  
4 dissectors, scissors, scalpels, forceps --

5 DR. MOLL: The added new instruments are forceps,  
6 scissors, scalpels, clip appliers, needle holders and  
7 electrocautery.

8 DR. WALKER: But scalpels here excludes harmonic  
9 scalpels.

10 DR. MOLL: Yes, absolutely.

11 MS. DUBLER: And it is my understanding that you  
12 are asking for approval for those instruments for the  
13 purposes discussed today.

14 DR. MOLL: I believe that is correct. It is  
15 approval of the instruments for the stated purpose of using  
16 them in laparoscopic surgery.

17 MS. DUBLER: So, if you wanted to use these same  
18 instruments for cardiac surgery you would feel compelled to  
19 come back to the FDA.

20 DR. MOLL: Correct.

21 DR. FERGUSON: I am impressed with the restoration  
22 of the surgical skills that this device might provide  
23 surgeons after we have been struggling with chopsticks, as  
24 he says, for a number of years. But my question is why does  
25 it take so much longer since all surgeons have learned those



1 skills when they were learning to be surgeons. Why does it  
2 take that much longer to do the cases when we have now  
3 forces applied that have been returned to us, and the  
4 ability to use scissors and so on?

5 DR. MOLL: Well, I think Dr. Gardiner attempted to  
6 address that, and we believe that the surgeon is given more  
7 dexterity but he is also given a new surgical tool, and with  
8 that tool and that system comes other types of new things  
9 that they need to learn about and be comfortable with. Like  
10 with any surgical device, there is training involved and a  
11 learning curve.

12 DR. FERGUSON: That doesn't answer my question. I  
13 am sorry.

14 DR. GARDINER: I think that we see on both sides  
15 of the study learning going on. Having been one of the four  
16 people that did the surgery, I felt I was, to some degree,  
17 learning again. I was using a new assistant, and I got  
18 better as I went along in doing the conventional  
19 laparoscopic that I had done thousands of times. I was  
20 using a new assistant who was unfamiliar to me, and I was  
21 using different instruments in a different location. So,  
22 there is learning that was going on, actually, on both arms  
23 of this study.

24 With regard specifically to the Intuitive device,  
25 you are learning a couple of things. There is learning that

1 is going on in terms of how to drive this system, how to  
2 operate with it, how does the surgeon actually use the  
3 console. That is one aspect of the learning. There is also  
4 learning that is going on in terms of how the assistant  
5 interacts with the device and the patient, and how the  
6 assistant and the surgeon interact together. And, I think  
7 that those factors together make this a little bit slower,  
8 and I think it is going to get better and faster as we get  
9 more experience. It certainly has with every other  
10 laparoscopic procedure I have done.

11 DR. FERGUSON: I guess what I am looking for -- I  
12 am not trying to be pejorative about it but I am looking for  
13 the positive fact that if you give surgeons back their  
14 surgical skills that is a good thing.

15 DR. GARDINER: Well, having sat down at this  
16 console and operated, and then going to the control patients  
17 and operating, there is no question, all four of us by the  
18 end of the study, when the card would get pulled up and we  
19 saw we had to do a control case, every one of us said, "oh,  
20 gee, I'd rather do this with the Intuitive system." Every  
21 one of us felt that way.

22 I think if you look at it and say what would I  
23 rather have? Would I rather have a more dexterous  
24 instrument or a less dexterous instrument; a more flexible  
25 instrument or a less flexible instrument; a more accurate

1 and precise one or a less accurate and precise one, as Dr.  
2 Talamini raised, you know, in the videos, I think we would  
3 all take the more capable equipment every time.

4 DR. TALAMINI: I really don't want to sound like  
5 an advocate but I don't think speed is always the end-all  
6 and be-all. One of Halstead's great contributions was to  
7 say, look, we can slow down and do things better. So, there  
8 certainly can be instances where a different technology that  
9 takes more time is superior to a quicker technology. I am  
10 not saying that is the case here but we need to have our  
11 minds open to that possibility I think.

12 DR. GARDINER: Having been responsible for doing a  
13 lot of these operations, I can tell you that was exactly  
14 going on. There is no question about it. We certainly  
15 didn't feel in a race and we weren't in a race. We were  
16 taking our time.

17 DR. CRITTENDEN: I am still troubled though with  
18 the fact that this should be better because it has an  
19 articulated wrist but if you look at some of the data that  
20 the FDA compiled, the person who did the most study  
21 laparoscopic cholecystectomies never really approached the  
22 conventional time, and these are people, by your own  
23 admission, who were the most advanced laparoscopic surgeons  
24 available. So, I just wonder if you turn this over to  
25 someone who is coming out of their residency how they are

1 going to be able to adapt to this technology.

2           What is more, we have only looked at laparoscopic  
3 cholecystectomies which, by the panel's acknowledgment and  
4 by your acknowledgment, are pretty easy procedures. But now  
5 we are talking about potentially doing laparoscopic  
6 splenectomies and colectomies as well, and I wonder whether  
7 we should not just limit it to laparoscopic  
8 cholecystectomies and Nissens because that is what we have  
9 data for.

10           DR. GARDINER: I don't believe that the intent of  
11 the FDA was to regulate this device procedure by procedure,  
12 and I think there is some element of physician judgment and  
13 surgeon judgment that comes into play here. Beyond that, I  
14 am not sure how much further to go.

15           This is certainly more capable than a  
16 straightforward conventional laparoscopic tool. Does a lap  
17 chole show that? Really not. But you can see it in the  
18 suturing part of the Nissen, for example. I guess the worry  
19 I would have -- and I am kind of editorializing here and  
20 maybe I ought to sit down but, you know, you take this  
21 device and hold it back from a surgeon that could use it  
22 perfectly and wonderfully to do a common-duct exploration,  
23 for example, and suture the duct closed but hold that back  
24 because we haven't demonstrated that. I think that the  
25 equipment is highly capable of that kind of procedure, and

1 whether this is going to be capable or would be something  
2 you would use for a splenectomy is questionable. I think  
3 that really kind of goes to surgical judgment.

4 DR. WHALEN: Correct me if I am wrong, Dr. Witten,  
5 but we are discussing efficacy and safety with lap choles  
6 and lap Nissens. That is what the panel is charged to do.  
7 Correct?

8 DR. WITTEN: No, actually -- can you put up your  
9 indication statement? Maybe we could ask the sponsor to put  
10 up their proposed indication statement.

11 DR. TALAMINI: While you are doing that, if I  
12 could just ask one other thing, Dr. Gardiner, in the control  
13 lap choles did you tie the cystic duct or not?

14 DR. GARDINER: I am sorry, I didn't hear you.

15 DR. TALAMINI: In the control lap choles did you  
16 clip and tie the cystic ducts or did you just clip them?

17 DR. GARDINER: No, we did them both the same way.

18 DR. TALAMINI: So you tied them as well with the  
19 control lap choles?

20 DR. GARDINER: Right.

21 DR. HANNAFORD: This is also for Dr. Gardiner.  
22 Sorry to make you jump up and down. For us, non-surgeons --  
23 actually, this is a question for any other surgeon in the  
24 room as well, can you give me an idea of the rate of trocar  
25 injuries in normal endoscopic lap chole? You had a certain

1 set number of them in the control arm, and what is the  
2 typical rate for those injuries in the normal operation?

3 DR. GARDINER: It would be significantly below one  
4 percent. I mean, we talked about that earlier this morning.  
5 It is hard for me to understand where those came from. I  
6 mean, we are not talking about surgeons that are  
7 inexperienced. They don't have that kind of experience in  
8 their own practice. We observed two of them.

9 DR. TALAMINI: There should be 1 in 500 or 1 in  
10 1000.

11 DR. GARDINER: Yes, right.

12 DR. CHANG: Dr. Gardiner, again with due respect,  
13 I want to revisit one other question since if I were to  
14 present this as a proposal to our surgical OR suite to say  
15 this is a wonderful new instrument, there would be a cost-  
16 benefit analysis. Based on your data, how could we present  
17 how this system benefits the patient undergoing lap chole or  
18 lap Nissen fundoplication?

19 DR. GARDINER: Well, I think that in the lap chole  
20 you have an operation which is well established to be able  
21 to be done now with conventional instruments, and I don't  
22 think it is going to benefit the lap chole. You know, that  
23 was put up there and put in the study primarily to evaluate  
24 the use of this device in a basic operation, and that is  
25 what that does. How much capability this system is going to

1 provide to the surgeon I think is going to be answered over  
2 time. We may well find out, for example, that you put this  
3 piece of equipment in the hands of that resident that is  
4 coming out and you may well see that the resident that is  
5 given back more capability can operate better, faster and  
6 more efficiently because he or she does have capability that  
7 they can use, whereas the highly experienced surgeon may  
8 have been able to get around a lot of the limitations of  
9 conventional laparoscopy.

10 DR. WHITE: A real quick response, the cost-  
11 benefit analysis is way too early to approach at this time.  
12 We went through that same argument with laparoscopic  
13 cholecystectomy, subsequently laparoscopic Nissen and,  
14 frankly, every procedure out there and there are still some  
15 of them where that is an argument. Not universally but very  
16 significantly, over time the cost-benefit analysis has  
17 become known in these procedures. For this technology it is  
18 way too early to enter into that.

19 But most importantly, we picked lap chole and lap  
20 Nissen because, one, they are commonly done by the  
21 individuals involved but allowed us to use the device to do  
22 all these things we have listed up here. We could have  
23 included other operations but these were the ones that we  
24 could universally use frequently and often, and demonstrate  
25 the capabilities of this device. How that is going to

1 figure into the cost-benefit analysis is way too early.

2 DR. WHALEN: The slide that the sponsor is  
3 projecting, requested by Dr. Witten, makes me stand  
4 corrected. It generically says laparoscopic procedures, just  
5 to point that out.

6 DR. FERGUSON: I really need a clarification about  
7 what we are voting on here because I came in with the idea  
8 that we are approving surgical instruments, the ones that  
9 are listed here, for any operation that they might wish to  
10 do with this new device. Yet, I hear them saying that they  
11 are going to bring back coronary bypass at another time  
12 where they can use these same instruments.

13 DR. GARDINER: No, what we are dealing with is  
14 laparoscopic indications, indications in the abdomen, not in  
15 the chest.

16 DR. FERGUSON: Could I ask then why does the  
17 request not limit itself to those two operations because you  
18 can use all of these surgical instruments to do a coronary  
19 bypass if you wished to do so.

20 DR. DILLARD: If I might be recognized, Dr.  
21 Whalen?

22 DR. WHALEN: Yes.

23 DR. DILLARD: This is Jim Dillard from the Food  
24 and Drug Administration. I think that is an issue that is  
25 important to talk a little bit about. Let me come from the



1 FDA perspective about how we handle indications for use and  
2 how we look at the various types of data.

3 Obviously from the FDA perspective, we regulate  
4 the medical device and we regulate the labeling and what the  
5 labeling says. Where we tend not to get too involved is in  
6 the practice of medicine. So, what we are asking you to do  
7 here today, and this is from our perspective and where we  
8 come from, is that we do not want to be looking at newer  
9 technology procedure by procedure, disease state by disease  
10 state to try to approve devices as they become validated in  
11 each individual type of operative procedure.

12 From this perspective, we tried to look at very  
13 representative general surgical kinds of procedures, picked  
14 a couple that would be representative of how the device  
15 would be used, would stress the system adequately so that we  
16 would get an understanding of the various types of  
17 instrumentation, and really from our perspective looking at  
18 this, we would be saying, since we will ultimately have to  
19 send a letter of whether we approve it or whether we not  
20 approve it, that if we approve it, it would be for these  
21 particular kinds of instruments, in this case the way the  
22 intended use is written, for a laparoscopic surgical kind of  
23 a procedure, but we would also be very specific in the  
24 labeling about the kinds of studies and the kinds of data  
25 and the two models that we had that went into the

1    approvability of the product in that labeling situation.

2                   We would not be saying, nor would we endorse if  
3    the manufacturer went out and started promoting the device  
4    for thoracoscopic, minimally invasive cardiac surgical  
5    procedures. We would say you are beyond the scope of your  
6    labeling; we did not approve that in the labeling. So, if  
7    the manufacturer wants to specifically state they can be  
8    used in other surgical procedures or in other surgical  
9    subspecialties, and they want that to be approved on the  
10   labeling the expectation is that they would present more  
11   data to us that we would then look at. Does that help?

12                  DR. FERGUSON: Very much.

13                  DR. DILLARD: Good.

14                  DR. WHALEN: Other questions?

15                  MS. DUBLER: I just want to follow-up on one  
16   discussion and the cost-benefit analysis is the basis for  
17   it. We have talked about individual benefit, and there is a  
18   school of ethical analysis these days that argues the  
19   following: that it is, in fact, unethical to approve new  
20   technologies that will add to the cost of medicine, given  
21   the number of people without health insurance and access to  
22   health care, unless there is a measurable benefit that  
23   proceeds from that technology.

24                  Now, I think we have been told that we can't  
25   assess that benefit at this point, and we can't produce the

1 cost-benefit analysis, but I would simply like to register  
2 my discomfort that we may be adding a substantial cost  
3 without commensurate benefit.

4 DR. WITTEN: I would like to make a comment from  
5 the FDA perspective. I think these are very important  
6 issues but from our perspective at the FDA in terms of what  
7 we are obligated to do to look at device approval for a  
8 device, it is to look at safety and effectiveness of the  
9 product. Although I think these questions about cost  
10 certainly are important, I don't think this is probably the  
11 place where we are going to hash them out since we don't  
12 really require the cost information as a part of the  
13 application from the sponsor for device approval.

14 **FDA Questions**

15 DR. WHALEN: Thank you. I would like to thank the  
16 sponsors for being in a position to answer those questions.  
17 We are about to embrace the questions that are put to the  
18 panel by the FDA. The questions have already been read into  
19 the record, and they will be projected as we encounter them  
20 so I will not reread them.

21 The process that we will follow is this, each  
22 member of the panel is going to be asked to comment upon the  
23 particular question. I will, in staggered fashion, start  
24 with a different individual each time so you won't be left  
25 as the twelfth person to try to say the same thing in a

1 different way every time. Please do not feel absolutely  
2 obligated to articulate words if they are only a rehashing  
3 of what has already been said. A simple "I agree" would be  
4 appreciated by anyone who has a flight in the next three  
5 hours, I am sure. But if you have something important to  
6 say, please, by all means, say it.

7           Following everyone's opportunity to comment, I  
8 will attempt to distill the consensus into a precise couple  
9 of sentences. The sponsor will then have the opportunity to  
10 make any comment upon that. I will then ask Dr. Witten, a  
11 as representative of the FDA, if what we have collectively  
12 stated satisfies the FDA.

13           That being said, we will now embrace the first  
14 question, which will shortly be projected, and, briefly, has  
15 to do with benefits and risks in lap chole and lap Nissen  
16 fundoplication. I will ask for comments of the panelists,  
17 and since we had started introductions with Dr. Burns, I  
18 will ask first Ms. Brinkman for any comments.

19           MS. BRINKMAN: Well, I certainly believe that  
20 these instruments will enhance our ability to do  
21 laparoscopic surgeries and, certainly, laparoscopic  
22 surgeries have demonstrated reduced patient morbidity. It  
23 is a new technology. I think it is exciting, and I feel  
24 that we are moving forward and I am very supportive of  
25 adding the addition of these instruments.

1 DR. WHALEN: Dr. DeMets?

2 DR. DEMETS: Well, I am waffling a bit on this  
3 question because I don't think that the primary question as  
4 stated was really adequately addressed because of the size  
5 of the study. If you have a study with failure rates as low  
6 as these are, there is just no way you can definitively sort  
7 that out with 50 patients an arm, even though it is very  
8 encouraging in terms of the estimates. But you would have  
9 known that in the design phase. I mean, you would have some  
10 sense of what your failure rate is for standard laparoscopy.

11 So you then go to the secondary questions that  
12 were listed, and some make the criterion and some don't.  
13 When you start looking at changes of size of 2 or larger,  
14 not being a surgeon, those sound large to me but they may  
15 not be clinically important, and I haven't heard any  
16 discussion around the panel as to what the minimum delta is.  
17 That is the key for me. Have you met that criteria? I  
18 don't know. I can't judge that surgically. But clearly  
19 from the criteria that have been bounced around, it is  
20 marginal. Some do; some don't, depending on where you draw  
21 the line. So, I am not convinced.

22 DR. WHALEN: Dr. Ferguson?

23 DR. FERGUSON: I respect those concerns and have  
24 those too, but I agree with Ms. Brinkman about the overall  
25 package.

1 DR. WHALEN: Dr. Hannaford?

2 DR. HANNAFORD: I just want to add to the  
3 discussion that when thinking about benefits, which are  
4 mentioned up there, I think there are benefits -- well, I  
5 think we should be able to consider potential benefits, and  
6 benefits that may accrue in the future when some extension  
7 to this is approved.

8 DR. WHALEN: Dr. Galandiuk?

9 DR. GALANDIUK: I think questions number one and  
10 two almost could be combined, and I agree with Dr. DeMets  
11 about the delta being defined. Depending on what variable  
12 they look at, the device clearly isn't identical because it  
13 does take longer. Is that clinically significant? I don't  
14 think so. So, in other words, it may be different but none  
15 of these differences are significant, and I think it has  
16 been shown both to be safe and effectiveness.

17 DR. WHALEN: Dr. Crittenden?

18 DR. CRITTENDEN: I am having a hard time with the  
19 question because I am not sure there has really been a  
20 benefit that has been established, but I do think they are  
21 equivalent technologies. So, as best I can answer the  
22 question, I kind of agree with what Dr. Ferguson and Ms.  
23 Brinkman said.

24 DR. WHALEN: Dr. Anderson?

25 DR. ANDERSON: I agree. I think this is a

1 valuable technology. The four variables were conversion  
2 rate, procedure duration, learning curve and hospital stay.  
3 Of those, the procedure duration of 51 minutes was a non-  
4 significant delta. I think that that is an acceptable  
5 difference, particularly for a new procedure. I was  
6 satisfied with the learning curve data which, because it  
7 looks like other learning curves that we see, such as in  
8 sentinel node mapping, and I think this is very promising.  
9 I am supportive.

10 DR. WHALEN: Dr. Chang?

11 DR. CHANG: In thinking about risks versus  
12 benefits, it appears that the major risks and the untoward  
13 events that occurred during the conduct of this clinical  
14 study was actually due to low tech instruments, such as the  
15 trocar and other instruments not related to the Intuitive  
16 system such as the harmonic scalpel. Those would be another  
17 issue related to just laparoscopic surgery in general and  
18 the risks for patient populations. So, I feel that this  
19 exciting technology is worthwhile in putting on the market.

20 DR. WHALEN: Dr. Talamini?

21 DR. TALAMINI: In short, I think the benefits to  
22 outweigh the risks for the device. I think the increased  
23 length in operative time is statistically significant but  
24 not important clinically, and I think the same in terms of  
25 blood loss. I think the benefits are that some day I will be

1 able to do my whipple using these types of tools.

2 DR. WHALEN: Ms. Dubler?

3 MS. DUBLER: I don't think we have heard  
4 sufficient discussion of benefits to be able to address  
5 number one. In regard to number two, my concern is the  
6 following: That in order to arrive at a positive risk-  
7 benefit ratio not only do we need more data on benefit, but  
8 I fear we need more data on risk, which would be related in  
9 my mind to the dimension of the training that will be  
10 required. People who are not sufficiently trained could, in  
11 fact, pose a risk using this procedure. Therefore, I think  
12 both of these questions have yet to be adequately addressed,  
13 whereas, I think that the notion of equivalence lets us  
14 address the safe and effective issue I don't think it lets  
15 us address the ratio issue.

16 DR. WHALEN: Dr. Walker?

17 DR. WALKER: Let me put in my two cents worth for  
18 what I see the benefit of this product to be. Clearly, the  
19 risks are about the same so we don't need to go into that.  
20 The numerator part, however, is that we are being asked to  
21 evaluate -- if we look at your ethics question, is there an  
22 ethical benefit? And my answer is, yes, there is a very  
23 real ethical benefit because this is an enabling technology.  
24 This is a first step to something that will have really far-  
25 reaching benefits for the patients. I agree that in this



1 particular application no benefit has been shown. But if we  
2 quash and don't approve, then that second step which really  
3 shows the promise will never be allowed to be taken by this  
4 company. So, I think we have a moral obligation, as long as  
5 they are safe and effective, to say go for it and show us  
6 the real benefit in step two now that we have given you  
7 permission in step one.

8 DR. WHALEN: Dr. Burns?

9 DR. BURNS: I pretty much concur with Dr. Walker.  
10 I feel this is an important first step. In an absolute  
11 sense, there perhaps hasn't been an enhanced benefit shown  
12 over standard laparoscopic procedures but, listening to the  
13 panel presentation as well as the investigators, it appears  
14 that it is at least equivalent, or can be equivalent, and  
15 potentially is an important first step into things that  
16 could potentially be much more beneficial in the future.

17 DR. WHALEN: I would summarize that the panel  
18 feels that the data has largely demonstrated, with some  
19 asterisks about its numeric versus clinical significance,  
20 that this is a safe and effective technology and,  
21 furthermore, that the preponderance of opinion would say  
22 that there is a net risk-benefit ratio in its favor, with an  
23 important disclaimer on the ethical side of things on our  
24 ethics expert.

25 Are there any comments in that regard from the

1 sponsor?

2 MR. DANIEL: No.

3 DR. WHALEN: Dr. Witten, with that consensus, has  
4 that successfully answered both questions one and two?

5 DR. WITTEN: Yes, thank you.

6 DR. WHALEN: Thank you. We will go to question  
7 number three which, as many may recall, is the lengthiest of  
8 the five questions, and has to do with some extensions into  
9 other arenas. Staggering again, we will start with Dr.  
10 DeMets.

11 DR. DEMETS: Not being a surgeon, I will pass.

12 DR. WHALEN: Dr. Ferguson?

13 DR. FERGUSON: Given some of the answers I got a  
14 minute ago, I don't think this falls within the purview of  
15 the group to concern itself with in the approval of this  
16 particular device. What is going to happen is that  
17 everybody in the room knows that this is going to be applied  
18 widely across all kinds of surgeries and all disciplines,  
19 and it will be inevitably misapplied, unfortunately, in some  
20 situations. I look on our job here, today, to regulate not  
21 number three but number five. So, I will hold off on that.

22 DR. WHALEN: Dr. Hannaford?

23 DR. HANNAFORD: I also will pass, not being a  
24 clinical person.

25 DR. WHALEN: Dr. Galandiuk, you can't pass because

1 you are a clinical person.

2 [Laughter]

3 DR. GALANDIUK: I don't really think it is  
4 relevant because it is the same issue as if you are doing an  
5 open operation. There is going to be a difference in tissue  
6 strength and fragility of tissues whether you are operating  
7 on a 90-year old woman or a 20-year old ma, and the same  
8 will be true for this technique. So, it is not any  
9 different than conventional surgery.

10 DR. WHALEN: Dr. Crittenden?

11 DR. CRITTENDEN: I agree with this. This is a  
12 purely clinical decision and I think really is beyond the  
13 purview of the panel to really talk about, and I think the  
14 surgeon on the scene has to make their own decision about  
15 this based on the clinical data they have at hand.

16 DR. WHALEN: Dr. Anderson?

17 DR. ANDERSON: I agree. This is a clinical  
18 decision. Bad decisions are made at all levels, not just  
19 with new technology, and this is not unique.

20 DR. WHALEN: Dr. Chang?

21 DR. CHANG: I would ditto. Once this is on the  
22 market, we really are dependent on clinicians' judgment in  
23 the proper use of the instrument.

24 DR. WHALEN: Dr. Talamini?

25 DR. TALAMINI: I agree with Dr. Ferguson and the

1 rest.

2 DR. WHALEN: Dr. Dubler?

3 MS. DUBLER: Pass.

4 DR. WHALEN: Dr. Walker obviously passes. Dr.  
5 Burns?

6 DR. BURNS: I agree in that any potential misuse  
7 can only be guarded against by the appropriate labeling and  
8 good judgment by the surgeons.

9 DR. WHALEN: Ms. Brinkman?

10 MS. BRINKMAN: I agree. Unfortunately, mistakes  
11 will be made but we will learn.

12 DR. WHALEN: So to be very concise, Dr. Witten, I  
13 would say that this falls under the practice of medicine in  
14 the opinion of the panel. Are there any comments by the  
15 sponsor?

16 MR. DANIEL: No.

17 DR. WHALEN: Does that satisfactorily answer the  
18 question?

19 DR. WITTEN: Yes.

20 DR. WHALEN: Thank you. Going to question number  
21 four then, it has to do with the device being labeled or  
22 claimed to be fail-safe. I will start this one with Dr.  
23 Ferguson.

24 DR. FERGUSON: I think this is one of the  
25 strongest aspects of the whole proposal, to me. I think the

1 way in which the company has approached the potential  
2 problems that could occur with complex machinery of this  
3 type is outstanding.

4 DR. WHALEN: I will ask Dr. Hannaford to  
5 disambiguate this issue for us a little bit --

6 [Laughter]

7 DR. HANNAFORD: Well, I generally share that  
8 opinion. My concerns about stability are more like  
9 something that, for the most part, are a potential cause of  
10 a failure, not so much a failure of fail-safe. In other  
11 words, there are these hardware mechanisms that are well  
12 documented in the filing which will catch this kind of  
13 failure, should it occur.

14 In the trial there were three instances of that  
15 happening. I was a little concerned at first when I read  
16 that the threshold for that failure detection was raised in  
17 response to those three failures. The instance in  
18 particular was a current check which, after a certain period  
19 of time, caused a fault if the current exceeded a certain  
20 value, just like a circuit breaker in your house. But, in  
21 detail, the time interval required to trip that safety  
22 feature was set to be extremely short, a few hundred  
23 nanoseconds, as I recall, which is a very, very, very  
24 conservative setting and one that is very likely to be  
25 triggered by non-dangerous, noise type of events. So, it

1 didn't bother me at all that they something like doubled  
2 that time. In fact, that almost is evidence of how well the  
3 safety features seemed to be working. So, I basically  
4 concur.

5 One other slight concern I would have is I wonder  
6 about if there are any effects of hardware failures on  
7 stability of the system. In other words, are there failure  
8 modes which could cause the control system to go unstable?

9 So, I am sort of a one-note person here, but this  
10 is where my expertise comes in. We built a similar system  
11 in the lab, much less sophisticated, of course, but we had  
12 an instance where if a cable was not tight enough -- if the  
13 cable became loose the system could go unstable. So, that  
14 kind of thing should at least be thought about, anyway.

15 DR. WHALEN: Dr. Galandiuk?

16 DR. GALANDIUK: Well, I have to defer to the  
17 experts.

18 DR. WHALEN: Dr. Crittenden?

19 DR. CRITTENDEN: I don't have much to add but just  
20 kind of wonder since this is an enabling technology whether  
21 or not some sort of objective performance criteria ought to  
22 be set in regards to the things that Dr. Hannaford talked  
23 about during his review. That is all I have.

24 DR. WHALEN: Dr. Anderson:

25 DR. ANDERSON: The safety concerns about the

1 software crashing was what I understood to be the major  
2 problem, and I don't put that in the same category as the  
3 device, you know, moving wildly and injuring a structure.  
4 There are lots of reasons that things can slow down in the  
5 OR, including the nurse not having the right equipment for  
6 you or the wrong cart being pulled. This happens all the  
7 time. I think that the safety mechanisms are fine for this.

8 DR. WHALEN: Dr. Chang?

9 DR. CHANG: I think the answer is yes.

10 DR. WHALEN: Dr. Talamini?

11 DR. TALAMINI: I am a little surprised at the word  
12 fail-safe. I have learned through painful experience to  
13 never say "never" or never say "always" in surgery because  
14 you are always being proven wrong. But I think the system,  
15 as I see it, is at least as fail-safe as my tired resident  
16 who is helping me do the case.

17 [Laughter]

18 I would feel better about it if I had used it five  
19 or ten times and felt how you can actually swing things out  
20 of the way quickly, but it seems from the videos and the  
21 data that that capability certainly exists.

22 DR. WHALEN: Dr. Dubler?

23 MS. DUBLER: Pass.

24 DR. WHALEN: Dr. Walker?

25 DR. WALKER: The term fail-safe bothers me as well

1 because it brings up visions of the Titanic and Dr.  
2 Strangelove.

3 [Laughter]

4 And, I would hope that the sponsor doesn't  
5 actually say their device is fail-safe but, rather, points  
6 to the adequate safeguards of it.

7 I share Dr. Hannaford's concern about the hardware  
8 failures and the myriad of possibilities for instability  
9 that that could lead to, but I feel reasonably comfortable  
10 that the existence of the product problem reporting system  
11 to FDA is going to mean that if those problems they can't be  
12 swept under the carpet and they will have to be dealt with.

13 DR. WHALEN: Dr. Burns?

14 DR. BURNS: I have to defer to Dr. Hannaford for  
15 some of his comments, but it would appear, based on what we  
16 have seen, that the system is fairly safe. I would look  
17 forward to the sponsor being able to respond to Dr.  
18 Hannaford's concerns.

19 DR. WHALEN: Ms. Brinkman?

20 MS. BRINKMAN: I just want to emphasize the fact  
21 that we remember to provide good technical training. I  
22 don't understand at all, but I am always sure that there is  
23 some good technician that does.

24 DR. WHALEN: Dr. DeMets?

25 DR. DEMETS: Computers fail, and they fail more



1 often because of software -- glitches in the system, than  
2 they do hardware, at least in my experience. And, I don't  
3 know how you would do it but I am not so confident that the  
4 software is fool-proof. I don't know how this division  
5 tests software, but I suspect that there are holes in that  
6 software that some day some realization will occur that we  
7 haven't anticipated, and goodness knows what the results  
8 would be.

9 DR. WHALEN: Dr. Witten, with the important caveat  
10 echoed by several members that nothing is a hundred percent  
11 and we wouldn't label this as totally fail-safe, I believe I  
12 can represent the consensus of the panel that we do feel  
13 that there is adequate safety built into the system.

14 Any specific questions from the sponsor or any  
15 other comments that someone would like to address?

16 DR. GOODHART: There were a couple of questions,  
17 most of them center around stability and on the issue  
18 stability analysis, we have done complete analysis in  
19 several different frames and, in fact, very much the terms  
20 that you presented in your summary. The FDA should have  
21 those. We are happy to provide additional detail in that  
22 sense. We have observed no instability in the system in any  
23 of its operative modes.

24 You brought up a second issue on stability, which  
25 is failure -- how did things fail. In fact, you brought up

1 a specific example of a cable. On our systems we have  
2 redundant sensing. Redundant sensing is on either side of  
3 cable drives. So, if we do see, for example, the thing that  
4 you mentioned -- a cable go loose -- we can detect it in  
5 advance of tip motion and transition to fail-safe.

6 So, we have done an analysis of how this failure  
7 affects stability as a whole, and that is a consistent set  
8 of performance evaluation tests that we do at every product  
9 release interval.

10 You had mentioned scaling, filtering and things  
11 like that, we do have performance criteria and performance  
12 tests available at the FDA, and in additional detail, if  
13 necessary.

14 On the software front, we do exhaustive testing.  
15 We also rely very heavily on following standard software  
16 procedures. One of the biggest things we want to make sure  
17 we do with our software is exercise every single mode that  
18 the system can get into, and we do that exhaustively. You  
19 had mentioned a specific concern of kind of an unexercised  
20 raffle that you never quite get into, and one of the things  
21 that we do is make sure that we hit every raffle that the  
22 system can head for. So, to the best of our ability, we  
23 addressed that concern.

24 DR. WHALEN: Do any panelists have any follow-up  
25 comments or questions?

1 DR. HANNAFORD: Yes, I will follow-up. I guess  
2 that pretty much satisfies me because my main desire was  
3 that the control design is well documented in the FDA just  
4 because it is really introducing a new thing and something  
5 that I think is going to show up in other devices, other  
6 areas, other companies -- and, just so that it is done  
7 right. That is my main thing. It just isn't in the  
8 material that I saw in advance of the meeting.

9 DR. WHALEN: Dr. Witten, does FDA feel that  
10 question four has been adequately addressed?

11 DR. WITTEN: Yes, thanks.

12 DR. WHALEN: Thank you. The fifth and final  
13 question has to do with training issues, and I will first  
14 ask that Dr. Hannaford make comment.

15 DR. HANNAFORD: My one point on the training --  
16 actually, I have two points on the training. One is that  
17 there is some very recent data from collaborators of mine at  
18 UW on training using some simple simulators, and then  
19 evaluating performance of those students in their further  
20 clinical training with animal surgery and human surgery.  
21 They found that their initial study showed no effect. Then  
22 they realized that the learning curve was, in fact, longer  
23 than they expected. Once they extended the number of  
24 training repetitions out to around nine repetitions, then  
25 they were able to measure, in a statistically valid way, the

1 effects of the training.

2           So, I think the number of operations required is  
3 on the order of ten or more to see that kind of training  
4 effect. Of course, this was a little different. This was  
5 residents or beginning surgeons. The authors of that study  
6 would be Dieter Pohl and Mika Sinanan.

7           The other thing is I think a training program for  
8 this device ought to cover the exceptional cases. So, the  
9 surgeons ought to be explicitly trained about failure modes  
10 or system fail-safe modes -- what the surgeon should do when  
11 X happens, even if that is a very unlikely thing. That is  
12 one of the benefits. In fact, just for the future I see a  
13 tremendous benefit related to training because ultimately  
14 you will unplug the console from the robot and plug it into  
15 a computer, and the surgeon will practice on a computer  
16 simulation. Again, that is not what we are seeing today but  
17 that is where this can go. So, I think there is a lot of  
18 benefit for training ultimately to this technology.

19           DR. WHALEN: Dr. Galandiuk?

20           DR. GALANDIUK: I agree that training has to be  
21 very clearly specified, not only for the operating surgeon  
22 but also for the nurse or technician in the operating room  
23 that will help. I was looking at the user manual here, and  
24 the control panel has five pages of terms and pictures on  
25 it. To me, that looks intimidating, and I think it would be

1 very important to ensure that anyone who purchased this  
2 device in terms of hospital -- whoever was going to be the  
3 technical advice person in the operating room go through all  
4 this and be very familiar with it.

5 Similarly, I think there should be a requirement  
6 in terms of hours that a surgeon trains for this. The  
7 surgeons are all, you know, "we know everything all the  
8 time, you know, just look at this and browse through this  
9 and be comfortable with it," but I think that with this you  
10 should require a certain amount of both didactic as well as  
11 animal lab or simulation training to "certify" surgeons for  
12 this because it still is different from regular laparoscopic  
13 surgery and I think it would be important to ensure safety  
14 on the part of the surgeon doing this rather than just the  
15 device.

16 DR. WHALEN: Dr. Crittenden?

17 DR. CRITTENDEN: I would echo those remarks. I  
18 don't have anything to add.

19 DR. WHALEN: Dr. Anderson?

20 DR. ANDERSON: The training I think is the most  
21 important part of what was discussed today, and I just want  
22 to make the point to the sponsors, echoing what Dr. Talamini  
23 said, with laparoscopic cholecystectomy we were lucky. We  
24 went up to a one percent common duct injury rate but because  
25 the operation was so much better we survived that. But in

1 this setting, this isn't that much better than standard  
2 laparoscopy and I think to protect your device for the  
3 future we really need to make sure that this is adequate,  
4 and very specific non-human practicing needs to be done  
5 because this is really a different technology.

6 DR. WHALEN: Dr. Chang?

7 DR. CHANG: Laparoscopic cholecystectomy is  
8 patient driven. There is no doubt about that. I mean, it  
9 just exploded because patients asked their surgeons for the  
10 mini-procedure, and so I would echo all the comments that  
11 physician education would be key in terms of success in  
12 advancing this instrument.

13 DR. WHALEN: Dr. Talamini?

14 DR. TALAMINI: Training and credentialing were the  
15 sticky wickets of laparoscopic cholecystectomy, and they are  
16 still being wrestled with for laparoscopic procedures. The  
17 issues are complex -- who pays for the training? How are  
18 surgeons going to take time off for the training? It is not  
19 simple and straightforward. Whose responsibility is it?  
20 Who decides when somebody can do this and when not? So, I  
21 don't think it is a simple issue but it is an awfully  
22 important issue, and I think it is important for this device  
23 because this is the first device where the surgeon affecting  
24 tissue manipulation is distant from the patient's tissue.

25 So, I would even propose that we ask the sponsor

1 to development some sort of a formal plan or proposal as  
2 part of the approval process. Being a novice, I don't know  
3 how that works but I think it is an important enough and a  
4 difficult enough issue that it warrants that much attention.

5 DR. WHALEN: Dr. Dubler?

6 MS. DUBLER: I think these are very helpful  
7 comments, and I too would be inclined to approve the  
8 technology subsequent to a training program or certification  
9 program being suggested by the sponsor. I don't know what  
10 analogous situations we might turn to for guidance. Perhaps  
11 there are some; perhaps it really is a matter of first  
12 impression. But, I think if we take the patient in this  
13 complex seriously and the patient is being asked to balance  
14 and choose between different technologies, the training of  
15 the person who is performing and using this new technology  
16 will determine the risk. If patients are to address for  
17 themselves the question of their own values and a risk-  
18 benefit ratio, then I think in fact there has to be a level  
19 of certification and training that becomes the base for  
20 patient choice.

21 So, I think I would like the sponsor to be able to  
22 tell us how they think this could be done, and I would be  
23 inclined to build that into the approval.

24 DR. WHALEN: Dr. Walker?

25 DR. WALKER: I am very uncomfortable with the

1 notion of an FDA regulatory mandate for a certain number of  
2 hours or a certain type of training, and would argue that  
3 that should not be a condition of approval. The reason is  
4 that if there are patient injuries from untrained surgeons  
5 using this device, then the sponsor is going to be the deep-  
6 pocket co-defendant sitting in court with the surgeon, and I  
7 think that probably is going to have the effect of ensuring  
8 adequate training before the device is sold, probably a  
9 better program of adequate training than we, sitting here as  
10 non-training experts today, could possibly come up with.

11 DR. WHALEN: Dr. Burns?

12 DR. BURNS: I agree that training for the surgeon  
13 as well as the surgical team is going to be important for  
14 this. I don't know if I would make that a requisite for  
15 approval, other than to note that it should be worked out  
16 between the sponsor and the FDA upon approval.

17 DR. WHALEN: Ms. Brinkman?

18 MS. BRINKMAN: If I were you, I would develop the  
19 hottest video game and get it in every doctor's lounge, and  
20 I bet you there is a high correlation between people who are  
21 good at video games and people who are good at this kind of  
22 surgery --

23 [Laughter]

24 -- and all of us who don't have that eye-hand  
25 coordination will probably never be very good at it. So,



1 make it fun and they will be doing it just for the fun of it  
2 and they will learn well.

3 DR. WHALEN: Dr. DeMets?

4 DR. DEMETS: Well, as I pointed out in my  
5 comments, I think it does take more than two or three cases  
6 to get the learning curve over with, and the data that has  
7 been submitted show that, I think. How you translate that  
8 into regulatory language, if at all, I don't know.

9 DR. WHALEN: Dr. Ferguson?

10 DR. FERGUSON: I agree basically with what has  
11 been said, particularly with Dr. Burns. I don't think it  
12 ought to be part of our charge to tell them how to train. I  
13 think the company will be responsible for that. I do think  
14 that in that analysis you have to carry it a long way past  
15 just the surgeon who is being trained but think about  
16 resident training and all of the problems that we have had  
17 with VATs technology and trying to teach residents in  
18 hospitals how to do it when there is sort of only one person  
19 who can be doing something on time. So, there are a lot of  
20 those issues that are going to come up for the future.

21 The other comment I would like to make is that,  
22 again looking at the future, I don't know how you do this  
23 but I would like to be assured somehow that the surgeon who  
24 is responsible for the case is in the room, and not removed  
25 to another tele-situation, and also that he was able to get

1 to the patient and get in, in case of a severe emergency.  
2 You are not going to see that in what we have been talking  
3 about today, but you can see it a lot in cardiovascular  
4 procedures. So, looking toward the future there, I think  
5 that it would be worthwhile thinking about those. Those are  
6 not on the table today, of course.

7 DR. WHALEN: Dr. Anderson has one more comment.

8 DR. ANDERSON: The question of should the FDA be  
9 telling this group how they should be teaching, I don't  
10 think that was the proposal that Dr. Talamini suggested.  
11 The question was to ask the group, you tell us how it is  
12 that you would have these people taught. The four surgeons  
13 who did this work know more about what this procedure is  
14 like than anyone else in the world, and they can design for  
15 you, I think, a protocol that looks reasonable. And, I  
16 think it is important to document that, that this is how it  
17 ought to be done, not to say you have to do eight of them  
18 and then you are done. We want to hear from you more than  
19 what was said today about how you would teach them.

20 DR. WHALEN: Dr. Witten, in answering question  
21 number five about the types of training that the panel feels  
22 are warranted, in summary, it should cover any failure modes  
23 that are inherent in the system. It should cover the entire  
24 operative team and not simply the operating surgeon. There  
25 should be didactic portions followed by either inanimate or

1 animal laboratory sessions prior to use in humans. There  
2 should be a formal plan or proposal as a component of the  
3 sponsor's application. And, there is discomfort in  
4 specifying any specifics as to the number of hours that  
5 would be entailed.

6 Are there any comments by the sponsor?

7 DR. MOLL: I just wanted to remark that Intuitive  
8 Surgical takes training very seriously. In fact, we believe  
9 it is one of the keys to both clinical and commercial  
10 success. I would also like to agree with a couple of very  
11 insightful comments that training and the methods of  
12 training will change over time, and we believe one of the  
13 really exciting parts of this technology is not only now it  
14 can impact the practice of surgery clinically but how it can  
15 affect clinical training.

16 MR. DANIEL: On behalf of Intuitive Surgical, I  
17 would like to thank the panel and FDA for a very thorough  
18 and thoughtful analysis. Thank you very much.

19 DR. WHALEN: Dr. Witten, does FDA feel that  
20 question five has been adequately addressed?

21 DR. WITTEN: Yes.

22 DR. WHALEN: Thank you. At this juncture, as  
23 mentioned this morning, we do have a second opportunity for  
24 any public comment that wishes to be made. Is there anyone  
25 in the audience who wishes to make comment at this time?

1 Seeing no hands being raised, we can proceed to summations.

2 Is there any final comment from FDA?

3 DR. WITTEN: No. No, I think we have said  
4 everything we have to say.

5 DR. WHALEN: Thank you. Is there any final  
6 comment from the sponsor?

7 MR. DANIEL: No.

8 **Concluding Panel Deliberations and Vote**

9 DR. WHALEN: Thank you. We will then proceed to  
10 voting, and I would remind everyone that the industry and  
11 consumer representatives do not participate in the voting.  
12 I will only vote in the case of a tie. Dr. Krause will read  
13 the voting instructions for the panel.

14 DR. KRAUSE: Everybody here in the audience and  
15 the panel, we all get to be the first ones for these new  
16 voting instructions. They used to be two pages; they are  
17 now one. So, we get to try it out here today.

18 The Medical Device Amendments to the Federal Food,  
19 Drug and Cosmetic Act, as amended by the Safe Medical  
20 Devices Act of 1990, allows the Food and Drug Administration  
21 to obtain a recommendation from an expert advisory panel on  
22 designated medical device premarket approval applications  
23 that are filed with the agency. The PMA must stand on its  
24 own merits, and your recommendation must be supported by  
25 safety and effectiveness data in the application or by

1 applicable publicly available information.

2           Safety is defined in the Act as reasonable  
3 assurance, based on valid scientific evidence that the  
4 probable benefits to health, under conditions on intended  
5 use, outweigh any probable risks.

6           Effectiveness is defined as reasonable assurance  
7 that in a significant portion of the population the use of  
8 the device for its intended uses and conditions of use, when  
9 labeled, will provide clinically significant results.

10           Your recommendation options for the vote are as  
11 follows: Option number one, approval -- you can vote  
12 approval with no conditions attached.

13           Option number two, approvable with conditions --  
14 the panel may recommend that the PMA be found approvable  
15 subject to specified conditions, such as physician or  
16 patient education, labeling changes or a further analysis of  
17 existing data. Prior to voting, all of the conditions  
18 should be discussed by the panel.

19           Number three, not approvable -- the panel may  
20 recommend that the PMA is not approvable if the data do not  
21 provide a reasonable assurance that the device is safe, or  
22 if a reasonable assurance has not been given that a device  
23 is effective under the conditions of use prescribed,  
24 recommended or suggested in the proposed labeling.

25           Following the voting, the chair will ask each

1 panel member to present a brief statement outlining the  
2 reasons for their vote.

3 DR. WHALEN: Thank you, Dr. Krause. The chair  
4 will entertain a motion. Dr. Anderson?

5 DR. ANDERSON: I have a motion that conditions for  
6 approval, if one were to vote for conditions, would be that  
7 the sponsors supply a training description -- what is  
8 considered standard training by your group.

9 DR. WHALEN: The motion has been made to approve  
10 with conditions as specified. Is there a second?

11 [Seconded]

12 The motion has been made and seconded. Since we  
13 have brought up conditions, we can now discuss those  
14 conditions before proceeding to a vote. Please raise your  
15 hand if you wish to discuss.

16 DR. CRITTENDEN: I would like to add to the  
17 motion, and that goes with product labeling, that we ought  
18 to label this as being an equivalent product, not one that  
19 has clinical benefit, and that the clip applicator was not  
20 clinically tested, and that the surgeon ought to be  
21 appropriately gowned.

22 DR. WHALEN: Restated then as approval with  
23 conditions, the conditions would be that the manufacturer  
24 have a training program created; that it be labeled as an  
25 equivalent program; that the clip applicator not be within the

1 approval as it has not been tested; and that the surgeon be  
2 sterile and gowned at the time of the operation of the  
3 instrument. Further discussion of those conditions in the  
4 approval?

5 DR. TALAMINI: Could you just state that again for  
6 me?

7 DR. WHALEN: Probably not but I will give it a  
8 whirl. It is to approve with conditions. The conditions  
9 include the outline of a training program as proposed by Dr.  
10 Anderson; that there be labeling of it as an equivalent  
11 product; that the clip applier not be a part of the approved  
12 instruments; and that the surgeon be sterile and gowned when  
13 operating at the console.

14 DR. CHANG: Is this an amendment to Dr. Anderson's  
15 original motion?

16 DR. WHALEN: Well, if we are within Roberts Rules  
17 of Order we would have to regard it as such, so it would  
18 require a second. Is there a second to the further  
19 additional conditions of approval?

20 DR. HANNAFORD: Second.

21 DR. WHALEN: It has been made and seconded. It is  
22 open for discussion.

23 DR. TALAMINI: I feel as if we may need to discuss  
24 this training. I think we are asking for what would be a  
25 best-case training protocol on the part of the sponsor.

1 Would you agree with that, Dr. Anderson?

2 DR. ANDERSON: Yes.

3 DR. WHALEN: Further comments?

4 DR. HANNAFORD: Yes, can you explain what a best-  
5 case training protocol is to me? My questions about  
6 training are is it within the scope of the FDA's actions --  
7 Can the FDA require a certain amount of training before this  
8 device can be used, or could one option be that the company  
9 can design a training program but that any user is required  
10 to follow that training program? What are the regulatory  
11 options for training?

12 DR. WITTEN: Maybe I can comment on that and if I  
13 am wrong maybe, Jim, you can jump in. We can certainly  
14 require that a sponsor design and provide a training  
15 program. But, as Jim already said, we don't regulate the  
16 practice of medicine so I think it would be unusual for us  
17 to require specific training for the individual user. That  
18 is, to have the training available is one thing but to  
19 require the training from the user would be, I think, beyond  
20 the scope of what we do.

21 DR. KRAUSE: Can I just read the approvable with  
22 conditions again? A panel may recommend that the PMA be  
23 found approvable subject to specified conditions, such as  
24 physician or patient education. So, it is definitely within  
25 the purview of the panel to recommend physician education.



1 MS. DUBLER: I think it is important to state that  
2 even though the FDA doesn't regulate practice, it is very  
3 powerful legally and morally for the company to say this is  
4 the program that the team has to undergo before the team is  
5 qualified to use this technology. So, I think that is very  
6 important and I hope the company would be comfortable  
7 developing it, given the fact that they are the only ones  
8 that have the expertise.

9 DR. WALKER: When we say with this addition or  
10 amendment that it be labeled as an equivalent device,  
11 legally what does that mean?

12 DR. WITTEN: Well, maybe I will give a more  
13 general answer, which is when we look at the panel  
14 recommendations we interpret what the panel is recommending  
15 in terms of what we can do. So, that might translate into,  
16 for example, providing the results of a clinical study on  
17 the label, which we would typically do in any case for a  
18 PMA. In describing the results in the label, you know, we  
19 might make the point that we think could be demonstrated by  
20 the study.

21 DR. FERGUSON: I have a little problem with that  
22 aspect. I thought we were headed on a track to talk about  
23 training but the words "equivalent device," to me, could be  
24 misinterpreted. There is no other device like this that I  
25 know of. So, when we say "equivalent device" what does that

1 mean? That means to us that they didn't prove superiority  
2 over the standard procedure. Is that not correct?

3 DR. WHALEN: Dr. Crittenden, would you answer?

4 DR. CRITTENDEN: Yes, it is not better than the  
5 conventional technique. So, I guess my particular concern  
6 is just that someone will market themselves as being a  
7 robotics surgeon and, hence, that is going to be better  
8 surgery, where it has not been demonstrated to be better but  
9 certainly equivalent.

10 DR. FERGUSON: I couldn't agree with that more,  
11 but when you say "equivalent device" I think that could be  
12 misconstrued perhaps.

13 DR. WHALEN: Other comments or questions? Dr.  
14 Galandiuk?

15 DR. GALANDIUK: I think we are putting too many  
16 conditions on this. I think the physician education is  
17 important, but whether or not the surgeon is gowned and  
18 sterile when he is operating with these controls I think is  
19 clinical judgment and should be left to the surgeon's  
20 discretion. Again, I think the mention of equivalency --  
21 equivalence means different things to different people. So,  
22 I don't think I would make that a condition of approval.

23 DR. WHALEN: If my Roberts Rules are correct, we  
24 will first vote on the amendment, which are those additional  
25 stipulations which Dr. Crittenden has brought up, before we

1 will vote on the original recommendation which has the  
2 training within it. Further comments? Dr. Burns?

3 DR. BURNS: In regards to Dr. Crittenden's comment  
4 that it would be equivalent but not better, that is already  
5 implied in the existing label in the sense that it would be  
6 approved for laparoscopic surgery, period. It is not saying  
7 anything about it being better than any other procedure;  
8 just that it would be approved for that type of surgical  
9 procedure. So, I think it is implicit in the label as it  
10 stands.

11 DR. CRITTENDEN: I think that fair-minded  
12 individuals would agree with you but, given the marketplace  
13 and other incentives, I just wonder if there is some role  
14 for mischief.

15 DR. WHALEN: The chair would suggest that if  
16 somebody is going to put an ad in the "Boston Globe" and  
17 suggest that there is something or other, they are going to  
18 do it, not matter what the FDA and the sponsor does.

19 MS. BRINKMAN: Could not the dressing of the  
20 physician, to be sterile, be included in the education  
21 portion of it rather than just added on? It seems like if  
22 you get a good, comprehensive education package -- and my  
23 assumption is that it is to the company's benefit to do that  
24 anyhow, they could put that in as part of the educational  
25 process and say why.

1 DR. WHALEN: Further question or comment?

2 [No response]

3 We are first voting then upon the Crittenden  
4 amendment to the proposal, which has the three sub-elements  
5 within it of equivalence, clip applier, sterilely gowned and  
6 gloved.

7 All those in favor, please raise their hands.

8 [One hand raised]

9 Dr. Crittenden votes for. All those opposed,  
10 please raise their hands.

11 [Show of hands]

12 Dr. DeMets, Dr. Ferguson, Dr. Galandiuk, Dr.  
13 Hannaford, Dr. Anderson, Dr. Chang, Dr. Talamini, Dr. Dubler  
14 and Dr. Walker. The amendment is defeated. Is there any  
15 subcomponent of that amendment that wishes to be re-proposed  
16 as an amendment before we go to the original question? Dr.  
17 Walker?

18 DR. WALKER: I reintroduce the clip applier  
19 because that has not fully been shown to have been tested.

20 DR. WHALEN: Is there a second to that amendment?

21 [Seconded]

22 DR. WALKER: But perhaps that can be done between  
23 FDA and the manufacturer rather than involving the panel.

24 DR. WHALEN: The amendment has been made and  
25 seconded. Perhaps I jumped in prematurely, but we did make

1 and second it. Is there any further discussion of that?  
2 Dr. Dubler?

3 MS. DUBLER: I don't understand what the  
4 implications of the amendment versus the FDA doing it are.

5 DR. WHALEN: Well, as I understand it, we are  
6 approving a list of instruments. On that list of  
7 instruments was a clip applier. When the motion was  
8 originally amended it was pointed out that there had not  
9 been demonstration and testing of this, and that is why it  
10 was suggested that it be subtracted from the list. And,  
11 that I believe -- I don't want to speak for you, Dr. Walker,  
12 is why it has been reintroduced and seconded.

13 DR. CHANG: Could we have clarification from the  
14 FDA? That was my question earlier to the sponsor, that the  
15 clip was not tested. So, I asked if it were to be developed  
16 in the future would that be submitted to FDA as an amendment  
17 to their PMA.

18 DR. WITTEN: Any instruments that aren't approved  
19 as part of this PMA would need to be submitted as  
20 supplements to the PMA.

21 DR. CHANG: So, there is a protocol to be followed  
22 if the sponsor wants to introduce a clip applier.

23 DR. WITTEN: If there is anything that is not part  
24 of this approval, then they would need to supply the  
25 information separately.

1 DR. FERGUSON: Does that include the harmonic  
2 scalpel?

3 DR. WITTEN: Yes.

4 DR. WHALEN: Any further questions or comments to  
5 be made? Dr. Dubler?

6 MS. DUBLER: Is it appropriate to ask the sponsor  
7 a question?

8 DR. WHALEN: I don't know whether the protocol  
9 dictates it but it is fine with me. So go for it.

10 MS. DUBLER: Why is this on the list with other  
11 things that were tested?

12 DR. WHALEN: Any of the sponsor wish to answer?

13 MR. DANIEL: As I thought we indicated, but let me  
14 clarify, the clip applier that we had available was a small  
15 clip applier, not something that a general laparoscopic  
16 surgeon would generally use. We did use our small clip  
17 three or four times to clip a cystic artery.

18 DR. WHALEN: Seeing no further questions or  
19 comments, the panel is now asked to vote upon the amendment  
20 which is to subtract from approval with conditions approval  
21 of the clip applier. Will those in favor of subtracting  
22 that clip applier please signify by raising their hand and  
23 leaving their hand raised?

24 [Show of hands]

25 DR. WHALEN: That vote is unanimous. Is there any

1 further question or comment upon the motion to approve with  
2 conditions of training? Seeing no further question or  
3 comment --

4 DR. HANNAFORD: Question, sorry. Just the word  
5 "training?" What is the exact amendment or condition that  
6 we are voting on?

7 DR. WHALEN: Dr. Anderson, would you care to  
8 restate?

9 DR. ANDERSON: The proposal is that the condition  
10 for the sponsor is that they provide a detailed description  
11 of the training program which they recommend and endorse for  
12 the use of their product.

13 DR. WHALEN: Is there any further question or  
14 comment?

15 [No response]

16 The panel is then asked to vote upon the motion of  
17 approval with conditions, the conditions being that the  
18 sponsor provide a detailed training program which they  
19 endorse, and also, as amended, that the clip applier is not  
20 a part of the approval. Would those in favor of that please  
21 signify by raising their hands?

22 [Show of hands]

23 It is not unanimous so I will read the names: Dr.  
24 Ferguson, Dr. Hannaford, Dr. Galandiuk, Dr. Crittenden, Dr.  
25 Anderson, Dr. Chang, Dr. Talamini, Dr. Dubler, and Dr.

1 Walker.

2 Will those not in favor, those oppose, please  
3 raise their hands?

4 [One hand raised]

5 That is Dr. DeMets. It is approved. The  
6 recommendation of the panel is that the premarket approval  
7 application for Intuitive Surgical endoscopic instrument  
8 control system and Intuitive Surgical endoscopic instruments  
9 from Intuitive Surgical Inc. be recommended for approval  
10 with conditions, those conditions being that the sponsor  
11 will provide a detailed training program which they endorse,  
12 and that it not be labeled for use with the clip applier.

13 DR. WITTEN: Maybe you were about to do this but  
14 don't we need to go around and have them state --

15 DR. WHALEN: That is what I am about to do, yes,  
16 ma'am. With that recommendation being made, if each of the  
17 panel members who have voted please indicate why they have  
18 so voted. Just for giggles at this last part of the  
19 afternoon, we will go to Dr. Walker first.

20 DR. WALKER: I voted as I did because it seems to  
21 be a safe, well designed product and I am convinced that the  
22 sponsor has done its homework adequately.

23 DR. WHALEN: Dr. Dubler?

24 MS. DUBLER: I voted as I did because I think the  
25 safety and effectiveness of the product has been



1 demonstrated and because I am not permitted, according to  
2 the rules of the FDA, to take into account that it will add  
3 substantially to the cost of surgery.

4 DR. WHALEN: Dr. Talamini?

5 DR. TALAMINI: I voted affirmatively because I  
6 believe it to be safe and effective based on the studies  
7 done and the video case histories that were provided to me,  
8 but believe that the training issues are important.

9 DR. WHALEN: Dr. Chang?

10 DR. CHANG: I voted yes because safety and  
11 effectiveness was demonstrated. The statistical differences  
12 between the two study groups were not clinically  
13 significant. Again, as I know the sponsor has said,  
14 teaching and training for safe use is a primary concern of  
15 theirs?

16 DR. WHALEN: Dr. Anderson?

17 DR. ANDERSON: I voted yes because I think this is  
18 a great product, and I have nothing to add beyond what the  
19 other panelists said.

20 DR. WHALEN: Dr. Crittenden?

21 DR. CRITTENDEN: I voted in the affirmative  
22 because I think this is a safe and effective product in  
23 well-trained individuals.

24 DR. WHALEN: Dr. Galandiuk?

25 DR. GALANDIUK: I also believe it is safe and

1 effective, and the differences that were shown were not  
2 clinically significant, and agree with Dr. Talamini that I  
3 think this will add greatly to the agility of the surgeon  
4 performing laparoscopic procedures.

5 DR. WHALEN: Dr. Hannaford?

6 DR. HANNAFORD: I voted yes because I thought they  
7 showed it was safe and effective, and also demonstrated a  
8 lot of potential for future enhancements to surgery.

9 DR. WHALEN: Dr. Ferguson?

10 DR. FERGUSON: I voted yes for the same reasons,  
11 and I particularly like the concept of what is coming down  
12 the road.

13 DR. WHALEN: So, to keep them in their pews until  
14 the last hymn is sung, we want to go to you, Dr. DeMets, and  
15 find out why you voted against.

16 DR. DEMETS: Despite my enthusiasm and interest in  
17 the product, I don't feel that the data presented met the  
18 criteria which were set out -- safety perhaps; efficacy not  
19 even close; equivalence not always close. So, I think the  
20 study on its primary endpoint was way under-powered, as was  
21 alluded to by the presentation, and the secondary endpoints  
22 were mixed. So, despite my interest and enthusiasm, I don't  
23 feel I can support it as effective or even equivalent.

24 DR. WHALEN: I would like to thank everyone who  
25 presented today, and especially thank the panel members for

1 all of their efforts. The meeting is adjourned.


2 [Whereupon, at 4:40 p.m. the panel adjourned.]

3 - - -

4

## *C E R T I F I C A T E*

I, ALICE TOIGO, the Official Court Reporter for Miller Reporting Company, Inc., hereby certify that I recorded the foregoing proceedings; that the proceedings have been reduced to typewriting by me, or under my direction and that the foregoing transcript is a correct and accurate record of the proceedings to the best of my knowledge, ability and belief.

A handwritten signature in cursive script, reading "Alice Toigo", is written over a horizontal line.

ALICE TOIGO