

1 another, you just tell us that and we'll make sure we act  
2 appropriately.

3 DR. DRUMMOND: It was more for completeness since  
4 you have silicate glass ionomers in phosphates.

5 DR. TYLEND: Okay.

6 DR. DRUMMOND: They are two other cements.

7 DR. TYLEND: We can add those to the list, if you  
8 want to discuss the appropriate priority, polycarboxylate  
9 and resin cements.

10 DR. DRUMMOND: We would probably suggest that  
11 polycarboxylate be low and resin cements be high.

12 CHAIRMAN ROBERTSON: I see. Polycarboxylate, like  
13 zinc silicate cements, would be low.

14 DR. DUNCANSON: I wanted to ask the Panel about  
15 the nomenclature of zinc silicate cement. Are we talking  
16 about zinc silicophosphate cement, or are we talking about  
17 the aluminosilicate mixture with phosphoric acid?

18 CHAIRMAN ROBERTSON: Somebody from FDA?

19 DR. TYLEND: Actually, Dr. Duncanson, why don't  
20 you tell us what you would like? We can subdivide that into  
21 two.

22 DR. DUNCANSON: All right. Since I believe--and  
23 my colleagues can correct me, but I believe that the zinc  
24 silicophosphate cement is not a big issue as far as a

1 marketable product, it's off the market, and likewise  
2 silicate cement has fallen into disuse a great deal also, so  
3 I'm just questioning whether that needs really to be  
4 included since we are talking about including the  
5 polycarboxylate and the glass ionomer cements which are in  
6 use.

7 CHAIRMAN ROBERTSON: Well, if they are no longer  
8 in use--and I will listen to what advice you have there--  
9 then it would not seem reasonable to keep them on the  
10 ingredient labeling list. It seems like we ought to move  
11 them over to the list of things that is not required for  
12 labeling. But I will listen to whatever Panel advice there  
13 is.

14 DR. TYLEND: Consider that we will still accept  
15 submissions for those products, and when a submission comes  
16 in, we will act appropriately as to which list it is on.

17 CHAIRMAN ROBERTSON: I guess the question is: If  
18 somebody for some reason did submit such a cement or a  
19 variant of that cement for use, which list would you want it  
20 on? And given the principle that Deb articulated, from a  
21 perspective of an outsider here, since it would have long-  
22 term contact with mucosa, it would seem that it needs to be  
23 on this list and not the other, irrespective of whether it's  
24 in use or not.

1 DR. DUNCANSON: Well, then I would have two  
2 categories, the zinc silicophosphate cement and the silicate  
3 cement, and I would put them in the low category.

4 CHAIRMAN ROBERTSON: Speak up.

5 DR. NORMAN: I don't see that a silicate cement  
6 ought to be categorized because they are available in  
7 offices that bought them about 20 years ago. But I don't  
8 know that they'd be used. I don't know that they're  
9 manufactured anymore. The silicate cement powder is very  
10 similar to what the glass ionomer cement powder is.

11 CHAIRMAN ROBERTSON: Well, then, I think that is  
12 information that FDA needs to use to refine their list.

13 MR. ULATOWSKI: I was looking at the list. At  
14 first I thought that we had just simply listed the products  
15 out by classification name in the regulations, but it looks  
16 like as the subcommittee had discussed these issues, it had  
17 begun to split things out and divvy things up. So I think  
18 if we have things missing or you'd like to sub-categorize  
19 further, it's open to whatever you choose.

20 CHAIRMAN ROBERTSON: Well, I think we'll then  
21 follow your advice, split them out and leave to FDA the  
22 ability to eliminate one simply since it no longer is  
23 available.

24 All right. Are there any other changes then to

1 the filling materials?

2 DR. TYLEND: As far as the non-Eugenol-containing  
3 temporary filling materials, you're not interested at this  
4 point in making a recommendation for a priority. Is that  
5 correct?

6 CHAIRMAN ROBERTSON: I certainly am not. I think  
7 the subcommittee didn't consider it. I think we just need  
8 to put it aside.

9 DR. DUNCANSON: Well, these cements, I believe,  
10 have chemical homologues of Eugenol, and if some of them  
11 were to come on--well, it's also important to know what the  
12 tissue response would be to those homologues as well as to  
13 Eugenol. And that may not be known at this point in time,  
14 so perhaps zinc oxide Eugenol and non-Eugenol temporary  
15 filling materials could be a category.

16 CHAIRMAN ROBERTSON: The point here is, however,  
17 that the principle, as I understand it, driving these  
18 decisions is primarily contact with oral tissues. And those  
19 materials that have long-term contact with oral tissues got  
20 included on this list. Is that right, Deb?

21 DR. GREENSPAN: Yes.

22 CHAIRMAN ROBERTSON: And a temporary filling  
23 material would have reasonably long-term contact with oral  
24 tissues and, therefore, would be included on the list. And

1 FDA may choose here to either leave it as it is, zinc oxide  
2 Eugenol and other temporary filling materials, or split it  
3 out into two categories.

4 Good. Other comments about filling materials?

5 [No response.]

6 CHAIRMAN ROBERTSON: Deb, shall we move on?

7 DR. GREENSPAN: Moving on to crown and bridge  
8 alloys, resins and materials, we looked at the five that are  
9 listed, and starting with gold-based alloys, that was given  
10 a medium priority; other precious metal alloys, medium  
11 priority; base metal alloys, high priority; temporary crown  
12 and bridge resin, high priority; and ceramics, low priority.

13 CHAIRMAN ROBERTSON: Comments from the Panel?  
14 Changes? Comfortable?

15 [No response.]

16 CHAIRMAN ROBERTSON: Full and partial denture  
17 materials.

18 DR. GREENSPAN: We looked at those materials that  
19 were used with full and partial dentures, so we included  
20 full and partial dentures in one group. Precision  
21 attachments, low; preformed clasp, low; hydrophilic resin  
22 coatings, high; denture adhesives, low. And then the  
23 following OTC preparations: denture cleansers, high;  
24 denture pads and cushions, high; denture reliners, high; and

1 denture repair kits, high.

2           Preformed gold denture teeth, low; preformed  
3 plastic denture teeth, low; partially fabricated denture  
4 kit, low; relining, repairing, or rebasing resins, including  
5 denture acrylic, high; porcelain teeth, low.

6           CHAIRMAN ROBERTSON: Any additions, changes,  
7 comments?

8           [No response.]

9           CHAIRMAN ROBERTSON: Endodontic materials.

10          DR. GREENSPAN: Root canal post, low; root canal  
11 filling resin, high; silver points, low; Gutta percha, low;  
12 stabilizing splint, high.

13          CHAIRMAN ROBERTSON: Stabilizing splint?

14          DR. DUNCANSON: Because of the resin--

15          CHAIRMAN ROBERTSON: Wire and acrylic?

16          DR. GREENSPAN: Yes.

17          CHAIRMAN ROBERTSON: Okay. Comments, concerns,  
18 changes?

19          [No response.]

20          CHAIRMAN ROBERTSON: Orthodontic materials.

21          DR. GREENSPAN: Do we want to add as a category  
22 here endodontic sealers?

23          CHAIRMAN ROBERTSON: We can recommend to FDA that  
24 they add endodontic sealers. And its priority? High.

1 Based on your notions of zinc oxide Eugenol before. Good.

2 DR. GREENSPAN: Orthodontic materials. Bracket  
3 adhesive resin and tooth conditioner, high; band material,  
4 low; metal brackets, low; plastic brackets, low; spring  
5 tube, expansion screws, and wires, low; ceramic bracket,  
6 low.

7 CHAIRMAN ROBERTSON: Changes, additions, concerns,  
8 amendments?

9 [No response.]

10 DR. GREENSPAN: We then come to a group of items  
11 that were put under the category other materials. They were  
12 all considered because of length of time that they would be  
13 in contact with human tissues or else the possibility of  
14 damage and, therefore, contact with human tissue.

15 The first one was the elastomeric impression  
16 material, which was high; oral cavity polishing agent, low;  
17 abrasive disk, low; abrasive point, low; a caries detection  
18 device that comes into contact with the mucosa, low; resin  
19 impression tray material, medium; intraoral ligature and  
20 wire lock, high; extraoral dental headgear, low; teething  
21 ring, low; ultraviolet light for polymerization, low; bite  
22 registration material, high; rubber dam, high; silicate  
23 protector, high; and intraoral dental wax, low.

24 CHAIRMAN ROBERTSON: I guess I have to ask a

1 question about intraoral ligature and wire lock. Lou, help.

2 DR. SINGLETON: That is a maxillofacial type of  
3 device that would be used to bring two bony sections into  
4 apposition and fixate them. If one of our oral surgeons is  
5 here, is Sue Runner here, by any chance?

6 CHAIRMAN ROBERTSON: Well, Willie is here. Does  
7 that make sense to you?

8 DR. STEPHENS: Is that already covered in the  
9 first section in intraosseous fixation?

10 DR. GREENSPAN: I think it was included in other  
11 materials because of the way that it was listed. So if it  
12 is felt that that is more appropriately moved--although it  
13 strictly speaking isn't an implant, I think that's why we  
14 didn't leave it under implants and, rather, why it was left  
15 under other materials. The rationale for giving it a high  
16 priority was based on what we were doing with the implant  
17 materials.

18 DR. SINGLETON: It might be more appropriate to  
19 list it under the implant section, actually.

20 DR. GREENSPAN: All right.

21 CHAIRMAN ROBERTSON: Good. Any other question?  
22 Dr. Patters?

23 DR. PATTERS: Yes, the ultraviolet light for  
24 polymerization, it's not clear to me why we would not want



1 to know the materials that it's constructed from, or we  
2 don't want to know the materials that a saliva ejector is  
3 constructed from, that being in contact much longer.

4 CHAIRMAN ROBERTSON: Yes.

5 DR. NORMAN: I don't think the light would ever--  
6 should ever be in contact with tissue, and the part that  
7 might come in contact with tissue would be the stainless  
8 steel covering. I don't think that it would have any  
9 bearing on any tissues.

10 DR. PATTERS: Is there a reason to recommend  
11 labeling for such?

12 CHAIRMAN ROBERTSON: The members of the  
13 subcommittee need to dredge back into their memories.

14 DR. GREENSPAN: Into their collective memories,  
15 yes.

16 DR. DRUMMOND: Actually, I have another point.  
17 The one, if we're really going to talk about it, should be a  
18 visible light, because ultraviolet light is not used that  
19 much anymore, anyway.

20 CHAIRMAN ROBERTSON: Right.

21 DR. DRUMMOND: So the issue is really visible  
22 lights.

23 DR. GREENSPAN: I don't remember why it was there  
24 without looking back at some of the discussions that came

1 up, or whether there was the possibility for mucosal damage  
2 if it were misused. It may be one of the reasons.

3 DR. PATTERS: How does that differ, for instance,  
4 from a facebow?

5 DR. GREENSPAN: Yes, if it's misuse.

6 DR. PATTERS: A facebow is not recommended for  
7 labeling.

8 DR. GREENSPAN: Right.

9 CHAIRMAN ROBERTSON: And it's not an ingredient--

10 DR. GREENSPAN: No.

11 CHAIRMAN ROBERTSON: --in the strict sense, that  
12 is causing that damage. Unless the emissions are considered  
13 an ingredient.

14 DR. SINGLETON: And the emissions portion of the  
15 device would be covered in another part of the regulation,  
16 anyway.

17 CHAIRMAN ROBERTSON: So the emissions would not be  
18 considered an ingredient.

19 DR. SINGLETON: Would not be.

20 CHAIRMAN ROBERTSON: Under those conditions, I  
21 guess I agree that it should not be here.

22 DR. GREENSPAN: I have no problem with removing  
23 it.

24 DR. DRUMMOND: None.

1 DR. NORMAN: None.

2 CHAIRMAN ROBERTSON: Consensus. Actually, we'll  
3 get it off the list, but a note to FDA that ultraviolet as  
4 the designator here may not be appropriate anymore, and  
5 curing light or whatever, a more appropriate term might be  
6 better.

7 MR. HLAVINKA: We have the ultraviolet under a 21  
8 CFR regulation, but no one uses them anymore. Everybody  
9 uses filters for light, so everything is vis now, but we  
10 just don't have a regulation for it. So it might be an  
11 oxymoron, we find these things substantially equivalent to  
12 UV, even though there's none.

13 I just noticed here, could we return to full and  
14 partial denture material?

15 CHAIRMAN ROBERTSON: Sure.

16 MR. HLAVINKA: On the partially fabricated denture  
17 kit, and it is in low. Let me read you our regulation for  
18 that. A partially fabricated denture kit is a device  
19 composed of connected preformed teeth that is intended for  
20 use in construction of a denture. A denture base is  
21 constructed using the patient's mouth as a mold by partially  
22 polymerizing the resin denture base materials while the  
23 materials are in contact with the oral tissues. After the  
24 denture base is constructed, the connective preformed teeth

1 are chemically bonded to the base. Because of the term  
2 partial polymerization within the patient's mouth, instead  
3 of complete.

4 DR. GREENSPAN: Can the panel look at that again?  
5 I think that we might want to reconsider changing that  
6 priority to high because of the procedures that are carried  
7 on intraorally.

8 CHAIRMAN ROBERTSON: One hopes this is silicate  
9 revisited, but I think we used the same rules. Is there  
10 consensus that we change the priority to high, given what  
11 was just read to us?

12 DR. SINGLETON: Well, one thing you might want to  
13 consider also is that this particular device does not  
14 include, from what I can gather, relining, repairing, and  
15 rebasing resin. That would be separate and apart from this  
16 device. In other words, you'd have to incorporate that into  
17 the partial denture. And the relining, repairing, and  
18 rebasing resin already is in high priority.

19 DR. GREENSPAN: Yes.

20 DR. SINGLETON: So you might want to consider  
21 that.

22 DR. NORMAN: If the kit comes with its resin, with  
23 its--

24 DR. GREENSPAN: Yes.

1 DR. NORMAN: --auto-polymerizing resin system,  
2 then it ought to be in a high category.

3 CHAIRMAN ROBERTSON: And everybody is agreed?  
4 Good. We'll change that.

5 Any other comments on the other materials? Are we  
6 comfortable with those given the general principles that  
7 we're trying to follow here?

8 [No response.]

9 DR. PATTERS: Could I ask that the Panel just take  
10 one more look at denture adhesives and see if they're  
11 comfortable with low priority?

12 CHAIRMAN ROBERTSON: Where, Mark?

13 DR. PATTERS: Full and partial denture materials,  
14 fourth item, denture adhesives. In my mind, these are in  
15 contact with mucosa for a relatively long period of time; in  
16 some cases, misused could be all the time. I don't know  
17 what the incidence of sensitivity to these agents are.

18 DR. NORMAN: They're reasonably high molecular  
19 weight polymers, and most of them are--I can't think of any  
20 that are not this type of structure right now. The allergic  
21 response to this would be very, very low. Irritation can  
22 occur, but that's not due to the adhesive per se.

23 DR. PATTERS: If the Panel is comfortable, I'm  
24 comfortable.

1           CHAIRMAN ROBERTSON: Good. So I don't hear any  
2 suggestion to change that from its present low priority. Is  
3 that correct?

4           [No response.]

5           CHAIRMAN ROBERTSON: Good. Any other changes,  
6 additions?

7           [No response.]

8           CHAIRMAN ROBERTSON: Good. Then--

9           DR. GREENSPAN: Do we want to look at the things  
10 we didn't consider?

11           CHAIRMAN ROBERTSON: No. I think what we might--  
12 you mean because we might move some back to this list? Good  
13 point. Then we'll move on to dental devices not recommended  
14 for ingredient labeling. Maybe I'll just read them quickly  
15 one at a time, and the question everybody asks themselves  
16 is: Given the principles articulated by the subcommittee,  
17 are there any here that we should consider moving to our  
18 list of dental devices recommended for ingredient labeling?

19           Paper saliva ejector, cotton roll; dental  
20 amalgamator, AC powered; dental handpiece and accessories;  
21 gingival fluid measurer--by your silence here--jump in if  
22 there's one that you're concerned about--pulp tester;  
23 electrode gel for pulp tester; extraoral source X-ray  
24 system; intraoral source X-ray system; dental X-ray exposure

1 alignment device; lead-lined X-ray position indicator;  
2 dental X-ray film holder; mercury and alloy dispenser;  
3 dental amalgam capsule; resin applicator; articulate;  
4 facebow; dental bur, and there's a note that with regard to  
5 surgical burs, the subcommittee wants to seek to establish  
6 consistency of policy with the Orthopedic Devices Branch.  
7 Any update on that?

8 [No response.]

9 CHAIRMAN ROBERTSON: No.

10 Mechanical denture cleaner; pantograph; bone  
11 cutting instruments and accessories; powered bone drill;  
12 gas-powered jet injector; spring-powered jet injector;  
13 dental hand instruments; fiber optic dental light; dental  
14 operating unit; dental injecting needle; rotary scaler;  
15 ultrasonic scaler; dental electrosurgical device; airbrush;  
16 anesthetic warmer; articulate paper, dental chair; rubber  
17 dam accessories.

18 Where did rubber dam in itself go?

19 DR. NORMAN: It went in high because of latex.

20 CHAIRMAN ROBERTSON: Okay, good.

21 Paper points; dental floss--

22 MR. ULATOWSKI: Mr. Chairman?

23 CHAIRMAN ROBERTSON: Yes.

24 MR. ULATOWSKI: In regard to dental floss, I know

1 there was some discussion in the past. Just a caveat on  
2 dental floss. From time to time, we receive applications  
3 for dental floss with fluoride and perhaps in the future  
4 other therapeutic agents. We have particular concerns about  
5 those ingredients and labeling for therapeutic ingredients  
6 and particular concerns about clinical studies and other  
7 aspects.

8 CHAIRMAN ROBERTSON: Yes, but not the dental  
9 floss, which is--

10 MR. ULATOWSKI: Per se, right. Not the base  
11 floss.

12 CHAIRMAN ROBERTSON: The point here--right.

13 DR. GREENSPAN: I seem to remember that in our  
14 discussions those claims would be--there is another  
15 mechanisms for dealing with those claims and that those  
16 products would be considered in that appropriate category  
17 rather than dental floss just on its own.

18 MR. HLAVINKA: That's correct.

19 CHAIRMAN ROBERTSON: Dental floss was the vehicle.

20 MR. HLAVINKA: There was a recent publication in  
21 the 21 CFR whereby the agency redefined dental floss. It  
22 has to be totally inert. Any other additive, it's no longer  
23 a 510(k).

24 CHAIRMAN ROBERTSON: And that's the context in



1 which this dental floss appears.

2           Massaging pick; boiling water sterilizer--I  
3 remember those; I remember silicate, too. Endodontic dry  
4 heat sterilizer; manual toothbrush.

5           DR. GREENSPAN: Do you remember that?

6           [Laughter.]

7           CHAIRMAN ROBERTSON: Good. And now we've added  
8 the ultraviolet light. Good.

9           DR. FRAZIER: Mr. Chairman, could we just for a  
10 moment go back to denture adhesives?

11          CHAIRMAN ROBERTSON: Sure.

12          DR. FRAZIER: Under full and partial denture  
13 materials. I have just been thinking about the length of  
14 time that these adhesives could be in someone's mouth, and I  
15 would just like to hear a more full discussion about the low  
16 versus medium or high priority for labeling.

17          DR. NORMAN: I don't have anything to add. I  
18 think they're low priorities. I do not see them to be--

19          CHAIRMAN ROBERTSON: Could you get the mike?

20          DR. NORMAN: Yes.

21          DR. FRAZIER: What do they contain?

22          DR. NORMAN: As far as I know, all of them contain  
23 a petrolatum base, and FDA may want to talk about this,  
24 because petrolatum is now being looked at with a different

1 viewpoint. But the actual material that is used to hold the  
2 denture in place is hydrophilic resin. Most of them were  
3 initially naturally occurring gums. Those have been  
4 replaced. I don't know that there's anything on the market  
5 at the present time that contains a natural gum.

6 CHAIRMAN ROBERTSON: Is there anybody in the  
7 audience who can satisfy the concerns of our public  
8 representative about the long-term exposure of mucosa to  
9 dental adhesives?

10 Well, we've done the best we can.

11 DR. FRAZIER: What is petrolatum?

12 DR. SINGLETON: Maybe I could just add something.  
13 I think in this particular case perhaps you ought to  
14 consider the possibility of the reaction to ingestion, not  
15 just to reaction of local tissue. I think long-term  
16 ingestion of the material may be a factor.

17 CHAIRMAN ROBERTSON: My sense is that the  
18 expertise herein assembled is insufficient to answer the  
19 question that you addressed, with the exception that it has  
20 been noted in the literature--

21 DR. NORMAN: There are people who put in as much  
22 as a gram of material at a particular time in an  
23 application. The amount that in an ill-fitting denture may  
24 be sufficient would be half that. So the great majority of

1 the material that you're ingesting would be petrolatum, and  
2 the other's a resin that's nonreactive from the information  
3 I at least have been able to get. I have looked at a lot of  
4 these over the years, and I don't recall any of the  
5 biological information supplied to me that shows any  
6 toxicity.

7 CHAIRMAN ROBERTSON: Okay.

8 DR. GREENSPAN: I would raise just one thought for  
9 this. The fact that it's--perhaps we don't have sufficient  
10 data, but also should the possibility for abuse of this  
11 product concern us, the possibility for overuse or misuse?  
12 Or is that already taken care of in the way that the product  
13 is packaged? Should we move this from a low to a medium  
14 priority, or do we feel that there are no data to suggest  
15 that we should do that?

16 CHAIRMAN ROBERTSON: My only observation, which I  
17 probably should not make, is that I'm quite proud of the  
18 notion that the need for such devices decreases every year.

19 DR. GREENSPAN: I know we felt when we first  
20 looked at this that because of the type of material, we felt  
21 that the priority was low.

22 DR. FRAZIER: But at the same time, we have been  
23 talking about the principle that we were trying to use of  
24 being length of contact rather than--

1 DR. GREENSPAN: And materials that produce mucosal  
2 response, if it's mucosal contact. You know, plastic teeth  
3 are in contact for a long time, but we gave them a low  
4 priority.

5 DR. BOUWSMA: I'm not aware of any of the safety  
6 issues associated with the denture adhesives, but it's one  
7 thing that I can provide perhaps at the next meeting. We do  
8 accumulate information like that, and I can bring an answer  
9 to that question.

10 DR. FRAZIER: The reason I brought it up is simply  
11 because it looks so strange on a listing to have all of the  
12 other OTC denture-related products be high, but it's low. I  
13 don't have any particular knowledge that it should be medium  
14 or high; it is just that it looks kind of strange.

15 DR. BOUWSMA: I think that was based on the safety  
16 concerns more so with the others rather than with the  
17 adhesive that's...

18 MR. HLAVINKA: The adhesive is furnished in a  
19 finished configuration while all the other ones are--you  
20 make them yourself, so to speak. It is in a kit, all those  
21 other OTCs.

22 CHAIRMAN ROBERTSON: Other concerns?

23 [No response.]

24 CHAIRMAN ROBERTSON: Then may I have a motion to

1 recommend these dental devices to the FDA for use in  
2 whatever approach FDA chooses or chooses not to take in  
3 device labeling? This is a recommended list of materials  
4 based in general on the principle of contact with oral  
5 tissues in categories of implant materials, in filling  
6 materials with the modification of the definition of  
7 cements, in crown and bridge alloy resins and materials,  
8 full and partial denture materials, with a change of  
9 partially fabricated denture kits priority to high,  
10 endodontic materials with the additional of endodontic  
11 sealers with a priority of high, orthodontic materials and  
12 other materials, intraoral ligature and wire lock, it's  
13 moved to the category of implants, and ultraviolet light for  
14 polymerization moves off the list.

15 DR. DRUMMOND: I'll make the motion.

16 CHAIRMAN ROBERTSON: Moved.

17 DR. PATTERS: Second.

18 CHAIRMAN ROBERTSON: And seconded by Dr. Patters.

19 Now, any discussion? And in the discussion, I  
20 think it needs to be said that we're responding to a request  
21 by FDA for a list based on some kind of principle. This is  
22 in no way meant to either approve or disapprove the notion  
23 of labeling.

24 All in favor of the motion, please raise your

1 hand.

2 [A show of hands.]

3 CHAIRMAN ROBERTSON: Three, four, five. There is  
4 unanimous support for the motion, and the motion carries--  
5 with considerable excitement.

6 [Laughter.]

7 CHAIRMAN ROBERTSON: We now will adjourn until--  
8 captain?

9 DR. TYLEND: 1:45.

10 CHAIRMAN ROBERTSON: Until 1:45. Okay, that gives  
11 us an hour and 15 minutes.

12 [Whereupon a luncheon recess was taken to  
13 reconvene at 1:45 p.m., this same day.]

## 1 A F T E R N O O N S E S S I O N

2 DR. ROBERTSON: Welcome back. We will now move to  
3 a discussion of a Guidance Document for Dental Handpieces.  
4 We will start the open public hearing with Ms. Dawn Johnson  
5 from Midwest. Welcome.

6 MS. JOHNSON: Thank you. I am just going to pass  
7 out a copy of what I am saying.

8 I am Dawn Johnson. I am with Midwest Dental  
9 Products. We are a division of Dentsply and we are the  
10 largest U.S. manufacturer of dental handpieces, high and low  
11 speed.

12 We would first like to thank the FDA for becoming  
13 more involved in handpiece standards. There is a lot of  
14 confusion and misinformation around handpieces and the  
15 proper use and care of them. As the market share leader, we  
16 find we are often caught in the middle of a lot of this  
17 information. We welcome some attention to the handpiece  
18 product category and hope that FDA initiatives will help to  
19 reduce misinformation and improve the credibility of those  
20 companies making an effort to sell products that are safe,  
21 effective, and with claims that are honestly communicated to  
22 the dental profession.

23 We also appreciate the opportunity to share our  
24 issues regarding the proposed guidelines. We have several

1 concerns and questions relative to the current proposal and  
2 will share them in the hopes that these concerns will be  
3 considered in drafting the final guidelines.

4           There are four areas we will cover. They are the  
5 product performance standards, sterilization efficacy  
6 testing and assurances, labeling requirements, and the  
7 impact of sterilizers and accessory products, and also  
8 after-market remanufacture of handpieces have on handpiece  
9 safety.

10           We understand that the FDA's primary  
11 responsibility is to protect the consumer regarding the  
12 safety of products within its jurisdiction. It is our  
13 understanding that performance standards are required for  
14 Class II and Class III medical devices but not for Class I  
15 devices. Given the handpiece is a Class I device, our first  
16 concern with the proposed standard is with the performance  
17 requirements which are included in the document. We  
18 understand and concur with the FDA review of safety-related  
19 performance issues but we are concerned with review of  
20 product performance issues completely unrelated to either  
21 the product's safety or its ability to perform its intended  
22 function.

23           Based on 68 years of handpiece design and  
24 manufacturing, we believe that the primary safety concerns



1 relative to handpieces are bur retention, cooling of the  
2 tooth surface, overheating at the cap, use for unintended  
3 procedures, and potential for disease transmission.

4           The following are specific performance inclusions  
5 to the proposed guideline which we would like reconsidered.  
6 These are inclusions which we believe relate only to  
7 performance and not to performance relative to product  
8 safety.

9           The first is the use of the ISO as a part of the  
10 guideline, and I have several points on the ISO. One is,  
11 ISO standards are designed to harmonize testing and device  
12 interface compatibility, not safety. Included in the ISO  
13 are many performance criteria which do not provide value  
14 relative to product safety.

15           The ISO is a voluntary compliance standard; it is  
16 not mandatory. Depending on which criteria of an ISO  
17 standard, many of the handpieces used and sold in the U.S.  
18 may not be in compliance with the ISO. These product  
19 differences versus the ISO are market-driven preferences  
20 versus safety-related product issues.

21           The ISO standard is a dynamic document subject to  
22 change. In fact, standard 7785-1, which is referred to on  
23 page three of the proposed guideline, is a proposed standard  
24 which has not yet been approved. The U.S. can and may

1 choose not to adopt future proposals, but the way the  
2 current guideline or proposed guideline is written, the  
3 product may be challenged on non-conformance to future and  
4 unknown standards that would be passed on future ISOs.

5           We request that the FDA select those portions of  
6 the ISO which relate to safety issues and include them in  
7 the guideline but exclude non-safety-related data. We also  
8 request that the FDA either insert relevant text into the  
9 guideline or specify the date of the ISO revision to be  
10 referenced in the guideline to ensure that we are not forced  
11 to comply with future, unknown changes to the ISO standard.

12           Second on performance is proposed light output and  
13 light measurement at submission. This requirement appears  
14 to be performance versus safety driven. We are not aware of  
15 any safety issue related to the amount of light projected by  
16 the handpiece. Over half of the handpieces currently  
17 purchased and used in the U.S. do not even offer a light  
18 feature. Many dental schools teach dentistry using non-  
19 lighted handpieces and dentists continue to choose this  
20 avenue out of both familiarity and also to reduce costs.

21           The proposed light output standard references  
22 10,000 lux as required light output, but we are unaware of  
23 any study which defines light output relative to consumer  
24 safety or that would support 10,000 lux being a minimal

1 requirement to ensure safe use.

2           We believe light output measurement is not a  
3 safety issue and we would like to see it eliminated from the  
4 final handpiece guideline. If the FDA determines that there  
5 is a related safety issue and is determined to keep it in  
6 the guideline, we request clarification on light measurement  
7 relative to how and where the light should be measured.  
8 This is to ensure consistency among manufacturers.

9           Two other issues on performance, one is the  
10 handpiece angle of visibility, which is included as a  
11 criteria. Again, we are unaware of any safety issue  
12 relative to the stated angle of measurement, and we are also  
13 unaware of any study which would suggest a safe versus  
14 unsafe angle, what that angle would be, so we do not see the  
15 purpose of this requirement and suggest it be deleted from  
16 the final guideline.

17           There is an additional performance issue about 12-  
18 ounce force on a push-button cap. We agree that a  
19 manufacturer should ensure lack of heat generation at the  
20 cap of the handpiece but we believe the 12-ounce button  
21 force is one avenue to that means and the guideline should  
22 request test submission on heat generation data, not impose  
23 a design tolerance to be followed.

24           On the subject of sterilization testing and

1 assurances, we believe the testing requirement for a  
2 manufacturer to confirm the sterilizability of a handpiece  
3 is a prudent requirement for the FDA and for the  
4 manufacturers. We would like to clarify several issues,  
5 though, around validation protocols.

6 First, we would recommend that the guideline  
7 require three runs of ten handpieces versus three runs of  
8 three, only in that gives us statistical significance and we  
9 would be more comfortable with that data.

10 We need clarification as a handpiece manufacturer  
11 on the definition of permissible load.

12 Then finally, as a handpiece manufacturer, we will  
13 be able to test a limited number and types of sterilizers.  
14 We would choose standard, conventional sterilization  
15 devices. We would hope that any new sterilizer entering the  
16 market as a device which can sterilize handpieces would be  
17 required to perform the same sterilization protocols on  
18 handpieces as the handpiece manufacturer is asked to  
19 perform. We are just not sure where the liability is here  
20 between the two companies or between the two products. We  
21 want to make sure that our testing is mated at some point so  
22 that we are all saying the same thing and testing to the  
23 same degree.

24 On reprocessing and labeling requirements, we

1 believe submission of proof that our product can be rendered  
2 sterile under a given set of sterilization parameters is an  
3 important safety issue and should be a requirement to obtain  
4 a 510(k). We believe that the number of reprocessing cycles  
5 through sterilization is an economic and performance issue,  
6 not a safety issue. We would like the panel to consider the  
7 following.

8           First of all, from an effectiveness standpoint,  
9 the required number of reprocessing cycles is primarily an  
10 economic issue, not a performance or a safety issue. Should  
11 a manufacturer choose to develop a \$50 handpiece which could  
12 be reprocessed 100 times and then thrown away, the dentist  
13 may very well be pleased with the performance because the  
14 economics would be superior to higher-priced units  
15 reprocessed through more cycles. For instance, it would be  
16 economically better than a handpiece that he purchased for  
17 \$500 that gave him 600 cycles. So we would like the FDA to  
18 consider that reprocessing is primarily economic.

19           Secondly, sterilization may accelerate handpiece  
20 failure, but the mechanical degradation and the mode of  
21 failure is not different than it was pre-sterilization. The  
22 current concern with reprocessing information required in  
23 both the validation and labeling sections of the proposed  
24 guideline suggests that the failure mode is changed and is

1 unsafe as a result of the reprocessing.

2           The proposed guideline requires all manufacturers  
3 to publish their reprocessing information on the package  
4 label. We do not believe that an accurate or clinically-  
5 relevant measurement of reprocessing life expectancy is  
6 attainable and believe an attempt to communicate such  
7 information via product labeling is extremely misleading to  
8 the dentist.

9           The labeling of reprocessing cycles suggests that  
10 heat sterilization is the sole input to handpiece life  
11 expectancy and handpiece failure. In fact, there are  
12 significant variations in handpiece life expectancy based on  
13 other variables. These variables include but are not  
14 limited to the type of practice, whether it is a pedodontist  
15 or a crown and bridge practice; the degree of pressure  
16 applied by the individual practitioner; the staff  
17 maintenance procedures; the shank quality, shape, and length  
18 of burs and diamonds used in the handpiece, and the delivery  
19 unit air pressure. Each of these variables are outside of  
20 the handpiece manufacturer's control and are beyond any  
21 measurement via the testing of reprocessing cycles.

22           We strongly suggest that emphasis on the number of  
23 reprocessing cycles be reconsidered and excluded from the  
24 final guideline. We also believe including laboratory

1 testing data in our instructions which claims an expected  
2 number of reprocessing cycles is extremely misleading and  
3 will misrepresent the product to the dentist. We do not  
4 believe this should be required on the labeling. We do  
5 believe if a manufacturer chooses to make a reprocessing  
6 claim, that claim should be supported by test data and the  
7 testing required for such a claim should not be accelerated.

8           Finally, if the FDA decides to include  
9 reprocessing considerations in the final guideline and  
10 allows accelerated testing of reprocessing cycles, a  
11 definite limit should be set on how much acceleration is  
12 allowed to ensure parity testing among manufacturers.

13           One other thing that is a little bit outside of  
14 this but we would like to bring up, in the interest of  
15 safety, we believe the FDA should have equal if not more  
16 concern around the remanufacture of handpieces as it does  
17 around new products. We estimate that approximately 50  
18 percent of U.S. handpiece repairs are done by other than the  
19 manufacturer using non-manufacturer parts, processes, and  
20 specifications. We also believe the handpiece repair market  
21 is as large or larger a market than the new product market.

22           This means that roughly half of the handpieces in  
23 use today and those that will be in use in the future are  
24 remanufactured in some way and will not be positively

1 affected relative to either the safety, sterilization  
2 assurance, or labeling as a result of new guidelines. The  
3 primary reason to source remanufactured parts is economic.

4 To highlight the extent of this after-market  
5 activity, we believe there are multiple hundreds of repair  
6 sources available to the dentist that use parts, processes,  
7 and tooling which are not the manufacturer's, are not to  
8 manufacturer specification, and are not in compliance with  
9 GMPs.

10 Independent handpiece repair is so common that one  
11 organization has sold handpiece franchises throughout the  
12 country. These franchises cost \$20,000 and provide tools,  
13 parts, and three-day training to franchisees to repair  
14 handpieces. A device is currently on the market, being sold  
15 to dentists, which allows the dentist to repair his own  
16 handpiece by removing the failed bearings and pressing on  
17 new ones.

18 Going back to the subject of handpiece safety, the  
19 handpiece is generally designed to have the bearing be the  
20 failure mode. Bearing failure is generally a safe failure  
21 mode. As people remanufacture handpieces, they are  
22 depending on other handpiece components to continue to  
23 perform safely through other bearing lives and bearing  
24 failures. We believe this can be an unsafe condition,



1 moving the failure mode to a more dangerous mode, such as  
2 chuck failure and possible bur ejection.

3           We request that the FDA do anything within its  
4 authority to regulate the handpiece after-market. While we  
5 support a higher FDA interest in handpiece products, we also  
6 recognize that this adds expense to the handpiece  
7 manufacturer. The after-market repair source's competitive  
8 advantage is low cost. Further regulating the handpiece  
9 manufacturer gives the non-regulated after-market companies  
10 an even greater competitive advantage over the manufacturer.  
11 This will result in even larger numbers of handpieces being  
12 used on patients that are remanufactured, have not proved  
13 sterilizability, have not had FDA review, and are not  
14 manufactured in accordance with GMP standards.

15           Given that about one-half of the handpiece in use  
16 are remanufactured by non-manufacturers, we believe  
17 increased regulation of the handpiece manufacturer without  
18 regulation of the after-market puts the manufacturer at a  
19 disadvantage and at the same time fails to address what we  
20 believe is one of the product category's most significant  
21 safety issues.

22           In addition to the after-market issue, we have  
23 also seen a myriad of lubrication products, sterilizers,  
24 sterilization aids, et cetera marketed with claims that they

1 will improve handpiece life. Some of these products have  
2 been registered with the FDA. We have asked companies for  
3 testing to support these claims but they have not been able  
4 to provide them to us. Information, photos, and other  
5 material included in their literature clearly misrepresent  
6 the product.

7           Some of the claims that we have seen in print  
8 include sterilizers which extend handpiece life four to five  
9 times, lubricants which extend handpiece life ten times, et  
10 cetera.

11           As we conduct our sterilizability testing, we will  
12 be using our own maintenance products and cannot verify  
13 sterilizability when used with other than our own tested  
14 accessories. We request that the FDA require companies  
15 marketing accessory products making claims relative to  
16 handpiece sterilization to follow the same guidelines as the  
17 handpiece manufacturer follows with regard to  
18 sterilizability and infection control claims.

19           That summarizes our input to the new guidelines  
20 for handpieces. We thank the FDA for its interest in the  
21 safety of handpieces and hope that some of the  
22 considerations we have presented are of value to the FDA and  
23 will be comprehended in the final guideline. Thank you.

24           DR. ROBERTSON: Thank you very much. That was a

1 fascinating presentation. It's a more complicated issue  
2 than I thought.

3 I have a quick question. I probably didn't get  
4 all of them, but in terms of the performance standards with  
5 respect to safety that you outlined, I got retention of the  
6 bur and cooling of the tooth and heat generation of the cap  
7 and microbial transmission. There may have been a fifth.

8 MS. JOHNSON: The fifth was use for unintended  
9 procedures.

10 DR. ROBERTSON: Okay.

11 MS. JOHNSON: Such as an air handpiece for oral  
12 surgery or something like that.

13 DR. ROBERTSON: My question is, to what extent  
14 does the sterilization of that handpiece, continued  
15 sterilization, to what extent does it affect those safety  
16 issues? I mean, you said that it was not possible for you  
17 to estimate--

18 MS. JOHNSON: I understand. I guess our  
19 experience, or to the best of our knowledge, and I have two  
20 counterparts here, as well, who might be better qualified to  
21 answer that, but we don't believe that handpiece  
22 sterilization affects safety in any way other than positive,  
23 by eliminating the potential for disease transmission.

24 DR. ROBERTSON: So that--

1 MS. JOHNSON: We do believe that handpiece  
2 sterilization decreases on average, for the average, if you  
3 could find him, practitioner, the life of the handpiece  
4 between repairs. So we believe that the performance  
5 characteristics and the safety issues are unchanged, but the  
6 time line of actual use of the handpiece is generally  
7 shortened.

8 DR. ROBERTSON: I understand that, but generally  
9 shortened, meaning it doesn't work anymore, or generally  
10 shortened meaning there has been an effect because of the  
11 sterilization on retention of the bur or cooling or heat  
12 generation of the tip?

13 MS. JOHNSON: No, the general failure is a bearing  
14 failure, which was the same failure as anyone would  
15 generally have experienced before. I would ask, we have an  
16 engineering and a regulatory, but I would say that we have  
17 seen no increase in any safety issues such as bur retention,  
18 cooling of the tooth surface, use for unintended procedures,  
19 or overheating of the cap as a result of handpiece  
20 sterilization. What we see is a quicker, a shorter time  
21 line for the bearing to fail, and a bearing failure is a  
22 slowdown of the speed of the handpiece that ultimately gets  
23 to a point where the handpiece doesn't cut but it's not  
24 unsafe.

1 DR. ROBERTSON: So that a number which got pasted  
2 on the outside of the handpiece which gave you the number of  
3 cycles would not, in fact, relate to the safety issues?

4 MS. JOHNSON: It's an economic issue. The other  
5 thing with the numbering is that we don't know how to get to  
6 that number, because if you are a pedodontist, the number is  
7 different than if you are a prosthodontist. And if you are  
8 a heavy-handed dentist, the number is different than if you  
9 are a light-handed dentist, and there are significant  
10 differences in each of those areas that we point out. They  
11 are not minor differences, they are significant.

12 DR. ROBERTSON: As a dean, I have the sense that I  
13 probably have a laboratory which will give you the lower  
14 limit.

15 MS. JOHNSON: It will give me one or the other,  
16 I'm not sure.

17 DR. ROBERTSON: I mean, I have no doubt that my  
18 students will go through the handpiece as fast as any of  
19 those other examples you gave, and I could probably  
20 determine for you the lower limit. Thank you.

21 MS. JOHNSON: Any other questions?

22 DR. TYLEND: Dawn, you gave some handouts but  
23 they didn't reach this far. Do you have any more?

24 MS. JOHNSON: I think I do. I'll walk them up

1 right after. Thank you.

2 DR. ROBERTSON: Wait, wait, any other questions  
3 from the panel?

4 DR. STEPHENS: Does your company and other major  
5 manufacturers do your own remanufacturing of handpieces?

6 MS. JOHNSON: Yes, yes, and I'm being careful for  
7 this reason. There are certain products we choose to  
8 remanufacture, but most, we choose to put in new components.

9 DR. STEPHENS: I see.

10 MS. JOHNSON: You know, to put in a brand new  
11 assembly. But I believe that all major handpiece  
12 manufacturers selling in the U.S. have their own repair  
13 service, and their judgment on when to use a remanufactured  
14 part or a new part varies by product and by company. But  
15 service is always available from the handpiece manufacturer.  
16 It's just more expensive.

17 DR. STEPHENS: Thank you.

18 DR. DRUMMOND: This may not be to you. I just  
19 have a question of clarification in terms of the FDA's role  
20 in this guidance document. Is this document to be for just  
21 safety or safety and effectiveness, and would effectiveness  
22 include performance standards?

23 MR. ULATOWSKI: The answer is, it pertains to both  
24 aspects. In our evaluation of equivalence, we take a look

1 at equivalent performance in terms of safety and  
2 effectiveness.

3 DR. ROBERTSON: I think probably the presentations  
4 by FDA will maybe try to clarify that and it will be an  
5 opportunity for the panel in this new area to see if we  
6 can't figure out where we are.

7 Any other questions from the panel? Dr. Patters?

8 DR. PATTERS: Your company then has no problem  
9 with certification of sterilization and with the drop test,  
10 is that correct?

11 MS. JOHNSON: The certification of sterilization,  
12 we have no problem with and we believe it's in both of our  
13 interests that we do that testing.

14 The drop test, I guess I understood to be related  
15 to disposable handpieces. Am I misinterpreting?

16 DR. TYLEND: That's correct.

17 MS. JOHNSON: Is that correct? We don't make a  
18 disposable handpiece. Our own irrelevant opinion, since we  
19 don't make one, is if you drop it and it can't be  
20 sterilized, you throw it away, so to drop test it doesn't  
21 make any sense because you ought to toss it. But because  
22 that's not our product category, I shouldn't have said that,  
23 but that's our opinion. Thank you very much.

24 DR. ROBERTSON: Thank you very much.

1           Mr. Tom Fise from the American Dental Trade  
2 Association?

3           MR. FISE: I did not want to confuse this with the  
4 prior statement, which is why I delayed handing this out.

5           DR. ROBERTSON: It doesn't guarantee that we won't  
6 be confused.

7           [Pause.]

8           MR. FISE: Good afternoon. I'm Tom Fise and  
9 appear, again, on behalf of the American Dental Trade  
10 Association as their Special Counsel on Regulatory Affairs.

11           The concept of providing guidance documents to  
12 help manufacturers understand the issues of keenest interest  
13 to FDA in its review of certain product submissions is  
14 inherently a good one so long as the purpose of the guidance  
15 document is not misunderstood. We have some questions  
16 relating both to the specific content of this guidance  
17 document as well as to how the document is used by FDA.  
18 This guidance document should not be construed as a catalog  
19 or checklist of FDA requirements which must be met to secure  
20 510(k) approval for handpieces.

21           Dental handpieces and their accessories are so-  
22 called pre-amendment devices, meaning that they were legally  
23 marketed before the 1976 Medical Device Amendments. Dental  
24 handpieces on the market at that time, and including any



1 that have filed one or more subsequent 510(k) notifications  
2 of modifications that have been accepted by FDA, are legally  
3 marketed today without respect to whether they meet all of  
4 the criteria established in the FDA guidance document.

5           The prospective criteria for having an acceptable  
6 510(k) filing for a new or modified device continues to be  
7 whether the device can demonstrate substantial equivalence  
8 to a device that is currently legally marketed, and that  
9 case may be a pre-1976 device or handpiece. So the idea of  
10 an FDA guideline is not a substitute for that legal  
11 criteria, which is the limitation that allows someone to  
12 come onto the market.

13           A 510(k) submission that does not meet each and  
14 every criterion set out in the 11-page FDA guidance document  
15 may still meet the legal requirements for substantial  
16 equivalency, and so our first point is that it is important  
17 to place the FDA guidance document in the proper context.  
18 Complete compliance with the guidance document may  
19 constitute a "safe harbor" to acceptability, but compliance  
20 with the guidance document is not a sine qua non. It is not  
21 the sole pathway to demonstrate substantial equivalency.

22           We would also like to highlight a few specific  
23 items in the FDA guidance document that we believe may be  
24 troublesome, and to some extent, a couple of these have been

1 touched upon previously.

2           We have noted the sections there, II.B.9. We do  
3 not think that the manufacturer should be required to list  
4 all accessories or attachments that may be used with the  
5 device, as this may include products manufactured by others.  
6 A current listing of accessories or attachments that the  
7 specific manufacturer may intend to market with that device  
8 seems justified.

9           With respect to II.D.1, we are uncertain with  
10 regard to the requirement that hose connections remain  
11 intact up to 150 percent of normal operating air pressure.  
12 This tolerance level may relate to a voluntary  
13 specification, and probably does, to the ISO or ADA/ANSI  
14 specs, but there is not a regulatory obligation for  
15 handpieces as Class I devices to meet such a voluntary  
16 standard.

17           With respect to II.D.5, as was mentioned before,  
18 we are a bit uncertain about the source of the regulatory  
19 status of a requirement for a minimum load of 12 ounces on  
20 the push button before the release mechanism contacts any  
21 rotating parts.

22           On II.D.6, we are also uncertain of the source or  
23 validity of the three-foot drop test requirement. Again, we  
24 were not clear until this point that that was limited to

1 disposables, and so that may be either an error in our  
2 reading or might be something that could have some further  
3 clarification. Moreover, the requirement for 100 percent  
4 pass rate, we had some question about. It says there be a  
5 minimum of ten and that all must pass.

6           With respect to II.E.1, we think this illustrates  
7 the difference between guidelines and the legal requirements  
8 for substantial equivalency, and we note two specific  
9 provisions under this section. The first would require  
10 identification of a specific number of reprocessing cycles,  
11 and I want to come back to the term reprocessing cycles.

12           The second would state that any advice to follow  
13 the recommendations of the sterilizer manufacturer is  
14 inadequate. Now, it is our understanding that, again, in  
15 theory, a device can be legally marketed today without  
16 either identifying a number of reprocessing cycles and that  
17 it could today bear the legend simply to follow the  
18 recommendations of the sterilizer manufacturer. So if that  
19 device can be legally marketed, there certainly should be  
20 able to be a substantially equivalent device that has the  
21 same information.

22           Clearly, if a manufacturer chooses to state a  
23 number of reprocessing cycles or cleaning cycles or  
24 sterilization cycles, the product must meet that criteria.

1 And likewise, if the manufacturer does state any  
2 instructions beyond following the sterilizer manufacturer's  
3 instructions, then the instructions they give must be  
4 adequate. However, we don't think that this guideline  
5 document can operate to change the rules, as it were, or  
6 change the requirements for substantial equivalency.

7           Moving then to II.E.2.a--I'm sorry to be a little  
8 bit convoluted--we do not see the basis on which FDA would  
9 maintain with respect to 510(k) demonstrations of  
10 substantial equivalency that the manufacturer be required to  
11 maintain "a record of the data that show that the handpiece  
12 can withstand the number of reprocessing cycles claimed in  
13 labeling."

14           Likewise, we are uncertain about the source or  
15 validity and the regulatory basis of the requirement or  
16 criteria "with less than ten percent decrease in performance  
17 characteristics" that is listed in II.L. Actually, we had  
18 some trouble finding Section II.L in the document itself, at  
19 least as we received it in advance of the meeting.

20           Finally, I want to get back to the reprocessing  
21 issue and say that, again, we question this requirement or  
22 prospective requirement that manufacturers be required to  
23 state a number of use or reprocessing cycles that the  
24 handpiece can withstand before disposal or repair is

1 required. Particularly, we are troubled by the use of the  
2 term "reprocessing," when what we seem to be referring to is  
3 sterilization cycles.

4           In the FDA's recent publication on good  
5 manufacturing practices, the FDA has defined reprocessing  
6 this way. Reprocessing means all or part of a manufacturing  
7 operation which is intended to correct non-conformance in a  
8 component or finished device before distribution. So we  
9 think the term reprocessing has a specific meaning that  
10 relates to something that happens before the device is sold  
11 or redistributed or whatever and does not refer to what  
12 happens in the dental office during sterilization.

13           So we think a correction just in that term is  
14 important, and it is important because the term reprocessing  
15 is carried on throughout these regulations, so something  
16 that is a little clearer, we would be happier with.

17           In conclusion, we are concerned about potential  
18 misunderstandings as well as the potential that the document  
19 in its current form might be misapplied, either by  
20 manufacturers or by evaluation or compliance staff at FDA.  
21 Many of the recommendations in the document, we think, are  
22 very valuable.

23           We would, however, recommend that because of  
24 possible misunderstandings that the best thing would be to

1 add some explanatory language and indicate with absolute  
2 clarity that the document defines a safe harbor and  
3 establishes targets, that it is not to establish the minimum  
4 requirements either for a 510(k) substantial equivalency or  
5 for compliance requirements of current products that are on  
6 the market. If that cannot be done, then perhaps the  
7 document ought to be withdrawn and looked at again. We do  
8 believe that the specific areas we have highlighted would  
9 merit some further review and consideration.

10 We appreciate again the chance to present on  
11 behalf of ADTA and would be happy to answer any questions.

12 DR. ROBERTSON: Thank you very much.

13 Questions from the panel?

14 [No response.]

15 DR. ROBERTSON: That was very useful. Thank you.  
16 I am sorry?

17 DR. TYLEND: I just want to tell Dawn that we  
18 were able to get some copies of your presentation made, so  
19 you don't have to worry about rooting through your briefcase  
20 to find them.

21 DR. ROBERTSON: Mr. Jeffrey Peinhardt from Den-  
22 Tal-Ez?

23 MR. PEINHARDT: Good afternoon. My name is  
24 Jeffrey Peinhardt from Star Dental, which is actually a

1 subsidiary of Den-Tal-Ez, Incorporated.

2 I would like to first thank the FDA for generating  
3 this draft document. It has been long overdue and needed  
4 and we basically agree that we do need a draft document. I  
5 appreciate the panel allowing us to review changes that we  
6 feel should be necessary.

7 Dr. Tylanda has distributed copies to the panel  
8 dated July 28 of our discussions. We have some areas of  
9 concern and I would like to go over these in particular.

10 Sterilization validation, Section II.E.1, the  
11 draft implies the handpiece manufacturers must test every  
12 available method of sterilizer, and when combined with other  
13 sections, print individual instructions for each.  
14 Additionally, if one model differs from another, the draft  
15 implies publishing separate instructions for a potential  
16 lengthy list of caveats for each individual difference.  
17 This would prove overly burdensome to the manufacturer, as a  
18 small manufacturer of dental products.

19 We would suggest FDA recommend a guideline for  
20 performance ranges of steam autoclaves, unsaturated chemical  
21 vapor sterilizers, and dry heat for handpieces to be tested  
22 in. In this way, a generic testing standard can be  
23 developed which would enhance the purpose of this document.

24 What I wasn't aware of is the labeling for

1 reusable medical handpieces that I picked up this morning.  
2 This is section reference Part C, Section 6, which may  
3 better describe what the FDA really intends to do as far as  
4 the sterilization, what they want the manufacturers to do,  
5 so we might want to try to harmonize those two thoughts. I  
6 would recommend having the panel review Part C, Section 6,  
7 Paragraph G and compare it with Section II.E.1 of the  
8 proposed draft.

9           The next item is also on sterilization validation,  
10 II.E.2.a. The draft, in the first paragraph of II.E.2.a,  
11 calls upon the manufacturer to be significantly more precise  
12 than a dental environment requires. For example, a ten  
13 percent reduction in power is probably not subjectively  
14 discernible. Additionally, different operators will use  
15 their handpieces through varying levels of performance.

16           We believe the issue here should be safety and  
17 efficacy. We would recommend that manufacturers establish a  
18 level of performance that is unacceptable and consider that  
19 level a failure. Manufacturers can then indicate a mean  
20 time between failures specification.

21           In no circumstance should the failure level  
22 represent a level of performance that would compromise safe  
23 operation of the handpiece, and that is the real issue here,  
24 is safety and efficacy of the product. The mean time



1 between failures are consistently used in many manufacturing  
2 specifications.

3           The second paragraph in II.E.2.a leaves a great  
4 deal of latitude to manufacturers to develop a testing  
5 program that could very well compromise or overstate the  
6 life or performance expectations of a handpiece. It could  
7 also understate life and performance. It is not clear in  
8 the proposal.

9           We would recommend at least a minimum guideline be  
10 suggested within the guidance. Care should be taken not to  
11 make the suggested guideline in excess of what dental office  
12 environments represent.

13           Lastly, the concerns in the labeling section  
14 II.F.3.g, Part 4, basically applied to our first point. In  
15 the standard, they want the labeling to reflect all  
16 sterilizers be listed and to qualify each sterilizer for our  
17 particular product.

18           While I'm on the subject, I'd just like to address  
19 the FDA, if I could, and ask them a question. In Section  
20 II.E.2.a, where did the ten percent figure come about, the  
21 ten percent reduction in power?

22           DR. ROBERTSON: I think it is a good question and  
23 I think that is a question we will certainly pass along to  
24 FDA.

1 MR. PEINHARDT: Those are the only items. If the  
2 panel has any questions for me, I'll certainly try to answer  
3 them.

4 DR. ROBERTSON: There was a statement you made  
5 that somehow I missed. It had something to do with making  
6 the guidelines in excess of what dental offices represent,  
7 and I didn't understand the point you were making there.

8 MR. PEINHARDT: Okay. Let me go over that again.  
9 This is in the section that--were you relating to the  
10 discernible, that the dentist would not be able to make a  
11 discernible difference?

12 DR. ROBERTSON: No, you made a statement that  
13 these guidelines should not be in excess of what dental  
14 offices represent toward the end of your presentation.

15 MR. PEINHARDT: Okay.

16 DR. ROBERTSON: I didn't know where you were going  
17 with that.

18 MR. PEINHARDT: That was in the performance area  
19 under handpiece testing, the second paragraph of II.E.2.a.  
20 This is a section that calls out, in order to minimize the  
21 time needed to obtain reprocessing data, accelerated wear  
22 testing, so on and so forth, and it leads into the bottom,  
23 and load parameters that are typical in a single dental  
24 appointment. The emphasis here is this leaves a great deal

1 of latitude for the dental as a manufacturer.

2 We feel there should be some sort of a more  
3 specific guideline so that a testing program that could  
4 compromise or overstate the life of performance of a device,  
5 I think there needs to be--and we, Den-Tal-Ez--thinks that  
6 there should be some sort of at least a guideline for  
7 testing.

8 DR. ROBERTSON: I'm surprised, but all right.

9 Thank you.

10 Any other questions?

11 [No response.]

12 DR. ROBERTSON: Thank you very much.

13 MR. PEINHARDT: Thank you.

14 DR. ROBERTSON: Mr. Steve Jefferies from Dentsply?

15 MR. JEFFERIES: No, I won't be making comments.

16 DR. ROBERTSON: Oh, I'm sorry. Well, if Steve  
17 doesn't want to talk to us, then is there anyone else who  
18 would like to address the panel? Anyone else who would like  
19 to address the panel?

20 [No response.]

21 DR. ROBERTSON: We will now hear from Dr. Michael  
22 Mendelson from FDA.

23 PRESENTATIONS

24 DR. MENDELSON: The Dental Handpiece Guidance

1 Document is directed toward FDA personnel, so that the  
2 review of pre-market notifications for 510(k)s will continue  
3 to be consistent, and it is designed to help industry so  
4 that submissions will be complete. That is, there will be  
5 enough information provided initially to allow the reviewer  
6 to evaluate the handpiece submission quickly.

7           The document addresses the following aspects of  
8 the 510(k). One, the physical description of the handpiece.  
9 Two, identification and description of a predicate device.  
10 Three, performance characteristics. Four, labeling. And  
11 five, ensuring that the handpiece can be sterilized.

12           We would like this document to compliment two  
13 existing documents that also address infection control. The  
14 first document is the draft guideline entitled, "Labeling  
15 Reusable Medical Devices for Reprocessing in Health Care  
16 Facilities, FDA Reviewer Guidance," which was distributed by  
17 the Infection Control Devices Branch of the FDA's Office of  
18 Device Evaluation. It provides items such as an overview of  
19 device reprocessing steps, instructions, information on  
20 documentation of sterilization validation, a check list to  
21 encourage consistent reviews by FDA personnel and  
22 references.

23           The second document is Technical Information  
24 Report, or TIR, No. 12, entitled "Designing, Testing, and

1 Labeling Reusable Medical Devices for Reprocessing in Health  
2 Care Facilities, A Guide for Device Manufacturers," which  
3 was published by the Association for the Advancement of  
4 Medical Instrumentation, or AAMI, in 1994. According to  
5 AAMI, it is intended to, "assist medical device  
6 manufacturers in the design, testing, and labeling of  
7 devices intended for reuse and reprocessing in health care  
8 facilities. Manufacturers may wish to reassess the labeling  
9 of existing products in the light of the recommendations of  
10 the TIR."

11           The Centers for Disease Control and Prevention  
12 recommends routine heat sterilization between patients of  
13 the following three items: High-speed dental handpiece,  
14 intra-oral components of low-speed dental handpiece, and  
15 reusable prophylaxis angles.

16           To ensure that handpieces are actually sterilized  
17 in the clinical setting is the responsibility of both the  
18 manufacturer and the user. The 510(k) should address two  
19 basic aspects: One, that it is physically possible to  
20 sterilize a handpiece, and two, that the labeling is  
21 adequate to allow the user to sterilize properly.

22           Therefore, the manufacturer should provide  
23 instructions on the use, decontamination, sterilization, and  
24 state how many times this cycle of use and reprocessing

1 steps can be repeated before safety or performance  
2 deteriorates excessively. If periodic disassembly,  
3 lubrication, disposal of components, or other service is  
4 needed, the instruction should include this information,  
5 also.

6           Before marketing a handpiece, an applicant must  
7 determine that a handpiece can, in fact, be rendered  
8 sterile. In other words, microbiological techniques should  
9 be used to determine that there is an acceptable probability  
10 that all viable forms of microbial life can be removed or  
11 destroyed from the handpiece if the manufacturer's  
12 instructions are followed. Manufacturers should also  
13 perform testing to verify that durability in terms of the  
14 number of use and reprocessing cycles claimed in the  
15 labeling is accurate.

16           There are several corrections that I would like to  
17 make to this guidance document. One is item II.E.2.a, that  
18 addresses the performance characteristics addressed to in  
19 the sterilization section. The performance characteristics  
20 referred to are listed in item II.D, not L.

21           DR. TYLEND: It would be helpful if you gave the  
22 page number, also.

23           DR. MENDELSON: Page nine. Also, in the labeling  
24 section--

1 DR. GREENSPAN: Could you repeat that?

2 DR. MENDELSON: Sure. When it refers to  
3 performance characteristics, it should be making reference  
4 to Section II.D, not II.L.

5 DR. GREENSPAN: Line three.

6 DR. MENDELSON: I discovered that before this  
7 meeting, but it was too late to mail anything out.

8 Also in the labeling section, II.F.2.a.3, the  
9 phrase "maximum number of use reprocessing cycles before  
10 disposal or repair is required" should have the word  
11 "maximum" deleted.

12 One note on the labeling section. This is not a  
13 sterilization issue. In II.F, item 3, it would be helpful  
14 to the user if the instruction manual for a handpiece were  
15 supplemented with several items. Therefore, paragraph E,  
16 installation and connection instructions, should be modified  
17 with the addition of the following basic items. One, the  
18 maximum free-running operating speed. Two, the minimum  
19 shank length to be fitted inside the chuck. Three, the  
20 maximum overall bur length recommended--

21 DR. TYLEND: Could you give us the page number  
22 that you're on so we could follow along?

23 DR. MENDELSON: Page ten.

24 DR. TYLEND: Page ten, item E?

1 DR. MENDELSON: Item E, yes, installation and  
2 connection instructions. I'll go over them again.

3 Maximum free-running operating speed was the  
4 first. Second, the minimum shank length to be fitted inside  
5 the chuck. I believe these items are listed earlier under  
6 the section where the handpiece is described to the FDA, but  
7 these would be valuable for the user to prevent safety  
8 problems.

9 DR. TYLEND: So if you go back to page four under  
10 the capital letter B, you want to add some of those items--  
11 to copy some of those items into page ten under number E?

12 DR. MENDELSON: That's right.

13 DR. TYLEND: You want to copy item three, item  
14 four, item five, item six, right? So this information  
15 coming into FDA should also be included in the instructional  
16 material for the user?

17 DR. MENDELSON: I don't think the user would be  
18 helped with all of it, but there are certain items that most  
19 clinicians would appreciate to prevent bur tips from flying  
20 off or other such accidents.

21 DR. TYLEND: So those are items--just give us the  
22 numbers and we'll just circle them.

23 DR. MENDELSON: I don't believe maximum operating  
24 speed is--it's on page four. Shank length is item number



1 three. The maximum overall bur length is not in page four.  
2 Air pressure range would be on page five. Item seven, and  
3 that's it.

4 Finally, as mentioned, when other guidance  
5 documents are presented, the dental handpiece guidance is  
6 presented now in draft form because the Dental Devices  
7 Branch is interested in the advice offered by panel members  
8 and industry. Also, this document is subject to continued  
9 changes as knowledge is extended and new designs are  
10 presented.

11 What can the panel do? There are three basic  
12 questions. Would this document be helpful to industry?  
13 Two, does it provide a strong enough framework upon which an  
14 applicant can easily build an adequate pre-market  
15 notification? Three, what suggestions for improvement do  
16 you have?

17 We have provided a list of more specific questions  
18 to help you make your assessment. There is one correction  
19 that needs to be made to these questions. In question  
20 number five, the second sentence lost a phrase when it was  
21 copied into your packet. If you want to get that out, I'll  
22 wait.

23 [Pause.]

24 DR. MENDELSON: It should read, "Is it more

1 appropriate to provide labeling stating the number of cycles  
2 a particular model can withstand, subject to forces such as  
3 the price consumers are willing to pay and the maintenance  
4 steps they are willing to perform?"

5 DR. GREENSPAN: Could you repeat that, please?

6 DR. MENDELSON: Sure.

7 DR. ROBERTSON: Where are you?

8 DR. MENDELSON: I'm on question five. It's the  
9 bottom of the first page of questions.

10 DR. ROBERTSON: Okay.

11 DR. MENDELSON: There is a phrase that was  
12 missing. "Is it more appropriate to provide labeling  
13 stating the number of cycles a particular model can  
14 withstand, subject to forces such as the price consumers are  
15 willing to pay and the maintenance steps they are willing to  
16 perform?"

17 That's it.

18 DR. ROBERTSON: Thank you.

19 I have an initial question, and that was to try to  
20 bring you back to this term reprocessing. I was, in fact,  
21 digging through this paper and the "Labeling Reusable  
22 Medical Devices for Reprocessing in Health Care Facilities",  
23 an FDA reviewer guidance document from the Office of Device  
24 Evaluation in March 1995, in fact, does define reprocessing

1 as cleaning, disinfecting, and sterilizing. But a companion  
2 document, the Technical Information Report, AAMI, 1994, does  
3 define reprocessing in the context of repairing.

4 DR. MENDELSON: I have this document in front of  
5 me. It was my understanding that the Technical Information  
6 Report chiefly addresses the decontamination and  
7 sterilization of reusable devices in the health care site.

8 DR. ROBERTSON: Somebody raised, and I forget who,  
9 maybe it was Tom, a concern about the term. Was it?

10 MR. FISE: It was.

11 DR. ROBERTSON: Maybe you could restate that  
12 concern and make sure I'm not confused, but I was confused  
13 about the use of the term reprocessing.

14 MR. FISE: Yes. Again, Tom Fise with the American  
15 Dental Trade Association.

16 The point we made is that reprocessing is a term  
17 of art that is used with one definition in the GMP document  
18 and it relates strictly to activities that occur in the  
19 manufacturing plant. We are somewhat concerned that this is  
20 confusing enough, that using the same term, even albeit with  
21 a different definition, in this document might confuse  
22 people, to say that everything that applies to reprocessing  
23 in here must apply to what happens in the dental office.

24 We just think the choice of a different term that

1 isn't kind of pregnant with other meanings would probably  
2 make sense.

3 DR. ROBERTSON: Dr. Patters? No?

4 MR. HLAVINKA: We will probably change the  
5 terminology to recleaning and resterilization of devices.  
6 That way, there won't be any confusion with the reprocessing  
7 term.

8 DR. ROBERTSON: I thought it was an interesting  
9 point. I don't care about the terminology, but it would be  
10 useful to make it clear.

11 MR. ULATOWSKI: The term reprocessing has been  
12 used in the infection control community to mean  
13 sterilization and decontamination, and yes, there is overlap  
14 in regard to GMPs. That document on reprocessing is in  
15 draft form and we are collecting comments right now and will  
16 revise it accordingly. It is a good point and one we'll  
17 consider.

18 DR. ROBERTSON: Any other questions from the  
19 panel?

20 MR. HLAVINKA: One last clarification. This  
21 guidance doesn't have anything to do with repair. It's just  
22 the reprocessing aspect, recleaning, resterilization. I  
23 also want to clarify that.

24 DR. O'NEILL: I have a question in that regard.

1 When we are talking about reprocessing and you are talking  
2 about renaming that cleaning and resterilization, I think  
3 there is more implied in this term of reprocessing and that  
4 is cleaning, sterilization, and then general use again. In  
5 other words, there was some discussion in these documents  
6 about taking a handpiece and sterilizing it 14 times versus  
7 taking a handpiece, sterilizing it, using it in a standard  
8 situation, then sterilizing again, then using again.

9           So I think in that term reprocessing, it implies  
10 more than just cleaning and sterilizing. It implies  
11 standard use of the instrument in between, is that not  
12 correct?

13           MR. ULATOWSKI: That's correct, yes.

14           DR. O'NEILL: So if you just go back to a term  
15 "recleaning and resterilizing", it might lose that other  
16 component, is what I'm saying.

17           MR. HLAVINKA: We'll work on the terminology so it  
18 won't be ambiguous.

19           DR. O'NEILL: Okay.

20           DR. TYLEND: The AAMI Technical Information  
21 Report dated 1994 that Dr. Mendelson referred to can be  
22 purchased from AAMI at a cost of \$41 for AAMI members and  
23 \$62 for non-members.

24           DR. ROBERTSON: Any other questions?

1 [No response.]

2 DR. ROBERTSON: We now have a presentation to the  
3 panel by Dr. Chris Miller from the University of Indiana  
4 School of Dentistry. Chris?

5 DR. MILLER: Thank you very much. It's indeed a  
6 pleasure and an honor to be here and to talk to such a  
7 distinguished panel. It's good to see many of you again.

8 When the variety of regulations came out in regard  
9 to heat processing handpieces, I immediately remembered that  
10 many have been doing that for quite a long time, but maybe  
11 in the lack of scientific information. So a few years ago,  
12 I went to look for some of this scientific information to  
13 see if, indeed, there is documentation that we can kill high  
14 levels of bacterial spores that are placed inside of  
15 handpieces. That information at the time, which was about  
16 three years ago, was essentially unavailable except for one  
17 publication that I found in a Swedish journal.

18 So I thought, well, maybe it would be a good idea  
19 to begin to generate some of this information. So I would  
20 like to present to you some of the work that we have been  
21 doing at our lab and simply acknowledge my coworkers here,  
22 and also acknowledge the fact that most of this work has  
23 been funded, at least in part, from a grant from the  
24 American Fund for Dental Health.

1           The purpose of this initial presentation, or  
2 study, I should say, is to determine if the inside of high-  
3 speed handpieces can be actually sterilized by steam and  
4 unsaturated chemical vapor sterilizers. The inoculum that  
5 we have used is the indicated organism for these two methods  
6 of sterilization and that is *Bacillus stearothermophilus*,  
7 and as we have used an organic load of ten percent  
8 defibrillated sheep blood. Ten percent of the total volume  
9 of the spore inoculation contains the sheep blood. As much  
10 as possible, the level of the organisms placed into the  
11 handpiece units are at six logs, or a million spores per  
12 handpiece unit.

13           Method one, the spores in the blood were placed on  
14 turbine fins. The handpiece was assembled and then pre-  
15 flushed for five seconds, which simply means it's hooked up  
16 to the unit and operated with air and water for five  
17 seconds. This was originally used in these early studies to  
18 distribute the inoculum fully within the turbine chamber.

19           Then, in this particular presentation here, spores  
20 were then placed into the water line of that same handpiece  
21 to test both internal sites of the turbine and the water  
22 line. Then the whole handpiece was dried, in this  
23 particular situation, at 50 degrees Centigrade for one hour,  
24 and then it was individually packaged in paper/plastic peel

1 pouches and heat processed.

2           Again, let me kind of visualize again what we did  
3 here. Disassembled a handpiece, and here in these tests we  
4 were using three different brands of handpieces. We  
5 disassembled the handpiece, inoculated the fins of the  
6 turbine with a measured amount of our spores and blood  
7 suspension of a known concentration and the resistance and  
8 so forth, assembled the handpiece as it would be normally  
9 used, hooked it up to the water/air lines, flushed it for  
10 five seconds with the bur in place, again, simply to  
11 distribute the inoculum within the turbine chamber, and  
12 we'll see where this causes some methodology problems here,  
13 but this was the protocol that we were testing.

14           Then the handpiece was inoculated internally into  
15 the water line with the same inoculum of spores and blood,  
16 individually packaged, and then heat processed through a  
17 variety of time/temperature combinations. So a pretty  
18 straightforward protocol in this particular situation.

19           After heat processing, we made attempts to recover  
20 live spores. The turbine and the head and the end cap were  
21 flushed ten times with 4.0 mL of sterile water. In other  
22 words, it was all disassembled and placed in a beaker under  
23 a biosafety hood for prevention of contamination and then  
24 flushed ten times with a single 4.0 mL volume, to knock the



1 spores off of the turbine and so forth into the recovery  
2 fluid. The water line was then flushed ten times with a  
3 separate 4.0 mL volume of sterile water to again recover any  
4 live spores remaining.

5           The recovery fluid was then cultured in double-  
6 strength tripticase soy broth, TSB, and incubated at  
7 standard temperatures for seven days. A viable cell count  
8 on TSA, when that was necessary to verify the challenge  
9 level, was again performed at 55 degrees centigrade for  
10 seven days. So that's how we recovered the spores to see if  
11 we had any sterilization failures.

12           Again, just a visualization here. Disassemble the  
13 handpiece, remove the turbine, place it in the bottom of a  
14 beaker, and just flush this ten times, up and down, up and  
15 down, to knock the spores off. This is, again, in lieu of  
16 submerging the handpiece at the end for culturing. And  
17 again, flushing the water line into a separate beaker and  
18 then culturing this recovery fluid.

19           The results of this particular series of testing.  
20 Using spores and ten percent blood, a gravity steam  
21 sterilizer operating at 121 degrees Centigrade for 30  
22 minutes--now, we did not have any thermocouples inside of  
23 the sterilizers that we used. We didn't do that. It's a  
24 thing that we should do, but it hasn't been done. So we

1 can't verify the internal temperatures here. However, these  
2 sterilizers have been verified by routine spore testing as  
3 well as all of our control spore testing along with the  
4 actual experiments.

5           On the left side here, we see simply three brands  
6 of handpieces, labeled A, B, and C. We just happened to do  
7 B twice in this particular set of experiments. The mean log  
8 of the base ten challenge in the turbine, we could recover  
9 five to six logs of spores. And in the water line, about  
10 the same level. So we had a fairly good challenge, real  
11 close in almost all instances of a million spores per  
12 handpiece.

13           Over here, we have failures. A failure is growth  
14 of the spores that were confirmed as the test organism from  
15 any single handpiece.

16           So in the turbine area, for example, in this first  
17 line here, we did a total of 12 handpieces that were  
18 inoculated both in the turbine and the water line. In  
19 culturing the turbines after heat processing, six of those  
20 12 handpieces had live spores. In this particular case,  
21 three of the 12 handpieces had live spores remaining in the  
22 water line.

23           In this sample down here with brand B, one of the  
24 12 handpieces failed with still growth present in the

1 turbine, and here, three of the 12 handpieces failed with  
2 growth in the water line. When we repeat the same handpiece  
3 in a separate cycle, again, zero of 12 and two of 12  
4 failures. In type C, zero of 12 and two of 12.

5           So, clearly, there are some differences,  
6 obviously, between handpieces, and this is an important  
7 concept to understand, because the internal situation is  
8 different and the availability and access of a sterilizing  
9 agent to the inside is certainly going to be affected by the  
10 physical arrangement of the internal portions.

11           So if you look at the total number of handpieces  
12 in this particular set of information here, there are 48 of  
13 them tested and 29 percent of them failed.

14           DR. ROBERTSON: While you have the slide up, can I  
15 just ask a quick question?

16           DR. MILLER: Yes. Let me put it back here.

17           DR. ROBERTSON: That's all right. You can just  
18 answer the question. Are there handpieces in this run that  
19 were not contaminated with spores but went through the whole  
20 procedure that these did? Are there any control handpieces  
21 where there were no spores, there was no contamination of  
22 the handpiece?

23           DR. MILLER: No.

24           DR. ROBERTSON: No?

1 DR. MILLER: Well, we verified sterility of the  
2 handpieces beforehand, yes. In other words, what we do is  
3 in preparation for these handpieces to--

4 DR. ROBERTSON: Yes, I understand. So there is  
5 not a parallel set of handpieces that were sham contaminated  
6 and then carried all the way through this same--

7 DR. MILLER: Oh, yes, that were contaminated, yes.

8 DR. ROBERTSON: No, not contaminated.

9 DR. MILLER: Oh, no. Not contaminated, no.  
10 There's nothing to recover in an uncontaminated handpiece.

11 DR. ROBERTSON: Dr. Rosan?

12 DR. ROSAN: Yes, these are the presence or absence  
13 of any organisms in terms of failure, is that correct?

14 DR. MILLER: That's right.

15 DR. ROSAN: Whether there was growth or no growth?

16 DR. MILLER: Yes.

17 DR. ROSAN: And that was carried out in a broth  
18 rather than a plate? They were not counted, in other words?  
19 You just--

20 DR. MILLER: No, this was qualitative, growth or  
21 no growth. The entire recovery fluid was cultured to get  
22 maximum chance for organisms to grow in recovery.

23 Yes, sir, Mark?

24 DR. PATTERS: How can you be certain that the

1 bacteria didn't come from the water line and is not part of  
2 your inoculum at all?

3 DR. MILLER: That's a good question. These  
4 particular spores, the *Bacillus stearothermophilus*, are not  
5 common environmental spores. You could probably find them  
6 if you looked hard enough, but in all of the water that I  
7 personally have cultured from dental units, we have never  
8 found anything that grows at 55 degrees Centigrade, so--

9 DR. ROBERTSON: But you were using the spores as  
10 part of this experiment.

11 DR. MILLER: Oh, absolutely.

12 DR. ROBERTSON: You had the inoculum around, so  
13 there was a source of these spores that could have well been  
14 your experimental design, a mistake--

15 DR. MILLER: Sure. The spores are in our  
16 laboratory and we're using them.

17 DR. ROBERTSON: Yes.

18 DR. MILLER: You could have experimental error,  
19 yes, accidental contamination with the test organism.

20 DR. ROBERTSON: Who knows what.

21 DR. MILLER: Understand, thank you.

22 Let's look at the same set of conditions except  
23 under a different sterilization method. Here we have 134  
24 degrees Centigrade for 30 minutes. This is a real serious

1 challenge here, and this is a pretty serious approach to  
2 killing microbes, at 134 degrees Centigrade for 30 minutes.

3           So again, similar challenge failures. A few  
4 failures, still, even at this maximum temperature. In this  
5 case, three of the 12, and here, three of the 12, two of the  
6 12, one of the 12, and again, if you look at the total  
7 number of handpieces in this particular analysis of being  
8 48, 17 percent of them end up failing.

9           A third method of sterilization, or condition of  
10 sterilization, the unsaturated chemical vapor sterilizer  
11 operating at its normal cycle time of 20 minutes at 134  
12 degrees Centigrade. Again, with the three brands, we see a  
13 considerably larger number of failures here, some failures  
14 here, none here. In this particular case, the water lines  
15 were not tested because they became plugged with this  
16 massive amount of spores and blood placed into them and then  
17 drying. We couldn't even get samples out of them, so that  
18 caused us some procedural problems, too.

19           When we got this information, this was instilled  
20 into our thinking, some concerns that, well, maybe we can't  
21 kill organisms inside these handpieces very reliably, so not  
22 as an indication for sterilizability but for some general  
23 information, what if saliva gets back up inside these  
24 handpieces, as we presume it does in many instances. Can we

1 kill the organisms that are present in normal saliva?

2           So we thought we would do a test using the  
3 identical conditions, except instead of using spores, simply  
4 use raw saliva. The raw saliva that we were using had about  
5 a million CFUs per milliliter, a million organisms per  
6 milliliter, but, of course, by the time you put them into  
7 the handpieces with the small volume, the challenge drops  
8 way down.

9           But nevertheless, I think it's important to note  
10 that even in this rather standard steam sterilization cycle  
11 of 121 degrees Centigrade for 30 minutes, we were able to  
12 show at least to be able to kill salivary organisms, which,  
13 again, are not challenging, but I think this is very  
14 important information to have in regard to efficacy and the  
15 reasons for heat processing handpieces in the private  
16 practice. It is important to do this. While we may not be  
17 able to show right now, at least, that we can kill spores in  
18 all test handpieces, we certainly can apparently kill  
19 salivary organisms.

20           DR. ROBERTSON: Chris, what does controls mean in  
21 this context?

22           DR. MILLER: The controls were inoculated  
23 handpieces that were not heat processed so we could measure  
24 the level of the challenge.

1 DR. GREENSPAN: What type of organisms did you  
2 recover from your controls?

3 DR. MILLER: We didn't do an analysis of the  
4 colonies, being so there were a tremendously large number in  
5 raw saliva.

6 DR. GREENSPAN: So you essentially were just doing  
7 bacterial testing--

8 DR. MILLER: Total counts, total counts.

9 DR. GREENSPAN: You weren't doing viral testing at  
10 all, just bacterial counts?

11 DR. MILLER: No, no, no, strictly bacteria,  
12 salivary bacteria.

13 A second method of testing handpiece sterilization  
14 has been used in a couple of instances and we thought we  
15 would like to look at that. In this particular method, the  
16 turbine was completely removed from the handpiece and  
17 replaced with one spore strip cut into six pieces. Then the  
18 end cap was placed back on. The bur hole was then sealed  
19 with a rubber dam patch. The head was wrapped three times  
20 with autoclave tape. It was individually packed and heat  
21 processed, again, looking at another challenge of the  
22 handpiece system.

23 Again, the turbine was removed. A spore strip was  
24 aseptically cut up into six pieces and placed into the



1 chamber and then the end cap placed back on. A rubber patch  
2 placed over the bur hole and then wrapped with tape to  
3 secure the rubber patch in place, again, to supposedly  
4 enhance the challenge to this system. Individually packaged  
5 and heat processed.

6           Culturing of spores, again, here, we simply  
7 aseptically removed the cut-up spore strip and cultured it  
8 in broth, as one would normally culture a spore strip.  
9 This, we found, was very much less challenging than the  
10 previous method we described of using a suspension of spores  
11 and blood.

12           The gravity steam sterilizer, 121 degrees  
13 Centigrade for 15 minutes or 30 minutes, we found total kill  
14 in all of the test handpieces, 12 in each of the runs. So  
15 again, less challenging, which we kind of suspected, but  
16 again, this method has been used in a couple of instances  
17 and we wanted to look at it.

18           Using 134 degrees Centigrade and steamed 15 to 30  
19 minutes, again, no failures.

20           And in an unsaturated chemical vapor sterilizer at  
21 134, pretty good. A couple of failures in one instance here  
22 at their normal half-cycle of ten minutes.

23           A third method, which is very much like the first  
24 except a couple of variations here and there in the drying.

1 Here, we used spores and blood placed on the turbine fins.  
2 We either pre-flushed or not pre-flushed, which means,  
3 again, circulating the spores inside the chamber.

4           The drying here was at a different temperature.  
5 Drying spores at 50 degrees Centigrade, some people may feel  
6 could begin to enhance their germination in the drying  
7 process and make them less susceptible to killing when you  
8 put them in a sterilizer. So we were previously drying at  
9 50 degrees Centigrade, so we thought, well, we don't know  
10 what that effect is, but let's drop the drying temperature  
11 down to room temperature. We know the spores won't  
12 germinate there, and for 24 hours. Then we individually  
13 packaged and heat processed.

14           The effect of pre-flushing on killing, we had the  
15 assumption that if you put spores on the turbine fin and  
16 then you operate that handpiece for a little bit, you'll  
17 blow all the spores around essentially inside the turbine  
18 chamber and get them all into the nooks and crannies. It  
19 will probably make a pretty good challenge. And I think the  
20 data suggests that that's indeed the case.

21           Here's a handpiece, type B, for example, that was  
22 not pre-flushed, and of those five handpieces, we detected  
23 live spores on the turbine fins of one. The same handpieces  
24 that were pre-flushed for five seconds and then heat

1 processed, we found three of the five failures--not a  
2 tremendously large number of samples here, I understand  
3 that, but this is ongoing research.

4           Type D, or brand D, again, no difference. In  
5 fact, it looked a little bit just the opposite. But with  
6 brand E, one failure with no flushing but three failures  
7 with flushing.

8           But there are some other problems with flushing.  
9 When you put spores inside of a turbine chamber and then  
10 hook it up to an air/water system and you operate it, what  
11 happens to those spores? Some of them come out--

12           DR. ROBERTSON: A good question. What happens to  
13 those spores?

14           DR. MILLER: I don't know for sure, but we can't  
15 recover nearly as many. So what are their choices? One,  
16 they're killed by what some friends of mine at CDC have  
17 referred to as centrifugal sterilization. You just blast  
18 them against the side of the turbine chamber wall and smash  
19 them, but that's being facetious. Secondly, they're going  
20 to come out the end cap, and thirdly, out the exhaust air  
21 line, and maybe a little bit out the bur shaft.

22           But whatever happens to them, look here. Again,  
23 just testing here, recovery of spores, with and without pre-  
24 flushing--this is just inoculated and then pre-flushed and

1 then the spores were recovered. This also tells you how  
2 good our recovery system is of flushing the handpiece  
3 turbine with 4.0 mL ten times.

4 Brand B, no pre-flushing, we put in 3.8 times ten-  
5 to-the-sixth spores per handpiece and we recovered 0.6 times  
6 ten-to-the-sixth handpieces, which is a 17 percent recovery  
7 with our normal recovery system. If we pre-flushed that  
8 handpiece, though, we only recovered 4,000 spores, 0.1  
9 percent, so a tremendous loss of your challenge by pre-  
10 flushing.

11 Similar situations with brand D. We recovered  
12 without our pre-flushing system 29 percent of the spores we  
13 put on there and very little after pre-flushing.

14 And here we got, which is not too bad, 55 percent  
15 recovery of the spores. I know it would be nice if we could  
16 recover 90 percent of them, but in this particular instance,  
17 we didn't.

18 So while pre-flushing may tend to blow the spores  
19 around and create a pretty serious challenge, you've got  
20 this other problem that you're losing some and you don't  
21 know where they're going and you don't know what's happening  
22 to them, and that's not a good thing to have in a validation  
23 test.

24 Another question: Is it important to have blood

1 with spores? What effect does that have on sterilization,  
2 killing of the spores inside of the turbine chamber? Again,  
3 effective on killing spores on handpiece turbines, we put,  
4 again, a little over six logs of spores into each handpiece.

5 Brand B, no blood, processed at 121 degrees  
6 Centigrade for 15 minutes, one of five failures. With  
7 blood, same thing, no difference. Of brand D, there was a  
8 difference. One out of five failure without blood but three  
9 out of five failures with blood. The same thing for brand  
10 E.

11 So it looks like, at least with some handpieces,  
12 the presence of blood presents a more difficult challenge,  
13 and this is not new information. This information is  
14 available in many past kinds of studies. Blood presents--or  
15 an organic material presents a more difficult challenge.

16 So here is what we have using blood and no pre-  
17 flushing, to show you a comparison here of three brands of  
18 handpieces. Killing of spores in blood on handpiece  
19 turbines with no pre-flushing, the challenge in each  
20 instance was 3.8 million spores per handpiece, and this, of  
21 course, again, is the number of handpieces that failed  
22 versus the total number tested.

23 Type B, under the gravity steam sterilizer, 121  
24 degrees Centigrade for 15 minutes, one of five failures, one

1 of three failures, and three of five failures.

2 Same steam sterilizer but boosted up at 134  
3 degrees Centigrade, and if you'll notice here, these times  
4 are what we referred to as half-cycles, half the normal  
5 recommended cycle. At 134 degrees Centigrade at two-and-a-  
6 half minutes, one of five failures, none of five, none of  
7 five.

8 In the unsaturated chemical vapor sterilizer, in  
9 their half cycle, one of five, three of five failures, and  
10 zero of five failures.

11 Yes?

12 DR. GREENSPAN: Did you look, and it may be your  
13 next slide, but in case it isn't, did you look at when you  
14 go to full time as opposed to choosing this--

15 DR. MILLER: Not in this particular series of  
16 studies. In the previous ones where we went at 121 degrees  
17 Centigrade for 30 minutes, we found that we still couldn't  
18 get--we had 29 percent failures of the 48 handpieces.

19 That's the information I care to present on high  
20 speeds. Now I have some additional information on slow  
21 speeds. Would you like for me to go ahead and present that  
22 now?

23 DR. ROBERTSON: Surely.

24 DR. MILLER: Okay. Potential for cross

1 contamination with the slow-speed handpiece, again,  
2 acknowledge my coworkers.

3           The purpose of this series of studies was to  
4 determine if the inside of slow-speed handpiece/prophy angle  
5 systems can become contaminated during use, and if so,  
6 possibly indicate that the entire system should be heat  
7 processed rather than covered or surface disinfected.

8           Method: One type of slow-speed motor was used  
9 here with six types of disposable and two types of reusable  
10 prophy angles. The motors and the nose cones and the  
11 reusable prophy angles, anything that could be was steam  
12 sterilized prior to use, and obviously, the disposable  
13 prophy angles weren't. That steam sterilization here, by  
14 the way, was 134 degrees at 30 minutes, a serious steam  
15 challenge. And again, sterility was confirmed prior to  
16 testing.

17           The methods here were performed in replicates of  
18 20. In other words, each time we did a test, we did 20  
19 units, 20 motors with nose cones and prophy angles attached.  
20 That was one test system.

21           Two methods: The first method, we put a test  
22 contaminant in the laboratory in the prophy angel and we  
23 looked at its route of spread from the prophy angle up  
24 inside to the gears of the motor. In the second method, we

1 inoculated the test organism on the gears of the motor and  
2 looked for its spread during operation down through the  
3 handpiece to the prophy angle and even out. So we looked at  
4 contamination coming from the prophy angle up as well as  
5 from the gears of the motor back down.

6 In method one, when we went from the angle to the  
7 motor, we had a handpiece--HP stands for handpiece/prophy  
8 angle system--wrapped up in plastic. So we had the whole  
9 thing connected, ready to go, as if it were going to be used  
10 on a patient, and we wrapped the entire outside with Saran  
11 wrap to make sure that any internal contamination came from  
12 the inside, not the outside.

13 The head of the handpiece system, or the prophy  
14 angle, I should say, was submerged in a test bacterium,  
15 *Serratia marcescens*. This particular bacterium is a gram  
16 negative rod. We use it because when it grows at room  
17 temperature, it produces a very vivid red pigment and we can  
18 detect it from other kinds of possible contaminants. So  
19 it's a marker organism, but similar to many gram negative  
20 rods that are elsewhere in the human body and nature.

21 So after we submerged the head into the suspension  
22 of this bacterium, we turned it on and we pressed the prophy  
23 cup up against the side of the beaker 30 times in a minute  
24 to simulate actual stress of the prophy cup during use.



1           Again, this is a pictorial of what we did. This  
2 is the handpiece/prophy angle system all wrapped up in  
3 plastic. There's the culture of the test bacterium here.  
4 We submerged the prophy cup in there, operated it for a  
5 minute, pressing the prophy cup against the side of the wall  
6 to stress it 30 times in that minute, and we carefully  
7 blotted off the outside of it.

8           We aseptically disassembled the handpiece. We  
9 sampled the inside shaft of the prophy angle. We then  
10 sampled the inside of the nose cone. These were not  
11 quantitative recovery procedures. We sampled the gears of  
12 the nose cone and the gears of the motor to see if the test  
13 bacterium got in the prophy angle and traveled up all the  
14 way.

15           Again, a sample of what we did. After it was  
16 aseptically disassembled, simply taking a paper point and  
17 swabbing it around the shaft of that prophy angle, again,  
18 non-quantitative culturing.

19           Sampling the gears of the nose cone that actually  
20 connect to the gears of the motor, non-quantitative, and  
21 then sampling the gears of the actual handpiece motor to see  
22 if the organism could be detected there.

23           The samples were then placed into--the paper  
24 points were then swabbed onto an auger plate and then--or

1 the swabs, and then the swab and the point were again  
2 dropped totally into broth medium and vortex to make sure we  
3 could recover any organisms that may be present. Then any  
4 growth that was detected was confirmed to be that of the  
5 test organism. We also did sterilized sham, in this  
6 particular case, no inoculation, because this organism  
7 could, by chance, may be a contaminant in the environment,  
8 and found no presence of the organism except when we put it  
9 there.

10 DR. ROBERTSON: So you sampled the lab top or  
11 something?

12 DR. MILLER: That's right. We do this in a  
13 biosafety hood, a laminar air flow system.

14 That's just showing culturing of the paper point  
15 and dropping them in fluid and incubation.

16 Results: Presence of internal contamination when  
17 inoculated at the prophylaxis angle end, and these data represent  
18 a percent of the 20 systems that became contaminated, and  
19 this is the particular prophylaxis angle we used. You can see  
20 here we've used brand names because we're using just about  
21 anything we could find out there in the marketplace. "D"  
22 stands for the disposable prophylaxis angle. "R" stands for  
23 reusables.

24 Culturing the inside of the prophylaxis angle, for

1 example, in this first one, we found inside of the 30  
2 percent of those prophylaxis angles had the test organism inside.  
3 In 15 percent of them, we found the organism in the nose  
4 gears. And in 40 percent of them, we found organisms  
5 present in the motor gears that came all the way up through  
6 the inside system to contaminate the gears of the motor.

7 The same brand here, again, a separate study with  
8 these 20 handpieces. Twenty percent inside the angle, zero  
9 here in the nose cone, and again, this is not quantitative  
10 recovery in any of these cases, just sampling for the  
11 presence or absence, and we certainly could miss some. Five  
12 percent here, and on down the list.

13 Really, the most important, I think, concern is,  
14 do they get all the way up to those gears of that motor?  
15 And in many instances, they do. If you look at the total  
16 number here of positives, this one slide represents a total  
17 of 220 tests. And of those 220 tests, 32, or 15 percent,  
18 ended up with contamination in the gears of the motor.

19 Are there any questions on what we're showing  
20 here?

21 DR. ROBERTSON: What percent of the controls were  
22 positive?

23 DR. MILLER: None. Of the sham inoculated  
24 controls, none, and that was--yes?

1 DR. PATTERS: Just how did you sterilize the motor  
2 to begin with?

3 DR. MILLER: At 134 degrees Centigrade at 30  
4 minutes in steam, and then we also cultured those to verify  
5 that they did not have the test bacterium present.

6 DR. PATTERS: You did that?

7 DR. MILLER: Yes, absolutely.

8 Any other questions before we look at the reverse  
9 of this?

10 [No response.]

11 DR. MILLER: The thinking is, in using a  
12 handpiece, are the organisms present on saliva going to get  
13 up inside and contaminate the inside of the motor, the gears  
14 of the motor? And apparently they do in some instances  
15 here, as you can see.

16 Then the next question is, well, if they're up  
17 there, can they get back out into the mouth of the next  
18 patient? So we again put the test organism now on the gears  
19 of the motor and tested to see if it would travel back down  
20 through the handpiece system and come out the prophy angle.

21 In this particular case, again, the gears of the  
22 motor inoculated with *Serratia marcescens*. The  
23 handpiece/prophy angle system was wrapped in plastic. The  
24 head was submerged now in sterile water as a recovery fluid.

1 It was turned on for one minute and the prophy cup stressed  
2 against the side of the test tube 30 times in that one  
3 minute, and then that water was then cultured.

4 But again, here we're placing now the inoculum on  
5 the handpiece motor gears. That's where the nose cone  
6 attaches to, right there. We put it all together, wrap it  
7 up, and be careful that will eliminate any outside  
8 contamination, put it in sterile water, operate it, and then  
9 we'll culture this sterile water for the presence of the  
10 test organism as well as the inside of the prophy angle and  
11 the nose cone by the same systems that we described before.  
12 Again, not quantitative culturing, but the results.

13 So again, presence of internal and external  
14 contamination when inoculated at the motor gears, percent of  
15 the 20 systems that became contaminated. Here, the inside  
16 of the nose cone, 15 percent of those 20, or three of those  
17 20, became contaminated inside of the angle, none, none, and  
18 so forth down the line.

19 Exit through the prophy angle, this is the  
20 recovery of the organisms that are coming out of the prophy  
21 angle into the sterile water. This is what ended up inside  
22 of the prophy angle from the motor gears. Again,  
23 collectively, this slide represents 200 tests, of which 35  
24 of the tests were positive in this column here, coming out

1 of the prophy angle, which is 18 percent. Again, quite  
2 variable, but appears to occur in some degree in all  
3 systems.

4 Lots of other studies need to be done, and this  
5 basically concludes what I have to say.

6 DR. ROBERTSON: Thank you, Chris.

7 Maybe before we take the break, while Chris is  
8 still here and has his slides up, we'll take questions from  
9 the panel.

10 DR. ROSAN: Chris, was there some consistency in  
11 the contamination? In other words, if you found it in your  
12 nose, did you find it in the other parts, or would these  
13 represent different sorts of things?

14 DR. MILLER: There was some degree of consistency,  
15 but since we didn't have a quantitative flush of a  
16 particular site, you can't really make any inferences, and I  
17 don't know how to quantitatively recover from the inside of  
18 the nose cone. That's my problem. It's a real tough  
19 situation there.

20 But I can tell you, I don't know how many times we  
21 attempted to verify that these were false results. We  
22 looked at all possible routes of contamination. We covered  
23 our handpiece system and put the test organism on the  
24 outside of the Saran wrap to see if there's any way it could

1 get inside. We really--and now we're into trying to do  
2 quantitative tests. How many organisms are getting back  
3 inside and how many are coming back out? Again, you have to  
4 have a quantitative recovery system.

5           The two tests that we have done quantitatively, we  
6 have put the organisms on the gears and then tested  
7 quantitatively the sterile water at the other end, and we've  
8 come up with 100 CFUs of organism in three of the 20  
9 handpieces that were positive. And that's just one test  
10 now, so it doesn't appear that there's a tremendous amount  
11 of influx of these organisms, but at least enough to warrant  
12 some concern.

13           DR. ROBERTSON: Dr. Patters?

14           DR. PATTERS: Chris, based on the sum total of  
15 your experiments using the both high-speed and slow-speed  
16 handpieces that you handled, do you feel that using standard  
17 autoclave practice, 121 degrees, 30 minutes, that you could  
18 sterilize the inside of a handpiece?

19           DR. MILLER: It depends upon how we will define  
20 sterilization. If we define it as I understand it in  
21 accepting the sterilizability that FDA is concerned about,  
22 then I think the answer would be, no, we can't guarantee  
23 that we will sterilize it under those conditions, but maybe  
24 at 134 degrees Centigrade for 30 minutes.

1 DR. PATTERS: Sterile is like pregnant. You are  
2 or you're not.

3 DR. MILLER: And if you assume that the test  
4 organism is--we showed that we could sterilize salivary  
5 bacteria, okay, but not bacterial spores.

6 DR. GREENSPAN: And that's viruses.

7 DR. MILLER: You understand what I'm saying.

8 DR. PATTERS: I do. You're saying you cannot.

9 DR. MILLER: Not bacterial spores, right. But on  
10 the other hand, I think that what is being done is  
11 reasonably okay because at least we can kill some salivary  
12 organisms. I understand, the hallmark is the spore, and  
13 that's what you must gauge everything on. I understand that  
14 very well.

15 DR. PATTERS: What it tells you is that there are  
16 certain parts of the handpiece that do not reach 121 degrees  
17 Centigrade and one atmosphere of pressure for 30 minutes or  
18 the spores would be dead, because they'll be dead--if you  
19 put them on a plate, they'll be dead.

20 DR. MILLER: Right. And again, remember how I'm  
21 coming to you with one particular steam sterilizer that was  
22 in use and these particular brands of handpieces.

23 DR. ROBERTSON: That was very nice, and I think it  
24 does ask questions, anyway.



1 DR. MILLER: Yes, absolutely.

2 DR. ROBERTSON: When you did the sterilization  
3 experiments, did you run parallel spore strips outside the  
4 handpieces?

5 DR. MILLER: Oh, yes, in every cycle.

6 DR. ROBERTSON: And you then cultured that spore  
7 strip as well? I didn't see any of that data.

8 DR. MILLER: No, it's not up in the slides, but  
9 that's standard practice.

10 DR. ROBERTSON: It's kind of important, because--

11 DR. MILLER: Oh, absolutely, and--

12 DR. GREENSPAN: He said it at the beginning.

13 DR. MILLER: And in addition to that,  
14 periodically, we would also put the actual spore strip, yes,  
15 but in addition to that spore strip, we would also place the  
16 same volume of spores put into the handpiece on a piece of  
17 aluminum foil and place that in the chamber, as well, so at  
18 least we could test if we could kill that same exact volume  
19 of spores in blood when it had total access to the steam,  
20 and we always did, under all of those conditions.

21 DR. ROBERTSON: The data wasn't there, but you've  
22 done it, and in the future, that control data will appear?

23 DR. MILLER: When this is published, that control  
24 data will lead the list.

1 DR. ROBERTSON: Good. And while you're doing  
2 controls, you're going to, for either putting the spore  
3 strips that you've cut up and stuffing them into the  
4 handpiece, you're also going to stuff into the handpiece  
5 some strips which have no spores on them?

6 DR. MILLER: We probably won't ever be doing that  
7 kind of testing again. I don't think it's a real good--

8 DR. ROBERTSON: But if you did that--

9 DR. MILLER: But I understand what you're saying.

10 DR. ROBERTSON: If you did do that again, you'd  
11 want to be sure that you had a control in which there were  
12 no spores initially, and when you inoculated your turbine  
13 fins, you'd want to inoculate your turbine spins with  
14 oatmeal or something that was sterile that didn't have any  
15 spores in it and run that in parallel, as well.

16 DR. MILLER: As an environmental check, yes.

17 DR. ROBERTSON: Well, if you're--

18 DR. MILLER: I don't know where else they'd come  
19 from.

20 DR. ROBERTSON: If you're floating spores around--

21 DR. MILLER: Well, they're out there. You've got  
22 to have it on the outside, yes.

23 DR. ROBERTSON: You've got to have a control.  
24 Otherwise, I don't know what to do with your positive

1 results.

2 Yes, Deborah?

3 DR. GREENSPAN: You had no problem killing the  
4 spore strips in the handpiece, and presuming, then, that the  
5 center of the handpiece without the turbine reached 121  
6 degrees Centigrade for that period of time, can you  
7 speculate why you think the autoclave cycles, and, indeed,  
8 the chemclave cycles, are not killing the spores?

9 DR. MILLER: When spores are placed on strips or  
10 anything else, we refer to whatever they're placed on as a  
11 carrier, and I think the carrier can influence, really, the  
12 true resistance of a spore when exposed to some sterilizing  
13 agent. The spore strips, in actuality, when you put a spore  
14 strip in a sterilizer with nothing else in there or fully  
15 exposed to the steam, as it turns out, in many instances,  
16 those spores are killed before it even reaches temperature,  
17 okay?

18 That's the way it is in real life when you do  
19 these kinds of studies. By the time the sterilizer gets up  
20 to temperature, there has been enough heat there to kill  
21 most of the spores within a matter of a very few seconds  
22 afterwards. So the carriers are different.

23 And I'll be honest with you. I've never done what  
24 we call D-value testing or validation of the resistance of

1 spores in a suspension on handpieces, in other words, done  
2 the incremental time exposure and then culture the number of  
3 spores still alive to calculate a D-value, which is a  
4 measure of the resistance of the spores, and that's a stated  
5 fact when you buy a spore strip.

6           It comes with a D-value from one to two minutes,  
7 but it's a totally different carrier, totally different  
8 environment that we're dealing with here, between spore  
9 strips and spores actually dried onto a metal surface. The  
10 heat-up of the metal is going to change things, versus the  
11 heat-up of the paper, strips, lots of variability.

12           DR. ROBERTSON: Other questions?

13           [No response.]

14           DR. ROBERTSON: Chris, thank you very much.

15           DR. MILLER: Sure.

16           DR. ROBERTSON: We will take a break until 4:00.

17           [Break.]

18           DR. ROBERTSON: Let us begin. Carolyn, you are  
19 assembling a table down there.

20           DR. TYLEND: Dr. Mulry, if you would join Dr.  
21 Mendelson and Dr. Kuehne at the table, we would appreciate  
22 it. Kevin, could you pull up one more chair for Dr. Miller?  
23 Thank you.

24           [Pause.]

1 DR. ROBERTSON: I was kind of waiting for Dr.  
2 Miller, as well, because I thought we could tail off of his  
3 presentation.

4 [Pause.]

5 OPEN COMMITTEE DISCUSSION

6 DR. ROBERTSON: Thank you for joining us, Chris.

7 I thought we might actually continue on from your  
8 presentation and everybody can jump in here. Your  
9 preliminary data certainly suggests the possibility that  
10 standard sterilization at present pressure and temperature  
11 could--it is possible that that process will not sterilize,  
12 by our presently-accepted definitions, the interior surface  
13 of handpieces.

14 DR. MILLER: Some of the lower-cycle conditions,  
15 time and temperature. Some of the higher-cycle conditions  
16 looked pretty good in the tests, in other words, 134 degrees  
17 Centigrade for 30 minutes, versus 121 for 15.

18 DR. ROBERTSON: Well, Mark said to you, under  
19 present conditions of sterilization, do you think that those  
20 processes sterilize all handpieces?

21 DR. MILLER: If you're looking for a 100 percent  
22 positive answer, the answer is no.

23 DR. ROBERTSON: But that was based on spores.

24 DR. MILLER: That's correct.

1 DR. ROBERTSON: What that has done is, I think,  
2 set up for you a very nice hypothesis. There was  
3 preliminary data we saw, and you now have, I assume, gone  
4 like a bat based on that preliminary data to do definitive  
5 studies on the sterilizability of handpieces. But based on  
6 that data, the implications of that are that it is possible  
7 that more rigorous sterilization procedures may be necessary  
8 for sterilization of handpieces.

9 DR. MILLER: That's a possibility, which makes the  
10 authors--that's important to the authors, it seems to me, of  
11 this guidance document, because one of the critical issues  
12 here is sterilization of the handpiece, such that that  
13 handpiece is not, as was suggested, one of the primary  
14 health hazards, a source of contamination.

15 DR. ROSAN: I just have a question about one of  
16 the conditions. Were those handpieces lubricated?

17 DR. MILLER: Yes.

18 DR. ROSAN: So they were oiled?

19 DR. MILLER: Yes. The handpieces were processed  
20 according to the manufacturer's directions before they were  
21 inoculated. So if they had to be sprayed before and after  
22 the previous sterilization, they were, or whatever the  
23 recommendations were, with the cleaner or lubricant. So in  
24 those that were regularly oiled, they had oil in them at the

1 time we put the spores on.

2 MR. ULATOWSKI: Mr. Chairman, a comment. I think  
3 Dr. Miller's data highlights one of the dilemmas in  
4 evaluating 510(k)s and how the guidance document will be  
5 implemented by manufacturers in regard to validating  
6 sterilization processes. If there is not a canned process,  
7 a cook book procedure for validation, and one accepts and  
8 relies upon data submitted by the manufacturers, you can get  
9 almost any result you want to get, depending on how you do  
10 your validation study. We see enormous variability in the  
11 manner in which devices are subject to validation and  
12 conflicting data, in many cases.

13 It is very troublesome to us. I think it is an  
14 area that deserves standardization and additional work by  
15 the manufacturing community, and we would solicit their help  
16 in regard to that through some standard-setting organization  
17 or whatever.

18 DR. NORMAN: Paul, might I ask Chris a question?  
19 As you look at the processes that you have been through, is  
20 cleaning more of a necessary step in sterilization than has  
21 been advocated, do you believe?

22 DR. MILLER: I think cleaning is an extremely  
23 important step, and as I look at it, from somebody  
24 interested in infection control and disease prevention,

1 cleaning is to be done to reduce the total bioburden down as  
2 low as possible so that when you get to the sterilization  
3 step, there will be as few organisms left as possible to  
4 kill. So cleaning is very, very important.

5 DR. NORMAN: Do you intend to include things like  
6 ultrasonic cleaning prior to this and removing of oil and a  
7 few steps like this in your process to more fully define the  
8 problems of sterilization?

9 DR. MILLER: That's a very good point. It's  
10 difficult for me, for example, to automatically do  
11 ultrasonic cleaning on handpieces which the manufacturers do  
12 not recommend. So that kind of information, probably, the  
13 best way to clean things, probably has to come from  
14 manufacturers. But nevertheless, the process has to be  
15 done. It's very important to define how things have to be  
16 cleaned prior to sterilization.

17 DR. ROBERTSON: Following up on your point, the  
18 manufacturers' kind of standards, it seems to me that I  
19 would not have expected the result that Dr. Miller reported.

20 MR. ULATOWSKI: I have seen Dr. Miller's data and  
21 it was troublesome, yes.

22 DR. ROBERTSON: Well, no--

23 MR. ULATOWSKI: Disturbing.

24 DR. ROBERTSON: Yes, yes. I would not have



1 expected his answer to be less than 100 percent of the  
2 handpieces were rendered sterile by that method. I mean, I  
3 think that's absolutely wonderful because that's what  
4 science is all about, and I'm thrilled because it opens up  
5 an area to really ask good questions about. But if we are  
6 not even sure, I couldn't blame the manufacturers for not  
7 coming up with that appropriate methodology. I would have  
8 guessed that the standard method of sterilization would have  
9 been sufficient. I worry a little bit about writing a  
10 guidance document based on information you don't know.

11 Deborah?

12 DR. GREENSPAN: I'd like some clarification about  
13 the guidance document for the manufacturers requiring that  
14 handpieces can be sterilized. Are the manufacturers then  
15 required to produce their own documentation showing that  
16 they, in fact, can effectively sterilize a handpiece?

17 MR. ULATOWSKI: As part of the validation process  
18 of the sterilization instructions, they are required to  
19 conduct tests to show that the product can be sterilized and  
20 to maintain the records at their facilities indicating that  
21 that is the case.

22 DR. GREENSPAN: They would be required to do this  
23 sort of tests, along the line of what Dr. Miller described?

24 MR. ULATOWSKI: Yes, exactly, exactly, without

1 doubt.

2 DR. TYLEND: And we ask the manufacturers that  
3 they make sure that the instructions for sterilization are  
4 included in the instructional material that accompanies the  
5 handpiece and that those conditions are the same conditions  
6 used in the validation studies, which seems--I mean, you  
7 would think that would be obvious, but we have seen cases  
8 where that hasn't happened. So we specifically ask that  
9 that be put in the instructional manual.

10 DR. ROSAN: But that seems to create a real  
11 problem, because with what we see here, there could be  
12 enough variation in perhaps even the sterilizer, so that if  
13 you have to follow the manufacturer's recommendation, you  
14 have to buy the same sterilizer, the same capacity. It's  
15 really an enormous problem, I think, in terms of how that's  
16 going to be done. I think we do need to get some real  
17 standardization here so that we don't have all these kinds  
18 of things to discuss.

19 MR. ULATOWSKI: We have come a long way in  
20 infection control procedures even over the last few years,  
21 and with regard to steam, it's been somewhat of a naive  
22 approach to instructions for use, reprocessing with steam,  
23 by manufacturers saying resterilize with steam, when we know  
24 there's many types of steam sterilizers and cycle conditions

1 and all may not be effective. So we asked for more rigor in  
2 terms of testing under controlled conditions to specify the  
3 types of procedures and processes that are necessary.

4 DR. KUEHNE: I just wanted to comment. I think we  
5 ought to make a distinction between something which can be  
6 autoclaved versus something which can be sterilized versus  
7 something which is always going to be sterilized. Like you  
8 pointed out, it depends upon the testing methodology that  
9 you use. You can pretty much come up with whatever you want  
10 to.

11 Dr. Miller has presented evidence to show that  
12 under a broad range of clinical conditions, we can't be sure  
13 that 100 percent of handpieces are actually achieving  
14 sterility by our present definition. However, we can  
15 require, at the very--and the fallout of that is going to  
16 come over the next few years as we do more testing and  
17 research to show the different effects.

18 But what we can right now require, in fact, is to  
19 show that the handpieces that are made can, in fact, be  
20 autoclaved or chemclaved. In other words, the materials  
21 will withstand the conditions the manufacturer has  
22 recommended. The handpiece will withstand 134 degrees for  
23 30 minutes and still function properly. That doesn't assure  
24 that in every clinical condition, bioload, they will

1 actually be sterilized, but they can be autoclaved.

2           The second thing you can require is that under  
3 some test conditions which are reasonable, you can achieve  
4 sterility, again, without the assurance it's always going to  
5 happen 100 percent.

6           Does that help?

7           MR. ULATOWSKI: Yes, very much so.

8           DR. ROBERTSON: That may be helpful, but I guess I  
9 disagree with it.

10           [Laughter.]

11           DR. ROBERTSON: From my perspective, Dr. Miller's  
12 data says we need to ask the question. That's what the data  
13 says. If it turns out that, in fact, handpieces can't be  
14 sterilized at a certain standard pressure and temperature,  
15 then I'm not sure why you'd go through the motions of asking  
16 manufacturers to meet the standard specification. What you  
17 have to do is develop either a new methodology to sterilize  
18 handpieces so you're assured that they are sterilizable or  
19 you've got to build handpieces to a different standard which  
20 the science says will sterilize them.

21           So I think until we know, in fact, that a standard  
22 autoclave procedure will, in fact, sterilize handpieces, I'm  
23 not sure why you'd set that up as a guidance standard. I  
24 think you have to know that first. And all Chris's data