

III. FY 01: BPD Reports Submitted by Manufacturers of Biological Products Other Than Blood and Blood Components (Non-Blood)

The number of reports submitted by non-blood manufacturers increased by 54% from FY 00 to FY01, with the greatest increase being in reports involving precipitate in final allergenic products. Allergenic manufacturers submitted 28 reports in FY-00 and 176 reports in FY-01 related to precipitate found in the final product, events most often discovered by customers.

Total BPDs By Manufacturing System

TYPE OF BPD	ALLERGENIC	DERIVATIVE	IN-VITRO DIAGNOSTIC	THERAPEUTIC	VACCINE	TOTAL	
Incoming Material	2	8	0	1	1	12	3.0%
Process Controls	6	6	4	2	5	23	5.8%
Testing	1	1	1	0	5	8	2.0%
Labeling	27	7	10	1	6	51	12.8%
Product Specifications	183	23	6	10	32	254	63.8%
Quality Control & Distribution	5	11	5	3	1	25	6.3%
Miscellaneous	0	0	1	0	1	2	0.5%
Not Reportable	19	0	1	2	1	23	5.8%
<i>Total</i>	<i>243</i>	<i>56</i>	<i>28</i>	<i>19</i>	<i>52</i>	<i>398</i>	<i>100%</i>

Potential Recalls By Manufacturing System

TYPE OF BPD	ALLERGENIC	DERIVATIVE	IN-VITRO DIAGNOSTIC	THERAPEUTIC	VACCINE	TOTAL	
Incoming Material	0	1	0	1	0	2	4.7%
Process Controls	0	0	3	1	0	4	9.3%
Testing	1	0	0	0	0	1	2.3%
Labeling	9	0	6	0	0	15	34.9%
Product Specifications	3	3	6	5	1	17	41.9%
Quality Control & Distribution	0	0	3	0	0	3	7.0%
Miscellaneous	0	0	1	0	0	1	0.0%
<i>Total</i>	<i>13</i>	<i>4</i>	<i>19</i>	<i>7</i>	<i>1</i>	<i>43</i>	<i>100%</i>

ALLERGENIC

BIOLOGICAL PRODUCT DEVIATION	# REPORTS
<i>INCOMING MATERIAL SPECIFICATIONS</i>	<i>2</i>
Container	
Specifications not met	1
Defective	1
<i>PROCESS CONTROLS</i>	<i>6</i>
Treatment set manufactured with incorrect diluent	1
Manufacturing or processing performed using incorrect parameters	
Aseptic processing procedures not performed according to specifications	1
Other	1
Bulk material does not meet specifications or otherwise determined to be unsuitable	
Discolored	1
Cardboard in vial	1
Contains precipitate	1
<i>TESTING</i>	<i>1</i>
Sterility-performed incorrectly	1
<i>LABELING</i>	<i>27</i>
2 different product species were from same source material	1
Product label incorrect	2
Expiration date – Extended	5
Expiration date – Missing	11
Lot number incorrect	2
Concentration or volume incorrect	6
<i>PRODUCT SPECIFICATIONS</i>	<i>183</i>
Produce specification not met	
Sterility test on associated lot failed	2
Fill volume	2
Appearance	1
Dull needles from vendor	1
Contains precipitate	176
Stability testing failed-potency	1
<i>QUALITY CONTROL & DISTRIBUTION</i>	<i>5</i>
Miscellaneous	
Distributed to incorrect patient	1
Incorrect product type distributed	2
Product distributed inappropriately-failed gown test	1
Shipping and storage-product shipped at incorrect temperature	1
<i>NOT REPORTABLE</i>	<i>19</i>
Release of product other than that which was ordered, labeled appropriately	10
Shipment to incorrect facility	3
Event not under control of reporting establishment	1
Product made available for distribution, but not distributed	5
<i>TOTAL</i>	<i>243</i>

DERIVATIVES

BIOLOGICAL PRODUCT DEVIATION	# REPORTS
<i>INCOMING MATERIAL SPECIFICATIONS</i>	8
Container – defective	1
Closures - specifications not met	1
Closures – defective	2
Source or raw material does not meet specifications or otherwise found to be unsuitable	
Source plasma donor determined to be unsuitable	3
Testing deviation	1
<i>PROCESS CONTROLS</i>	6
Manufacturing or processing performed using incorrect parameters	
Aseptic processing procedures not performed according to specifications	1
Process/procedures - failure to follow prescribed documentation procedures	1
Process water - specifications not met - purified water	1
Bulk material does not meet specifications or otherwise determined to be unsuitable	
Contaminated with microorganisms	2
Stored for an excessive hold time	1
<i>TESTING</i>	1
Sterility-not performed or not documented	1
<i>LABELING</i>	7
Product label missing	3
Carton label incorrect	1
Carton label missing	2
Expiration date missing	1
<i>PRODUCT SPECIFICATIONS</i>	23
Miscellaneous	
Component packaged with final product failed to meet specifications	4
Elevated bilirubin levels	1
Produce specification not met - contains precipitate	4
Stability testing failed	1
Vacuum failure	1
Unacceptable "other protein" peak %	1
% Albumin	1
Non-uniform pellets	1
Stopper discoloration	2
Potency	4
Container closure integrity	1
Chemical analysis/purity	1
<i>QUALITY CONTROL & DISTRIBUTION</i>	11
Product distributed inappropriately-distributed prior to validation of process	1
Shipping and storage-product shipped at incorrect temperature	9
Product identified as unacceptable and not quarantined-unacceptable vacuum results	1
<i>TOTAL</i>	56

IN-VITRO DIAGNOSTICS

BIOLOGICAL PRODUCT DEVIATION	# REPORTS
<i>PROCESS CONTROLS</i>	<i>4</i>
Failure to test treated lots with appropriate antisera	1
Filling not performed according to specifications	1
Process water - specifications not met -reagent water	1
Bulk Material contaminated with microorganisms	1
<i>TESTING:</i> Identity testing not performed or not documented	<i>1</i>
<i>LABELING</i>	<i>10</i>
Product label incorrect	7
Product label missing	2
Lot number incorrect	1
<i>PRODUCT SPECIFICATIONS</i>	<i>6</i>
Produce specification not met	
Unexpected appearance to negative test reactions	1
Increased initial reactive rate	1
Contaminated with green particles	1
Contaminated with microorganisms	1
Potency	2
<i>QUALITY CONTROL & DISTRIBUTION</i>	<i>5</i>
Incomplete number of final packages inspected	1
Shipping and storage	
Product shipped at incorrect temperature	1
Product stored at incorrect temperature	3
<i>MISCELLANEOUS</i>	<i>1</i>
Computer software contained error in dilution calculation	1
<i>NOT REPORTABLE</i>	<i>1</i>
Automated Blood Bank instrument reported ABO or Rh mistypings (device not subject to BPD reporting)	1
<i>TOTAL</i>	<i>28</i>

THERAPEUTICS

BIOLOGICAL PRODUCT DEVIATION	# REPORTS
<i>INCOMING MATERIAL SPECIFICATIONS</i>	<i>1</i>
Container-defective	1
<i>PROCESS CONTROLS</i>	<i>2</i>
Manufacturing or processing performed using incorrect parameters	
Maximum column load exceeded	1
Bulk material does not meet specifications or otherwise determined to be unsuitable	
Contains precipitate	1
<i>LABELING</i>	<i>1</i>
Incorrect patient identification number	1
<i>PRODUCT SPECIFICATIONS</i>	<i>10</i>
Product specifications not met - contaminated with microorganism	5
Stability testing failed	
Percent single chain failure	1
Percent monomer	1
Potency	2
Chemical analysis/purity	1
<i>QUALITY CONTROL & DISTRIBUTION</i>	<i>3</i>
Product distributed inappropriately-distributed prior to validation of process	1
Shipping and storage	
Stored in liquid nitrogen instead of vapor phase	1
Shipping and storage-product shipped at incorrect temperature	1
<i>NOT REPORTABLE</i>	<i>2</i>
Product labeled with incorrect facility identifiers, product acceptable	1
Prior to distribution, product determined to be suitable	1
<i>TOTAL</i>	<i>19</i>

VACCINES

BIOLOGICAL PRODUCT DEVIATION	# REPORTS
<i>INCOMING MATERIAL SPECIFICATIONS</i>	<i>1</i>
Container-defective	1
<i>PROCESS CONTROLS</i>	<i>5</i>
Discrepancy between SOP batch records and product license application	1
Manufacturing or processing performed using incorrect parameters - incorrect column volume	1
Process/procedures – failed media fill	1
Bulk material does not meet specifications or otherwise determined to be unsuitable	
Impurities exceeded specification	1
Other	1
<i>TESTING</i>	<i>5</i>
Sodium Chloride not qualified	1
Safety - performed incorrectly	1
Sterility – not performed or not documented	1
Stability - not performed or not documented	2
<i>LABELING</i>	<i>6</i>
Ink smeared on vial labels	1
Package insert incorrect	1
Expiration date extended	2
Concentration or volume incorrect	2
<i>PRODUCT SPECIFICATIONS</i>	<i>32</i>
Component packaged with final product failed to meet specifications	2
Product specifications not met	
Appearance	1
Contains precipitate	3
Moisture	2
Preservative content	1
Potency	2
Stability testing failed	
pH	2
Fill volume	1
PT responder rate	1
Potency	13
Preservative content	1
Chemical analysis/purity	3
<i>QUALITY CONTROL & DISTRIBUTION</i>	<i>1</i>
Shipping and storage -vials broken or cracked during shipment	1
<i>MISCELLANEOUS</i>	<i>1</i>
Test animals exhibited tetanus toxicity	1
<i>NOT REPORTABLE</i>	<i>1</i>
Intermediate product affected, final product not affected	1
<i>TOTAL</i>	<i>52</i>

NON-BLOOD MANUFACTURES

Timeliness Of Reportable BPDs

Number Of Days From Date BPD Discovered To FDA Received


CUMULATIVE PERCENT OF REPORTS	Allergenic (Days)	Derivative (Days)	In-Vitro Diagnostic (Days)	Therapeutic (Days)	Vaccine (Days)	TOTAL (Days)
10%	18	19	10	32	26	19
25%	25	32	25	42	34	27
50%	36	44	42	47	45	43
75%	53	68	67	53	63	56
90%	88	92	96	140	77	140
# REPORTS	214	53	24	17	47	355
RANGE	3-899	10-848	9-160	32-581	9-207	3-899
AVERAGE	76	87	57	109	53	75
# Reports lacking date discovered	10	3	3	0	4	20

Adherence To 45 Day Required Time For Reporting

(Reporting Time = Date of FDA receipt – Date of discovery of BPD)

Reporting Time (days)	Allergenic		Derivatives		In-Vitro Diagnostics		Therapeutics		Vaccines		Total	
< or = 45	137	64.0%	24	45.3%	12	50.0%	5	29.4%	23	48.9%	201	56.6%
Between 45 and 90	44	20.6%	16	30.2%	7	29.1%	9	52.9%	21	44.7%	97	27.3%
> 90	33	15.4%	13	24.5%	5	20.8%	3	17.6%	3	6.4%	57	16.1%
Total	214	100%	53	100%	24	100%	17	100%	47	100%	355	100%

FY-2001 Non-Blood BPD Reporting Timeliness



Total Reports = 355 Allergenics = 214; Derivatives = 53; In-Vitro
Diagnostics = 24; Therapeutics = 17; Vaccines = 47