



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

Manufacturing standards for plasma for fractionation

Scientific relevance and regulatory requirements

Albert Farrugia

Senior Principal Research Scientist and Head,

Blood and Tissues Unit

Office of Devices, Blood and Tissues



Presentation outline

- Current standards
- Empirical observations
- Basic science
- Resulting tensions
- Possible approaches

WARNING :


THIS PRESENTATION IS PEPPERED WITH
NOSTALGIA AND OTHER MANIFESTATIONS
OF PERSONAL INDULGENCE



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

Available standards



2001: 0853

HUMAN PLASMA FOR FRACTIONATION

Plasma humanum ad separationem

2.2.4 Conditions of storage and transport of plasma.

See Annex IV and V.

Describe the conditions for freezing and storage of plasma for every establishment responsible for collecting blood/plasma including the following:

- Sites/organisations which are involved in the storage and indicate whether they have been inspected by a Competent Authority.
- Compliance with Ph. Eur. with respect to freezing and storage.
- Conditions of storage (temperature and maximum time).

Describe the conditions of transport of plasma including the following:

- Transport flows from centres of collection to interim storage sites, if relevant, and further to fractionation sites.
- Organisations which are involved in the transport (own and contractors) and indicate whether they have been inspected by a Competent Authority.
- Conditions of transport (maximum time and temperature).

⁹ ... for Fractionation and if applicable, with any Ph. Eur. requirements for particular products.



Council of Europe Guide for Blood Components



Guide to the preparation,
use and quality assurance of
blood components
introduction



- Includes chapters on FFP, cryo-poor plasma
- FFP standards sometimes at variance with EP monograph
- NOT APPLICABLE for fractionation - refers to EP monograph



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

Available standards



National Archives and
Records Administration

code of
federal regulations

TITLE 21

FOOD AND DRUGS

PART 640

ADDITIONAL STANDARDS FOR
HUMAN BLOOD AND BLOOD
PRODUCTS

Subpart G--Source Plasma



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

A contentious statement

*The regulatory requirements
underpinning blood and plasma storage,
freezing and frozen storage are
predicated on the needs of Factor VIII*



FVIII in 2004

- Plasma-derived FVIII production is becoming increasingly marginal in the developed blood economies
- Fractionators still ship plasma for FVIII manufacture in the hope of supplying the “emerging” markets
- Factor VIII is the most labile plasma therapeutic protein
- Conditions affecting FVIII may affect other proteins in ways which are still unknown
- Tailoring conditions to optimising FVIII preservation is therefore still a valid goal



Australian Government
Department of Health and Ageing
Therapeutic Goods Administration

European plasma standards

FVIII levels

Council of Europe

(for transfusion)

Requirement for $\geq 70\%$ of the “*average normal value*” controlled through measurement of FVIIIc every two months on a pool of six units of mixed blood groups during the first and last months of storage

European Pharmacopeia


(for fractionation)

On a pool of not fewer than ten units, measurement of factor VIII, using the EP reference method and a reference plasma calibrated against the International Standard for blood coagulation factor VIII in plasma. The activity is not less than 0.7 I.U. per millilitre.



Factors claimed to affect FVIII yield in fractionated concentrates

- Anticoagulant
- Collection method
- Time/Temperature to separation/freezing
- Freezing rate
- Storage conditions of frozen plasma
- Thawing conditions
- Purification chemistry
- Viral inactivation






Australian Government

Department of Health and Ageing
Therapeutic Goods Administration

Plasma quality

The role of Factor VIII

FACTOR VIII IN BLOOD BANK CPD DONATIONS



Farrugia et al 1984




Influence of blood cooling on FVIII in plasma

Cooling temperature	FVIII IU/ml	Protein g/l
0-4°C n=9	0.45 \pm 0.06	63.7 \pm 2.2
20°C n=9	0.84 \pm 0.1	65.6 \pm 2.7



Factors claimed to affect FVIII yield in fractionated concentrates

- Anticoagulant
- *Collection method*
- Time/Temperature to separation/freezing
- Freezing rate
- Storage conditions of frozen plasma
- Thawing conditions
- Purification chemistry
- Viral inactivation



Source vs recovered plasma


FVIII yield in low purity concentrates

	Cryo g/kg plasma	FVIII IU/kg plasma	
		Cryo extract	Pre-finish
Manual <18 h, CPDA n=13¹	9.1±0.5	391±31	290±25
Haemonetics CPD¹	10.9±1.1	461±50	319±33
Whole blood²		403	283
Apheresis²		450	317
Apheresis/low citrate²		507	362
¹ Smith et al (1985) - Low purity concentrate			
² Ribeiro et al 1997 - Intermediate purity concentrate			



Factors claimed to affect FVIII yield in fractionated concentrates

- Anticoagulant
- Collection method
- ***Time/Temperature to separation/freezing***
- Freezing rate
- Storage conditions of frozen plasma
- Thawing conditions
- Purification chemistry
- Viral inactivation





Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

DONATION-FREEZING INTERVAL EFFECT ON FVIII LEVELS ARCBS data mid 1990's

DELAY TO FREEZING

% < 0.7 IU/mL

APHERESIS < 12 HOURS	1%
WHOLE BLOOD < 12 HOURS	13%
WHOLE BLOOD < 18 HOURS	27%
WHOLE BLOOD < 24 HOURS	40%



Factor VIII content of 5 litre plasma packs of different ages

	Factor VIII (IU/kg) in plasma samples		
Time (h) between blood collection and plasma freezing	Before freezing pack	After thawing pack to 20°C	Core sample of frozen pack
3-4 - special collection	930	870	880
4-8 - routine fresh frozen	840	820	790
16-18 - routine overnight frozen	800	740	730



Effect of pack type and freezing rate on plasma FVIII

	Plasma IU/ml	
	6-8 hours	18 hours
5 L packs over \approx 2 h	0.84 \pm 0.13 (n=10)	0.77 \pm 0.08 (n=10)
SD packs over \approx 15 min	0.86 \pm 0.10 (n=10)	0.77 \pm 0.15 (n=10)



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**


DOES IT MATTER?

- There is no doubt that delayed blood processing to frozen plasma decreases FVIII levels in plasma for fractionation
- Does this affect the yields and quality of fractionated products?



Distribution of FVIII in small scale plasma cryoprecipitation

Effect of overnight storage



Hellings et al (1982)




Plasma Quality

Effect of delayed freezing on FVIII

Hughes et al 1989

“..... plasma intended for the recovery of proteins that are labile in plasma is frozen by cooling rapidly at – 30 °C or below as soon as possible and at the latest within 24 h of collection.” EP Monograph



Effect of separation/freezing interval on FVIII yield in low purity concentrates

	FVIII recovery IU/kg plasma			
	ACD		CPD	
	6-8 hours	18 hours	6-8 hours	18 hours
Cores of frozen packs (mean of 10 cores of five L packs) ¹	720	700	840	770
Freeze-dried concentrate (mean of 8 batches) ¹	214	208	255	257
Continuous thaw ²			293	268
Batch thaw ²			204	174

¹ Smith et al (1985) Dev Hem & Imm 13 :15-23

¹ Foster et al (1985) Dev Hem & Imm 13 :15-23



Australian Government

Department of Health and Ageing
Therapeutic Goods Administration

DOES IT MATTER?

- There is no doubt that delayed blood processing to frozen plasma decreases FVIII levels in plasma for fractionation
- Does this affect the yields and quality of fractionated products?


IT DEPENDS

- Cryo yield affected
- LP & IP sometimes affected
- No data for current generation of FVIII concs



Factors claimed to affect FVIII yield in fractionated concentrates

- Anticoagulant
- Collection method
- Time/Temperature to separation/freezing
- ***Freezing rate***
- Storage conditions of frozen plasma
- Thawing conditions
- Purification chemistry
- Viral inactivation





Australian Government

Department of Health and Ageing
Therapeutic Goods Administration

“Plasma should be frozen at.....

.....-18^oC, -20^oC, -23^oC,-30^oC?”

- Remarkably ambiguous language in standards
 - “*cooling rapidly at -30^oC, frozen at -20^oC*” (EP)
 - “*shall be stored at a temperature not warmer than -20^oC*” (CFR)
- Little recognition of the important - obvious - parameter

THE FREEZING ***RATE***




Figure 4. Freezing and supercooling in plasma. **: Plasma placed directly on -50°C cold shelf; ••: Plasma cooled to below 0°C on a -13°C shelf, then shelf temperature reduced to -50°C .

Plasma & freezing medium temperatures


Farrugia et al 1992

Plasma - -30°C cold room

Plasma - -30°C Instacool

Medium - -30°C Instacool


Medium - -30°C cold room



Plasma freezing rates and FVIII

PLASMA FREEZING RATES

Farrugia et al 1992



FVIII YIELDS

EFFECT OF PLASMA FREEZING RATE


<u>FREEZING MEDIUM</u>	<u>FVIII U/KG</u> <u>(Blood Bank</u> <u>CRYO)</u>
(1) INSTACOOL FREEZER	575±122
(2) - 30 ° C COLD ROOM	443±62
(3) -40 ° C CABINET FREEZER	303±130 *

* P<0.01 vs (1) & (2)

Farrugia et al 1992




Plasma freezing time to -25°C with different equipment





FVIII recovery in cryoprecipitate and cryosupernatant with different freezing methods




Carlebjörk et al (1986)

Transfusion 26:159-162




Plasma frozen in different media



FVIII related activities in plasma fractions

Thaw-siphon cryoprecipitation - effect of plasma freezing rate





Effect of freezing rate on cryoprecipitate quality

Freezing condition medium	Temp	Thickness of plasma layers		FVIII IU/L	Fibrinogen g/L	Total protein g/L
		4 cm	2 cm			
Alcohol dry ice	-70°C	11 min	7 min	467	10.5	23.2
Circulating N ₂ -gas	-100 °C	12 min	10 min - fast	513	11.2	23.7
Stationary air	-30 °C	4 hours	3 hours	490	13.0	29.5
Idem + insulation	-30 °C	19 hours	15 hours - slow	433	14.6	34.9



Plasma freezing


What is important?

- Rapid freezing - ca -30°C in 30 minutes - results in better FVIII yields in cryo relative to slower freezing - - ca -30°C in 3-4 hours
- The ice crystal structure and the physical nature of cryoprecipitate are affected by the plasma freezing rate.
- Slower freezing also increases fibrinogen in cryo; this has its pro's and con's
- The effect of freezing rates on FVIII yields in current concentrates is not well recorded



Factors claimed to affect FVIII yield in fractionated concentrates

- Anticoagulant
- Collection method
- Time/Temperature to separation/freezing
- Freezing rate
- *Storage conditions of frozen plasma*
- Thawing conditions
- Purification chemistry
- Viral inactivation





Plasma freezing and storage

Effect on thaw siphon blood bank cryoprecipitate

Freezing	Storage period	Storage temperature	FVIII cryo yield IU/kg plasma	Fibrinogen in cryo mg/kg plasma
Fast	16 h	- 20 ° C	426	605
	3 mo		500	609
	6 m		416	538
	16 h	- 40 ° C	493	607
	3 mo		522	577
	6 mo		449	542
Slow	16 h		318	522
	3 mo		306	502



Effect of variation in plasma cold storage on FVIII recovery and cryoprecipitate quality



	Cold Storage °C	
	-40 °C to -20 °C (1 month) to -40 °C	-40 °C
Number of batches	5	9
Batch size (L)	751 \pm 80	718 \pm 77
Plasma FVIII (IU/L)	795 \pm 84	790 \pm 77
Cryoprecipitate weight (g/L plasma)	11.63 \pm 0.48	10.98 \pm 0.57
Cryoprecipitate extract FVIII (IU/L plasma)	597 \pm 25	594 \pm 40
Cryoprecipitate quality (IU/L)	51.4 \pm 2.7	54.2 \pm 2.3



Plasma Quality

Effect of poor storage conditions

“.....*Source Plasma intended for manufacture into injectable products that is inadvertently exposed (i.e., an unforeseen occurrence in spite of compliance with good manufacturing practice) to a storage temperature warmer than -20 deg.C and colder than +10 deg.C may be issued only if labeled as ``Source Plasma Salvaged.'`*” CFR 21- 640





Effect of thawing and re-freezing CPD plasma on plasma FVIII and concentrate yield

	Stage yield IU/kg		
Stage	Once-frozen	Twice-frozen	Significance
Plasma cores (plasma standard)	593 \pm 99	301 \pm 116	Significant
Cryo extract (conc. Standard)	330 \pm 35	291 \pm 35	Not significant
Dried low purity concentrate	194 \pm 19	192 \pm 19	Not significant




Figure 1. Schematic representation of the present Scottish procedure for thawing plasma.




Australian Government

Department of Health and Ageing
Therapeutic Goods Administration

Plasma conditioning

Effect on FVIII concentrate





Effect of plasma conditioning on cryoprecipitate FVIII

Conditioning regimens				
■ PCU temperature	-10 °C	-15 °C	+4 °C	-16 °C
■ Conditioning time	6 hours	6 hours	2 hours	5 hours
■ Plasma temp at crushing	-8 °C	-12 °C	-10 \pm 2 °C	-10.5 \pm 0.5 °C
■ Plasma storage temp	-20 °C	-20 °C	-40 °C	-40 °C
Number of batches	18	5	36	13
Batch size L	750 \pm 54	789 \pm 84	761 \pm 104	785 \pm 129
Plasma FVIII (IU/L)	770 \pm 68	789 \pm 84	764 \pm 130	769 \pm 72
Cryoprecipitate weight g/L plasma)	12.21 \pm 1.19	11.63 \pm 0.48	10.7 \pm 0.64	10.05 \pm 0.34
Cryoprecipitate extract FVIII (IU/L plasma)	567 \pm 47	597 \pm 24.8	509 \pm 42	551 \pm 50.3
Cryoprecipitate quality (IU/g)	46.4 \pm 4.1	51.4 \pm 2.7	47.8 \pm 4.6	54.4 \pm 5.5



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

Plasma conditioning

Effect on blood bank cryoprecipitate

Plasma processing conditions	FVIII IU	Fibrinogen mg	Fibronectin mg	Adhesive strength g	VWF U
Stored -30°C, immediately thawed in 4 °C WB	132 \pm 24	213 \pm 82	53 \pm 21	7.6 \pm 1.8	189 \pm 31
Stored -30 °C, then 4 °C cold room (in polystyrene box) for 18 h, then thawed in 4 °C WB	110 \pm 22	457 \pm 125 p<0.001	63 \pm 17	23.5 \pm 4.8 p<0.05	221 \pm 25

Farrugia et al (1992) Transfusion 32:755-759



Plasma Storage

What is important?

- As long as freezing is optimised, storage requirements appear to be flexible in the range -20°C to -40°C
- Maintaining a steady storage temperature is more important than the absolute storage temperature, within this range
- While temperature changes can affect the quality of cryoprecipitate, this can be exploited to improve both blood bank and industrial cryo

The “theory” of plasma freezing - *CBBS e-Network Forum*

A transfusion medicine physician in the Netherlands reports that the optimal storage temperature of fresh frozen plasma is **minus 30 degrees** centigrade or colder as is generally agreed upon in Europe and as can be read in "The guide to the preparation, use and quality assurance of blood components", Council of Europe Publishing, ISBN 92-871-3530-4. The scientific background of this choice is the following:

*"Plasma is a solution of proteins and salt in water. Freezing of such a solution results in the formation of pure water ice until the **eutectic** freezing point at minus 23 degrees centigrade is reached and then also the solutes like proteins and salts start to form crystals. Finally, the total plasma mixture has frozen solid after the eutectic point temperature has been reached. During this freezing process the remaining liquid gradually becomes an ever more concentrated solution of the salts and proteins.....Labile proteins like factor VIII and other clotting enzymes will denature when exposed for a long time to these highly concentrated acid salt solutions. During the frozen storage, the freezer with plasma will be opened and closed regularly in a blood bank, and increases in temperature inside a freezer of up to 5 degrees centigrade can easily occur. At any temperature higher than minus 23 degrees centigrade the original eutectic mixture of salts, proteins and water will become liquid again and the deterioration of the labile proteins will resume. If you keep the average storage temperature of the frozen plasma below minus 30 degrees centigrade, you avoid the critical eutectic temperature and the plasma will maintain its original quality."*




Australian Government

Department of Health and Ageing
Therapeutic Goods Administration

Oh yeah?



Eutectic point of plasma?



Transition	Temperature of transition (°C)	
	Slow freezing (2°C/min)	Rapid freezing (200°C/min)
1 Glass transition (onset of motion of water molecules)	- 80	- 85
2 Antemelting (onset of molecular mobility of proteins)	- 42 to -38	- 38 to -35
3 Incipient melting (beginning of thermodynamic melting of ice)	- 27	- 27
4 Melting point (final melting of ice)	- 0.5	- 0.5
MacKenzieAP 1980		



Phase transitions in frozen plasma

Resistivity of normal plasma & saline
measurements on slow thawing after fast freezing
McIntosh 1990

Freezing of 700 ml plasma

Energy consumption at different stages of freezing

Carlebjörk et al (1986) Transfusion 26:159-162





Plasma freezing and storage

- Conventional eutectics offer no guidance
- Freezing so that phase change is as rapid as possible
- Storage so that this is maintained - -20°C is adequate

⇒ AND WHY SHOULD THIS BE AN ISSUE FOR
REGULATORS ANYWAY?

⇒ IS THERE ANY EVIDENCE THAT
BLOOD/PLASMA PROCESSING AFFECTS SAFETY
AND QUALITY (AS OPPOSED TO YIELD)?



FVIII and activation of coagulation in plasma freezing

Split 300 ml pairs from 600 ml source
plasma units, n=12

	-30°C	-80°C	% difference	p
FVIII:C (u/dL)	103 (84-16)	121 (108-149)	16	0.0005
FVIII:Chr (U/dL)	96 (85-107)	105 (90-117)	8.9	0.0005
C:Chr ratio	1.05 (0.88-1.21)	1.15 (1.05-1.43)	9.0	0.0005
F1+2 (nmol/L)	0.60 (0.36-1.74)	0.84 (0.43-2.11)	33.3	0.001


NB - *FAST* FREEZING RESULTS
IN *HIGHER* ACTIVATION



FVIII from concentrates made from plasma pools with evidence of coagulation activation

	Normal coagulation markers	Elevated coagulation markers
FPA in plasma pools	$1.7 \leq \text{FPA} \leq 4.5 \text{ } \mu\text{g/ml}$	$11.2 \leq \text{FPA} \leq 13.4 \text{ } \mu\text{g/ml}$
TAT in plasma pools	$2.1 \leq \text{TAT} \leq 3.0 \text{ ng/ml}$	$11.4 \leq \text{TAT} \leq 15.9 \text{ ng/ml}$
40 kDa fragment FVIII heavy chain fragment	Absent	Present
Inhibitor development	None of the batches resulted in inhibitors	All inhibitor patients received these batches
Binding of FVIII in product to PS/PC	156	30
Binding of FVIII in product to mcab to C2 domain overlapping with VWF/PL	156	30
Binding of FVIII in product to mcab to C2 domain not overlapping with VWF/PL	13	15

Effect of mixing and donation time on plasma FpA



Blood processing and activation of coagulation

- Standard blood collection measures minimise activation
- No special measures are needed
- Reported FpA levels leading to problems are abnormally high



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**


And of course, there are other things
one can get out of plasma.....



Plasma Quality

Effect on IMiG fragmentation during storage

“..... When obtained from whole blood, plasma intended solely for the recovery of proteins that are not labile in plasma is separated from cellular elements and frozen at – 20 °C or below as soon as possible and at the latest within 72 h of collection..” EP Monograph



McIver 1980



Plasma quality

Effect on albumin solutions

- Albumin made from plasma recovered from outdated blood shows higher PKA levels (BoB W/S 1977)
- Albumin made from haemolytic plasma was unstable at 25% concentration (Boros et al 1974)



Australian Government

**Department of Health and Ageing
Therapeutic Goods Administration**

.....are these issues mainly of historical interest.....or can other plasma proteins be affected by poor storage conditions?.....

.....is this part of the great unknown.....and therefore subject to regulatory precautionism?.....




Plasma Manufacture

What is a quality product?

The characteristics of an item or process that indicate its conformance to designated parameters, and its degree of perceived customer acceptance or satisfaction. Quality characteristics often include ***reliability, consistency and the ability to continue performance in stress or volume situations***, but are critical only in relation to the value placed on them by the user or customer.

<http://www.bridgefieldgroup.com/glos7.htm>

How can plasma be assured to a high level of

- 
- **Reliability**
 - **Consistency**
 - **Ability to continue performance in stress or volume situations**

⇒ A defined manufacturing process

⇒ Specified freezing and storage conditions

⇒ Robustness to volume and temperature changes



Tentative conclusions and possible approaches

- There is a need for clear and unambiguous standards for plasma freezing and storage
- A process which results in a consistent product, irrespective of scale and location, should form the basis of any standard
- Empirical observations appear to support greater flexibility than some current requirements
- There is little evidence that any of these requirements have a bearing on product safety
- Basic conditions for minimising microbial contamination and preserving product integrity should be defined
- Other requirements reflecting product yield eg FVIII levels should be left to be negotiated between the manufacturer and plasma supplier

Thank you for reminding me of

When We Were
Very Young

