
FDA PLASMA STANDARDS WORKSHOP

**August 31st, 2004
Bethesda, MD**

**Barbara Glantschnig, MSc.
Head Quality Assurance Plasma
Octapharma**

Plasma Specifications

- Octapharma uses different plasma types for
 1. Production of SD-Plasma
 2. Fractionation of coagulation factors and factor complex
 3. Fractionation of IVIG and Albumin

Use of Plasma Qualities

Type	SD Plasma	F IX and Complex	F VIII	IVIG, Albumin
RP < 8h	X	X	X	X
SP	X	X	X	X
RP < 24h	-	-	X	X
RP < 72h	-	-	-	X

Plasma Specification Details

Type	Separate	Freeze	Store	Ship	Expires
RP < 8h	Whole Blood/ High Spin	< -30°C core T/1h	<-20°C	<-20°C	1 year
SP	Apheresis	< -30°C	<-20°C	<-20°C	2 years
RP < 24h	Whole Blood/ High Spin	< -30°C	<-20°C	<-20°C	2 years
RP < 72h	Whole Blood/ High Spin	< -30°C	<-20°C	<-20°C	2 years

Yield Influencing Parameters – Coagulation Factors

- Time between blood collection and plasma separation
 - RP 8h and SP are preferred for labile products
 - RP 24h still acceptable for FVIII
- Method of plasma separation
 - Low cell contamination in plasma
 - Prevention of hemolysis
- Method of plasma freezing
 - SP → normal freezing process sufficient
 - RP 8h, RP 24h → shock freezing at $< -30^{\circ}\text{C}$ improves yield

Conclusion

- Octapharma plasma specifications are based on current Ph.Eur. regulations.
- For SD plasma, stricter requirements apply as a minimum level of defined coagulation factors must be achieved in the final product.
- For storage and transport of plasma, $<-20^{\circ}\text{C}$ are sufficient based on experience with the given expiration dates.