

FDA Plasma Standards Workshop

August 2004

Mary Ann Lamb, Ph.D.
Senior Director, Regulatory Affairs

Bayer Plasma Summary

Raw Material	pdFVIII	IGIV IGIM	Alpha₁-PI	Albumin PPF
Source Plasma	Yes	Yes	Yes <i>and</i> Fr IV-1 Paste Intermediate	Yes
Recovered Plasma	No	No	Fr IV-1 Paste Intermediate	No

Bayer Plasma Summary

Raw Material	Freeze Temp.	Store Temp.	Ship Temp.	Time to Expiry
Source Plasma (Plasmapheresis)	-20°C or Colder	-20°C or Colder	-20°C or Colder	3 Years
Plasma Used to Produce Fr IV-1 Paste Intermediate	SP -20°C or Colder RP -18°C or Colder	SP -20°C or Colder RP -18°C or Colder	SP and RP : -20°C or Colder	As Specified by Fractionator

Bayer Source Plasma Specifications

- Plasma is considered suitable for fractionation for the US market if the following conditions are fulfilled
 - The total period of time during which the temperature exceeds -20°C does not exceed 72 hours and the plasma remains frozen solid during storage
 - The temperature at no time exceeds $+10^{\circ}\text{C}$
 - 21 CFR §640.76
- Plasma is considered suitable for fractionation for Europe if the following conditions are fulfilled
 - The total period of time during which the temperature exceeds -20°C does not exceed 72 hours
 - The temperature does not exceed -15°C on more than one occasion
 - The temperature at no time exceeds -5°C
 - EP Monograph *Human Plasma for Fractionation*

On the Use of Plasma Derived from Whole Blood (Recovered Plasma)

- Bayer does not currently fractionate recovered plasma for commercial product manufacture and distribution
- Bayer manufactures Alpha₁-Proteinase Inhibitor from Fraction IV-1 paste intermediates
 - Recovered Plasma stored at -18°C or colder and shipped at -20°C or colder
 - Source Plasma stored and shipped at -20°C or colder
 - Fraction IV-1 paste stored and shipped at -20°C or colder
- Alpha₁-Proteinase Inhibitor manufactured from intermediate paste has the same final product specifications and shelf-life with no differences in the final product stability as compared to Alpha₁-Proteinase Inhibitor manufactured from Source Plasma

Bayer Supporting Data Summary

- Bayer Plasma Stability Study (10 year study)
 - Data available for 3 years of storage for Source Plasma collected by plasmapheresis
 - Stability indicating parameters monitored: FVIII potency, anti-HBs Ag potency
- Source Plasma collected from 52 random donors was thawed, pooled and transferred to polyethylene plasma collection bottles
- Polyethylene bottles stored at -20°C or colder and thawed at 5°C for 19 hours prior to testing
- Evaluation of the pooled data trend for the parameters monitored indicate:
 - FVIII potency - no significant change through 36 months
 - Anti-HBs Ag potency - no significant change through 36 months

Conclusions

- The existing US regulations regarding freezing, storage and transport of plasma for fractionation have been in place for decades and have served the consumers and industry well
- The decreased demand of pdFVIII brings into question the need to invest millions into increasing factor VIII yield through plasma collection activities
- In the absence of a demonstrated improvement in the quality of derivatives manufactured from plasma frozen rapidly after collection, manufacturers should be afforded the flexibility to improve yield through processing innovation and optimization