

**FDA TSE Advisory Committee
Rockville, Maryland
June 26-27, 2002**

Perspective on TSE Clearance Studies

Robert G. Rohwer, Ph.D.

**VA Medical Center
Baltimore, MD**

file: FDA TSEAC Clearance Studies 6-26,27-02.ppt

Clearance Studies 1:

- | | |
|------------------------------------|---------------|
| ■ Bovine albumin | brain spike |
| ■ Bovine aprotinin | brain spike |
| ■ Cohn Fractionation | |
| ● Fukuoka GSS | endogenous |
| ● Hamster 263K scrapie | brain cells |
| ● Hamster 263K scrapie | endogenous |
| ● Factor VIII | hamster brain |
| ● IVIG | hamster brain |
| ■ Kistler-Nitschmann Fractionation | |
| ● Human albumin | hamster brain |
| ● Human albumin | fibrils |
| ● IVIG | hamster brain |
| ● IVIG | fibrils |

file: FDA TSEAC Clearance Studies 6-26,27-02.ppt

Clearance Studies 2:

- | | |
|------------------------------------|-------------|
| ■ Bovine collagen 1 | brain spike |
| ■ Bovine collagen 2 | brain spike |
| ■ Bovine gelatin | brain spike |
| ■ Filtration Asahi Planova filters | brain spike |
| ■ Filtration Asahi Planova filters | fibrils |
| ■ Irradiation Clearant Process | brain spike |
| ■ Others | |

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Process steps:

- Depth filtrations
- Membrane filtrations
- Phase separations
- Extractions
- Precipitations
- Column chromatography
- Thermal inactivation
- Irradiation
- Chemical inactivation
- Others

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Key Elements:

- **Scale down**
- **Agent**
- **Host**
- **Spike**
- **Processing**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Scale down:

- **Usually left to client**
 - Plug into existing viral validation protocol
- **Smaller scale is seldom better**
 - Surface affects dominate
 - More idiosyncratic
 - Less flexibility in sampling and assay
 - Proportionately larger sampling losses

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Agent:

- **Scrapie – many strains**
- **CJD – many strains**
- **BSE/v-CJD – one strain**

- **BSE/vCJD is the relevant strain for this agent**
- **Otherwise the choice is arbitrary**
 - **It is not clear that the variation within a disease class is any less than that between disease classes**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Host 1:

- **Mouse**
- **Tg Mice**
 - **If PrPSc is the agent – it is model of choice, but**
 - **The most effective tg mice are chimeras; carry the doppel gene; have random and multiple insertions, aberrant expression, in an abnormal context.**
 - **Tg gene replacement mice by Jean Manson – normal context but not necessarily fast or most susceptible.**
- **Hamster – convenient, well characterized**
- **Sheep – endemic scrapie, disappearance of susceptible genotypes**
- **Cow – should receive more direct study**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Host 2:

- BSE presents very differently in cattle and humans
- Hamster 263K scrapie is clinically more similar to BSE in cattle than BSE in the VM mouse.
- BSE in the VM mouse is clinically more similar to vCJD in humans than BSE in cattle.
- Even the selection of a host is somewhat arbitrary.

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Spike 1: Source of Problem

- TSE infectivity is polydisperse in its physical and chemical properties.
- Subpopulations with differing properties fractionate differently
- Steps that reliably remove most of the infectivity may never be able to remove all of the infectivity
- High titer infectivity is limited to CNS tissues.
- Monodisperse virus preparations like PPV can be introduced anywhere in a process and display the same properties as virus introduced early and surviving to that point.

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Spike 2: Brain derived spikes

- **High titer but questionable relevance**
- **Brain fractions:**
 - Microsomes
 - Liposomes
 - DLPCs
 - Caveolar domains
 - Fibrils
 - May behave better but are not necessarily more relevant
- **Appropriate if source tissue is brain or cross-contaminated by brain.**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Spike 3: Endogenous Infectivity

- **High relevance**
- **Low titer**
 - Have developed methods for accurate titer
- **Rodent tissues may not scale well**
 - Dura mater is fragile
 - Hearts, kidneys, pituitaries, livers, bones, tendons, etc. are small

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Spike 3: Introduction of the Spike

- **Intrinsic infectivity**
 - Blood - Infectivity in proportion to blood volume of tissue
 - Other intrinsic infectivity
 - Difficult to mimic intrinsic infectivity in solid tissues
 - Dura mater, tendons, hides, bone
 - Failure reduces to testing the spike in the presence of the tissue rather than in the tissue itself.
- **Extrinsic infectivity**
 - Cross-contamination
 - Brain may be relevant
 - May be the most significant source of infectivity
 - Spiking is the same as contamination

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Other issues

- **Sampling**
- **Sample preparation for titration**
- **Assay methods**
- **Control of cross-contamination**
- **Husbandry logistics and other issues**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Experimental Design 1

- **Wherever possible test endogenous infectivity first**
 - Highest possible relevance for most tissues
 - Tests removal of intrinsic infectivity
 - Low titer
 - Highly accurate methods for measurement
 - Carry the process forward as far as possible or until there is not possibility of remaining infectivity
 - The claimed efficiency of removal from exogenous spikes can not exceed that obtained for endogenous infectivity

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Example: Cohn Fractionation


- **Infectivity from endogenous infected blood could be detected as far as the Fraction $IV_1 + IV_4$ pellet**

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt


file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Experimental Design 2: Condition the Spike




file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Experimental Design 3: Condition the Spike



file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Experimental Design 4: Measure the product stream directly



file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Evaluation of Total Clearance 1.

- Endogenous studies take precedence over spikes
- Continuous processing takes precedence over step-wise values
- Great caution should in interpreting the cumulative removal from very similar steps tested step-wise.

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Evaluation of Total Clearance 2.

- The exercise is worthwhile
- High values are better than low ones
- However, the actual values should not be interpreted too literally

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

Convergence = Confidence

As data accumulate for multiple agents, spikes, assays, scale downs and animal models, convergence of diverse approaches on the same result provides the greatest security for the outcome.

file:FDA TSEAC Clearance Studies 6-26,27-02.ppt

