

FOOD AND DRUG ADMINISTRATION

WORKSHOP ON PRECLINICAL TESTING FOR ENDOVASCULAR

1076 GRAFTS 18 12 50

THURSDAY, JULY 29, 2004

The workshop came to order at 9:00 a.m. in the Grand Ballroom of the Hilton Washington, DC North, 620 Perry Parkway, Gaithersburg, MD. Dorothy B. Abel presiding.

Steering Committee:

Dorothy B. Abel
Marianne Grunwaldt
Angela C. Smith

Scientific Advisory Committee:

Mark M. Dehdashtian
Stuart T. Rodger
Louis J. Smith
Matthew S. Waninger, Ph.D.

PRESENTERS:

Timothy A.M. Chuter, M.D.
Louis J. Smith
Michael Brown
Kurt Liffman
Roy K. Greenberg, M.D.

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

2001 N-0463

TR4

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

I-N-D-E-X

Device Integrity Fatigue and Durability

Clinical Perspective,

Dr. Timothy A.M. Chuter

Scientific Perspective,

Dr. Louis J. Smith 20

Cook, Inc. 29

Kurt Liffman 32

Participant and Audience Participation 53

Clinical and Preclinical Performance:**Past, Present and Future**

Clinical Perspective,

Dr. Roy K. Greenberg 186

Regulatory/Scientific Perspective

Dorothy B. Abel 199

Participant and Audience Participation 211

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

P-R-O-C-E-E-D-I-N-G-S

9:08 a.m.

MS. ABEL: We'll do our usual formal introductions.

DR. CHUTER: Well, I was just going to leave this title as fatigue. Many of you were here last night. I thought I'd narrow this subject a bit and talk about the device.

Let's see where we go. All right. Okay. So there was a time when just about every stent graft that went in was destined to come out in pieces, and it was a very instructive time. There were examples. If you ever want to find an example of any form of failure of the stent graft, you just have to look at the old MinTec Vanguard experience and there will be plenty of pictures with everything you want to know.

I'm just going to plug through examples of failure of every part of the stent graft and every part of the stent graft has failed at some time in some patient with some device. Then I'm going to go through some of the factors that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 predispose towards failure. Obviously, a lot of the
2 progress that we've made has been in identifying the
3 failure modes, characterizing them, looking for the
4 causes, finding out what works, what doesn't work
5 and then doing a little bit of evolution to winnow
6 out the solutions.

7 So here's fabric failure. And you can
8 see this is a micrograph of fabric from a couple of
9 pieces that have been subject to repetitive impact
10 with an adjacent stent. This is what happens if the
11 flat surface of the stent is impinging on the
12 fabric. You get some fiber flattening, perhaps a
13 little bit of a defect, but nothing like the
14 horrendous defect that you get if the apex of the
15 stent is impinging on the graft. And it's just the
16 kind of recipe for disaster that you would imagine
17 where you have this soft pliable object moving under
18 hemodynamic forces against a far more resilient
19 object. And the answer that most people have found
20 to that particular problem is either to get more of
21 the stents in between the orificial sealing and
22 attachment stents, get rid those altogether, or to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 strap the stents down so that the potential for
2 movement between the fabric and the stents is
3 eliminated. It's kind of like, you know, putting on
4 seat belts. The problem is that the seat belts
5 themselves can do some damage. And those, too, can
6 hurt the fibers.

7 The underlying problem in the MinTec and
8 Vangard experience was the malalignment of the
9 stents that was caused by breakages of the sutures
10 that held them together. But they're not the only
11 ones where the sutures have broken. And if you look
12 at the AneuRx explant data you'll see that suture
13 breakages are very common. It's just that there's
14 sufficient redundancy in that system that the suture
15 breakages don't seem to have caused much of a
16 problem.

17 One of the problems they do seem to have
18 caused, at least early in the experience, is with a
19 more loosely woven version of the fabric is to
20 produce some suture hole leaks. The job of holding
21 the stent to the fabric can also injure that fabric,
22 as I said, and you don't always see a manifestation

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of it like this where you have sort of a water
2 sprinkler opening into the aneurysm. But what you
3 do see if you compare the rates of shrinkage between
4 a device like the AneuRex and some of the others as
5 has been recorded, is a profound difference in the
6 rates of aneurysm shrinkage. And you wonder what
7 that represents.

8 The excluder, obviously, is something of
9 a different phenomenon, perhaps on a more
10 microscopic level. But you wonder what that
11 represents in terms of sac pressures. Of course
12 both of those fabrics are being replaced. And in
13 light of the discussions that we had yesterday
14 relating sac pressure to migration force, it will be
15 interesting to see how those changes in fabric
16 impact the migration rate of those particular
17 devices. That's a very nice experiment that's going
18 to be going on right now.

19 One of the big sources of fatigue
20 failure with these devices has been in the skeletal
21 elements either the stents or longitudinal struts,
22 and all of the different materials have been prone

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 to breakage. But Nitinol in its early forms seems
2 to have been particularly prone.

3 If you looked at the old Nitinol basic
4 systems, the wires were black. And what that meant
5 is that they were covered with oxide. And it meant
6 that they hadn't been properly electropolished. And
7 if you looked at these microscopically there were
8 multiple surface defects which either became the
9 focus for a stress strain propagation of fractures
10 or you've got this kind of funky corrosive
11 phenomenon going on, perhaps also aided and abetted
12 by some of the repetitive stress and strain. So you
13 can go from this to this with an improperly treated
14 Nitinol wire.

15 And this, you would see fractures in all
16 parts of this system, particularly susceptible
17 elements in some of the systems just to call out a
18 few examples of the longitudinal struts. These are
19 in Talent, was particularly prone to breakage,
20 especially if it was on the inside of the curve,
21 longitudinal struts and the excluder. But not only
22 Nitinol stents, the stents of the bifurcation of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Zenith device, the original stents were fairly prone
2 to breakage.

3 And of course barbs. The prototype of
4 that was the barb fractures that we saw in the
5 Yanker device. And really, those barbs were a very
6 active securing mechanism. They were all that held
7 that device in place. So that was a bit of a
8 disaster for those. And the problem seems to have
9 been a manufacturing one where the radius of
10 curvature was just a little bit too tight on the
11 butt.

12 These barbs are attached in a totally
13 different way. They're passively deployed. They
14 point down. There is no acute curvature but still
15 they will come off.

16 The answer the Cook people seemed to
17 have found is redundancy, just a multiplicity of
18 these barbs. And if you look at one of the devices
19 when you first implant it in an angulated neck, it's
20 very easy to imagine how if you tilt this stent,
21 let's say you tilt it that way, the aorta is
22 angulated in the opposite direction, how you would

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 be loading up a couple of these. And, in fact, Bob
2 Skelter is very big on this concept where
3 implantation of this device into an angulated aorta
4 is often accompanied by some sort of settling over
5 the course of the next year or two as the load
6 shifts and the device becomes reoriented and spreads
7 that load out.

8 These sutures were also the site of
9 breakage, probably as much as anything evidenced by
10 the fact that the stent was doing something to keep
11 the device in place. Because once all those sutures
12 fractured, then the remainder of the device would
13 migrate distally. The answer to that, of course,
14 again was redundancy; more sutures.

15 And that's probably a pretty good
16 example
17 of the kind of testing that we were talking about
18 yesterday where, you know, you have a basis for
19 comparison in the old device, you can see how it
20 relates to the new device and you can extrapolate
21 that into a change in clinical performance.

22 So what is hurting these devices? It's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 not steady state forces, the kind of things that
2 have been in models and measured, and used to do
3 device testing. You know, it's pulsatile forces.
4 It's this kind of little dance that goes on time and
5 time again through the life of the device.

6 And we followed a bunch of these devices
7 at various stages from implantation through about 4½
8 years. And you can see some very interesting
9 changes going on.

10 As you look at the device the Bedford is
11 most prone to movement at all. There's two kinds of
12 movement. There's the translational movement where
13 the things moves up and down and from side to side,
14 and then there's a pulsatile movement where it
15 changes its diameter.

16 I don't know if we could get the room
17 lights down a bit so you can see this, because it's
18 quite subtle. Could that be possible? See, it's
19 magnified as much as I can.

20 If you look at the pulsatile expansion
21 contraction of this graft, when you first put these
22 devices in, the bit inside the aneurysm is pulsating

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the most. You'll see that wall flopping back and
2 forth, back and forth. And you put your hand on the
3 patient's belly and it's easy to understand why
4 there's a massive pulse there because everything is
5 moving. You can see here everything is moving.

6 This patient actually is a little longer
7 than that, I think it's a month or six months out,
8 and that is all quieting down. This is the location
9 of the stent that used to fracture. There used to
10 be one stent there and when you fluor the examples,
11 I've done a couple of those, you can see this stent
12 really is taking the heat right at the bifurcation.
13 Less so now that there are two.

14 But if you look at this segment now,
15 we're a month or six months out, the pulsatile
16 aspect of this, the in and out, has disappeared.
17 Why would that be the case? Well, I think it
18 relates to the pressure environment of the aneurysm,
19 and I think we're going to find more and more that
20 that's a very important factor in the durability of
21 these stent grafts is the behavior of the sac and
22 the pressure differentials between inside and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 outside.

2 If you imagine a situation where -- if
3 you measure these pressures, they don't have much of
4 a pulse pressure unless you got a Fran type 1 or
5 type 3 endoleak. The pulse pressure inside there,
6 it's fairly flat wave form. If that is above
7 diastolic pressure, you're going to have a phase in
8 the cardiac cycle where the pressure is bigger on
9 the outside than it is on the inside. And as it
10 cycles through there, you're going to have this
11 flapping back and forth. So systolally that's going
12 to got and distally it's going to come back. Now,
13 it's not going to move that much, because obviously
14 those pressure changes in the sac would be abolished
15 by any movement that they can generate in that graft
16 wall. That's the sort of the capacities chamber for
17 this whole thing.

18 And I suspect that's why these early
19 stent grafts are flapping around, because clinically
20 you put your hand on the patient's belly you will
21 see a corresponding decline in the pulse of the
22 aneurysm with a decline in the movement of that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 stent graft. So that disappears first.

2 The next one to disappear is here. If
3 you look at this neck very carefully, let's say you
4 look up there, you'll see that there's pulsatile
5 movement of this. It's somewhat obscured by the
6 translational movement, but it is going in and out.
7 And that is because compare the diameters here and
8 here. This is partially constrained. It can go in
9 and out. Obviously, it has the additional problem
10 of it's facing -- it's compliance at the stent graft
11 added to the -- well, it's compliance with the stent
12 graft aorta component there is what is influencing
13 the movement, but obviously there's some potential
14 for movement because it's not fully expanded.

15 When you look at here, there's no graft,
16 this one can move quite freely.

17 You follow these, this one is the next
18 to go. So this disappears first, then this
19 disappears. And the reason is this: Just about all
20 of these stent grafts start out looking like a
21 bottle, sometimes even more so than that. And we
22 followed a lot of these patients just on serial x-

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 rays comparing the diameter here to the diameter
2 here, which I call the oversizing index. And it's
3 sometimes shocking when I look at my numbers of how
4 I calculated the stent graft to find how much these
5 things were oversized. Because is the unconstrained
6 diameter and this is the constrained diameter.
7 Well, you'll see over time this is eliminated, and
8 it actually happens fairly quickly. By four years
9 these all look at like Coke cans. They're as
10 straight as can be. They are completely fully
11 dilated. And if you plop them out, you'll see that
12 the diameter rises and rises, and the gradient seems
13 to be fairly constant and then it just hits a
14 plateau. And where they hit the plateau; sometimes
15 they'll get to one-to-one. It's usually in the sort
16 of .9 to 1 ratio between this diameter and this
17 diameter, or this and this.

18 So there's not necessarily full
19 expansion, and I'm not entirely clear why, but
20 something in that stent graft is stopping the
21 further dilatation. Well, whatever it is, it's also
22 stops the pulsation. These things stop pulsating

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 once they get to look like that.

2 This, of course, is unconstrained. It
3 continues to pulsate.

4 And then what you end up this guy is
5 four years out. I apologize for the graininess, but
6 he weighs about 350 pounds. You've just got a
7 straight sided tube here. This thing is completely
8 dilated, there's no longer any pulse there. There's
9 not been any pulse there for a long time. And the
10 blood comes rushing straight down, impinges on the
11 bifurcation and this is the kind of movement that it
12 generates. And if you look at the changes in
13 translational movement over time, they don't go
14 away. In fact, there's a trend towards them
15 increasing with time. And probably the increase
16 relates to the increase in diameter of the inlet.
17 Those of you who model the forces on here will know
18 that that inlet diameter is a fairly strong
19 component in influencing the force on the stent
20 graft.

21 If you look down here these stent grafts
22 are implanted pretty close to the iliac bifurcation,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 which is a fairly fixed point in the arterial tree.
2 You've got the internal iliac coming off there, it
3 carries probably about here. It doesn't move very
4 much. This is moving a lot, this is moving not very
5 much. You can see what it's doing to that limb
6 there. And the only thing that gives me any comfort
7 there is that that limb is sitting inside the iliac
8 artery because I can see those sutures fairly
9 quickly wearing their way through the fabric.

10 So this is how we made the measurement.
11 From those semi-loops we would generate systolic and
12 diastolic images just from the extremes of the range
13 of movement, and we would plot the positions of
14 certain points on the stent graft and obtain
15 estimates of the translational movement and also the
16 pulsatile in terms of a percentage diameter change.
17 So you can see down here the pulsatile movement is
18 pretty much gone, the translational movement is
19 really quite striking.

20 And you just go back to some of the
21 analyses that have been done. You'll hear more about
22 them later. And the force is determined by the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 diameter up at the top there, the pressure, the
2 angle and I think sac pressure because I think what
3 we're looking at here is the difference between the
4 two.

5 Interestingly, if you look at the
6 migration rates, obviously this was not a pulsatile
7 analysis, but it correlated pretty well with the
8 findings from the Urista database in terms of
9 mitigation risk. Same factors: diameter pressure,
10 angle and interestingly type 2 endoleak seemed to
11 provide some protection against rupture and I
12 suspect that it's providing protection by the same
13 mechanism, sac pressure force, migration rights.
14 Migration being the primary predictor of rupture.

15 So how can we look into the forces a
16 little bit more? Well, what we've been doing is
17 taking CT scans from these people, segmenting them
18 out and doing some computational flow dynamics on
19 these. Go back to that one. And we can compare
20 them with the fluoroscopic movements. It doesn't
21 seem we can compare them in real time. But believe
22 me, they both move.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 This is just flow, but it's relatively
2 easy to map out both temporally and spacially
3 surface pressure, surface sheer, other factors in
4 the generation force, although pressure of course
5 predominates in all of that, and compare them with
6 the movements of the stent graft, the corresponding
7 movements of the stent graft. And interestingly,
8 you can do the same kind of analysis where you're
9 looking at flow and compare it to the CT scans and
10 you'll find that these points in recirculation tend
11 to end up as points in mirror thrombosis. There is
12 correlation between this. What we haven't done yet
13 is take that luminal profile and model that and see
14 if the formation of the thrombus has eliminated the
15 eddies and the spaces and so on that where
16 generating a thrombus.

17 So, to conclude, the pulsatile diameter
18 changes by location and time from implantation. The
19 bit in the aneurysm pulsates the most, but that
20 doesn't last very long. The bit in the neck seems
21 to pulsate probably for a couple of years, depending
22 upon the degree of oversizing and how long that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 oversizing persists, and then the pulsation of the
2 top stent just seems to go on forever.

3 Translational movement varies by
4 location in that the aortic elements seem to have
5 more movement than the more distal elements. But it
6 doesn't change with time. In fact, it seems to
7 increase if anything.

8 The only consequence that we could find
9 in these analyses of oversizing is that it perhaps
10 correlates with the extent of neck dilatation, and
11 we're not the first people to notice that. And by
12 that the mechanism may be increasing the hemodynamic
13 forces because there is some relationship between
14 neck dilatation and migration risk. And it seems
15 that based upon the temporal mapping of these forces
16 that they are primarily acting on the bifurcation,
17 which may become an issue for stent graft redesigns,
18 and I think perhaps accounts for some of the
19 astonishingly low rates of renal loss in the
20 fenestrated experience. With the fenestrated Zenith
21 where they're using a two component device. One
22 device is attached to the aorta and doesn't seem to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 move at all, and the other device is -- the other
2 part of that device is the bifurcation and that's
3 the bit that's getting all the heat. And it seems
4 that that sort of separation between the two
5 components functional and mechanical seems to have
6 some beneficial effects in terms of proximal
7 migration.

8 Thank you.

9 MS. ABEL: Well now I don't know what to
10 call Lou's talk, because that was pretty scientific.
11 You've got a challenge.

12 DR. SMITH: Well, thanks. Thanks, Dr.
13 Chuter. I hope I do justice in following you, but I
14 wanted to talk about the scientific perspective, the
15 device integrity, fatigue and durability. And
16 really it's all about how do you test for all of the
17 stuff that Dr. Chuter just talked about.

18 There is some testing in the ISO, I
19 wanted to just kind of go over that briefly. There
20 is characteristics in there that are addressed in
21 the testing, and then there are some limitations and
22 some other special considerations.

1 This is the lists of tests that you see
2 in the section of ISO 25539. The ones I've
3 highlighted in yellow are the ones the workshop
4 folks were asked to respond to, but there are
5 several others. And I just wanted to point out that
6 they're very interactive; you know, the strength of
7 the material, the corrosion, factory and anastomosis
8 strength is very similar to the strength of stent
9 and attachment systems in the graft. You know,
10 we've talked about that a little bit. A lot of
11 these things are interactive, so no one test -- my
12 point here is no one test covers everything. And I
13 think that's really important.

14 So I just want to jump into the big
15 ones. Fatigue testing. Obviously, Dr. Chuter
16 showed us the kind of forces in motion and movement
17 that we can see. There's a lot of failure modes that
18 can be identified by fatigue testing. Primarily
19 everyone's doing it to look at whether their stent
20 is going to fracture due to pulsatile motion. But
21 there's a lot of other stuff that goes on in vivo
22 and that we're trying to replicate on the bench,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 especially when the stents detached from fabric
2 where there's tearing or wearing of the fabric and
3 abrasion between overlapping components.

4 Typical fatigue tester. Just wanted to
5 have one up there in case folks don't know what they
6 look like.

7 There's a lot of types of fatigue
8 testing. In the standard we addressed pulsatile
9 fatigue, but there's also you know attachment method
10 fatigue; whether you're using anchors or barbs or
11 hooks, or whatever, and that could be a separate
12 fatigue test that could be done.

13 There is wear and migration. Often in
14 pulsatile fatigue testing we will generate forces
15 that are worse case in terms of pulsatility, but
16 those aren't always the worst case situation for
17 overlap component movement or abrasion or wear. So
18 you may have to do the test more than once. You may
19 have to do it at a highly oversized condition or
20 then maybe at a low oversized condition.

21 There's bending going on. You saw some
22 of that. And there are other compressive forces. You

1 saw Dr. Chuter's translational movement up and down.
2 That causes comprehension in some cases.

3 This is just an equation that's been
4 rearranged a little bit, but it's out of the ISO.
5 It talks about diametric deflection. It basically
6 speaks to how you can calculate that based on
7 compliance, percent compliance. And pulse pressure
8 you might see, the delta p. And it's a very
9 important thing to familiarize yourself with if
10 you're going to set up a pulsatile fatigue test.

11 I want to jump right into all the other
12 things. Corrosion. Versions 1 we saw some pictures
13 there. There are several types of corrosion that
14 the standard addresses. There's pitting corrosion,
15 galvanic corrosion, crevice corrosion. We can keep
16 going on with the list.

17 The examples that were shown are mostly
18 a surface finish issue and it can create pitting.
19 Galvanic, obviously, dissimilar metals. If you have
20 gold markers in touch with your stent material,
21 something you need to be concerned about.

22 And crevice corrosion where you get a

1 micro-environment being created, say, underneath
2 sutures or underneath bonding tape and how that's
3 going to effect your metallic components.

4 This is just a typical result of ASTM
5 F2129 where it's a potential dynamic curve. You can
6 see current density as a function of the potential
7 placed upon the sample. And these are like the
8 typical curves that you should be used to seeing.

9 Obviously, Dr. Chuter kind of showed you
10 the MinTec example, but corrosion can lead to the
11 picture here where you really should be looking at
12 devices after explant that look like the picture on
13 the other side of the screen, depending which way
14 you're facing.

15 In all of this, of course we're looking
16 at stresses and strains trying to determine what are
17 the loads. In the standard it talks about
18 calculating loads due to manufacturing, deployment
19 and in vivo conditions not just the in vivo
20 conditions. And what in vivo conditions are we
21 talking about? We can talk about pulsatile pressure
22 and flow, that's what we have been talking about.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1 There's also bending and translational motion. And
2 all of these analyses can be used to feed into how
3 you're going to do your durability testing; whether
4 it's just a diameter kind of test in a pulsatile
5 fatigue testing or bending, or whatever.

6 Obviously, the strength of how your
7 stent is attached to the graft is important.
8 Failure modes that are evaluated by this kind of
9 testing are pretty obvious, material tears or
10 sutures breaks. Basically the separation of the two
11 components which leads to many of the issues we've
12 seen.

13 There's standard characteristics that
14 have not been addressed in our testing, and that's
15 probably the primary reason why we're all here to
16 see, you know, about changes in the shape and
17 diameters, tortuosity, disease vessels and
18 angulation.

19 Tortuosity and angulation. That's seems
20 to be our mantra this weekend. I think it's real
21 important, it leads to bending, it leads to a whole
22 bunch of other forces. I do have a video here. It's

1 not as good as Dr. Chuter's. But, you know as it
2 runs you'll see the motion. This is mostly the
3 thoracic aorta, but you can see that motion being
4 translated down into the abdomen and it continue.
5 If you just focus on these bends. I guess you can't
6 see it. That's good for me. So we'll just stop that
7 little work.

8 But you can do fatigue testing out of
9 the pulsatile machine where you can see here we've
10 got some samples set up to be bent and moved and
11 under bending fatigue versus just your normal
12 pulsatile fatigue. Tortuosity and angulation leads
13 to this kind of motion. There's no way to get
14 around it and Dr. Chuter had some eloquent video
15 that I couldn't reproduce.

16 So, again, we got the same old monster
17 here. Limits of our testing are, you know, trying
18 to incorporate tortuosity, neck angulation, changes
19 after you implant. All these things are very
20 difficult to reproduce in testing.

21 A lot of the limitations are based on
22 what we've learned. It's difficult to put everything

1 into one test. I think we're getting better and
2 better at defining the clinical forces and what
3 tests to put together to at least study those, but
4 it's common sense after a while that it's not really
5 possible to assimilate all the known forces, at
6 least in one test. Trying to do flexing and bending
7 is difficult while you're also trying to do
8 pulsatile motion on the stent.

9 Of course, we've talked about the other
10 things, calcification and thrombus, changes in the
11 sac and the differential pressure across the wall.
12 All of these change, you know, effect the way we
13 want to determine our durability testing.

14 To try to determine worse case,
15 engineers are always looking at worst case. You
16 know, worst case delta p is where the sac pressure
17 is zero and these kinds of things to set up a test
18 that subjects devices to what they would consider to
19 be some of the worst conditions that they're going
20 to see in the clinic.

21 And I just wanted to touch real briefly
22 on special considerations for extenders and cuffs

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1 and stuff like that. It's not possible,
2 necessarily, to put all your components into the
3 same test. And like I said earlier, you may really
4 have to do different types of testing. Setting up a
5 fatigue test to understand how many times your stent
6 is going to fracture under just basic pulsatile
7 motion is a different game when you start
8 overlapping extenders or cuffs, or whatever you want
9 to call them.

10 At different oversizing conditions
11 there's different kind of abrasion that occurs. In
12 some conditions no abrasion can occur. But in the
13 same two devices together in another condition you
14 can get wear and abrasion and fabric pulls and what
15 have you.

16 When someone throws a Palmaz stent in on
17 top of a Nitinol stent or whatever, you could have
18 issue with galvanic corrosion. It hasn't really
19 shown up that much, but it's a potential that needs
20 to be examined.

21 Basically, in conclusion, there's
22 guidance for these things. If you look at the ISO

1 standards it not only in and of itself is a
2 guidance, but it points to a lot of other guidance
3 documents, whether they're ASTM methods or ISO
4 methods. And it's really an ongoing process.
5 There's never going to be one document that a
6 manufacturer can pick up and say this is all I have
7 to do. A manufacturer has really got to look at
8 their design, look at the target that they're going
9 for in terms of the clinical use and you may have to
10 come up with an array of tests. I mean, there may be
11 30 tests in that standard, or 34, or whatever the
12 number is but coming into submission, you've
13 probably got 50 different tests that you've
14 performed. So it's really not, you know, one
15 durability test or one fatigue test or one material
16 test that can cover everything.

17 Thanks.

18 MS. ABEL: We have a special guest
19 speaker.

20 PARTICIPANT (Cook, Inc.): I'd like to
21 introduce Kurt, who I work with. Kurt works with
22 CSIRO, which is an Australian research organization.

1 And we had a grant from the government,
2 a bit like your NIH grant because we wanted to look
3 at the effect of pulsatile flow. So we pulled a
4 pulsatile rig to see how these things work. And our
5 concern was how we build things inside an aneurysm.
6 So the simplest is the connection of the short leg.
7 But if we go for fenestration like Tim mentioned, we
8 separate the top from the bottom to separate the
9 forces out, but also to separate out the
10 orientation.

11 And when we come to the thoracic we
12 don't have that bifurcation situation at the bottom,
13 but we have a lot of force on the curve of the arch.
14 So now it's a different force mechanism that tends
15 to separate the modules. And the reason for having
16 modules is it's very hard to get a length
17 assessment. So to get a proper length assessment use
18 the trombone principle. And most of the times
19 people will underestimate the length.

20 So we wanted to know what was the stable
21 position of two pieces on a curve or what was the
22 force that was going to drive these things apart.

1 So we built the pulsatile module, the pulsatiler, a
2 flow situation which wasn't that easy to do in the
3 beginning.

4 We acquired an artificial heart to
5 discover that the artificial heart really doesn't
6 produce anything like the forces and pressures that
7 we need. So Kurt had to build a pressure pump. So I
8 just want to show you that and the effect.

9 And it's interesting that you always
10 think you're going to solve it and you find someone
11 else has worked it all out before you get there.
12 And So Tim's worked this out, and we discussed it
13 yesterday, and we have the new term that Lou said
14 was a differential pressure. And we thought of
15 pressure differential. So we really ought to use
16 that term from the inside of the systolic to the
17 sac.

18 And we found that the movement was quite
19 a lot, as Tim did, and especially in the model where
20 all the pressures are equal because he explained it:
21 You have a pulsatile pressure inside.

22 And then when the pressure differential

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 is such that it's great enough so the graft in fact
2 inflexs and it can come apart at that stage, but the
3 movement ceases.

4 And the point there is that if you look
5 on the -- now we can study this. If you look at it
6 you find that it's working all the time. So when
7 it's moving it's working. And that's a real
8 fatiguing motion.

9 It also I think -- I'd like to put
10 forward the concept that while we have to have tests
11 for many different things, we can set up the
12 condition and we can say well if there's only one
13 test, is it going to work in a patient? So we set
14 up the worst scenario. So we say, okay, now 10
15 centimeter thoracic aneurysm on the arch with a
16 pressure differential of this, and we set what we
17 ought to be, let's put it in there and see you meet
18 the standard.

19 So we'll see if our video flows. Give
20 it to Kurt.

21 MR. LIFFMAN: Right. I'd like to not
22 take too much time because I know time of the

1 essence.

2 When I first started talking to Michael
3 or working with Michael -- Michael you'll have to
4 come back here, you know. But we're looking at
5 forces on grafts. And so Michael asked me to work
6 out the equations for that, which we did and we
7 published a number of years ago. And as Michael
8 said, we want to -- we went through quite a few
9 machinations to produce the pulsatile flow rig, and
10 this is our current evolution. It's basically -- I
11 don't have it with me. But basically it's a
12 computer controlled unit and we can dial in any
13 particular pulsatile profile that we'd like.

14 Our motor is -- the pump is over there.
15 And you just basically put in your wave form and
16 you're fine. But we had problems with the pressure
17 and so -- because when you put in a motor which can
18 generate any wave form that you like, you still have
19 to worry how pressure propagates through the system.
20 And so we had to put a damper in. It's called a
21 windkessel. It's that device over there. And so
22 the pressure wave forms still nearly could work.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We want to measure the forces on the
2 grafts. Thank you very much. That's exactly what I
3 need.

4 So we want to work out the pressures on
5 the grafts. And this is our main units over here.
6 And what we constructed was a system like this.
7 This is an acrylic model of a symmetric bifurcated
8 graft. We have a load cell, which is connected up
9 to a computer. We just measure the deflections of
10 the wave, the load cell and that gets translated
11 into voltages which can then be translated into
12 pressures by calibration.

13 And we have flexible rubber membranes
14 just from rubber gloves. So attached you got
15 pressurized fluid, in this case water or water
16 glycerine mixture going through here. And we can
17 measure with a pulsatile motion what the forces
18 should be.

19 So the first case we looked at was the
20 steady state flow case. I haven't got the equations
21 here, but when I first arrived at them and published
22 them, there was some concern and still is some

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 concern in the medical community that we developed
2 equations which are for steady state, that is
3 continuous flow. But of course in the body it is
4 pulsatile flow and so then people are saying well
5 how realistic is this?

6 Well, first of all we looked at the
7 steady state flow. And the dots here, the measured
8 forces and the line here basically is from our
9 equation. This is assuming the initial -- in the
10 inlet diameter of the actual acrylic graft that I
11 showed you before, but when you actually pressurize
12 the system this joint goes out about a millimeter,
13 it expands. And so when you put that in, that
14 expansion in, you get a pretty good fit between
15 theory and what you measure.

16 So in this case we've validated the
17 pressure flow or the force flow equations that we
18 derived. But then we had to look at the pulsatile
19 flow equation. And when you start going to
20 pulsatile flow you get into quite complex
21 mathematics. And so I tried to simplify it as much
22 as possible by using a standard momentum equation,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 which were shown yesterday.

2 And so I took the bifurcated symmetric
3 graft and just looked at a pressure wave propagating
4 through the system. And you have this symmetric
5 graft. We have a body about -- LB is the length
6 there and you have the legs. And they're all at an
7 angle, alpha is the half angle. And when you go
8 through and see, I don't know if you want to see,
9 but here's the restraining force and here's what you
10 get to the steady state flow. If you just had
11 continuous steady state flow you get out this
12 equation. And as was shown yesterday, the dominate
13 term is this entrance term when you have the
14 pressure in the area. That's what really determines
15 the forces on the graft, because all the other terms
16 tend to be quite small. But when you put in
17 pulsatile flow and you start looking at it, you get
18 this extra term. And this is the density here of
19 your fluid. In this case it's going to be close to
20 water. And then you have actually the length of the
21 graft becomes important. All right.

22 So when you start worrying about

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 pulsatile forces in this very simple approximation,
2 the length of the graft starts to become important
3 and also the flow rate. This is the flow rate. So
4 this is the pressure as the flow is a function of
5 time.

6 But it turns it out that this term can
7 be neglected. It's very small. It's about the one
8 percent level compared to this term. So you're when
9 you're looking pulsatile flow, at least in this
10 theoretical sense, the steady state flow equation
11 which people have been so concerned, is actually the
12 appropriate one.

13 And these are some forces or PSB from
14 our initial experiments.

15 So these are experiments we've done.
16 Here is the pulse time just in one second. And this
17 is what we get from the experiment from our load
18 cell. This is a bit of hysteresis and time lag in
19 our load cell. And this is what we get from our
20 pressure equation.

21 Just basically by looking at this
22 particular terms. You could even neglect most of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 these and just look for P times A. And you get a
2 moderately good fit. We're trying to improve the
3 fit. So -- but that's basically where we are at the
4 moment with that, and we're doing some more
5 experiments.

6 Okay. Stability of modular grafts is
7 what Michael was talking about. Michael is very
8 concerned with modular graft system, about how
9 they're coming apart and what sort of pressures. And
10 we're still doing some ongoing work on that.

11 And we have our flow rate. Basically
12 this is our aneurysm. It's made out of acrylic and
13 we put in two modular grafts. And we have the
14 systemic pressure going through here. Typically we
15 go 130 on 80, pulsatile time of one second or so.
16 And we can change the pressure inside the aneurysm.
17 So we can look at how the behavior -- we can put the
18 grafts at different angles, different overlaps. So
19 there's some variability in the system but we can
20 tighten all the variables and so we can make sure
21 everything is consistent.

22 And this gives you an idea of the

1 movement that Michael is talking about. When you
2 have a main pressure difference, and I'll get him to
3 come up here soon and discuss that, but when you
4 have a main pressure difference of zero between the
5 graft -- here two grafts are basically joined
6 together. And we call two stents -- this is our
7 stent, unit of weight stint. So we have two stent
8 overlap. I think hopefully most people understands
9 what that means.

10 And there's a pressure difference of
11 zero between the aneurysm and the systemic pressure.
12 You get this movement. This shows -- this redline
13 shows the diastolic and when it's on systolic it's
14 up here. So it's going backwards and forwards.

15 When you have a pressure difference
16 which is equal to the pulse pressure, maybe you've
17 got 130 on 80, so the pulse pressure is 50 mil,
18 everything gets pressurized. So when the main
19 pressure difference is there, it just becomes
20 tightened up and it has actual structure and it
21 stops pulsating.

22 So hopefully this will work. WE shall

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 see and keep our figures crossed. It doesn't come
2 up.

3 MS. ABEL: All of you have to gather
4 around the laptop.

5 PARTICIPANT: Copy off that computer,
6 and just put it on the video.

7 MR. LIFFMAN: Sorry. Trouble with the
8 computer. I don't understand. Function F7.

9 You can come up later on and have a
10 movie. That's fine.

11 All right. But anyway, this is on the
12 side view. It's the same sort of thing. You get
13 this pulsatile system and there's zero millimeters
14 difference. And when it's fully pressurized, you
15 get the idea. Otherwise it's just pulsating
16 backwards and forwards and you can see the movement
17 between the stent systems. And same sort of thing
18 over here as well.

19 Okay. And that's the end of the
20 presentation.

21 PARTICIPANT (Cook, Inc.): You got
22 bitten by the technology a little bit. Need a bit

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of preclinical testing there, I think.

2 What the video showed was that the graft
3 did come apart. And one stent overlap was very
4 unstable. And at two stents -- two stents overlap
5 it would part would you got to the pulse pressure.
6 And you needed three stents overlap to sustain a
7 pulse -- a pressure differential of 50 millimeters.
8 So that transcribes to IFUs and a testing situation.
9 You know that if you haven't got at least two stents
10 overlap than you have a situation where it could
11 come apart.

12 It won't always come apart and people
13 will say in a clinical situation how come it doesn't
14 come apart? It won't come apart if it gets to the
15 wall before it comes apart. So that's why in a big
16 aneurysm it's more likely to come apart in a smaller
17 aneurysm because of the distance of travel. So the
18 reason for having a long overlap if you are using
19 modular grafts that can move is to allow for a
20 certain amount of travel as well.

21 MR. CARDELLA: While they're setting
22 that up, can I ask a question again? When you say

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the pressure difference is zero --

2 PARTICIPANT (Cook, Inc.): Well, you can
3 just see it moving. This is when the systemic
4 pressure inside is the same as outside. And this is
5 when the pulse -- when it's down to 18, the pressure
6 differential is 18. So the pressure differential is
7 18 and there's still some movement. And then it
8 goes up to 50 and it's essentially fully inflated,
9 and that's why you don't see the movement.

10 If we go back and you watch the join --
11 will it go back? But if you watch this part here
12 when it goes around again, you see here. You can
13 see this is actually working. So you just wonder
14 sometimes whether you start to get a leak if you
15 don't get a proper seal.

16 MR. CARDELLA: So a question again, in
17 the upper left hand panel before you take that down,
18 when you say the pressure difference is zero does
19 that mean that you are helping 51 millimeters of
20 pressure inside the graft and 51 millimeters of
21 pressure in the sac chamber in that aneurysm sac
22 chamber. That would be the difference of zero.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 That looks to me like you'd have a high pressure
2 inside the graft and zero outside the graft. So it's
3 not a pressure difference. You're just saying the
4 pressure in the sac on this picture -- well now on
5 this picture. The pressure in the sac is 18.

6 COOK, INC. The pressure -- the pressure
7 differential was zero between the systolic and the
8 pressure in the sac in the first one. The pressure
9 differential here is the pressure differential
10 between the systolic inside and the pressure in the
11 sac. And what Tim said there was very little pulse
12 pressure within the sac. So we have a -- and may
13 have a pressure differential between systolic and
14 the pressure in the sac of 50. And it's essentially
15 inflated because it's above diastolic. The pressure
16 differential is above the pulse pressure. Pressure
17 differential is above the pulse pressure it's
18 stabilized.

19 MR. ELLER: Just a question here. Was
20 the stent graft sealed in this experiment?

21 MR. LIFFMAN: No, they weren't sealed.
22 We just put them el natural and we just had water

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 going through here. So there is some leakage going
2 through the grafts. We just wanted to make sure we
3 didn't interfere with the friction and we just had
4 water and maybe next time we'll put glycerol and
5 that apparently stops leakage.

6 PARTICIPANT (Cook, Inc.): It seems the
7 same question because you have outflow from the
8 bottom. So there's differences, and like if you
9 have a lot of outflow form the bottom, a little
10 outflow through the graft, you've still got outflow.
11 So it's really the -- it's the pressure
12 differentials that matter.

13 MR. LIFFMAN: So Michael, what was the
14 pulse pressure on the inside?

15 PARTICIPANT (Cook, Inc.): It's 130 to
16 80. So as someone mentioned, maybe that's a bit
17 low. If we're going to test, maybe we should test
18 when someone's running upstairs or something and
19 their blood pressure is 220 or something like that
20 so we can have a pulse pressure of 200 on 100, which
21 gives a pulse pressure of 100. And as Michael said
22 yesterday, maybe we should test for a pressure

1 differential of a 100.

2 DR. CHUTER: Could I just ask Kurt a
3 question?

4 Kurt, you said that the -- when you ran
5 the equation with the pulsatile flow that the
6 pulsatile element really was quite trivial compared
7 to the pressure elements. Did you run the equations
8 with pulsatile pressure as opposed to pulsatile
9 flow?

10 MR. LIFFMAN: What I'm talking about
11 when I'm saying pulsatile flows is I mean the flow
12 and the pressure. So the steady state flow
13 equations you can also just -- the pressure can be a
14 function of time. So it can be pulsatile as well.

15 DR. CHUTER: I see.

16 MR. LIFFMAN: All I'm saying is that
17 when you look at the equation that's been shown
18 yesterday, the momentum equation fits into two
19 terms. One it's pressure time dependent, it's a
20 volume integral. And the other one is a surface
21 integral. But the one with the steady state flow is
22 a surface effect, the surface forces that we're

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 looking at. And that's where the pressure and the
2 area comes in. When you worry about the volume and
3 the time part, that's what we're looking at in that
4 extra bit of that equation. That's where that
5 pulsatile flow comes in, the flow through the graft.
6 It turns out it's very small.

7 If that term were to be of the same
8 magnitude as the pressure area term, you'd need a
9 pulse rate of 100 times per second. That's how
10 small it is. And as a pulse rate is about one per
11 second, you can basically neglect.

12 And so the steady state force equations
13 that we've been using appear to be okay, at least
14 for the approximation we've looked at for if you
15 just make the pressure change with time.

16 DR. CHUTER: All right. I'm still
17 puzzled. Let me ask you another question.

18 MR. LIFFMAN: Okay.

19 DR. CHUTER: Okay. If you neglect the
20 flow related effects is there a temporal change in
21 the pressure related effect according to the
22 pressure at that particular instance? So, as the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 pressure varies through the cardiac cycle, does the
2 force vary through the cardiac cycle?

3 MR. LIFFMAN: Yes, it does. It's
4 directly proportionate.

5 DR. CHUTER: So there is pulsatile
6 variation in the force, it's just that the flow
7 related element of that is small?

8 MR. LIFFMAN: That's right.

9 DR. CHUTER: Okay.

10 PARTICIPANT (Cook, Inc.): There were a
11 couple of things on that -- on a chemical side if
12 you -- if you do a fluoroscopy and you see your
13 graft is moving, then you suspect you have a
14 significant endoleak just by inference.

15 And the second thing is that very early
16 on we put a spine along the graft. And we took it
17 out for -- I can't remember the exact reasons. But
18 I just wonder whether -- the reason I never
19 understood proximal migration at the top end was
20 because there is a certain amount of -- continuing
21 there. And if you don't allow it to do that, maybe
22 you can push it upwards. I'm not sure about that.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 So maybe the rigidity of the graft is a factor.

2 MR. SMITH: I have a question for Kurt
3 as well. You know, you're putting this graft into
4 an acrylic tube or an acrylic aneurysm. How much of
5 that do you think dampens the effect of what's going
6 on? Is there any issue with wave reflection and
7 things that we're not considering here?

8 PARTICIPANT (Cook, Inc.): There may be
9 a very important issue with wave reflection. In all
10 the analysis I've done I haven't looked at the
11 reflection of waves. And what you can have happen
12 in the graft is you can have the waves interacting
13 such they either reenforce or they negate. And so
14 you might have greater pressure differentials.

15 The way I've tried to look at the
16 problem is I go step from step. And so we do
17 something simple first, do steady state flow. The
18 next one is well we don't worry about the pressure
19 waves. Let's just look at the momentum equation,
20 look at this extra time then in turn. And then the
21 next one is going we worry about the pressure waves.
22 And there's sure to be important facts there, in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 particular with restenosis probably, you know, at
2 the ends of the graft. Because the graft's
3 structure is different from the artery and you've
4 got all of this sudden change within the artery. So
5 there are going to be pressure changes there, so
6 that's definitely going to effect the environment of
7 the artery, so it would be important.

8 DR. CHUTER: I got another question for
9 Kurt. As you got this thing pulsating and it's
10 dilating and contracting as it goes through the
11 cardiac cycle, obviously that's driving fluid in and
12 out of the space around the graft. Presumably your
13 pressure control mechanisms were capable of
14 eliminating the pressure variations that that would
15 have put on that space --

16 MR. LIFFMAN: No, but it's an important
17 question because what's going to happen is the graft
18 expands and contracts, it's going to change the
19 pressure within the sac.

20 DR. CHUTER: Exactly.

21 MR. LIFFMAN: But all our analyses
22 assumes a rigid tube. I mean, that's standard

1 analyses. That's an approximation. So the next
2 step is also we have to look at the -- in the graft.
3 And that will change the pressure in the sac, as
4 you're explaining. If the graft were perfect like a
5 concrete tube, there'd be no pressure transfer into
6 the sac.

7 DR. CHUTER: Right. But it's not.

8 MR. LIFFMAN: But because it expands it
9 contracts, that's where you get the pressure
10 transfer.

11 DR. CHUTER: Were you controlling the
12 pressure in the sac or the volume in the sac,
13 though?

14 MR. LIFFMAN: In the sac we were
15 controlling the pressure. We just had its --

16 DR. CHUTER: So it was an infinitely
17 compliant chamber?

18 MR. LIFFMAN: Basically, yes.

19 DR. CHUTER: Okay. So if you were to
20 mess the compliance of that chamber, i.e., the way
21 that pressure and volume were related in that
22 chamber, you could probably try to mimic the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 compliance of an aneurysm, and that would give you a
2 closer approximation to sustain --

3 MR. LIFFMAN: Absolutely correct. So
4 one of the next things we're going to do is the
5 latex aneurysm if we can get the funding. Hint.
6 Hint. We're going to do a latex aneurysm and just
7 measure everything in there.

8 You're absolutely right. Again, it's
9 the first approximation.

10 DR. CHUTER: Okay.

11 DR. FILLINGER: If you have -- in your
12 aneurysm sac though, in a sense you've got -- I mean
13 you have a place for that fluid to go when it
14 expands, so it's really not very dissimilar to what
15 you're already doing unless, of course, you have the
16 instance where all of them are thrombosed. But
17 that's not that common.

18 MR. LIFFMAN: You're very kind. But in
19 your particular case that you're talking about, you
20 don't guarantee that the pressure in the sac is
21 always the same. We've sort of guaranteed that the
22 main pressure is always the same.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MS. ABEL: Very well. I think we better
2 move on and start talking about the other things we
3 have on the side.

4 I think there are very good
5 presentations this morning that really helped to
6 demonstrate that this is a very complex situation
7 that we're trying to deal with in terms of all the
8 forces that are actually on the grafts.

9 I just went to mention a couple of
10 things with respect to the compiled work assignment.
11 I keep losing my point. Sorry.

12 In this side I just wanted to point out,
13 again we really can't calculate any two compile
14 rates based on the information that we provided.
15 But it is interesting to see that people who did
16 observations with respect to loss of graft integrity
17 or suture integrity where anywhere from 1 to 27.9
18 percent. So there are some devices that have had
19 some pretty significant issues with respect to these
20 things.

21 And also if you look at the identified
22 fractures, and these are from the clinical studies

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 that people reported to us, here I think it's
2 interesting to see when they happened. It looks
3 like most of what's going on is between 3 and 12
4 months, which gives credence to FDA's requiring a 12
5 month data. At least we're capturing that
6 information within our clinical studies.

7 We also asked people to tell us about
8 their explant analyses, and we had a total of 329
9 explants reported. In this table it's a bit busy,
10 but again if you just look at kind of the grouping
11 of where most of the observations occurred for an
12 awful lot of stent fractures are occurring in the 12
13 to 30 month time frame -- now we can't be sure that
14 that isn't going to continue in terms of those
15 numbers because we don't know the number of patients
16 at risk farther out. But there's a pretty
17 significant number of factors that have been
18 identified. And this is there were five respondents
19 that had observations of fractures and only two that
20 said that they didn't have any at all. And it's
21 interesting, too, to see that the hooks and barbs
22 only have one case identified.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 And so, you know, this obviously isn't
2 totally reflective of what we know is happening in
3 the clinic, but this from the explants.

4 The additional explant observations
5 we've got graft material there. Again, five people
6 with observations, two don't. And a lot of that is
7 going on in the 12 to 30 months time frame.

8 So I'm not going to spend any time to
9 talk about the information that people sent back to
10 us because I think we've got a good basis based on
11 the prior presentations. Yes. So I think we'd be
12 best to not break yet but go on to the slides that
13 we want to be discussing.

14 MS. SMITH: And I think, like yesterday,
15 we want to look at whether the -- that have been
16 seen and that were identified by respondent can
17 actually be evaluated in a fatigue and variability
18 test. And the other thing that we wanted to
19 evaluate was if there are things that are seen in
20 the testing but not in clinical use and can be
21 defined as testing artifacts.

22 MS. ABEL: So ISO says we're looking for

1 stent fractures and most people who responded said
2 they did, but not everyone which was kind of
3 interesting. Because I would assume that the focus
4 of a fatigue durability test would be to look at
5 device integrity. Would everyone agree? The answer
6 should be yes.

7 And we wanted to talk briefly about what
8 you do in terms of testing artifacts. And we have a
9 lot of books, quite honestly, who say we always saw
10 these things falling apart but it's because the
11 testing is too severe or, you know, it's worst case,
12 etcetera, etcetera. So how likely is it truly that
13 you're developing a test that is so rigorous that if
14 you see a stent fracture that it wouldn't be a
15 realistic observation? You now, if you think you
16 would see fractures in your bench test in your
17 fatigue test and that would not be an indicator or
18 potential to see that in the clinic?

19 MR. VASUTEK: Well, you test to
20 destruction. Surely you're going to push, so you're
21 going to see failures but a higher level of
22 compliance or whatever than you would expect to see

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in vivo.

2 MS. ABEL: So I think that's a very good
3 point. You know, I think most of the testing that's
4 been done so far hasn't been necessarily, you know,
5 to destruction hasn't been to failure. It's been,
6 you know, trying to simulate reality. And,
7 obviously, we're not very good at that right now.
8 But even in an attempt to simulate reality would you
9 expect to see failures or be able to explain them
10 away?

11 MR. VASUTEK: Well, if you're simulating
12 reality, then you would not expect to see failures.

13 MS. ABEL: But you're attempting to
14 simulate.

15 MR. VASUTEK: Okay. If you're
16 attempting to simulate, you wouldn't see failures.
17 But if you're taking it beyond what you think is
18 reality, then you would expect to see failures. So
19 for example, if we try and simulate in bench testing
20 or screening, we've doubled the loads when we get a
21 failure in half the time or if we double the motion
22 from the testing, will we get a failure in half the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 time as well.

2 So when we use, I guess, conditions that
3 are accelerative, then we produce failures in a
4 relatively short period of time. Are they
5 clinically relevant and how do we interpret those?

6 DR. WHITE: There's a problem with that
7 accelerated testing which is probably hierarchy, but
8 I mean the 10 ten cycles and accelerated are for
9 straight segment stent tests and it's not predicted
10 to failures which have been at angles or transition
11 zones clinically. So my observation would be is
12 that accelerated testing would give you artifacts if
13 you do get it in a straight model, and it ought to
14 physiologic in the angles if you're trying to
15 predict the event.

16 MS. ABEL: Dan, did you have something
17 you wanted to say?

18 MR. CHWIRUT: Yes. I think there's two
19 items that you can look at in trying to determine if
20 the failure in the bench test is artifactual. One
21 is where it occurs and second is when it occurs.
22 Obviously, if you put something in, put it on test

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 and then the first, you know, simulated six months
2 or something like that something breaks, probably
3 something was a amiss. It was deployed improperly,
4 the device was not representative of clinical
5 quality or something like that.

6 The other thing is you look at the
7 comparison between your analysis and your test. And
8 if you've done an FEA or some other stress analysis
9 of the test conditions that you're trying to
10 produce, and you're getting failures that are
11 totally at odds with what your analysis tells you,
12 again something is amiss and you might look at that
13 failure as being artifactual.

14 MS. ABEL: Thank you.

15 Okay. We can move on to the failure
16 mode state, graft fatigue. Have not shown what the
17 difference is between fatigue and fracture other
18 than maybe it just gets tired, it doesn't break all
19 the way. So we're just going to kind of say -- yes,
20 more people actually looked at fatigue an fracture.

21 Well, that's actually separate. The
22 fabric we look at separately. This is with respect

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to the stent, but maybe that's what people -- maybe
2 they grouped it together.

3 Barb and hook fractures, are these --you
4 know, again, just the straight durability tests that
5 we're talking about described in the ISO standard.
6 Is this test really something that you can use to
7 evaluate barb or hook fractures or is that a
8 separate test that you need to design?

9 PARTICIPANT: Are you talking about just
10 a radial fatigue test or are you talking about
11 longitudinal fatigue or multiple possibilities?
12 This table is geared just towards a pulsatile
13 fatigue, correct?

14 MS. ABEL: Yes. So does the pulsatile
15 test tell you anything about your hooks and barbs?

16 PARTICIPANT (Cook, Inc.): I think it
17 depends on the test set up. We have test -- fatigue
18 test which will test longitudinal fatigue as well.

19 MS. ABEL: So you incorporate
20 longitudinal forces?

21 MR. SMITH: My experience, we've had to
22 develop a separate test to specifically apply worse

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 case forces to anchors.

2 MR. RODGER: Yes, we also have developed
3 a separate test.

4 MR. CARDELLA: Whether it's a separate
5 test or not, I think it should be tested if that's
6 the issue.

7 MS. ABEL: We're just talking about in
8 this test. And so if you were to say that it should
9 be evaluated in this test, obviously the test
10 grammars would need to be modified so that you're
11 adequately doing that. But if it needs to be
12 evaluated in a separate test, we don't have to try
13 to incorporate anything that would effect those
14 hooks and barbs. Okay.

15 And this is where we're talking about
16 the tearing, other failures of the fabric. It is
17 mentioned in the ISO standard as something that you
18 would be looking for. Only half the respondents were
19 looking for it in their tests.

20 Is this a reasonable test to use to look
21 at graft wear and tear and that sort of thing or are
22 the conditions being, you know, inside of the mock

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 artery and just not simulating enough in terms of
2 that particular failure mode?

3 Where's Frank? Frank, get to a
4 microphone.

5 PARTICIPANT: I would say you don't use
6 that test to determine that because you don't have
7 the movement that is -- you should have a second
8 test to determine friction between the two different
9 components. You can get misled.

10 MS. ABEL: For those of you who don't
11 know, Frank -- can I say it?

12 PARTICIPANT: Sure. Why not.

13 MS. ABEL: Worked on the Vanguard
14 project.

15 So I think it's certainly something that
16 you would document if you did see an observation.
17 But the test isn't specifically designed and
18 adequate to look for that type of a problem?

19 PARTICIPANT: I think I would say again
20 it depends on the design of your set up.

21 MS. ABEL: Yes. And you must have a
22 respectable one that's different -- you should have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 filled out your little homework and then maybe we
2 could have --

3 PARTICIPANT: The dog ate it
4 unfortunately,

5 DR. MARIN: Dorothy, I would add that
6 I'm not clear exactly what you mean by other
7 failures of the fabric. The fabrics, as you know,
8 can restore before failure. And I guess where I
9 think this test can be useful is identifying areas
10 of distortion or movement of yarns that can create
11 openings and so on that aren't necessarily failure
12 points but suggest modificational changes in the
13 structure during exposure to this type of durability
14 test.

15 MS. ABEL: I think that's a very fair
16 observation. But again, I think what we're saying
17 is you need to design a separate test to look at,
18 like you say, the changes in the fabric or the
19 material. Tailoring it or wear or you know.

20 DR. MARIN: Distortion.

21 The point about fabrics is that they're
22 usually, you know, the yarns are laying at the right

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 angle to each other and often under particularly
2 nonuniform loading, you get distortion of the yarns
3 -- no longer at right angles to others. So you can
4 identify -- you can use the structure to identify
5 places of stress concentration

6 MS. ABEL: Sure.

7 DR. MARIN: So I don't know if you want
8 to add this to a specific others at the bottom.
9 This is not particularly considered failure mode,
10 but it is an observation of changes in the
11 structure.

12 MS. ABEL: I think we'll just put it in
13 with that and acknowledge again that that's
14 something that should be observed in a test.
15 Unfortunately, you can't really -- this test isn't
16 designed to look at the -- or to challenge the
17 fabric appropriately for the interaction of the
18 fabric with the stent.

19 Okay. Detachment of the stent from the
20 fabric. So the suture breaks. Is this the sort of
21 test that helps you to identify that or do you need
22 a separate one?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MR. WANINGER: I think that's a
2 situation where you're going to document your
3 observations, but I think you're going to want a
4 separate test. At least in our experience we
5 developed another one.

6 DR. MATSUMURA: I think when you're
7 talking about that testing, too, it's important not
8 just to test the device itself, but to test it after
9 it's been loaded and the delivery system deployed in
10 anatomy, you might expect, and then treat it
11 sometimes by the physician in ways that you might
12 not expect.

13 MS. ABEL: So you should bring in some
14 physicians to deploy the devices in your testing in
15 your marked arteries.

16 In the standard it does specify that the
17 device should be manufactured and loaded and
18 everything in accordance with the IFU and then
19 deployed in accordance with the IFU.

20 DR. MATSUMURA: Well, just to emphasize
21 the point that Rod and Tim made, on curved anatomy.
22 I think the big thing that we miss on this

1 preclinical testing was on curved anatomy so when
2 you deploy a graft encurved anatomy and then if you
3 apply a balloon on the inside of that, that may do
4 something different because it translates the load
5 to one area of the graft radially as opposed to the
6 whole circumstance which on curved anatomy is
7 testing.

8 MS. ABEL: Yes. So we'll need to get
9 into that when we talk about the specific
10 modifications and things to consider in this test.

11 Migration, is this something that you
12 can really look -- evaluate in this type of a test?
13 I know you guys can evaluate everything. Anyone
14 else, do you think you really can look for
15 migration?

16 MR. SMITH: We did a very similar test
17 but modified to make it worse case for potential
18 migration.

19 MS. ABEL: So other than what we just
20 talked about was pretty much you're looking for some
21 sort of a fatigue or a fracture of the stent
22 material or the attachment system, whatever you want

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to call it. Is that really the only thing that you
2 can look at in this test other than documenting
3 observations?

4 Dan? Sorry. Okay. Jim. He's giving
5 you the floor.

6 PARTICIPANT: Well, we have seen testing
7 in the past that has shown --we have seen testing
8 results in the past that have shown fabric abrasions
9 similar to clinical results, the suture breaking and
10 certain stent migration. I think a lot of it has to
11 do with the protocols that one uses in their
12 testing. And over the years we've made a lot of
13 strides in the standards committees trying to focus
14 in on that.

15 Sometimes we forget that the earlier
16 data that's been prepared was not prepared with the
17 most recent protocols.

18 And in those same meetings we have had
19 statements from manufacturers who have indicated
20 that there are protocols that do predict the number
21 of cycles and the location of failure once those
22 products go into the clinical setting.

1 So we need to remember that we're
2 evolving in our testing and can't get stuck on the
3 interpretation of all the data that seems to be
4 unpredicted because there were certain protocol
5 flaws associated with them.

6 MS. ABEL: Well, what's we're talking
7 about right now is, like you say, the evolution.
8 Should we be changing these tests so that we can
9 better evaluate these different types of failure
10 modes and what the majority of the people are
11 saying, at least what I've heard so far, is that in
12 one group because of the way that the tests are set
13 up you aren't really capable of truly looking at
14 whether you're going to have wear and abrasion. So
15 it's interesting that --

16 PARTICIPANT: It's historical. I mean,
17 that's true. That's historical information.

18 I think we have protocols at the various
19 standards committees that we work on together that
20 have addressed these issues. And it's possible that
21 if the more modern protocols that we're coming up
22 with we used throughout history, then we wouldn't

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 have so much of the concern and the unpredictability
2 of the past testing.

3 MS. ABEL: Could some of the other folks
4 who have been in the standards meetings respond?

5 PARTICIPANT (Medtronic): One of the
6 things we discussed last time is that some of those
7 things we cannot study because the grafts not
8 pressurized. And also, there's a lot of
9 nonphysiological graft material movement in that
10 test at the accelerated testing speed.

11 So, just want to maybe ask the question
12 whether the audience thinks that we should do
13 something about trying to pressurize the graft from
14 inside?

15 PARTICIPANT (Lombard Medical): We think
16 it's a good idea to do that. We don't want to be
17 too coy, but there are a number of really quite
18 radical problems with traditional fatigue testing.
19 And principally you tend to use superphysiological
20 pressures because you have to have to have a very
21 stiff rubber tube in order to get the durability,
22 which means that you get the correct movement but

1 the pressures and forces involved inside the stent
2 graft tend to be different. Very often those sorts
3 of models emit an aneurysm which is quite a major
4 omission for a stent graft and of course are very
5 difficult to put angulated necks in. And we've done
6 that both to try to come up with a system which
7 addresses those factors and is quite a bit more
8 representative of the physiological situations. I
9 would encourage people to move in that direction.

10 MS. ABEL: I know we're in the midst of
11 a hot discussion here, but I think maybe we'd better
12 take a break because I see quite a few are needing
13 to get up, and I hate to be rude and not let Angie
14 have a break, too. So if we can come back in ten
15 minutes, we'll resume this discussion.

16 (Whereupon, at 10:21 a.m. a recess until
17 10:45 a.m.)

18 MS. ABEL: Welcome back, everybody.

19 I think if I remember right, you were
20 standing there. Jim or Joe -- Jim, you're standing
21 there.

22 We need to start with a little

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 announcement with respect to transportation.

2 MS. SMITH: I think the hotel actually
3 had a sign-up sheet out on our registration table
4 for those who need transportation from the hotel
5 this afternoon or this evening. So if you do need
6 transportation, sign up, I guess as soon as
7 possible. Lunch, something like that. I think
8 that's about it.

9 MS. ABEL: All right. Now that we have
10 our housekeeping out of the way, Dan, do you
11 remember what it was you wanted to say?

12 MR. CHWIRUT: I believe so.

13 Under the last row there specify others.
14 I know you've got a totally separate section on
15 corrosion, but I want to ask a question now is this
16 test, the radial dilation pulsatile whole device
17 durability test an appropriate way to test threading
18 and galvanic corrosion for these devices if the
19 appropriate environment is specified? So could they
20 be put in as additional failure modes that can be
21 assessed by this test? My opinion is yes.

22 MS. ABEL: The question that you missed,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Lou, whether corrosion could be a part of this test.

2 MR. SMITH: Threading and galvanic
3 corrosion.

4 MR. CHWIRUT: If everybody agrees yes on
5 this test, we can cross those off the discussion on
6 corrosion.

7 MS. ABEL: Yes, right.

8 PARTICIPANT: I think it's still
9 replicated to the test, corrosion.

10 MS. ABEL: Anyone else?

11 PARTICIPANT (Medtronic): I have a
12 concern with respect to corrosion because this
13 testing often induces non-physiological type of
14 contact due to accelerated testing speed which may
15 give you some artifacts.

16 DR. FOGARTY: We could slow the testing
17 down.

18 MS. ABEL: Yes. Yes.

19 MR. SMITH: I think corrosion
20 evaluations take combinations of tests, though. I
21 mean, in this you can look at and say okay, yes, did
22 I get any or whatever. But that's no different than

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 really flowing the metal in a bucket of saline for a
2 certain amount of time, which was the old way of
3 doing it. So, you know, just the interactions with
4 the FDA and all other, and ASTM and all, there's
5 been a proof toward this potential dynamic testing
6 to kind of get a relative measure compared to other
7 materials.

8 So I don't know -- even though you can
9 look at it in this test, I don't know if the other
10 types of evaluations go away.

11 MS. ABEL: So we'll get into that when
12 we get to corrosion segment of the workshop. Thanks.

13 Jim, you had something else you wanted
14 to say.

15 DR. CONTI: No, I'm okay.

16 MS. ABEL: You're okay now? Actually,
17 he's not okay because he was asking me what it is
18 we're trying to accomplish here. And I was saying
19 what we're trying to figure out is, you know,
20 exactly what does this test tell us now, are there
21 things that we can do to make it tell us more or,
22 you know, or do we just live with the fact that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 there are different tests. And from what I've heard
2 this morning, I think there are two ways of going
3 about looking at what you need to look at. And I
4 don't know that we know yet that when you come up
5 with a more complicated model that our friends have
6 done across the pond, whether that's going to give
7 you all the information that you need. Because, you
8 know, can you just tell us how much clinical
9 experience you have, or you don't want to tell us?

10 PARTICIPANT: We have 36 cases implanted
11 with a maximum plantation time of two years.

12 MS. ABEL: So you have at least some
13 patients out longer term.

14 PARTICIPANT: Yes.

15 MS. ABEL: But a relatively small number
16 of patients. So I guess, you know, time will tell
17 if you've been able to predict. Although may you
18 just have the best device in the world and it
19 doesn't matter, you didn't even need to test it.
20 I'm not sure the new test has been validated
21 necessarily, but there are two ways of going about
22 it. You can come up with new testing that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 incorporates a lot of the things and, hopefully,
2 you'll be willing to share with us again everything
3 that you've incorporated into your tests or you can
4 do the individual testing tests to look at the
5 various aspects, the various failures that we're
6 talking about.

7 PARTICIPANT: The other thing we've been
8 doing is testing competitive devices and seeing how
9 they check out in comparison with the clinical
10 experiences.

11 MS. ABEL: Sure. So you'd better do
12 that to validate your test.

13 And I had asked Jim if he wanted to tell
14 us what he believes the state of the art is. You
15 know, he said that he thinks that you can now
16 evaluate things like wear and those sorts of things.
17 And I just want to know why people sitting here
18 don't think that you can look at fabric issues with
19 the type of testing that we're talking about and Jim
20 says that's possible, and he says that the test has
21 evolved. And so I want to know what the evolutions
22 are.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 PARTICIPANT (Medtronic): I think one of
2 the things why we cannot evaluate the wear at the
3 present time is because of the extremely high
4 testing speeds. So I think we should separate the
5 metals from the fabric. And in metal fatigue in
6 some development industries, you know, there are
7 testing methodologies according to which you cannot
8 -- you can test metals provided that you know the
9 real service conditions. But however in fabric
10 fatigue due to extremely accelerated speeds, you
11 know, we can get a lot of artifacts because a lot of
12 the nonphysiological graft material movement.

13 My experience is that we can better
14 failure modes by slowing things down. Like fabric,
15 usually one to hertz -- I can replicate the clinical
16 failing ones. But at like 100 or 200 hertz, I don't
17 think you can do that.

18 MS. ABEL: Okay. So it has to do with
19 you're talking about different tests than what Dr.
20 Chwirut uses?

21 Tom?

22 MR. GREENAN: Even the conduit material

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 could have very different results on fabric ware. So
2 they're concerned about testing under these tests
3 something that has the same displacement running at
4 different speeds and running with different conduit
5 material can have significantly different results.
6 And I don't think we know how these may relate to
7 the clinical conditions. Those are some of the
8 concerns.

9 MS. ABEL: Yes. Tom?

10 DR. FOGARTY: Yes. Part of the issue,
11 I'm sure you would test for anything. But in
12 reality, the points of stress and strain between the
13 metal and the fabric continually changes as the
14 aneurysm reconfirms. So I don't know if it's
15 possible to accommodate that in any testing because
16 you don't know in what direction it's going to
17 reconfirm. I mean, you can test all day in many,
18 many cycles, but it really doesn't test what in fact
19 happens in most patients.

20 MS. ABEL: Well, you're getting -- you
21 know, it may not be predictive of exactly how the
22 device is being challenged in a clinical situation.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 But it certainly can give you some information.

2 I mean, we know that there's motion, so
3 that if you test with motion, you can look at
4 whether you've got some weak points. I mean, I
5 would agree it's not -- I mean --

6 MR. GREENAN: Yes, no --

7 MS. ABEL: -- very clearly it's not
8 going to tell you everything you need to know.

9 MR. GREENAN: Yes. I agree with that.
10 But if you interpolate that to you're going to
11 prevent erosion or help prevent erosion, it may or
12 may not.

13 DR. MARIN: Could I just clarify my
14 understanding of what we're trying to measure here.
15 Because let's be clear that bending fatigue is a
16 very different phenomenon from abrasion. And so the
17 test that you would design for measuring bending
18 fatigue, for example, of the components whether it
19 be the fabric or the stent components would be a
20 separate test from looking at the surface abrasion
21 between those two components. And bending fatigue,
22 for example, you would want to look at a rate at

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 which you can accelerate the bending of the
2 components and which would not necessarily be the
3 same rate or the same displacement for the stent as
4 for the fabric. Because their bending module is so
5 very, very different. So that's one issue that is
6 -- and, again, totally different from the question
7 of ware. And surface ware and abrasion has to be
8 taken as a separate issue where there is micromotion
9 between the two surfaces moving against each other.
10 And clearly that is a condition that one needs to
11 try and reproduce, but does require control both of
12 the area of contact and the pressure between the two
13 as well as the speed and frequency of the cycling.

14 So yet there are very specific issues
15 here and one ware fatigue, abrasion test will not do
16 it all.

17 MS. ABEL: Right. Agreed.

18 Okay. Jim, did you want to talk about
19 some of the improvements in the methodology over
20 time?

21 DR.CONTI: Over the past several years
22 two different groups have convened experts from

1 around the world to evaluate possible modifications
2 to the durability testing done on stents and stent
3 grafts. One is an ASTM committee and one is the
4 AIME that reports to -- or represents us as a nation
5 for the ISO.

6 MS. ABEL: It's actually the ISO
7 Committee.

8 DR. CONTI: A lot of work has been done,
9 hundreds and hundreds of hours of just committee
10 alone, to say nothing of what's on the outside. But
11 to summarize where we have come, we're trying to be
12 sensitive to developing a replicate system that
13 exposes the implantables to the kinds of chemistry
14 and mechanical loading that they'll experience in
15 vivo. In addition, we're trying to be sensitive to
16 the fact that we need to generate enough information
17 to predict safety of these products without
18 burdening all of us that do testing with the
19 enormously long lead times in the testing areas.
20 We've been very sensitive to that.

21 How can we do the best testing that we
22 can and do it in as quick a time as possible?

1 Now the things that we've agreed upon.
2 I don't necessarily agree upon all the issues, but
3 I'm pretty happy with in general where we've ended
4 up. And it's really a very logical approach.

5 If you take a material, if you take what
6 we call a mock vessel that has the appropriate
7 biologically relevant properties and you put your
8 product inside of that vessel and expose it to
9 pulsatile loading, and you monitor with whatever
10 techniques you can; high speed photography, an
11 array of sensors, whatever it is that you're an
12 expert at or have available and you determine the
13 kind of motions that that vessel is applying to that
14 product, then you are free to go ahead and test that
15 faster and faster as long as you can verify that all
16 the things that you saw at biologically relevant
17 frequencies are being replicated at the higher
18 frequencies. And that's within the general
19 guidelines I think is a lot of ability for
20 individual scientists and engineers to perhaps be a
21 little bit creative, simply they have to validate
22 what they're doing. But that's really the idea.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We don't just ahead and limit testing
2 frequency just because we are philosophically
3 against it and we don't go wild with the testing
4 frequency. We use some rational measuring points to
5 determine how fast can we test.

6 Now, in general there's a lot of
7 pressure on all of us to try and test a little bit
8 faster. We try to encourage individuals to try and
9 test early and test often. Because if you get into a
10 development project and you need to have something
11 on the table in a year, and it's going to take you
12 14 months to test, well you're just under a terrible
13 pressure situation. Earlier testing will help a lot
14 to pick out things.

15 And in the data that you presented this
16 morning I find one of the most encouraging things
17 I've seen in a long time, and that is that we might
18 be able to do a pretty good job deciding whether or
19 not we go to the next step with the design based
20 upon 36 months worth of testing, maybe 48 months
21 worth of testing. A lot of stuff goes bad in that
22 short period of time.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 MS. ABEL: I'd be careful to make that
2 assumption, because this is obviously just
3 information that was pulled together. We don't know
4 how many patients were out further. It may be that
5 there are more failures going on, but certainly
6 there are observations in the relatively shorter
7 term.

8 DR. CONTI: Yes. And I think if it's
9 true, it's very exciting because now that will give
10 everyone -- everyone wants to do the best job they
11 possibly can. And nobody wants to design, sell a bad
12 medical product, particularly one that's so
13 critically important to life. But if we can, you
14 know, perhaps make certain judgments about things
15 within a couple years worth of testing, then
16 everybody can relax a little bit and maybe slow
17 their testing down some and end up with more
18 informative test results than we've gotten in
19 general in the past.

20 So it's sort of where we've come after
21 all these years.

22 MS. ABEL: Thank you.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 So I think we've heard some of the
2 problems with respect to trying to evaluate the
3 additional failure modes is that people do the
4 faster accelerated testing. And so there are
5 different ways you can possibly address the problem.
6 You can slow it down or you can do additional tests,
7 or you can come with the more complicated model that
8 lets you do everything in one single test. I can't
9 wait to hear about it.

10 MS. WOODS: I have one more question.
11 This is Terry Woods from the FDA.

12 I would just like know how many of you
13 manufacturers think you get useful information of
14 this test you're doing right now, this pulsatile
15 fatigue test or are you just doing it because
16 Dorothy and I saw you have to do it?

17 MS. ABEL: That a good question. Okay.
18 How many think they get useful information? If
19 you'll raise your hand? Get them up there, guys.

20 MR. SMITH: I think there's a lot of
21 benefit in the tests. I mean, even if it's not--
22 for instance just a general fatigue tests we're

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 talking about. Even if it doesn't specifically
2 address, say, fabric wear, it is an observation that
3 can be made. And then if you see an issue in this
4 test, it should lead you to potentially other tests
5 or resolution of the issue through whatever means.

6 Just my colleague from Medtronic talked
7 about 100 to 200 hertz in terms of maybe too fast
8 for fabricware, and that's true. But generally
9 these tests are run, you know, in the 50/40/30 hertz
10 range.

11 PARTICIPANT (Medtronic): I would just
12 like to caution here that I think the results of
13 this test are quite useful, but I do believe that we
14 have to be careful how we interpret data that we got
15 from that test. And that we have to be familiar
16 with, you know, metal fatigue and the fabric fatigue
17 and what that test can give us. And, like, you know
18 if we tested the highly elevated testing conditions,
19 constant temperature test, it may give us a wrong
20 failure mode. So we have to know at what level we
21 have to be and at what testing we have to be to
22 attack certain failure modes.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 MS. ABEL: That's fair.

2 PARTICIPANT: This is a fine test, and I
3 wouldn't want to do away with it. But you have to be
4 aware of the limitations of it and if you want to
5 play around with the edges where these different
6 components are going to fail, you have to stress
7 those components in the worse case cardinal
8 scenario. And unless you're extremely lucky, it is
9 unlikely that this test that will stress the
10 skeletal elements to its maximal conditions is also
11 going to stress some of the components. So you also
12 have to do some individually specific individual
13 test for a particularly designed feature that
14 maximally stresses that particular design feature.

15 MS. ABEL: Thank you.

16 Is there anyone in the room that's
17 willing to talk about having designed a product,
18 found out in the clinic that it broke? Did anyone
19 try to retest it, see if they could duplicate the
20 failures? Anything like that? Okay, we have a
21 couple who will admit it. Will you talk about it?

22 MR. SMITH: I'll talk about. It goes to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 thoracic experience, though, it's not abdominal
2 experience. Still want to talk about it?

3 I think the one thing and a lot of the
4 stuff you've done today shows how early days have
5 underestimated the types of forces. I mean we're
6 talking forces, forces, forces. But, you know,
7 there's motion and then bending independent of how
8 much force it takes to do that. And I think in our
9 thoracic device we had longitudinally oriented wires
10 to provide columnar strength for deployment and
11 short term antimigration, basically. And over time
12 those fractures -- those spines were seen to
13 fracture in a small number of cases.

14 And I think the first step when you have
15 such an incident is to understand why something like
16 that is occurring and then immediately try to
17 reproduce it. And in trying to reproduce it, that's
18 where you end up developing other tests. Okay.

19 And I showed a picture of a bending
20 fatigue test setup that we have. Specifically we've
21 been using it for the thoracic device. And the goal
22 of that test was to: (1) reproduce the fractures

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 that we've seen clinically in the mode that we would
2 see clinically. Not some benchtop mode. But, you
3 know, to replicate what the explants were. And I
4 think then now we have a baseline to improve upon.

5 MS. ABEL: What I was curious about, you
6 know, did you test it in your old testing system?

7 MR. SMITH: Okay. Yes. You know, this
8 is where I think Dr. Conti's been talking about how
9 fatigue testing it has improved. It depends on how
10 far back you go, you know, to say whether it's
11 improved or not. It's obviously improved from the
12 early '90s.

13 In the types of forces that were causing
14 these longitudinal spines to fracture are not
15 generated in the pulsatile fatigue test in and of
16 itself. They can do some longitudinal stresses, you
17 have some pressure but there's not actual bending.
18 You can try to incorporate it in there. You had to
19 develop a completely different test to replicate the
20 failure mode.

21 MS. ABEL: Right. So is there anyone
22 that had failures, went back to do a pulsatile

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 fatigue test with different parameters and then
2 could duplicate it? Testing the old design?

3 MR. DEHDASHTIAN: Well, we did that,
4 Dorothy in the Life Generation I. We had
5 experienced some fracture, went back and revisited
6 the FEA analysis as well as the testing how come it
7 did not predict. And we found out that there is a
8 lot more loads and not uniform loading on the graft
9 that really causes miscalculation from our end.
10 Forced us to rigorously test the device again on the
11 generation II was duplicating the fractures that
12 essentially happened in the clinical environment.

13 MS. ABEL: And what was the change to
14 your test to make it a more rigorous test?

15 MR. DEHDASHTIAN: Just additional
16 forces. The test didn't essentially change per se
17 even though we knew we need to induce not uniform
18 load. We weren't able, that comes back to Terry's
19 question. We were not able to duplicate, even
20 though we know what we want to do, were unable to do
21 it in the best testing. We just increased the
22 loading and the displacement essentially. And --

1 MS. ABEL: Did you change the components
2 in your mock artery or anything?

3 MR. DEHDASHTIAN: Yes.

4 MR. SMITH: I can speak to that a little
5 bit, too. I think that when you do the test, I
6 don't know, seven years ago you have certain
7 assumptions and if you're going to make a new device
8 or modify, you have new assumptions today based on
9 all the information that we get in workshops like
10 this and other meetings. And I think the way that
11 the standard pulsatile fatigue test has evolved, at
12 least in my mind is, you know there are some
13 critical things when you're running that test to
14 ensure you're still getting the diametric
15 deflections that you set the test up for and ensure
16 that your device is following the mock artery if
17 that's what you're going to rely on to create those
18 diametric deflections. And also determining, you
19 know, from either your own clinical data or
20 literature or whatever, what types of compliance you
21 really should be using, what types of pressure pulse
22 you really should be using. The literature says 5

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 to 7 percent on triple A compliance. It's generally
2 accepted in the ISO standard. We've written that
3 down. But there's a whole range of pulse pressure
4 that goes with that.

5 And you know, it's one thing to assume
6 120 over 80, it's another to assume 160 over 80.
7 And then it's another to look at these patients and
8 determine, gees, those with low compliance generally
9 have high delta ps, and those with high compliance
10 in their aorta tend to have low delta p issues. I
11 mean, that's something that we've discovered out of
12 this.

13 So how these tests have evolved in my
14 mind is not necessarily by the equipment or the mock
15 artery, but all the finer details. Okay. How do I
16 make sure my graft is still following that mock
17 artery? How do I make sure that I'm still getting
18 the pressure pulse and the diameter deflections that
19 I thought I had at day one at, you know, at the
20 single aided 100, 200, 300 or 400 million cycle
21 points. So that's how the test has evolved.

22 And in doing so, you know, you put in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the whole device instead of pieces and components
2 and you can see the interaction. I think that's
3 real important.

4 MR. DEHDASHTIAN: I think in
5 continuation of what Lou said, the current test or
6 the tests that we did may not duplicate the
7 clinical. But what helped us, it's kind of evolving.
8 It teaches us, it takes us to the next step -- it's
9 like a stepping stone. Teaches us to go to the next
10 step.

11 MS. ABEL: Medtronic, did you have
12 something you wanted to say?

13 PARTICIPANT (Medtronic): No.

14 MS. ABEL: You had your hands up
15 earlier.

16 MR. VASUTEK: And following on from Dr.
17 Chuter's beautiful picture, it's quite clear that
18 the motions change over the course of the implant as
19 well. I think that's something that should be
20 considered -- might be that you consider oversizing
21 at worse case for the full duration of ten years or
22 whatever. But it may be that the stent itself is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 going to expand both in the mock artery and in the
2 patient. So you may have a different boundary
3 condition from the start of the test until the end
4 of the test.

5 MS. ABEL: All right. Well, I think
6 we're ready to move on to our next table then. So
7 considering I think what I'd like to do here, I
8 think what we talked about as far as pulsatile
9 fatigue testing, the way that it's done by most
10 people right now, we're really looking at primary
11 the stent fracture of the stent fatigue. So if we
12 think just with respect to that, are there other
13 characteristics that should be incorporated within
14 this testing to make sure that we've got as rigorous
15 a test as possible to evaluate that parameter. And
16 you can see that the majority of the people that
17 responded to our survey said that they used straight
18 devices, they didn't use bifurcated devices in their
19 fatigue testing.

20 Is it we're at the point now where all
21 these devices should be tested in a bifurcated
22 system? Does that give you more information?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 PARTICIPANT (Cook, Inc.): It gives you
2 more information. On the mathematical calculations
3 and in the test model we know that the strength
4 comes on the stent immediately above the
5 bifurcation. That's why it fractures there. And it
6 also depends on the angle which strut will break.

7 So, but whether it makes any difference,
8 I'm not sure other than that the broken strut can
9 penetrate the fabric. But it's certainly a stress
10 point. And it's a stress point that you would only
11 see in the bifurcated model.

12 MS. ABEL: Yes.

13 DR. CONTI: I agree completely. I think
14 the shape of the vessel that you put the device into
15 has a huge amount of influence on where loading
16 points actually occur. And if we're aware of
17 certain risky scenarios in the shape of the
18 recipient vessel, then we should try and in
19 corporate that into the more quality of the mock
20 vessel, because it does make a very big difference.

21 MS. ABEL: Thank you.

22 Robert?

1 DR. WHIRLEY: I think I very much agree
2 that the important underlying premise is to identify
3 and test whatever is worse case for a particular
4 device. And in some cases that may well be in a
5 bifurcated configuration. But I wouldn't want to
6 see us move to mandating a bifurcated configuration
7 if that required us to compromise on some other test
8 parameters which could no longer be as challenging
9 as they would be, say, in a bent tube configuration
10 that might have a lot more angulation.

11 So I think the question may be a little
12 more complex than just straight versus bifurcated.

13 MS. ABEL: So for the basic tests that
14 people are doing right now, though, you know if you
15 didn't incorporate any sort of a curvature or
16 whatever, certainly bifurcated makes sense? I think
17 your point's well take, though, if you make other
18 changes.

19 PARTICIPANT (Cook, Inc.): I think
20 they're different. I think that the stress on a
21 curve is on the center of the curve. And the stress
22 on a bifurcated system everything is different. So

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 it depends what graft you're testing. If you're
2 testing a bifurcated graft, then you test that. And
3 if it's a straight curve graft, you test that. Even
4 if it's a tube curved graft.

5 MS. ABEL: Right. I guess what I
6 understood, and maybe I misunderstood you, Robert,
7 but I thought you meant that, you know, it could be
8 that what you really want to look at is the curved
9 vessel and to do that in a bifurcated model may be
10 not possible or it may make your results difficult
11 to interpret, or something.

12 DR. WHIRLEY: That's right. My point
13 was just I would caution us from thinking that
14 because it's a bifurcated model well everything is
15 good and we've automatically incorporate the worst
16 case configuration, whereas a worse case
17 configuration might could be replicated in a much
18 more complex angulated test setup but in a straight
19 tube. And I think that has to be evaluated on a
20 case-by-case basis to come up with what's worse
21 case.

22 MS. ABEL: Mark?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. FILLINGER: I was pretty much going
2 to say the same thing. I think guess the only thing
3 I would sort of add to reinforce that is that when
4 we're sort of setting IFUs for a certain degree of
5 angulation and other sorts of things for the device,
6 that the testing should be sort of directed to that.
7 So you've got a bifurcated graft, you should test
8 the bifurcated graft because if you say the IFU is
9 going to have a certain degree of angulation
10 allowable, the you should try and incorporate in
11 your testing something that will replicate the
12 stresses that would be induced by that degree of
13 angulation.

14 MS. ABEL: Now is that something you can
15 do in this test or is that something for a separate
16 test to work out the angulation? What do people
17 think about that?

18 DR. WHIRLEY: I think you could do it in
19 this test. I think it depends on what you're looking
20 at. If you're looking at the effect of angulation
21 on dilatational fatigue, then you can probably do
22 that in this test. But if you're looking at the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 effect of angulation on longitudinal fatigue or some
2 of the other aspects that we talked about this
3 morning, that may well be better addressed in a
4 different test.

5 MS. ABEL: And so what's most critical
6 to address in this? I mean, you know, would you
7 have value added by incorporating angulated mock
8 artery or a --

9 MR. SMITH: I can say from experience
10 it's somewhat difficult in this standard test setup
11 to put in a certain type of angulation and still try
12 to maintain everything else, especially when you're
13 deciding what type of mock artery to use.

14 So, you know, yes it can be done because
15 that's what engineers say they can do. But actually
16 making that happen is a different thing.

17 In terms of what Dr. Fillinger said,
18 which is pretty good, you can come up with ways to
19 determine what the stresses are due to angulation
20 and try to replicate them in different ways, either
21 by more compression or more deflection or other
22 things to say, okay, if I do see these additional

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 stresses, although it's not the exact type of thing
2 that would happen with the material, it can happen
3 easily with the metallic components and therefore
4 you can say, okay, if I test at a higher stress even
5 mean or alternating stresses, do I see more
6 fractures. And then how does that relate to the
7 stresses caused by angulation.

8 MS. ABEL: So I guess, you know, if we
9 just look at this list here are there any
10 characteristics that should be incorporated in the
11 standard fatigue tests or should we leave that test
12 alone and just note that other tests should be
13 designed to evaluate these failure or these
14 characteristics, or the effects of these
15 characteristics?

16 DR. WHITE: From the clinical data that
17 we've got, two, three, four year stuff there is a
18 pattern in the devices where you can predict where
19 the fractures are. And it distributed, at least as
20 far as I've seen, predictably in each device. So
21 why not model for that device or what we know from
22 that two or three year data and forget all this term

1 stuff with bifurcations and 400 million cycle stuff
2 which hasn't been predictive? I mean, I think you
3 can refine the testing to the clinical scenario that
4 we know exists rather than what's turned out to be a
5 theoretical that didn't predict it.

6 MS. ABEL: Well, I don't know that you
7 can say it didn't predict, because we don't know
8 about the devices that have failed the test and
9 having then gone on to be developed. So you don't
10 know. I mean, it could have shown exactly what it
11 should have shown.

12 DR. WHITE: And I understand that. But
13 it didn't predict in the devices that are clinically
14 used that failure mode.

15 MS. ABEL: But that doesn't mean that
16 you should not do this test?

17 DR. WHITE: Yes, it does to me.

18 MS. ABEL: Because what I'm saying is
19 there may have been ten devices. People were
20 developing devices, they did test and they saw
21 fractures and so then they never went on to develop
22 them or they modified their devices. And then they

1 retested, they didn't have the fractures --

2 DR. WHITE: But there's an assumption in
3 there that that test meant something to start with
4 in this area, which it doesn't. It came from stent
5 technology transferred to this field and it turns
6 out, it doesn't predict anything.

7 MS. ABEL: Actually, it was more graft
8 testing I would say it evolved from.

9 DR. WHITE: Well, okay, that's right.
10 Yes, conventional vascular grafts for stent
11 technologies which in this case in a straight tube
12 400 million cycle ten year thing isn't predicting
13 anything. And it's not predicting --

14 MS. ABEL: How can you say it doesn't
15 predict anything if you've not seen the test results
16 from everyone?

17 DR. WHITE: I think we have seen the
18 test results from everybody.

19 MS. ABEL: Do you think everyone tells
20 you all their testing and all their --

21 DR. WHITE: Well, what they tell you is
22 published in public record. I mean, I don't --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701