



# Rh Complexity

## Serology and DNA genotyping

---

**Connie M. Westhoff, SBB, PhD**  
**Scientific Director, American Red Cross, and**  
**Adjunct, University of Pennsylvania,**  
**Department of Pathology and Laboratory Medicine**



**Molecular Blood Group  
and Platelet Testing Laboratory**

# Rh typing

## U.S.A.- American Red Cross referral experience

- Why serologic D typing is not always straightforward
  - C and e typing also
- Summarize problems
  - How they present
  - Clinical implications
  - Examples
- Perspective – role for genotyping in transfusion medicine
  - Current
  - Future



American  
Red Cross

# Why is D typing sometimes problematic ?

In U.S.- Large number of variables in serologic testing

1. **Multiple Methods** – Hospitals - 91% tube test, 1% Solid phase, 8% Gel  
Some perform AHG test for weak D, others do not  
Donor Centers - Automated analyzers (Olympus PK), tube tests
2. **Different Reagents** – Contain different clones  
**Can react differently with weak or variant D antigens**  
FDA - only reactivity with DIV, DVa, and DVI need be specified  
Results in D typing discrepancies
3. **Variability in expression of RhD protein** ( ~120 different genes = variations)  
All due to changes at the DNA level from “conventional” sequence

Weak D	-	53	different	mutations
Partial D	~	45	“	“
D <sub>el</sub>	~	8	“	“
“others”	~	18		



American  
Red Cross

Molecular Blood Group  
and Platelet Testing Laboratory

# Variation in expression of RhD



**D Positive** - Majority are “conventional”


---

1. Weak D - previously “Du”
  - incidence- 0.2-1% - wide population differences
  - requires indirect antiglobulin test to detect (depends on reagent)
2. D<sub>el</sub> - very weak expression of D
  - “el” because adsorb and elute anti-D
3. Partial D - previously called “mosaic”
  - “missing” part of RhD
  - type as D-positive
  - make anti-D
4. D epitopes - expressed on Rhce proteins
  - cause D typing discrepancies
    - D<sup>HAR</sup>, Crawford, ceRT, ceSL

# Weak D

(Wagner et al. Blood 93:385, 1999)

- single gene mutations – intracellular or cytoplasmic
- many different weak D (Type 1 thru 53 as circles)
- effect quantity of protein, but not D epitopes  
usually do not make anti-D



**Reactivity of weak D is variable with different monoclonal antibodies and with different techniques**

many 3+ , but some very weak +/- or missed in IAT  
Our experience. - weak D type 2, but by Gel testing stronger



**D<sub>el</sub>**


– type as D negative, (including IAT)

- adsorb and elute anti-D

8 different mutations

1/3 of Asians who type D negative; are also C+

In Caucasians - are C+ or E+



Recently “in the news” - have stimulated anti-D in recipients


- most would agree should be D positive as donors

- C or E serologic screening would be eliminated from donor pool

# Partial D

- type as D positive, make anti-D
- don't detect until make antibody
- **mutations located on the extracellular surface**  
altered D epitopes -

some are due to single mutations



DMH - L54P

DVII - L110P

DFW - H166P

DHR - R229K

DHO - K235T

#weak D type 15 - G282D

D<sup>HMI</sup> - T283I


DIM - C285Y

DNU - G353R

D<sup>II</sup> - A354R

\*DNB – G355S- more common in Europe

# Most Partial D's - *RHD* replaced with *RHCE*



## *RHD* - Gene Conversion

D IIIa	
D IIIb	
D IIIc	
D IVa	
D IVb	
D IVbIII	
D IVbIV	
D V (6)	
D VI	
D VI	
D VI	
DFR	
DFR	
DBT	
DBT	


## Normal *RHCE*

DAK	
G-	
Go <sup>a</sup>	
BARC	
BARC	
FPTT	
FPTT	
Rh32	
Rh32	

Alter D epitopes and create new antigens



# D epitopes expressed on Rhce



**D-specific amino acid (s) in the Rhce protein**  
strong reactivity with some monoclonal anti-D

In the US these are a major cause of D typing discrepancies that are referred for *RHD* gene investigation



# FDA licensed reagents in use in U.S.


**4 reagents for Tube testing, 1 for gel**

REAGENT	IgM monoclonal	IgG
Gammaclone	<b>GAMA401</b>	F8D8 monoclonal
Immucor Series 4	<b>MS201</b>	MS26 monoclonal
Immucor Series 5	<b>Th28</b>	MS26 monoclonal
Ortho BioClone	<b>MAD2</b>	Human-polyclonal
Ortho Gel (ID-MTS)	<b>MS201</b>	

**Only two contain same IgM clone**

**Clones can differ in reactivity with variant D antigens**

# Difference in reactivity of D<sup>Har</sup> and Crawford+ RBCs with FDA licensed reagents



REAGENT	IgM monoclonal	IgG	D <sup>Har</sup>	Crawford
Gammaclone	GAMA401	F8D8 monoclonal	POS	POS
Immucor Series 4	MS201	MS26 monoclonal	POS	NEG
Immucor Series 5	Th28	MS26 monoclonal	POS	NEG
Ortho BioClone	MAD2	Human-polyclonal	NEG	NEG
Ortho Gel (ID-MTS)	MS201		POS	NEG

D<sup>Har</sup> and Crawford+ RBCs are non-reactive with human source polyclonal anti-D

Has contributed to perception in U.S. that D typing discrepancies are greater with monoclonal reagents



# Case

---

81 year old African-American woman  
2003- typed as D positive (3+ IS) - received 3 units

6/22/2006- admission for anemia, GI bleed  
D positive, antibody screen negative  
received 3 units D positive blood

7/4/2006 - strongly +DAT, low hgb, elevated bilirubin

**Delayed transfusion reaction**  
Eluate and serum: Anti-D



American  
Red Cross

Molecular Blood Group  
and Platelet Testing Laboratory



# Case

<u>Anti-D</u>	IS	RT	AHG
Gamma ( <i>DMB40-2</i> )	3+	4+	NT
Gamma ( <i>DMB36-1</i> )	3+	3+ <sup>s</sup>	NT
Gamma ( <i>D139-2</i> )	3+	3+	0
Immucor Series 4	0	0	0
Immucor Series 5	0	0	0
Ortho Bioclone	0	0	0

Gamma IgM clone – strongly D positive

Ortho and Immucor clones - D negative, weak D neg



American  
Red Cross

Molecular Blood Group  
and Platelet Testing Laboratory

# D typing “recommendations”

As patients/recipients – should be D negative

		D <sup>HAR</sup>	Crawford	
Gammaclone	GAMA401	POS	POS	→ Donor reagent
Immucor Series	MS201	POS	NEG	
Immucor Series	TH28	POS (varies)	NEG	
Ortho BioClone	MAD2	NEG	NEG	→ Patient reagent
Ortho Gel Card	MS201	POS	NEG	

How often are these encountered ?

Estimated that Crawford- 1:900 AA in Southern US

D<sup>HAR</sup> ? German background - higher in Midwest

Our experience – Crawford represent large number samples referred for D typing discrepancies



American  
Red Cross

Molecular Blood Group  
and Platelet Testing Laboratory

# Summary of *RHD* referrals

2006 (Jan - Sept ) 9 months/ 210 samples

<i>RHD</i>	num
<b><u>D discrepancies</u></b>	
Crawford	15
weak D type 2	7
type 18	3
type 1	2
type 15	1
D <sub>el</sub> (M295I)	1
Partial DAU-2	<u>1</u>
	<b>30</b>
<b><u>D+ with anti-D</u></b>	
DAR	3
DIIIa	2
Crawford	2
DIVa type 1	3
DIVb	1
DVI	1
weak D type 2	<u>2</u>
	<b>14</b>
<b><i>RHD</i> zygosity</b>	16



# Serologic D typing For Donor Testing

**Goal:** to label all donor RBCs with D antigen as D positive

## Problem:

Weak D - some are missed - even with IAT testing

D<sub>el</sub> - all are typed as D negative

**These can stimulate anti-D in D negative patients**

**Important:** all D<sub>el</sub> (to date) and majority of weak D are inherited with C+ or E+, so can be removed from the D negative donor pool by serologic typing for C and E.

**Questions not yet answered:**

- Will we accept no anti-D ?
- In specific group - girls and women of child bearing age
- Accept anti-K (HDFN) and anti-c in U.S., so inconsistent ?





# Serologic D typing

## For patient and OB testing

---

**Goal:** to detect those at risk for anti-D

**“Most” weak D-** are not at risk for anti-D (there are exceptions)

**Partial D** - at risk for anti-D, but type as D positive so are not detected

- female children & women of child-bearing age better served treated as D negative for transfusion and RhIG candidates

### Problems:

- Serologic tests cannot distinguish weak D from partial D
- Weak D mutations may also alter D epitopes



# DNA genotyping for *RHD*

---

Would this resolve the “D problem” ?

- **More complex than change in methodology**
- **How to act on the results ?**

**Variant D (weak and partial)- if treat all as D negative- would be significant burden to D negative donor pool**

- **Additional data needed -which *RHD* change/generate new epitopes ?**

What is the Goal ? To have no anti-D produce in any patient ?

**High throughput platforms needed:**

- Many regions of gene must be sampled
- Complex algorithm for interpretation

# Summary of case referrals

2006 (Jan - Sept ) 9 months/ 210 samples

<b><i>RHD</i></b>	num
<b><u>D discrepancies</u></b>	
Crawford	15
weak D type 2	7
type 18	3
type 1	2
type 15	1
D <sub>el</sub> (M295I)	1
Partial DAU-2	1
	<b>30</b>
<b><u>D+ with anti-D</u></b>	
DAR	3
DIIIa	1
DIVa type 1	4
DIVb	1
DVI	1
weak D type 2	2
	<b>14</b>
<b><i>RHD</i> zygosity</b>	16

<b><i>RHCE</i></b>	num
<b><u>e+ variants</u></b> (most with anti-e)	<b>43</b>
<b><u>C+ variants</u></b> (with anti-C or -Ce)	<b>8</b>
<u>e discrepancies</u>	6
<u>E discrepancies</u>	2




**Many RHCE problems  
also**

Others	Of 210 samples
ABO discrepancies	14
Dombrock screening	30
Typing multiply transfused pts.	14
Duffy typing discrepancies	3
McLeod	4
MNSs, U-	13
Misc. Nulls/new polymorphisms	6

# When is serologic C typing not straightforward ?

In African Americans, Hispanic, and mixed ethnic groups with a specific *RHD-CE (3-8)-D* hybrid gene



Type as C positive; make anti-C

Hybrid *D-CE-D* gene is linked to variant e; make anti-e

This RH haplotype is prevalent in sickle cell patients  
estimated 22% African Americans



American  
Red Cross


Molecular Blood Group  
and Platelet Testing Laboratory


# When is serologic e typing not straightforward ?


In African Americans, Hispanic, and mixed ethnic groups


 “conventional” ce


Many different genes, all encode altered expression of e antigen

 **ce<sup>s</sup>**

 **ceMO**

 **ceEK**

 **ceAR**

 **ceBI**

Type as e positive

Make anti-e and/or anti-ce

Prevalent in sickle cell patients



# When is serologic e typing not straightforward ?

In African Americans, Hispanic, and mixed ethnic groups

associated with altered C





D/Ce/D hybrid

Type as C positive

Make anti-C






associated with altered D



Type as D positive

Make anti-D

altered expression of e



Type as e positive

Make anti-e and/or anti-ce



American  
Red Cross

Molecular Blood Group  
and Platelet Testing Laboratory

Prevalent in sickle cell patients

# Role for genotyping in transfusion medicine

## Current:

- Type multiply transfused patients
- Screen donor units when reagents are not available
- Resolve reagent typing discrepancies
  - Reveal shortcomings of current reagents- important to design of future reagents
- Predict fetal risk for hemolytic disease
  - Paternal RHD zygosity
  - Genotype fetus- amniotic fluid or maternal plasma
- Resolve antibody identification (is it allo or auto)
- Provide compatible donor units for sensitization sickle cell patients (variant D, C, and e antigens)



# Role for genotyping in transfusion medicine

## Going Forward:

Donors: Screen for antigen negative units with high throughput

- 2 types
  - One serologic
  - One DNA based

**Will provide the opportunity to validate genotyping and to detect any exceptions in different populations**

Patients: Sickle cell “antigen matching” programs

- Transfusion predicted to increase with STOP trial outcome
  - Variant C, e, and D makes “genetic match” for Rh superior
  - Detect those at risk for production of antibodies to high-incidence Rh antigens.





# Role for genotyping in transfusion medicine

## Future:

Integration of DNA-based assays into blood bank

- as genomics is applied to diagnosis and treatment

- Potential for “antigen” matching
  - How many antigens (15 or 50 ?)
  - Which patient population
    - Those needing long term transfusion support
    - For all
- At what cost? Or.....at what savings?
  - Impact workload
    - Walk-away systems
    - Decrease sensitization=decreased ABID

**Require major changes in management of donor inventory**



# Blood group “typing”


power in combination of both technologies

Phenotype



agglutination

Genotype



DNA

Focus should not be

to prove or decide if one method is superior OR  
to replace all serologic typing with DNA-based testing

Focus

use each to strengthen the other

Shortcoming of D typing reagent revealed by DNA testing

– use knowledge to improve serologic reagents

Shortcomings of SNP testing revealed by serology

- use knowledge to design additional SNP



# Acknowledgements

---

- **Molecular Blood Group and Platelet Testing Laboratory**
  - Sunitha Vege
  - Wynn Meyer
  - Michelle Henry
  - Trina Whorley
  - Cindy Marks

## **Support - National Reference Laboratory Blood Group Serology**

- Sandra Nance
- Staff

