

**FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH**

INSTRUCTIONS FOR COMPLETING FORM 3500A
(Specific to Medical Device Reporting)

With

CODING MANUAL FOR FORM 3500A

December 14, 1995

**FOR FURTHER INFORMATION OR QUESTIONS CONTACT
THE REPORTING SYSTEMS MONITORING BRANCH
DIVISION OF SURVEILLANCE SYSTEMS
OFFICE OF SURVEILLANCE AND BIOMETRICS
1350 PICCARD DRIVE, HFZ-533
ROCKVILLE, MD 20850
(301) 594-2735**

FOREWORD

The Food and Drug Administration (FDA), has established the final reporting requirements for Device User Facilities and Manufacturers under the revised Medical Device Reporting Regulation, 21 CFR, Part 803. The final reporting requirements for Distributors is contained in 21 CFR Part 804. These regulations enact a uniform system for the submission of reports of deaths, serious injuries and certain malfunctions associated with the use of medical devices. These reports must be submitted to the FDA and/or the manufacturer on FDA Form 3500A (MedWatch form).

The FDA has developed this guidance document to assist User Facilities, Distributors and Manufacturers in complying with the requirements established in Parts 803 and 804 when completing applicable sections of Form 3500A. The Coding Manual for Form 3500A is also included, and it contains instructions for proper identification and utilization of the adverse event and manufacturer evaluation codes which must be provided when filling out certain sections of Form 3500A. Various references, including telephone and facsimile numbers, addresses for obtaining copies of this document and Form 3500A and other sources of information are also included in this document.

These final instructions for completing FDA Form 3500A with the coding manual replaces all versions previously in

TABLE OF CONTENTS

Instructions for Completing FDA Form 3500A	1-11
Instructions for Using the Device Coding Manual for FDA Form 3500A (Parts I and II).....	12-13
Important Reference Information and Mailing Addresses	14
Sample of FDA Form 3500A	15a-15b
 PART I - EVENT PROBLEM CODES AND TERMS (For Section F, Block F10 of Form 3500A)	
Part I, Subpart A - Alphabetic Listing of Event Terms with Their Corresponding Numeric Event Code	
<i>Patient-Related Event Terms</i>	16-26
<i>Device-Related Event Terms</i>	27-40
Part I, Subpart B - Numeric Listing of Event Codes with Their Corresponding Event Terms.....	41-54
Part I, Subpart C - Definitions for Event Terms Listed Numerically by Event Code	55-68
 PART II - DEVICE MANUFACTURER CODES (For Section H, Blocks H3 and H6 of Form 3500A)	
Part II, Subpart A - Device Evaluated by Manufacturer Codes (Block H3).....	69
Part II, Subpart B - Evaluation Codes (Block H6)	69-77
 Method of Evaluation Codes:	
<i>Source of Device Evaluated</i>	69
<i>Type of Evaluation Performed</i>	69
 Results of Evaluation Codes:	
<i>Category A - Device</i>	69-70
<i>Category B - Use of Device</i>	70-71
<i>Category C - Physiological/Procedural Factors</i>	71
<i>Category D - Device Component/Subassembly Failures</i>	71-76
<i>Category E - Computer-, Imaging System-, Microprocessor-Controlled Device Problems</i> ..	76
Conclusions Codes	76-77

INSTRUCTIONS FOR COMPLETING FDA FORM 3500A

These instructions and FDA Form 3500A should be used by user facilities, distributors, and manufacturers for **MANDATORY** reporting of adverse events and product problems associated with the use of medical devices, as designated in the applicable statutes and FDA regulations. Form 3500A should also be used by User Facilities to report, voluntarily, malfunctions and other product problems which did not cause or contribute to a death or serious injury.

- All entries should be typed.
- Complete all sections that apply. When information is not readily available, select one of the following:
 - **NA** for not applicable
 - **NI** for no information at this time (but may be available at a later date)
 - **UNK** for unknown
- Dates should be entered as month/day/4-digit calendar year (e.g. June 3, 1995 = 06/03/1995). If exact dates are unknown, or only a portion of the date is known (e.g. month/year), provide the best estimate.
- A computer-generated facsimile of the form may be submitted in lieu of the preprinted form if the submitter has received written preapproval from the MEDWATCH office. (see address on page 14). It is acceptable to reproduce each side of the pre-printed form on two separate, front-only pages.
- At the top-right corner, enter the manufacturer report number, user facility number or distributor number in the correspondingly labeled box. Complete both report numbers, if applicable, in order to cross-reference this report with a report from another source on the same event. For follow-up (i.e. supplemental) reports, the number should be identical to the number assigned to the original report.
- The **distributor report number** consists of three components: The FDA-assigned registration or identification number for the distributor of the device, the four-digit calendar year in which the report is submitted, and a consecutive five-digit sequence number for each report filed during the year by the distributor (e.g. 1234567-1995-00001, 1234567-1995-00002). If a distributor does not have an assigned identification number, it should use all zeros in the appropriate space on the initial report, and continue to use zeros on subsequent reports until the FDA-assigned number is received. The distributor report number is also entered in block F2.

The **manufacturer report number** consists of three components: The manufacturer's FDA registration number for the manufacturing site of the reported device, the four-digit calendar year in which the report is submitted, and a consecutive five-digit sequence number for each report filed during the year by the manufacturer (e.g. 1234567-1995-00001, 1234567-1995-00002). If the manufacturing site does not have a registration number, then FDA will assign a temporary one, to be used until the site is officially registered. The manufacturer report number is also entered in block G9.

MANUFACTURER AND DISTRIBUTOR REPORT NUMBERS (continued)

In cases where a reporting site is registered as both a manufacturer and a distributor, and the registration and/or FDA-assigned identification numbers are identical for both, then the five-digit sequence number for reports submitted during the year by either one, may not be duplicated. For example, in the role of manufacturer, the report number would consist of the registration number, a 4-digit calendar year, and a five-digit sequence number, (e.g. 1234567-1995-00001, 1234567-1995-00002, etc.). In the role of distributor, the registration number and year would remain the same, but the five-digit sequence number must be different, (e.g. 1234567-1995-00003, 1234567-1995-00004, etc.).

- The **user facility report number** consists of three components: The facility's 10-digit Health Care Financing Administration (HCFA) number, the 4-digit calendar year in which the report is submitted, and a consecutive 4-digit sequence number for each report filed during the year by the facility (e.g. 1234567890-1995-0001, 1234567890-1995-0002). If the HCFA number has fewer than 10 digits, fill in any remaining digits with zeros. If a facility does not have a HCFA number, the first report and any subsequent reports, should be submitted with all zeros in the HCFA space, and FDA will assign a number to be used in future reports (e.g. 0000000000-1995-0001). If a facility has more than one HCFA number, the facility must select one of those numbers, and use the same number for subsequent submissions. The user facility report number is also entered in block F2.

If a user facility has multiple sites, the primary site can report centrally and use one reporting number for all sites if the primary site provides the name, address and HCFA number for each respective site.

- If the fields do not provide adequate space for narrative entries, attach an additional page(s) and indicate the appropriate section and block number next to the narrative continuation.
- At the top of both the front and back page, enter the page number and the total number of pages submitted (include attachments in the total). All attached pages should be identified as page _of_ and should display the user facility, distributor, or manufacturer report number in the upper right corner as applicable.

SECTION A: PATIENT INFORMATION

If more than one patient was involved in the same event, complete all applicable sections of the form for the first patient and all of Section A and blocks B2, B5, B6, B7, D10, and F10 for each additional patient. For each additional block F10 enter the corresponding patient identifier in block F3.

- A1: **Patient Identifier** - Provide the patient's initials or some other type of identifier that will allow both the submitter and the initial reporter (if different) to identify the report if contacted for follow-up. **DO NOT** use the patient's name or social security number. (The patient's identity is held in strict confidence by FDA and protected to the fullest extent of the law).

SECTION A: PATIENT INFORMATION (continued)

A2: **Age** - Enter the patient's birthdate, if known, or the patient's age at the time of the event onset. For age indicate time units used, e.g. years, months, days. Provide the best estimate if exact age is unknown.

- if the patient is **3 years or older**, use **years** (e.g., 3 years).
- if the patient is **less than 3 years old**, use **months** (e.g. 24 months).
- if the patient is **less than 1 month old**, use **days** (e.g., 5 days).

If the adverse event is a congenital anomaly, use the age or birthdate of the child or the date pregnancy is terminated. If information is available as to the time during pregnancy when exposure occurred, provide that information in narrative block B5.

A3: **Sex** - Enter the patient's gender. If the adverse event is a congenital anomaly, report the sex of the child.

A4: **Weight** - Indicate whether the weight is in pounds (lbs.) or kilograms (kgs). Make a best estimate if exact weight is unknown. If the adverse event is a congenital anomaly, use the weight of the child.

SECTION B - ADVERSE EVENT OR PRODUCT PROBLEM

B1: **Adverse event and/or product problem** - Choose box 1 to report an adverse event and/or box 2 to report a product problem. An adverse event is selected when reporting a death or serious injury. Product problem is selected when reporting malfunctions that could lead to a death or serious injury if it were to recur. Both boxes should be checked if a malfunction or product problem caused a death, or serious injury.

B2: **Outcomes attributed to adverse event** - Check all that apply to the reported event.

Death - Check if a device has or may have caused or contributed to the death of a patient (this includes employees).

Life-threatening - Check if the patient was at risk of dying at the time of the adverse event or if the continued use of the device might have resulted in the death of the patient.

Disability - Check if a device has or may have caused or contributed to a permanent injury or impairment to the patient.

Congenital anomaly - Check if exposure to a medical device prior to conception or during pregnancy has or may have resulted in an adverse outcome in the child.

Required intervention to prevent permanent impairment or damage - Check if medical or surgical intervention was necessary to preclude permanent impairment of a body function or to prevent permanent damage to a body structure that was or may have been due to the use of a medical device.

SECTION B - ADVERSE EVENT OR PRODUCT PROBLEM - BLOCK B2 (continued)

Other - Check only if the other categories are not applicable to the event. Briefly describe the patient outcome in the space provided. The actual narrative of the event will be entered in block B5.

- B3: **Date of event** - Provide the actual or best estimate of the date of first onset of the adverse event. For congenital anomalies, the date of birth or the date pregnancy is terminated should be used. If day is unknown, month and year are acceptable. If day and month are unknown, year is acceptable.
- B4: **Date of report** - The date the initial reporter provided the information about the event (i.e. the first person or entity who initially provided the information to the user facility, manufacturer or distributor).
- B5: **Describe event or problem:**

Adverse event - Describe the event in detail. Include a description of the nature of the event, how the device was involved, any environmental conditions that may have influenced the event, patient follow-up or required treatment (e.g. relevant clinical information, medical status prior to the event, signs, symptoms, diagnoses, clinical course, treatment, outcome, etc.). If available, and if relevant, include synopses of any office visit notes or the hospital discharge summary. If permitted by the institution, copies of these records may be attached with any confidential information deleted. **DO NOT identify any patient, physician or institution by name. The initial reporter's identity should be provided in full in section E.** Enter results of relevant tests and laboratory data in block B6 and preexisting medical conditions and other relevant history in block B7.

Product Problem - Describe the problem and the circumstances surrounding the defect or malfunction of the medical device. The results of any evaluation of the malfunctioning device and, if known, any relevant maintenance/service information should be included in this section.

- B6: **Relevant tests/laboratory data, including dates** - Include any relevant baseline laboratory data prior to the administration or use of the medical device, all laboratory data used in diagnosing the event and any available laboratory data/engineering analyses that provide further information on the course of the event. Include any available pre and post-event medication levels and dates if applicable. Include synopses of any relevant autopsy, pathology, engineering or lab reports, if available. If preferred, copies of any reports may be submitted as attachments with all confidential information deleted. **DO NOT identify any patient, physician or institution by name. The initial reporter's identity should be provided in full in section E.**
- B7: **Other relevant history, including preexisting medical conditions** - If relevant, provide information on other known conditions in the patient (e.g. hypertension, diabetes, renal/hepatic dysfunction, etc.) and significant history (e.g. allergies, race or ethnic origin, pregnancy, smoking and alcohol use, drug abuse, etc.).

SECTION C: SUSPECT MEDICATION(s)

This section is intended to be used when reporting adverse events associated with the use of suspect medications.

C1 - C10: **Skip - Not intended for medical device use.**

SECTION D: SUSPECT MEDICAL DEVICE

The suspect medical device is the device that may have caused or contributed to the adverse event or the device that malfunctioned. If more than one suspect medical device was involved in the event, complete all of sections D and F for the first device and a separate section D and blocks F9, F10, F13, and F14 for each additional device. Pair each section D with it's corresponding section F by marking each as follows: "Device 1", "Device 2", etc.

D1: **Brand name** - The trade or proprietary name of the suspect medical device as used in product labeling or in the catalog. (e.g. Flo-Easy Catheter, Reliable Heart Pacemaker, etc.). This information may be on a label attached to a durable device, may be on a package of a disposable device, or may appear in labeling materials of an implantable device.

D2: **Type of device** - The generic or common name of the suspect medical device or a generally descriptive name (e.g. urological catheter, heart pacemaker, patient restraint, etc.).

D3: **Manufacturer name & address** - If available, list the full name and mailing address of the manufacturer of the suspect medical device.

D4: **Operator of device** - Indicate the type (**not the name**) of person operating or using the suspect medical device on the patient at the time of the event.

Health professional = physician, nurse, respiratory therapist, etc.

Lay user/patient = person being treated, parent/spouse/friend of patient.

Other = Nurse's aide, orderly, etc.

D5: **Expiration date** - If applicable. This date can often be found on the device itself or printed on the accompanying package.

D6: **Product identification numbers** - If available. Provide any or all identification numbers associated with the suspect medical device exactly as they appear on the device or device labeling.

Model # - The exact model number found on device label or accompanying packaging.

Catalog # - The exact catalog number as it appears in the manufacturer's catalog, device labeling or accompanying packaging.

Serial # - Can be found on device label or accompanying packaging. This number, assigned by the manufacturer should be specific to each device.

SECTION D: SUSPECT MEDICAL DEVICE - BLOCK D6 (continued)

Lot # - Can be found on the label or packaging material.

Other # - Any other applicable identification number (e.g. component number, product number, part number, barcoded product ID, etc.).

- D7: **If implanted, give date** - For medical devices that are implanted in the patient, provide the implant date or best estimate. If day is unknown, month and year are acceptable. If month and day are unknown, year is acceptable.
- D8: **If explanted, give date** - If an implanted device was removed from the patient, provide the explant date or best estimate. If day is unknown, month and year are acceptable. If month and day are unknown, year is acceptable.
- D9: **Device available for evaluation?** - Indicate whether the device is available for evaluation. Indicate if the device was returned to the manufacturer and if so, the date of the return. (**DO NOT send the device to FDA**).
- D10: **Concomitant medical products and therapy dates** - List and provide device product names and therapy dates for any other medical devices that were being used on a patient at the time of the event. (**DO NOT include products used to treat the event**).

SECTION E: INITIAL REPORTER

Indicate the person who initially provided information about the event to the user facility, manufacturer or distributor. For user facilities, this person may or may not be the designated medical device reporting (MDR) contact.

- E1: **Name, address & phone #** - Please provide the name, mailing address and phone number of the first person or entity that provided the information to the user facility, manufacturer or distributor.
- E2: **Health professional?** - Indicate whether the initial reporter is a health professional (e.g. physician, pharmacist, nurse, etc.) or not.
- E3: **Occupation** - Indicate the initial reporter's occupation, and include specialty if appropriate.
- E4: **Initial reporter also sent report to FDA** - Indicate whether the initial reporter also notified or submitted a copy of this report to FDA.

SECTION F: FOR USE BY USER FACILITY/DISTRIBUTOR - DEVICES ONLY

- F1: **Check one** - Indicate whether the report is from a user facility or distributor.
- F2: **UF/Dist report number** - Enter the complete number of the report exactly as entered in the upper right corner of the front page.

SECTION F: FOR USE BY USER FACILITY/DISTRIBUTOR - DEVICES ONLY

- F3: **User facility or distributor name/address** - Enter the full name and address of the user facility or distributor reporting site.
- F4: **Contact person** - Enter the full name of the medical device reporting (MDR) contact person. This is the person designated by the facility's most responsible person as the device user facility/distributor contact for this requirement. FDA will conduct its MDR correspondence with this individual. The contact person may or may not be an employee of the facility. However, the facility and its responsible officials will remain the parties ultimately responsible for compliance with MDR reporting requirements.
- F5: **Phone number** - Enter the phone number of the medical device reporting (MDR) contact person.
- F6: **Date user facility or distributor became aware of event** - Enter the date that the user facility's medical personnel or the distributor became aware that the device may have caused or contributed to the reported event.
- F7: **Type of report** - Check the appropriate box to identify the type of report being filed.
- If a follow-up report, make sure that the UF/Dist report number for the previously submitted initial report is recorded in block F2. In the blank provided in block F7, record the appropriate sequence of follow-up to that particular initial report (e.g. first follow-up report = follow-up #1, second follow-up report = follow-up # 2, etc.). Follow-up reports should not repeat material that was submitted in the initial report but should only provide additional or corrected information.
- F8: **Date of this report** - Enter the date that the report was forwarded to the manufacturer and/or the FDA.
- F9: **Approximate age of device** - Enter the age of the device or a best estimate. (Include unit of time used, e.g. month, year).
- F10: **Event problem codes** - Enter up to 3 "patient" and 3 "device" codes from the coding manual (See Part I) that most accurately describe the event. Patient codes describe what happened to the patient as a result of the event and device codes describe device failures or problems during the event. If more than 3 "patient" codes or more than 3 "device" codes are needed, record them on a separate sheet, mark it "F10" and provide the report number and page number.
- F11: **Report sent to FDA?** - Check "yes" or "no" and indicate the date sent, if applicable.
- F12: **Location where event occurred** - Check the location of the actual occurrence of the event. If none of the designated location options apply, check the "other" box and provide the location.
- F13: **Report sent to manufacturer?** - Check "yes" or "no" and indicate the date sent, if applicable. Reports must be sent to the manufacturer if the manufacturer is known.
- F14: **Manufacturer name/address** - Enter full name and address of the device manufacturer, if available.

SECTION G: ALL MANUFACTURERS

- G1: **Contact office - name/address(& mfring site for devices)** - Enter the full name and address of the manufacturer reporting site [contact office] including contact name. If the manufacturing site is different than the contact office, enter "Site:" and the name and address of that site.
- G2: **Phone number** - Enter the telephone number of the contact office.
- G3: **Report source** - Check the box(s) that most accurately describes how the manufacturer [contact office] became aware of the reported adverse event or where the information about the adverse event came from.

Foreign - Foreign sources include foreign governments, foreign medical facilities, etc.

Study - Postmarketing, clinical trial, surveillance, or other study which involves a systematic collection of adverse events from a protocol designed specifically to investigate product safety. Check this box when reporting on required postmarket and discretionary postmarket studies under Section 522 of the FD&C Act.

Literature - If the report source is the scientific literature or an unpublished manuscript, then attach a copy of the article or manuscript, in English. Record the date of the article as the date of the event (block B3), and provide a full literature citation in block H10.

Consumer - (including attorneys) - Generally, additional information should be sought from the treating health care provider. A determined effort should be made to obtain additional detailed information from health professionals for all events initially reported by consumers. When this additional information is obtained, the follow-up report should check health professional rather than consumer.

Health professional - self explanatory.

User facility - self explanatory.

Company representative - Check if the company representative reported the event to the contact office based on information from a health professional. The health professional should be listed as the initial reporter in section E.

Distributor - Distributor should be checked if the manufacturer received the report from a distributor of the suspect medical device.

Other - Check if none of the available options is appropriate and specify source.

- G4: **Date received by manufacturer** - This is the date when the manufacturer became aware of information that an adverse event or medical device malfunction occurred. This would apply anywhere in the world. For follow-up (i.e. supplemental) reports, use the date that the follow-up (supplemental) information was received.
- G5: **Skip - Not intended for medical device reports.**

SECTION G: ALL MANUFACTURERS (continued)

G6: **Skip - Not intended for medical device reports.**

G7: **Type of report:**

"5-day" - Check if submitting a 5-day report.

"Initial" - Check if first or initial submission of a 5-day or 30-day report.

"Periodic" - Not intended for medical device reports.

"Follow-up (i.e. supplemental) #" - Check if follow-up (supplemental) submission. The manufacturer report number for the previously submitted initial report is recorded in block G9. In the blank provided in block G7 after "follow-up", record the appropriate sequence of follow-up (supplement) to that particular initial report (e.g. first follow-up (supplemental) report = follow-up #1, second follow-up (supplemental) report = follow-up # 2, etc). Follow-up (supplemental) reports should not repeat any information that was submitted in the initial report, but should only provide additional or corrected information.

G8: **Skip - Not intended for medical device reports.**

G9: **Mfr. report number** - Enter the manufacturer report number exactly as it appears in the upper right corner of the front page. For a follow-up (supplemental) report the manufacturer report number is to be identical to the number assigned to the initial report.

SECTION H: DEVICE MANUFACTURERS ONLY

H1: **Type of reportable event** - Check the appropriate box. These choices represent the categories of events that device manufacturers are required to report. ("See guidelines" refers to the applicable sections in 21 CFR Part 803 reporting guidelines associated with device malfunctions).

Other - Specify the type of report in the space provided. This option is intended to capture reports that a manufacturer believes the agency should be aware of that are not covered by death, serious injury, or malfunction as these terms are defined by the statute, regulation or guidelines.

H2: **If follow-up (supplemental) report, what type?** - Check the box(s) that most accurately describes the nature of the follow-up (supplemental) report.

Correction - Changes to previously submitted information.

Additional information - Information concerning the event that was not provided in the initial report because it was not known/available when the report was originally submitted.

Response to FDA request - Additional information requested by FDA concerning the device/event.

Device evaluation - Evaluation/analysis of device.

SECTION H: DEVICE MANUFACTURERS ONLY (continued)

- H3: **Device evaluated by mfr?** - Check the box marked "not returned to mfr." if an evaluation could not be made because the device was not returned to or made available to the manufacturer. Check the box marked "yes" if an evaluation was made of the suspect or related medical device. If an evaluation was conducted attach a summary of the evaluation and check the box marked "evaluation summary attached". If an evaluation of a returned suspect or related medical device was not conducted, check the box marked "no" and attach a page to explain why not or provide the appropriate code from the coding manual (Part II, Subpart A) in the space provided.
- H4: **Device manufacture date** - Enter the month and year of manufacture of the suspect medical device using a MM/YYYY date format.
- H5: **Labeled for single use?** - Indicate whether the device was labeled for single use or not. If the question is not relevant to the device being reported, (e.g. an x-ray machine), check "no".
- H6: **Evaluation codes** - Enter the applicable codes from the coding manual (See Part II, Subpart B) for one or more of the categories listed. Conclusion codes must be entered even if the device was not evaluated.
- H7: **If remedial action initiated, check type** - Indicate the applicable actions(s). If other, specify the type of action in the space provided. Most of these terms are defined or further explained in the Act or in the FDA regulations concerning remedial action (see 21 U.S.C. 360h and 21 CFR part 803).
- H8: **Usage of device** - Indicate whether the use of the suspect medical device was the initial use, reuse or unknown.
- H9: **If action reported to FDA under 21 U.S.C. 360i(f), list correction/removal reporting number**
Enter the number that FDA assigned to the corrective action. If a number has not yet been assigned by FDA, the internal number assigned by the firm for the action may be used.
- H10: **Additional manufacturer narrative** - Enter any additional information, evaluation, or clarification of data presented in previous sections. Do not duplicate information that has already been provided elsewhere.

SECTION H: DEVICE MANUFACTURERS ONLY (continued)

H11: **Corrected data** - Provide the following additional, corrected or missing information, identifying each data item by the applicable section and block number:

(1) any information missing on the user facility or distributor report, including any missing or incomplete event codes required by block F10,

(2) information corrected on the user facility or distributor report form after verification, including any corrected event codes required by block F10,

(3) for each event code provided in block F10, a statement of whether the type of event represented by the code is addressed in the device labeling, e.g., code # 1738 - labeled, code # 1701 - not labeled, and

(4) an explanation of why any required information was not provided and the steps taken to obtain such information.

INSTRUCTIONS FOR USING THE DEVICE CODING MANUAL FOR FORM 3500A

The Device Coding Manual is to be used to fill out certain sections of Form 3500A when reporting adverse events and product problems associated with the use of medical devices. It has been divided into two parts. Part 1 is to be used by User Facilities when filling out Section F, block F10 of the form, and Part 2 is to be used by Device Manufacturers when filling out Section H, blocks H3 and H6 of the form. When entering the applicable codes, place one entire code per box, not one digit per box. Instructions for using each part are detailed below.

PART I - EVENT PROBLEM CODES AND TERMS (For Section F, Block F10 of Form 3500A)

This part contains patient- and device-related event terms and corresponding event codes that, when selected, will best describe what happened to the patient as a result of the event or the device failure or problem at the time of the event. At least one code should be selected for patient and/or device. If more than 3 codes are needed for either category, a supplemental page containing those additional codes, may be attached following the general instructions for completing form 3500A. In order to make this part more flexible, it has been divided into three subparts, A, B and C.

SUBPART A - ALPHABETIC LISTING OF EVENT TERMS WITH THEIR CORRESPONDING EVENT CODE

This subpart is divided into two lists. One list contains patient-related event terms and the other device-related event terms. Each list is alphabetically ordered by event term with its corresponding numeric event code.

SUBPART B - NUMERIC LISTING OF EVENT CODES WITH THEIR CORRESPONDING EVENT TERM

This subpart contains a list of all patient- and device-related event terms, numerically ordered by event code.

SUBPART C - DEFINITIONS FOR EVENT TERMS LISTED NUMERICALLY BY EVENT CODE

This subpart contains a glossary of all patient- and device-related event terms that are not self-explanatory, and is intended to be used as a reference guide for Subpart A. Only one list is provided in this part with both patient- and device-related event term definitions ordered numerically by event code.

PART II - DEVICE MANUFACTURER CODES (For Section H, Blocks H3 and H6 of Form 3500A)

This part contains the codes that are to be used by Device Manufacturers when reporting the results of their evaluation of the suspect or related medical device. The evaluation codes for this part have been divided into two subparts, A and B. The sections within each subpart, are alphabetically ordered by evaluation term with its corresponding numeric evaluation code.

PART II - DEVICE MANUFACTURER CODES (For Section H, Blocks H3 and H6 of Form 3500 A)
(continued)

SUBPART A - DEVICE EVALUATED BY MANUFACTURER (Block H3)

The codes from this section are designed to be used when the device was not evaluated. The codes selected should indicate why no evaluation was performed.

SUBPART B - EVALUATION CODES (Block H6)

The codes from this section are divided into three groups, Method of Evaluation, Results of Evaluation and Conclusions. When entering the applicable codes, place one entire code per box, not one digit per box. Select only those evaluation codes that best identify the primary contributing factors. Four boxes for each group have been provided for your codes, however, if more are needed, a supplemental page may be attached following the general instructions for completing form 3500A. Indicate the appropriate section and block number and the additional codes.

METHOD OF EVALUATION CODES

There are two lists in this group, one to indicate the source of the device evaluated and the other to indicate the type of evaluation performed. At least one code from each of the two lists should be selected.

RESULTS OF EVALUATION CODES

There are five categories in this group. The first three, A, B and C are broad characterizations of what the evaluation found. The second two, D and E are very specific in nature.

Select one or more codes from at least one of the first three categories. You may select codes from more than one category, if applicable.

- Category A - Device
- Category B - Use of Device
- Category C - Physiological/Procedural Factors

When the result involves component failures and computer related problems, select as many of the applicable codes as needed from the following categories and enter them on the same line:

- Category D - Device Component/Subassembly Failures
- Category E - Computer-, Imaging System-, Microprocessor-Controlled Device Problems

CONCLUSIONS CODES

In this group, select the code(s) that best describe your conclusions.

IMPORTANT REFERENCE INFORMATION AND MAILING ADDRESSES

HOW TO OBTAIN COPIES OF FDA FORM 3500A AND INSTRUCTIONS FOR COMPLETING FORM 3500A WITH THE DEVICE CODING MANUAL:

Bulk copies of FDA Form 3500A can be obtained from:

Consolidated Forms and Publications Distribution Center
Washington Commerce Center
3222 Hubbard Road
Landover, Maryland 20785

Ten copies or less of the form and/or copies of the "Instructions for Completing Form 3500A with the Device Coding Manual for Form FDA 3500A" can be obtained from:

Division of Small Manufacturers Assistance (HFZ-220)
Center for Devices and Radiological Health
Food and Drug Administration
5600 Fishers Lane
Rockville, MD 20857

Fax number: 301-443-8818

DSMA FACTS-ON-DEMAND fax number: 1-800-899-0381 or 301-827-0111.

Copies of a blank FDA Form 3500A may be duplicated by the applicant.

The "Instructions for Completing Form 3500A with the Device Coding Manual for Form 3500A" are also available through the Center for Devices and Radiological Health's (CDRH) Electronic Docket/Home Page at:

Via modem: 1-800-252-1366, 301-594-2741 or 1-800-222-0185

Telnet Access Via Internet = fedworld.gov (192.239.92.3)

World Wide Web (Home Page) = www.fda.gov/cdrh/cdrhhome.html (this is the URL)

FOR APPROVAL OF COMPUTER-GENERATED FACSIMILES OF 3500A, CONTACT:

MEDWATCH: The FDA Medical Products Reporting Program
Office of the Commissioner HF-2
Food and Drug Administration, Rm. 9-57
5600 Fishers Lane
Rockville, MD 20857

SEND COMPLETED FORM 3500A TO:

Food and Drug Administration
Center for Devices and Radiological Health
Medical Device Reporting
P.O. Box 3002
Rockville, MD 20847-3002

Each envelope containing a completed Form 3500A should be specifically identified as a "User Facility Report", "Distributor Report", "Manufacturer Report", "Five-Day Report", etc.

MEDWATCH

THE FDA MEDICAL PRODUCTS REPORTING PROGRAM

For use by user-facilities,
distributors and manufacturers for
MANDATORY reporting

Page ____ of ____

Form Approved: OMB No. 0910-0291 Expires: 12/31/94
See OMB statement on reverse

Mfr report #
UR/Dist report #
FDA Use Only

A. Patient information

1. Patient identifier In confidence	2. Age at time of event: or _____ Date of birth: _____	3. Sex <input type="checkbox"/> female <input type="checkbox"/> male	4. Weight ____ lbs or ____ kgs
--	--	--	---

B. Adverse event or product problem

1. <input type="checkbox"/> Adverse event and/or <input type="checkbox"/> Product problem (e.g., defects/malfunctions)	
2. Outcomes attributed to adverse event (check all that apply)	
<input type="checkbox"/> death _____ (mo/day/yr)	<input type="checkbox"/> disability
<input type="checkbox"/> life-threatening	<input type="checkbox"/> congenital anomaly
<input type="checkbox"/> hospitalization – initial or prolonged	<input type="checkbox"/> required intervention to prevent permanent impairment/damage
	<input type="checkbox"/> other: _____
3. Date of event (mo/day/yr)	4. Date of this report (mo/day/yr)

5. Describe event or problem

6. Relevant tests/laboratory data, including dates

7. Other relevant history, including preexisting medical conditions (e.g., allergies, race, pregnancy, smoking and alcohol use, hepatic/renal dysfunction, etc.)

C. Suspect medication(s)

1. Name (give labeled strength & mfr/labeler, if known)	
#1	
#2	
2. Dose, frequency & route used	3. Therapy dates (if known, give duration) from/to (or best estimate)
#1	#1
#2	#2
4. Diagnosis for use (indication)	5. Event abated after use stopped or dose reduced
#1	#1 <input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> doesn't apply
#2	#2 <input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> doesn't apply
6. Lot # (if known)	7. Exp. date (if known)
#1	#1
#2	#2
9. NDC # – for product problems only (if known)	8. Event reappeared after reintroduction
-	#1 <input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> doesn't apply
	#2 <input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> doesn't apply
10. Concomitant medical products and therapy dates (exclude treatment of event)	

D. Suspect medical device

1. Brand name	
2. Type of device	
3. Manufacturer name & address	4. Operator of device
	<input type="checkbox"/> health professional <input type="checkbox"/> lay user/patient <input type="checkbox"/> other: _____
6. model #	5. Expiration date (mo/day/yr)
catalog #	7. If implanted, give date (mo/day/yr)
serial #	8. If explanted, give date (mo/day/yr)
lot #	
other #	
9. Device available for evaluation? (Do not send to FDA)	
<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> returned to manufacturer on _____ (mo/day/yr)	
10. Concomitant medical products and therapy dates (exclude treatment of event)	

E. Initial reporter

1. Name, address & phone #		
2. Health professional?	3. Occupation	4. Initial reporter also sent report to FDA
<input type="checkbox"/> yes <input type="checkbox"/> no		<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> unk



FDA Form 3500A (6/93)

Submission of a report does not constitute an admission that medical personnel, user facility, distributor, manufacturer or product caused or contributed to the event.

Medication and Device Experience Report

(continued)

Refer to guidelines for specific instructions

Submission of a report does not constitute an admission that medical personnel, user facility, distributor, manufacturer or product caused or contributed to the event.

U.S. DEPARTMENT OF HEALTH AND HUMAN SERVICES
Public Health Service • Food and Drug Administration

Page ____ of ____

FDA Use Only

F. For use by user facility/distributor—devices only

1. Check one <input type="checkbox"/> user facility <input type="checkbox"/> distributor		2. UF/Dist report number	
3. User facility or distributor name/address			
4. Contact person		5. Phone Number	
6. Date user facility or distributor became aware of event (mo/day/yr)		8. Date of this report (mo/day/yr)	
7. Type of report <input type="checkbox"/> initial <input type="checkbox"/> follow-up # _____			
9. Approximate age of device		10. Event problem codes (refer to coding manual)	
		patient code _____ - _____ - _____ device code _____ - _____ - _____	
11. Report sent to FDA? <input type="checkbox"/> yes _____ (mo/day/yr) <input type="checkbox"/> no		12. Location where event occurred <input type="checkbox"/> hospital <input type="checkbox"/> outpatient diagnostic facility <input type="checkbox"/> home <input type="checkbox"/> ambulatory surgical facility <input type="checkbox"/> nursing home <input type="checkbox"/> outpatient treatment facility <input type="checkbox"/> other: _____ specify _____	
13. Report sent to manufacturer? <input type="checkbox"/> yes _____ (mo/day/yr) <input type="checkbox"/> no			
14. Manufacturer name/address			

G. All manufacturers

1. Contact office — name/address (& mfring site for devices)		2. Phone number	
4. Date received by manufacturer (mo/day/yr)		3. Report source (check all that apply) <input type="checkbox"/> foreign <input type="checkbox"/> study <input type="checkbox"/> literature <input type="checkbox"/> consumer <input type="checkbox"/> health professional <input type="checkbox"/> user facility <input type="checkbox"/> company representative <input type="checkbox"/> distributor <input type="checkbox"/> other: _____	
6. If IND, protocol #		5. (A)NDA # _____ IND # _____ PLA # _____ pre-1938 <input type="checkbox"/> yes OTC product <input type="checkbox"/> yes	
7. Type of report (check all that apply) <input type="checkbox"/> 5-day <input type="checkbox"/> 15-day <input type="checkbox"/> 10-day <input type="checkbox"/> periodic <input type="checkbox"/> Initial <input type="checkbox"/> follow-up # _____		8. Adverse event term(s)	
9. Mfr. report number			

H. Device manufacturers only

1. Type of reportable event <input type="checkbox"/> death <input type="checkbox"/> serious injury <input type="checkbox"/> malfunction (see guidelines) <input type="checkbox"/> other: _____		2. If follow-up, what type? <input type="checkbox"/> correction <input type="checkbox"/> additional information <input type="checkbox"/> response to FDA request <input type="checkbox"/> device evaluation	
3. Device evaluated by mfr? <input type="checkbox"/> not returned to mfr. <input type="checkbox"/> yes <input type="checkbox"/> evaluation summary attached <input type="checkbox"/> no (attach page to explain why not) or provide code: _____		4. Device manufacture date (mo/yr)	
		5. Labeled for single use? <input type="checkbox"/> yes <input type="checkbox"/> no	
6. Evaluation codes (refer to coding manual)			
method _____ - _____ - _____ - _____			
results _____ - _____ - _____ - _____			
conclusions _____ - _____ - _____ - _____			
7. If remedial action initiated, check type <input type="checkbox"/> recall <input type="checkbox"/> notification <input type="checkbox"/> repair <input type="checkbox"/> inspection <input type="checkbox"/> replace <input type="checkbox"/> patient monitoring <input type="checkbox"/> relabeling <input type="checkbox"/> modification/adjustment <input type="checkbox"/> other: _____		8. Usage of device <input type="checkbox"/> initial use of device <input type="checkbox"/> reuse <input type="checkbox"/> unknown	
9. If action reported to FDA under 21 USC 360(i)(f), list correction/removal reporting number:			
10. <input type="checkbox"/> Additional manufacturer narrative and/or 11. <input type="checkbox"/> Corrected data			

The public reporting burden for this collection of information has been estimated to average one-hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send your comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Reports Clearance Officer, PHS
Hubert H. Humphrey Building, Room 721-B
200 Independence Avenue, S.W.
Washington, DC 20201
ATTN: PRA

and to:
Office of Management and Budget
Paperwork Reduction Project (0910-0291)
Washington, DC 20503

Please do NOT return this form to either of these addresses.

PART I, SUBPART A - ALPHABETIC LISTING OF EVENT TERMS WITH THEIR CORRESPONDING NUMERIC EVENT CODE

PATIENT-RELATED TERMS

1685	Abdominal pain	2472	Aortic valve replacement
1034	Abnormal blood gas measurements	1805	Apgar, decreased
2370	Abnormal mammogram	1718	Apgar score, decreased or low
2078	Abnormality of ventricle	1719	Apheresis
1688	Abortion	2372	Apicectomy
2211	Abortion, artificial	1720	Apnea
2212	Abortion, complete	2390	Arachnoiditis, spinal
2213	Abortion, incomplete	1696	ARDS, adult respiratory distress syndrome
2214	Abortion, induced	2351	Arm fracture
2215	Abortion, missed	2353	Arm/hand hypoesthesia
2216	Abortion, therapeutic	1762	Arrest, cardiac
1689	Abrasion	1765	Arrest, cardiopulmonary
1789	Abrasion, corneal	2044	Arrest, respiratory
1690	Abscess	1721	Arrhythmia
2261	Absence seizures	2480	Arterial wedge pressure, high pulmonary
1945	Acanthamoeba keratitis	2479	Arterial wedge pressure, low pulmonary
1770	Accident (CVA), cerebrovascular	2481	Arterial wedge pressure, normal pulmonary
1692	Achalasia	1722	Arteriosclerosis
2482	Acidosis, respiratory	2355	Arthralgia
1693	Acoustic shock	1723	Arthritis
1694	Acoustic trauma	1724	Arthritis, rheumatoid
1684	Acquired immunodeficiency syndrome, (AIDS)	2211	Artificial abortion
2269	Acquired toxoplasmosis	1851	Asphyxia
1695	Adhesion(s)	2209	Asphyxia, fetal
1696	Adult respiratory distress syndrome, (ARDS)	1725	Aspiration
1684	AIDS, acquired immunodeficiency syndrome	1726	Asthma
1697	Air embolism	1727	Asystole
1698	Airflow, restricted	1728	Atherosclerosis
1699	Airway obstruction	1729	Atrial fibrillation
2483	Alkalosis, respiratory	1730	Atrial flutter
1701	Allergic reaction	1731	Atrial tachycardia
2217	Amnionitis	2109	Attack, transient ischemic
1702	Amputation	2058	Aureus, staphylococcus
2015	ANA, positive antinuclear antibodies	1732	Autoimmune disease
1703	Anaphylactic shock	1733	Autoimmune reaction
2218	Anaphylactoid	1735	Bacterial infection
1704	Anaphylaxis	1736	Biliary cirrhosis
1028	Anastomose, failure to	1738	Bleeding
1705	Anastomosis	1739	Bleeding, cerebral
1706	Anemia	1741	Bleeding, intracranial
2279	Anemia, hemolytic	1742	Blindness
1707	Anesthesia, insufficient, light or patchy	2229	Blinking, excessive
1708	Aneurysm	1743	Blister(s)
1710	Angina	1744	Bloating
1711	Anoxia	2026	Block, pupillary
1713	Antibiotics, reaction to	2332	Blockage
2015	Antibodies (ANA), positive antinuclear	1034	Blood gas measurements, abnormal
1714	Anticoagulation	2451	Blood glucose, high
2015	Antinuclear antibodies (ANA), positive	2450	Blood glucose, low
2328	Anxiety	2143	Blood in vitreous fluid
1715	Aortic insufficiency	1747	Blood pooling
1716	Aortic regurgitation	2488	Blood pressure, high
1717	Aortic stenosis	2489	Blood pressure, low

PATIENT-RELATED TERMS (continued)

2337	Blood products, transfusion of	1773	Changes, cervical
2025	Blood vessels, punctured	1831	Changes, emotional
1749	Blood, transfusion with incompatible	2455	Chemical contamination
1745	Blood/fluids, contamination from	1775	Chemosis
2137	Blurred vision	1776	Chest pain
2230	Blurring	2463	Chest tightness/pressure
2461	Body temperature, decreased	2098	Children, toxins in
2096	Body temperature, elevated	2191	Chills
1803	Bone shedding debris	2464	Choking
1756	Bowel burn	1777	Chorioamnionitis
1751	Bradycardia	2237	Chronic obstructive pulmonary disease, (COPD)
2394	Brain damage	2060	Cicatrix
2219	Brain injury	1736	Cirrhosis, biliary
1852	Brain injury, fetal	2222	Clonic convulsion
1388	Breakage of membrane	2228	Clouding, central corneal
2106	Breakdown of optical tissue	2220	CMV, cytomegalovirus
1759	Breast cancer	1778	Coagulation
2439	Breast lumps	1813	Coagulation (DIC), disseminated intravascular
2438	Breast neoplasm	1779	Coagulopathy
2401	Breathing difficulties	2189	Collagen disease
1593	Breaths, stacking	2416	Collapse
1752	Bronchitis	2417	Coma
2437	Bronchopneumonia	1781	Comatose
1754	Bruise	2331	Complaint, ill-defined
1757	Burn(s)	2212	Complete abortion
1756	Burn, bowel	2192	Concussion
1755	Burn, radiation	2155	Condition of patient, improved
2146	Burning sensation	1782	Congenital defect/deformity
1774	C-section (cesarean section), delivery by	2270	Congenital toxoplasmosis
1758	Calcification	1783	Congestive heart failure
2147	Calcium deposit(s)	1784	Conjunctivitis
1759	Cancer, breast	1785	Conjunctivitis, giant papillary
1760	Cancer, other	1786	Connective tissue disease
1761	Capsular contracture	2418	Consciousness, loss of
1551	Capsular contracture, rupture due to stress from	2199	Consequences or impact to patient, none
1762	Cardiac arrest	1745	Contamination from blood/fluids
1763	Cardiac contusion	2455	Contamination, chemical
1838	Cardiac enzyme elevation	1761	Contracture, capsular
2334	Cardiac insufficiency	1125	Contraindicated patient
2226	Cardiac tamponade	1787	Contusion
2262	Cardiogenic shock	1763	Contusion, cardiac
1764	Cardiomyopathy	1953	Contusion, liver
1765	Cardiopulmonary arrest	1763	Contusion, myocardial
1766	Cataract	2221	Convulsion
1767	Cataract, induced	2222	Convulsion, clonic
1768	Cellulitis	2223	Convulsion, tonic
2228	Central corneal clouding	2237	COPD, chronic obstructive pulmonary disease
1739	Cerebral bleeding	1788	Core temperature rise
1889	Cerebral hemorrhage	1854	Core temperature rise, fetal
1771	Cerebral, infarction	1792	Cornea, perforation of
1772	Cerebrospinal fluid leakage	1789	Corneal abrasion
1770	Cerebrovascular accident, (CVA)	2228	Corneal clouding, central
1773	Cervical changes	1790	Corneal decompensation
1927	Cervix, incompetent	1791	Corneal edema
1774	Cesarean section (c-section)	2173	Corneal implant
1817	Changes in EKG/ECG	2231	Corneal infiltrates
2162	Changes in/loss of nipple sensation		

PATIENT-RELATED TERMS (continued)

1793	Corneal scar	2047	Detached retina
2233	Corneal sensitivity	2445	Detachment of vitreous
1794	Corneal touch	1204	Device embedded in vessel or plaque
1795	Corneal transplant	1810	Device induced
1796	Corneal ulcer	2364	Diabetic ketoacidosis
2193	Cramp(s)	2406	Diagnosis, identified
1797	Crushing injury	2224	Dialysis dementia
2453	Cut(s)	2452	Diaphoretic
1417	Cuts, nicks or tears of dura or other tissues by device	1811	Diarrhea
1770	CVA, cerebrovascular accident	1813	DIC, disseminated intravascular coagulation
1798	Cyanosis	2401	Difficulty breathing
1799	Cyclitis	2070	Digits (finger or toe), severed
1800	Cyst(s), formation of	2244	Direct infection
1957	Cytomegaloviral retinitis	2371	Disability
2220	Cytomegalovirus (CVM)	2225	Discharge
1952	Damage to ligament(s)	1812	Discharge, purulent
1979	Damage to nerve(s)	2123	Discharge, vaginal
1981	Damage to neural tissue	2074	Discoloration of skin
1986	Damage to optical nerve(s)	2330	Discomfort
2048	Damage to retina	2237	Disease (COPD), chronic obstructive pulmonary
2394	Damage, brain	2039	Disease (ESRD), end stage renal
2104	Damage, tissue	2000	Disease (PID), pelvic inflammatory
2124	Damage, vaginal mucosa	2071	Disease (VD), venereal
1801	Deafness	1732	Disease, autoimmune
1855	Death, intrauterine fetal	2189	Disease, collagen
1802	Death/expired	1786	Disease, connective tissue
1803	Debris, bone shedding	2002	Disease, peripheral vascular
1790	Decompensation, corneal	2053	Disease, rheumatic heart
1805	Decreased apgar	2360	Disfigurement
2461	Decreased body temperature	2374	Dislocated joint
2430	Decreased forced expiratory volume (FEV)	2373	Disorder, joint
1718	Decreased or low apgar score	1353	Disruption, leaflet
2435	Decreased peak expiratory flowrate	1333	Dissection, intimal
2485	Decreased respiratory rate	1813	Disseminated intravascular coagulation (DIC)
2271	Decreased therapeutic response	1826	Dissociation, electro-mechanical
1782	Defect/deformity, congenital	2329	Distress
1982	Deficit/dysfunction, neurological	1856	Distress, fetal
1974	Deformities, neonatal	2045	Distress, respiratory
1782	Deformity/defect, congenital	2140	Disturbances, visual
1806	Degeneration	2194	Dizziness
2049	Degeneration of retina	1814	Dry eye(s)
1154	Dehiscence wound	1417	Dura nicks, cuts or tears by device
1807	Dehydrated	1947	Dysfunction, left ventricular
1156	Delamination	1954	Dysfunction, liver
2369	Delayed union fracture	2019	Dysfunction, pulmonary
1936	Delayed, uncontrolled intraocular pressure (IOP)	2054	Dysfunction, right ventricular
1774	Delivery by cesarean section (c-section)	1982	Dysfunction/deficit, neurological
1808	Dementia	2363	Dyskinesia
2224	Dementia, dialysis	1815	Dysphagia
2395	Dependent, ventilator	2195	Dysphasia
2147	Deposit(s), calcium	1816	Dyspnea
1809	Deposits	2103	Ears, ringing in
2361	Depression	1818	Ecchymosis
2487	Depression, ST segment	1817	ECG/EKG changes
2253	Dermatomyositis	1819	Ectopic pregnancy
		1820	Edema

PATIENT-RELATED TERMS (continued)

1791	Edema, corneal	1783	Failure, congestive heart
1822	Edema, macular	2206	Failure, heart
1823	Edema, microcytic	1948	Failure, left ventricular
2020	Edema, pulmonary	2041	Failure, renal
1824	Edema, stromal	2484	Failure, respiratory
2273	Effects, teratogenic	2055	Failure, right ventricular
2099	Effects, unexpected therapeutic	2105	Failure, tissue
2010	Effusion, pleural	1610	Faint
1817	EKG/ECG changes	1847	Fainting
1826	Electro-mechanical dissociation	1848	Fall
1827	Electrocution	2375	Fasciitis
2196	Electrolyte imbalance	1849	Fatigue
2096	Elevated body temperature	1850	Feeding problems
1838	Elevation of cardiac enzymes	2209	Fetal asphyxia
2059	Elevation, ST segment	1852	Fetal brain injury
1829	Embolism	1854	Fetal core temperature rise
1697	Embolism, air	1855	Fetal death, intrauterine
1498	Embolism, pulmonary	1856	Fetal distress
1830	Embolus	2210	Fetal hypoxia
1831	Emotional changes	1857	Fetal scalp laceration(s)
1832	Emphysema, pulmonary	2430	FEV (forced expiratory volume), decreased
2429	Encephalitis	2431	FEV, (forced expiratory volume), increased
1833	Encephalopathy	1858	Fever
2039	End stage renal disease (ESRD)	1859	Fiberoptic fragments, unretrieved in body
1834	Endocarditis	1729	Fibrillation, atrial
1835	Endophthalmitis	2130	Fibrillation, ventricular
2236	Endotoxin	1860	Fibromyitis
2327	Entrapment	1861	Fibrosis
1838	Enzyme elevation, cardiac	2070	Finger or toe, severed
2150	Epithelial marsupialization	2239	First use syndrome
2232	Epithelial microcysts	1862	Fistula
1750	Erosion	1972	Flap tissue, necrosis of
2013	Erosion, pocket	1864	Flashers
2075	Erosion, skin	2153	Flashes, hot
1839	Eructate	1865	Flatus
1840	Erythema	1866	Floaters, vitreous
2398	Esophagus, laceration(s) of	2052	Flow, retrograde
2399	Esophagus, perforation of	2435	Flowrate, decreased peak expiratory
2039	ESRD, end stage renal disease	2436	Flowrate, increased peak expiratory
2229	Excessive blinking	2143	Fluid, blood in vitreous
2235	Excessive tearing	1772	Fluid, leakage of cerebrospinal
2358	Excision, scar	2142	Fluid, loss of vitreous
1843	Exhaustion, extreme	1745	Fluids/blood, contamination from
1802	Expired/death	1730	Flutter, atrial
1841	Exsanguination	2131	Flutter, ventricular
1890	Extradural hemorrhage	2064	Focal motor seizures
1842	Extravasation	2260	Focal seizures
1843	Extreme exhaustion	2354	Foot/leg hypoaesthesia
1844	Extrusion	2430	Forced expiratory volume (FEV), decreased
2402	Extubate	2431	Forced expiratory volume (FEV), increased
1845	Eye injury	1868	Foreign body reaction
1814	Eye(s), dry	1869	Foreign body sensation
2038	Eye(s), red	2365	Foreign body, removal of
1901	Face, rupture of hyaloid	1800	Formation of cyst(s)
1846	Facial nerve paralysis	1447	Formation, pannus
1924	Failure of implant	1870	Fracture(s)
1028	Failure to anastomose	2351	Fracture, arm

PATIENT-RELATED TERMS (continued)

2369	Fracture, delayed union	2478	High oxygen saturation
2349	Fracture, hip	2480	High pulmonary arterial wedge pressure
2077	Fracture, skull	2349	Hip fracture
2428	Fracture, tooth	2197	HIV, human immunodeficiency virus
1859	Fragments unretrieved in body, fiberoptic	1900	Hives
2275	Frequency, urinary	1297	Host-tissue reaction
2419	Fungal infection	2153	Hot flashes
1872	Fungus	2197	Human immunodeficiency virus (HIV)
1873	Gangrene	1901	Hyaloid face, rupture of
1874	Gastritis	1903	Hyperbilirubinemia
1785	Giant papillary conjunctivitis	1904	Hyperemia
2092	Glands, swollen	1905	Hyperglycemia
1875	Glaucoma	2242	Hypernatremia
2451	Glucose, high blood	1906	Hyperplasia
2450	Glucose, low blood	1907	Hypersensitivity
2168	Grand-mal seizures	1908	Hypertension
1876	Granuloma	1909	Hyperthermia
2152	Great vessel perforation	1949	Hypertrophy, left ventricular
1877	Hair loss	2056	Hypertrophy, right ventricular
2227	Halo	1910	Hyperventilation
2353	Hand/arm hypoesthesia	1911	Hyphema
1878	Haze	2352	Hypoesthesia
1879	Head injury	2353	Hypoesthesia, arm/hand
1880	Headache	2354	Hypoesthesia, foot/leg
2186	Headache, lumbar puncture	1912	Hypoglycemia
2186	Headache, post spinal	1913	Hypopyon
2378	Healing, impaired	1914	Hypotension
1881	Hearing impairment	1915	Hypothermia
1975	Hearing impairment, neonatal	1916	Hypoventilation
1882	Hearing loss	2243	Hypovolemia
1976	Hearing loss, neonatal	1917	Hypovolemic shock
2053	Heart disease, rheumatic	1918	Hypoxia
2206	Heart failure	2210	Hypoxia, fetal
1783	Heart failure, congestive	1922	Iatrogenic lesion
2476	Heart valve repair	2406	Identified diagnosis
1883	Heartburn	1923	Idioventricular rhythm
1884	Hematoma	2331	Ill-defined complaint
1288	Hemoconcentration	2196	Imbalance, electrolyte
1885	Hemodialysis	2156	Immuno-deficiency
1886	Hemolysis	2199	Impact or consequences to patient, none
2279	Hemolytic anemia	2378	Impaired healing
1887	Hemoptysis	2138	Impaired vision
1888	Hemorrhage	1975	Impairment to hearing, neonatal
1889	Hemorrhage, cerebral	1881	Impairment, hearing
1890	Hemorrhage, extradural	2173	Implant, corneal
1891	Hemorrhage, intracranial	1924	Implant, failure of
1892	Hemorrhage, intraventricular	1925	Impotence
1893	Hemorrhage, subarachnoid	2155	Improved patient condition
1894	Hemorrhage, subdural	2388	Inadequate pain relief
1895	Hemostasis	1749	Incompatible blood, transfusion with
1896	Hemothorax	1927	Incompetent cervix
1897	Hepatitis	2213	Incomplete abortion
2240	Hernia	1928	Incontinence
1898	Herpes	2431	Increased forced expiratory volume (FEV)
1899	Hiccups	2436	Increased peak expiratory flowrate
2451	High blood glucose	2486	Increased respiratory rate
2488	High blood pressure	2272	Increased therapeutic response

PATIENT-RELATED TERMS (continued)

2245	Indirect infection	1937	IOPR, intraocular pressure rise
2214	Induced abortion	1940	Iritis
1810	Induced by device	2469	Irregular pulse
1767	Induced cataract	1959	Irregularities, menstrual
1969	Infarction (MI), myocardial	2421	Irritability
1771	Infarction, cerebral	1941	Irritation
2021	Infarction, pulmonary	2076	Irritation of skin
1930	Infection	1942	Ischemia
1735	Infection, bacterial	2109	Ischemic attack, transient
2244	Infection, direct	1943	Itching
2419	Infection, fungal	2187	Jaundice
2245	Infection, indirect	2374	Joint dislocation
1933	Infection, intraocular	2373	Joint disorder
2446	Infection, post-operative wound	2355	Joint pain
2447	Infection, post-traumatic wound	2356	Joint swelling
2246	Infection, pyrogenic	1944	Keratitis
2247	Infection, subclinical	1945	Keratitis, acanthamoeba
2420	Infection, upper respiratory tract	2364	Ketoacidosis, diabetic
2120	Infection, urinary tract	2465	Labor, premature
2248	Infection, viral	1946	Laceration(s)
2174	Infiltrates	2398	Laceration(s) of esophagus
2231	Infiltrates, corneal	1857	Laceration(s), fetal scalp
1931	Infiltration	1955	Laceration(s), liver
2087	Infiltration, subepithelial	2003	Laceration(s), peritoneal
1932	Inflammation	1341	LAL test (limulus amebocyte lysate)
2443	Inflammation of skin	2366	Laparotomy
2000	Inflammatory disease (PID), pelvic	1353	Leaflet disruption
2442	Injection site reaction	2201	Leak(s), perivalvular
2348	Injury	1772	Leakage of cerebrospinal fluid
2219	Injury to the brain	1947	Left ventricular dysfunction
1797	Injury, crushing	1948	Left ventricular failure
1845	Injury, eye	1949	Left ventricular hypertrophy
1852	Injury, fetal brain	2354	Leg/foot hypoaesthesia
1879	Injury, head	1950	Lesion
2081	Injury, spinal	1922	Lesion, iatrogenic
2432	Injury, spinal cord	1951	Lesions in the target points
1715	Insufficiency, aortic	1952	Ligament(s), damage to
2334	Insufficiency, cardiac	1707	Light, patchy or insufficient anesthesia
1963	Insufficiency, mitral	1341	Limulus amebocyte lysate (LAL) test
2022	Insufficiency, pulmonary	1953	Liver contusion
2111	Insufficiency, tricuspid	1954	Liver dysfunction
1926	Insufficiency, valvular	1955	Liver laceration(s)
1707	Insufficient, light or patchy anesthesia	2035	Local reaction
2264	Insulin shock	2418	Loss of consciousness
1333	Intimal dissection	1877	Loss of hair
1741	Intracranial bleeding	1882	Loss of hearing
1891	Intracranial hemorrhage	1976	Loss of hearing, neonatal
1933	Intraocular infection	1958	Loss of memory
1936	Intraocular pressure (IOP), delayed, uncontrolled	2032	Loss of range of motion
1937	Intraocular pressure rise, (IOPR)	2139	Loss of vision
1855	Intrauterine fetal death	2142	Loss of vitreous fluid
1938	Intravasation	2162	Loss of/changes in nipple sensation
1813	Intravascular coagulation (DIC), disseminated	2450	Low blood glucose
1892	Intraventricular hemorrhage	2489	Low blood pressure
1936	IOP (intraocular pressure), delayed, uncontrolled	1718	Low or decreased apgar score
		2477	Low oxygen saturation
		2479	Low pulmonary arterial wedge pressure

PATIENT-RELATED TERMS (continued)

2186	Lumbar puncture headache	1986	Nerve damage, optical
2439	Lumps, breast	2362	Nerve palsy, peroneal
2397	Lung, overinflation of	1846	Nerve paralysis, facial
1956	Lupus	1980	Nerve stimulation, undesired
2093	Lymph nodes, swollen	1981	Neural tissue damage
1822	Macular edema	2265	Neurogenic shock
2359	Malaise	1982	Neurological deficit/dysfunction
2370	Mammogram, abnormal	1983	Neuropathy
2150	Marsupialization, epithelial	2046	Newborns (RDS), respiratory distress syndrome of
2276	Materials (nonresorbable), unretrieved in body	1417	Nicks, cuts or tears of dura or other tissues by device
1034	Measurements, abnormal blood gas	2162	Nipple sensation, changes in/loss of
2198	Mediastinal shift	2441	Nipple ulceration
2117	Medication, not under	2199	No consequences or impact to patient
2014	Megophthalmos	2093	Nodes, swollen lymph
1388	Membrane, breakage of	2250	Nonpyrogenic
1958	Memory loss	2276	Nonresorbable materials, unretrieved in body
2389	Meningitis	2411	Nonspecific
1959	Menstrual irregularities	2481	Normal pulmonary arterial wedge pressure
1969	MI, myocardial infarction	2415	Numbness
1960	Microcyst(s)	2422	Obstruction
2232	Microcysts, epithelial	1699	Obstruction of airway
1823	Microcytic edema	2466	Obstruction, nasal
1961	Microperforation(s)	2237	Obstructive pulmonary disease (COPD), chronic
2249	Migration	1984	Occlusion
1962	Miscarriage	1986	Optical nerve damage
2159	Misdiagnosis	2106	Optical tissue, breakdown of
2215	Missed abortion	1987	Organ(s), perforation of
1963	Mitral insufficiency	2110	Organ/tissue transplant
1964	Mitral regurgitation	1428	Ossification
1965	Mitral stenosis	2377	Osteolysis
2471	Mitral valve prolapse	2200	Other (for use when an appropriate patient code cannot be identified)
2376	Mitral valve replacement	1436	Overcorrection
2032	Motion, loss of range of	1988	Overdose
2124	Mucosa damage, vaginal	1510	Overdose of radiation
1966	Muscle spasm(s)	2397	Overinflation of lung
1412	Muscle stimulation	1989	Overmedicated
1967	Muscle weakness	1990	Oversedated
1968	Muscular rigidity	1991	Overstimulation
2161	Muscular tics	2234	Overwear syndrome
2238	Myalgia	2478	Oxygen saturation, high
1763	Myocardial contusion	2477	Oxygen saturation, low
1969	Myocardial infarction, (MI)	1994	Pain
2470	Myocarditis	1776	Pain in chest
2466	Nasal obstruction	2388	Pain relief, inadequate
1970	Nausea	1685	Pain, abdominal
2433	Neck pain	2355	Pain, joint
2434	Neck stiffness	2433	Pain, neck
1971	Necrosis	1995	Painful stimulation
1972	Necrosis of flap tissue	2468	Pallor
1413	Necrosis, pressure	2467	Palpitations
2462	Needle stick/puncture	2362	Palsy, peroneal nerve
1974	Neonatal deformities	1447	Pannus formation
1975	Neonatal hearing impairment	1785	Papillary conjunctivitis, giant
1976	Neonatal hearing loss		
2438	Neoplasm, breast		
1978	Neovascularization		
1979	Nerve damage		

PATIENT-RELATED TERMS (continued)

1997	Paralysis	2463	Pressure/tightness in chest
1846	Paralysis of facial nerve	1850	Problems, feeding
2448	Paraplegia	2102	Problems, thyroid
1998	Paresis	2042	Procedure repeated, surgical
1707	Patchy, light or insufficient anesthesia	2357	Procedure, surgical
2155	Patient condition, improved	2475	Prolapse
1125	Patient, contraindicated	2471	Prolapse, mitral valve
2199	Patient, no consequences or impact to	2177	Prolonged surgery
2435	Peak expiratory flowrate, decreased	2480	Pulmonary arterial wedge pressure, high
2436	Peak expiratory flowrate, increased	2479	Pulmonary arterial wedge pressure, low
1999	Peeling	2481	Pulmonary arterial wedge pressure, normal
2000	Pelvic inflammatory disease, (PID)	2237	Pulmonary disease (COPD), chronic obstructive
2001	Perforation	2019	Pulmonary dysfunction
1792	Perforation of cornea	2020	Pulmonary edema
2399	Perforation of esophagus	1498	Pulmonary embolism
1987	Perforation of organ(s)	1832	Pulmonary emphysema
2277	Perforation of sinus	2021	Pulmonary infarction
2121	Perforation of uterus	2022	Pulmonary insufficiency
2135	Perforation of vessel(s)	2023	Pulmonary regurgitation
2152	Perforation, great vessel	2024	Pulmonary stenosis
2002	Peripheral vascular disease	2474	Pulmonary valve replacement
2003	Peritoneal laceration(s)	2469	Pulse, irregular
2252	Peritonitis	2462	Puncture/stick from needle
2201	Perivalvular leak(s)	2025	Punctured blood vessels
2362	Peroneal nerve palsy	2026	Pupillary block
2169	Petit-mal seizures	1812	Purulent discharge
2367	Pharyngitis	2027	Pus
2034	Phenomenon, Raynaud's	2028	Pyrogenic
2004	Phlebitis	2246	Pyrogenic infection
2164	Phosphene visualization	2254	Pyrogenic reaction
2165	Phototoxicity	2449	Quadriplegia
2000	PID, pelvic inflammatory disease	1755	Radiation burn
1204	Plaque or vessel, device embedded in	1510	Radiation overdose
2009	Pleomorphism	2256	Radiation sickness syndrome
2010	Pleural effusion	2166	Radiation underdose
2011	Pneumonia	2255	Radiodermatitis
2012	Pneumothorax	2032	Range of motion, loss of
2013	Pocket erosion	2033	Rash
1992	Polymyositis	2485	Rate of respiration, decreased
1747	Pooling of blood	2486	Rate of respiration, increased
2015	Positive antinuclear antibodies (ANA)	2034	Raynaud's phenomenon
2186	Post spinal headache	2046	RDS, respiratory distress syndrome of newborns
2446	Post-operative wound infection	2414	Reaction
2447	Post-traumatic wound infection	1297	Reaction of host-tissue
2266	Postoperative shock	1713	Reaction to antibiotics
1819	Pregnancy, ectopic	1701	Reaction, allergic
2465	Premature labor	1733	Reaction, autoimmune
1936	Pressure (IOP), delayed uncontrolled intraocular	1868	Reaction, foreign body
1937	Pressure (IOPR), rise in intraocular	2442	Reaction, injection site
2326	Pressure sores/ulcers	2035	Reaction, local
2488	Pressure, high blood	2254	Reaction, pyrogenic
2480	Pressure, high pulmonary arterial wedge	2036	Reaction, systemic
2489	Pressure, low blood	2038	Red eye(s)
2479	Pressure, low pulmonary arterial wedge	2259	Regurgitation
1413	Pressure, necrosis	1716	Regurgitation, aortic
2481	Pressure, normal pulmonary arterial wedge		

PATIENT-RELATED TERMS (continued)

1964	Regurgitation, mitral	1857	Scalp laceration(s), fetal
2023	Regurgitation, pulmonary	2358	Scar excision
2112	Regurgitation, tricuspid	2060	Scar tissue
2335	Regurgitation, valvular	1793	Scar, corneal
2403	Reinfusion	2061	Scarring
2392	Reintubate	2062	Scleroderma
2365	Removal of foreign body	1718	Score, decreased or low apgar
2039	Renal disease (ESRD), end stage	2368	Sedation
2041	Renal failure	2118	Sedation, not under
1985	Reocclusion	2063	Seizures
2476	Repair, heart valve	2261	Seizures, absence
2042	Repeated surgical procedure	2260	Seizures, focal
2472	Replacement of aortic valve	2064	Seizures, focal motor
2474	Replacement of pulmonary valve	2168	Seizures, grand-mal
2473	Replacement of tricuspid valve	2169	Seizures, petit-mal
2376	Replacement, mitral valve	2146	Sensation, burning
2482	Respiratory acidosis	2162	Sensation, changes in/loss of nipple
2483	Respiratory alkalosis	1869	Sensation, foreign body
2044	Respiratory arrest	2065	Sensitivity
2045	Respiratory distress	2427	Sensitivity of teeth
1696	Respiratory distress syndrome of adults (ARDS)	2233	Sensitivity, corneal
2046	Respiratory distress syndrome of newborns, (RDS)	2066	Sensitization
2484	Respiratory failure	2067	Sepsis
2485	Respiratory rate, decreased	2068	Septic shock
2486	Respiratory rate, increased	2069	Seroma
2420	Respiratory tract infection, upper	2070	Severed digit(s), (finger or toe)
2271	Response, decreased therapeutic	1803	Shedding of bone debris
2272	Response, increased therapeutic	2198	Shift, mediastinal
1698	Restricted airflow	2072	Shock
2400	Resuscitation	2108	Shock syndrome (TSS), toxic
2119	Retention, urinary	1693	Shock, acoustic
2048	Retina, damage to	1703	Shock, anaphylactic
2049	Retina, degeneration of	2262	Shock, cardiogenic
2047	Retina, detached	1917	Shock, hypovolemic
2050	Retina, tear(s) in	2264	Shock, insulin
1957	Retinitis, cytomegaloviral	2265	Shock, neurogenic
2051	Retrograde	2266	Shock, postoperative
2052	Retrograde flow	2068	Shock, septic
2053	Rheumatic heart disease	2267	Shock, surgical
1724	Rheumatoid arthritis	2268	Shock, traumatic
1923	Rhythm, idioventricular	2256	Sickness syndrome, radiation
2054	Right ventricular dysfunction	2277	Sinus, perforation of
2055	Right ventricular failure	2073	Sjogren's syndrome
2056	Right ventricular hypertrophy	2074	Skin discoloration
1968	Rigidity, muscular	2075	Skin erosion
2103	Ring in ears	2443	Skin inflammation
1788	Rise in core temperature	2076	Skin irritation
1854	Rise in fetal core temperature	2077	Skull fracture
1937	Rise in intraocular pressure, (IOPR)	2251	Sneezing
2208	Rupture	2396	Sore throat
1551	Rupture due to stress from capsular contracture	2079	Soreness
1901	Rupture of hyaloid face	2326	Sores/ulcers, pressure
2058	S. aureus	1966	Spasm(s), muscle
2478	Saturation, high oxygen	2390	Spinal arachnoiditis
2477	Saturation, low oxygen	2432	Spinal cord injury
		2081	Spinal injury
		2082	Spotting

PATIENT-RELATED TERMS (continued)

2083	Sprain	1951	Target points, lesions in
2487	ST segment depression	2050	Tear(s) in retina
2059	ST segment elevation	2235	Tearing, excessive
1593	Stacking breaths	1417	Tears, nicks or cuts of dura or other tissues by device
2178	Staining	2097	Tears, tentorial
2058	Staphylococcus aureus	2427	Teeth, sensitivity of
2263	Stenosis	1788	Temperature rise, core
1717	Stenosis, aortic	1854	Temperature rise, fetal core
1965	Stenosis, mitral	2461	Temperature, decreased body
2024	Stenosis, pulmonary	2096	Temperature, elevated body
2113	Stenosis, tricuspid	2097	Tentorial tears
2462	Stick/puncture from needle	2273	Teratogenic effects
2434	Stiffness, neck	1341	Test, LAL (limulus amebocyte lysate)
1412	Stimulation of muscle	2216	Therapeutic abortion
1980	Stimulation of nerve, undesired	2099	Therapeutic effects, unexpected
1995	Stimulation, painful	2271	Therapeutic response, decreased
2084	Strangulation	2272	Therapeutic response, increased
1551	Stress from capsular contracture, rupture due to	2396	Throat, sore
2085	Stricture	2100	Thrombosis
2086	Stroke	2101	Thrombus
2336	Stroke syndrome	2102	Thyroid problems
1824	Stromal edema	2161	Tics, muscular
1893	Subarachnoid hemorrhage	2463	Tightness/pressure in chest
2247	Subclinical infection	2171	Tingling
1894	Subdural hemorrhage	2103	Tinnitus
2087	Subepithelial infiltration	2104	Tissue damage
2088	Suffocation	1981	Tissue damage, neural
2089	Sunset syndrome	1786	Tissue disease, connective
2177	Surgery, prolonged	2105	Tissue failure
2357	Surgical procedure	1417	Tissue nicks, cuts, or tears by device
2042	Surgical procedure, repeated	2106	Tissue, breakdown of optical
2267	Surgical shock	1972	Tissue, necrosis of flap
2444	Sweating	1297	Tissue, reaction of host
2091	Swelling	2060	Tissue, scar
2356	Swelling, joint	2110	Tissue/organ transplant
2092	Swollen glands	2070	Toe or finger, severed
2093	Swollen lymph nodes	2223	Tonic convulsion
1610	Syncope	2428	Tooth fracture
1684	Syndrome (AIDS), acquired immunodeficiency	2107	Torsades-de-Pointes
1696	Syndrome (ARDS), adult respiratory distress	1794	Touch, corneal
2108	Syndrome (TSS), toxic shock	2207	Toxemia
2046	Syndrome of newborns (RDS), respiratory distress	2108	Toxic shock syndrome, (TSS)
2239	Syndrome, first use	2333	Toxicity
2234	Syndrome, overwear	2098	Toxins in children
2256	Syndrome, radiation sickness	2269	Toxoplasmosis, acquired
2073	Syndrome, sjogren's	2270	Toxoplasmosis, congenital
2336	Syndrome, stroke	2393	Tracheostomy
2089	Syndrome, sunset	2120	Tract, infection of urinary
2114	Syndrome, twiddler's	2337	Transfusion of blood products
2115	Syndrome, uveitis-glaucoma-hyphema (UGH)	1749	Transfusion with incompatible blood
2094	Synovitis	2109	Transient ischemic attack
2036	Systemic reaction	2110	Transplant of organ/tissue
2095	Tachycardia	1795	Transplant, corneal
1731	Tachycardia, atrial	1694	Trauma, acoustic
2132	Tachycardia, ventricular	2268	Traumatic shock
2226	Tamponade, cardiac	2111	Tricuspid insufficiency

PATIENT-RELATED TERMS (continued)

2112	Tricuspid regurgitation	2055	Ventricular failure, right
2113	Tricuspid stenosis	2130	Ventricular fibrillation
2473	Tricuspid valve replacement	2131	Ventricular flutter
2108	TSS, toxic shock syndrome	1949	Ventricular hypertrophy, left
2114	Twiddler's syndrome	2056	Ventricular hypertrophy, right
2172	Twitching	2132	Ventricular tachycardia
2115	UGH (uveitis-glaucoma-hyphema) syndrome	2133	Ventriculomeglia
2274	Ulcer	2134	Vertigo
1796	Ulcer, corneal	1204	Vessel or plaque, device embedded in
2116	Ulceration	2135	Vessels, perforation of
2441	Ulceration, nipple	2025	Vessels, punctured blood
2326	Ulcers/sores, pressure	2248	Viral infection
1936	Uncontrolled/delayed intraocular pressure (IOP)	2136	Virus
2117	Under no medication	2197	Virus (HIV), human immunodeficiency
2118	Under no sedation	2137	Vision, blurring of
2166	Underdose, radiation	2138	Vision, impaired
1980	Undesired nerve stimulation	2139	Vision, loss of
2099	Unexpected therapeutic effects	2140	Visual disturbances
2369	Union fracture, delayed	2164	Visualization, phosphene
2202	Unknown (for use when the patient's condition is not known)	1866	Vitreous floaters
2490	Unresponsive	2143	Vitreous fluid, blood in
1859	Unretrieved fiberoptic fragments in body	2142	Vitreous fluid, loss of
2276	Unretrieved nonresorbable materials in body	2445	Vitreous, detachment of
2420	Upper respiratory tract infection	2181	Vitritis
2188	Uremia	2430	Volume (FEV), decreased forced expiratory
1871	Urgency	2431	Volume (FEV), increased forced expiratory
2275	Urinary frequency	2144	Vomiting
2119	Urinary retention	2145	Weakness
2120	Urinary tract infection	1967	Weakness, muscular
2278	Urticaria	2480	Wedge pressure, high pulmonary arterial
2121	Uterine perforation	2479	Wedge pressure, low pulmonary arterial
2122	Uveitis	2481	Wedge pressure, normal pulmonary arterial
2115	Uveitis-glaucoma-hyphema (UGH) syndrome	2241	Welt(s)
2123	Vaginal discharge	2241	Wheal(s)
2124	Vaginal mucosa damage	1154	Wound dehiscence
2376	Valve replacement, mitral	2446	Wound infection, post-operative
2471	Valve, prolapse of mitral	2447	Wound infection, post-traumatic
2476	Valve, repair of heart		
2472	Valve, replacement of aortic		
2474	Valve, replacement of pulmonary		
2473	Valve, replacement of tricuspid		
1926	Valvular insufficiency		
2335	Valvular regurgitation		
2002	Vascular disease, peripheral		
2126	Vasoconstriction		
2127	Vasodilatation		
2128	Vasospasm		
2071	VD, venereal disease		
2071	Venereal disease (VD)		
2129	Venipuncture		
2395	Ventilator dependent		
2078	Ventricle, abnormality of		
1947	Ventricular dysfunction, left		
2054	Ventricular dysfunction, right		
1948	Ventricular failure, left		

DEVICE-RELATED TERMS

2158	Aberration distortion, lens	1009	Arterial pressure alarm, failure of
1002	Abnormal	1036	Artifact
2424	Abnormal patch test results	1039	Artificial
2288	Aborted charge	1040	Aspherical lens
1387	Abrasion from instrument/object	1041	Aspiration, excessive
1003	Absorption	1042	Aspiration, incomplete
1082	Absorption, carbon	1050	Asymmetrical balloon
1001	AC/DC power, failure to run on	1441	Asynchronous pacing
1004	Accessories, incompatible	1010	Audible alarm
1591	Accidental spillage	1016	Audible alarm, low
1512	Accuracy rate	1019	Audible alarm, none
1675	Accuracy, volume	1043	Augmentation, loss of
1619	Activate system, failure to	1044	Autofill, unable to
1525	Activated, device remains	1046	Automatic injection system overinfusion
1486	Activation of system, premature	2307	Automatic injection system underinfusion
1557	Activation/keying, self	1734	Automatic injection system, failure to infuse
1006	Adaptor, failure of	1045	Automatic injection system, incompatible
1572	Advise shock, failure to	1440	Back-up mode, pacer found in
1007	Agglutinate, failure to	1047	Back-up, failure to
2300	Air eliminator, defective	1048	Back-up, failure to convert to
1008	Air leak(s)	2303	Bacterial contamination
1022	Alarm not visible	1050	Balloon asymmetrical
1620	Alarm system, failure of check-catheter	1051	Balloon burst
1621	Alarm system, failure of gas-leak/loss	1052	Balloon leak(s)
1615	Alarm system, failure of message-battery status	1053	Balloon mushroomed
1616	Alarm system, failure of message-check electrode	1049	Balloon rupture
1617	Alarm system, failure of message-leads off	1054	Battery charger, defective
1618	Alarm system, failure of message-service	1055	Battery failure
1010	Alarm, audible	1057	Battery, premature discharge of
1014	Alarm, defective	1355	Bed system, leak(s) from hydraulic
1011	Alarm, delayed	2310	Beds, rips/tears in specialty
1024	Alarm, error of warning	1059	Bend
1009	Alarm, failure of arterial pressure	1278	Bent haptic(s)
1020	Alarm, failure of fast heart rate	1060	Bicarbonate
2386	Alarm, failure of high inspiratory pressure	1061	Bifocal lens
1017	Alarm, failure of low flow	1062	Biofilm coating
2387	Alarm, failure of low inspiratory pressure	1063	Blank screen
1406	Alarm, failure of monitor display	1064	Bleed back
1025	Alarm, failure of warning	1065	Blockage
1012	Alarm, failure to	2317	Blood contaminated device
1013	Alarm, false	2290	Blood gas measurements, erroneous
1015	Alarm, intermittent	1746	Blood in tubing
1016	Alarm, low audible	2318	Blood pooling
1021	Alarm, no lead	1067	Bolus mechanism failure
1019	Alarm, not audible	1068	Bond, failure to
2286	Algorithms, inconsistent	1132	Break of counterbalance
1026	Altitude variations	1238	Break(s), filter
1027	Ambient temperature, changes in	1076	Break, cable
1193	Analysis, failure to perform EKG/ECG	1219	Break, external fiberoptic
1539	Analyze rhythm, failure to	1257	Break, footrest
1622	Angulation of table, unintended	1274	Break, handpiece
1029	Antenna, failure of	1277	Break, handrest
1031	Application problems	1445	Break, pallet
2289	Arcing at electrodes	1565	Break, shaft
1032	Arcing at paddles	1625	Break, table
1033	Arm motion, unintended	1229	Break/separation, fiberoptic
		1073	Breakage of burr

DEVICE-RELATED TERMS (continued)

1069	Breakage of device	1359	Cloudy lens
1347	Breakage of lead(s)	1095	Clumping
1638	Breakage of tip	1096	Coagulation
2347	Breakage of wire(s)	1062	Coating, biofilm
1088	Breaker, tripped circuit	2291	Code, failure to override semiautomatic
1103	Broken component(s)	1097	Coefficient error
1279	Broken haptic(s)	1098	Coiled
2297	Bubble detector, failure of	1099	Collapse
1070	Bubble(s)	1626	Collapse of table
2346	Buckling of tube(s)	1188	Collapse, dome
2305	Burn hole(s)	1100	Collimator design problem
1348	Burn(s) from lead(s)	1429	Collision, unintended
1071	Burned	1101	Colony forming units
1073	Burrs, breakage of	1330	Communication, interlumen
1074	Burst	1102	Compatibility
1051	Burst, balloon	2345	Compatibility/incompatibility of electro-magnetic interference
1206	Button or switch, failure of emergency stop	1105	Component(s) falling
1076	Cable break	1103	Component(s), broken
2302	Cable, defective	1104	Component(s), detachment of
1077	Calcification	1108	Component(s), incompatible
1189	Calculation error, dose	1106	Component(s), overheating of
1449	Calculations of parameter, incorrect	1107	Component(s), worn
1495	Calculations, incorrect programming of	2292	Components, defective
2440	Calibrate, failure to	2306	Components, missing
1078	Calibration error	1110	Computer failure
1079	Capacitive coupling	1111	Computer hardware error
1150	Capsulotomy, deflation due to	1112	Computer software error
1547	Capsulotomy, rupture due to	1113	Concentrate
1081	Capture, failure to	1114	Conduct, failure to
1080	Capture, intermittent	1115	Conductivity
1082	Carbon absorption	1116	Connection error
2319	Cardiac enzyme evaluation, erroneous	1371	Connection(s), loose
1670	Care/use of device, incorrect	1117	Connector pin failure
1234	Cartridge, unintended ejection of film	1118	Contact lens, problem with trial set
1920	Catheter, pull out of IV	1200	Contact(s), problem with electrode
1075	Cathode ray tube (CRT) failure	1120	Contamination
1148	Cause of deflation, unknown	2317	Contamination of device by blood
1548	Cause of rupture, unknown	1187	Contamination of distilled water
1083	Cautery	2303	Contamination, bacterial
1027	Changes in ambient temperature	1121	Continuity, intermittent
1084	Channeling	1122	Continuous
2288	Charge, aborted	1123	Continuous firing
1085	Charge, failure to	1124	Continuous mode failure
1054	Charger, defective battery	1128	Contraindicated coolant
1086	Charred	1126	Control settings incorrect
1620	Check-catheter alarm system, failure of	1540	Convert rhythm, failure to
1087	Chemical reaction	1048	Convert to back-up, failure to
1088	Circuit breaker tripped	1127	Cool, failure to
1089	Circuit failure	1128	Coolant, contraindicated
1231	Cladding material, separation of fiberoptic	1129	Cooling system, failure of
1224	Claim, false	1130	Cooling, inadequate
1090	Clean, failure to	1131	Corrode
1213	Cleaner, failure to remove enzymatic	1132	Counterbalance break
1091	Cleaning, inadequate	1571	Counters, resetting of shock
1092	Clearance	1133	Countershock, failure to deliver
1093	Cleaver, failure of fiberoptic	1134	Couple, failure to
1094	Clogged		

DEVICE-RELATED TERMS (continued)

1079	Coupling, capacitive	1508	Delivery of radiation to incorrect body area
1135	Crack(s)	2339	Delivery, inaccurate
1446	Crack(s) in pallet	1157	Deploy, difficult to
1639	Crack(s) in tomographic pallet	1158	Deploy, failure to
1549	Crease or fold due to rupture	1484	Deployment, premature
1136	Cross connection	1232	Deposits on lens
1137	Cross reactivity	1159	Deprimed
1075	CRT (cathode ray tube) failure	1222	Design failsafe, failure of
1138	Crush	1100	Design problem, collimator
1139	Cuff degeneration	1161	Design/structure problem
2454	Cut(s)	1280	Detached haptic(s)
1140	Cutoff point	1321	Detached insulation
1141	Cutter-torque-cable twisted or tangled	2043	Detachment from respirator
1142	Cycle, failure to	1586	Detachment from source
1405	Damage from moisture	1104	Detachment of component(s)
1570	Damage from shipping	1256	Detachment of footrest
1151	Damage from surgical instrument, deflation due to	1275	Detachment of handpiece
1550	Damage from surgical instrument, rupture due to	1276	Detachment of handrest
2284	Damage, internal/external	1302	Detachment of parts from IV pole
1143	Damaged screw tapper	1576	Detachment of shoulder rest
1216	Date of expiration exceeded	1577	Detachment of shoulder support
1804	Debris, metal shedding	2297	Detector, failure of bubble
1360	Decentration of lens	2160	Detector, failure of motion
1144	Declotting	1069	Device breakage
1145	Decoupling	2317	Device contaminated by blood
1500	Decrease of pump speed	1162	Device energization, unintended
1490	Decreased pressure	2379	Device failure
1146	Decreased suction	1628	Device material, tears, rips, holes in
2300	Defective air eliminator	1525	Device remains activated
1014	Defective alarm	1526	Device remains implanted
1054	Defective battery charger	1670	Device, incorrect care/use of
2302	Defective cable	1663	Device, inoperable
2292	Defective components	1328	Device, interference with monitoring
2294	Defective electrical wires	1527	Device, or device fragments remain in patient
1147	Defective photomultiplier (PM) tube	1628	Device, tears, rips holes in
2301	Defective seal	1260	Device/material fracture(s)
1612	Defective syringe	1541	Diagnose rhythm, failure to
2299	Defective tube(s)	1165	Dialysate
2298	Defective valve(s)	1166	Dialyzer, failure of
1149	Deflation difficulties	1149	Difficult to deflate
1150	Deflation due to capsulotomy	1157	Difficult to deploy
1151	Deflation due to damage from surgical instrument	1251	Difficult to flush
1152	Deflation due to trauma	1254	Difficult to fold
1148	Deflation, cause unknown	1310	Difficult to inflate
1139	Degeneration of cuff	1316	Difficult to insert
1201	Degeneration of electrode(s)	1331	Difficult to interrogate
1153	Degradation	1467	Difficult to position
1320	Degradation of insulation	1487	Difficult to prep
1155	Deionizer, failure of	1496	Difficult to program
1011	Delayed alarm	1528	Difficult to remove
1133	Deliver countershock, failure to	1668	Difficult to unwrap
1211	Deliver energy, failure to	1680	Difficult to wedge
1573	Deliver shock, failure to	1168	Disassembly
2338	Deliver, failure to	1352	Disc escape, leaflet
		1361	Disc lens
		1309	Disc, incomplete opening of
		2149	Discharge, electro-static

DEVICE-RELATED TERMS (continued)

1169	Discharge, failure to	1206	Emergency stop button or switch failure
1057	Discharge, premature battery	1207	Emergency table stop failure
1362	Discoloration of lens	2345	EMI, (electro-magnetic interference), compatibility/incompatibility
1170	Discolored	1194	EMI, electro-magnetic interference
1322	Discolored insulation	1480	End-of-life (EOL) indicator, premature
1171	Disconnect	1208	Endoscopic accessory fire or melt
1629	Discrepancies, telemetry	1162	Energization of device, unintended
1376	Discrepancy of magnet mode	1431	Energy output incorrect
1172	Disengaged	1209	Energy output to patient tissue incorrect
1175	Disinfect, failure to	1210	Energy spectrum incorrect
1176	Disinfection, inadequate/improper	1211	Energy, failure to deliver
1177	Disintegrate	1212	Entrapment
1178	Dislocated	1213	Enzymatic cleaner, failure to remove
1281	Dislocated haptic(s)	2319	Enzymes, erroneous evaluation of cardiac
1934	Dislocated intraocular lens (IOL)	1480	EOL (end-of-life) indicator, premature
1179	Dislodged	1483	ERI (elective replacement indicator), premature
1180	Displacement	1214	Erosion
1181	Display misread	1182	Erratic display
1182	Display, erratic	2290	Erroneous blood gas measurements
1183	Display, failure to	2319	Erroneous cardiac enzyme evaluation
1184	Display, incorrect	1345	Error of laser pulse timing
1304	Displayed image error	1024	Error of warning alarm
1185	Disposable	1215	Error or warning message, failure to produce
1607	Disruption of suture line	1078	Error, calibration
2412	Dissatisfaction	1097	Error, coefficient
1186	Dissection	1111	Error, computer hardware
1187	Distilled water, contaminated	1112	Error, computer software
2158	Distortion of lens aberration	1116	Error, connection
1188	Dome collapse	1189	Error, dose calculation
1189	Dose calculation error	1304	Error, image display
1191	Drift	1317	Error, installation
1611	Drive problem with syringe	2319	Evaluation of cardiac enzymes, erroneous
1410	Drive unit, motor stalled or jammed	1216	Exceeded expiration date
1192	Dry, failure to	1567	Exceeded shelf life
2407	Dull	1041	Excessive aspiration
1414	Dull needle	1464	Exchange of gas, poor
1234	Ejection of film cartridge, unintended	1216	Expiration date exceeded
1193	EKG/ECG analysis, failure to perform	1603	Expiratory or inspiratory phase, stuck in
1195	Elective removal	1485	Explantation, premature
1196	Elective replacement	1217	Explanted
1483	Elective replacement indicator (ERI), premature	1218	Explode
1198	Electrical failure	1651	Exploding tube(s)
1197	Electrical shock	1301	Explosion during hydrogen peroxide sterilization
2294	Electrical wires, defective	1509	Exposure to radiation, unintended
1194	Electro-magnetic interference (EMI)	1219	External fiberoptic break
2345	Electro-magnetic interference (EMI), compatibility/incompatibility	2284	External/internal damage
2149	Electro-static discharge	1220	Extra needle
1349	Electrocution from lead(s)	1221	Extraneous radiofrequency wave transmission
1200	Electrode contact(s), problem with	2154	Extrusion of implant
1201	Electrode(s), degeneration of	1222	Failsafe design failure
1199	Electrode(s), failure of	1006	Failure of adaptor
1202	Electrode(s), fracture of	1029	Failure of antenna
1203	Electrode(s), migration of	1009	Failure of arterial pressure alarm
2289	Electrodes, arcing	1734	Failure of automatic injection system to infuse
1205	Emergency power failure	1055	Failure of battery

DEVICE-RELATED TERMS (continued)

1067	Failure of bolus mechanism	1646	Failure of transducer
2297	Failure of bubble detector	1647	Failure of transmitter
1075	Failure of cathode ray tube (CRT)	1648	Failure of trocar
1620	Failure of check-catheter alarm system	1672	Failure of valve(s)
1089	Failure of circuit	1025	Failure of warning alarm
1117	Failure of connector pin	2315	Failure of water purification system
1129	Failure of cooling system	1677	Failure of water softener process
1155	Failure of deionizer	1572	Failure to advise shock
2379	Failure of device	1007	Failure to agglutinate
1166	Failure of dialyzer	1012	Failure to alarm
1199	Failure of electrode(s)	1539	Failure to analyze rhythm
1205	Failure of emergency power	1047	Failure to back-up
1206	Failure of emergency stop button or switch	1068	Failure to bond
1207	Failure of emergency table stop	2440	Failure to calibrate
1222	Failure of failsafe design	1081	Failure to capture
1020	Failure of fast heart rate alarm	1085	Failure to charge
1093	Failure of fiberoptic cleaver	1090	Failure to clean
2151	Failure of film processor	1114	Failure to conduct
1258	Failure of footswitch	1540	Failure to convert rhythm
1263	Failure of function indicator light(s)	1048	Failure to convert to back-up
1265	Failure of gas delivery system	1127	Failure to cool
1621	Failure of gas-leak/loss alarm system	1134	Failure to couple
1283	Failure of head immobilizer	1142	Failure to cycle
2386	Failure of high inspiratory pressure alarm	2338	Failure to deliver
1300	Failure of hydraulic system	1133	Failure to deliver countershock
2157	Failure of interlock(s)	1211	Failure to deliver energy
1367	Failure of lens washers	1573	Failure to deliver shock
1369	Failure of lockout mechanism	1158	Failure to deploy
1017	Failure of low flow alarm	1541	Failure to diagnose rhythm
2387	Failure of low inspiratory pressure alarm	1169	Failure to discharge
1375	Failure of magnet	1175	Failure to disinfect
1615	Failure of message-battery status alarm system	1183	Failure to display
1616	Failure of message-check electrode alarm system	1192	Failure to dry
1617	Failure of message-leads off alarm system	1252	Failure to flush
1618	Failure of message-service alarm system	1255	Failure to fold
2382	Failure of meter	1286	Failure to heat
1407	Failure of monitor	2340	Failure to infuse
1406	Failure of monitor display alarm	1332	Failure to interrogate
2160	Failure of motion detector	1379	Failure to maintain
1409	Failure of motor	1863	Failure to osseointegrate
1489	Failure of pressure sensor	1542	Failure to override rhythm
1493	Failure of probe	2291	Failure to override semiautomatic code
1513	Failure of rate modulated pacing sensor	1439	Failure to pace
1519	Failure of receiver	1193	Failure to perform EKG/ECG analysis
1520	Failure of receiver stimulator unit	1476	Failure to power-up
1568	Failure of shielding	1488	Failure to prep
1578	Failure of shunts	1492	Failure to prime
1579	Failure of shutter	1215	Failure to produce error or warning message
1580	Failure of side rails	1497	Failure to program
2324	Failure of solder joint	1502	Failure to pump
1590	Failure of speech processor	1581	Failure to read input signal
2304	Failure of stopcock valve	1517	Failure to recalibrate
2323	Failure of surgical graft	1213	Failure to remove enzymatic cleaner
1619	Failure of system to activate	1532	Failure to reset
1632	Failure of threader	1536	Failure to retract
1635	Failure of timer	1001	Failure to run on AC/DC power
		1466	Failure to run on portable mode

DEVICE-RELATED TERMS (continued)

1582	Failure to select signal	1634	Fluctuations of tidal volume
1563	Failure to service	1250	Fluid leak(s)
2170	Failure to suction	1251	Flush, difficult to
1521	Failure to transmit record	1252	Flush, failure to
1669	Failure to unwrap	1290	Flux rate, high
1683	Failure to zero	1253	Fogging
1110	Failure, computer	1549	Fold or crease due to rupture
1124	Failure, continuous mode	1254	Fold, difficult to
1198	Failure, electrical	1255	Fold, failure to
1299	Failure, hybrid	1257	Footrest break
1323	Failure, insulation	1256	Footrest detachment
1223	Failure, intermittent	1258	Footswitch failure
1384	Failure, mechanical	1259	Foreign material
1312	Failure, no infusion of injector system	1202	Fracture of electrode(s)
2311	Failure, out-of-box	1282	Fracture of haptic(s)
1313	Failure, overinfusion of injector system	1350	Fracture of lead(s)
1314	Failure, underinfusion of injector system	2413	Fracture of leaflet
1105	Falling component(s)	1260	Fracture(s) of device/material
1013	False alarm	1602	Fracture, strut
1224	False claim	1261	Fragmentation
1226	False output	1334	Fragmentation or fiberoptic break during intra-cavity procedure
1228	False readings	1527	Fragments of device remain in patient
1225	False-negative test result	1262	Frayed
1227	False-positive test result	1263	Function indicator light(s), failure of
1020	Fast heart rate alarm, failure of	2316	Fungus
1649	Fiber/trocar incompatibility	1264	Gain, inadequate
1230	Fiberguide incompatibility	1265	Gas delivery system failure
1334	Fiberoptic break or fragmentation during intra-cavity procedure	1464	Gas exchange, poor
1219	Fiberoptic break, external	1621	Gas leak/loss alarm system, failure of
1229	Fiberoptic break/separation	1266	Gas output, improper
1231	Fiberoptic cladding material separation	1267	Gel leakage
1093	Fiberoptic cleaver, failure of	1268	Germicide
1233	Filling problem	1269	Glass, shattered
1234	Film cartridge ejection, unintended	1270	Gradient increase
1235	Film processing, inadequate	2323	Graft, failure of surgical
2151	Film processor failure	1271	Grounding malfunctions
1236	Filter	1272	Growth, will not support
1238	Filter break(s)	1274	Handpiece break
1237	Filter leak(s)	1275	Handpiece detachment
1242	Filter tear(s)	1277	Handrest break
1243	Filter, assembly	1276	Handrest detachment
1244	Filter, inadequate	1278	Haptic(s), bent
1679	Filter, problem with wedge	1279	Haptic(s), broken
1556	Filter, problems with sediment	1280	Haptic(s), detached
2308	Filtration process, inadequate	1281	Haptic(s), dislocated
1245	Fire	1282	Haptic(s), fractured
1335	Fire or flash during intraprocedure	1111	Hardware error, computer
1208	Fire or melting of endoscopic accessories	1283	Head immobilizer failure
1123	Firing, continuous	1284	Head motion, unintended
2183	Fitting problems	1681	Head(s) on screw, incorrect
2383	Fixation, revision of internal	1285	Heat
1246	Flaking	1286	Heat, failure to
1247	Flammable	1287	Heating, inadequate
1335	Flash or fire during intraprocedure	1289	Hemofiltration
1249	Flowrate, inaccurate	1290	High flux rate
1248	Flowrate, restricted	1291	High impedance

DEVICE-RELATED TERMS (continued)

2386	High inspiratory pressure alarm, failure of	1287	Inadequate heating
1433	High output	1318	Inadequate instructions
2426	High pH	1442	Inadequate pacing
1470	High potentiometer readings	1533	Inadequate resistance
1491	High pressure	2167	Inadequate safety interlock(s)
2459	High readings	1564	Inadequate service
2457	High test results	1569	Inadequate shielding
1293	Hole(s)	1596	Inadequate sterilization
1325	Hole(s) in insulation	1600	Inadequate storage
1613	Hole(s) in syringe	1643	Inadequate training
2305	Hole(s), burn	1176	Inadequate/improper disinfection
1628	Holes, rips, tears in device, device material	2280	Inappropriate prompts
2385	Holes, rips, tears in packaging	1574	Inappropriate shock
1553	Homemade saline, use of	2283	Inaudible voice prompts
1294	Hose line occlusion	2345	Incompatibility/compatibility of electro-magnetic interference (EMI)
1295	Hose line rupture	1004	Incompatible accessories
1296	Hose mismatch	1045	Incompatible automatic injection system
1298	Hot oil leak	1108	Incompatible component(s)
1299	Hybrid failure	1230	Incompatible fiberguide
1355	Hydraulic bed system, leak(s) from	1630	Incompatible test strips
1300	Hydraulic system failure	1649	Incompatible trocar/fiber
1301	Hydrogen peroxide sterilization, explosion during	1650	Incompatible trocar/instrument
1303	Ignited	1042	Incomplete aspiration
1304	Image display error	1309	Incomplete opening of disc
1408	Image on monitor of poor quality	2312	Incomplete/missing packaging
1305	Image orientation incorrect	2286	Inconsistent algorithms
1306	Image resolution poor	1126	Incorrect control settings
1358	Image reversal, left or right	1508	Incorrect delivery of radiation to a body area
1291	Impedance, high	1184	Incorrect display
2285	Impedance, low	1431	Incorrect energy output
2154	Implant extrusion	1209	Incorrect energy output to patient tissue
2384	Implant, lens	1210	Incorrect energy spectrum
2320	Implant, removal of	1305	Incorrect image orientation
2321	Implant, repositioning of	1319	Incorrect instructions
2322	Implant, reprogramming of	1383	Incorrect measurements
1526	Implanted, device remains	1424	Incorrect occlusion
1307	Imprecision	1449	Incorrect parameter calculations
1266	Improper gas output	1654	Incorrect placement of tubing
2309	Improper rinsing	1494	Incorrect product
1176	Improper/inadequate disinfection	1495	Incorrect programming of calculations
1393	Improper/incorrect method	1535	Incorrect results
2017	Improper/incorrect procedure	1543	Incorrect rhythm interpretation
1337	Inability to irrigate	1681	Incorrect screw head(s)
2339	Inaccurate delivery	1555	Incorrect seal
1249	Inaccurate flowrate	1583	Incorrect size
1614	Inaccurate markings on syringe	1588	Incorrect source
1005	Inaccurate measurements	1623	Incorrect table incrementation
1609	Inaccurate synchronization	1670	Incorrect use/care of device
2456	Inaccurate test results	1676	Incorrect warning light
1662	Inactivation of unit	1682	Incorrect wavelength
1091	Inadequate cleaning	1393	Incorrect/improper method
1130	Inadequate cooling	2017	Incorrect/improper procedure
1235	Inadequate film processing	2405	Incorrectly placed syringe markings
1244	Inadequate filter	1374	Increase milliamps (MA), unable to
2308	Inadequate filtration process	1270	Increase of gradient
1264	Inadequate gain	1501	Increase of pump speed

DEVICE-RELATED TERMS (continued)

1491	Increased pressure	1331	Interrogate, difficult to
1604	Increased suction	1332	Interrogate, failure to
1623	Incrementation of table incorrect	1334	Intra-cavity procedure, fiberoptic break or fragmentation during
1263	Indicator of function light(s), failure of	1335	Intraprocedure, fire or flash during
1310	Inflation difficulties	1919	Intrauterine contraceptive device (IUD), removal of
1734	Infuse, failure of automatic injection system to	1336	Invagination
1312	Infuse, failure of injection system to	2293	Invalid sensing
2340	Infuse, failure to	1934	IOL (intraocular lens), dislocated
2341	Infusion, intermittent	1935	IOL (intraocular lens), migration of
1734	Injection system (automatic), failure to infuse	1337	Irrigate, inability to
1045	Injection system (automatic), incompatible	1919	IUD (intrauterine contraceptive device), removal of
1046	Injection system (automatic), overinfusion of	1920	IV catheter pull out
2307	Injection system (automatic), underinfusion of	1302	IV pole, detachment of parts from
1312	Injector system failure, no infusion	1410	Jammed or stalled motor drive unit (MDU)
1313	Injector system failure, overinfusion of	1557	Keying/activation, self
1314	Injector system failure, underinfusion of	1338	Kinetic
1663	Inoperable device	1339	Kink
1581	Input signal, failure to read	1340	Knot
1316	Insertion difficulties	2381	Laparoscopic sterilization
1603	Inspiratory or expiratory phase, stuck in	1342	Laser nozzle separation
2386	Inspiratory pressure alarm, failure of high	1343	Laser output above specifications
2387	Inspiratory pressure alarm, failure of low	1344	Laser output, unintended
1317	Installation error	1345	Laser pulse timing error
1318	Instructions, inadequate	1346	Laser tip separation
1319	Instructions, incorrect	1021	Lead alarm, none
1387	Instrument/object, abrasion from	1347	Lead(s), breakage of
1650	Instrument/trocar incompatibility	1348	Lead(s), burn(s) from
2391	Insufficient pressure	1349	Lead(s), electrocution from
1320	Insulation degradation	1350	Lead(s), fracture of
1323	Insulation failure	1351	Lead(s), shock from
1321	Insulation, detached	1352	Leaflet disc escape
1322	Insulation, discolored	2413	Leaflet fracture
1325	Insulation, hole(s) in	1298	Leak of hot oil
1326	Insulation, none	1354	Leak(s)
1327	Interference	1355	Leak(s) from hydraulic bed system
2345	Interference (EMI), compatibility, incompatibility of electro-magnetic	1587	Leak(s) from source
1194	Interference (EMI), electro-magnetic	1008	Leak(s), air
2314	Interference (RFI), radiofrequency	1052	Leak(s), balloon
1328	Interference with monitoring device	1237	Leak(s), filter
1329	Interference with pacemaker	1250	Leak(s), fluid
1378	Interference, magnetic	1390	Leak(s), membrane
2157	Interlock(s), failure of	1450	Leak(s), paravalvular
2167	Interlock(s), inadequate safety	1457	Leak(s), perivalvular
1330	Interlumen communication	1465	Leak(s), port
1015	Intermittent alarm	1357	Leak, radiation
1080	Intermittent capture	1267	Leakage of gel
1121	Intermittent continuity	1358	Left or right image reversal
1223	Intermittent failure	1934	Lens (IOL), dislocated intraocular
2341	Intermittent infusion	1935	Lens (IOL), migration of intraocular
1443	Intermittent pacing	2158	Lens aberration, distortion of
1558	Intermittent sensing	2384	Lens implant
2287	Intermittent shock	2184	Lens replacement
1599	Intermittently stops	1367	Lens washers, failure of
2383	Internal fixation, revision of	1040	Lens, aspherical
2284	Internal/external damage		
1543	Interpretation of rhythm, incorrect		

DEVICE-RELATED TERMS (continued)

1061	Lens, bifocal	2290	Measurements, erroneous blood gas
1359	Lens, cloudy	1005	Measurements, inaccurate
1360	Lens, decentration of	1383	Measurements, incorrect
1232	Lens, deposits on	1384	Mechanical failure
1361	Lens, disc	1067	Mechanism, failure of bolus
1362	Lens, discoloration of	1369	Mechanism, failure of lockout
1363	Lens, malposition of	1208	Melting or fire of endoscopic accessories
1411	Lens, multifocal	1385	Melts
1364	Lens, opacification of	1390	Membrane leak(s)
1118	Lens, problems with trial set of contact	1391	Membrane tear(s)
1365	Lens, repositioning of	1215	Message, failure to produce error or warning
1640	Lens, toric	1615	Message-battery status alarm system, failure of
1366	Lens, vaulting	1616	Message-check electrode alarm system, failure of
1263	Light(s), failure of function indicator	1617	Message-leads off alarm system, failure of
1676	Light, incorrect warning	1618	Message-service alarm system, failure of
1294	Line, occlusion of hose	1804	Metal shedding debris
1295	Line, rupture of hose	1392	Meter
1368	Linear	2382	Meter failure
1369	Lockout mechanism failure	1393	Method, improper/incorrect
1370	Looping	1395	Migration
1372	Loose	1203	Migration of electrode(s)
1371	Loose connection(s)	1935	Migration of intraocular lens (IOL)
2296	Loose temperature probe	1374	Milliamps (MA) rate, unable to increase
1043	Loss of augmentation	1396	Mineralization
2408	Loss of osseointegration	1397	Misapplication
1475	Loss of power	1398	Misasassembly
1534	Loss of resistance	1473	Miscalculation of power
1633	Loss of threshold	2410	Miscalibration
1671	Loss of vacuum	1399	Misconnection
1016	Low audible alarm	1401	Misfocusing
1017	Low flow alarm, failure of	1402	Mislabeled
2285	Low impedance	1403	Mislocation
2387	Low inspiratory pressure alarm, failure of	1296	Mismatched hose
1434	Low output	1404	Misplacement
2425	Low pH	1181	Misreading of display
1471	Low potentiometer readings	2306	Missing components
1490	Low pressure	2312	Missing/incomplete packaging
2460	Low readings	2295	Mistiming of water softener regeneration cycle
2458	Low test results	1124	Mode failure, continuous
1373	Lubrication	1466	Mode, failure to run on portable
1374	MA (milliamps), unable to increase	1440	Mode, pacer found in back-up
1375	Magnet failure	1405	Moisture damage
1376	Magnet mode discrepancy	1406	Monitor display alarm, failure of
1377	Magnet quench, unintended	1407	Monitor failure
1378	Magnetic interference	1408	Monitor, image quality poor
1379	Maintain, failure to	1328	Monitoring device, interference with
2409	Malfunction	2160	Motion detector failure
1271	Malfunction of grounding	1430	Motion of system, unintended
1637	Malfunction of timer	1624	Motion of table top, unintended
1363	Malposition of lens	1033	Motion, unintended arm
1614	Markings on syringe, inaccurate	1410	Motor drive unit (MDU) stalled or jammed
2405	Markings on syringe, incorrectly placed	1409	Motor failure
1381	Markings unclear	1411	Multifocal lens
1382	Markings, none	1053	Mushroomed balloon
1259	Material, foreign	1414	Needle, dull
1260	Material/device fracture(s)	1220	Needle, extra
1410	MDU (motor drive unit), stalled or jammed		

DEVICE-RELATED TERMS (continued)

1415	Needle, separation	1439	Pace, failure to
1416	Needle, unsheathed	1329	Pacemaker, interference with
1019	No audible alarm	1440	Pacer found in back-up mode
1326	No insulation	1441	Pacing asynchronously
1021	No lead alarm	1442	Pacing inadequately
1382	No markings	1443	Pacing intermittently
1435	No output	1513	Pacing sensor, failure of rate modulated
1559	No sensing	2312	Packaging, incomplete/missing
1022	No visible alarm	2385	Packaging, tears, rips, holes in
2282	No voice prompts	1444	Packaging, unsealed
2163	Noise	1032	Paddles, arcing at
1418	Non-linear	1445	Pallet break
1419	Nonpyrogenic	1446	Pallet crack(s)
1420	Nonstandard	1639	Pallet crack(s), tomographic
1421	Nonsterility	1448	Paramagnetic
1422	Normal	1449	Parameter calculations, incorrect
1342	Nozzle separation from laser	1450	Paravalvular leak(s)
2423	Obstruction	1451	Particulates
1516	Obtain readings, unable to	2424	Patch test, abnormal results of
1423	Occlusion	1527	Patient, device or device fragments remain in
1294	Occlusion of hose line	1454	Peeling
1424	Occlusion, incorrect	1455	Percutaneous
1425	Odor	2205	Perforation
1298	Oil leak, hot	1456	Performance
1364	Opacification of lens	1457	Perivalvular leak(s)
1426	Opaque	1301	Peroxide sterilization, explosion during hydrogen
1309	Opening of disc incomplete	2426	pH, high
1305	Orientation of image incorrect	2425	pH, low
1538	Osmosis, reversed	1458	Phantom problem
1863	Osseointegrate, failure to	1507	Phenomenon, R on T
2408	Osseointegration, loss of	1147	Photomultiplier (PM) tube, defective
2203	Other (for use when an appropriate device code cannot be identified)	1459	Pierce
2311	Out-of-box failure	1117	Pin connector, failure of
1432	Output above specifications	1460	Pitting
1431	Output energy incorrect	1654	Placement of tubing, incorrect
1343	Output from laser above specifications	1461	Plugged
1344	Output from laser, unintended	1462	Plunge
1209	Output of energy to patient tissue incorrect	1147	PM (photomultiplier) tube, defective
1266	Output of gas, improper	1463	Pocket stimulation
1226	Output, false	1302	Pole, detachment of parts from IV
1433	Output, high	2318	Pooling of blood
1434	Output, low	1464	Poor gas exchange
1435	Output, none	1306	Poor image resolution
2342	Overdelivery	1506	Poor or unsatisfactory quality
2313	Overexposure, radiation	1408	Poor quality of monitor image
2404	Overfill	1465	Port leak(s)
1437	Overheat	1466	Portable mode, failure to run on
1106	Overheating of component(s)	1467	Positioning difficulties
1644	Overheating of transducer	1469	Potency
1645	Overheating of transducer probe	1470	Potentiometer readings, high
1311	Overinfusion	1471	Potentiometer readings, low
1046	Overinfusion of automatic injection system	1474	Power surge
1313	Overinfusion of injector system	1205	Power, failure of emergency
1542	Override rhythm, failure to	1001	Power, failure to run on AC/DC
2291	Override semiautomatic code, failure to	1475	Power, loss of
1438	Oversensing	1473	Power, miscalculation of

DEVICE-RELATED TERMS (continued)

1476	Power-up, failure to	1502	Pump, failure to
1477	Pre/post-pumping problems	1477	Pumping problems, pre/post
1478	Precipitate	1503	Pumping, stopped
1479	Precision	1504	Puncture
1484	Premature deployment	1505	Pyrogenic
1057	Premature discharge of battery	1408	Quality of monitor image, poor
1480	Premature end-of-life (EOL) indicator	1506	Quality, unsatisfactory or poor
1483	Premature ERI (elective replacement indicator)	1377	Quench of magnet, unintended
1485	Premature explantation	1507	R on T phenomenon
1486	Premature system activation	1508	Radiation delivered to incorrect body area
1487	Prep, difficult to	1509	Radiation exposure, unintended
1488	Prep, failure to	1357	Radiation leak
1009	Pressure alarm, failure of arterial	2313	Radiation overexposure
2386	Pressure alarm, failure of high inspiratory	2031	Radiation underexposure
2387	Pressure alarm, failure of low inspiratory	1511	Radio signal problems
1489	Pressure sensor failure	2314	Radiofrequency interference (RFI)
1490	Pressure, decreased	1221	Radiofrequency wave transmission, extraneous
1491	Pressure, high	1580	Rails, failure of side
1491	Pressure, increased	1513	Rate modulated pacing sensor, failure of
2391	Pressure, insufficient	1512	Rate, accuracy
1490	Pressure, low	1514	Reaction
1492	Prime, failure to	1087	Reaction, chemical
1493	Probe failure	1137	Reactivity, cross
1645	Probe overheating, transducer	1581	Read input signal, failure to
2296	Probe, loose temperature	1470	Readings on potentiometer, high
1100	Problem with collimator design	1471	Readings on potentiometer, low
1118	Problem with contact lens trial set	1228	Readings, false
1161	Problem with design/structure	2459	Readings, high
1200	Problem with electrode contact(s)	2460	Readings, low
1511	Problem with radio signals	1516	Readings, unable to obtain
1561	Problem with sensor(s)	1517	Recalibrate, failure to
1570	Problem with shipping	1518	Recannulation
1592	Problem with spring loading mechanism	1519	Receiver failure
1594	Problem with steering wire	1520	Receiver stimulator unit, failure of
1611	Problem with syringe drive	1521	Record, failure to transmit
1679	Problem with wedge filter	1522	Reflux
1233	Problem, filling	1523	Regenerate
1458	Problem, phantom	2295	Regeneration cycle, mistiming of water softener
1477	Problems with pre/post pumping	1524	Rejection
1556	Problems with sediment filter	1527	Remains in patient, device or device fragments
1031	Problems, application	1528	Removal difficulties
2183	Problems, fitting	2320	Removal of implant
2017	Procedure, improper/incorrect	1919	Removal of intrauterine contraceptive device (IUD)
1677	Process, failure of water softener	1195	Removal, elective
1235	Processing of film, inadequate	1213	Remove enzymatic cleaner, failure to
2151	Processor, failure of film	1529	Repair
1494	Product, incorrect	1530	Replace
1496	Program, difficult to	2184	Replacement of lens
1497	Program, failure to	1196	Replacement, elective
1495	Programming calculations, incorrect	2321	Repositioning of implant
2281	Prompts will not clear	1365	Repositioning of lens
2280	Prompts, inappropriate	2322	Reprogramming of implant
2283	Prompts, inaudible voice	1531	Reservoir
2282	Prompts, no voice	1532	Reset, failure to
1920	Pull out of IV catheter	1571	Resetting of shock counters
1500	Pump speed, decreased		
1501	Pump speed, increased		

DEVICE-RELATED TERMS (continued)

2325	Residue	1561	Sensor problems
1533	Resistance, inadequate	1489	Sensor, failure of pressure
1534	Resistance, loss of	1513	Sensor, failure of rate modulated pacing
1306	Resolution, poor image	1562	Separates
2043	Respirator, detachment from	1342	Separation from laser nozzle
1248	Restricted flowrate	1346	Separation from laser tip
1225	Result, false-negative test	1231	Separation of fiberoptic cladding material
1227	Result, false-positive test	1415	Separation of needle
1535	Results incorrect	1608	Separation of suture line
2424	Results of patch test, abnormal	1229	Separation/break, fiberoptic
2457	Results of tests, high	1563	Service, failure to
2456	Results of tests, inaccurate	1564	Service, inadequate
2458	Results of tests, low	1118	Set, problems with trial contact lens
1631	Results, unexpected therapeutic	1126	Setting of controls, incorrect
1536	Retract, failure to	1565	Shaft break
1537	Reuse	1566	Shaft, split
1358	Reversal of image, left or right	1269	Shattered glass
1538	Reverse osmosis	1652	Shattered tube(s)
2383	Revision of internal fixation	1804	Shedding of metal debris
2314	RFI, radiofrequency interference	1567	Shelf life exceeded
1539	Rhythm, failure to analyze	1568	Shielding failure
1540	Rhythm, failure to convert	1569	Shielding, inadequate
1541	Rhythm, failure to diagnose	1570	Shipping damage or problem
1542	Rhythm, failure to override	1571	Shock counters, resetting of
1543	Rhythm, incorrect interpretation of	1351	Shock from lead(s)
1358	Right or left image reversal	1197	Shock, electrical
1544	Rigid	1572	Shock, failure to advise
2309	Rinsing, improper	1573	Shock, failure to deliver
1628	Rips, tears, holes in device, device material	1574	Shock, inappropriate
2385	Rips, tears, holes in packaging	2287	Shock, intermittent
2310	Rips/tears in specialty beds	1575	Short fill
1545	Runaway	1576	Shoulder rest, detachment of
1546	Rupture	1577	Shoulder support, detachment of
1547	Rupture due to capsulotomy	1578	Shunts, failure of
1550	Rupture due to damage from surgical instrument	1579	Shutter failure
1552	Rupture due to trauma	1580	Side rails, failure of
1049	Rupture of balloon	1511	Signal problems, radio
1295	Rupture of hose line	1581	Signal, failure to read input
1548	Rupture, cause unknown	1582	Signal, failure to select
1549	Rupture, fold or crease due to	1583	Size, incorrect
2167	Safety interlock(s) inadequate	1584	Slippage
1553	Saline, use of homemade	1585	Smoke
1554	Salt tablet(s), use of	1112	Software error, computer
1063	Screen, blank	2324	Solder joint failure
1681	Screw head(s) incorrect	1586	Source, detachment from
1143	Screw tapper, damaged	1588	Source, incorrect
2301	Seal, defective	1587	Source, leak(s) from
1555	Seal, incorrect	1665	Source, unretracted
1556	Sediment filter problems	2310	Specialty beds, rips/tears in
1582	Select signal, failure to	1343	Specifications, laser output above
1557	Self-activation/keying	1432	Specifications, output above
2291	Semiautomatic code, failure to override	1589	Specificity
1558	Sensing intermittently	1210	Spectrum, incorrect energy
2293	Sensing, invalid	1590	Speech processor failure
1559	Sensing, none	1500	Speed of pump, decreased
1560	Sensitivity	1501	Speed of pump, increased
		1591	Spillage, accidental

DEVICE-RELATED TERMS (continued)

1566	Split, shaft	1486	System, premature activation of
1653	Splitting of tube(s)	2307	System, underinfusion of automatic injection
1592	Spring loading mechanism problem	1622	Table angulation, unintended
1410	Stalled or jammed motor drive unit (MDU)	1625	Table break
2149	Static discharge, electro-	1626	Table collapse
1594	Steering wire problem	1623	Table incrementation incorrect
1595	Sterility	1207	Table stop, failure of emergency
2381	Sterilization of laparoscope	1624	Table top motion, unintended
1301	Sterilization using hydrogen peroxide, explosion during	1554	Tablet(s), use of salt
1596	Sterilization, inadequate	1141	Tangled or twisted cutter-torque-cable
1597	Sticking	1627	Tap water, use of
1673	Sticking valve(s)	1242	Tear(s), filter
1463	Stimulation, pocket	1391	Tear(s), membrane
1520	Stimulator receiver unit, failure of	1628	Tears, rips, holes in device, device material
1206	Stop button or switch, failure of emergency	2385	Tears, rips, holes in packaging
2304	Stopcock valve, failure of	2310	Tears/rips in specialty beds
1503	Stopped pumping	1629	Telemetry discrepancies
1599	Stops intermittently	2296	Temperature probe, loose
1600	Storage, inadequate	1027	Temperature, changes in ambient
1601	Stretched	1225	Test result, false-negative
1630	Strips used for testing, incompatible	1227	Test result, false-positive
1161	Structure/design problem	2424	Test results, abnormal patch
1602	Strut fracture	2457	Test results, high
1603	Stuck in inspiratory or expiratory phase	2456	Test results, inaccurate
2170	Suction failure	2458	Test results, low
1146	Suction, decreased	1630	Test strips, incompatible
1604	Suction, increased	1631	Therapeutic results, unexpected
1605	Superheat	1632	Threader failure
1606	Supplier	1633	Threshold, loss of
1474	Surge of power	1634	Tidal volume fluctuations
2323	Surgical graft, failure of	1635	Timer failure
1151	Surgical instrument, deflation due to damage from	1637	Timer malfunction
1550	Surgical instrument, rupture damage due to	1345	Timing error, laser pulse
1607	Suture line disruption	1638	Tip breakage
1608	Suture line separation	1346	Tip of laser, separation from
1206	Switch or button stop, failure of emergency	1209	Tissue, incorrect energy output to patient's
1609	Synchronization, inaccurate	1639	Tomographic pallet crack(s)
1611	Syringe drive problem	1640	Toric lens
1614	Syringe markings inaccurate	1641	Torqued
2405	Syringe markings, incorrectly placed	1642	Tracking
1612	Syringe, defective	1643	Training, inadequate
1613	Syringe, hole(s) in	1646	Transducer failure
1619	System fails to activate	1644	Transducer overheating
1620	System failure of check-catheter alarm	1645	Transducer probe overheating
1621	System failure of gas-leak/loss alarm	1221	Transmission, extraneous radiofrequency wave
1615	System failure of message-battery status alarm	1521	Transmit record, failure to
1616	System failure of message-check electrode alarm	1647	Transmitter failure
1617	System failure of message-leads off alarm	1152	Trauma, deflation due to
1618	System failure of message-service alarm	1552	Trauma, rupture due to
1430	System motion, unintended	1678	Treatment, water
1129	System, failure of cooling	1118	Trial contact lens set, problems with
1265	System, failure of gas delivery	1088	Tripped circuit breaker
1300	System, failure of hydraulic	1648	Trocar failure
2315	System, failure of water purification	1649	Trocar/fiber incompatibility
		1650	Trocar/instrument incompatibility
		2346	Tube(s), buckling of
		2299	Tube(s), defective

DEVICE-RELATED TERMS (continued)

1651	Tube(s), exploding of	1366	Vaulting of lens
1652	Tube(s), shattering of	1674	Vibration
1653	Tube(s), splitting of	1022	Visible alarm, none
1746	Tubing, blood in	2283	Voice prompts, inaudible
1654	Tubing, incorrect placement of	2282	Voice prompts, none
1141	Twisted or tangled cutter-torque-cable	1675	Volume accuracy
1655	Twisting	1634	Volume fluctuations, tidal
1656	Ultrafiltration	1024	Warning alarm, error of
1657	Ultraviolet	1025	Warning alarm, failure of
1658	Ultraviolet absorbing	1676	Warning light, incorrect
1044	Unable to autofill	1215	Warning or error message, failure to produce
1374	Unable to increase milliamps (MA)	1367	Washers, failure of lens
1516	Unable to obtain readings	2315	Water purification system, failure of
1381	Unclear markings	1677	Water softener process, failure of
1659	Uncoiled	2295	Water softener regeneration cycle, mistiming of
1660	Undercorrection	1678	Water treatment
2343	Underdelivery	1187	Water, contamination of distilled
2031	Underexposure, radiation	1682	Wavelength, incorrect
2182	Underinfusion	1679	Wedge filter problem
2307	Underinfusion of automatic injection system	1680	Wedge, difficult to
1314	Underinfusion of injector system	1272	Will not support growth
1661	Undersensing	2347	Wire(s), breakage of
1631	Unexpected therapeutic results	2294	Wires, defective electrical
1033	Unintended arm motion	1107	Worn component(s)
1429	Unintended collision	1683	Zero, failure to
1162	Unintended device energization		
1234	Unintended film cartridge ejection		
1284	Unintended head motion		
1344	Unintended laser output		
1377	Unintended magnet quench		
1509	Unintended radiation exposure		
1430	Unintended system motion		
1622	Unintended table angulation		
1624	Unintended table top motion		
1662	Unit inactivated		
1520	Unit, failure of receiver stimulator		
1101	Units, colony forming		
2204	Unknown (for use when the device problem is not known)		
1148	Unknown cause of deflation		
1548	Unknown cause of rupture		
1664	Unravel		
1665	Unretracted source		
1506	Unsatisfactory or poor quality		
1444	Unsealed packaging		
1416	Unsheathed needle		
1667	Unstable		
1668	Unwrap, difficult to		
1669	Unwrap, failure to		
1627	Use of tap water		
1670	Use/care of device, incorrect		
1671	Vacuum, loss of		
2298	Valve(s), defective		
1672	Valve(s), failure of		
1673	Valve(s), sticking		
2304	Valve, failure of stopcock		
1026	Variations in altitude		

PART I, SUBPART B - NUMERIC LISTING OF EVENT CODES WITH THEIR CORRESPONDING EVENT TERM

CODE	TERM	CODE	TERM
1001	AC/DC, failure to run on	1062	Biofilm coating
1002	Abnormal	1063	Blank screen
1003	Absorption	1064	Bleed back
1004	Accessories, incompatible	1065	Blockage
1005	Measurements, inaccurate	1067	Bolus mechanism failure
1006	Adaptor, failure of	1068	Bond, failure to
1007	Agglutinate, failure to	1069	Device breakage
1008	Air leak(s)	1070	Bubble(s)
1009	Alarm, failure of arterial pressure	1071	Burned
1010	Alarm, audible	1073	Burrs, breakage of
1011	Alarm, delayed	1074	Burst
1012	Alarm, failure to	1075	Cathode ray tube (CRT) failure
1013	Alarm, false	1076	Cable break
1014	Alarm, defective	1077	Calcification
1015	Alarm, intermittent	1078	Calibration error
1016	Alarm, low audible	1079	Capacitive coupling
1017	Alarm, failure of low flow	1080	Capture, intermittent
1019	Alarm, not audible	1081	Capture, failure to
1020	Alarm, failure of fast heart rate	1082	Carbon absorption
1021	Alarm, no lead	1083	Cautery
1022	Alarm not visible	1084	Channeling
1024	Alarm, error of warning	1085	Charge, failure to
1025	Alarm, failure of warning	1086	Charred
1026	Altitude variations	1087	Chemical reaction
1027	Ambient temperature, changes in	1088	Circuit breaker tripped
1028	Anastomose, failure to	1089	Circuit failure
1029	Antenna, failure of	1090	Clean, failure to
1031	Application problems	1091	Cleaning, inadequate
1032	Arcing at paddles	1092	Clearance
1033	Arm motion, unintended	1093	Fiberoptic cleaver, failure of
1034	Blood gas measurements, abnormal	1094	Clogged
1036	Artifact	1095	Clumping
1039	Artificial	1096	Coagulation
1040	Lens, aspherical	1097	Coefficient error
1041	Aspiration, excessive	1098	Coiled
1042	Aspiration, incomplete	1099	Collapse
1043	Augmentation, loss of	1100	Collimator design problem
1044	Autofill, unable to	1101	Colony forming units
1045	Automatic injection system, incompatible	1102	Compatibility
1046	Automatic injection system overinfusion	1103	Component(s), broken
1047	Back-up, failure to	1104	Component(s), detachment of
1048	Back-up, failure to convert to	1105	Component(s) falling
1049	Balloon rupture	1106	Component(s), overheating of
1050	Balloon asymmetrical	1107	Component(s), worn
1051	Balloon burst	1108	Component(s), incompatible
1052	Balloon leak(s)	1110	Computer failure
1053	Balloon mushroomed	1111	Computer hardware error
1054	Battery charger, defective	1112	Computer software error
1055	Battery failure	1113	Concentrate
1057	Battery, premature discharge of	1114	Conduct, failure to
1059	Bend	1115	Conductivity
1060	Bicarbonate	1116	Connection error
1061	Lens, bifocal	1117	Connector pin failure

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1118	Contact lens, problem with trial set	1178	Dislocated
1120	Contamination	1179	Dislodged
1121	Continuity, intermittent	1180	Displacement
1122	Continuous	1181	Display misread
1123	Continuous firing	1182	Display, erratic
1124	Continuous mode failure	1183	Display, failure to
1125	Contraindicated patient	1184	Display, incorrect
1126	Control settings incorrect	1185	Disposable
1127	Cool, failure to	1186	Dissection
1128	Coolant, contraindicated	1187	Distilled water, contaminated
1129	Cooling system, failure of	1188	Dome collapse
1130	Cooling, inadequate	1189	Dose calculation error
1131	Corrode	1191	Drift
1132	Counterbalance break	1192	Dry, failure to
1133	Countershock, failure to deliver	1193	EKG/ECG analysis, failure to perform
1134	Couple, failure to	1194	Electro-magnetic interference (EMI)
1135	Crack(s)	1195	Elective removal
1136	Cross connection	1196	Elective replacement
1137	Cross reactivity	1197	Shock, electrical
1138	Crush	1198	Electrical failure
1139	Cuff degeneration	1199	Electrode(s), failure of
1140	Cutoff point	1200	Electrode contact(s), problem with
1141	Cutter-torque-cable twisted or tangled	1201	Electrode(s), degeneration of
1142	Cycle, failure to	1202	Electrode(s), fracture of
1143	Screw tapper, damaged	1203	Electrode(s), migration of
1144	Declotting	1204	Vessel or plaque, device embedded in
1145	Decoupling	1205	Emergency power failure
1146	Suction, decreased	1206	Emergency stop button or switch failure
1147	Photomultiplier (PM) tube, defective	1207	Emergency table stop failure
1148	Deflation, cause unknown	1208	Endoscopic accessory fire or melt
1149	Deflation difficulties	1209	Energy output to patient tissue incorrect
1150	Deflation due to capsulotomy	1210	Energy spectrum incorrect
1151	Deflation due to damage from surgical instrument	1211	Energy, failure to deliver
1152	Deflation due to trauma	1212	Entrapment
1153	Degradation	1213	Enzymatic cleaner, failure to remove
1154	Wound dehiscence	1214	Erosion
1155	Deionizer, failure of	1215	Error or warning message, failure to produce
1156	Delamination	1216	Expiration date exceeded
1157	Deploy, difficult to	1217	Explanted
1158	Deploy, failure to	1218	Explode
1159	Deprimed	1219	External fiberoptic break
1161	Design/structure problem	1220	Needle, extra
1162	Device energization, unintended	1221	Extraneous radiofrequency wave transmission
1165	Dialysate	1222	Failsafe design failure
1166	Dialyzer, failure of	1223	Failure, intermittent
1168	Disassembly	1224	False claim
1169	Discharge, failure to	1225	False-negative test result
1170	Discolored	1226	False output
1171	Disconnect	1227	False-positive test result
1172	Disengaged	1228	False readings
1175	Disinfect, failure to	1229	Fiberoptic break/separation
1176	Disinfection, inadequate/improper	1230	Fiberguide incompatibility
1177	Disintegrate	1231	Fiberoptic cladding material separation

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1232	Lens, deposits on	1290	High flux rate
1233	Filling problem	1291	Impedance, high
1234	Film cartridge ejection, unintended	1293	Hole(s)
1235	Film processing, inadequate	1294	Hose line occlusion
1236	Filter	1295	Hose line rupture
1237	Filter leak(s)	1296	Hose mismatch
1238	Filter break(s)	1297	Host-tissue reaction
1242	Filter tear(s)	1298	Hot oil leak
1243	Filter, assembly	1299	Hybrid failure
1244	Filter, inadequate	1300	Hydraulic system failure
1245	Fire	1301	Hydrogen peroxide sterilization, explosion during
1246	Flaking	1302	IV pole, detachment of parts from
1247	Flammable	1303	Ignited
1248	Flowrate, restricted	1304	Image display error
1249	Flowrate, inaccurate	1305	Image orientation incorrect
1250	Fluid leak(s)	1306	Image resolution poor
1251	Flush, difficult to	1307	Imprecision
1252	Flush, failure to	1309	Disc, incomplete opening of
1253	Fogging	1310	Inflation difficulties
1254	Fold, difficult to	1311	Overinfusion
1255	Fold, failure to	1312	Injector system failure, no infusion
1256	Footrest detachment	1313	Injector system failure, overinfusion of
1257	Footrest break	1314	Injector system failure, underinfusion of
1258	Footswitch failure	1316	Insertion difficulties
1259	Foreign material	1317	Installation error
1260	Fracture(s) of device/material	1318	Instructions, inadequate
1261	Fragmentation	1319	Instructions, incorrect
1262	Frayed	1320	Insulation degradation
1263	Function indicator light(s), failure of	1321	Insulation, detached
1264	Gain, inadequate	1322	Insulation, discolored
1265	Gas delivery system failure	1323	Insulation failure
1266	Gas output, improper	1325	Insulation, hole(s) in
1267	Gel leakage	1326	Insulation, none
1268	Germicide	1327	Interference
1269	Glass, shattered	1328	Interference with monitoring device
1270	Gradient increase	1329	Interference with pacemaker
1271	Grounding malfunctions	1330	Interlumen communication
1272	Growth, will not support	1331	Interrogate, difficult to
1274	Handpiece break	1332	Interrogate, failure to
1275	Handpiece detachment	1333	Intimal dissection
1276	Handrest detachment	1334	Intra-cavity procedure, fiberoptic break or fragmentation during
1277	Handrest break	1335	Intraprocedure, fire or flash during
1278	Haptic(s), bent	1336	Invagination
1279	Haptic(s), broken	1337	Irrigate, inability to
1280	Haptic(s), detached	1338	Kinetic
1281	Haptic(s), dislocated	1339	Kink
1282	Haptic(s), fractured	1340	Knot
1283	Head immobilizer failure	1341	LAL test (limulus ameobocyte lysate)
1284	Head motion, unintended	1342	Laser nozzle separation
1285	Heat	1343	Laser output above specifications
1286	Heat, failure to	1344	Laser output, unintended
1287	Heating, inadequate	1345	Laser pulse timing error
1288	Hemoconcentration		
1289	Hemofiltration		

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1346	Laser tip separation	1406	Monitor display alarm, failure of
1347	Lead(s), breakage of	1407	Monitor failure
1348	Lead(s), burn(s) from	1408	Monitor, image quality poor
1349	Lead(s), electrocution from	1409	Motor failure
1350	Lead(s), fracture of	1410	Motor drive unit (MDU) stalled or jammed
1351	Lead(s), shock from	1411	Lens, multifocal
1352	Leaflet disc escape	1412	Muscle stimulation
1353	Leaflet disruption	1413	Necrosis, pressure
1354	Leak(s)	1414	Needle, dull
1355	Leak(s) from hydraulic bed system	1415	Needle, separation
1357	Leak, radiation	1416	Needle, unsheathed
1358	Image reversal, left or right	1417	Nicks, cuts or tears of dura or other tissues by device
1359	Lens, cloudy	1418	Non-linear
1360	Lens, decentration of	1419	Nonpyrogenic
1361	Lens, disc	1420	Nonstandard
1362	Lens, discoloration of	1421	Nonsterility
1363	Lens, malposition of	1422	Normal
1364	Lens, opacification of	1423	Occlusion
1365	Lens, repositioning of	1424	Occlusion, incorrect
1366	Lens, vaulting	1425	Odor
1367	Lens washers, failure of	1426	Opaque
1368	Linear	1428	Ossification
1369	Lockout mechanism failure	1429	Collision, unintended
1370	Looping	1430	System motion, unintended
1371	Connection(s), loose	1431	Output energy incorrect
1372	Loose	1432	Output above specifications
1373	Lubrication	1433	Output, high
1374	Milliamps (MA) rate, unable to increase	1434	Output, low
1375	Magnet failure	1435	Output, none
1376	Magnet mode discrepancy	1436	Overcorrection
1377	Magnet quench, unintended	1437	Overheat
1378	Magnetic interference	1438	Oversensing
1379	Maintain, failure to	1439	Pace, failure to
1381	Markings unclear	1440	Pacer found in back-up mode
1382	Markings, none	1441	Pacing asynchronously
1383	Measurements, incorrect	1442	Pacing inadequately
1384	Mechanical failure	1443	Pacing intermittently
1385	Melts	1444	Packaging, unsealed
1387	Abrasion from instrument/object	1445	Pallet break
1388	Membrane, breakage of	1446	Pallet crack(s)
1390	Membrane leak(s)	1447	Pannus formation
1391	Membrane tear(s)	1448	Paramagnetic
1392	Meter	1449	Parameter calculations, incorrect
1393	Method, improper/incorrect	1450	Paravalvular leak(s)
1395	Migration	1451	Particulates
1396	Mineralization	1454	Peeling
1397	Misapplication	1455	Percutaneous
1398	Misassembly	1456	Performance
1399	Misconnection	1457	Leak(s), perivalvular
1401	Misfocusing	1458	Phantom problem
1402	Mislabeled	1459	Pierce
1403	Mislocation	1460	Pitting
1404	Misplacement	1461	Plugged
1405	Moisture damage		

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1462	Plunge	1521	Record, failure to transmit
1463	Pocket stimulation	1522	Reflux
1464	Poor gas exchange	1523	Regenerate
1465	Port leak(s)	1524	Rejection
1466	Portable mode, failure to run on	1525	Device remains activated
1467	Positioning difficulties	1526	Device remains implanted
1469	Potency	1527	Device, or device fragments remain in patient
1470	Potentiometer readings, high	1528	Removal difficulties
1471	Potentiometer readings, low	1529	Repair
1473	Power, miscalculation of	1530	Replace
1474	Power surge	1531	Reservoir
1475	Power, loss of	1532	Reset, failure to
1476	Power-up, failure to	1533	Resistance, inadequate
1477	Pre/post-pumping problems	1534	Resistance, loss of
1478	Precipitate	1535	Results incorrect
1479	Precision	1536	Retract, failure to
1480	Premature end-of-life (EOL) indicator	1537	Reuse
1483	Premature ERI (elective replacement indicator)	1538	Reverse osmosis
1484	Premature deployment	1539	Rhythm, failure to analyze
1485	Premature explantation	1540	Rhythm, failure to convert
1486	Premature system activation	1541	Rhythm, failure to diagnose
1487	Prep, difficult to	1542	Rhythm, failure to override
1488	Prep, failure to	1543	Rhythm, incorrect interpretation of
1489	Pressure sensor failure	1544	Rigid
1490	Pressure, decreased	1545	Runaway
1491	Pressure, increased	1546	Rupture
1492	Prime, failure to	1547	Rupture due to capsulotomy
1493	Probe failure	1548	Rupture, cause unknown
1494	Product, incorrect	1549	Rupture, fold or crease due to
1495	Programming calculations, incorrect	1550	Rupture due to damage from surgical instrument
1496	Program, difficult to	1551	Rupture due to stress from capsular contracture
1497	Program, failure to	1552	Rupture due to trauma
1498	Pulmonary embolism	1553	Saline, use of homemade
1500	Pump speed, decreased	1554	Salt tablet(s), use of
1501	Pump speed, increased	1555	Seal, incorrect
1502	Pump, failure to	1556	Sediment filter problems
1503	Pumping, stopped	1557	Self-activation/keying
1504	Puncture	1558	Sensing intermittently
1505	Pyrogenic	1559	Sensing, none
1506	Quality, unsatisfactory or poor	1560	Sensitivity
1507	R on T phenomenon	1561	Sensor problems
1508	Radiation delivered to incorrect body area	1562	Separates
1509	Radiation exposure, unintended	1563	Service, failure to
1510	Radiation overdose	1564	Service, inadequate
1511	Radio signal problems	1565	Shaft break
1512	Accuracy rate	1566	Shaft, split
1513	Rate modulated pacing sensor, failure of	1567	Shelf life exceeded
1514	Reaction	1568	Shielding failure
1516	Readings, unable to obtain	1569	Shielding, inadequate
1517	Recalibrate, failure to	1570	Shipping damage or problem
1518	Recannulation	1571	Shock counters, resetting of
1519	Receiver failure		
1520	Receiver stimulator unit, failure of		

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1572	Shock, failure to advise	1625	Table break
1573	Shock, failure to deliver	1626	Table collapse
1574	Shock, inappropriate	1627	Tap water, use of
1575	Short fill	1628	Tears, rips, holes in device, device material
1576	Shoulder rest, detachment of	1629	Telemetry discrepancies
1577	Shoulder support, detachment of	1630	Test strips, incompatible
1578	Shunts, failure of	1631	Therapeutic results, unexpected
1579	Shutter failure	1632	Threader failure
1580	Side rails, failure of	1633	Threshold, loss of
1581	Input signal, failure to read	1634	Tidal volume fluctuations
1582	Signal, failure to select	1635	Timer failure
1583	Size, incorrect	1637	Timer malfunction
1584	Slippage	1638	Tip breakage
1585	Smoke	1639	Tomographic pallet crack(s)
1586	Source, detachment from	1640	Lens, toric
1587	Source, leak(s) from	1641	Torqued
1588	Source, incorrect	1642	Tracking
1589	Specificity	1643	Training, inadequate
1590	Speech processor failure	1644	Transducer overheating
1591	Spillage, accidental	1645	Transducer probe overheating
1592	Spring loading mechanism problem	1646	Transducer failure
1593	Stacking breaths	1647	Transmitter failure
1594	Steering wire problem	1648	Trocar failure
1595	Sterility	1649	Trocar/fiber incompatibility
1596	Sterilization, inadequate	1650	Trocar/instrument incompatibility
1597	Sticking	1651	Tube(s), exploding of
1599	Stops intermittently	1652	Tube(s), shattering of
1600	Storage, inadequate	1653	Tube(s), splitting of
1601	Stretched	1654	Tubing, incorrect placement of
1602	Strut fracture	1655	Twisting
1603	Stuck in inspiratory or expiratory phase	1656	Ultrafiltration
1604	Suction, increased	1657	Ultraviolet
1605	Superheat	1658	Ultraviolet absorbing
1606	Supplier	1659	Uncoiled
1607	Suture line disruption	1660	Undercorrection
1608	Suture line separation	1661	Undersensing
1609	Synchronization, inaccurate	1662	Unit inactivated
1610	Syncope	1663	Device, inoperable
1611	Syringe drive problem	1664	Unravel
1612	Syringe, defective	1665	Source, unretracted
1613	Syringe, hole(s) in	1667	Unstable
1614	Syringe markings inaccurate	1668	Unwrap, difficult to
1615	Alarm system, failure of message-battery status	1669	Unwrap, failure to
1616	Alarm system, failure of message-check electrode	1670	Device, incorrect care/use of
1617	Alarm system, failure of message-leads off	1671	Vacuum, loss of
1618	Alarm system, failure of message-service	1672	Valve(s), failure of
1619	System fails to activate	1673	Valve(s), sticking
1620	Alarm system, failure of check-catheter	1674	Vibration
1621	Alarm system, failure of gas-leak/loss	1675	Volume accuracy
1622	Table angulation, unintended	1676	Warning light, incorrect
1623	Table incrementation incorrect	1677	Water softener process, failure of
1624	Table top motion, unintended	1678	Water treatment
		1679	Wedge filter problem
		1680	Wedge, difficult to

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1681	Screw head(s), incorrect	1741	Bleeding, intracranial
1682	Wavelength, incorrect	1742	Blindness
1683	Zero, failure to	1743	Blister(s)
1684	AIDS, acquired immunodeficiency syndrome	1744	Bloating
1685	Pain, abdominal	1745	Contamination from blood/fluids
1688	Abortion	1746	Blood in tubing
1689	Abrasion	1747	Blood pooling
1690	Abscess	1749	Transfusion with incompatible blood
1692	Achalasia	1750	Erosion
1693	Shock, acoustic	1751	Bradycardia
1694	Acoustic trauma	1752	Bronchitis
1695	Adhesion(s)	1754	Bruise
1696	Adult respiratory distress syndrome, (ARDS)	1755	Burn, radiation
1697	Air embolism	1756	Burn, bowel
1698	Airflow, restricted	1757	Burn(s)
1699	Airway obstruction	1758	Calcification
1701	Allergic reaction	1759	Cancer, breast
1702	Amputation	1760	Cancer, other
1703	Shock, anaphylactic	1761	Capsular contracture
1704	Anaphylaxis	1762	Cardiac arrest
1705	Anastomosis	1763	Myocardial contusion
1706	Anemia	1764	Cardiomyopathy
1707	Anesthesia, insufficient, light or patchy	1765	Cardiopulmonary arrest
1708	Aneurysm	1766	Cataract
1710	Angina	1767	Cataract, induced
1711	Anoxia	1768	Cellulitis
1713	Antibiotics, reaction to	1770	Cerebrovascular accident, (CVA)
1714	Anticoagulation	1771	Cerebral, infarction
1715	Aortic insufficiency	1772	Cerebrospinal fluid leakage
1716	Aortic regurgitation	1773	Cervical changes
1717	Aortic stenosis	1774	Cesarean section (c-section)
1718	Apgar score, decreased or low	1775	Chemosis
1719	Apheresis	1776	Chest pain
1720	Apnea	1777	Chorioamnionitis
1721	Arrhythmia	1778	Coagulation
1722	Arteriosclerosis	1779	Coagulopathy
1723	Arthritis	1781	Comatose
1724	Arthritis, rheumatoid	1782	Congenital defect/deformity
1725	Aspiration	1783	Congestive heart failure
1726	Asthma	1784	Conjunctivitis
1727	Asystole	1785	Conjunctivitis, giant papillary
1728	Atherosclerosis	1786	Connective tissue disease
1729	Atrial fibrillation	1787	Contusion
1730	Atrial flutter	1788	Core temperature rise
1731	Atrial tachycardia	1789	Corneal abrasion
1732	Autoimmune disease	1790	Corneal decompensation
1733	Autoimmune reaction	1791	Corneal edema
1734	Automatic injection system, failure to infuse	1792	Cornea, perforation of
1735	Infection, bacterial	1793	Corneal scar
1736	Biliary cirrhosis	1794	Corneal touch
1738	Bleeding	1795	Corneal transplant
1739	Bleeding, cerebral	1796	Corneal ulcer
		1797	Crushing injury
		1798	Cyanosis

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1799	Cyclitis	1858	Fever
1800	Cyst(s), formation of	1859	Fiberoptic fragments, unretrieved in body
1801	Deafness	1860	Fibromyitis
1802	Death/expired	1861	Fibrosis
1803	Debris, bone shedding	1862	Fistula
1804	Debris, metal shedding	1863	Osseointegrate, failure to
1805	Decreased apgar	1864	Flashers
1806	Degeneration	1865	Flatus
1807	Dehydrated	1866	Floater, vitreous
1808	Dementia	1868	Foreign body reaction
1809	Deposits	1869	Foreign body sensation
1810	Device induced	1870	Fracture(s)
1811	Diarrhea	1871	Urgency
1812	Purulent discharge	1872	Fungus
1813	Disseminated intravascular coagulation (DIC)	1873	Gangrene
1814	Dry eye(s)	1874	Gastritis
1815	Dysphagia	1875	Glaucoma
1816	Dyspnea	1876	Granuloma
1817	EKG/ECG changes	1877	Hair loss
1818	Ecchymosis	1878	Haze
1819	Ectopic pregnancy	1879	Head injury
1820	Edema	1880	Headache
1822	Edema, macular	1881	Hearing impairment
1823	Edema, microcytic	1882	Hearing loss
1824	Edema, stromal	1883	Heartburn
1826	Electro-mechanical dissociation	1884	Hematoma
1827	Electrocution	1885	Hemodialysis
1829	Embolism	1886	Hemolysis
1830	Embolus	1887	Hemoptysis
1831	Emotional changes	1888	Hemorrhage
1832	Emphysema, pulmonary	1889	Hemorrhage, cerebral
1833	Encephalopathy	1890	Hemorrhage, extradural
1834	Endocarditis	1891	Hemorrhage, intracranial
1835	Endophthalmitis	1892	Hemorrhage, intraventricular
1838	Enzyme elevation, cardiac	1893	Hemorrhage, subarachnoid
1839	Eructate	1894	Hemorrhage, subdural
1840	Erythema	1895	Hemostasis
1841	Exsanguination	1896	Hemothorax
1842	Extravasation	1897	Hepatitis
1843	Extreme exhaustion	1898	Herpes
1844	Extrusion	1899	Hiccups
1845	Eye injury	1900	Hives
1846	Facial nerve paralysis	1901	Hyaloid face, rupture of
1847	Fainting	1903	Hyperbilirubinemia
1848	Fall	1904	Hyperemia
1849	Fatigue	1905	Hyperglycemia
1850	Feeding problems	1906	Hyperplasia
1851	Asphyxia	1907	Hypersensitivity
1852	Fetal brain injury	1908	Hypertension
1854	Fetal core temperature rise	1909	Hyperthermia
1855	Death, intrauterine fetal	1910	Hyperventilation
1856	Fetal distress	1911	Hyphema
1857	Fetal scalp laceration(s)	1912	Hypoglycemia
		1913	Hypopyon

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
1914	Hypotension	1969	Myocardial infarction, (MI)
1915	Hypothermia	1970	Nausea
1916	Hypoventilation	1971	Necrosis
1917	Shock, hypovolemic	1972	Necrosis of flap tissue
1918	Hypoxia	1974	Neonatal deformities
1919	IUD (intrauterine contraceptive device), removal of	1975	Neonatal hearing impairment
1920	IV catheter pull out	1976	Neonatal hearing loss
1922	Iatrogenic lesion	1978	Neovascularization
1923	Idioventricular rhythm	1979	Nerve damage
1924	Implant, failure of	1980	Nerve stimulation, undesired
1925	Impotence	1981	Neural tissue damage
1926	Insufficiency, valvular	1982	Neurological deficit/dysfunction
1927	Incompetent cervix	1983	Neuropathy
1928	Incontinence	1984	Occlusion
1930	Infection	1985	Reocclusion
1931	Infiltration	1986	Optical nerve damage
1932	Inflammation	1987	Organ(s), perforation of
1933	Infection, intraocular	1988	Overdose
1934	Lens (IOL), dislocated intraocular	1989	Overmedicated
1935	Lens (IOL), migration of intraocular	1990	Oversedated
1936	Intraocular pressure (IOP), delayed, uncontrolled	1991	Overstimulation
1937	Intraocular pressure rise, (IOPR)	1992	Polymyositis
1938	Intravasation	1994	Pain
1940	Iritis	1995	Painful stimulation
1941	Irritation	1997	Paralysis
1942	Ischemia	1998	Paresis
1943	Itching	1999	Peeling
1944	Keratitis	2000	Pelvic inflammatory disease, (PID)
1945	Keratitis, acanthamoeba	2001	Perforation
1946	Laceration(s)	2002	Peripheral vascular disease
1947	Left ventricular dysfunction	2003	Peritoneal laceration(s)
1948	Left ventricular failure	2004	Phlebitis
1949	Left ventricular hypertrophy	2009	Pleomorphism
1950	Lesion	2010	Pleural effusion
1951	Lesions in the target points	2011	Pneumonia
1952	Ligament(s), damage to	2012	Pneumothorax
1953	Liver contusion	2013	Pocket erosion
1954	Liver dysfunction	2014	Megophthalmos
1955	Liver laceration(s)	2015	Positive antinuclear antibodies (ANA)
1956	Lupus	2017	Procedure, improper/incorrect
1957	Cytomegaloviral retinitis	2019	Pulmonary dysfunction
1958	Memory loss	2020	Pulmonary edema
1959	Menstrual irregularities	2021	Pulmonary infarction
1960	Microcyst(s)	2022	Pulmonary insufficiency
1961	Microperforation(s)	2023	Pulmonary regurgitation
1962	Miscarriage	2024	Pulmonary stenosis
1963	Mitral insufficiency	2025	Punctured blood vessels
1964	Mitral regurgitation	2026	Pupillary block
1965	Mitral stenosis	2027	Pus
1966	Muscle spasm(s)	2028	Pyrogenic
1967	Muscle weakness	2031	Radiation underexposure
1968	Muscular rigidity	2032	Range of motion, loss of
		2033	Rash
		2034	Raynaud's phenomenon

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
2035	Reaction, local	2093	Swollen lymph nodes
2036	Reaction, systemic	2094	Synovitis
2038	Red eye(s)	2095	Tachycardia
2039	Renal disease (ESRD), end stage	2096	Body temperature, elevated
2041	Renal failure	2097	Tentorial tears
2042	Surgical procedure, repeated	2098	Toxins in children
2043	Respirator, detachment from	2099	Therapeutic effects, unexpected
2044	Respiratory arrest	2100	Thrombosis
2045	Respiratory distress	2101	Thrombus
2046	Respiratory distress syndrome (RDS) of newborns	2102	Thyroid problems
2047	Retina, detached	2103	Tinnitus
2048	Retina, damage to	2104	Tissue damage
2049	Retina, degeneration of	2105	Tissue failure
2050	Retina, tear(s) in	2106	Tissue, breakdown of optical
2051	Retrograde	2107	Torsades-de-Pointes
2052	Retrograde flow	2108	Toxic shock syndrome, (TSS)
2053	Rheumatic heart disease	2109	Transient ischemic attack
2054	Right ventricular dysfunction	2110	Transplant of organ/tissue
2055	Right ventricular failure	2111	Tricuspid insufficiency
2056	Right ventricular hypertrophy	2112	Tricuspid regurgitation
2058	Staphylococcus aureus	2113	Tricuspid stenosis
2059	ST segment elevation	2114	Twiddler's syndrome
2060	Scar tissue	2115	UGH (uveitis-glaucoma-hyphema) syndrome
2061	Scarring	2116	Ulceration
2062	Scleroderma	2117	Under no medication
2063	Seizures	2118	Under no sedation
2064	Seizures, focal motor	2119	Urinary retention
2065	Sensitivity	2120	Infection, urinary tract
2066	Sensitization	2121	Uterine perforation
2067	Sepsis	2122	Uveitis
2068	Shock, septic	2123	Vaginal discharge
2069	Seroma	2124	Vaginal mucosa damage
2070	Severed digit(s), (finger or toe)	2126	Vasoconstriction
2071	Venereal disease (VD)	2127	Vasodilatation
2072	Shock	2128	Vasospasm
2073	Sjogren's syndrome	2129	Venipuncture
2074	Skin discoloration	2130	Ventricular fibrillation
2075	Skin erosion	2131	Ventricular flutter
2076	Skin irritation	2132	Ventricular tachycardia
2077	Fracture, skull	2133	Ventriculomeglia
2078	Ventricle, abnormality of	2134	Vertigo
2079	Soreness	2135	Vessels, perforation of
2081	Spinal injury	2136	Virus
2082	Spotting	2137	Vision, blurring of
2083	Sprain	2138	Vision, impaired
2084	Strangulation	2139	Vision, loss of
2085	Stricture	2140	Visual disturbances
2086	Stroke	2142	Vitreous fluid, loss of
2087	Subepithelial infiltration	2143	Vitreous fluid, blood in
2088	Suffocation	2144	Vomiting
2089	Sunset syndrome	2145	Weakness
2091	Swelling	2146	Burning sensation
2092	Swollen glands	2147	Calcium deposit(s)
		2149	Electro-static discharge

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
2150	Epithelial marsupialization	2206	Heart failure
2151	Film processor failure	2207	Toxemia
2152	Great vessel perforation	2208	Rupture
2153	Hot flashes	2209	Asphyxia, fetal
2154	Implant extrusion	2210	Hypoxia, fetal
2155	Improved patient condition	2211	Abortion, artificial
2156	Immuno-deficiency	2212	Abortion, complete
2157	Interlock(s), failure of	2213	Abortion, incomplete
2158	Lens aberration, distortion of	2214	Abortion, induced
2159	Misdiagnosis	2215	Abortion, missed
2160	Motion detector failure	2216	Abortion, therapeutic
2161	Muscular tics	2217	Amnionitis
2162	Nipple sensation, changes in/loss of	2218	Anaphylactoid
2163	Noise	2219	Brain injury
2164	Phosphene visualization	2220	Cytomegalovirus (CVM)
2165	Phototoxicity	2221	Convulsion
2166	Radiation underdose	2222	Convulsion, clonic
2167	Safety interlock(s) inadequate	2223	Convulsion, tonic
2168	Seizures, grand-mal	2224	Dementia, dialysis
2169	Seizures, petit-mal	2225	Discharge
2170	Suction failure	2226	Cardiac tamponade
2171	Tingling	2227	Halo
2172	Twitching	2228	Clouding, central corneal
2173	Corneal implant	2229	Blinking, excessive
2174	Infiltrates	2230	Blurring
2177	Surgery, prolonged	2231	Corneal infiltrates
2178	Staining	2232	Microcysts, epithelial
2181	Vitritis	2233	Corneal sensitivity
2182	Underinfusion	2234	Overwear syndrome
2183	Fitting problems	2235	Tearing, excessive
2184	Lens replacement	2236	Endotoxin
2186	Headache, lumbar puncture	2237	Chronic obstructive pulmonary disease, (COPD)
2187	Jaundice	2238	Myalgia
2188	Uremia	2239	First use syndrome
2189	Collagen disease	2240	Hernia
2191	Chills	2241	Wheal(s)
2192	Concussion	2242	Hypernatremia
2193	Cramp(s)	2243	Hypovolemia
2194	Dizziness	2244	Infection, direct
2195	Dysphasia	2245	Infection, indirect
2196	Electrolyte imbalance	2246	Infection, pyrogenic
2197	HIV, human immunodeficiency virus	2247	Infection, subclinical
2198	Mediastinal shift	2248	Infection, viral
2199	No consequences or impact to patient	2249	Migration
2200	Other (for use when an appropriate patient code cannot be identified)	2250	Nonpyrogenic
2201	Leak(s), perivalvular	2251	Sneezing
2202	Unknown (for use when the patient's condition is not known)	2252	Peritonitis
2203	Other (for use when an appropriate device code cannot be identified)	2253	Dermatomyositis
2204	Unknown (for use when the device problem is not known)	2254	Reaction, pyrogenic
2205	Perforation	2255	Radiodermatitis
		2256	Radiation sickness syndrome
		2259	Regurgitation
		2260	Seizures, focal

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
2261	Seizures, absence	2314	Radiofrequency interference (RFI)
2262	Shock, cardiogenic	2315	Water purification system, failure of
2263	Stenosis	2316	Fungus
2264	Shock, insulin	2317	Blood contaminated device
2265	Shock, neurogenic	2318	Blood pooling
2266	Shock, postoperative	2319	Cardiac enzyme evaluation, erroneous
2267	Shock, surgical	2320	Implant, removal of
2268	Shock, traumatic	2321	Implant, repositioning of
2269	Toxoplasmosis, acquired	2322	Implant, reprogramming of
2270	Toxoplasmosis, congenital	2323	Surgical graft, failure of
2271	Therapeutic response, decreased	2324	Solder joint failure
2272	Therapeutic response, increased	2325	Residue
2273	Teratogenic effects	2326	Pressure sores/ulcers
2274	Ulcer	2327	Entrapment
2275	Urinary frequency	2328	Anxiety
2276	Nonresorbable materials, unretrieved in body	2329	Distress
2277	Sinus, perforation of	2330	Discomfort
2278	Urticaria	2331	Complaint, ill-defined
2279	Hemolytic anemia	2332	Blockage
2280	Prompts, inappropriate	2333	Toxicity
2281	Prompts will not clear	2334	Cardiac insufficiency
2282	Prompts, no voice	2335	Regurgitation, valvular
2283	Prompts, inaudible voice	2336	Stroke syndrome
2284	Damage, internal/external	2337	Transfusion of blood products
2285	Impedance, low	2338	Deliver, failure to
2286	Algorithms, inconsistent	2339	Delivery, inaccurate
2287	Shock, intermittent	2340	Infuse, failure to
2288	Charge, aborted	2341	Infusion, intermittent
2289	Arcing at electrodes	2342	Overdelivery
2290	Blood gas measurements, erroneous	2343	Underdelivery
2291	Semiautomatic code, failure to override	2345	Electro-magnetic interference (EMI), compatibility/incompatibility
2292	Components, defective	2346	Tube(s), buckling of
2293	Sensing, invalid	2347	Wire(s), breakage of
2294	Electrical wires, defective	2348	Injury
2295	Water softener regeneration cycle, mistiming of	2349	Fracture, hip
2296	Temperature probe, loose	2351	Fracture, arm
2297	Bubble detector, failure of	2352	Hypoesthesia
2298	Valve(s), defective	2353	Hypoesthesia, arm/hand
2299	Tube(s), defective	2354	Hypoesthesia, foot/leg
2300	Air eliminator, defective	2355	Arthralgia
2301	Seal, defective	2356	Joint swelling
2302	Cable, defective	2357	Surgical procedure
2303	Contamination, bacterial	2358	Scar excision
2304	Stopcock valve, failure of	2359	Malaise
2305	Burn hole(s)	2360	Disfigurement
2306	Components, missing	2361	Depression
2307	Automatic injection system underinfusion	2362	Peroneal nerve palsy
2308	Filtration process, inadequate	2363	Dyskinesia
2309	Rinsing, improper	2364	Diabetic ketoacidosis
2310	Specialty beds, rips/tears in	2365	Foreign body, removal of
2311	Out-of-box failure	2366	Laparotomy
2312	Packaging, incomplete/missing	2367	Pharyngitis
2313	Radiation overexposure	2368	Sedation

NUMERIC LISTING (continued)

CODE	TERM	CODE	TERM
2369	Fracture, delayed union	2424	Patch test, abnormal results of
2370	Mammogram, abnormal	2425	pH, low
2371	Disability	2426	pH, high
2372	Apicectomy	2427	Teeth, sensitivity of
2373	Joint disorder	2428	Fracture, tooth
2374	Joint dislocation	2429	Encephalitis
2375	Fasciitis	2430	Forced expiratory volume (FEV), decreased
2376	Mitral valve replacement	2431	Forced expiratory volume (FEV), increased
2377	Osteolysis	2432	Spinal cord injury
2378	Healing, impaired	2433	Pain, neck
2379	Device failure	2434	Stiffness, neck
2381	Laparoscopic sterilization	2435	Peak expiratory flowrate, decreased
2382	Meter failure	2436	Peak expiratory flowrate, increased
2383	Internal fixation, revision of	2437	Bronchopneumonia
2384	Lens implant	2438	Breast neoplasm
2385	Packaging, tears, rips, holes in	2439	Breast lumps
2386	Alarm, failure of high inspiratory pressure	2440	Calibrate, failure to
2387	Alarm, failure of low inspiratory pressure	2441	Nipple ulceration
2388	Pain relief, inadequate	2442	Reaction, injection site
2389	Meningitis	2443	Skin inflammation
2390	Arachnoiditis, spinal	2444	Sweating
2391	Pressure, insufficient	2445	Vitreous, detachment of
2392	Reintubate	2446	Wound infection, post-operative
2393	Tracheostomy	2447	Wound infection, post-traumatic
2394	Brain damage	2448	Paraplegia
2395	Ventilator dependent	2449	Quadriplegia
2396	Sore throat	2450	Blood glucose, low
2397	Lung, overinflation of	2451	Blood glucose, high
2398	Esophagus, laceration(s) of	2452	Diaphoretic
2399	Esophagus, perforation of	2453	Cut(s)
2400	Resuscitation	2454	Cut(s)
2401	Breathing difficulties	2455	Contamination, chemical
2402	Extubate	2456	Test results, inaccurate
2403	Reinfusion	2457	Test results, high
2404	Overfill	2458	Test results, low
2405	Syringe markings, incorrectly placed	2459	Readings, high
2406	Identified diagnosis	2460	Readings, low
2407	Dull	2461	Body temperature, decreased
2408	Osseointegration, loss of	2462	Needle stick/puncture
2409	Malfunction	2463	Chest tightness/pressure
2410	Miscalibration	2464	Choking
2411	Nonspecific	2465	Labor, premature
2412	Dissatisfaction	2466	Nasal obstruction
2413	Leaflet fracture	2467	Palpitations
2414	Reaction	2468	Pallor
2415	Numbness	2469	Pulse, irregular
2416	Collapse	2470	Myocarditis
2417	Coma	2471	Mitral valve prolapse
2418	Consciousness, loss of	2472	Aortic valve replacement
2419	Infection, fungal	2473	Tricuspid valve replacement
2420	Infection, upper respiratory tract	2474	Pulmonary valve replacement
2421	Irritability	2475	Prolapse
2422	Obstruction	2476	Heart valve repair
2423	Obstruction	2477	Oxygen saturation, low

NUMERIC LISTING (continued)

CODE	TERM
2478	Oxygen saturation, high
2479	Pulmonary arterial wedge pressure, low
2480	Pulmonary arterial wedge pressure, high
2481	Pulmonary arterial wedge pressure, normal
2482	Respiratory acidosis
2483	Respiratory alkalosis
2484	Respiratory failure
2485	Respiratory rate, decreased
2486	Respiratory rate, increased
2487	ST segment depression
2488	Blood pressure, high
2489	Blood pressure, low
2490	Unresponsive

PART I, SUBPART C - DEFINITIONS FOR EVENT TERMS LISTED NUMERICALLY BY EVENT CODE

CODE	DEFINITION
1003	The uptake of substances into or across tissues; e.g. skin, intestine, and kidney tubules. In radiology, the taking up of energy by matter with which the radiation interacts.
1007	The clumping together in suspension of antigen-bearing cells, microorganisms, or particles in the presence of specific antibodies.
1027	The temperature of the surrounding medium, such as gas or liquid, which comes into contact with the apparatus.
1028	Failure to create an opening by surgical, traumatic, or pathological means between two normally separate spaces or organs.
1036	Anything in a graphic record or medical image that is caused by the technique used to produce the record and is not a natural occurrence but merely incidental.
1040	A lens in which one or both surfaces are not spherical, so designed to minimize certain optical aberrations.
1041	The removal by suction of excess fluid or gas from a body cavity.
1042	The removal by suction of excess fluid or gas from a body cavity.
1043	Loss of size, degree or effect.
1050	Lack of proportion; not symmetrical.
1060	A salt of carbonic acid; an acid carbonate, as sodium bicarbonate.
1061	Lens having two portions of different focal power.
1064	Backflow of blood and/or fluid/solution.
1073	A thin, ragged fin left on the edge of a piece of metal by a cutting or punching tool.
1077	A calcified substance or part.
1079	Use of a capacitor to transfer energy from one circuit to another.
1080	Ineffective and inconsistent depolarization of the heart resulting from the electrical stimulus of the pacemaker.
1081	Inability to achieve effective and consistent depolarization of the heart resulting from the electrical stimulus of the pacemaker.
1082	Use of activated carbon to remove a variety of organic chemicals in the water used in hemodialysis.
1083	The application of a caustic substance, a hot instrument, an electric current or other agent to destroy tissue.
1084	Blood leakage around the catheter upon insertion into the patient.
1086	To scorch or become scorched.
1095	The aggregation of particles, such as bacteria, into irregular masses.
1096	To congeal, solidify, thicken, curdle.
1097	Measurement value associated with dialyzer clearance.
1101	A collection or group of bacteria in a culture derived from the increase of an isolated single organism or group of organisms.
1113	Chemicals, either liquid or powder, used to make dialysate in hemodialysis.
1115	Measurement of the ionic concentration of a solution used in dialysis.
1121	Continuous effective contact of all components of an electric circuit to give it high conductance by providing low resistance.
1123	Fire conducted at a normal rate without interruption, for application of adjustment corrections or for other causes.
1124	Continuous failure of one of several alternate conditions or methods of operation.
1125	Any special symptom or circumstance that renders the use of a remedy or the carrying out of a procedure inadvisable, usually because of the risk.
1133	A high intensity direct current shock delivered to the heart to interrupt ventricular fibrillation and restore synchronous electrical activity.
1137	The degree to which an antibody or antigen participates in cross reactions.
1139	Chemical and/or material (heart valve or tracheal) deterioration.
1140	The point at which there is a transition from spiral flow in the housing of a centrifugal fan to a straight line flow in the connected duct.
1145	Preventing transfer or feedback of energy from one circuit to another.
1147	A photo tube with one or more dynodes between its photocathode and the output electrode.
1153	The reduction of a chemical compound to one less complex, as if by splitting off one or more groups.
1154	Separation of the layers of a surgical wound; it may be partial and superficial only, or complete with disruption of all layers.

DEFINITIONS (continued)

- 1155** Removes cations and anions used in water treatment systems for hemodialysis.
- 1156** Division into separate layers.
- 1157** Spread out, utilize or arrange.
- 1165** That part of a mixture that passes through a dialyzing membrane.
- 1166** The apparatus for performing dialysis; a membrane used in dialysis.
- 1191** Slow movement away from the normal or original position.
- 1194** Interference, generally at radio frequencies, that is generated inside systems, as contrasted to radio-frequency interference coming from sources outside a system.
- 1197** The immediate effects produced by the passage of an electric current through any part of the body, e.g., painful stimulation of nerves or tetanic contractions of muscles.
- 1204** Device and/or fragments of device are embedded in patient's vessel and/or plaque.
- 1212** Device and/or accessories caught within patient vasculature or tissue.
- 1219** Coated with glass or plastic fibers having special optical properties.
- 1225** A test result which erroneously excludes an individual from a specific diagnostic or reference group, due particularly to insufficiently exact methods of testing.
- 1227** A test result which erroneously assigns an individual to a specific diagnostic or reference group, due particularly to insufficiently exact methods of testing.
- 1268** An agent that kills pathogenic microorganisms.
- 1270** Rate of change of temperature, pressure, or other variable as a function of distance, time, etc.
- 1278** A haptic is a loop or foot of an intraocular lens implant that supports the lens against the iris.
- 1288** An increase in the concentration of blood cells resulting from the loss of plasma or water from the blood stream.
- 1289** A device that corrects abnormalities of the blood by bulk convection of solutes and water through a semipermeable membrane.
- 1290** The rate per unit of area at which water passes through a semi-permeable membrane, such as those used for ultrafiltration or reverse osmosis.
- 1291** The opposition to the flow of an alternating current, which is the vector sum of ohmic resistance plus additional resistance, if any, due to induction, to capacity or to both.
- 1297** Growth of tissue in or around a foreign body as the body's antibody response to the foreign body.
- 1301** Hydrogen peroxide sterilization is a method of sterilizing soft contact lenses in a peroxide solution.
- 1320** Loss of continuity of insulation via wear, split, crack and/or breach.
- 1328** See also: Electro-magnetic interference (EMI) codes.
- 1329** See also: Electro-magnetic interference (EMI) codes.
- 1330** Two or more lumens intercrossed.
- 1331** To give or send out a signal to a transponder or computer for triggering an appropriate response.
- 1333** Pertaining to the inner layer of the blood vessels.
- 1336** The ensheathing, infolding or insertion of a structure within itself or another.
- 1338** Relating to motion or movement.
- 1341** An extract of blood cells from the horseshoe crab is exposed to a blood sample from a patient to test for gram negative endotoxin.
- 1352** Leakage of blood around leaflet disc.
- 1353** A separation of the cusp of a heart valve.
- 1364** The development of opacity as of the cornea or lens.
- 1374** A milliamps is a measurement of the flow of electricity.
- 1376** Alternation in the activation of a pulse generator reed switch by a magnet.
- 1377** Any type of interference that reduces the intensity of fluorescence.
- 1411** A lens with various dioptric powers, such as bifocal, trifocal or a progressive lens.
- 1413** Death of tissue due to external pressure.
- 1419** Free of fever producing substances.
- 1423** The act of closing or the state of being closed; an obstruction.
- 1426** Impervious to light; not translucent or only slightly so.
- 1428** The formation of bone or of a bony substance; the conversion of fibrous tissue or of cartilage into bone or a bony substance.

DEFINITIONS (continued)

- 1429** Collision occurring with the patient or other room components caused by device failure or unintended system motion.
- 1430** Any motion of the system or components that was not initiated by the user.
- 1438** The pacemaker does not release a stimulus when it is needed because it senses intracorporeal and extracorporeal signals as intrinsic cardiac activity.
- 1441** The occurrence at distinct times of events normally synchronous; disturbance of coordination.
- 1447** Superficial vascularization of the cornea with infiltration of granulation tissue; an inflammatory exudate overlying the lining layer of synovial cells on the inside of a joint.
- 1448** The property of being attracted by a magnet, and of assuming a position parallel to that of a magnetic force, but not of becoming permanently magnetized.
- 1449** An arbitrary constant in a mathematical expression that distinguishes specific cases; a parameter has a definite fixed value in one case but will have different values in another case.
- 1451** Particles introduced by the device during the procedure.
- 1455** Effected or performed through the skin.
- 1457** The escape of blood around a heart valve, particularly around its leaflets.
- 1458** In radiology, a device that stimulates the conditions encountered when radiation or radioactive material is deposited in vivo and permits a quantitative estimation of its effects.
- 1463** Electrical stimulation of the patient's pocket of skin in which the pulse generator is housed.
- 1470** A potentiometer is a device used for the measurement of an electromotive force by comparison with a known potential difference.
- 1471** A potentiometer is an device used for the measurement of an electromotive force by comparison with a known potential difference.
- 1478** A substance separating, in solid particles, from a liquid as the result of a chemical or physical change; to form a precipitate.
- 1498** The closure of the pulmonary artery or one of its branches by an embolus, sometimes associated with infarction of the lung.
- 1505** A fever producing substance.
- 1507** The occurrence of a premature ventricular complex near the peak of the T wave in electrocardiography; it may lead to ventricular tachycardia or fibrillation.
- 1518** A tube for insertion into a duct or cavity.
- 1522** A backward or return flow.
- 1538** Membrane separation process for removing solvent from a solution; typically used as a water purification process in hemodialysis.
- 1556** The sediment filters used to remove particulates from the water supply, in order to protect other water treatment equipment or dialysate delivery devices, become clogged or damaged.
- 1578** To provide with or divert by means of an electrical shunt; to divert blood from one part to another by a surgical shunt.
- 1593** Condition where there are a greater number of inhalations relative to the number of exhalations.
- 1602** Breakage of supporting pieces of a heart valve.
- 1610** A temporary suspension of consciousness due to generalized cerebral ischemia; a faint or swoon.
- 1629** Discrepancies occurring during automatic data measurement and transmission, as by wire or radio, from remote sources, to a receiving station for recording and analysis.
- 1631** Use of a medical device for therapeutic purposes.
- 1633** Loss of the minimum amount of energy, voltage, or current needed to consistently stimulate the heart muscle.
- 1634** The amount of gas that is inspired and expired during one respiratory cycle.
- 1640** Lens with a cylindrical component; used for correcting an eye's astigmatic refractive error.
- 1648** A sharp-pointed instrument equipped with a cannula, used to puncture the wall of a body cavity and withdraw fluid.
- 1656** The transfer of fluid between the blood and dialysate through the dialysis membrane due to a pressure gradient (transmembrane pressure) existing between the blood and dialysate compartments.
- 1677** Ionic exchange process used in hemodialysis to treat water.
- 1684** A disease characterized by opportunistic infections.
- 1688** The premature expulsion from the uterus of the products of conception, of the embryo or of a nonviable fetus.

DEFINITIONS (continued)

- 1689** Abraded wound; excoriation or circumscribed removal of the superficial layers of the skin or mucous membrane.
- 1690** Circumscribed collection of pus appearing in an acute or chronic, localized infection and associated with tissue destruction and frequent swelling.
- 1692** Failure of the smooth muscle fibers of the gastrointestinal tract to relax at any one point of junction of one part with another.
- 1695** The process of adhering or uniting of two surfaces or parts, especially the union of the opposing surfaces of a wound.
- 1696** Fulminant pulmonary interstitial and alveolar edema, which usually develops a few days after the initiating trauma.
- 1697** The presence of bubbles of gas in the vascular system; occurrence is related to the entry of air into the venous circulation following trauma or surgery.
- 1701** Hypersensitivity; a local or general reaction of an organism following contact with a specific allergen to which it has been previously exposed and to which it has become sensitized.
- 1703** Pertaining to anaphylaxis.
- 1704** Systemic anaphylaxis is the most dramatic example of an immediate hypersensitivity reaction; is uncommon and unexpected in nature and occasionally results in a fatal outcome.
- 1705** An opening created by surgical, traumatic or pathological means between two normally distinct spaces or organs.
- 1706** Equilibrium between blood loss and blood production is disturbed, resulting in a reduction in the number of erythrocytes.
- 1707** The amount of anesthesia administered to the patient was inadequate.
- 1708** A sac formed by the dilation of the wall of an artery, a vein or the heart.
- 1710** Severe constricting pain, usually relating to the heart muscle.
- 1711** Absence or lack of oxygen; reduction of oxygen in body tissues below normal levels.
- 1714** Serving to prevent the coagulation of blood. Any substance that delays, or nullifies coagulation of the blood.
- 1715** Defective functioning of the aortic valve, with incomplete closure resulting in aortic regurgitation.
- 1716** The backward flow of blood from the aorta into the left ventricle, owing to insufficiency of the aortic semilunar valve; it may be chronic or acute.
- 1717** Narrowing of the orifice of the aortic valve or of the supra- or subvalvular regions.
- 1718** A numerical expression of a newborn infant's, physical condition, usually determined at 60 seconds after birth.
- 1719** Any procedure in which blood is withdrawn from a donor, a portion is separated and retained, and the remainder is retransfused into the donor.
- 1720** Cessation of breathing.
- 1721** Any variation from the normal sinus rhythm of the heart.
- 1722** A group of diseases characterized by thickening and loss of elasticity of arterial walls.
- 1724** A chronic, systemic disease primarily of the joints, usually polyarticular, marked by inflammatory changes in the synovial membranes and articular structures and by muscle atrophy and rarefaction of the bones.
- 1725** The inspiratory sucking into the airways of fluid or a foreign body.
- 1727** Cardiac standstill or arrest; absence of a heartbeat.
- 1728** An extremely common form of arteriosclerosis in which deposits of yellowish plaques containing cholesterol, lipid material and lipophages are formed within the intima and inner media of large and medium sized arteries.
- 1729** An arrhythmia in which minute areas of the atrial myocardium are in various uncoordinated stages of depolarization and repolarization; instead of intermittently contracting, the atria quiver continuously in a chaotic pattern, causing a totally irregular, often rapid ventricular rate.
- 1730** A condition of cardiac arrhythmia in which the atrial contractions are rapid (250 to 350 per minute) but regular.
- 1731** A rapid cardiac rate usually between 160 and 190 beats per minute.
- 1732** Any disorder in which destruction of normal tissue arises from humoral or cellular immune responses of the individual to their own tissue constituents; it may be systematic or organ specific.
- 1733** The immune response in which antibodies or immune lymphoid cells are produced against the body's own tissues.
- 1736** Cirrhosis of the liver due to obstruction or infection of the major extra- or intrahepatic bile ducts.
- 1739** Bleeding within the cerebrum.
- 1741** Bleeding within the cranium.

DEFINITIONS (continued)

- 1742** Loss of the sense of sight.
- 1745** The patient's blood becomes contaminated from contact with blood/fluids on a device.
- 1747** Unintended collection/stasis of blood in patient.
- 1750** An eating away; destruction of the surface of a tissue, material or structure.
- 1751** Slowness of the heart beat, as evidenced by slowing of the pulse rate to less than 60.
- 1754** An injury of a part without a break in the skin.
- 1755** A burn caused by exposure to x-ray, radium, sunlight, atomic or any other type of radiant energy.
- 1757** Injury to tissues caused by contact with dry heat, moist heat, chemicals, electricity, friction or radiant and electromagnetic energy. A first degree burn is associated with redness, a second degree burn with vesication and a third degree burn with necrosis through the entire skin.
- 1758** The process by which organic tissue becomes hardened by a deposit of calcium salts within its substance.
- 1761** The tightening of scar tissue that forms around the implant.
- 1762** Sudden cessation of the pumping function of the heart, with disappearance of arterial blood pressure, connoting either ventricular fibrillation or ventricular standstill.
- 1763** A bruise to the heart.
- 1764** A general diagnostic term designating primary noninflammatory disease of the heart muscle, often of obscure or unknown etiology.
- 1765** Cessation of breathing and/or cardiac function.
- 1766** A loss of the transparency of the lens of the eye or its capsule.
- 1767** A cataract that has been produced artificially or by induction, e.g. as a result of device use, medication, trauma, tears, falls, accidental injury, etc.
- 1768** An acute, diffuse, spreading, edematous, suppurative inflammation of the deep subcutaneous tissues and sometimes muscle, which may be associated with abscess formation.
- 1770** See also Stroke Syndrome.
- 1771** An ischemic condition of the brain, producing a persistent focal neurological deficit in the area of distribution of the cerebral arteries.
- 1773** Changes such as abnormal, pathological, benign, malignant, etc.
- 1775** Excessive edema of the ocular conjunctiva.
- 1777** Inflammation of the chorion and amnion.
- 1778** Process of clot formation.
- 1779** Any disorder of blood coagulation.
- 1783** A clinical syndrome due to heart disease, characterized by breathlessness and abnormal sodium and water retention, often resulting in edema.
- 1784** Inflammation of the conjunctiva.
- 1785** Allergic type of conjunctival inflammation associated with continuous wearing of soft contact lenses. Hard, flat papillae form a cobblestone pattern on undersurface of upper eyelid.
- 1786** A group of generalized diseases affecting the connective tissue.
- 1787** A bruise; an injury of a part without a break in the skin.
- 1788** Elevation of total body temperature.
- 1789** Scraped area of corneal surface, accompanied by a loss of superficial tissue.
- 1790** Inability to maintain corneal integrity; or corneal edema resulting from failure of the corneal endothelium to maintain detumescence.
- 1791** Hazy, swollen cornea.
- 1794** Complications of intraocular lens surgery; usually refers to intraocular lens contact with cornea, can be intermittent or chronic.
- 1796** Area of epithelial tissue loss from corneal surface; associated with inflammatory cells in the cornea and anterior chamber.
- 1797** Any type of crushing injury associated with device use or operation.
- 1798** A bluish discoloration, applied especially to such discoloration of skin and mucous membranes due to excessive concentration of reduced hemoglobin in the blood.
- 1799** Inflammation of the ciliary body.
- 1803** Shedding of accumulated bone fragments.
- 1806** Deterioration; change from a higher to a lower form; especially change of tissue to a lower or less functionally active form.

DEFINITIONS (continued)

- 1808** An organic mental syndrome characterized by a general loss of intellectual abilities involving impairment of memory, judgement, and abstract thinking as well as changes in personality.
- 1812** Consisting of or containing pus; associated with the formation of or caused by pus.
- 1813** A disorder characterized by reduction in the elements involved in blood coagulation due to their utilization in widespread blood clotting within the vessels.
- 1815** Difficulty in swallowing.
- 1816** Difficult or labored breathing.
- 1817** Changes in cardiac electrical activity.
- 1818** A small hemorrhagic spot in the skin or mucous membrane forming a nonelevated, rounded, or irregular blue or purplish spot.
- 1819** Development of the fertilized ovum outside of the uterine cavity.
- 1820** Accumulation of an excessive amount of fluid in or around cells, tissues or body cavities.
- 1822** Swelling of the macula.
- 1823** Swelling of abnormally small red blood cells.
- 1826** Continued electrical rhythmicity of the heart in the absence of effective mechanical function.
- 1827** The passage of electrical current through the body.
- 1829** Sudden blocking of an artery by a clot or foreign material which has been brought to its site of lodgement by the blood current.
- 1830** A mass of clotted blood or other formed elements, such as bubbles of air, calcium fragments, etc. brought by the blood from another vessel and forced into a smaller one, thus obstructing the circulation.
- 1832** A condition of the lung characterized by increase beyond normal in the size of air spaces distal to the terminal bronchioles, either from dilatation of the alveoli or from destruction of their walls.
- 1833** Any degenerative disease of the brain.
- 1834** Inflammation of the endocardium (the endothelial lining membrane of the heart and the connective tissue bed on which it lies).
- 1835** Inflammation involving the ocular cavities and their adjacent structures.
- 1839** To belch, the casting of upwind from the stomach.
- 1840** Redness of the skin produced by congestion of the capillaries.
- 1841** Extensive loss of blood due to internal or external hemorrhage.
- 1842** A discharge or escape, as of blood, from a vessel into the tissues.
- 1843** Extreme fatigue; inability to respond to stimuli.
- 1844** Thrusting or pushing out; expulsion by force.
- 1846** Weakening or paralysis of the facial nerve.
- 1847** Extremely weak; threatened with syncope.
- 1851** Stopping of the pulse. A condition due to lack of oxygen in respired air, resulting in impending or actual cessation of life.
- 1852** Injury occurring as a result of the delivery process.
- 1854** Elevation of total body temperature.
- 1855** Death in utero; failure of the product of conception to show evidence of respiration, heart beat, or definite movement of a voluntary muscle after expulsion from the uterus, with no possibility of resuscitation.
- 1858** Elevation of body temperature above normal.
- 1860** Inflammation and fibrous degeneration of a muscle.
- 1861** The formation of fibrous tissue; fibroid or fibrous degeneration.
- 1862** An abnormal passage or communication, usually between two internal organs, or leading from an internal organ to the surface of the body.
- 1863** Failure of direct anchorage of an implant by the formation of bony tissue around the implant without the growth of fibrous tissue at the bone-implant interface.
- 1864** A sudden or brief burst of light.
- 1865** Gas or air in the gastrointestinal tract.
- 1866** "Spots before the eyes"; deposits in the vitreous of the eye, usually moving about and probably representing fine aggregates of vitreous protein occurring as a benign degenerative change.
- 1868** A granulomatous inflammatory reaction evoked by the presence of an exogenous material in the tissues, a characteristic feature of which is the formation of foreign body giant cells.

DEFINITIONS (continued)

- 1869** Feeling of grittiness or having something in the eye; frequently caused by a foreign body. Other possible causes include corneal abrasion, corneal ulcer, inturned eye lash or acute conjunctivitis.
- 1870** The breaking of a part; e.g. a bone, tooth, etc.
- 1871** A sudden compelling urge to urinate.
- 1873** Death of tissue, usually in considerable mass and generally associated with loss of vascular (nutritive) supply and followed by bacterial invasion and putrefaction.
- 1874** Inflammation of the stomach.
- 1875** A group of eye diseases characterized by an increase in intraocular pressure which causes pathological changes in the optic disk and typical defects in the field of vision.
- 1876** Any small nodular delimited aggregation of mononuclear inflammatory cells.
- 1884** A localized collection of blood, usually clotted, in an organ, space, or tissue, due to a break in the wall of a blood vessel.
- 1885** The removal of certain elements from the blood by virtue of the difference in the rates of their diffusion through a semipermeable membrane.
- 1886** The destruction of red blood cells and the subsequent release of hemoglobin into the plasma.
- 1887** The expectoration of blood or of blood-stained sputum.
- 1888** The escape of blood from the vessels; bleeding.
- 1889** Hemorrhage into the cerebrum. See also Stroke Syndrome.
- 1890** Intracranial hemorrhage into the epidural space.
- 1891** Bleeding within the cranium, which may extradural, subdural, subarachnoid, or cerebral.
- 1892** Extravasation of blood into the ventricular system of the brain.
- 1893** Intracranial hemorrhage into the subarachnoid space.
- 1894** Cerebral hemorrhage into the subdural space. See also Stroke Syndrome.
- 1895** The arrest of bleeding, either by the physiological properties of vasoconstriction and coagulation or by surgical means.
- 1896** A collection of blood in the pleural cavity.
- 1897** Inflammation of the liver.
- 1898** Any inflammatory skin disease caused by a herpesvirus and characterized by the formation of small vesicles in clusters.
- 1900** Urticaria.
- 1903** Excessive concentrations of bilirubin in the blood, which may lead to jaundice.
- 1904** The presence of an increased amount of blood in a part or organ; engorgement.
- 1905** Abnormally high concentration of glucose in the circulating blood.
- 1906** The abnormal multiplication or increase in the number of normal cells in normal arrangement in a tissue.
- 1907** A state of altered reactivity in which the body reacts with an exaggerated immune response to a foreign substance.
- 1908** Persistently high arterial blood pressure.
- 1909** Abnormally high body temperature; especially that induced by therapeutic purposes.
- 1910** Abnormally prolonged, rapid, and deep breathing.
- 1911** Hemorrhage within the anterior chamber of the eye; bloodshot.
- 1912** An abnormally small concentration of glucose in the blood.
- 1913** An accumulation of pus in the anterior chamber of the eye.
- 1914** Abnormally low blood pressure.
- 1915** Abnormally low body temperature.
- 1916** A state in which there is a reduced amount of air entering the pulmonary alveoli.
- 1917** Shock resulting from insufficient blood volume for the maintenance of adequate cardiac output, blood pressure and tissue perfusion.
- 1918** Decrease below normal levels of oxygen in inspired gases, arterial blood or tissue.
- 1922** Denoting any adverse condition in a patient occurring as the result of treatment by a physician or surgeon, especially infections acquired by the patient during the course of the treatment.
- 1923** Relating to or affecting the cardiac ventricles alone.
- 1924** An object or material, such as an alloplastic or radioactive material or tissue, partially or totally inserted or grafted into the body for prosthetic, therapeutic, diagnostic, or experimental purposes.
- 1925** Lack of power.

DEFINITIONS (continued)

- 1926** Dysfunction of one of the cardiac valves, with incomplete valve closure resulting in valvular regurgitation.
- 1927** One that is abnormally prone to dilate in the second trimester of pregnancy, resulting in premature expulsion of the fetus.
- 1928** Inability to control excretory functions.
- 1931** The diffusion or accumulation in a tissue or cells of substances not normal to it or in amounts in excess of the normal.
- 1933** Infection within the eye.
- 1936** Changes in the fluid pressure inside the eye.
- 1937** Increase of pressure of the intraocular fluid in the eye.
- 1938** The entrance of foreign matter into a blood vessel.
- 1940** Inflammation of the iris.
- 1942** Deficiency of blood in a part due to functional constriction or actual obstruction of a blood vessel.
- 1944** Inflammation of the cornea.
- 1945** Keratitis due to infection by *acanthamoeba*; it is usually associated with soft contact lens wear, particularly overnight wear.
- 1948** Failure of adequate output by the left ventricle despite an increase in distending pressure and in end-diastolic volume, with dyspnea, orthopnea, and other signs and symptoms of pulmonary congestion and edema.
- 1949** Enlargement or overgrowth of the myocardium of the left ventricle, due to chronic pressure overload.
- 1950** Any pathological or traumatic discontinuity of tissue or loss of function of a part.
- 1953** Bruising of the liver.
- 1956** A localized destruction or degeneration of the skin caused by various cutaneous diseases.
- 1957** Opportunistic infection of the retina by cytomegalovirus, seen in patients with immunodeficiency; symptoms include retinal necrosis and hemorrhage, leading to blindness.
- 1959** Deviations from the normal process; e.g. delayed, difficult, profuse, scanty, unusual bleeding, etc.
- 1960** A tiny cyst, frequently of such dimensions that a magnifying lens or microscope is required for observation.
- 1961** Small holes made through a part or a substance.
- 1962** Loss of the products of conception from the uterus before the fetus is viable; spontaneous abortion.
- 1963** Defective functioning of the mitral valve, with incomplete closure resulting in mitral regurgitation.
- 1964** The backward flow of blood from the left ventricle into the left atrium, owing to insufficiency of the mitral valve; it may be acute or chronic, usually due to mitral valve prolapse, rheumatic heart disease or a complication of cardiac dilatation.
- 1965** Narrowing of the left atrioventricular mitral orifice.
- 1966** A sudden, violent, involuntary contraction of a muscle or group of muscles.
- 1968** Stiffness or inflexibility.
- 1969** Gross necrosis of the myocardium, as a result of interruption of the blood supply to the area, as in coronary thrombosis.
- 1970** An unpleasant sensation; symptoms resulting from an inclination to vomit.
- 1971** Pathologic death of one or more cells, or a portion of tissue or organ, resulting from irreversible damage.
- 1978** New blood vessel formation in abnormal tissue or in abnormal positions.
- 1983** Functional disturbances or pathological changes in the peripheral nervous system.
- 1984** The act of closure or the state of being closed; an obstruction.
- 1992** A chronic inflammatory disease of the skeletal muscle, characterized by symmetrical weakness of the limb girdles, neck and pharynx, usually associated with pain and tenderness, and sometimes preceded or followed by manifestations typical of scleroderma, arthritis, systemic lupus, erythematosus or Sjorgen's Syndrome.
- 1997** Loss or impairment of motor function in a part due to lesion of the neural or muscular mechanism.
- 1998** A slight or incomplete paralysis.
- 1999** A peeling off or loss of epidermis, as in sunburn, postscarlatinal peeling, or toxic epidermal necrolysis.
- 2000** Any pelvic infection involving the upper female genital tract beyond the cervix.
- 2004** Inflammation of a vein.
- 2009** The assumption of various distinct forms by a single organism or species; also the property of crystallizing in 2 or more forms.
- 2010** The presence of fluid in the pleural space.
- 2011** Inflammation of the lungs with consolidation.

DEFINITIONS (continued)

- 2012** Accumulation of air or gas in the pleural space, which may occur spontaneously or as a result of trauma or a pathological process or to be induced deliberately.
- 2013** Erosion of patient's skin pocket which houses a device.
- 2014** Abnormally large eyeball.
- 2015** Antibodies directed against nuclear antigens; almost invariably found in systemic lupus erythematosus and are frequently found in rheumatoid arthritis, scleroderma, Sjogren's Syndrome and mixed connective tissue disease.
- 2020** Abnormal, diffuse, extravascular accumulation of fluid in the pulmonary tissues and air spaces due to changes in hydrostatic forces in the capillaries or to increased capillary permeability.
- 2021** Localized necrosis of lung tissue caused by obstruction of the arterial blood supply, most often due to pulmonary embolism.
- 2022** Defective functioning of the pulmonary valve, with incomplete closure resulting in pulmonic regurgitation.
- 2023** The backflow of blood from the pulmonary artery into the right ventricle, owing to insufficiency of the pulmonic semilunar valve.
- 2024** Narrowing of the opening between the pulmonary artery and the right ventricle, usually at the level of the valve leaflets.
- 2026** An obstruction of the pupil.
- 2027** A liquid inflammation product consisting of leukocytes and the debris of dead cells and tissue elements.
- 2028** Inducing fever.
- 2034** Intermittent bilateral attacks of ischemia of the fingers or toes and sometimes of the ears or nose, marked by severe pallor, and often accompanied by paresthesia and pain.
- 2035** A reaction that occurs at or about the site of infection or the point of injection.
- 2036** A reaction pertaining to or affecting the body as a whole.
- 2038** Lay term applied to any condition with dilation of conjunctival or ciliary blood vessels; innumerable causes, especially irritation and infection.
- 2039** Chronic, irreversible renal failure.
- 2041** The inability of a kidney to excrete metabolites at normal plasma levels under conditions of normal loading or the inability to retain electrolytes under conditions of normal intake.
- 2044** Cessation of breathing.
- 2045** Difficulty breathing.
- 2046** A condition of the newborn marked by dyspnea with cyanosis, most frequently occurring in premature infants, children of diabetic mothers and infants delivered by cesarean section, and sometimes with no predisposing cause.
- 2051** Moving backward; degenerating; reversing the normal order of growth and development.
- 2052** Moving backward or against the usual direction of flow.
- 2053** The most important manifestation of and sequel to rheumatic fever, consisting chiefly of valvular deformities.
- 2055** Failure of proper functioning of the right ventricle, with venous engorgement, hepatic enlargement, and subcutaneous edema.
- 2056** Enlargement or overgrowth of the myocardium of the right ventricle, due chronic pressure overload.
- 2058** A species comprising the yellow-pigmented, coagulase-positive pathogenic forms of the genus, causing serious suppurative infections and systemic disease; they produce toxins that cause food poisoning and toxic shock syndrome.
- 2059** The interval from the end of ventricular depolarization to the onset of the T wave; it is usually isoelectric in normal subjects.
- 2060** Formation of new tissue formed in the healing of a wound.
- 2062** Chronic hardening and thickening of the skin, which may be a finding in several different diseases, occurring in a localized or focal form and as a systemic disease.
- 2063** An attack or recurrence of a disease; a single episode of epilepsy.
- 2064** A simple partial seizure consisting of clonus or spasm of a muscle or muscle group; it may be single or in a continuous and repetitive series or may spread to adjacent muscles.
- 2066** Administration of antigen to induce a primary immune response; exposure to allergen that results in the development of hypersensitivity.
- 2067** The presence in the blood or other tissues of pathogenic microorganisms or their toxins.
- 2068** Shock associated with overwhelming infection, most commonly infection with gram-negative bacteria.
- 2069** A tumor-like collection of serum in the tissues.

DEFINITIONS (continued)

- 2071** Any contagious disease acquired during sexual contact; eg. syphilis, gonorrhea, chancroid.
- 2072** A sudden disturbance of mental equilibrium; a condition of profound hemodynamic and metabolic disturbance characterized by failure of the circulatory system to maintain adequate perfusion of vital organs.
- 2073** A symptom complex of unknown etiology, marked by the triad of keratoconjunctivitis sicca, xerostomia, and the presence of a connective tissue disease, usually rheumatoid arthritis but sometimes systemic lupus erythematosus, scleroderma or polymyositis.
- 2075** A gradual breakdown or very shallow ulceration of the skin which involves only the epidermis and heals without scarring.
- 2084** Chocking or arrest of respiration due to occlusion of the air passageway; arrest of circulation in a part due to compression.
- 2085** Decrease in the caliber of a canal, duct, or other passage, as a result of cicatricial contraction or the deposition of abnormal tissue.
- 2086** A sudden and severe attack.
- 2088** The stopping of respiration or the asphyxia that results from it.
- 2089** Decentration, or malposition of intraocular lens; zonular dehiscence or lens malposition caused by zonular or pars plana fixation.
- 2094** Inflammation of a synovial membrane.
- 2095** Rapid beating of the heart, usually applies to a heart rate above 100 beats per minute.
- 2097** Pertaining to tears involving the tentorium of the cerebellum.
- 2098** Children affected by toxins may be due to silicone breast implants and their effects on unborn children and from breastfeeding.
- 2100** Clotting within a blood vessel which may cause infarction of tissues supplied by the vessel.
- 2101** An aggregation of blood factors, primarily platelets and fibrin with entrapment of cellular elements, frequently causing vascular obstruction at the point of its formation.
- 2103** A noise in the ears, such as ringing, buzzing, roaring, clicking.
- 2107** "Fringe of pointed tips"; An atypical rapid ventricular tachycardia with periodic waxing and waning of amplitude of the QRS complexes on the electrogram as well as rotation of the complexes about the isoelectric line.
- 2108** A severe illness caused by infection with staphylococcus aureus and characterized by high fever of sudden onset, vomiting, diarrhea, and myalgia, followed by hypotension and in severe cases, shock; a sunburn-like rash with peeling of the skin, especially of the palms and soles, occurs during the acute phase.
- 2109** A brief attack (from a few minutes to an hour) of cerebral dysfunction of vascular origin, with no persistent neurological deficit.
- 2111** Incomplete closure of the tricuspid valve, resulting in tricuspid regurgitation; it is usually secondary to systolic overload in the right ventricle.
- 2112** The backflow of blood from the right ventricle into the right atrium, owing to imperfect functioning/insufficiency of the tricuspid valve.
- 2113** Narrowing or stricture of the tricuspid orifice of the heart.
- 2114** Dislodgement, breakdown, or other malfunction of an artificial cardiac pacemaker, chemotherapy port, drip infusion valve, or similar implanted diagnostic or therapeutic device as a result of unconscious or habitual manipulation by the patient.
- 2115** (Uveitis, Glaucoma, Hyphema) For complications that occur secondary to intraocular lens implantation in the anterior chamber; caused by rubbing of the lens loop (haptics).
- 2116** The formation or development of an ulcer.
- 2117** Denotes that a patient is resting and not under the influence of any medications.
- 2118** Denotes that the patient is resting and not under the influence of any sedatives.
- 2119** Accumulation of urine within the bladder because of the inability to urinate.
- 2122** An inflammation of part or all of the uvea, the middle (vascular) tunic of the eye, and commonly involving the other tunics (the sclera, cornea and retina).
- 2126** The diminution of the caliber of vessels, especially constriction of arterioles leading to decreased blood flow to a part.
- 2127** A state of increased caliber of the blood vessels.
- 2128** Spasm of the blood vessels, resulting in decrease in their caliber.
- 2129** Puncture of a vein, usually to draw blood or to inject a solution.

DEFINITIONS (continued)

- 2130** Arrhythmia characterized by fibrillary contractions of the ventricular muscle due to rapid repetitive excitation of myocardial fibers without coordinated contraction of the ventricle.
- 2131** A ventricular tachyarrhythmia characterized electrocardiographically by smooth undulating waves with QRS complexes merged with T waves, a rate of approximately 250 per minute.
- 2132** An abnormally rapid ventricular rhythm with aberrant ventricular excitation, usually in excess of 150 beats per minute.
- 2134** An illusory sense that either the environment or one's own body is revolving.
- 2150** Establishing a pouch of what was formerly an enclosed cyst.
- 2152** The large vessels entering the heart, including the aorta, the pulmonary arteries and veins, and the venae cavae.
- 2154** Thrusting or pushing out; expulsion by force.
- 2156** A deficiency of immune response or a disorder characterized by deficient immune response.
- 2158** Unequal refraction or focalization of light rays by a lens, resulting in degradation of the image they produce.
- 2159** A wrong or mistaken diagnosis.
- 2161** An involuntary, compulsive, repetitive stereotyped movement, usually of the face or shoulders.
- 2164** An objective visual sensation that appears with the eyes closed and in the absence of visual light.
- 2165** A nonimmunologic, chemically induced type of photosensitivity.
- 2168** A symptomatic form of epilepsy often preceded by an aura; characterized by loss of consciousness with generalized tonic-clonic seizures.
- 2169** Epilepsy characterized by absence seizures.
- 2171** A sensation as of repetitive pin pricks, caused by cold or by striking a nerve, or as a result of various diseases of the central or peripheral nervous system.
- 2172** The occurrence of a single contraction or a series of contractions of a muscle.
- 2174** To penetrate a tissue or substance.
- 2178** Associated with contact lens wear, due to an inadequate spreading of the tear film over the cornea as a result of incomplete and/or infrequent blinking.
- 2181** Inflammatory intraocular reaction with clouding and cells in vitreous; often accompanies inflammation of ciliary body, iris, choroid, or retina.
- 2186** Headache in erect position, after lumbar puncture; due to lowering of intracranial pressure by leakage of cerebrospinal fluid through the needle tract.
- 2187** A syndrome characterized by hyperbilirubinemia and deposition of bile pigment in the skin and mucous membranes with a resulting yellow appearance of the patient.
- 2188** An excess in the blood of urea, creatinine, and other nitrogenous end products of protein and amino acid metabolism.
- 2189** Any of a group of diseases that although clinically distinct and not necessarily related etiologically, have in common widespread pathologic changes in the connective tissue; they include lupus erythematosus, dermatomyositis, scleroderma, polyarteritis nodosa, thrombotic purpura, rheumatic fever, and rheumatoid arthritis.
- 2191** Shivering or shaking, accompanied by a sense of cold and pallor of the skin.
- 2192** A violent jar or shock, or the condition which results from such an injury.
- 2193** Painful, spasmodic muscular contractions.
- 2194** A sensation of unsteadiness with a feeling of movement within the head.
- 2195** Impairment of speech.
- 2196** Higher or lower than normal values for the serum electrolytes; usually affecting NA, K, CHL, CO₂, glucose, bun.
- 2197** A cytopathic retrovirus. It is the etiological agent of AIDS.
- 2198** Normal location of the mediastinum changes to the right or left, depending upon the underlying cause.
- 2201** The escape of blood around a heart valve, particularly around its leaflets.
- 2206** Inability of the heart to pump blood at an adequate rate to fill tissue metabolic requirements or the ability to do so only at an elevated filling pressure.
- 2207** The condition resulting from the spread of bacterial products (toxins) by the bloodstream.
- 2209** Asphyxia in utero due to hypoxia. See also fetal hypoxia.
- 2210** Hypoxia in utero, caused by conditions such as inadequate placental function (often abruptio placentae), preeclamptic toxicity, prolapse of the umbilical cord, or complications from anesthetic administration.
- 2211** Induced abortion.

DEFINITIONS (continued)

- 2212 All products of conception are expelled and identified.
- 2213 The uterus is not entirely emptied of its contents.
- 2214 Abortion brought on intentionally.
- 2215 Retention in uterus of an abortus that has been dead at least 4 weeks.
- 2216 Abortion induced to save the life or health of a pregnant woman.
- 2217 Inflammation of the amnion.
- 2218 Resembling anaphylaxis.
- 2220 Appearance of enlarged infected cells altering cell formation and possible loss of fetus.
- 2221 A violent involuntary contraction or series of contractions of the voluntary muscles; seizure.
- 2222 A convulsion marked by alternating contracting and relaxing of the muscles.
- 2223 Prolonged contraction of the muscles, as the result of an epileptic discharge.
- 2224 A severe, often fatal encephalopathy which has been attributed to accumulation in the brain of aluminum from dialysate prepared with inadequately purified water.
- 2225 An excretion or substance evacuated.
- 2226 Acute compression of the heart caused by increased intrapericardial pressure due to the collection of blood or fluid in the pericardium from rupture of the heart, penetrating trauma, or progressive effusion.
- 2227 Hazy ring around bright lights seen by some patients with refractive error or optical defects, e.g. cataracts, or corneal swelling.
- 2228 Diffuse edema of the central region of the cornea, usually associated with the wearing of hard contact lenses, but may also occur in keratoconus.
- 2231 Discrete, small lesions present in the cornea as a result of corneal inflammation and, in some cases, after soft contact lens wear especially extended-wear lenses.
- 2232 Very small, round vesicles containing fluid and cellular debris observed on the surface of the cornea under slit-lamp examination in some types of corneal dystrophy and in wearers of extended-wear lenses.
- 2233 The capability of the cornea to respond to stimulation.
- 2234 Ocular pain which may be very intense, accompanied by corneal epithelium damage, conjunctival injection, lacrimation, blepharospasm, photophobia and hazy vision following corneal edema caused by overwear of contact lenses, principally the PMMA type.
- 2236 Toxic substance from gram-negative bacteria that has a broad spectrum of biological activities, including pyrogenicity.
- 2237 A group of lung disorders characterized by obstruction of air flow, particularly on exhalation. The most common COPD disorders include chronic bronchitis, emphysema, asthma, and bronchiectasis.
- 2238 Pain in a muscle or muscles.
- 2239 A symptom complex characterized by nervousness, chest pain, back pain, palpations, pruritus, and other usually mild symptoms occurring minutes following the initiation of dialysis with a new dialyzer.
- 2240 A protrusion of a loop or knuckle or an organ or tissue through an abnormal opening.
- 2241 A smooth, slightly elevated area on the body surface, which is redder or paler than the surrounding skin. It is the typical lesion of urticaria, the dermal evidence of allergy, and in sensitive persons may be provoked by mechanical irritation of the skin.
- 2242 Excessive amount of sodium in the blood. A hazard associated with the use of softeners in hemodialysis.
- 2243 Abnormally decreased volume of circulating fluid (plasma) in body.
- 2244 Infection produced by direct contact with another person.
- 2245 Infection transmitted by water, food or other means of conveyance.
- 2246 An infection caused by pus-producing organisms.
- 2247 Infection associated with no detectable symptoms but caused by microorganisms capable of producing easily recognizable diseases, such as poliomyelitis or mumps.
- 2249 Pass from one part of the body or organ to another.
- 2250 Free from fever producing substances.
- 2252 Inflammation of the peritoneum.
- 2253 Polymyositis occurring in association with characteristic inflammatory skin changes.
- 2254 A physical response to the presence of endotoxins in the bloodstream, which is characterized by fever and occasionally chills and shaking rigors.
- 2255 A cutaneous inflammatory reaction occurring as a result of exposure to biologically effective levels of ionizing radiation.

DEFINITIONS (continued)

- 2256 Radiation induced skin desquamation.
- 2259 Flow in the opposite direction from normal, as the casting up of undigested food or gas from the stomach, or the backward flowing of blood into the heart, or between the chambers of the heart when a valve is incompetent.
- 2260 Partial seizure.
- 2261 The seizure seen in absence epilepsy, consisting of a sudden momentary break in consciousness of thought or activity, often accompanied by automatisms or clonic movements, especially of the eyelids.
- 2262 Shock resulting from primary failure of the heart in its pumping function, as in myocardial infarction, severe cardiomyopathy, or mechanical obstruction or compression of the heart.
- 2263 Narrowing or stricture of a duct or canal.
- 2264 A hypoglycemic reaction to overdosage of insulin, a skipped meal, or strenuous exercise in an insulin-dependent diabetic.
- 2265 Shock resulting from neurogenic vasodilation, which can be produced by cerebral trauma or hemorrhage, spinal cord injury, deep general or spinal anesthesia, or toxic central nervous system depression.
- 2266 A state of shock following a surgical operation.
- 2267 Shock that occurs during or after a surgical operation.
- 2268 Any shock produced by trauma.
- 2269 Acquired as an adult.
- 2270 Passed from mother to fetus.
- 2273 An agent or factor that induces or increases the incidence process that affects the production of a malformed fetus.
- 2274 A local defect, or excavation of the surface of an organ or tissue, which is produced by the sloughing of inflammatory necrotic tissue.
- 2275 Urination at short intervals without increase in daily volume of urinary output, due to reduced bladder capacity.
- 2278 A vascular reaction, usually transient, involving the upper dermis, representing localized edema caused by dilatation and increased permeability of the capillaries, and marked by the development of wheals.
- 2279 An anemia resulting from the destruction of red blood cells.
- 2285 Impedance is the resistance met by alternating electrical currents passing through a conductor.
- 2295 Process which involves the exchange of ions of sodium for the calcium and magnesium in the water supply, thus lowering to an acceptable range the otherwise high levels of these elements typically found in feed water.
- 2297 Failure to indicate the presence of bubbles in the extracorporeal circuit during cardiopulmonary bypass.
- 2300 The air eliminator vents and eliminates air from the cardioplegia delivery system during cardiopulmonary bypass.
- 2304 Failure of the stopcock valve to regulate or stop the flow of a fluid through a tubule structure.
- 2308 Process used in water purification systems in hemodialysis.
- 2314 Interference from sources of energy outside a system or systems, as contrasted to electromagnetic interference generated inside systems.
- 2315 System used to purify water in hemodialysis.
- 2317 Equipment becomes contaminated from blood/fluids coming from the patient, and must be discarded.
- 2318 Unintended collection/stasis of blood in device.
- 2327 Patient becomes entangled or caught in side rails of bed.
- 2336 A condition with sudden onset caused by acute vascular lesions of the brain, such as infarction from hemorrhage, embolism, thrombosis or rupturing aneurysm. See also: cerebrovascular accident (CVA), cerebral hemorrhage and subdural hemorrhage.
- 2352 Consisting of abnormally decreased sensitivity, particularly to touch.
- 2353 Consisting of abnormally decreased sensitivity, particularly to touch.
- 2354 Consisting of abnormally decreased sensitivity, particularly to touch.
- 2355 Pain in a joint.
- 2358 Removal of scar tissue.
- 2359 A vague feeling of bodily discomfort and fatigue.
- 2362 Paralysis of the nerves located in the legs.
- 2363 Difficulty moving; distortion or impairment of voluntary movement, as in tic, spasm, or myoclonus.
- 2364 A type of metabolic acidosis produced by accumulation of ketone bodies resulting from uncontrolled diabetes mellitus.

DEFINITIONS (continued)

- 2367** Inflammation of the pharynx.
- 2372** Excision of the apex of the petrous portion of the temporal bone.
- 2375** Inflammation of fascia.
- 2389** Inflammation of the meninges, usually by either a bacterium or a virus.
- 2390** A chronic adhesive arachnoiditis in the spinal arachnoid, with root and spinal cord symptoms similar to those caused by pressure from a tumor.
- 2393** The creation of an opening in the anterior trachea for insertion of a tube to relieve upper airway obstruction and to facilitate ventilation; a surgical creation of an opening into the trachea through the neck.
- 2430** A decrease in the fraction of the forced vital capacity that is exhaled in a specific number of seconds.
- 2431** An increase in the fraction of the forced vital capacity that is exhaled in a specific number of seconds.
- 2435** A decrease of the greatest rate of flow that can be achieved during forced expiration, beginning with the lungs fully inflated.
- 2436** An increase of the greatest rate of flow that can be achieved during forced expiration, beginning with the lungs fully inflated.
- 2479** The measurement of the mean left arterial pressure, as measured by a catheter introduced into the distal pulmonary artery, is low.
- 2480** The measurement of the mean left arterial pressure, as measured by a catheter introduced into the distal pulmonary artery, is high.
- 2481** The measurement of the mean left arterial pressure, as measured by a catheter introduced into the distal pulmonary artery, is normal.

PART II - DEVICE MANUFACTURER CODES
(For Section H, Blocks H3 and H6 of Form 3500A)

SUBPART A - DEVICE EVALUATED BY MANUFACTURER (Block H3)

CODES USED FOR A "NO" RESPONSE (ordered alphabetically)

02	Device evaluation anticipated, but not yet begun	01	Device received in a condition which made analysis impossible
03	Device not made by company	81	Other (code unspecified, describe in H10)
04	Device problem already known, no evaluation necessary		

SUBPART B - EVALUATION CODES (Block H6)

METHOD OF EVALUATION CODES (ordered alphabetically within each section)

SOURCE OF DEVICE EVALUATED

11	A device from same lot of the actual device involved in incident was evaluated	13	Device from controlled/non-released inventory evaluated
10	Actual device involved in incident was evaluated	12	Device from reserve sample evaluated

TYPE OF EVALUATION PERFORMED

36	Analysis of labeling performed	28	Mechanical tests, static - tension or compression failure
39	Chemical tests; e.g. corrosion, reaction	84	Optical tests of all specifications performed
21	Computer hardware performance tests conducted	86	Other (code unspecified, describe in H10)
22	Computer software performance tests conducted	83	Pathological evaluation of returned device (i.e. vascular grafts, heart valves, etc.)
20	Device evaluated with respect to operational environment	33	Performance tests of all specifications performed
25	Electrical tests of all specifications performed	31	Performance tests performed
23	Electrical tests performed	34	Performed test to determine if incident was the result of interaction with another device(s)
82	Environmental tests; temperature, humidity and vibration	37	Photographic images made during evaluation
30	Mechanical tests of all specifications performed	35	User-device interface test performed
26	Mechanical tests performed	38	Visual examination
27	Mechanical tests, dynamic - fatigue test		

RESULTS OF EVALUATION CODES (ordered alphabetically within each category)

CATEGORY A - DEVICE

703	Alarms inadequate or absent	704	Critical information not displayed
101	Component/subassembly failure (select specific code from Category D)	110	Design (unspecified)
103	Computer hardware problem	115	Design - maintenance difficult
104	Computer software problem	112	Design - not fail safe
105	Computer/human interface problem	114	Design - operational context
705	Controls, switches, keypads difficult to read or use	111	Design - packaging deficiencies
		113	Design - user interface
		706	Device difficult to assemble or set up

RESULTS OF EVALUATION CODES

CATEGORY A - DEVICE (continued)

708	Device performed according to specifications	174	Manufacturing - materials
702	Display hard to read or misleading	171	Manufacturing - packaging
120	Electrical problem (unspecified)	172	Manufacturing - process
122	Electrical problem - open circuit	173	Manufacturing - quality control
123	Electrical problem - operating outside specifications	175	Manufacturing - user interface
121	Electrical problem - short circuit	135	Material degradation/deterioration, anticipated or expected
197	Electro-magnetic interference problem	137	Material degradation/deterioration, unanticipated or unexpected
198	Electro-static interference	136	Material problem - post production
133	End of life - expected	180	Mechanical problem (i.e. pump, motor, wiring, cable, battery, parts, friction, etc.)
134	End of life - premature	710	Negative results of device testing
138	Exceeded expected life of device	709	None
140	Expected wear/deterioration	707	Operating steps confusing
141	Failure to cycle	100	Other (code unspecified, describe in H10)
142	Foreign material contamination	190	Out of specification
196	Inadequate optical transmission	160	Package insert information inadequate/incorrect
143	Inadequate quality assurance/control	167	Package insert missing
102	Incompatible components	162	Package labeling inadequate/incorrect
168	Instructions not available	163	Package mislabeled
144	Insulation degradation/deterioration	131	Premature power source depletion
154	Labeling contains inadequate instructions for use/maintenance	132	Premature RRT indicator (recommended replacement time)
157	Labeling contains incorrect instructions for use/maintenance	191	Quality assurance
151	Labeling difficult to read/see	164	Service manual (use in conjunction with another labeling code)
152	Labeling difficult to understand	192	Shelf life exceeded
153	Labeling does not provide necessary information	195	Sterilization
156	Labeling is incomplete	193	Storage/shipment
159	Labeling is not available	194	Telemetry failure
150	Labeling problems (unspecified)	176	Tolerance stack-up (hardware/software)
145	Lubrication excessive	165	User instruction manual (use in conjunction with another labeling code)
146	Lubrication insufficient	701	User/device interface
147	Lubrication, loss of		
148	Lubrication, wrong SAE (specifications)		
149	Lubrication, wrong type used - oil vs grease		
170	Manufacturing		

CATEGORY B - USE OF DEVICE

201	Cleaning error	209	Misapplication/misuse of device
213	Device used according to labeled indications	223	Modification of device - by authorized service organization
214	Device used in appropriate environment	225	Modification of device - by distributor
202	Device used with inappropriate material	224	Modification of device - by other service organization
215	Device used with incompatible equipment	221	Modification of device - by user/user facility biomedical engineering department
203	Device used with incompatible medium/material	220	Modification of device - unspecified
204	Disinfection error	236	Negative results of device testing
205	Failure to follow instructions	235	None
206	Failure to remove packaging per use instructions	200	Other (code unspecified, describe in H10)
207	Failure to service/maintain according to manufacturer recommendations	230	Reuse of device - unspecified
208	Incorrect technique/procedure	232	Reuse of device beyond labeled specifications

RESULTS OF EVALUATION CODES

CATEGORY B - USE OF DEVICE (continued)

234	Reuse of device without following disinfection/sterilization instructions	211	Timing error
231	Reuse of disposable device	212	Unapproved use of device
233	Reuse of single use device		

CATEGORY C - PHYSIOLOGICAL / PROCEDURAL FACTORS

311	Anticipated adverse reaction - long term	316	Patient's condition - predisposed event
312	Anticipated adverse reaction - short term	317	Patient's condition affected effectiveness of device
310	Anticipated or known (physiological/ procedural related, but code unspecified)	318	Patient's condition contraindicated use of device
321	Caused by another drug/device	322	Related to another drug/device
323	Caused by magnetic field - MRI	326	Related to operational context
324	Caused by metal - microwave	334	Support system problem (e.g., facility oxygen flow system, electrical Failure of entire facility)
325	Caused by operational context	340	Unanticipated (physiological / procedural related, but code unspecified)
331	Environmental factors	341	Unanticipated adverse reaction - long term
332	High/low temperature or humidity	342	Unanticipated adverse reaction - short term
313	Inherent risk of procedure	343	Unanticipated long term complication of procedure
314	Known long term complication of procedure	344	Unanticipated short term complication of procedure
315	Known short term complication of procedure		
346	Negative results of device testing		
345	None		
300	Other (code unspecified, describe in H10)		
333	Patient diagnosis contraindicated use of device		

CATEGORY D - DEVICE COMPONENT / SUBASSEMBLY FAILURES

401	Absorber (CO2)	597	Attenuator
548	Accumulator	598	Backpanel
402	Actuator	599	Backplane
431	Adapter	417	Bacterial filter
552	Air cleaner	730	Baffle
549	Air pump assembly	553	Bags
403	Alarm	418	Ball
404	Alarm, assembly	731	Ballast
405	Alarm, audible	419	Balloon
407	Alarm, high inspiratory pressure	732	Barrier
408	Alarm, LED	420	Battery
409	Alarm, low inspiratory pressure	733	Beam splitter
410	Alarm, oxygen pressure	734	Bearings
411	Alarm, power	735	Bell
412	Alarm, pressure	737	Belt
550	Alarm, self	738	Blower
413	Alarm, visual	739	Bobbin
414	Alarm, volume	740	Bolometer
590	Alternator	741	Bolt
591	Altimeter	554	Bottom dead center sensor bracket
592	Amplifier	742	Brake
593	Analyzer	422	Breathing circuit
594	Annunciator	743	Brush
595	Antenna	745	Bus
596	Arrester	746	Bushing

RESULTS OF EVALUATION CODES

CATEGORY D - DEVICE COMPONENT/SUBASSEMBLY FAILURES (continued)

736	Buzzer	440	Cylinder
423	Cable	441	Cylinder valve
747	Calibrator	778	Damper
748	Camera	779	Data acquisition unit
749	Cannula, inner	780	Decoder
750	Cannula, outer	442	Defibrillator paddles
424	Cap	443	Defibrillator subassembly (only use when part of another device)
751	Capacitor	781	Dehumidifier
555	Capacitor (ceramic chip)	782	Dehydrator
752	Carrier	783	Delay line
556	Cascade	784	Demagnetizer
425	Cassette	785	Demodulator
753	Caster	786	Densitometer
754	Cell	787	Detector
426	Charger	788	Device deployer
755	Chassis	789	Dial
756	Choke	445	Diaphragm
427	Circuit board	790	Digitizer
757	Clamp	447	Diode
758	Clip	448	Discrete component/device (resistor, capacitor, diode)
759	Clock	791	Disk
428	CO2 monitor	416	Display
429	CO2 monitor subassembly (only use when part of another device)	792	Divider
760	Coder	793	Dome
761	Coil	794	Doubler
762	Collimator	795	Driver
763	Comparator	796	Duct
557	Compressor (air pump)	797	Dummy load
558	Concentrator	798	Duplexer
559	Conductor coil	450	EKG/ECG monitor
560	Connecting rod	449	EKG/ECG subassembly (only use when part of another device)
435	Connector	452	Electrical lead
764	Contact lens	451	Electrode
436	Control switches	799	Emitter
765	Controller	800	Enclosure
766	Converter	801	Encoder
767	Cooling module	802	Endoscope
768	Cord	454	EPROM (erasable programmable read only memory)
769	Core	803	Equalizer
770	Counter	804	Evaporator
771	Coupler	455	Exhalation filter
772	Cover	805	Expander
437	CPU (central processing unit) board	806	Extender
561	Crank arm assembly	807	Extractor
438	CRT (cathode ray tube)	808	Eyelet
773	Crystal	456	Fail-safe system
562	Cuff	809	Fan
563	Cup seal	810	Fastener
774	Cups	811	Feedthru
775	Current limiter		
776	Current source		
439	Cusp		
777	Cutter		

RESULTS OF EVALUATION CODES

CATEGORY D - DEVICE COMPONENT/SUBASSEMBLY FAILURES (continued)

812	Ferrule	845	Inductor
813	FET (field effect transistor)	473	Insulation
814	Fiber	846	Integrator
815	Film	847	Intercom
816	Filter	569	Interconnection board intubate
457	Flange	848	Interface
817	Flasher	849	Interrupter
458	Flowmeter	850	Inverter
818	Foil	851	IOL (intraocular lens) implant
819	Foot pedal	852	Isolator
459	Foot switch	853	Jack
820	Frame	854	Joint
460	Function indicator	855	Joystick
564	Function lid	856	Jumper
821	Furnace	857	Junction
822	Fuse	475	Keyboard
461	Gas scavenging	858	Keyer
823	Gasket	859	Keypad
462	Gauges/meters	860	Kit
824	Gear	861	Knob
565	Gears	862	Label
825	Generator	863	Laminate
826	Grille	864	Lamp
827	Grommet	865	Laser
828	Ground strap	476	Lead conductor
829	Guide	477	LED (light emitting diode)
463	Guidewire	866	Lenses
830	Handset	478	Limit switch
464	Header	867	Limiter
831	Headphone	868	Line conditioner
832	Heads	570	Line cord
833	Headset	479	Lithium iodide battery
465	Heart valve cage	869	Load
466	Heart valve leaflet	870	Lock
834	Heart valve sewing ring	571	Locking sleeve
835	Heat exchanger	480	Logic board
836	Heat sink	572	Lower piston bearing
467	Heater	481	Luer valve
566	Heater rod	871	Magazine
567	Helifix coil	482	Magnet
837	Hinge	872	Magnetizer
838	Holder	873	Magnetometer
468	Hollow fiber	874	Magnetoresistor
839	Housing	875	Magnifier
469	Hub	483	Manifold
470	Humidifier	876	Matrices
471	Hybrid circuit	484	Membrane
474	IC (integrated circuit)	877	Memory
472	IC (integrated circuit) chip	878	Meter
840	Igniter	573	Meter needle
841	Imager	879	Microcircuit
842	Impedance converter	880	Microphone
843	Indexer	881	Microscope
844	Indicator	882	Mirror

RESULTS OF EVALUATION CODES

CATEGORY D - DEVICE COMPONENT/SUBASSEMBLY FAILURES (continued)

883	Mixer	580	Power lamp
884	Modem	913	Power module
885	Modulator	498	Power supply
886	Monitor	499	Power switch
485	Mother board	914	Preamplifier
486	Motor	915	Prescaler
574	Motor roll pin	916	Preselector
887	Mount	500	Pressure sensor
888	Multiplier	581	Pressure signal
889	Nebulizer (only use when part of another device)	501	Pressure tubing
890	Network	917	Printer
1023	None	918	Probe
891	Nuts	919	Processor
487	O2 monitor subassembly (only use when part of another device)	920	Programmer
892	O2 sensor	503	PROM (programmable read only memory)
488	Obturator	504	Prong
894	Optical cable	921	Protector
893	Optical disk	922	Proximity switch
546	Optical fiber	923	Pulley
895	Optocoupler	924	Pulser
896	Oscillator	925	Pump
897	Oscilloscope	926	Pyrometer
400	Other (code unspecified, describe in H10)	927	Rail
898	Outlet	505	RAM (random access memory)
899	Oven	928	Reactor
575	Overlay	929	Reader
489	Oximeter	930	Receiver
490	Oximeter (only use when part of another device)	507	Recorder (tape, stripchart, etc.)
491	Oxygen analyzer	931	Rectifier
900	Pad	582	Reed (electrically open)
901	Panel	932	Regulator
492	PC (printed circuit) board	508	Relay
902	Photodetector	933	Repeater
903	Photomultiplier	934	Reservoir
904	Phototransistor	935	Resistor
905	Pickup	936	Resolver
493	Pilot balloon valve	937	Resonator
906	Pin	938	Retainer
576	Piston guides	583	Retaining pin
577	Piston heads	939	Reticle
578	Piston Pads	940	RFI (radio-frequency interference)
494	Pivot	941	Rheostat
907	Plate	942	Ribbon
908	Plotter	943	Ring
579	Plug	944	Rivet
910	Pointer	945	Robot
911	Polarizer	946	Rod
495	Port	509	ROM (read only memory)
912	Post	947	Safety interlock
496	Potentiometer	984	SAW (surface acoustic wave device)
497	Power cord	949	Scale
		948	Scaler
		950	Scanner
		951	Scintillometer

RESULTS OF EVALUATION CODES

CATEGORY D - DEVICE COMPONENT/SUBASSEMBLY FAILURES (continued)

952	Scrambler	993	Teleprinter
953	Screen	994	Teletypewriter
568	Screw	995	Television unit
954	Scriber	996	Temperature compensator
432	Seal	997	Terminal
584	Seal tracheal	585	Terminal assembly
510	Sensor	998	Terminator
551	Servo	999	Thermistor
955	Shaft	1000	Thermocouple
956	Shield	547	Thermocouple wire
957	Shift register	1001	Thermometer
958	Shifter	1002	Thermostat
959	Shunt	1003	Thyristor
960	Shutter	1004	Timer
961	Signal conditioner	586	Tip conductor coil
962	Simulator	1005	Trackball
963	Slide	1006	Transceiver
964	Slip ring	522	Transducer
965	Slotter line	523	Transformer
967	Socket	524	Transistor
966	Socket adaptor	1007	Translator
511	Solder joint	1008	Transmission line
512	Solenoid	1009	Transmitter
968	Sound absorber	1010	Transponder
969	Spacer	1011	Transport
970	Spark gap	1012	TRIAC
971	Speaker	1013	Trimmer
972	Spectrometer	525	Tube
514	Spirometer	587	Tube capillary
973	Splicer	1014	Tuner
974	Spooler	444	User interface
975	Spring	527	Valve
976	Stabilizer	528	Valve (one-way, hemo, stopcock, etc.)
977	Stand	529	Valve, control
515	Stent	530	Valve, directional
516	Stepper motor	589	Valve, dump
978	Stiffener	531	Valve, exhalation
517	Stopcock	532	Valve, flap
979	Strain relief	533	Valve, flow
980	Strip line	534	Valve, inhalation
518	Stylet	535	Valve, inlet port
981	Substrate	536	Valve, inspiratory
982	Sunlamp	537	Valve, outlet port
983	Suppressor	538	Valve, PEEP (positive end expiratory pressure)
519	Switches	539	Valve, pressure limit
985	Synchronizer	540	Valve, relief
986	Synthesizer	541	Valve, safety
987	Syringe	542	Valve, selector
988	Table	543	Vaporizer
989	Tachometer	1015	Varistor
990	Tape	588	VDT (video display terminal)
991	Tee	1016	Vibrator
520	Telemetry equipment	1017	Viewer
992	Telephone	1018	Voltmeter

RESULTS OF EVALUATION CODES

CATEGORY D - DEVICE COMPONENT/SUBASSEMBLY FAILURES (continued)

1019	Washer	544	Wiring harness
1020	Waveguide	545	Y-piece connector
1021	Window	1022	Yoke
430	Wire		

CATEGORY E - COMPUTER-, IMAGING SYSTEM-, MICROPROCESSOR-CONTROLLED DEVICE PROBLEMS

617	Algorithm analysis	606	Magnetic interference
618	Algorithm problem	607	Medical diagnosis software
657	Bus conflicts	654	Memory access errors
633	Bus corruption	629	Metastable state
635	Calculation error (pentium)	608	Modem/communications failure
653	Cache controller	665	None
652	Cache memory	600	Other (code unspecified, describe in H10)
601	Central processing unit/operating system - error/design	609	Poor image resolution
643	Cold	659	Power fluctuation
628	Communication	610	Power loss/surge related complication
663	Compression errors	662	Power sequencing
603	Computer - user interface	661	Power spike
602	Computer hardware problem	660	Power transient
605	Computer software problem	640	Power-up errors
616	Computer software problem - error/design	641	Radiation effects
658	Conducted interference	622	Real-time operating system
626	Data corruption	639	Reset errors
631	Deadlock	611	RF interference
647	Dirt	612	Safety interlock/feature failure
655	DMA controller (direct memory access)	625	Sequence processing
646	Dust	650	Shock
664	Electromagnetic interference	636	Stack errors
648	Electrostatic discharge	619	State machine
630	Glitch	651	Temperature fluctuation
642	Heat	634	Timing violation
645	Humidity	620	Token
656	I/O (input/output) controller	632	Transmission effects
627	I/O (input/output) integrity	613	User installation error
637	Interrupt errors	614	User maintenance error
621	Kernel	624	Vector processing
638	Latch-up	649	Vibration
644	Liquid ingress	623	Watchdog

CONCLUSION CODES (ordered alphabetically)

40	Another device caused failure	64	Device failed during pre-test/pre-trial
70	Device discarded - unable to follow-up	42	Device failed just prior to use
71	Device evaluated and alleged failure could not be duplicated	43	Device failure directly caused event
72	Device evaluated and alleged failure could not be duplicated - cause of event unknown	44	Device failure directly contributed to event
41	Device failed during assembly	45	Device failure indirectly caused event
		46	Device failure indirectly contributed to event

CONCLUSION CODES(CONTINUED)

47	Device failure occurred and was related to event	74	No failure detected and product within specification
48	Device failure occurred but not related to event	75	No failure detected but product out of specification
49	Device failure related to maintenance	76	Operational context caused event
61	Device failure related to user handling	77	Operational context contributed to event
50	Device failure/lack of effectiveness related to patient condition	68	Other (code unspecified, describe in H10)
51	Device maintenance contributed to event	58	Software/firmware caused event
65	Device operated according to specifications	59	Software/firmware contributed to event
63	Device repaired and returned	88	This is a report of an accidental radiation occurrence (ARO), submitted pursuant to 21 CFR 1002.20
52	Device was out of calibration	66	Unusual event
53	Device was out of specification but this does not relate to event	79	User error caused event
54	Device was out of specification in a manner that relates to event	80	User error contributed to event
55	Intermittent failure directly caused event	62	User interface contributed to event
56	Intermittent failure directly contributed to event	60	User-interface caused event
57	Labeling related		
67	No conclusion can be drawn		
78	No device failure		