Instruction Sheet Owner's Certification of Lasers as PMA Approved Devices

Attached is a sample form for your certification that the laser you own is identical in all relevant aspects to a PMA approved medical device. This sample form contains information adequate for such a certification. Also attached is a blank form which you may use for your certification.

Please submit the completed certification form to FDA for a determination of the adequacy of the certification. FDA will notify you in writing either that the information is determined to be adequate or inadequate. If the information is inadequate, an investigational device exemption (IDE) application is required.

The adequacy of the certification does not depend on any proprietary information. All of the information is available in the approval letter, the operator's manual, the Summary of Safety and Effectiveness Data, and the laser owner's (or agent's) experience with the laser.

Please send your completed certification, signed and dated, to:

Division of Physical Sciences HFZ-130 Office of Science and Technology Center for Devices and Radiological Health 12721 Twinbrook Parkway Rockville, MD 20852

SAMPLE

Owner's Certification of Lasers as PMA Approved Devices [The information adequate for certification is in parentheses.]

	-	-		-	_	
1.	Name of manufact	urer:	Summit Tech	nology	VISX	
2.	Model name and n				_	
3.	Serial number(s): [Submission of this number is to ensure that the owner identifies a unique laser. This number is needed in case there are device failures, recalls or other actions required.]					
4.	Ablation zone si	ze(s): [≤6.5 mm]			
5.	Maximum energy per pulse: [180mJ/cm², Summit; 160mJ/cm², VISX]					
6.	Calibration techniques and schedule: [The applicant states that the calibration techniques and schedule are the same as described in the appropriate operator's manual for the Summit or VISX laser.]					
7.	Is the computer hardware, firmware, and software of your device identical to that distributed in the PMA approved device? [Yes. This device is enabled only for those indications and conditions for which the model received PMA approval. Please obtain certification from a qualified engineering company and submit their certification as an attachment to your certification. This certification specifically excludes any computer hardware, firmware or software used to obtain fees for use of the device, e.g., Pillar Point counters.]					
8.	Do you have a consumpropriate oper brochure if prover this information Act	rator's mar rided, and n is avai	nual, a pation other releva lable throu	ent information ant material?	n [Yes.	
the mak Section	make this Owner's Certificing or submission of false s 1001. Penalties include up ure of Owner (Certifie	to \$10,000 in fi	sent violations of U	nited States Code title rears imprisonment.	_	
Typed	Name of Owner (Cert	ifier)	Address of Ow	vner (Certifier)		

Owner	's Certification of Lasers as PMA Approved Devices					
1.	Name of manufacturer: Summit Technology VISX					
2.	Model name and number:					
3.	Serial number(s):					
4.	Ablation zone size(s):					
5.	Maximum energy per pulse:					
6.	Calibration techniques and schedule:					
7.	Is the computer hardware, firmware, and software of your device identical to that distributed in the PMA approved device?					
8.	Do you have a copy of the approved labeling, e.g., the appropriate operator's manual, a patient information brochure if provided, and other relevant material?					
full kn United	by make this Owner's Certification of Lasers as PMA Approved Devices to FDA with owledge that the making or submission of false statements represent violations of States Code title 18, Chapter 47, Section 1001. Penalties include up to \$10,000 in and up to five years imprisonment.					
	Signature of Owner (Certifier)					
	Typed Name of Owner (Certifier) Date					
	Address of Owner (Certifier)					