

TERANG NUSA Sdn Bhd

510(k) Summary for NUZONE Nitrile Surgical Glove Powderfree

MAY 1 5 2000

K000179

510(k) Summary

Submitter Name	Terang Nusa Sdn Bhd
Submitter Address	1, Jalan 8
	Pengkalan Chepa 2 Industrial Zone
	16100 Kota Bharu,
	Kelantan , Malaysia.
Submitter Telephone	+60 9 7735133
Submitter Fax	+60 9 7737755
Contact Person	LOW, Chin Guan
Date of preparation	15 Dec 1999
Trade Name	NUZONE
Common Name	Surgical Glove
Classification	Surgeon's Glove
Legally marketed device to which	The NUZONE nitrile surgical powderfree glove
substantial equivalence is being	described in this 510(k) is substantially equivalent
claimed.	to the Pure Advantage Powderfree Nitrile Surgical
	Gloves that is currently marketed.
Description of device	NUZONE meets the requirements for surgical
	gloves described by the American Standard for
	Testing and Material ASTM D3577. Type II



TERANG NUSA Sdn Bhd

510(k) Summary for NUZONE Nitrile Surgical Glove Powderfree

Intended Use of the device	NUZONE Powderfree Nitrile surgical gloves are disposable and sterile devices intended to be worn by healthcare personnel to prevent cross contamination between the user and the patient during procedures.
Brief description of non-clinical tests	Test conducted per ASTM D3577, ASTM D512 indicates that the product meet the requirements. Primary Skin Irritation test ASTM F 719-81 and Dermal Sensitization Test ASTM F 720-81 (86) indicates no sensitization or irritation. Final product is iodine tested for starch free status.
Brief description of clinical tests	Not required
Conclusion drawn from clinical and non clinical tests	It can be concluded that NUZONE Powderfree Nitrile Surgical Gloves will perform according to the performance standards referenced and therefore meets ASTM standards., FDA requirements and labelling claims. This device is substatially equivalent to the currently marketed devices.
Additional information deemed necessary by the FDA	None



MAY 1 5 2000

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Chin-Guan Low M. Director Terang NUSA Sdn. Bhd. 1 Jalan 8, Pengkalan Chepa 2 Industrial Zone, 16100 Kota Bharu, Kelantan, Malaysia

Re: K000179

Trade Name: Nuzone Nitrile Surgical Gloves Powder-Free,

Green

Regulatory Class: I Product Code: KGO Dated: April 18, 2000 Received: April 21 2000

Dear Mr. Low:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of

the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Timothy A. Ulatowski

Director

Division of Dental, Infection Control, and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure



TERANG NUSA Sdn Bhd

510(k) Submission for NUZONE Nitrile Surgical Glove Powderfree

3. Indication for use Statement

OR

Prescription Use

Per 21 CFR 801.109

Over the counter