K000295

FFR 1 6 2000

Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim)

**Ansell Perry** 

1875 Harsh Avenue SE Massillon, Ohio 44646 Telephone: 330-833-2811

Fax: 330-833-6213

[1] Summary

[2] Ansell Perry Inc.

1875 Harsh Avenue SE Massillon, Ohio 44646

Contact:

James R. Chatterton

Telephone: Fax:

330-833-2811 330-833-6213

January 28, 2000

[3] Trade Name: Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label

Claim)

Common Name:

Surgical Gloves Classification Name: Surgeon's Glove

Primary Skin Irritation in Rabbits

Guinea Pig Sensitization

- [4] Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim) meet all of the requirements of ASTM D 3577-99, Type 1.
- [5] Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim) meet all the current specifications for ASTM D 3577-99 Rubber Surgical Gloves.
- [6] Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim) are sterile disposable devices intended to be worn by operating room personnel to protect a surgical wound from contamination.
- [7] Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim) are summarized with the following technological characteristics compared to ASTM or equivalent standards.

| Characteristics     | Standard   |  |  |
|---------------------|--|--|--|
| Dimensions          | Meets ASTM D 3577-99   |  |  |
| Physical Properties | Meets ASTM D 3577-99, Type 1   |  |  |
| Freedom from holes  | Meets ASTM D 3577-99<br>Meets ASTM D 5151-92   |  |  |
| Protein Label Claim | This latex glove contains 50 micrograms or less of total water extractable protein per gram.  Meets ASTM D 5712-95 Standard Test Method for Analysis of Protein in Natural Rubber and Its Products |  |  |
| Biocompatability    |  |  |  |

**Passes** 

**Passes** 

Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim)
Ansell Perry
1875 Harsh Avenue SE
Massillon, Ohio 44646
Telephone: 330-833-2811

Fax: 330-833-6213

- [8] The performance test data of the non clinical tests are the same as mentioned immediately above.
- [9] Clinical data is not needed for medical gloves or for most devices cleared by the 510(k) process.
- [10] It is concluded that the Encore® MicrOptic® Powder Free Latex Surgical Gloves (Protein Label Claim) are as safe, as effective, and perform as well as the glove performance standards referenced above and therefore meet:

ASTM listed standards, FDA hole requirements, and labeling claims for the product.

[11] This summary will include any other information reasonably deemed necessary by the FDA.

## DEPARTMENT OF HEALTH & HUMAN SERVICES



FEB 1 6 2000

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. James R. Chatterton Vice President Regulatory Ansell Healthcare Products, Incorporated 1875 Harsh Avenue S.E. Massillon, Ohio 44646

Re: K000295

Trade Name: Encore MicrOptic Powder Free Latex

Surgical Gloves With Protein Content Labeling Claim (50

Micrograms or Less) Regulatory Class: I Product Code: KGO

Dated: January 28, 2000 Received: January 31, 2000

Dear Mr. Chatterton:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the <u>Code of Federal Regulations</u>, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in

the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Timothy A. Ulatowski

Ditector

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

## Attachment 2

## **Indications for Use Statement**

| 510(k)<br>Number<br>(if known)          | <del>_ к932665-</del> К   | 00020         | 15                   |             |
|---|---|---------------|----------------------|-------------|
| Device Name                             | Encore® Powder Free Surgical Gloves, MicrOptic® WITH PROTEIN CONTENT LABELING CLAIM (50 MICROGRAMS OR LESS)   |               |                      |             |
| Indications for Use                     | Encore® MicrOptic® Powder Free Latex Surgical Gloves intended use is to be to by operating room personnel to protect a surgical wound from contamination. |               |                      |             |
|   |   |               |                      |             |
|   |   |               |                      |             |
|   |   |               |                      |             |
| PLEASE DO NO                            | OT WRITE BELOW T  | HIS LINE - CO | ONTINUE ON ANOTH     | HER PAGE IF |
| Co                                      | oncurrence of CDRH  | Office of Dev | ice Evaluation (ODE) |             |
|   |   |               |                      |             |
| Prescription Use _<br>Per 21 CFR 801.10 | 9   | OR            | Over-The-Counter U   | se <u>X</u> |
| Bary A JW Cler<br>Division Sign-Off)    |   |               |                      |             |
| Division of Dental, Infection Co        | ontrol,   | 12            |                      |             |

and General Hospital Devices
510(4) Number <u>£ 000 2 95</u>