

K 000300



FEB 15 2000

GE Medical Systems

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

This 510(k) summary of Safety and Effectiveness information is submitted in accordance with the requirement of 21 CFR Part 807.87(h).

Submitter: Larry A. Kroger, Ph.D.
Senior Regulatory Programs Manager
TEL: 262-544-3894, FAX: 262-544-3863
Summary prepared 17 January, 2000

PRODUCT IDENTIFICATION

Name: LightSpeed 2.0 CT Scanner System
Classification Name: Computed Tomography X-ray System

Manufacturer: General Electric Medical Systems
16800 W. Ryerson Road
New Berlin, WI 53151

Distributor: Same as Manufacturer

Marketed Devices:

The LightSpeed 2.0 CT Scanner System is of comparable type and substantially equivalent to currently marketed Computed Tomography X-ray Systems that comply with the same or equivalent standards and have the same intended uses.

DEVICE DESCRIPTION

The LightSpeed 2.0 CT Scanner System is composed of a gantry, patient table, console, computer, and associated accessories.

Materials: Materials and construction are equivalent to the HiSpeed QX/i CT Scanner System (K980163) are compliant with UL 2601, IEC 601-1, and 21CFR Subchapter J.

Design: The System is designed to be a head and whole body CT scanner utilizing a solid state detector, and an intuitive Operator Console with the same tube and similar features to the HiSpeed QX/i CT Scanner System (K980163).

Indications for Use:

The LightSpeed 2.0 CT Scanner System is indicated for head and whole body x-ray computed tomography applications.

Comparison with Predicate:

It is the opinion of GE Medical Systems that the LightSpeed 2.0 CT Scanner System is of a type and substantially equivalent to currently marketed head and whole body ix-ray computed topography systems with respect to design, material composition, energy source, and radiation characteristics. It will comply with the x-ray requirements of 21 CFR as well as the safety requirements of UL2601, IEC 601-1 and collateral standards.

Adverse Effects on Health:

Potential electrical, mechanical and radiation hazards are identified in a hazard analysis and controlled by:

- System verification and validation to ensure performance to specifications, Federal Regulations, and user requirements.
- Adherence to industry and international standards. (UL and IEC).

CONCLUSIONS

The LightSpeed 2.0 CT Scanner Systems does not result in any new potential safety risks and performs as well as or better than devices currently on the market.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB 15 2000

GE Medical Systems, Inc.
C/O Reiner Krumme
TUV Rheinland of North America
12 Commerce Road
Newton, CT 06470

Re: K000300
LightSpeed 2.0 CT Scanner System
Dated: January 27, 2000
Received: January 31, 2000
Regulatory class: II
21 CFR 892.1750/Procode: 90 JAK

Dear Mr. Krumme:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,


Capt. Daniel G. Schultz, M.D.
Acting Director, Division of Reproductive,
Abdominal, Ear, Nose and Throat,
and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K000300

Device Name: LightSpeed 2.0CT Scanner System

Indications For Use:

The LightSpeed 2.0CT Scanner System is indicated for head and whole body x-ray computed tomography applications.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

David A. Syron
(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K000300