

MAR 3 0 2000

K000588

II. 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

Submitter

Name:.....ESPE Dental AG  
Street:.....ESPE Platz  
ZIP-Code, City:.....D-82229 Seefeld  
Federal State:.....Bavaria  
Country:.....Germany  
Establishment Registration Number:.....9611385  
Contact:.....Dr. Andreas Petermann, Manager U.S.  
Regulatory Affairs  
Phone:.....011-49-8152-7001395  
Fax:.....011-49-8152-7001869  
E-mail.....Andreas\_Petermann@ESPE.de  
Date:.....February 16, 2000

Name of Devices

Proprietary Name: .....DIMENSION® GARANT® L  
DIMENSION® GARANT® L QUICK  
Classification Name: .....Impression material  
Common Name: .....Polyvinyl siloxane based impression material

Predicate Devices

DIMENSION® PENTA®  
DIMENSION® GARANT® L by ESPE.....K 960547

### Description for the Premarket Notification

DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are classified as impression materials (21 C.F.R. § 872.3660) because they are devices intended to reproduce the structure of a patient's teeth.

DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are low consistency impression materials designed to be used as wash materials in dual phase impression techniques. Therefore, they are similar in intended use and substantially equivalent to ESPE's polyvinylsiloxane based impression material DIMENSION® GARANT® L (old version).

DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are two-component materials (base paste and catalyst) packaged in cartridges which are intended to be mixed in ESPE's mixing, dosing and dispensing system GARANT®.

In recent years ESPE was already marketing low consistency impression materials tradenamed DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK in the U.S.A. The marketing of this material was based on the 510(k) for DIMENSION® PENTA®/DIMENSION® GARANT® L (K 960547). The development of a fast setting low consistency material in (DIMENSION® GARANT® L QUICK) was considered to not requiring the submission of a new 510(k) premarket notification because the quantitative chemical composition was changed only slightly and the indications for use were not affected.

However, in this 510(k) premarket notification submission the latest product versions DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are described in terms of chemical composition and physical and mechanical properties. Furthermore, the range of indications will no longer be limited to impressions for inlay, onlay, crown and bridge preparations and impressions of edentulous jaws. The range of indications will be expanded to allow all types of dual phase impressions, e.g. also for orthodontics.

DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK have the following similarities to the 510(k)ed DIMENSION® GARANT® L:

- DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are used by the same operating principle

- DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK incorporate the same basic chemical design
- DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK have the same shelf life
- DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are manufactured and packaged using the same materials and processes

The physical and mechanical properties of the new DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK have been compared to those of the already 510(k)ed DIMENSION® GARANT® L.

The compositions of DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK contain the same ingredients as the 510(k)ed material, variations only occur in quantitative terms. Therefore, additional biocompatibility testing is not necessary in our point of view.

ESPE's impression materials DIMENSION® PENTA® and DIMENSION® GARANT® L are well established and considered to be safe and effective. Comparison of chemical constitution and physical and mechanical properties show that the new materials DIMENSION® GARANT® L and DIMENSION® GARANT® L QUICK are substantially equivalent to the well established material.



MAR 3 0 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Dr. Andreas Petermann  
Manager U.S. Regulatory Affairs  
ESPE Dental AG  
ESPE Platz  
D-82229 Seefeld, Bavaria  
GERMANY

Re: K000588  
Trade Name: Dimension® Garant® L, Dimension® Garant® L  
Quick  
Regulatory Class: II  
Product Code: ELW  
Dated: February 16, 2000  
Received: February 22, 2000

Dear Dr. Petermann:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any

obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4690. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director

Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

III. STATEMENT OF INDICATIONS FOR USEDevice Name:

DIMENSION® GARANT® L

DIMENSION® GARANT® L QUICK

Indications for use:

Dental impression material for mixing and dispensing in a GARANT® mixing device:

Wash material for all kinds of dual phase impression techniques



(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K000588