



MAY 26 2000

GE Medical Systems

page 1 of 2

PO Box 414, W-709
Milwaukee, WI 53201
USA**SUMMARY OF SAFETY AND EFFECTIVENESS**

- This 510(k) summary of safety and effectiveness information is submitted in accordance with the requirements of 21 CFR Part 807.87(h).

- Identification of Submitter

Larry A. Kroger, Ph.D., 414-544-3894, March 10, 2000

- Identification of the Product

Spiral Imaging Option

Manufactured by: GE Medical Systems
3200 N. Grandview Blvd.
Waukesha, WI 53188

- Device Description

Spiral Imaging is a data acquisition method in which K-space is filled in a spiral fashion, as opposed to being filled in a uniform rectilinear grid, as in other imaging techniques.

- Indications for Use

The Spiral Imaging Option is intended for whole body use and is capable of producing images with high spatial and temporal resolution. This acquisition sequence can be useful for imaging fine structures in organs where motion is a problem, such as coronary arteries in the heart.

- Comparison with Predicate

It is the opinion of GE medical Systems that the Spiral Imaging Option is substantially equivalent to the Fast Gradient Echo image acquisition sequence in the Signa CV/i MRI System (K980114). Waveform changes permit more data to be acquired with each excitation without generating artifacts associated with extended rectilinear acquisition sequences. This in turn shortens the overall time of the acquisition.

**GE Medical Systems**

page 2 of 2

PO Box 414, W-709
Milwaukee, WI 53201
USA

SUMMARY OF SAFETY AND EFFECTIVENESS

- **Summary of Studies**

The Spiral Imaging Option was evaluated to the IEC 601-2-33 International medical equipment safety standard for Magnetic Resonance Systems. Evaluation testing was done to verify the performance of the option as well as to verify the present dB/dt limits of the Signa CV/i are maintained.

- **Conclusions**

It is the opinion of GE that the Spiral Imaging Option does not result in any new potential hazards.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 26 2000

Larry A. Kroger Ph.D.
Senior Regulatory Programs Manager
GE Medical Systems
P.O. Box 414, W-709
Milwaukee, WI 53201

Re: K000802
Spiral Imaging Option for MRI
Dated: March 10, 2000
Received: March 13, 2000
Regulatory Class: II
21 CFR §892.1000/Procode: 90 LNH

Dear Dr. Kroger:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure(s)

510(k) Number (if known): K000802

Device Name: Spiral Imaging Option

Indications For Use:

The Spiral Imaging Option is intended for whole body use, and is capable of producing images with high spatial and temporal resolution. This acquisition sequence can be useful for imaging fine structures in organs where motion is a problem, such as coronary arteries in the heart.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
 (Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K000802