



JUN - 9 2000

**Corporate**

595 Commerce Drive  
Buffalo, NY 14228  
(800) 343-2324  
(716) 835-7000  
(716) 835-3414 fax  
www.buffalofilter.com  
email to:  
info@buffalofilter.com

**West Coast**

1235 Francis Avenue  
Suite F  
Ontario, CA 91761  
(888) 311-6881  
(909) 923-2346  
(909) 923-2114 fax

**European Service  
And Repair Center**

Gruvgatan 37  
S-421 30 Vastra  
Frolunda, Goteborg  
Sweden  
+46-31-498165  
+46-31-476016 fax

**World Leader in  
Surgical Smoke  
Evacuation™**

K000904

**510(k) Summary  
PenAdapt 10**

**Manufacturer**

Medtek Devices, Inc.  
(d/b/a Buffalo Filter)  
595 Commerce Drive  
Buffalo, NY 14228

**Manufacturing Location**

Same

**Contract Sterilizer**

For Gamma/E-beam Radiation Sterilization:  
SteriGenics International  
305 Enterprise Drive  
Westerville, OH 43081

For Ethylene Oxide Sterilization:  
Not yet determined

**Telephone/FAX/Email**

Tel (716)-835-7000  
Fax: (716)-835-3414  
Email: kathy.selover@buffalofilter.com

**Contact Person**

Kathleen H. Selover

**Device Trade Name**

PenAdapt 10

**Common Name**

Smoke Aspiration Tip

**Classification Name**

Accessory to surgical exhaust apparatus

**ProCode**

79FYD

**Predicate Device**

SE-100 Smoke Aspiration Tip  
manufactured by Implants Associates.  
Originally submitted by Inamed under  
510(k) K931548.

KE00904

<b>Description</b>	An elastomer sheath with a built-in air passage designed to fit co-axially over an electrosurgical pencil or pencil-like device. The distal end is connected via surgical tubing and/or hoses to a vacuum source. Supplied both as a stand- alone device and configured with a variety of lengths of surgical tubing/hoses.
<b>Intended Use</b>	PENADAPT 10 is intended for removing smoke, particles and body fluids from the point of surgical activity during medical procedures that use an electrosurgical (ESU) pencil or pencil-like device for cutting and/or cauterizing.
<b>Physical/Technical Comparison</b>	Both the device and the predicate are made from a silicone elastomer, are of the same general size and shape and have the same intended use. Both devices require attachment to a vacuum or suction apparatus in order to aspirate smoke and fluid from the surgical site.
<b>Performance Summary</b>	PenApapt 10 is capable of maintaining adequate air flow to allow for adequate aspiration of smoke and body fluids found at the site of surgical activity.
<b>Biocompatibility Testing</b>	The material used for the manufacture of PenAdapt 10 has undergone biocompatibility testing and has been found to be non-irritating, non-cytotoxic and non-sensitizing.
<b>Sterility and Shelf Life</b>	The device is sterile processed using gamma (or e-beam) radiation or ethylene oxide. Method of sterilization is dependent upon the device's final configuration.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN - 9 2000

Ms. Kathleen H. Selover  
Regulatory Affairs Manager  
Buffalo Filter, Company, Incorporated  
595 Commerce Drive  
Buffalo, New York 14228

Re: K000904  
Trade Name: PenAdapt 10  
Regulatory Class: II  
Product Code: FYD  
Dated: May 11, 2000  
Received: May 16, 2000

Dear Ms. Selover:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Selover

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number, if known: K000904

Device Name: Medtek Devices, Inc. (d/b/a Buffalo Filter) PenAdapt 10

Indications for Use:

**Indications for Use:**

This device is an aspiration sheath that fits over an electrosurgical pencil body and leaves the tip or blade exposed. This device is an accessory to an electrosurgical unit (ESU).


The device is used to remove smoke, particles and body and other casual fluids from the point of surgical activity during surgical procedures that use ESU for cutting and cauterizing. This device is used in conjunction with a suction (vacuum) source.

**Contraindications**

This device should not be used for microsurgery.

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices

510(k) Number K000904

Prescription Use \_\_\_\_\_ OR Over-the-Counter-Use ☒ (Per 801.109)