

JUN - 8 2000

K000931 P. 1 of 3

510(k)
AU5 Ultrasound Imaging System
Biosound Esaote

510(k) Summary

The following safety and effectiveness summary has been prepared pursuant to requirement for 510(k) summaries specified in 21CFR§807.92(a).

807.92(a)(1)

Submitter Information

Colleen Hittle, Official Correspondent
8000 Castleway Drive
Indianapolis, IN 46250
Phone: (317) 849-1916
Facsimile: (317) 577-9070

Contact Person: Colleen Hittle

Date: March 15, 2000

807.92(a)(2)

Trade Name: AU5 Ultrasound Imaging System
(Addition of 3D Imaging Mode)

Common Name: Ultrasound Imaging System

Classification Name(s): Ultrasonic pulsed doppler imaging system 892.1550

Classification Number: 90IYN
90IYO

807.92(a)(3)

Predicate Device(s)

Esaote	AU5	K980468
--------	-----	---------

Additional Substantial Equivalence Information is provided in the following substantial Equivalence Comparison Table.

510(k)
AU5 Ultrasound Imaging System
Biosound Esaote

807.92(a)(5)

Intended Use(s)

The AU5 ultrasound imaging system is intended to be used by a physician for diagnostic imaging in cardiac, abdominal, peripheral vessel and fetal applications.

510(k) Summary
 AU5 Ultrasound Imaging System
 Biosound Esaote

Comparison Chart for Substantial Equivalence

General Characteristics	<u>Esaote</u>	<u>Esaote</u>
	<u>AU5</u>	<u>AU5</u>
	K980468 (predicate)	(New mode)
Transducer Type	Annular Array	
	Mechanical Sector	
	Linear	Linear
	Convex	Convex
	Phased Array	
2D Freq MHz	2.5/15	3.5/10
PW Freq MHz	2.25/10	2.3/7
CW Freq MHz	2.25/5.0	
Imaging Modes	Real-time/2D	Real-time/2D
	M Mode	M Mode
	PW Doppler	PW Doppler
	CW Doppler	
	CFM Doppler	CFM Doppler
	Power Doppler	Power Doppler
		3D
Probes MHz		
Annular array	2.5, 3.5, 7.5, 10, 13	
Mechanical Sector	13	
Linear	3.5-7.5	7.5-10
Convex	3.5-5.0	3.5-7.5
Multifrequency probes	Yes	Yes
Special probes	IVT transvaginal	IVT transvaginal
	TRT transrectal	
Biopsy attachments	Linear Array	
	Convex	
Monitor size (inches)	14	12
Programmability	6 presets	6 presets
Pulsed/CW Doppler	Yes	Yes
HIPRF	No	No
2D Updating	Yes	Yes
CW steerable	Yes	
Audio stereo	Yes	Yes
Color doppler upgrade	Yes	Yes
ECG	Option	Option
Computer interface	Centronics output	Centronics output
DSM (Dicom storage module)	Yes	No



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN - 8 2000

Colleen Hittle,
Official Correspondent
Biosound Esaote, Inc.
8000 Castleway Drive
Indianapolis, IN 46250

Re: K000931
AU5 Ultrasound Imaging Systems with 3-D Imaging
Dated: March 15, 2000
Received: March 22, 2000
Regulatory Class: II
21 CFR 892.1550/Procode: 90 IYN
21 CFR 892.1560/Procode: 90 IYO
21 CFR 892.1570/Procode: 90 ITX

Dear Ms. Hittle:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the AU5 Ultrasound Imaging Systems, as described in your premarket notification:

Transducer Model Number

CA11
LA13A
IVT22


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic QS inspections, the FDA will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the Federal Register. *Please note:* this response to your premarket notification does not affect any obligation you may have under sections 531 and 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Paul M. Gammell, Ph.D. at (301) 594-1212.

Sincerely yours,

for 
Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosures

AU5 System

Appendix F

Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:


Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler CFM	Amplitude Doppler PD	Color Velocity Imaging	Combined (specify)	Other (specify) 3D
Ophthalmic										
Fetal		P	P	P		P	P		See additional comments	N
Abdominal		P	P	P		P	P		See additional comments	N
Intraoperative (specify) Abdominal		P	P	P		P	P		See additional comments	N
Intraoperative										
Pediatric		P	P	P		P	P		See additional comments	N
Small Organ (specify)		P	P	P		P	P		See additional comments	N
Neonatal Cephalic		P	P	P		P	P		See additional comments	N
Adult Cephalic		P	P	P	P	P	P		See additional comments	N
Cardiac		P	P	P	P	P	P		See additional comments	N
Traneseophageal										
Transrectal		P	P	P		P	P		See additional comments	N
Transvaginal		P	P	P		P	P		See additional comments	N
Transurethral										
Intravascular										
Peripheral Vascular		P	P	P	P	P	P		See additional comments	N
Laparoscopic		P	P	P		P	P		See additional comments	N
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication: P= previously cleared by FDA; E= added under Appendix E


Additional Comments: Small Organs (specifically: thyroid, testicles, and breast); PV to include vein mappingApplicable combined modes: B+PW+CFM+M+PD

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices510(k) Number K000931Prescription Use 
(Per 21 CFR 801.109)

Diagnostic Ultrasound Indications for Use Form


Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler CFM	Amplitude Doppler PD	Color Velocity Imaging	Combined (specify)	Other 3D
Ophthalmic										
Fetal		P	P	P		P	P		See comments	N
Abdominal		P	P	P		P	P		See comments	N
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		E	E	E		E	E		See comments	N
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		P	P	P		P	P		See comments	N
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: TEI=3D Imaging Mode; Applicable combined modes: B+PW+CFM+M+PD


 (Division Sign-Off) **DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)**
 Division of Reproductive, Abdominal, ENT,
 and Radiological Devices
 510(k) Number K000931

Concurrence of CDRH, Office of Device Evaluation (ODE)

 Prescription Use ☒
 (Per 21 CFR 801.109)

Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler CFM	Amplitude Doppler PD	Color Velocity Imaging	Combined (specify)	Other (specify) 3D
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Peripheral vascular										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		See comments	N
Small Organ (specify)		P	P	P		P	P		See comments	N
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Tranosophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		P	P	P		P	P		See comments	N
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Small Organs (specifically, thyroid, testicles, and breast); Peripheral Vascular to include Vein Mapping & Sclerotherapy; TEI=3D Imaging Mode

Applicable combined modes: B+PW+CFM+M+PD

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

(Division Sign-Off)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Division of Reproductive, Abdominal, ENT,
and Radiological DevicesPrescription Use ☒
(Per 21 CFR 801.109)

510(k) Number

K000931

Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler CFM	Amplitude Doppler PD	Color Velocity Imaging	Combined (specify)	Other (specify) 3D
Ophthalmic										
Fetal		E	E	E		E	E		See comments	N
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Traneseophageal										
Transrectal		E	E	E		E	E		See comments	N
Transvaginal		E	E	E		E	E		See comments	N
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: TEI=3D Imaging Mode; Applicable combined modes: B+PW+CFM+M+PD

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K000931

Prescription Use ☒
(Per 21 CFR 801.109)