

MAY 26 2000

K 000960 P142

510(K) SUMMARY
(as required by 807.92(c))

Submitter of 510(k):	AJW Technology Consultants, Inc. 962 Allegro Lane Apollo Beach, FL 33572 Phone: 813-645-2855 Fax: 813-645-2856
Contact Person:	Art Ward
Date of Summary:	March 9, 2000
Trade Name:	Barzell-Whitmore Maroon Bells, Inc. 1121 Lewis Avenue Sarasota, FL 34237
Classification Name:	Image Processing System
Predicate Device:	Life Imaging Systems LIS 8030 A K961403
Device Description/ Comparison:	Colo-Rectal Mover

ABBREVIATED COMPARATIVE STUDY for COLO-RECTAL MOVER

**Life Imaging Systems
Translational Mover
LIS 8030**

**Barzell-Whitmore
Colo-Rectal Mover**

Indications for Use	Linear Endocavity Probe Positioner	Same
Adjustment Capabilities	Longitudinal Movement of at least 60mm	Longitudinal 5mm Movement of at least 160mm
Used with	Trans-Rectal Ultrasound System	Same
Energy	6V DC	6V DC
Sterility	Non-Sterile	Non-Sterile
Intended Use:	Intended to move a standard endocavity probe while it is positioned in its normal imaging position within a patient's body.	



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 26 2000

Barzell-Whitmore Maroon Bells, Inc.
c/o Mr. Arthur Ward, Consultant
AJW Technology Consultants, Inc.
962 Allegro Lane
Apollo Beach, FL 33572

Re: K000960
Colo-Rectal Mover (pullback device for 3D imaging)
Dated: March 9, 2000
Received: March 24, 2000
Regulatory Class: II
21 CFR §892.1560/Procode: 90 IYO

Dear Mr. Ward:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure(s)

510(k) Number (if known): K000960

Device Name: BARZELL-WHITMORE BELLS INC. Colo-Rectal Mover

Indications For Use:

The BARZELL-WHITMORE BELLS INC. Colo-Rectal Mover is intended to move a standard ultrasound endocavity probe while it is positioned in its normal imaging position within a patient's body.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

David A. Segman
(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K000960

[Optional Format 1-2-96]