

JUN - 8 2000

Page -1- of -2

510(k) SUMMARY

This summary of 510(k) safety and effectiveness is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92

The assigned 510(k) number is: K 000 989

Submitted by: InVitroCare, Inc.
11408 Sorrento Valley Rd.
Suite 202
San Diego, CA 92121
Telephone: (858) 452-1986
Facsimile: (858) 452-1828
E-mail: Invitroc@aol.com
Contact: Robert E. Lovins, PhD
Date Submitted: October 13, 1999

Device Identification:

Trade Name: SPERM WASH Medium
Common Name: Sperm processing and washing medium,
Human Tubal Fluid medium
Classification Name: Reproductive Media (21CFR, 886.6180)

Predicate Device:

Notice of Final Rule, 63 FR 48428, Docket number 97N-0335 and
510(k) Reference Number K874668

Description:

SPERM WASH Medium is a synthetic, defined culture medium supplemented with 5mg/ml human serum albumin and is intended for use in assisted reproductive technology procedures. It has been formulated to mimic the composition of the fluid found in fallopian tubes as defined by Quinn et al (Quinn P, Kerin JF, Warnes GM: Fertil Steril 1985;44:493-498). SPERM WASH Medium uses a combined sodium bicarbonate/HEPES ([4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid]) buffering system and is appropriate for those procedures that do not use a carbon dioxide atmosphere.

SPERM WASH Medium is intended for use in assisted reproductive technology procedures that involve the manipulation of sperm. Specifically, SPERM WASH Medium is intended for use as a sperm-processing medium in washing procedures and for transport of sperm during reproductive technology procedures.

Design Characteristics:

SPERM WASH Medium is intended for use as a culture medium for a variety of assisted reproductive procedures. It has been used for a number of years as a sperm-washing medium. During sperm wash procedures, viable sperm cells are separated from the other constituents of seminal fluid in an effort to concentrate the viable sperm and increase the number of sperm available for fertilization. A culture medium such as SPERM WASH Medium is used to suspend the semen, the sample is centrifuged to pellet the viable sperm and after the supernatant is decanted, the pellet is resuspended in fresh medium. After a brief incubation period during which the motile sperm "swim up" into the fresh medium, the sperm are aspirated and used for the desired fertilization procedure. SPERM WASH Medium is therefore intended for use as a sperm-washing and transport medium in assisted reproductive technology procedures.

Performance Data:

SPERM WASH Medium is subjected to cytotoxicity testing and sperm motility/hyperactivation analysis. Each lot of SPERM WASH Medium is also assayed by a mouse embryo assay prior to its release to market. These assays assure that the product is both functional for its intended use, and that no toxic components are present in the formulation. Sperm washing media have been used in a variety of clinical settings, for the intended use for a number of years. In that time the product has become the standard media used for sperm isolation, washing and transport.

Additional Information:

Mouse embryo testing will be performed as a condition of release for SPERM WASH Medium as well as endotoxin and sterility testing. Results of all release assays will be reported on a lot-specific certificate of analysis and will be indicated on the labeling.

Conclusion:

The conclusion from performance testing, as well as a review of published historical information contained in the professional literature shows that SPERM WASH Medium is suitable for its intended use and meets the criteria outlined in the Final Rule, 63 FR48428, Docket number 97N-0335.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN - 8 2000

Robert E. Lovins, Ph.D.
President
In VitroCare, Inc.
11408 Sorrento Valley Road
Suite 202
San Diego, CA 92121

Re: K000989
SPERM WASH Medium
Dated: March 27, 2000
Received: March 28, 2000
Regulatory Class: II
21 CFR §884.6180/Procode: 85 MQL

Dear Dr. Lovins:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure(s)

InVitroCare, Inc.

March 22, 2000

Page -10-

INDICATIONS FOR USE STATEMENT (Page 1 of 1)510(k) number: K000989

Device Names: SPERM WASH Medium

Indications for Use:

SPERM WASH Medium is intended for use in assisted reproductive technology procedures that involve the manipulation of sperm. Specifically, SPERM WASH Medium is intended for use as a sperm processing medium in washing procedures and for transport of sperm during reproductive technology procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use
(per 21 CFR 801.109)

David A. Segura
(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K000989