

AUG 7 2000

510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of the SDMA 1990 and 21 CFR Part 807.92.

The assigned 510(k) number is: K 001126

1. Date of summary: March 31, 2000
2. Submitted by: aLATEX Scientific Inc.
6355 Westheimer Rd.
Suite 105
Houston, Texas 77057
1-888-513-TEST
Contact: Mr. Stephen Gilson CEO
3. Device Name: Peace of Mind Multiple Drugs of Abuse Test.
4. Device Classification: CFR 21 862.3870, 862.3650, 862.3250, 8610.3610, Class II, Panel 91 Toxicology
5. Device description: The Peace of Mind Multiple Drugs of Abuse Test is an immunochromatographic based one step *in vitro* test.
6. Intended Use: The Peace of Mind Multiple Drugs of Abuse Test is a qualitative, one step immunochromatographic competitive assay used to screen human urine for the presence of the five most commonly abused drugs at cut off concentrations of 50 ng/mL THC, 300ng/mL benzoylecgonine/cocaine, 2,000ng/mL opiates, 1,000ng/mL methamphetamine and 25ng/mLPCP. The test is qualitative and provides only a preliminary analytical result, which must be confirmed by an alternate methodology preferably, GC/MS.
7. Substantial Equivalence: The Peace of Mind Multiple Drugs of Abuse Test was found substantially equivalent to the Instacheck Drug Screen Panel. Both products are immunoassays and use specific antibodies to detect various drug compounds. Both assays are preliminary screens for human urine and require confirmation with alternate methods such as GC/MS. The sensitivity of the tests is similar. The tests demonstrated 98% correlation overall when approximately 100 clinical specimens (50 negative and 50 positive) were compared for each of the five drugs. The tests are similar in sensitivity, specificity, accuracy and precision.

Conclusion:

The Peace of Mind Multiple Drugs of Abuse Test and the Instacheck Drug Screen Panel are substantially equivalent in performance characteristics. The correlation of the two tests was 98%.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

AUG 7 2000

Ms. Janis Freestone
Director, Regulatory Affairs
Advantage Diagnostics Corporation
1201 Douglas Avenue
Redwood City, California 94063

Re: K001126
Trade Name: Peace of Mind Multiple Drugs of Abuse Test
Regulatory Class: II
Product Code: LDJ, DJG, DJC, LCM, DIX
Dated: June 26, 2000
Received: June 27, 2000

Dear Ms. Freestone:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in dark ink, appearing to read "Steven Gutman". The signature is fluid and cursive, with the first name "Steven" and last name "Gutman" clearly distinguishable.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

K001126


510k Number: K001126

Device Name:

Peace of Mind Multiple Drugs of Abuse Test.

Indications for Use:

The Peace of Mind Multiple Drugs of Abuse Test is a qualitative, one step, immunochromatographic competitive assay used to screen human urine for the presence of five Drugs of Abuse at the following cut offs: THC 50ng/mL, phencyclidine, 25ng/mL Opiates 2000 ng/mL, Cocaine 300ng/mL and Methamphetamine 1,000 ng/mL. The test is qualitative and provides only a preliminary analytical result, which must be confirmed by an alternate methodology preferably, GC/MS.


(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K001126

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over the counter use _____