

AUG 31 2000

ATTACHMENT H

SUMMARY OF SAFETY AND EFFECTIVENESS

Pursuant to §513(i)(3)(A) of the Food, Drug, and Cosmetic Act, Boston Scientific Corporation is required to submit with this Premarket Notification either an "... adequate summary of any information respecting safety and effectiveness or state that such information will be made available upon request of any person." Boston Scientific Corporation chooses to submit a summary of information respecting safety and effectiveness. According to §513(i)(3)(B), "Any summary under subparagraph (A) respecting a device shall contain detailed information regarding data concerning adverse health effects..."

The summary regarding the adverse health effects of the proposed GAZELLE Balloon Dilatation Catheter is as follows:

Trade Name: GAZELLE Balloon Dilatation Catheter

Manufacturer: Boston Scientific Corporation
Ballybrit Business Park
Galway, Ireland

Device Generic Name: Balloon Dilatation Catheter

Classification: According to Section 13 of the Federal Food, Drug and Cosmetic Act, the device classification is Class II, Performance Standards.

Predicate Devices: The following device is referenced in this premarket notification as the predicate device for the GAZELLE Balloon Dilatation Catheter:

Boston Scientific Corporation -- Symmetry™ Balloon Dilatation Catheter (K953602)

The device mentioned above has been determined substantially equivalent by FDA.

Device Description:

The GAZELLE Balloon Dilatation Catheter is an over-the-wire catheter designed for single operator exchange over guidewires which have outer diameters of 0.018" or smaller. The catheter is constructed with a double lumen shaft at the distal end and a single lumen shaft at the proximal end (single operator exchange/monorail catheter). The double lumen ends 22cm from the distal tip and is used for the advancement of the balloon segment over a guidewire. The single lumen shaft is joined to the double lumen shaft and is used for the inflation and deflation of the balloon.

Indications for Use:

The GAZELLE Balloon Dilatation Catheter is indicated for PTA of the iliac, femoral, ilio-femoral, popliteal, renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

Safety and Performance:

Bench testing and biocompatibility testing were performed to support a determination of substantial equivalence. The results of these tests provide reasonable assurance that the proposed device has been designed and tested to assure conformance to the requirements for its intended use. No new safety or performance issues were raised during the testing regimen.

Conclusion:

Based on the Indication for Use, technological characteristics and safety and performance testing, the GAZELLE Balloon Dilatation Catheter has been shown to be safe and effective for its intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 31 2000

Ms. Jennifer Bolton
Regulatory Affairs Specialist
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537

Re: K001134
GAZELLE Balloon Dilatation Catheter
Regulatory Class: II (two)
Product Code: 74 LIT
Dated: July 31, 2000
Received: August 2, 2000

Dear Ms. Bolton:

We have reviewed your Section 510(k) notifications of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submissions does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Jennifer Bolton

This letter will allow you to begin marketing your device as described in your 510(k) premarket notifications. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

James E. Dillard III
Director

Division of Cardiovascular and
Respiratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): New Application


Device Name: GAZELLE Balloon Dilatation Catheter

Indications for Use:

The GAZELLE Balloon Dilatation Catheter is recommended for the percutaneous transluminal angioplasty of the iliac, femoral, ilio-femoral, popliteal, renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED.)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Signatory)
Deputy Director
and _____
510(k) Number K001134

☒ Prescription Use

OR

☐ Over-The-Counter Use

(Per 21 CFR 801.109)

(Optional Format 1-2-96)