

SEP 21 2000

K001176



---

510-528 S. Vermont Avenue Glendora, CA 91741 (626) 914-2891 FAX (626) 914-2285

**Ref:** 510(k) Premarket Notification Summary

**To:** Document Control Clerk:

This is to notify you of the intention of OASIS Medical, Inc. to manufacture and market the following device:

**Microkeratome Blades - PE**

**Establishment Registration Number:** 2083373

This 510(k) summary of safety and Effectiveness for the OASIS Microkeratome Blades is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92, and follows the Office of Device Evaluation guidance concerning the presentation and content of a 510(k) summary.

1. Submitter's name, address, telephone number, contact person, and date the summary was prepared:
  - a. **Applicant:** OASIS Medical, Inc.  
514 South Vermont Avenue  
Glendora, CA 91741
  - b. **Telephone Number:** (626) 914-2891  
**Facsimile Number:** (626) 914-9372
  - c. **Contact Person:** Yvonne Fernandez- RA/QA Director
  - d. **Date Summary Prepared:** August 18, 2000
2. Name of the Device, including trade name, the common or usual name, and the classification:
  - a. **Trade/Proprietary Name:** Disposable Microkeratome Blade – PE
  - b. **Common/Usual Name:** Keratome Blade
  - c. **Classification Name:** Keratome (Blade Only) - 21CFR §886.4370
  - d. **Classification:** Class I
  - e. **Product Code:** 86 HNO
  - f. **Classification Panel:** Ophthalmic



---

510-528 S. Vermont Avenue    Glendora, CA 91741    (626) 914-2891    FAX (626) 914-2285

**3. Identification of legally marketed devices to which equivalence is being claimed:**

The OASIS Medical, Inc. disposable Microkeratome Blades are substantially equivalent in design, material and function to the devices as marketed by:

| <u>Company</u>              | <u>Device</u>  | <u>510(k) Number</u> |
|-----------------------------|----------------|----------------------|
| Chiron Vision Corporation   | ACS Blade Only | #K941550             |
| Plancon Instruments (Moria) | LSK Blade Only | #K970377             |

**4. Description of the Device:**

The OASIS Microkeratome-PE blades are replacement stainless steel blades for the Chiron Automatic Corneal Shaper blade and the Moria LSK-One blade. The two stainless steel blade styles have very slight differences in dimension to serve as keratome blade replacements. The catalog number 0410 blade is designed for use with the Chiron Automated Corneal Shaper (ACS) Keratome, while the 0414 blade is used for the Moria LSK-1 Keratome. Both styles are manufactured out of the same materials (Stainless Steel), packaged and sterilized using the same methods. The OASIS Microkeratome Blades are single-use, disposable blades.

**Certification of Safety and Effectiveness:**

When used according to the keratome manufacturer's instructions, there are no adverse safety indications for either the 0410 (ACS) or 0414 (LKS) blade.

**Sterilization Methodology:**

All blades are sterilized by exposure to gamma radiation to a Sterility Assurance Level (SAL) of  $10^{-6}$  according a validated process in compliance with AAMI/ISO 11137:1994.

**Labeling:**

The pouch will indicate OASIS name, address, product identification, lot number, sterilization notes, single use, and federal law statements.

**5. Intended Use for the Device:**

The OASIS 0410 ACS-PE and OASIS 0414 LSK-PE disposable microkeratome blades are designed as replacement blades for the Chiron Automated Corneal Shaper and Moria LSK-1 microkeratomes, respectively, for lamellar resection of the cornea.



510-528 S. Vermont Avenue    Glendora, CA 91741    (626) 914-2891    FAX (626) 914-2285

**6. Summary of the technological characteristics of the submitted device compared to predicate devices:**

**ACS Blade – Summary of Technological Characteristics of Device Compared to predicate device [Section 807.92(a)(6)]**

| <b>Characteristics</b>     | <b>PD* - Chiron ACS Blade</b>                | <b>Oasis 0410 ACS-PE</b> |
|----------------------------|--|--------------------------|
| Intended Use               | Replacement blades for Chiron's ACS keratome | Same                     |
| Portion to Contact Patient | Blade  | Same                     |
| Materials                  | Stainless steel (400 series)                 | Same                     |
| Sterilization Method       | Gamma Radiation                              | Gamma Radiation          |

**LSK Blade – Summary of Technological Characteristics of Device Compared to predicate device [Section 807.92(a)(6)]**

| <b>Characteristics</b>     | <b>PD* - Moria LSK Blade</b>                | <b>Oasis 0414 LSK-PE</b> |
|----------------------------|---|--------------------------|
| Intended Use               | Replacement blades for Moria LSK-1 keratome | Same                     |
| Portion to Contact Patient | Blade                                       | Same                     |
| Materials                  | Stainless steel (400 series)                | Same                     |
| Sterilization Method       | Gamma Radiation                             | Gamma Radiation          |

**Performance Tests and Conclusions:**

1. **Dimensional Equivalency** – Physical measurements of the predicate device are consistent with those of OASIS Medical, Inc. Fit into the respective Microkeratome has been tested and shown to be acceptable when used according to the keratome labeling.
2. **Sharpness Tests** – Non-clinical testing on porcine eyes resulted in corneal lamellar sections equivalent to the predicate devices. Corneal tissue photographs show equivalence in terms of quality of cut when used according to the keratome labeling.
3. **100% inspection** ensures zero defects. SEM photographs show equivalence in terms of surface finish and edge quality.



SEP 21 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Yvonne Fernandez  
Director  
OASIS Medical, Inc.  
514 S. Vermont Ave.  
Glendora, CA 91741

Re: K001176  
Trade Name: Disposable Microkeratome Blades  
Regulatory Class: I  
Product Code: 86 HNO  
Regulation: 886.4370  
Dated: August 18, 2000  
Received: August 21, 2000

Dear Ms. Fernandez:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-6413. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in cursive script that reads "A. Ralph Rosenthal".

A. Ralph Rosenthal, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health



510-528 S. Vermont Avenue Glendora, CA 91741 (626) 914-2891 FAX (626) 914-2285

**OASIS Medical, Inc.**  
**Disposable Microkeratome Blades**  
**Indications For Use**

510(k) Number (if known): K001176

Device Name: Disposable Microkeratome Blades

The OASIS 0410 ACS-PE and OASIS 0414 LSK-PE disposable microkeratome blades are designed as replacement blades for the Chiron Automated Corneal Shaper and Moria LSK-1 microkeratomes, respectively, for lamellar resection of the cornea.

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

-----  
Concurrence of CDRH, Office of Device Evaluation (ODE)

*Bruce Drum*

(Division Sign-Off) and ENT  
Division of Ophthalmic Devices

510(k) Number K001176

Prescription Use: X OR Over The Counter Use: \_\_\_\_\_  
(Per 21 CFR 801.109) (Optional Format 1-2-96)