## MAY 1 8 2000

## Premarket Notification (510(k)) Summary

510(k) Number:

KØØ1257

Product Name:

IntraCoil™ Stent

Common Name:

Tracheal prosthesis

Class:

per 21 CFR 878.3720 (tracheal prosthesis)

Submitter's Name:

Official Contact:

IntraTherapeutics, Inc.

Maria Brittle

651 Campus Drive St. Paul, MN 55112 Sr. Reg/Clin/Training Associate Telephone: 651-697-2018

Fax: 651-697-4808

Summary Preparation Date:

April 18, 2000

This summary is provided in compliance with section 513(I)(3)(A) of the Act and summarizes the safety and effectiveness information contained in this premarket notification submission. Substantial equivalence is claimed to the IntraTherapeutics, Inc. IntraCoil™ Stent, K990221.

The IntraCoil™ Stent is a self-expanding nickel-titanium (Nitinol) coil premounted on a delivery catheter. The stent is provided in diameters 4 to 8-mm, and lengths of 40 and 60 mm. The intended use is "in the treatment of bronchial strictures produced by malignant neoplasms or in benign strictures after all alternative therapies have been exhausted." Upon deployment the stent expands to conform to the bronchial lumen surface.

This 510(k) covers addition of the longer length, 60 mm, in diameters 4 to 7 mm. Otherwise, the device is identical to the IntraCoil™ Stent as cleared under K990221. A subset of the *in vitro* performance tests conducted for K990221, and relevant to the stent length, were repeated for design verification and product validation.



MAY 1 8 2000

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Maria E. Brittle Senior Regulatory/Clinical/Training Associate IntraTherapeutics, Inc. 651 Campus Drive St. Paul, Minnesota 55112

Re:

K001257

Trade Name: IntraCoil<sup>TM</sup> Stent

Regulatory Class: II Product Code: JCT Dated: April 18, 2000 Received: April 19, 2000

Dear Ms. Brittle:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

| 510(k) Number (if                                      | known):                 |  |
|--|-------------------------|--|
| Device Name:   | IntraCoil™ Stent        |  |
| Indications for Use                                    | <b>:</b> :              |  |
| The IntraCoil™ S<br>produced by mal<br>have been exhau | ignant neoplasms, or in | in the treatment of bronchial strictures<br>benign strictures after all alternative therapie |
|  |                         |  |
|  |                         |  |
|  |                         |  |
|  |                         |  |
|  |                         |  |
|  |                         |  |
| (PLEASE DO NO<br>NEEDED)                               | T WRITE BELOW THIS      | LINE-CONTINUE ON ANOTHER PAGE IF   |
|  | Concurrence of CDRH, O  | ffice of Device Evaluation (ODE)   |
|  |                         | (Division Sign-Off) Division of General Restorative Devices  510(k) Number KOO1257           |
|  |                         |  |
| _  |                         |  |
| Prescription Use 2<br>(per 21 CFR 801.1                | OR (109)                | Over-The-Counter Use   |