

JAN 12 2001

Attachment I
510(K) Summary

K00 1259

B.A.S.I.C. Dental Implant System Post and Core Attachment

This 510(K) Summary of safety and effectiveness for the B.A.S.I.C. Dental Implant System Post and Core Attachment is submitted in accordance with the requirements of the SMDA 1990 and following guidance concerning the organization and content of a 510(K) summary.

Applicant:	B.A.S.I.C. Dental Implant Systems, Inc
Address:	3321 Columbia NE Albuquerque, New Mexico 87107 USA
Contact Person:	Dan Blacklock, Vice-President
Telephone / Fax / Email	505.881.1376 – Phone 505.884.1923 - Fax
Preparation Date:	March 26, 2000
Device Trade Name:	B.A.S.I.C. Dental Implant System Post and Core Attachment
Common Name:	Accessory to a Dental Implant
Classification:	DZE
Legally Marketed Predicate Device:	B.A.S.I.C. Dental Implant System Post and Core Attachment K number K960868
Description of the Post and Core Attachment:	The Post and Core Attachment is a three-piece attachment. The top is attached to the base with a hex screw. The base is cemented into a dental implant. Artificial teeth are then attached to the Post and Core using conventional techniques.
Intended use:	The Post and Core Attachment is intended to attach artificial teeth to a dental implant.
Performance Data:	None
Results of Clinical Study:	None
Conclusion:	The Post and Core attachment is substantially equivalent to other existing Post and Core attachments in commercial distribution.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 12 2001

Mr. Dan Blacklock
Vice President
Basic Dental Implant Systems, Incorporated
3321 Columbia Northeast
Albuquerque, New Mexico 87107-2001

Re: K001259
Trade Name: B.A.S.I.C Dental Implant System Post and
Core Attachment
Regulatory Class: III
Product Code: DZE
Dated: October 16, 2000
Received: October 17, 2000

Dear Mr. Blacklock:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531


Page 2 - Mr. Blacklock

through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act ~~may be~~ obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATION FOR USE STATEMENT

510(k) Number: Pending KCC 1259

Device Name: B.A.S.I.C. Dental Implant System Post and Core Attachment

Indications for Use:

The B.A.S.I.C. Dental Implant System Post and Core Attachment is indicated for attaching artificial teeth to a dental implant.

(Please do not write below this line - Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(per 21 CFR 801.109)

OR

Over-the-Counter Use _____

Susan Rimmer

(Division Sign-Off)

Division of Dental, Infection Control,

General Hospital Devices

Number

K001259