

JUN 20 2000

K001311

14. **Summary of Safety and Effectiveness Information:**

510(k) SUMMARY

Submitter	Synthes (USA) 1690 Russell Road Paoli, PA 19301
Company Contact	Bonnie Smith (610) 647-9700
Name of the Device	Synthes (USA) Orbital Mesh Plate
Common or Usual Name	Single/multiple component metallic bone fixation appliance.
Predicate Device	Synthes (USA) Midfacial System
Device Description	Synthes Orbital Mesh Plates have a semi-circular shape with a radially designed mesh pattern. Orbital Mesh Plates are available in 0.2, 0.3 and 0.4 mm profile thickness. Standard 1.0 mm screw holes positioned along the outer arc of the Orbital Mesh Plate accept 1.0 mm self-tapping bone screws and 1.2 mm emergency screws. Synthes Orbital Mesh Plates for the Midfacial System are provided nonsterile.
Intended Use	Synthes Orbital Mesh Plates for the Midfacial System are indicated for selective trauma of the midface and craniofacial skeleton; craniofacial surgery, reconstructive procedures; and selective orthognathic surgery of the maxilla and chin.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 20 2000

Ms. Bonnie J. Smith, RAC
Senior Regulatory Affairs Associate
SYNTHES (USA)
1609 Russell Road
Post Office Box 1766
Paoli, Pennsylvania 19301

Re: K001311

Trade Name: Orbital Mesh Plates for Synthes Midfacial System
Regulatory Class: II
Product Code: HRS
Dated: April 24, 2000
Received: April 25, 2000

Dear Ms. Smith:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Bonnie J. Smith, RAC

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,


Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative and
Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

2. Indications for Use

Premarket Notification [510(k)]

INTENDED USE STATEMENT

510(k) Number (if known):

K001311

Device Name:

Synthes (USA) Orbital Mesh Plates

Indications

Synthes Orbital Mesh Plates for the Midfacial System are indicated for selective trauma of the midface and craniofacial skeleton; craniofacial surgery, reconstructive procedures; and selective orthognathic surgery of the maxilla and chin.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

Premarket Notification 510(k):

Synthes (USA) Orbital Mesh Plates for the Midfacial System
CONFIDENTIAL

Donna R. Lechner

(Division Sign-Off)

000004

Division of General Restorative Devices

Number K001311