

Attachment 4**510 (k) SUMMARY**

Submitter's name: Ann Marie Pahlman MPGR A-2E
Address: 1620 Waukegan Rd.
McGaw Park, IL 60085

Phone: 847-473-6078
Fax: 847-473-6952

Contact: Ann Marie Pahlman or Robert L. Wilkinson
Date Prepared: May 4, 2000

Trade name: Polysynthane (PSN™) Hemodialyzer
Common name: Hemodialyzer
Classification name: High Permeability Hemodialysis System
per 21 CFR 876.5860

Equivalent predicate: PSN™ Polysynthane Dialyzers

Device Description: Model PSN 170 Hemodialyzer
Model PSN 210 Hemodialyzer

Intended Use: Intended specifically for use in patients with acute or chronic renal failure when conservative therapy is judged to be inadequate. It may also be indicated in the treatment of patients intoxicated with poisons or drugs.

Summary of the technological predicate device: The general function and materials of the subject PSN™ Hemodialyzers are the same as the Baxter PSN™ Dialyzers cleared under K963933.

Clinical data: N/A

Conclusions drawn: Components of the subject PSN™ Hemodialyzers have met the biological requirements of ISO 10993-1: Biological Evaluation of Medical Devices – Part: Guidance on selection of tests. The validation of the sterilization cycle for the PSN™ Hemodialyzer is based upon the Association for the Advancement of Medical Instrumentation (AAMI) Guidelines (ST-27-Industrial Ethylene Oxide (EO))

K001422

510(k) SUMMARY
May 4, 2000
Polysynthane (PSN™) Hemodialyzer
Page 2 of 2

**Conclusions
drawn continued:**

sterilization of Medical Devices) to ensure a sterility assurance level (SAL) of 1×10^{-6} . Prior to release, sterilant residues of ethylene oxide and ethylene chlorhydrin are in compliance with the prolonged exposure category of ANSI/AAMI/ISO10993-7:1995 Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.

Pyrogen testing of the subject dialyzers meets the requirements of Chapter 161, Transfusion and Infusion Assemblies and Similar Medical Devices of Supplement 2 of the USP.


Particles are counted per USP XXIII <788>.

Functional testing for blood side integrity and conformance to manufacturing specifications are performed as in-process and/or final inspections prior to product release ensuring a quality product.

**Additional
information
requested by FDA:**

none to date.


Official Correspondent:



Robert L. Wilkinson
Director, Regulatory Affairs

5/4/00
Date

Prepared by:



Ann Marie Pahlman
Manager, Regulatory Affairs

5-4-2000
Date



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN - 2 2000

Ms. Ann Marie Pahlman
Manager, Regulatory Affairs
Baxter Healthcare Corporation
1620 Waukegan Road
McGaw Park, IL 60085

Re: K001422
PSN™ Hemodialyzer,
Models PSN 170 and PSN 210
Dated: May 4, 2000
Received: May 5, 2000
Regulatory Class: II
21 CFR§ 876.5860/Procode: 78 KDI

Dear Ms. Pahlman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure(s)

Attachment 2

Indications for Use Statement

510 (k)
Number

K001422

Device Name

PSN™ Hemodialyzer

Indications for
Use

Hemodialysis with the PSN™ Hemodialyzer is indicated for patients with acute or chronic renal failure when conservative therapy is judged to be inadequate. It may also be indicated in the treatment of patients intoxicated with poisons or drugs.


PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF
NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K001422