

JUN 2 0 2000



P.O. Box 591 Edendale 1610
Republic of South Africa
Telephone: +27 (11) 699-8640
Telefax: +27 (11) 452-8918

X061478

510 (k) Summary

Trade name: Aurecast Super Inlay
Common name: Dental casting alloy
Classification name: Gold based alloys and precious metal alloys for clinical use
Classification number: EJT

Legally marketed device: Jensen JRVT
Description of the device: High gold casting alloy
Intended use of the device: Type II restoration

Summary of the technological characteristics

Test methods applied: as in ANSI/ADA 5 and ISO 1562

Comparison of composition:

ALLOY		COMPOSITION (WEIGHT%)								
	Name	Au	Ag	Cu	Pt	Pd	Zn	In	Ir	Ru
Legal	Jensen JRVT	77.0	13.0	8.5	-	1.0	x	x	x	
New	Aurecast Super Inlay	75.1	14.8	6.9	1.5	-	1.4	-	-	0.3

x is less than 1%

Comparison of physical and mechanical properties

ALLOY		Melting point range (°C)		Hardness (Vickers 5/30)		Yield strength (Mpa)		Elongation (%)		CTE (x10-6/°C)	Density (g/cm3)
Name		solid	liquid	hard	soft	hard	soft	hard	soft		
Jensen JRVT		900	955	120		245		55			15.4
Aurecast Super Inlay		888	935	140	115	205	185	40	50		15.7

Discussion

Aurecast Super Inlay has a noble metal content of 76.9% which guarantees a high corrosion resistance. The compositional difference is in the base metals and it is less than 2% in all cases. Physical and mechanical properties of the two alloys are also very similar.

Conclusion

The main elements and their concentration are almost identical.

SK-29

RECEIVED

MAY 08 11 48

FDA/CDRH/ODE/CMO



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 20 2000

Mr. J.D. Davis
Chief Executive Officer
Aurex SA Pty Ltd.
24 Plantation Road, Eastleigh
Edenvale, Gauteng, Rep South Africa 1610

Re: K001478
Trade Name: Aurecast Super Inlay
Regulatory Class: II
Product Code: EJT
Dated: April 17, 2000
Received: May 11, 2000

Dear Mr. Davis:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Davis

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K001478

Device Name: AURECAST SUPER INLAY

Indications For Use:

Dental casting alloy for metallic restorations

It cannot be used in combination with dental porcelains

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Rumer
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K001478

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐