

JUL 11 2000



K001524 Page 1 of 1

## SECTION 10

### 510(k) SUMMARY

**Submitted By:** Richard L. Miller  
Manager – Regulatory Compliance  
Forma Scientific, Inc.  
Mill Creek Road  
Marietta, Ohio 45750

May 11, 2000

**Names of Devices:**

Trade Name: Forma Scientific Incubator  
Common/Usual Name: Embryo Incubator  
Classification Name: Embryo Incubator  
21 CFR 884.6120

**Predicate Device:** 63 FR 48428, September 10, 1998

**Device Description:**

The Forma Scientific incubators are bench top, or floor standing units. They control carbon dioxide (CO<sub>2</sub>), temperature and provide elevated humidity. Controlled parameters and alarm functions are microprocessor controlled. The volume of each chamber is approximately 6.5 cubic feet (184 liters).

**Intended Use:**

The Forma Scientific Direct Heat Incubators are intended to be used to store and preserve gametes and/or embryos at or near body temperature.

**Substantial Equivalence:**

In accordance with the Final Rule on reclassification of Medical Devices Used for In Vitro Fertilization, Forma Scientific cites the Final Rule as support for substantial equivalence.

**Discussion of Tests and Test Results:**

The Forma Scientific Direct Heat Incubators were subjected to electrical safety, electromagnetic compatibility acceptability and operating performance. The incubators passed all these tests.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL 11 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Richard L. Miller  
Manager of Regulatory Compliance  
ThermoQuest Corporation  
Forma Scientific Division  
Mill Creek Road  
Marietta, OH 45750

Re: K001524  
Forma Scientific Direct Heat Incubator  
Dated: May 11, 2000  
Received: May 16, 2000  
Regulatory Class: II  
21 CFR §884.6120/Procode: 85-MQG

Dear Mr. Miller:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.  
Captain, USPHS  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure(s)

510(k) Number (if known): K001524

Device Name: Embryo Incubator

Indications For Use:

The intended use of these incubators is to provide an environment with controlled temperature, CO<sub>2</sub> and an elevated humidity, for the development of ova or embryos.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Optional Format 3-10-98)

Prescription Use ✓  
(Per 21 CFR 801.109)

Christ A. Segerson  
(Division Sign-Off)  
Division of Reproductive, Abdominal, ENT,  
and Radiological Devices  
510(k) Number K001524