

NOV 17 2000

K001544

Summary of Safety and Effectiveness
Compliance with 513 (i) of the Federal Food, Drug and Cosmetic Act

August 21, 2000

1. General Provisions

Common/Usual Name: Applicator Manual Radionuclide Applicator System

Proprietary Name: Vaginal Applicator: Shielded

Applicant Name and Address:
 Mick Radio-Nuclear Instruments, Inc.
 1470 Outlook Avenue
 Bronx, New York 10465

2. Name of Predicate Devices:

(1)

Manufacturer	K Number
Mick Radio-Nuclear Instruments, Inc. Hilaris-Nori Endometrial Applicator	K871216/A
Mick Radio-Nuclear Segmented Vaginal Applicator	K993472
Aktina Medical Physics Corporation Beam Blocks	K974239

¹ Any statement made in conjunction with this submission regarding substantial equivalence to any other product only relates to whether the product can be lawfully marketed without pre-market approval or reclassification and is not to be interpreted as an admission or used as evidence in patent infringement litigation. As the Commissioner of the FDA has indicated, "...a determination of substantial equivalence under the Federal Food, Drug, and Cosmetic Act relates to the fact that the product can be lawfully marketed without pre-market approval or reclassification. This determination is not intended to have any bearing whatsoever on the resolution of patent infringement suits." 42 Fed. Reg. 42,520 et seq. (1977).

3. **Classification**

14061544

This device is classified as a class II (90-JAQ) device according to 21 CFR 892.5700.

4. **Performance Standards**

Performance standards for applicators for remote controlled afterloading brachytherapy have not been established by the FDA under Section 514 of the Food, Drug and Cosmetic Act.

5. **Intended Use and Device Description**

The Mick Radio-Nuclear Instruments, Inc. Shielded Vaginal Applicators are intended for use in Brachytherapy. The delivery of intracavitary radiation therapy requires not only proper visualization and localization of the applicator within the treatment volume, but precise dosimetry and a stable delivery system from which treatment can be administered. The Mick Radio-Nuclear Shielded Vaginal Applicator meets these requirements by providing a system that can be adjusted in diameter to meet the dimensions of the treatment volume and providing shields for blocking radiation exposure to non-involved tissue (similar to common practice in external beam radiation therapy with beam blocks (Aktina Medical Physics Corporation K974239). The system described in this 510(k) has been developed to function as an applicator for the positioning of sealed sources in the intracavitary treatment cancers of the vagina.

6. **Biocompatibility**

No new issues of biocompatibility are raised with regard to this device.

7. **Summary of Substantial Equivalence**

This device is similar in design and construction, utilizes the identical materials, and has the same intended use and performance characteristics to the predicate devices. No new issues of safety or effectiveness are introduced by using this device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 17 2000

Felix Mick
President
Mick Radio-Nuclear Instruments, Inc.
P.O. Box 99
Bronx, NY 10465

Re: K001544
Vaginal Applicator Set: Shielded (Catalog 8818)
Dated: October 3, 1999
Received: October 5, 1999
Regulatory class: II
21 CFR 892.5700/Procode: 90 JAQ

Dear Mr. Mick:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4639. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Daniel G. Schultz, M.D.
Captain, USPHS
Acting Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure (s)

Indications for Use

510(k) Number: *K001544*


Device Name: Vaginal Applicator: Shielded

Indications for Use:

The system described in this 510(k) has been developed to function as an applicator for the positioning of sealed sources in the intracavitary treatment cancers of the vagina.

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: ✓ or Over-The Counter Use: (Per 21 CFR 801.109)


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number *K001544*