

JUN 30 2000

**Section 510(k) Premarket Notification Summary
(as required by 807.92 (j))**

*K001682
Page 1 of 3*

Submitter: Vital Images, Inc.
3300 Fernbrook Lane North - Suite 200
Plymouth, MN 55447-5341
Phone: (612) 915-8001
Fax: (612) 915-8030

Date Prepared: May 30, 2000

Contact Person: Robert C. Samec

Device Trade Name: VScore™

Device Common Name: Image Processing Software for CT Scanners

Classification Name: ^{JAK CT} 90LLZ - System, Image Processing

Substantially Equivalent To:
CT Coronary Artery Calcification Scoring (CACS) (K982004)
General Electric Co.

Indications for Use: Cardiac scoring from whole body computed tomography derived measurements. For non-invasive detection and quantification of atherosclerotic plaque.

Device Description: The VScore EKG Gating Option for Coronary Artery Calcification Scoring (CACS) is an additional image processing option for K990442, Cardiac Scoring, which was subsequently marketed as VScore™ by Vital Images, Inc. This image processing option allows the operator to select images with reduced motion artifacts when processing data from General Electric CT Scanners for Coronary Artery Calcification Scoring.

A recording EKG machine is used to aid in the selection and review of images used for calcium scoring. The EKG device (In vivo research, Millennia Model 3500 CT-K950688) is used to record the patients EKG during scanning. A time stamp is placed on the EKG data to show when the imaging started. This allows the operator to select images with reduced motion artifacts.

Software Development: The software utilized was designed, developed, tested and validated according to written procedures. These procedures specify individuals within the organization responsible for developing and approving product specifications, coding and testing, validation and maintenance.

Performance Testing: The EKG Gating Option has successfully completed Integration testing/verification.

Clinical Evaluation: Software Beta testing will be successfully completed validating the EKG Gating Option integration into VScore prior to market release.

Substantial Equivalence Comparison Chart

System:	<u>General Electric</u> CT Coronary Artery Calcification Scoring (CACS) (K982004)	<u>Vital Images, Inc.</u> EKG Gating Option for VScore™ (K990442)
Intended Use:	Cardiac scoring from whole body computed tomography derived measurements.	Cardiac scoring from whole body computed tomography derived measurements.
Data Source:	CT Scanner (General Electric) (EKG DATA) In Vivo Millennia 3500CT (K950688)	CT Scanner (General Electric) (EKG DATA) In Vivo Millennia 3500CT (K950688)
Physical Characteristics: (Workstation)	DICOM 3.0 compatible Archive capability Manual segmentation/contour	DICOM 3.0 compatible Archive capability Manual segmentation/contour
Performance Measurement Testing:	Clinical Comparison Summary/Data	
Safety:	Physician review of data/scoring integral to use of feature.	Physician review of data/scoring integral to use of feature.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 30 2000

Robert C. Samec
Vice President QA/RA
Vital Images, Inc.
3300 Fernbrook Lane North, Suite 200
Plymouth, Minnesota 55447-5341

Re: K001682
VScore Image Processing Software with EKG Signal
Gating Option for Cardiac Scoring
Dated: May 30, 2000
Received: June 1, 2000
Regulatory class: II
21 CFR 892.1750/Procode: 90 JAK

Dear Mr. Samec:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.
Captain, USPHS
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure(s)

510(k) Number (if known): K 001682

Device Name: EKG Signal Gating Option for Cardiac Scoring (K990442)

INDICATIONS FOR USE:

Intended Use:

Indications for Use: Cardiac Scoring from whole body computed tomography derived measurements. For non-invasive detection and quantification of atherosclerotic plaque.


(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
Per 21 CFR 801.109

OR

Over-The-Counter Use) _____

(Optional Format 1-2-96)



(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K001682